Medtronic

EverCross™

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
1 Device description

The EverCross PTA Catheter is an over-the-wire (OTW) dual lumen catheter (1). The catheter is compatible with 0.035 in (0.89 mm) guidewires. A semi-compliant inflatable balloon is mounted at the distal end (2). The distal atraumatic tip is tapered (3). The distal portion of the catheter has a hydrophilic lubricious coating (4).

The luer marked THRU (6) on the manifold (5) is the proximal opening for the central lumen. The central lumen ends at the distal tip of the catheter. This lumen is used to pass the catheter over a guidewire with a maximum diameter of 0.035 in (0.89 mm). The luer marked BALLOON (7) is used to inflate and deflate the dilatation balloon, using a mixture of contrast medium and saline solution. The balloon has 2 radiopaque marker bands (8) for positioning the balloon relative to the stenosis. The radiopaque marker bands indicate the dilating or working section of the balloon (9).

Figure 1. EverCross PTA Catheter

The EverCross PTA Catheter is available in multiple balloon sizes. Nominal balloon diameter and length are printed on the manifold. Refer to the package label for information about catheter length and sheath compatibility.

All balloons expand to sizes above the nominal size at pressures greater than the nominal pressure. Consult the balloon compliance card that is packaged with the device for balloon diameters at given pressures.

2 Intended purpose

The intended purpose of the EverCross PTA Catheter is to restore lumen patency and blood flow to peripheral arteries, or through native arteriovenous dialysis fistulae, or to facilitate stent post-dilatation after placement in peripheral arteries.

2.1 Intended patient population

The EverCross PTA Catheter is intended for adult patients in need of percutaneous transluminal angioplasty to improve and maintain luminal diameter in peripheral arteries or native arteriovenous dialysis fistulae.

2.2 Indications for use

The EverCross PTA Catheter is indicated to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, and infrapopliteal arteries, and for the treatment of obstructive lesions of native arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral arteries.

2.3 Contraindications

The EverCross PTA Catheter is contraindicated for use in the coronary arteries or in the neurovasculature, or when unable to cross the target lesion with a guidewire.

2.4 Clinical benefits

The clinical benefits of restoring lumen patency and blood flow with the EverCross PTA Catheter in the treatment of peripheral arterial occlusions or lesions, or native arteriovenous dialysis fistulae include: resolution of arterial disease related symptoms, improvement in quality of life, increased mobility, and preservation of hemodialysis access sites.

2.5 Intended users

Interventional procedures using the EverCross PTA Catheter should only be performed by physicians who have experience with interventional techniques in the vascular system.

This device should only be used in a clinical setting in a sterile environment.

2.6 Performance characteristics

The EverCross PTA Catheter is a dual lumen over-the-wire catheter with a balloon that is bonded onto the distal tip. The overall design and material selection allow for treatment of stenotic lesions in the peripheral arteries. The materials provide flexibility to move through tortuous anatomy and the tip design provides ease in crossing lesions. The balloon is inflated to its predetermined diameter and exerts a radial force to establish patency.

3 Warnings

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

- To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.
- Never use air or any gas to inflate the balloon.
- Use direct fluoroscopic observation while manipulating the catheter in the vascular system.
- Do not manipulate the balloon when it is inflated.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. If resistance is felt, carefully withdraw the dilatation catheter.
- Consider precautions to prevent or reduce clotting. Physician experience and discretion determines the appropriate anticoagulation therapy for each patient.
- Use the catheter with caution for procedures involving calcified lesions or vascular grafts.
- Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.

4 Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the sterile packaging has been unintentionally opened or the device is damaged.
- Balloons should not be inflated in excess of the rated burst pressure.
- Do not use contrast media that are contraindicated for intravascular use with this device.
- Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to the product labeling for sizing information.
- Ensure that the guidewire is in place before changing the position of the balloon catheter.
- Larger and longer models of the catheter may exhibit longer deflation times, particularly on long catheter shafts.
- Do not retract the balloon catheter after the stent has been fully dilated, unless the balloon is free from the stent and fully deflated under vacuum. Prematurely retracting the catheter can cause vessel trauma or product damage.
- To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Do not wipe the catheter surface with dry gauze.
- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

5 Potential adverse events

The potential adverse events (or complications) that may occur or require intervention with the use of this device include, but are not limited to, the following conditions:

- Allergic reaction
- Aneurysm
- Arrhythmia
- Arteriovenous fistula
- Bleeding requiring transfusion
- Contrast medium reaction or renal failure
- Death
- Embolism
- Endocarditis
- Fever
- Hematoma
- Hemorrhage
- Hypertension or hypotension
- Infection
- Inflammation
- Intraluminal thrombus
- Myocardial infarction
- Pseudoaneurysm
- Regional ischemia
- Sepsis
- Shock
- Stenosis or restenosis
- Stroke
- Transient ischemic attack
- Venous thromboembolism
- Vessel injury (such as dissection, perforation, or rupture)
- Vessel spasm

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority.

6 How supplied

The device is packaged with a protective sheath positioned over the balloon.

The device is supplied sterile and is sterilized using ethylene oxide.

Warning: This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

7 Storage

Store the device in a dry area away from sunlight.

Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

8 Device selection

Consider the following details when selecting a catheter:

- The diameter of the expanded balloon should not exceed the diameter of the artery immediately distal to or proximal to the stenosis.

Warning: To reduce the potential for vessel damage, select a balloon size with an inflation diameter that approximates the diameter of the lumen at the intended inflation site.

- Verify that the selected accessories accommodate the balloon catheter based on sizing information on the packaging label. Refer to the product labeling for size information.

Caution: Carefully inspect the sterile package and the device before use. Do not use the device if the sterile packaging has been unintentionally opened or the device is damaged.

9 Recommended items

Prepare the following items using a sterile technique:

- 10 mL syringe filled with sterile heparinized saline
- Three-way stopcock
- Contrast medium (The standard inflation medium is a 1:1 mixture of contrast medium and sterile saline.)

Caution: Do not use contrast media that are contraindicated for intravascular use with this device.

- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)
- Inflation device with manometer

10 Preparation

Note: As supplied, the balloon and the guidewire lumen of the EverCross PTA Catheter contains air. This air must be displaced to make certain that only liquid fills the balloon while the catheter is in the bloodstream. Refer to *Section 10.2* for air removal instructions.

- Carefully remove the catheter from the inner package.
- Remove and discard the protective tubing.

10.1 Activating the hydrophilic coating

Activate the hydrophilic coating before inserting the catheter. To activate the hydrophilic coating, immerse the catheter in normal saline for approximately 30 to 60 seconds, or wipe the catheter shaft with a saturated gauze sponge.

Caution: Do not wipe the catheter surface with dry gauze.

10.2 Displacing air

- Flush out the lumen for the guidewire with heparinized saline solution through the wire luer marked THRU.
- Mix equal volumes of contrast media and normal saline.
- Fill a 10 mL syringe with approximately 4 mL of the mixed contrast medium solution.

Warning: Before beginning a percutaneous transluminal angioplasty procedure, identify allergic reactions to contrast media.

Caution: Do not use contrast media that are contraindicated for intravascular use with this device.

- Evacuate air from the balloon segment:
 - Attach the syringe to the luer marked BALLOON. Hold the syringe pointing down with the plunger pointed upward.
 - Apply negative pressure and aspirate the lumen for 15 seconds. Slowly release the pressure to neutral so that the contrast medium fills the shaft of the dilatation catheter.
 - Disconnect the syringe from the BALLOON luer of the dilatation catheter.
 - Remove all air from the syringe barrel. Reconnect the syringe to the BALLOON luer of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
 - Slowly release the device pressure to neutral.

Note: On first use, avoid any positive pressure to prevent the balloon from unfolding and partially inflating outside of the vessel.

- Repeat steps 4b to 4e as necessary to remove all air from the system.
- Replace the syringe with a manometer-controlled dilatation system. Do not introduce air into the system.

11 Insertion and dilatation

- If applicable, use an appropriately sized introducer sheath to introduce the catheter percutaneously. Refer to the product labeling for sizing information.

Caution: Do not attempt to pass the catheter through an introducer sheath that is smaller than the size indicated on the product label. Refer to the product labeling for sizing information.

- After the guidewire is in place through the lesion, place the prepared catheter over the guidewire and advance the catheter tip to the introduction site. Use a guidewire that is of suitable length to maintain control and position of the guidewire at all times.

Warning: Use direct fluoroscopic observation while manipulating the catheter in the vascular system.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. If resistance is felt, carefully withdraw the dilatation catheter.

Caution: To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.

- Position the catheter with the center of the balloon in the stenotic area of the vessel. The radiopaque marker bands indicate the working length of the balloon.

Caution: Ensure that the guidewire is in place before changing the position of the balloon catheter.

- When the radiopaque marker bands are appropriately positioned, inflate the balloon to dilate the target area.

Warning: Never use air or any gas to inflate the balloon.

Warning: Do not manipulate the balloon when it is inflated.

Caution: Balloons should not be inflated in excess of the rated burst pressure.

- Using the inflation device, deflate the balloon by aspirating the catheter.

Caution: Larger and longer models of the catheter may exhibit longer deflation times, particularly on long catheter shafts.

- Withdraw the catheter while maintaining a vacuum in the balloon.

Caution: Do not retract the catheter unless the balloon is free and fully deflated under vacuum. Prematurely retracting the catheter can cause vessel trauma or product damage.

Note: Multiple balloon inflations and deflations may cause resistance during device withdrawal. If resistance is encountered, use a gentle counterclockwise twisting motion to ease withdrawal through the sheath. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath together, while maintaining vessel access with the guidewire.

- Use angiography to verify dilatation results.

12 Disposal

Caution: Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

13 Disclaimer of warranty

Important: This disclaimer of warranty does not apply in any countries where such a disclaimer is not permitted by law.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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© 2023 Medtronic
Medtronic and Medtronic logo are trademarks of Medtronic. ™** Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic, Inc.
Medtronic B.V.
Medtronic Limited