



BENCHMARK 071

INTRACRANIAL ACCESS SYSTEM



Penumbra 

BENCHMARK 071

Key Design Features

.071" (1.80 mm) Inner Diameter
Optimised for Delivering
Embolisation Devices

**Soft 9 cm Platinum Coil Wind
Flexible Segment**
Designed for distal trackability

**Advanced Polymer
Construction**
Enables navigation and visibility

**Hybrid Flat/Round
Stainless Steel Reinforcement**
Enables arch support and
navigation in tortuous vessels

Polished, Atraumatic Tip
Designed for distal placement

10 Material Transitions
Designed for efficient energy
transfer and trackability


Proximal Shaft Construction
Engineered to provide stability
and support

**2 mm Outer Diameter
Low Profile System**
Compatible with 6 F
femoral and radial sheaths

BENCHMARK 071 packaged with 5 F Penumbra Select Catheter


Advantages of Rapid Primary Access

- Faster procedures – no need for an exchange
- Allows easy selection off of arch into desired vessel
- Facilitates atraumatic placement into distal vasculature
- Can be used for diagnostic angiogram – .040" (1.02 mm) lumen
- Compatible with .035"–.038" (.89 mm–.97 mm) guidewires




Performance Testing

Designed to Maximise Arch Support



Buckling Test: Higher values correspond with more proximal support^a

Advanced Level of Tip Softness




Deflection Test: Lower values correspond to softer distal tip^b


Femoral Access, Defined

Designed for stability and compatibility with atraumatic navigation


Standalone or compatible with 6 F long sheath




Rapid Primary Access — Typical Approach via Femoral Artery




1 Introduce the preloaded system over a guidewire and advance to a straight section of the abdominal aorta




2 Advance the tip of 5 F Select over the guidewire into the ascending aorta while maintaining the position of BENCHMARK 071




3 Advance 5 F Select and BENCHMARK 071 over the guidewire into the internal carotid artery



4 With sufficient wire purchase, torque 5 F Select towards the first major turn








5 While maintaining position of 5 F Select and guidewire, advance BENCHMARK 071 into desired position



6 Remove 5 F Select while holding BENCHMARK 071 in position



BENCHMARK 071 Intracranial Access via Femoral Approach




Images courtesy Drs. Yasha Kayan and Josser Delgado
Abbott Northwestern Hospital, Minnesota, USA

Radial Advantage


Low profile 2 mm OD compatible with the typical radial artery



Stent-Assisted Coiling
Dr. Levansri Makalanda
The Royal London Hospital
London, UK




Primary Coiling
Dr. Jean Delbrune
Northside Hospital
Florida, USA



MicroVention® WEB® Embolisation
Dr. Justin Singer
Spectrum Health
Michigan, USA

Rapid Primary Access — Typical Approach via Radial Artery



Procedural and operative techniques and considerations are illustrative examples from physician experience. Physicians' treatment and technique decisions will vary based on their medical judgment.

Data Presented at SNIS 2018

Radial Access for Cerebrovascular Intervention Using Penumbra BENCHMARK 071 Guiding Catheter^c

36 Patients	
10	Primary Coiling
8	Stent-Assisted Coiling
5	Flow-Diverter Embolisation
2	Balloon-Assisted Coiling
2	Wingspan™ Stent-Assisted Coiling
2	AVM or Dural AV Fistula
2	Vessel Sacrifice
1	Subclavian Stent
1	Vasospasm Treatment
3	Other

Key Results

No catheter-related complications

No major radial access site complications

Images used with permission. Consents on file at Penumbra, Inc. Case examples presented are for informational purposes only. Results may not be predictive for all patients and may vary based on patient-specific attributes and other factors.

c. Satti S, Rastogi S, Eden T, et al. E-031 Radial access for cerebrovascular intervention using penumbra benchmark 071 guiding catheter. J NeuroInterv Surg. 2018;10:A63-A64. doi: 10.1136/neurintsurg-2018-SNIS.107






BENCHMARK™ 071 Kits

BENCHMARK 071				5 F Select™ Catheter	
Catalog Number	Description	Length (cm)	Shape	Length (cm)	Shape
BMK6F95BER120	BENCHMARK 071 KIT	95	Straight	120	BER
BMK6F95MBER120	BENCHMARK 071 KIT	95	MP	120	BER
BMK6F105BER130	BENCHMARK 071 KIT	105	Straight	130	BER
BMK6F105MBER130	BENCHMARK 071 KIT	105	MP	130	BER

BENCHMARK 071

BENCHMARK 071			
Catalog Number	Description	Length (cm)	Shape
BMK6F95	BENCHMARK 071	95	Straight
BMK6F95M	BENCHMARK 071	95	MP
BMK6F105	BENCHMARK 071	105	Straight
BMK6F105M	BENCHMARK 071	105	MP
BMK6F115	BENCHMARK 071	115	Straight
BMK6F115M	BENCHMARK 071	115	MP

Tip Shapes

BENCHMARK 071		Select Catheter		
				
Straight	MP	H1	BER	SIM

5 F Select Catheters

Catalog Number	Description	Working Length (cm)	Inner Diameter (in / mm)	Wire Compatibility (in / mm)	Shape
PNS5F120BER	5F Select Catheter	120	.040 (1.02)	.035-.038 (.89-.97)	BER
PNS5F130BER	5F Select Catheter	130	.040 (1.02)	.035-.038 (.89-.97)	BER
PNS5F130SIM	5F Select Catheter	130	.040 (1.02)	.035-.038 (.89-.97)	SIM
PNS5F120H1	5F Select Catheter	120	.040 (1.02)	.035-.038 (.89-.97)	H1
PNS5F130H1	5F Select Catheter	130	.040 (1.02)	.035-.038 (.89-.97)	H1

Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use.

BENCHMARK Intracranial Access System – Intended Use

The BENCHMARK Intracranial Access System is intended for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

Potential Adverse Events

Possible complications include, but are not limited to, the following: acute occlusion; air embolism;

death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

NEURON MAX System – Intended Use

The NEURON MAX System is intended for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

Potential Adverse Events

Possible complications include, but are not limited to, the following: acute occlusion; air embolism;

death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

NEURON Intracranial Access System – Intended Use

The NEURON Intracranial Access System is intended for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

Potential Adverse Events

Possible complications include, but are not limited to, the following: acute occlusion; air embolism;

death; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at puncture site; infection; intracranial hemorrhage; ischemia; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation.

Penumbra 

Penumbra, Inc. USA
One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

Penumbra Europe GmbH
Am Borsigturm 4413507 Berlin
Germany
T +49 30 2005 676-0
F +49 30 2005 676-10
de-order@penumbrainc.com
info@penumbrainc.de

Penumbra Neuro Australia Pty Ltd
Suite 3, Level 5, 1 Oxford Street
Darlinghurst NSW 2010
Australia
T +61-1300 817 025
F +61-1300 817 026
order.anz@penumbrainc.com

Penumbra Latin America Distribuidora de Equipamentos e Produtos Médicos Ltda
Avenida Brigadeiro Luís Antônio
3421 cj 201 CEP 01401-001
São Paulo, Brazil
T 5511.2883.5825
order.la@penumbrainc.com

Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Renderings for illustrative purposes only. Images used with permission. Consents on file at Penumbra, Inc. Please contact your local Penumbra representative for more information.

Copyright ©2020 Penumbra, Inc. All rights reserved. The Penumbra P logo, BENCHMARK, Neuron, MAX, Neuron MAX, and Select are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. All other trademarks are the property of their respective owners. 16720, Rev. A 02/20 OUS