Reacor D

External pacemaker
en • Technical Manual

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en • English

© BIOTRONIK SE & Co. KG
All rights reserved. Specifications
subject to modification, revision
and improvement.

® All product names in use may
be trademarks or registered
trademarks held by BIOTRONIK or
the respective owner.

CE0482

Revision: Y (2025-03-17)

480533--Y

BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin / Germany
Tel +49 (0) 30 68905-0
sales@biotronik.com
www.biotronik.com

 **BIOTRONIK**
excellence for life

en • English	2
bg • Български	35
cs • Český	70
de • Deutsch	103
el • Ελληνικά	137
es • Español	171
fr • Français	204
hu • Magyar	238
it • Italiano	271
pl • Polski	304
pt • Português	338

Table of Contents

1	About this Technical Manual	4
	Technical Manuals	4
	Conventions	4
2	Product Description	6
	Intended Purpose	6
	Indications	6
	Contraindications	7
	Intended Clinical Benefit	7
	Patient Target Group	7
	Intended Part of the Body or Type of Tissue	7
	Expected Lifetime	7
	Product Identification	7
	Intended Conditions of Use	8
3	Safety	9
	Reporting of Serious Incidents	9
	Residual Risks	9
	Technical Safety	9
	Electromagnetic Interference	10
4	First Steps	12
	Checking the Package Contents	12
	Inserting the Battery	12
	Self-Test	13
5	Operating Elements and LEDs	14
6	Operating the Reacor D	17
	Connecting the Redel Adapter	17
	Patient Cable with Redel Plug	18
	Patient Cable Without Redel Plug	18
	Setting Parameters	19
	Connecting the Leads	20
	Closing the Protective Cover	21
	Securing	21
7	Pacing Modes	22

8	Handling	24
	Changing the Battery	24
	Cleaning	24
	Disinfection	24
	Sterilization	25
	Maintenance and Inspection	25
	Reuse	25
	Disposal	25
9	Appendix	26
	Technical Data	26
	Parameters	26
	Connection Variants	29
	Accessories and Spare Parts	31
	Adapters for PK-67-L and PK-67-S	32
	Symbols	33
	Legend for the Label	34

1 About this Technical Manual

Technical Manuals

Technical manuals are either included in hard copy form in the storage package or available in digital form on the internet: <https://manuals.biotronik.com>.

1. Consult the technical manual.
2. Keep the technical manual for future reference.

To ensure safe operation, in addition to this technical manual, please also consult the following technical manuals:

- Technical manuals for cables, adapters, accessories

Conventions

Identification of Safety Warnings

The following symbol draws attention to potential safety risks:

Follow all safety warnings indicated by this symbol to avoid possible serious or fatal injury or damage to the system.

Safety warnings are also indicated by a signal word.

- **DANGER:** Non-compliance may immediately lead to severe injury or death.
- **WARNING:** Non-compliance leads to a potentially dangerous situation that can cause severe injuries or death.
- **Caution:** Non-compliance leads to a potentially dangerous situation that can cause moderate injuries.
- **Attention:** Non-compliance leads to a potentially dangerous situation that can cause minor injuries and/or material damage.

Typographical Conventions

The following typographical conventions are used in this technical manual:

Instructions

The individual steps of an instruction are numbered. Prerequisites, intermediate results, and results may be specified.

Prerequisite

- This is a prerequisite.

1. Step
2. Step
 - ▶ Intermediate result
3. Step

Result

This is the result.

Cross references

Cross references are indicated using "see" or "see also".

Highlights

Text that needs to be highlighted is shown in **bold**.

Figures

Figures that show the product or the user interface are used for illustration. The details shown in the figure may differ from that of the delivered product or your software version.

Gender

If personal designations are used only in the male or female form for reasons of better readability, this expressly includes all other gender identities.

2 Product Description

Intended Purpose

Reocor D is a battery-powered, external dual-chamber pacemaker. In conjunction with compatible patient cables and temporary pacemaker leads (including myocardial heart wires and transvenous implantable catheters) the device is intended for dual-chamber temporary pacing and termination of bradyarrhythmias in a clinical environment.

The primary functions and associated objectives of the Reocor D are:

- Sensing of the intrinsic heart rhythm and heart rate to automatically detect bradycardia
- Pacing by electrical pulses of low energy to ensure a stable heart rate or to support the intrinsic heart rate when needed
- High-rate pacing (burst function) to treat atrial tachycardia

External temporary pacemakers ensure temporary symptomatic therapy for patients requiring pacing or pacing backup for safety reasons in emergency situations, during and after heart surgery or during treatment of infections that require an explantation of an implanted device. In these situations, the temporary pacemaker serves as a bridge until patient's recovery or permanent pacemaker implantation.

Indications

Temporary pacing with Reocor D is suitable for the following applications for patients of any age:

- Heart block
- Acute myocardial infarction
- Sinus node disease
- Bradycardia (sinus bradycardia, sinus pause, sinus arrest, AV nodal escape rhythm)
- Permanent pacemaker failure
- Atrial and ventricular arrhythmias
- Asystole or cardiac arrest
- Prophylactic or periprocedural temporary pacing (prophylactic, diagnostic, required for procedure)
- Bridge to permanent pacemaker
- Pre-, intra- and postoperative pacing of patients with heart surgery
- Atrial tachyarrhythmias requiring high-rate pacing for treatment
- Treatment of supra-ventricular tachycardia (SVT) (atrial flutter, atrial fibrillation, AV reentrant tachycardia)
- Drug toxicity-induced bradycardia arrhythmias

Contraindications

For Reocor D, the following contraindications exist:

- The device cannot be sterilized. If no sterile cover is applied, it is not suitable for use within the sterile field.
- The atrially triggered pacing modes (DDD and VDD) are contraindicated for atrial fibrillation, atrial flutter and other fast atrial rhythms.
- When high ventricular rates are not well tolerated by the patient (e.g., in the presence of angina pectoris), atrium-controlled modes can be contraindicated.
- Atrial single-chamber pacing is contraindicated for patients with existing AV conduction disturbances.
- The use of an external pacemaker is contraindicated in the presence of an active, implanted pacemaker.
- MRI scan procedures are contraindicated in the presence of an external pacemaker.

Intended Clinical Benefit

The clinical benefits for the patients inherent to the use of the temporary pacemaker are the detection of an unphysiological low heart rate and unphysiological rhythm (bradycardia arrhythmia) and, subsequently, the restoration of a physiological heart rate.

Atrial therapy provides the additional clinical benefit of terminating stable atrial tachyarrhythmia after detection.

Patient Target Group

The external pacemaker is intended for patients undergoing cardiac surgery and pacemaker-dependent patients who are temporarily not supported by an implanted device. The patient population includes pediatric patients including neonates and infants.

Intended Part of the Body or Type of Tissue

The external pacemaker is attached to an infusion stand, the patient's bed or the patient via an armband fastened over the clothing. Therefore, it has no direct contact with the patient's body.

The external pacemaker interacts with the body of the patient via a compatible lead system.

Expected Lifetime

The lifetime of the external pacemaker is 12 years.

Product Identification

Each product is identified by a so-called unique device identification (UDI), which contains product-specific information. In addition, a basic identifier called B-UDI-DI (basic – unique device identification – device identifier) is assigned to several products.

The B-UDI-DI for the external pacemaker is: 4035479BUDI00072QB.

Using this B-UDI-DI, it will be possible to search the European Database on Medical Devices (EUDAMED) for additional information on the product.

Intended Conditions of Use

Application Environment

Reocor D is used in health care facilities with equipment according to international guidelines (ECG monitor, defibrillator):

- In an operating theatre or a cath lab of a hospital outside of the sterile field
- At the patient's bedside inside the hospital
- Mobile use inside the hospital attached to the patient or an infusion stand

Intended User

Reocor D is used for temporary pacing by trained technical and medical personnel of a hospital and physicians. In addition to having basic medical and cardiological knowledge, the user must be thoroughly familiar with the operation and the operating conditions of the device system.

Only qualified medical specialists having this required special knowledge are permitted to use the devices. If users do not possess this knowledge, they must be trained accordingly.

BIOTRONIK offers user trainings for specific target groups. Current information on training and education opportunities can be requested from: education.training@biotronik.com.

Frequency and Duration of Application

Duration and frequency of usage is strongly dependent on the application manner of the device. The duration of usage can range from 1 hour to 7 days. The frequency of usage can range between temporary manner, such as in an emergency situation, to 3-5 applications per month for intensive care.

Compatibility

For compatible and accessories, see Accessories and Spare Parts [Page 31].

3 Safety

Reporting of Serious Incidents

Report any serious incident that has occurred with the device to the manufacturer and to the competent authority.

The competent authority can be found at: <https://ec.europa.eu>

Residual Risks

- Infection in blood circulation
- Ongoing ventricular tachycardia
- Pain
- Toxic / allergic reaction (acute and chronic)
- Nausea / sickness / slight dizziness
- Syncope
- Prolonged undesired medical condition

Technical Safety

The external pacemaker Reocor D meets the international standards for the safety of electromedical devices according to IEC 60601-1 and IEC 60601-1-2, as well as the international standard IEC 60601-2-31 for temporary, external pacemakers.

The following special features offer safety for the patient:

- No metal parts that can be touched by the operator in accordance with the requirements of the IEC 60601-1 standard.
- The external pacemaker has an applied part (part that touches or is electrically connected to the patient) of type CF (cardiac floating) and is approved for direct treatment of the heart. The external pacemaker complies with the requirements for defibrillation protection stipulated in the international standards.
- Due to the closed protective cover, the parameters for the external pacemaker cannot be adjusted.

WARNING

When handling already implanted leads, their connector pins and metal contact surfaces must not touch or come into contact with electrically conductive or damp surfaces.

WARNING

Before using the pacemaker, the patient cable, or the leads, the user should touch the patient to balance differences in electrostatic potential.

Electromagnetic Interference

Possible Electromagnetic Interference

WARNING

Risk of electromagnetic interference

The use of this device close to or stacked with other devices should be avoided, as this may lead to the device operating incorrectly.

- Where usage in such a manner is unavoidable, you should monitor this device and the other device(s) to ensure correct operation.

WARNING

Risk of electromagnetic interference through the use of portable RF communication equipment

If portable RF communication devices (including peripheral devices such as antenna cables and external antennae) are operated closer than 30 cm (12 inches) from this device, this can result in a reduction in its performance. This applies even when using the specified cables.

- When operating portable RF communication devices (including peripheral devices such as antenna cables and external antennae), keep such devices at a distance of at least 30 cm (12 inches) from this device.

WARNING

Risk of electromagnetic interference through the use of unauthorized accessories

The use of accessories, transducers or cables not listed by BIOTRONIK or of accessories other than those specified by BIOTRONIK, can produce elevated electromagnetic emissions or cause degradation in the device's resistance to electromagnetic interference. Such effects can lead to the faulty operation of the device.

- Only use accessories authorized by BIOTRONIK.

This device is protected against electromagnetic interference and electrostatic discharges in the specialized environment in the vicinity of high-frequency surgical instruments. At the same time, the emitted interference is reduced to a minimum. The device thus meets all requirements of IEC 60601-1-2.

Section of IEC 60601-1-2	Test	Test level
7.1	EN 55011 (CISPR 11) Radiated emission	- Group 1 - Class B
8.9	IEC 61000-4-2 Electrostatic discharge (ESD)	- Air discharge ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV - Contact discharge ± 8 kV
	IEC 61000-4-3 Electromagnetic fields	- Modulation: 2 Hz 3 V/m - Limits for RF communication equipment per Table 9 in IEC 60601-1-2
8.9	IEC 61000-4-6 Conducted high-frequency interference	- Modulation: 2 Hz 3 V 6 V ISM
	IEC 61000-4-8 AC frequency magnetic fields	30 A/m 50/60 Hz
8.11	IEC 61000-4-39 Proximity magnetic fields	Limits according to Table 11 in IEC 60601-1-2 (7.5–65 A/m)

- A value set to < 1.0 mV for the device's atrial sensitivity may result in interferences from electromagnetic fields. If clinically possible, a sensitivity value ≥ 1.0 mV should be set. Programming sensitivity values to < 1.0 mV requires specific clinical indication. Such values may only be set and retained under physician supervision.
- In the case of interference with modulation frequencies in the range of heart rates (approx. 30–300 bpm) or single pulse events, the pacemaker inputs can detect false cardiac activity. In this case, selecting a pacing mode without sensing (x00) is recommended.

4 First Steps

Checking the Package Contents

The Reacor D is delivered in a padded carrying case.

Included in the delivery are the Reacor D, a battery, the Redel adapter, the protective cover, the armband, and the technical manual.

Attention

Before each use, visually inspect the external pacemaker and its accessories for detectable damage, deformation, loose parts, dirt, and moisture.

Never use a damaged device, but return it to BIOTRONIK.

WARNING

Only the manufacturer BIOTRONIK or a party expressly authorized by the manufacturer may perform preventative or corrective maintenance, enhancements, or modifications to the device.

Inserting the Battery

The external pacemaker only runs on a battery.

Attention

If you use a battery other than the one supplied, do not use rechargeable batteries.

The service time of these batteries is difficult to estimate.

Use a recommended battery, see Changing the Battery [Page 24].

1. Open the battery magazine by reaching into the blue recess on the right-hand side of the device and pushing the slider upwards.
2. Pull the battery magazine all the way out.
3. If necessary, remove the rubber plug that may be stuck on the battery poles.
4. Put the battery in the battery magazine.
5. Close the battery magazine.
6. Push the slider downwards, until it snaps in place with an audible click.

WARNING

If you do not use the external pacemaker for a long time, take out the battery to prevent leakage.

Fig. 1: Open the battery magazine

Attention

Insert battery correctly

The preferred pole orientation is marked in the battery magazine, the position of the plus and minus poles can be freely selected.

- When inserting the battery, ensure that the battery poles point to the middle of the housing.

Self-Test

Power on the Reocor D by selecting any pacing mode on the Mode control dial. After powering on the device, a self-test is always carried out, which only lasts a few seconds.

The self-test includes:


- Checking the program code and the microprocessor
- Memory test
- Function test of the LED indicators and the acoustic signals
- Testing of the pacing and sensing capability
- Testing of the efficacy of high rate protection

During the self-test, all LED indicators light up and brief, acoustic signals can be heard. When the self-test is complete, the acoustic and visual signals turn off after a few seconds.

If the self-test finds a defect, all LED indicators flash permanently and acoustic warning signals sound.

Turn the external pacemaker off and, if necessary, repeat the self-test. If the self-test finds a defect again, return the device to BIOTRONIK.

5 Operating Elements and LEDs

	Name	Function
1a	Atrial channel connection	For patient cables and leads with 2-mm connectors or for the Redel adapter (red = plus; blue = minus)
1b	Ventricular channel connection	
4	Control dial Ventricle Sens.	Setting the ventricular sensitivity
6	Control dial Atrium Sens.	Setting the atrial sensitivity
7	Control dial AV delay	Setting the AV delay
8	Control dial Mode	Setting the pacing mode and turning off
9	Control dial Rate	Setting the pacing rate
11	Control dial Burst rate atrium	Setting the atrial burst rate
13	Control dial Atrium Ampl.	Setting the atrial pulse amplitude (cannot be used in VDD pacing mode)
15	Control dial Ventricle Ampl.	Setting the ventricular pulse amplitude
17	Battery magazine	For 9 V block battery

	Name	LED behavior
2	Low battery	LED indicator flashes red if the battery is almost depleted. Additional information can be found under: Changing the Battery [Page 24]
3	Ventricle Sense	LED indicator flashes green briefly for each sensed R wave
5	Atrium Sense	LED indicator flashes green briefly for each sensed P wave
10	Select burst	Selection of the atrial burst function
12	Start burst	Starts the atrial burst function
14	Atrium Pace	LED indicator flashes yellow briefly for each atrial event
16	Ventricle Pace	LED indicator flashes yellow briefly for each ventricular paced event

The LED indicators and acoustic signals also provide the following warnings during operation:

Note	Meaning	Error correction
Acoustic signal for 2 s	A pulse amplitude of < 1 V or a rate of > 180 bpm is set	Verify whether the set values are suitable for the patient.
Fast sequence of sounds	Impedance outside the permissible range	Verify that all connectors are plugged in securely. Verify that the leads are in the required position.
Acoustic signal and flashing of the Pace and Sense LEDs	High rate protection has been triggered, self-test failed	Turn off the device and return it to BIOTRONIK.

If the self-test did not find any error, the LED indicators and warning signals will cease, and the device will start to deliver pacing pulses in accordance with the set parameters. The negative electrode (cathode) should, therefore, only be connected when it is ensured that the pacing mode, pacing rate, pulse amplitude, and sensitivity have been set correctly. Setting the Mode rotary switch to OFF prevents pacing pulses from being delivered to the patient immediately after connecting the leads.

6 Operating the Reocor D

This device is robust, however, it must still be treated with care. Mechanical impact, such as dropping the device, can impair its function.

Attention


Technical failure due to condensation

If the external pacemaker has been transported at low temperatures, it must be stored for at least 2 hours under the ambient conditions given for operation.


Connecting the Redel Adapter

The external pacemaker is delivered with a Redel adapter, which is needed when using patient cables with Redel plugs, see Patient Cable with Redel Plug [Page 18].

1. Insert the connector pins of the Redel adapter into the Reocor connector ports.
2. Gently push the Redel adapter onto the Reocor.
3. Screw the Redel adapter to the Reocor on both sides.

	Name
1	Touch-proof 2-mm connector pins
2	Redel port
3	Screw connections

	Name
1	2-mm connector ports

Patient Cable with Redel Plug

The Redel adapter is needed for the following patient cables, as these patient cables have a Redel plug:

- PK-175 with 2 touch-proof screw terminals each for temporary leads on the patient side
- PK-83-B for single-chamber pacing with 2 touch-proof screw terminals for temporary leads on the patient side
- PK-83-C for single-chamber pacing with 2 touch-proof 2 mm connector ports for temporary catheters on the patient side
- PK-141 with 4 insulated alligator clips on the patient side
- PK-67-L
- PK-67-S

- The following adapters fit on the PK-67-L and PK-67-S patient cables:
- PA-1-B and PA-1-C for the connection of touch-proof 2-mm plugs or MHW adapters (adapters for heart wires)
 - PA-2 for connection to IS-1

The PK-67-L and PK-67-S patient cables are connected to the patient via alligator clips with the patient cable PK-155 or PK-155-B. The PK-155 and PK-155-B patient cables are intended for **single use** and cannot be resterilized.

Patient Cable Without Redel Plug


The patient cables below can be connected directly to the Reocor using the touch-proof 2-mm plugs:

- PK-82 patient cable with 2 insulated alligator clips for temporary leads on the patient side
- PK-83 patient cable for single-chamber pacing with 2 touch-proof screw terminals for temporary leads on the patient side

Setting Parameters

If the protective cover is closed, you cannot use the control dial and keys. Open the protective cover and remember to close it after setting the parameters.

- Push the release lever up with one hand.
- Slide the cover downwards.

All names in the descriptions refer to the drawing of the device, see Operating Elements and LEDs [Page 14]. The parameters are set using the control dials.

The sensitivities can be set from 1 mV to 20 mV with the Ventricle Sens. control dial and from 0.2 mV to 10 mV with the Atrium Sens. control dial.

The pulse amplitudes can be set from 0.1 V to 17 V with the Ventricle Ampl. and Atrium Ampl. control dials.

If a value less than 1 V is set, the device will sound a warning signal for 2 seconds. The pulse width is 1 ms.

The amplitudes and sensitivities should be checked regularly to ensure that correct sensing takes place and that a sufficient safety margin has been set.

The atrial rate can be set from 60 bpm to 1000 bpm with the Burst rate atrium control dial.

1. Use the control dial to set the rate.
2. Press the Select burst key.
3. After 2 seconds at the latest, press the Start burst key.

The pulse delivery then lasts as long as this key is pressed. The ventricular channel will continue to pace at the set rate, which can also be adjusted during this process. If an inhibiting pacing mode has been set, pacing will be inhibited in the ventricle.

The AV delay can be set from 15 ms to 400 ms with the AV delay control dial, but the AV delay is limited to half the basic interval. A short AV delay may be selected for special indications, e.g., in case of reentry tachycardias.

The mode can be set using the Mode control dial. For more information on the modes DDD, D00, VDD, VVI, V00, and VVT, see Pacing Modes [Page 22].

The external pacemaker is turned off using the Mode control dial, by turning it to the OFF position.

The rate can be set from 30 bpm to 250 bpm with the Rate control dial. If a value greater than 180 bpm is set, a warning signal will sound for 2 seconds.

Connecting the Leads

- Prepare the external pacemaker:
 - Set the pacing rate with the Rate control dial.
 - Set the AV delay with the AV delay control dial.
 - Set the pulse amplitudes for the atrium¹ and ventricle with the Atrium Ampl.¹ and Ventricle Ampl. control dials.
- Prepare the patient by placing the leads but not yet connecting them to the external pacemaker.
- Either connect the patient cable directly to the external pacemaker or connect the patient cable via the Redel adapter.
- Connect the lead to the patient cable.
- Set the mode with the Mode control dial and thus switch on the device, see Self-Test [Page 13].
- Verify the setting of the pulse amplitude.

The yellow Atrium Pace¹ and Ventricle Pace LEDs flash synchronously with the atrial and ventricular pacing pulses.

¹ Only for dual-chamber pacing.

Attention

Reocor D has 4 connector ports for direct connection of leads with the touch-proof 2-mm connectors. Temporary leads and heart wires can be connected with the PK-82 and PK-83 patient cables to the A+, A- and V+, V- connector ports.

Attention

Keep external devices ready

As a precaution, keep an external defibrillator and, for pacemaker-dependent patients, an emergency pacemaker ready.

WARNING

Connect the patient cable to the external device first before connecting intracardiac leads or heart wires. Power on the external device afterwards.

Caution

Monitor the heart rate

During use, monitor the patient's heart rate using an ECG monitor with alarm function.

Caution

Regular check of the connections

The connections of the external pacemaker and the temporary pacing leads must be secured and checked regularly.


Closing the Protective Cover

Caution

During use, the protective cover must be locked to prevent inadvertent adjustment of the control dials, causing a possible change in the program parameters.

1. Slide the transparent cover for the operating panel upwards until the release lever is stopped at the stop.
2. Push the release lever up with one hand.
3. Slide the cover upwards a little further until the release lever snaps into place behind the stop.

The protective cover can be removed completely for cleaning. To remove, push the cover down to the stop and then beyond the stop.

Securing

The external pacemaker should be operated lying on a horizontal, non-slip surface. It can also be attached to the patient using the armband, or operated hanging from an infusion stand using the hanger on the back of the device.

Caution

If an armband is used, the device or the armband must not be worn directly on the skin.

7 Pacing Modes

Reocor D operates in one of the following 6 modes:

DDD	Synchronous AV pacing with sensing and pacing in the atrium as well as in the ventricle
VDD	Synchronous, ventricular pacing with atrial tracking
D00	Asynchronous AV pacing without sensing
VVI	Sensing and pacing in the ventricle
V00	Asynchronous pacing in the ventricle
VVT	Like VVI, but immediate pulse delivery upon sensing of a ventricular event outside the refractory period

In case of disturbances caused by electromagnetic interference (EMI), Reocor D will switch to asynchronous pacing when certain limits are exceeded.
Depending on whether the interference was sensed in the atrium or in the ventricle, the following modes will result for the duration of the interference:

Interference-free mode	In case of interference by EMI
VDD	VAT, VVI, V00
DDD	DAD, DVI, D00

Cross channel blanking

After a stimulus has been delivered, sensing is suppressed for 80 ms ± 5 ms in the opposite channel to prevent far-field sensing.

Interference interval

The interference interval is used to protect against erroneous sensing caused by noise. Sensing events due to interference or noise are shown by short interval times. If they fall within the 80 ms interference interval, they are discarded and the interference interval is reset.
The result is asynchronous pacing at the programmed rate as long as the interference persists.
For instance, in DDD mode, an atrial interference that does not affect the ventricular channel leads to DVI pacing. Sensing of noise in the ventricular channel leads to DAD pacing. Interference in both channels results in D00 pacing.

Refractory periods

The total atrial refractory period (TARP) of Reocor D cannot be programmed, but is automatically derived from the sum of 175 ms and the set AV delay. Its minimum value is 400 ms below a pacing rate of 121 bpm.

Above this rate, its minimum value is 240 ms.

The atrial refractory period (ARP) is started by atrial sensed and paced events, as well as sensed premature ventricular contractions.

Reocor D's ventricular refractory period (VRP) is:

Pacing rate	Refractory period VRP
Below 150 bpm	225ms
150 bpm to 200 bpm	200ms
Above 200 bpm	175ms

If the upper rate is exceeded, every second atrial pulse will fall into the ARP, will not be sensed, and will not trigger a ventricular pulse. The ventricular pacing rate will continue at a ratio of 2:1.

8 Handling

Changing the Battery

When the Low battery LED flashes, it means that less than 36 hours of service time remain. The battery should be replaced immediately. A battery is replaced in the same way as it is inserted, see Inserting the Battery [Page 12].

For safety reasons, the patient should be paced by another source during the battery replacement.

Attention

If you use a battery other than the one supplied, do not use rechargeable batteries. The service time of these batteries is difficult to estimate.

Attention

Replacement during operation

Battery replacement during operation is possible. When the battery is removed, the external pacemaker remains ready for use at an ambient temperature of $20 \pm 2^{\circ}\text{C}$ for approx. 30 seconds.

Attention

Recommended batteries

The service time may vary depending on the battery used.
The specified service time is achieved with the following batteries:
Duracell MN1604 (6LR61)
Duracell Ultra Power MX 1604 (6LR61)
Energizer EN22 (6LR61)
Energizer 522 (6LR61)

Cleaning

To clean, use a cloth slightly moistened with water and if necessary, a mild cleanser can also be used. However, avoid bringing the device into direct contact with water or other solvents.

Attention

Strong cleaning agents or organic solvents (such as ether or gasoline) must not be used, as these can corrode the plastic housing.

Disinfection

For disinfection, wipe the device with a cloth soaked with a disinfectant solution containing alcohol (e.g., Aerodesin 2000). When mixing the solution, follow the dilution measure stated by the manufacturer.

Sterilization

The device cannot be sterilized. If the device needs to be used in a sterile environment, a sterile cover can be used.

Maintenance and Inspection

The only required maintenance action is replacement of the battery, see Inserting the Battery [Page 12].

An inspection should be performed once a year. Furthermore, after use in conjunction with high-frequency surgical instrument or cardiac defibrillators and if malfunctions are suspected.

The inspection must conform with the manufacturer's specifications, which are available upon request. They list all necessary test steps and the necessary equipment.

Reuse

The external pacemaker is intended for permanent use. After the end of use, it is not intended to reprocess and reuse the device.

Disposal

Used batteries must be treated as hazardous waste and disposed of by the user in an environmentally sound manner.

Patient cables that have been in contact with blood must be disposed of as contaminated medical waste.

Attention

Environmentally sound disposal

Reocor D contains materials that must be correctly disposed of in accordance with environmental protection regulations.

Old devices and uncontaminated accessories can be disposed of as electronic waste according to the applicable regulations or returned to BIOTRONIK.

BIOTRONIK ensures disposal in accordance with the national versions of the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE 2).

9 Appendix

Technical Data

General Data on Reacor

- Operating mode: continuous operation
- Service time: 500 hours at 70 bpm, 5 V, DDD mode, 500 ohms
- IP 31
- Protection class: internal power supply
- Type CF defibrillation protected applied part
- Operating temperature: +10°C to +40°C
- Storage and transport temperature: 0°C to +50°C
- Relative humidity for storage and transport 20% to 75% (non-condensing)
- Relative humidity during operation: 20% to 95% (non-condensing)
- Atmospheric pressure: 700 hPa ... 1060 hPa
- Information according to Section 33 of REACH, Regulation (EC) No. 1907/2006, see <https://www.biotronik.com/material-compliance>
- Reacor dimensions (WxHxD): 160 x 75 x 35 mm (without Redel adapter)
- Reacor weight:
 - 325 g with battery and with Redel adapter
 - 280 g without battery and with Redel adapter
 - 240 g without battery and without Redel adapter
- Redel adapter dimensions (WxHxD) 76 x 36 x 30 mm
- Redel adapter weight: 40 g

Parameters

Adjustable parameters	
Pacing mode (Mode)	OFF, DDD, D00, VDD, WI, V00, VVT
Basic rate (Rate)	30 ... 250 ppm At a rate of > 180 bpm a warning signal sounds.
AV delay (AV delay)	15 ... 400 ms
Burst rate (Burst rate atrium)	60 ... 1000 bpm
Pulse amplitude (Atrium Ampl.) and (Ventricle Ampl.)	0.1 ... 17 V At a pulse amplitude of < 1 V a warning signal sounds.
Sensitivity (Atrium Sens.)	0.2 ... 10 mV With respect to 15 ms sin ² pulse.
Sensitivity (Ventricle Sens.)	1 ... 20 mV With respect to 40 ms sin ² pulse.

Fixed parameters	
Pulse width	1 ms
Auto short after pace	< 20 ms
Interference interval	80 ms
Cross channel blanking	80 ms
Total atrial refractory period (TARP)	AVD +175 ms
TARP minimum	
30 ... 120 bpm	400 ms
121 ... 250 bpm	240 ms
Refractory period (V)	
30 ... 150 bpm	225 ms
151 ... 200 bpm	200 ms
201 ... 250 bpm	175 ms
Upper rate	260 bpm
High Rate Protection	
1 ... 180 bpm	286 ms
181 ... 250 bpm	286 ms = 210 bpm, does not apply for burst
	214 ms
	214 ms = 280 bpm, does not apply for burst
Pulse waveform	Asymmetric, biphasic

Battery data	
Battery	Alkaline manganese type: IEC 6LR61/ANSI 1604A 9 V, leakproof
Polarity	Cathodic
Reverse polarity protection	None: polarity is irrelevant
Power consumption	Typically 1 mA (70 bpm, 5 V, DDD mode, 500 ohms)
Service time with new battery	500 hours at 20°C At 70 bpm, 5 V, DDD mode, 500 ohms until ERI signal (EOS warning)
Remaining service time after ERI signal	36 hours At 70 bpm, 5 V, DDD mode, 500 ohms
End of service (EOS)	Flashing Low battery LED
Behavior during battery replacement	Device remains ready for use for at least 30 Seconds when the battery is removed. The mode set is retained.

Connection Variants

Temporary Catheters with 2-mm Connectors or Heart wire with 2-mm Adapter

You can connect Reocor D directly to a temporary catheter with touch-proof 2-mm plug or to a heart wire with 2-mm adapter, without any other patient cables or adapters. All other options are listed in the following table.

Patient-side connection	Patient cable	Device-side connection	Reocor D connection
Recommended connections			
Direct connection			2-mm connector port
2-mm adapter	PK-67-L with PA-1-B	Redel plug	Redel adapter
2-mm adapter	PK-83 with TC Adapt	2-mm plug	2-mm connector port
2-mm connector port	PK-83-C	Redel plug	Redel adapter
Possible connections			
Alligator clip	PK-141	Redel plug	Redel adapter
Alligator clip	PK-67-L with PK-155/PK-155-B (1x or 2x each) PK-67-S with PK-155/PK-155-B (1x or 2x each)	Redel plug	Redel adapter
Alligator clip	PK-82 (1x or 2x)	2-mm plug	2-mm connector port

Heart Wire with Break-Off Needle or with Flexible End

Patient-side connection	Patient cable	Device-side connection	Reocor D connection
Recommended connections			
Screw terminal	PK-175	Redel plug	Redel adapter
Screw terminal	PK-83	2-mm plug	2-mm connector port
Screw terminal	PK-83-B	2-mm plug	Redel adapter
Possible connections			
Alligator clip	PK-141	Redel plug	Redel adapter
Alligator clip	PK-67-L with PK-155/PK-155-B (1x or 2x each) PK-67-S with PK-155/PK-155-B (1x or 2x each)	Redel plug	Redel adapter
Alligator clip	PK-82 (1x or 2x)	2-mm plug	2-mm connector port

Implanted Lead with IS-1 Connector

Patient-side connection	Patient cable	Device-side connection	Reocor D connection
Recommended connections			
IS-1 connector port	PK-67-L with PA-2	Redel plug	Redel adapter
Possible connections			
Alligator clip	PK-141	Redel plug	Redel adapter
Alligator clip	PK-67-L with PK-155/PK-155-B (1x or 2x each) PK-67-S with PK-155/PK-155-B (1x or 2x each)	Redel plug	Redel adapter
Alligator clip	PK-82 (1x or 2x)	2-mm connector	2-mm connector port

Accessories and Spare Parts

Not all accessory products are available in every country.

Model	Order number	Description
Standard armband	103704	Standard armband
Short armband	391843	Short armband of smaller size suitable for small arms
BAT-1 9V	103688	Battery
Softcase	369603	Padded case for storing and carrying
Protective cover	378007	Protection against inadvertent adjustment of parameters
Redel adapter	371262	For connecting the PK-67-L, PK-67-S, PK-83-B, PK-83-C, PK-141, and PK-175 patient cables

Model	Order number	Description	Connection
PK-67-S (0.8 m)	128085	Patient cable for combining with the PK-155 and PK-155-B patient cables, PA-1-B, PA-1-C, and PA-2 patient adapters, can be resterilized	Redel adapter
PK-67-L (2.6)	123672		
PK-82 (1.5 m)	128564	Patient cable, with 2 insulated alligator clips, can be resterilized	Direct connection
PK-83 (1.5 m) PK-83 (2.5 m)	128562 128563	Patient cable, with 2 touch-proof screw terminals, can be resterilized	Direct connection
PK-83-B (1.5 m) PK-83-B (2.5 m)	347606 347485	Patient cable, with 2 touch-proof screw terminals, can be resterilized	Redel adapter
PK-83-C (2.5m)	478437	Patient cable, with 2 touch-proof 2-mm connector ports, can be resterilized	Redel adapter
PK-141 (2.8 m)	353181	Patient cables, each with 2 insulated alligator clips, can be resterilized	Redel adapter
PK-175 (2.3 m)	333959	Patient cables, each with 2 touch-proof screw terminals, can be resterilized	Redel adapter


Adapters for PK-67-L and PK-67-S

Not all accessory products are available in every country.

Item	Order no.	Description
PA-1-B	123751	For connection to 2-mm adapter or MHW adapter (adapter for heart wires), can be re-sterilized
PA-1-C	349723	For connection to 2-mm adapter or MHW adapter (adapter for heart wires), can be re-sterilized
PA-2	123157	For connection to IS-1 connector, can be re-sterilized
PK-155 (2 m) (set with 2 cables)	337358	Sterile patient cable, 2-wire, with alligator clips for single use
PK-155-B (2 m) (set with 2 cables)	404027	Sterile patient cable, 2-wire, with alligator clips for single use

Symbols

The label icons on the device symbolize the following:

Symbol	Meaning
	Comply with the technical manual (white image on blue background)
	Type CF applied part, defibrillation-protected
	Manufacturer
	Medical device
	Unique device identification
	Serial number
	Manufacturing date
	CE mark
	Protection against solid foreign objects with a diameter of 2.5 mm or more Protection against vertical dripping water
	Device contains materials that must be correctly disposed of in accordance with environmental protection regulations. The European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE 2) applies. Return devices that are no longer used to BIOTRONIK.

Legend for the Label

The label icons symbolize the following:

Symbol	Meaning
	BIOTRONIK order number
	Serial number
	Manufacturing date
	Unique device identification
	Reocor D
	Patient with implanted lead
	Contents
	Consult the instructions for use!
	Device contains materials that must be correctly disposed of in accordance with environmental protection regulations. The European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE 2) applies. Return devices that are no longer used to BIOTRONIK.
	Temperature limit
	Air pressure limit
	Humidity limit
	CE mark

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