Certificate

The SQS herewith attests that the organisation named below has a management system that meets the requirements of the normative basis mentioned.

Socorex ISBA SA
Chemin Champ-Colomb 7A
1024 Ecublens VD
Switzerland

Scope

**Development, production and distribution
of components and assemblies for laboratory
applications**

Normative base

ISO 9001:2015

Quality Management System

Reg. no. 10557

Validity 09.08.2021 – 19.07.2024
Issue 09.08.2021

A. Grisard, President SQS

F. Müller, CEO SQS

Swiss Association for Quality
and Management Systems (SQS)
Bernstrasse 103, 3052 Zollikofen, Switzerland

Certificate

The SQS herewith attests that the organisation named below has a management system that meets the requirements of the normative basis mentioned.

Socorex ISBA SA
Chemin Champ-Colomb 7A
1024 Ecublens VD
Switzerland

Scope

**Development, production and distribution
of components and assemblies for laboratory
and medical applications**

Normative base

**EN ISO 13485:2016 Medical devices –
Quality Management System**

Reg. no. 43173

Validity 09.08.2021 – 19.07.2024
Issue 09.08.2021

A. Grisard, President SQS

F. Müller, CEO SQS

Swiss Association for Quality
and Management Systems (SQS)
Bernstrasse 103, 3052 Zollikofen, Switzerland

Schweizerische Eidgenossenschaft
Confédération suisse
Confederazione Svizzera
Confederaziun svizra

Federal Department of Economic Affairs,
Education and Research EAER
State Secretariat for Economic Affairs SECO
Swiss Accreditation Service SAS

Swiss Confederation

Based on the Accreditation and Designation Ordinance dated 17 June 1996 and on the advice of the Federal Accreditation Commission, the Swiss Accreditation Service (SAS) grants to

SOCOREX ISBA S.A.
Laboratoire d'étalonnage
Chemin Champ-Colomb 7A
1024 Ecublens/Lausanne

Period of accreditation:
18.11.2020 until 17.11.2025
(1st accreditation: 18.11.2005)

the accreditation as

Calibration laboratory for volumetric instruments (liquid handling)

International standard: ISO/IEC 17025:2017
Swiss standard: SN EN ISO/IEC 17025:2018

3003 Berne, 18.11.2020
Swiss Accreditation Service SAS

Head of SAS
Konrad Flück

SAS is a signatory of the multilateral agreements of the European co-operation for Accreditation (EA) for the fields of testing, calibration, inspection and certification of management systems, certification of personnel and certification of products, processes and services, of the International Accreditation Forum (IAF) for the fields of certification of management systems and certification of products, processes and services and of the International Laboratory Accreditation Cooperation (ILAC) for the fields of testing and calibration.

To whom it may concern

CONFORMITY STATEMENT

We, SOCOREX ISBA SA, herewith certify that the products,

**Bottle-top dispensers Calibrex™ 525, 530,
and relevant accessories**

- laboratory dedicated instruments intended for safe and reliable liquid handling - are manufactured in our plant in Ecublens, Switzerland, under controlled conditions in accordance with ISO 9001 and ISO 13485 standards.

The manufacturing process includes functional checks and performance testing.

An individual calibration certificate, bearing the instrument's serial number and showing volumetric results, is supplied with each dispenser, to confirm the conformity with the ISO 8655 requirements.

SOCOREX ISBA S.A.

Jean-Marc Ammann
Quality System & Regulatory Affairs Manager

quality-system@socorex.com

Ecublens, November 21st, 2019

Chemical resistance of Socorex® dispensers

Calibrex™ models 525 / 530

Bottle-top dispensers are used daily for dispensing a wide range of chemicals. Therefore, instruments have to meet various requirements assuring safety of the laboratory staff and their work. Dispensers shall not release any substances which may interfere with trace analysis, have cytotoxic properties, distort optical tests or influence chromatographic methods and residue analysis.

Materials

Special attention is given to component materials (see charts below). All parts of the Calibrex™ dispensers in contact with the liquid are made of robust and chemically inert materials providing for long instrument life.

Parts	Calibrex™ 525	Calibrex™ 530
Feed tube	FEP	
Valve body	Ceramic - Aluminum oxide	
Valve balls	Ceramic - Aluminum oxide	
Valve springs	Platinum Iridium	
Plate	PTFE	
Barrel	Borosilicate glass	
Plunger	Ground Borosilicate glass	PFA coated glass
Body	ETFE	
Delivery jet	FEP / PCTFE	
Cap	ETFE	

Chemicals from A to Z

The following list includes most frequently used chemicals. It provides useful information for the safe and adequate use of Calibrex™ 525/530 dispensers. However, safety precautions and recommendations in operating instructions must be followed carefully.

Code explanations

- ++ = Good resistance
- + = Acceptable

Chemicals A - Z	Calibrex™ 525	Calibrex™ 530
A		
Acetaldehyde (Ethanal)	++	++
Acetic acid 96%	++	+
Acetic acid 100% (Glacial)	+	+
Acetic anhydride	+	+
Acetone (Propanone)	++	++
Acetonitrile (MECN)	+	+
Acetophenone	+	+
Acetyl Chloride	+	+
Acetylacetone	++	++
Acrylic acid	++	++
Acrylonitrile	+	+
Adipic acid		++
Allyl alcohol	++	++
Aluminum chloride		++
Amino acids		++
Ammonia <20%	+	+
Ammonia 20-30%	+	+
Ammonium chloride		++
Ammonium fluoride		++
Ammonium hydroxide	+	+
Ammonium molybdate		++
Ammonium sulfate		++
Amyl alcohol (Pentanol)	++	++
Amyl chloride (Chloropentane)	+	+
Aniline	++	++
Antimony trichloride		++
Ascorbic acid		++
n-Amyl acetate	+	+
B		
Barium chloride		++
Benzaldehyde	++	++
Benzene	+	+
Benzine	++	++
Benzoyl chloride	+	+
Benzyl alcohol	++	++
Benzyl chloride	+	+
Benzylamine	+	+
Bis(2-ethylhexyl) phthalate	+	+
Boric acid 10%	+	++
Bromine		
Bromobenzene	+	+
Bromonaphtalene	++	++
Butanediol	+	++
Butanol	++	++
Butanone (MEK)		
Butyl acetate	+	+
Butyl acrylate	+	+
Butyl methyl ether	+	+
Butylamine	+	+

Chemicals A - Z	Calibrex™ 525	Calibrex™ 530
B		
Butyric acid	+	+
C		
Calcium carbonate		
Calcium chloride		++
Calcium hydroxide		+
Calcium hypochlorite		+
Carbon disulfide	+	+
Carbon tetrachloride Thertracholomethane	+	+
Chlorine dioxide	+	+
Chloronaphthalene	+	+
Chloroacetaldehyde 45%	+	++
Chloroacetic acid	+	++
Chloroacetone	+	+
Chlorobenzene	+	+
Chlorobutane	+	+
Chloroethanol	+	+
Chloroform	+	+
Chloronitric acid 100%	+	
Chlorosulfuric acid	+	+
Chlorosulfuric acid 100%	+	+
Chromic acid 100%	+	+
Chromosulfuric acid 100%		+
Citric acid	+	++
Copper fluoride		+
Copper sulfate		++
Cresol	+	++
Cumene (Isopropylbenzene)	+	+
Cyanoacrylate		
Cyclohexane	+	+
Cyclohexanone	+	+
Cyclopentane	+	+
D		
1,2-Diethylbenzene	+	+
1,4-Dioxane (Diethylene dioxide)	+	+
1-Decanol	++	++
Decane	++	++
Di-(2-ethylhexyl) peroxydicarbonate	+	+
Dibenzyl ether	+	+
Dichloroacetic acid	++	++
Dichlorobenzene	++	++
Dichloroethane (DCE)	++	++
Dichloromethane (DCM)	+	+
Dichloroethylene	+	+
Diesel oil (Heating oil)	++	++
Diethanolamine	++	++
Diethylamine	+	+
Diethylene glycol	++	++
Diethyl ether	+	+

Chemicals A - Z	Calibrex™ 525	Calibrex™ 530
D		
Dimethyl sulfoxide (DMSO)	+	+
Dimethylaniline	++	++
Dimethylformamide (DMF)	+	+
Dimethylglycol / Dimethoxyethane (DME)	+	+
Dioxide chlorine	+	+
Diphenyl ether	+	+
E		
Essentials oils	+	+
Ethanol	++	++
Ethanolamine	+	+
Ether	+	+
Ethyl acetate	+	+
Ethylbenzene	+	+
Ethylene chloride	+	+
Ethylenediamine	++	++
Ethylene glycol	++	++
F		
Fluoroacetic acid	+	+
Formaldehyde (Formalin)	++	++
Formamide	++	++
Formic acid	++	++
G		
Gamma-butyrolactone	++	++
Gasoline	+	+
Glycerin <40%	++	++
Glycolic acid <50%	+	++
H		
Heating oil (Diesel oil)	++	++
Heptane	++	++
Hexane	++	++
Hexanoic acid	+	++
Hexanol	++	++
Hydriodic acid	+	+
Hydrobromic acid	++	++
Hydrochloric acid <20% (HCL) 10 to 100 mL	++	++
Hydrochloric acid <20% (HCL) 1 to 5 mL	+	+
Hydrochloric acid 20 to 37% (HCL) 10 to 100 mL	+	+
Hydrochloric acid 20 to 37% (HCL) 1 to 5 mL	+	+
Hydrofluoric acid (HF)		
Hydrogen peroxide	++	+
I		
Iodine		+
Iodine bromide		
Iodine chloride		

Chemicals A - Z	Calibrex™ 525	Calibrex™ 530
I		
Isoamyl alcohol	++	++
Isobutanol	++	++
Isooctane	++	++
Isopropanol	++	++
Isopropyl ether	+	+
Isopropylamine	+	+
K		
Kerosene	+	+
L		
Lactic acid		++
M		
2-Methoxyethanol	++	++
Methanol	++	++
Methoxybenzene (Anisol)	+	+
Methyl benzoate	+	+
Methyl chloride (Chloromethane)	+	+
Methyl ethyl ketone peroxide (MEKP)		+
Methyl formate	++	++
Methyl iodine (Iodomethane)	+	+
Methyl methacrylate (MMA)	+	+
Methyl n-butyl keton (MEK)		
Methyl propyl ketone (2-Pentanone)	++	++
Methyl tert-butyl ether (MTBE)	+	+
Methylene chloride (Dichloromethane) (DCM)	+	+
Methylpentanone	++	++
Mineral oil (engine oil)	++	++
N		
N-Butylamine	+	+
Nitric acid <30% - 10 to 100mL	++	++
Nitric acid <30% - 1 to 5mL	+	+
Nitric acid 30-70% - 10 to 100mL	+	+
Nitric acid 30-70% - 1 to 5mL		
Nitric acid >70% - 10 to 100mL		
Nitric acid >70% - 1 to 5mL		
Nitro-hydrochloric acid (Aqua regia)	+	+
Nitrobenzene	+	+
Nitromethane	+	+
N-methyl-2-pyrrolidone (NMP)	++	++
O		
Octane	++	++
Octanol	++	++
Oil (vegetable, animal)	+	+
Oil of turpentine	+	+
Oleic acid	+	++
Oxalic acid		++

Chemicals A - Z	Calibrex™ 525	Calibrex™ 530
P		
Pentane	+	+
Peracetic acid	++	++
Perchloric acid 100%	+	+
Perchloric acid diluted	++	++
Perchloroethylene	+	+
Petroleum	+	+
Petroleum ether / spirit	+	+
Phenol	++	++
Phenylethanol (2-phenylethanol)	+	+
Phenylhydrazine	+	+
Phosphoric acid <100%	++	++
Phosphoric acid <85%	++	++
Piperidine	+	+
Potassium chloride		++
Potassium dichromate		+
Potassium fluoride		
Potassium hydroxide		++
Potassium iodide		++
Potassium permanganate		+
Potassium peroxydisulfate (Potassium persulfate)		+
Potassium sulfate		+
Propionic acid (Propanoic acid)	++	++
Propylene glycol (Propane-1,2-diol)	++	++
Propylene oxide	++	++
Picric acid (Trinitrophenol)	+	+
Pyridine	+	+
Pyruvic acid	+	++
R		
Resorcin		++
S		
Salicylaldehyde	++	++
Scintillation fluid	++	++
Silver acetate		
Silver nitrate		++
Sodium acetate		++
Sodium chloride (Kitchen salt)		++
Sodium dichromate		++
Sodium fluoride		+
Sodium hydroxide 30%		+
Sodium hypochlorite		+
Sodium thiosulfate		++
Sulfonitric acid 100%	+	+
Sulfochromic acid 100%	+	+
Sulfur dioxide	+	+
Sulfuric acid <60% - 10 to 100mL	++	++
Sulfuric acid <60% - 1 to 5mL	+	+
Sulfuric acid >60% - 10 to 100mL	+	+
Sulfuric acid >60% - 1 to 5mL	+	

Chemicals A - Z	Calibrex™ 525	Calibrex™ 530
T		
Trichlorotrifluoroethane	+	+
Terebentine oil	++	++
Tartaric acid		++
Tetrachloroethane	+	+
Tetrachloroethylene / methylene	+	+
Tetrahydrofuran (THF)	+	+
Tetramethylammonium hydroxide		+
Tetramin	++	++
TKN Digest		+
Toluene	++	++
Trichlorethylene	+	+
Trichloroacetic acid	+	+
Trichlorobenzene	+	+
Trichloroethane / Methane	+	+
Trichloromethane (Chloroform)	+	+
Triethanolamine	++	++
Triethylene glycol	++	++
Trifluoroacetic anhydride (TFAA)	+	+
Trifluoroacetic acid (TFA)	+	+
Trifluoromethane (Fluoroform)	+	+
U		
Urea		++
X		
Xylene	+	+
Z		
Zinc chloride 10%		++
Zinc sulfate 10%		++

The above guidelines have been carefully reviewed prior to publication. Should you require information on chemicals not listed or contribute to some comments, please feel free to contact us.

CALIBREX™ organo 525

CALIBREX™ solutae 530

BOTTLE-TOP DISPENSERS


FLASCHEN-DISPENSER

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OPERATING INSTRUCTIONS

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MODE D'EMPLOI

SAFETY PRECAUTIONS

- Read operating instructions carefully, observe manufacturer's recommendations and chemical compatibility limits.
- Refer to and follow regulations about handling of potentially hazardous reagents.
- Always control proper working, tightness and bottle stability.
- Delivery jet should never point towards a person.
- Do not activate plunger when screw-locking stopper is attached and/or flow control stopcock on STOP position.
- Seized parts should not be separated by applying force.
- When carrying instrument, set volume on LOCK position; attach locking stopper. If available, set flow control stopcock on STOP position, hold bottle firmly.
- Prime and rinse dispenser carefully before disassembling any part or prior to storing.
- Take adequate measures to avoid electrostatic discharge when dosing flammable liquids.
- When using small reagent bottles or extension tubing, use workstation to secure dispenser stability.
- Operating temperature: 15°C - 40°C (60°F - 105°F).
- Instrument should only be used for its intended purpose.

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DESCRIPTION

The line includes robust dosing instruments with excellent chemical resistance, high volumetry performance and simplified maintenance. Volume range stretches between 0.1 and 100 mL. Each instrument fits most laboratory bottles either directly or using adapters supplied. Flow control stopcock, when available, enables liquid priming and recycling without reagent loss or contamination.

Calibrex™ organo 525 includes ground glass or ceramic plunger, best suited for organics, non-crystallizing acid and base solutions.

Calibrex™ solutae 530 has a PFA coated plunger adequate for weak, strong acids, bases, salt and/or crystallizing solutions.

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OPERATION

Installation

Screw by hand delivery jet adapter nut No. 15 (Fig. 1.1) without forcing. Optional longer delivery jet available, see page 7.

Cut feed tube No.19 according to bottle size (Fig. 1.2) and insert into body connection (Fig. 1.3). Install instrument on bottle using appropriate adapter if needed. Bottle tightness is guaranteed when neck is properly sitting in adapter. Optional telescopic tubing is available on request, see page 7.

Notes: overtightening dispenser on bottle reduces instrument rotation. Improve dispenser stability by using appropriate workstation.

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Delivery jet position (Fig. 2.1)

Rotation possible over 360°. Hold connecting body only (and not delivery jet) to reach optimal position.


Note: flexible connection of delivery jet absorbs shocks and reduces risks of breakage (Fig. 2.2).

Extended delivery jet

Longer delivery jet installation.

- ① Screw locking stopper (No. 17)
- ② Unscrew delivery jet adapter (No.15)
- ③ Replace with extended delivery jet (Fig. 1.1)

Note: using stopcock adds another 20 mm to overall delivery jet size.

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

DOSING INSTRUCTIONS

Volume setting

Press with thumb and move sliding cursor (No. 8) up and down to desired volume (Fig. 3). Release to lock and prevents unwanted change. The dual graduation scale enables accurate reading from a wide angle.

Screw-type button (No. 9) available as accessory. If preferred, change button according to Fig. 9.1 and 9.2. When using screw type button, carefully align index to desired volume scale before locking. Screw tightly to avoid unwanted changes (Fig. 4).

Priming


Note: before first use, carefully rinse instrument.

Model without stopcock

Remove stopper No.17 and slowly activate dispenser several times over short distance (≈ 30 mm) to eliminate any air bubble in barrel chamber and delivery jet. Avoid hitting upper and lower stops. Collect liquid in a clean vessel for recycling (Fig. 5.1).

Model with flow control stopcock

Place stopcock selector in recycle position ➡ (Fig. 5.2) and smoothly activate over a short distance (≈ 30 mm) to eliminate any air bubble in barrel chamber. The liquid is returned to bottle. Avoid hitting upper and lower stops. Thereafter, select distribution position ▶ to start dispensing.

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
Dispensing

Remove locking stopper (No.17). Raising plunger head fills the dispenser barrel. Depress to dispense liquid. Best results are obtained by moving the plunger gently and regularly between upper and lower stops.

Liquid density or viscosity may change dosing speed and instrument performance. Keep working without applying force.

Warning: never try to distribute liquid if stopper is screwed on delivery jet or if stopcock selector is in STOP position. Do not leave instrument unattended with liquid in barrel.

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End of dosing and storing

Purge and rinse instrument carefully. Replace stopper (No. 17) on delivery jet. If available, move stopcock selector to STOP position to prevent liquid leakage (Fig.5.3).


Fig. 5.2

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
Fig. 6

MAINTENANCE

Regular maintenance substantially contributes to optimal performance and long life span. Rinse for instance when bottle is emptied or before distribution of another liquid.

Disassembling (Fig. 6)

Purge properly before removing instrument from the bottle. Pull out feed tube No. 19.

- ① Unscrew plunger head (No. 4) and remove plunger. Never apply force to avoid breakage.
- ② Unscrew connecting body nut to access glass barrel, valve body (No. 13) and plate (No. 11 – not available on 1 and 2.5 mL models). Do not misplace valve ball (No. 12).
- ③ Unscrew delivery jet nut (No. 15) and separate from dispenser.

Cleaning

Soak parts in a mild laboratory detergent. Take particular care of plunger, glass barrel and valve body. Remove dirt from valve to avoid instrument blocking. Eliminate residues on valve seat or ball with ultra-sonic bath to prevent tightness problem. Rinse with water.

Reassembling

Damaged parts must be replaced by original spare parts only. To reassemble, proceed step by step in the disassembly reverse order. The exploded view on page 8 helps positioning each part. Make sure plunger and all connections are correctly tightened. Check for leakfree aspiration and distribution functions.

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TROUBLE SHOOTING

Observation	Possible cause	Action
Liquid aspiration impaired	Plunger blocked by crystallization or dried residues.	Never apply force to resolve. Soak dispenser in warm water or mild laboratory detergent. Move parts gently.
	Aspiration valve ball blocked or missing.	Dismantle to access valve body. Use warm water or ultrasonic bath to free valve ball. Replace ball if missing.
	Feed tubing blocked.	Remove and clean feed tubing.
Erratic plunger movement	PFA coating damaged due to lack or improper cleaning.	Replace plunger. Follow cleaning instructions.
	PFA coating swolled by inappropriate chemicals.	Replace plunger. Check chemical compatibility.
	Dirty surface.	Clean plunger.
	Distortion / damage on parts resulting from improper autoclaving.	Replace parts. Respect autoclaving conditions.
Air bubbles in system	Dispenser not properly primed.	Prime according to "Dosing instructions" paragraph.
	Feed tubing improperly connected or damaged.	Connect properly or change.
	Insufficient liquid level in bottle or inadequate tubing length.	Refill bottle, replace feed tubing.
	Connecting body nut is loose	Tighten connecting body nut
No liquid distribution	Stopcock set on wrong position.	Check for correct stopcock position.
	Stopper locking delivery jet.	Remove locking stopper.
	Dispensing valve ball blocked.	Dismantle to access valve body. Use warm water or ultrasonic bath to free valve ball.
Liquid leakage	Loose connecting screws.	Tighten all connexions without excessive force.
	Glass barrel cracked or broken.	Replace part.
	Stopcock set on STOP position.	Check for correct stopcock position.
Unstable volume setting	Sliding cursor assembled upside down.	Follow assembly instructions in "Volume setting button" paragraph.
	Screw type button loose.	Tightly screw button.
Unstable performance	Uneven dosing movements.	Move plunger gently and regularly between upper and lower stops without shocks.
	Valves untight due to residues on ball or seat.	Clean valves.
	Liquid leakage.	See above.
	Liquid viscosity, density or temperature different from original calibration conditions.	Perform new calibration with liquid used or at actual temperature.

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
Fig. 8

Calibration (Fig. 7.1 and 7.2)

Instrument is factory calibrated and delivered with an individual QC certificate. A volumetric control is recommended at least after every 12 month of use.

Should performance no longer correspond to specifications, or if conditions of use (liquid density, temperature, etc.) require recalibration, this can be easily performed using an analytical balance.

A new calibration is also recommended when replacing either barrel/plunger assembly (No. 5), main body (No. 7) or valve body (No. 13).

- ① Remove calibration safety seal (No. 1)
- ② Remove cap (No. 2)
- ③ Lift calibration key (No. 3)
- ④ Adjust according to the +/- arrows engraved on calibration key. Each notch corresponds to a volume correction of:

Range	Correction
0.1 - 1 mL	± 0.9 µL
0.25 - 2.5 mL	± 2.3 µL
0.5 - 5 mL	± 4.6 µL
1 - 10 mL	± 9.3 µL
2.5 - 25 mL	± 17.4 µL
5 - 50 mL	± 34.7 µL
10 - 100 mL	± 52.1 µL

Replace cap (No. 2) and apply a new calibration safety seal (No. 1).

Sterilization (Fig. 6)

Instrument designed for sterilization in the autoclave at 121°C / 250°F (20 min), fully assembled. Instrument must be purged and rinsed prior autoclaving. Delivery jet parts (Nos. 15, 16, 18) must be properly mounted.

Slightly untighten connecting body nut ②, delivery jet (or stopcock) nut ③, and stopper ④. Set volume on minimal value. Place instrument horizontally in the autoclave, avoiding any direct contact with metal.

Let instrument dry and cool down to room temperature (min 2 hours) before securely re-tightening all connections (without excessive force).

Verify inaccuracy after first few autoclaving cycles to guarantee best instrument performance. Otherwise inaccuracy value may drift to higher values, particularly on small volumes.

Notes: repeated autoclaving may affect some parts and material colors. Correct autoclaving and resulting sterility are the responsibility of the user.

Materials

Parts in contact with the liquid are chemically inert:

- Feed tube: PTFE (up to 10 mL) FEP (as of 25 mL)
- Valve body: Ceramic
- Valve ball: glass (up to 10 mL) or ceramic (as of 25 mL)
- Valve spring: Platinum-iridium
- Plate: PTFE (as of 5 mL)
- Barrel: Borosilicate glass
- Plunger: Glass or ceramic (Calibrex™ 525), PFA coated ceramic or glass (Calibrex™ 530)
- Connecting body: ETFE
- Delivery jet tubing: FEP/PTFE
- Stopper: ETFE

Warning: Vapors of concentrated strong acids and some solvents may damage your instrument after some time and reduce its life span. It is recommended, not to leave the dispenser permanently on a bottle containing concentrated strong acids and to rinse it after use. To prevent long term loss of product by evaporation, the same recommendation is valid when using volatile solvents. Damages due to incompatible reagents are not covered by the warranty. The QR code printed on the instrument provides access to the latest chemical compatibility chart. In case of doubt, check for chemical compatibility or refer to manufacturer before use.

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USE OF ACCESSORIES

Note: exclusive use of Socorex® accessories.
Verify performances after the addition of any accessories.
purge dispenser prior to any change of accessory.

Sterile air filter (Fig. 8)

Air inlet can be enlarged to receive a membrane filter, if airborne contamination is a concern.

- ① Cut out air inlet protection cover using sharp blade
- ② Insert Luer side of filter in the opening

Remove filter prior to autoclave instrument

Volume setting button (Fig. 9.1 and 9.2)

To replace sliding cursor with screw type button:

- ① Place cursor to the bottom
- ② Unscrew plunger head (No. 4)
- ③ Remove plunger
- ④ Push body in the upper position. Remove upper apron (No. 6) by rocking toward the front side.
- ⑤ Slide cursor with thumb toward the top and remove it from the body.
- ⑥ Insert screw type button, slide it toward the bottom, and tighten the cursor to the desired volume.

Replace upper apron (No. 6), respecting the index position until it "clicks" (Fig. 9.3). Insert plunger completely and securely tighten.

Note: upon reinstalling sliding cursor, make sure the arrow is pointing up (Fig. 9.4).


Fig. 9.1


Fig. 9.2


Fig. 9.3


Fig. 9.4

Flow control stopcock (Fig. 10)

Accessory can be mounted and removed at anytime.

- ① Remove delivery jet (Fig. 10.1), then delivery jet adapter (Fig. 10.2) and keep for further use
- ② Mount stopcock on dispenser (Fig. 10.3)
- ③ Attach delivery jet to stopcock (Fig. 10.4)

Note: see "Dosing instructions" (Fig. 5.2 and 5.3) for proper use of flow control stopcock.


Fig. 10.1


Fig. 10.2


Fig. 10.3


Fig. 10.4

Extension tubing (Fig. 11 and 12)

Spiral tubing and Jet pen™ can replace standard delivery jet. Fully compatible with stopcock features.

- ① Remove delivery jet (Fig. 11.1 and 12.1)
 - ② Install extension tubing (Fig. 11.2 and 12.2)
 - ③ Fix Jetpen™ holder (Fig. 12.3)
- If required, empty Jet pen™ holder by unscrewing the reservoir holder (cap is maintained on the holder)

Materials in contact with the liquid:

- Spiral tubing: FEP
- Adapter: PCTFE
- Jet pen™ holder: PP coated glass, ETFE

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Remote aspiration kit (Fig.13)

Kit facilitates liquid intake from a drum or an other remote container (<10 m distance, <2 m elevation). It requires workstation # 320BC.100 and feed tubing set # 1.525.581 (to be ordered separately)

- ① Position ring spacer and fit with adapter (Fig. 13.1)
- ② Attach feed tubing
- ③ Place dispenser on workstation and screw the two clamps. (Fig. 13.2)

Note: accessory not compatible with recycling position of stopcock.

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
Fig. 13.1


Fig. 13.2

ORDERING INFORMATION - ACCESSORIES

BESTELLINFORMATIONEN - ZUBEHÖR

CATALOGUE - ACCESSOIRES

Description	Packaging	Cat. No.
Work station for dispenser		
Work station for dispenser stability	1 / pk	320.SB100
Remote aspiration work station and feed tubing set		
Work station for remote aspiration*	1 / pk	320.BC100
Feed tubing set, 2.5 m with connector (as of 25 mL)	1 / pk	1.525.581
* Dispenser and feed tubing / connector to be ordered separately		
Flow control stopcock		
Up to 10 mL models	1 / pk	1.525.544
As of 25 mL models	1 / pk	1.525.546
Delivery jet assembly, 90 mm		
Fits 1 mL model	1 / pk	1.525.090
Fits 2.5, 5 and 10 mL models	1 / pk	1.525.091
Delivery jet assembly, 120 mm		
Fits 25 mL model	1 / pk	1.525.120
Fits 50, 100 mL models	1 / pk	1.525.121
Delivery jet assembly, extended, 120 mm		
Fits 1 mL model	1 / pk	1.525.123
Fits 2.5, 5 and 10 mL models	1 / pk	1.525.125
Delivery jet assembly, extended, 150 mm		
Fits 25 mL model	1 / pk	1.525.150
Fits 50 and 100 mL models	1 / pk	1.525.151
Safety seal sticker		
Safety seal sticker	25 / pk	1.525.525
Safety seal sticker	100 / pk	1.525.526

Description	Packaging	Cat. No.
Extension tubing and Jet-Pen™ 600mm, int. ø 4mm		
Fits models up to 10 mL	1 / pk	1.525.610
Fits 25 mL model	1 / pk	1.525.625
Fits 50 and 100 mL models	1 / pk	1.525.650
Replacement feed tubing, int. ø 5mm		
Fits models up to 10 mL, cut, 300mm, PTFE	1/pack	511.707
Fits models up to 10 mL, by the meter, uncut, PTFE	meter	511.709
Fits models as of 25 mL, cut, 350mm, PEP	1/pack	525.350
Fits models as of 25 mL, by the meter, uncut, FEP	meter	525.706
Telescopic feed tubing, FEP		
Fits models up to 10 mL, 150 – 255 mm, ø ext 6.5 mm, ø int 5 mm	1 / pk	1.525.352
Fits models as of 25 mL, 195 – 345 mm, ø ext 8.5 mm, ø int 7 mm	1 / pk	1.530.355

Volume adjustment button		
Screw type setting button for organo model, yellow	1 / pk	1.525.918
Screw type setting button for solutae model, red	1 / pk	1.530.918

Adapter and bottles		
Threaded adapter (PP) ø 25, 28, 32, 38, 40 mm	1/pack	1.525.GLxx
Nested threaded adapter (PP) to be used with 1.525.GL32 ø 22, 30, 34, 36	1/pack	GLPxx
Amber glass reservoir (ø 32/45 mm) 250, 500, 1000, 2500 mL	1/pack	314.xxxx
Amber glass reservoir (ø 32/45 mm) 500, 1000, 2500 mL with anti-spill PE coating	1/pack	314.xxxxPE
Pyrex glass reservoir (ø 32 mm) 500 mL with side neck for filter	1/pack	314.0500P
Polyethylene reservoir (ø 32/45 mm) 250, 500, 1000, 2500 mL	1/pack	315.xxxx

Replace xx by adapter diameter and xxxx by reservoir volume

WARRANTY

Dispensers are guaranteed for a period of two years against any factory or material defect. Glass breakage, cases due to non respect of manufacturer's and safety instructions and dam-ages due to incompatible reagents are excluded from the warranty. Repairs and changing parts do not extend the warranty period.

Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased, after obtaining return authorization.

Purge, dry and decontaminate the instrument prior to returning it for service.

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SUPPLY CONTENT

Instrument with 45 mm base thread. Supplied with feed tubing, delivery jet with stopper, three bottle neck adapters, QC certificate and operating instructions.

Up to 2.5 mL: 300 mm feed tubing, 90 mm delivery jet, 25, 28 and 32 mm adapters.
From 5 mL to 10 mL: 300 mm feed tubing, 90 mm delivery jet, 28, 32 and 40 mm adapters.
As of 25 mL: 350 mm feed tubing, 120 mm delivery jet, 32, 38 and 40 mm adapters

PERFORMANCE DATA
AND ORDERING INFORMATION

Volume mL	Increment mL	Inaccuracy (E %)			Imprecision (CV %)			Calibrex organo 525		Calibrex solutae 530	
		Min. vol.	Mid. vol.	Max. vol.	Min. vol.	Mid. vol.	Max. vol.	Cat. No.	With stopcock Cat. No.	Cat. No.	With stopcock Cat. No.
0.1 – 1	0.02	<+/- 3.0 %	<+/- 1.8 %	<+/- 0.6 %	< 1.2 %	< 0.7 %	< 0.17 %	525.001	525.001FC	530.001	530.001FC
0.25 – 2.5	0.05	<+/- 2.7 %	<+/- 1.6 %	<+/- 0.6 %	< 0.9 %	< 0.55 %	< 0.17 %	525.002.5	525.002.5FC	530.002.5	530.002.5FC
0.5 – 5	0.1	<+/- 2.0 %	<+/- 1.3 %	<+/- 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %	525.005	525.005FC	530.005	530.005FC
1 – 10	0.2	<+/- 1.5 %	<+/- 1.2 %	<+/- 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %	525.010	525.010FC	530.010	530.101FC
2.5 – 25	0.5	<+/- 1.5 %	<+/- 1.1 %	<+/- 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %	525.025	525.025FC	530.025	530.025FC
5 – 50	1.0	<+/- 1.5 %	<+/- 1.1 %	<+/- 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %	525.050	525.050FC	530.050	530.050FC
10 – 100	1.0	<+/- 1.5 %	<+/- 1.1 %	<+/- 0.6 %	< 0.5 %	< 0.35 %	< 0.1 %	525.100	525.100FC	530.100	530.100FC

- Performance values obtained by a smooth and steady movement, with bidest. water at constant temperature (±0.5°C) comprised between 20 and 25°C, in accordance with ISO 8655.
- Refer to chapter "Sterilization" when autoclaving the instrument
- Omission to untighten connecting body ring before autoclaving, and/or over-tightening ring when dispensing, may reduce performance.

Vol. (mL)	Max. cinematic viscosity (cst)	Example of liquid
1	1500	Silicone oil
2.5	1200	Silicone oil
5	1000	Pure Glycerol
10	250	SAE 20 oil

Cinematic viscosity is indicated in centistokes (cst). The dynamic viscosity in centipoises (cps) is obtained by multiplying the cinematic viscosity by the density of the liquid (kg/m³).

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
Vol. (mL)	Max. cinematic viscosity (cst)	Example of liquid
25	400	Silicone oil
50	250	Silicone oil
100	150	Silicone oil

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
Test conditions for all models (without stopcock)

Dispensing time < 10s.

Dispensing force < 3kg for 5 to 100mL,
<0.5kg for 1 and 2.5mL

Pressure in barrel < 3bar

Liste de pieces

1. Pastille de protection

2. Capuchon

3. Clé de calibration

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5. Piston/cylindre sous-ensemble

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7. Manchon

8. Curseur cranté

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11. Platine (sauf 1 et 2.5 ml)

12. Bille de soupape

13. Corps de soupape

14. Embase

15. Raccord de bec

16. Bec verseur

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18. Robinet (accessoire)

19. Tuyau d'aspiration

Chemical resistance
Scan QR code to access
chemical resistance chart.

Calibrex™

your choice for flexibility and performance

bottle-top dispensers

*additional
volumes*

The Calibrex™ line includes a large selection of robust dosing instruments with excellent chemical resistance, high performance and simplified maintenance. Intended for the safe and reproducible liquid distribution in volumes ranging from 0.1 mL to 100 mL.

Choose from three different models ; whatever the application, there will be a Calibrex™ fitting your needs.

The Calibrex™ line features:

- Selection between three different models
- Colour coding identification
- Integrated calibration key
- Long lasting performance stability
- QR coded chemical compatibility
- Quick disassembling, no tool needed
- Autoclavable at 121° C fully assembled

The Calibrex™ models

520 **525** **530**
universal organo solutae

 **SOCOREX**
SWISS


Calibrex™ *organo / solutae*

your choice for flexibility and performance

► Ease of use ①

Smooth aspiration and distribution for trouble-free routine dispensing. Dual scale enables volume reading from a wide angle.

► Volume setting

Spring loaded sliding cursor ② softly moves up and down dual scale and precisely stops at desired graduation. Easy exchange for the alternative classical screw button ③, available as accessory.

► Permanent visibility ④

Safety transparent thick sleeve around barrel and delivery jet window allow liquid flow monitoring.

► Optimal working position ⑤

The instrument rotates 360° for adequate positioning on the bottle. Freely rotating dispenser body makes volume clearly visible from any position.

► Robustness

Glass barrel with thick wall contributes to overall dispenser robustness.

Calibrex™ *organo*

Calibrex™ *organo* 525 includes a ground glass or ceramic plunger, both suited for organics and non-crystallizing acid and base solutions.

525

organo

0.1 – 1 mL
0.25 – 2.5 mL
0.5 – 5 mL
1 – 10 mL
2.5 – 25 mL
5 – 50 mL
10 – 100 mL

Calibration safety seal sticker

Integrated calibration key under cap

Dual graduation

Volume adjustment sliding cursor

Ground glass plunger

Disassembling nut, no tool needed

Air inlet for filter connection

Shock absorbing flexible connexion

Screw-locking stopper

Superior material selection

Parts in contact with liquid flow are chemically inert, providing for stability and long instrument life.

Parts	525 <i>organo</i>	530 <i>solutae</i>
Feed tubing 1 to 10 mL	PTFE	
Feed tubing 25 to 100 mL	FEP	
Valve body	Ceramic	
Valve balls 1 to 10 mL	Pyrex glass	
Valve balls 25 to 100 mL	Ceramic	
Valve spring	Platinum-iridium	
Valve plate	PTFE	
Barrel	Borosilicate glass	
Plunger 1 to 5 mL	Ceramic	PFA coated ceramic
Plunger 10 to 100 mL	Ground glass	PFA coated glass
Connecting body	ETFE	
Delivery jet	FEP/PCTFE	
Stopper	ETFE	


Chemical resistance ⑦

Printed QR code for instant access to chemical resistance chart.

Calibrex™ solutae


Calibrex™ solutae 530 has a glass or ceramic PFA coated plunger preventing the crystallization of chemicals. It enables trouble free distribution of salt solutions, weak and strong acids, as well as bases.

530


solutae

0.1 – 1 mL
0.25 – 2.5 mL
0.5 – 5 mL
1 – 10 mL
2.5 – 25 mL
5 – 50 mL
10 – 100 mL

► Flow control stopcock ⑧

Both 525 and 530 models are available with or without stopcock, adding flexibility to safety. Next to its dispensing position (b), device enables liquid priming and recycling (a) without reagent loss or contamination. Locking position (c) for safe dispenser transportation.

► Easy In-lab calibration ⑨

Access to mechanism protected by safety seal sticker. Integrated key located under plunger cap. Engraved +/- arrows for easy and precise setting. Spare seal stickers are available.

► Maintenance ⑩

Disassembling/reassembling facilitated by limited number of elements – no tool needed. Key parts such as valves, plunger, barrel and delivery jet are removable in seconds for cleaning. Fully autoclavable at 121°C / 250°F.

► Air filter ⑪

Air inlet can be enlarged to receive a Luer membrane filter, if airborne contamination is a concern.

Performance – Calibrex™ 525 and 530

Volume mL	Division mL	Inaccuracy (E%)			Imprecision (CV%)		
		Min. vol.	Mid. vol.	Max. vol.	Min. vol.	Mid. vol.	Max. vol.
0.1 – 1	0.02	<+/- 3.0%	<+/- 1.8%	<+/- 0.6%	< 1.2%	< 0.7%	< 0.17%
0.25 – 2.5	0.05	<+/- 2.7%	<+/- 1.6%	<+/- 0.6%	< 0.9%	< 0.55%	< 0.17%
0.5 – 5	0.1	<+/- 2.0%	<+/- 1.3%	<+/- 0.6%	< 0.5%	< 0.35%	< 0.1%
1 – 10	0.2	<+/- 1.5%	<+/- 1.2%	<+/- 0.6%	< 0.5%	< 0.35%	< 0.1%
2.5 – 25	0.5	<+/- 1.5%	<+/- 1.1%	<+/- 0.6%	< 0.5%	< 0.35%	< 0.1%
5 – 50	1.0	<+/- 1.5%	<+/- 1.1%	<+/- 0.6%	< 0.5%	< 0.35%	< 0.1%
10 – 100	1.0	<+/- 1.5%	<+/- 1.1%	<+/- 0.6%	< 0.5%	< 0.35%	< 0.1%

Performance values obtained by a smooth and steady pace movement, with bidest water at constant temperature ($\pm 0.5^\circ\text{C}$) comprised between 20 and 25°C , according to EN ISO 8655.

Warning: omission to untighten connecting body ring before autoclaving, and/or over-tightening ring when dispensing, may reduce performance.

Calibrex™ universal

your choice for flexibility and performance

► Friendly volume setting ①

Fast and precise adjustment. Efficient click-stop mechanism prevents unwanted alteration. Large display is easy to read; the window adjusts to desired side of instrument body. Each step on the cylindrical cam is pre-calibrated and corresponds to one division on the volume graduation.

► Optimized working position ②

The instrument rotates 360° for adequate positioning on top of the bottle.

► Air filter ③

Air inlet can be enlarged to receive a Luer membrane filter, if airborne contamination is a concern.

► In-lab calibration ④

Dispensers are factory calibrated and can be easily recalibrated. Integrated adjustment screw bears clear setting indications.

► Easy maintenance ⑤

Disassembling/reassembling is facilitated by the limited number of elements – no tool needed. Fully autoclavable at 121°C / 250°F.


Calibrex™ universal

Calibrex™ universal 520 has a PFA coated plunger preventing the crystallization of chemicals. Made of high-tech materials, it provides for broad chemical resistance and is universally adapted for most laboratory reagents.

520

universal

0.25 – 2 mL
1 – 5 mL
1 – 10 mL

Superior material selection

Special attention is paid to component materials, providing for long instrument life. Parts coming in contact with the liquid flow are chemically inert.

Chemical resistance

Printed QR code for instant access to chemical resistance chart.

Parts	Materials
Feed tube	PTFE
Intake valve	Ceramic
Valve balls	Pyrex glass
Valve spring	Platinum-iridium
Barrel	Borosilicate glass
Barrel plate / base	PTFE
Plunger	Glass with PFA coating
Outlet valve	Ceramic
Body	ETFE
Delivery jet assembly	PTFE/ETFE

Performance – Calibrex™ 520

Volume mL	Division mL	Inaccuracy (E%)			Imprecision (CV%)		
		Min. vol.	Mid. vol.	Max. vol.	Min. vol.	Mid. vol.	Max. vol.
0.25 – 2	0.05	< ± 3.0%	< ± 1.8%	< ± 0.6%	< 0.5%	< 0.35%	< 0.1%
1 – 5	0.1	< ± 2.0%	< ± 1.3%	< ± 0.6%	< 0.5%	< 0.35%	< 0.1%
1 – 10	0.25	< ± 1.5%	< ± 1.1%	< ± 0.6%	< 0.5%	< 0.35%	< 0.1%

Performance values obtained with bidest water at constant temperature ($\pm 0.5^\circ\text{C}$) comprised between 20 and 25°C in accordance with ISO 8655.

Accessories for Calibrex™ dispensers

Ordering information - Accessories

Description	Packaging	Cat. No.
Work station for dispenser stability ①		
Fits Calibrex™ 520 dispenser	1/pk	320.SB050
Fits Calibrex™ 525/530 dispensers	1/pk	320.SB100

Remote aspiration work station and feed tubing set ②		
Fits Calibrex™ 520 dispenser*	1/pk	320.BC050
Fits Calibrex™ 525/530 dispensers* as of 25 mL	1/pk	320.BC100
Feed tubing set, 2.5 m with connector for remote aspiration, fits Calibrex™ 525/530 as of 25 mL	1/pk	1.525.581

* Dispenser and feed tubing / connector to be ordered separately

Description	Material	Length	int. Ø	Cat. No.
Extension tubing and Jet-Pen™ ③				
Fits Calibrex™ 520	PTFE	600 mm	2.2 mm	1.524
Fits Calibrex™ 525/530, up to 10 mL	FEP/PCTFE	600 mm	4 mm	1.525.610
Fits Calibrex™ 525/530, 25 mL	FEP/PCTFE	600 mm	4 mm	1.525.625
Fits Calibrex™ 525/530, 50 and 100 mL	FEP/PCTFE	600 mm	4 mm	1.525.650

Description	Packaging	Cat. No.
Flow control stopcock for Calibrex™ 525/530 ④		
Fits models up to 10 mL	1/pk	1.525.544
Fits models as of 25 mL	1/pk	1.525.546

Delivery jet assembly, 90 mm, for Calibrex™ 525/530 ⑤		
Fits 1 mL model	1/pk	1.525.090
Fits 2.5, 5 and 10 mL models	1/pk	1.525.091

Delivery jet assembly, 120 mm, for Calibrex™ 525/530 ⑤		
Fits 25 mL model	1/pk	1.525.120
Fits 50 and 100 mL models	1/pk	1.525.121

Delivery jet assembly, extended, 120 mm, for Calibrex™ 525/530 ⑤		
Fits 1 mL models	1/pk	1.525.123
Fits 2.5, 5 and 10 mL models	1/pk	1.525.125

Delivery jet assembly, extended, 150 mm, for Calibrex™ 525/530 ⑤		
Fits 25 mL models	1/pk	1.525.150
Fits 50 and 100 mL models	1/pk	1.525.151

Description	Colour	Packaging	Cat. No.
Screw type volume setting button ⑥			
Fits Calibrex™ 525 dispensers	Yellow	1/pk	1.525.918
Fits Calibrex™ 530 dispensers	Red	1/pk	1.530.918

Description	Length	int. Ø	Cat. No.
Replacement feed tubing, PTFE ⑦			
Fits Calibrex™ 520/525/530 up to 10 mL, cut	300 mm	5 mm	511.707
Fits Calibrex™ 520/525/530 up to 10 mL, by the meter, uncut	specify	5 mm	511.709

Replacement feed tubing, FEP ⑦			
For Calibrex™ 525/530 as of 25 mL, cut	350 mm	7 mm	525.350
For Calibrex™ 525/530 as of 25 mL, by the meter, uncut	specify	7 mm	525.706

Telescopic feed tubing, FEP ⑦			
Fits Calibrex™ 520/525/530 up to 10 mL	150 - 255 mm	1/pk	1.525.352
Fits Calibrex™ 525/530 as of 25 mL	195 - 345 mm	1/pk	1.525.355


Description	Packaging	Cat. No.
Safety seal sticker		
Fits Calibrex™ 525/530 dispensers	25/pk	1.525.525
Fits Calibrex™ 525/530 dispensers	100/pk	1.525.526

► Dispenser stability ①

Stand holds dispenser when screwed on small-size bottle.

► Remote aspiration ②

Work station facilitates liquid intake from drum or other remote container (<10 m distance, <2 m elevation).

► Extension tubing ③

Spiral tubing and delivery Jet-Pen™ help dispense into vessels with maximum comfort within a 60 cm distance. No tool needed to assemble or to remove for cleaning. Autoclavable.

► Flow control stopcock ④

Fitting Calibrex™ 525/530 models. No tool needed to assemble or to remove for cleaning. Autoclavable.

► Delivery jet ⑤

Extended sizes exchangeable without any tool on Calibrex™ 525/530 models. Autoclavable.

► Alternative setting button ⑥

Classical screw-type button easily replaces the original sliding cursor, according to user preference.

► Standard and telescopic feed tubings ⑦

Supplied either cut to standard size or uncut, by the meter. Optional telescopic tubing for instant fitting to bottle size.

Calibrex™ *universal / organo / solutae*

your choice for flexibility and performance

Chemical resistance chart

Chemicals A – H	Calibrex™		
	520	525	530
Acetaldehyde (Ethanal)	A	A	A
Acetic acid 96%	A	A	B/2
Acetic acid 100% (Glacial)	A	B/4	B/2/4
Acetone (Propanone)	B/4	B/4	B/4
Acetonitrile (MECN)	A	B/4	B/4
Amino acids	A	C/1	A
Ammonium hydroxide (amonia)	A	B/4	B/4
Amyl alcohol (Pentanol)	A	A	A
Aniline	A	A	A
Ascorbic acid	A	C/1	A
Benzaldehyde	A	A	A
Benzene	B/4	B/4	B/4
Boric acid 10%	A	B/1	A
Bromine	B/2	C/4	C/2/4
Butanol	A	A	A
Butanone (MEK)	B/4	B/4	B/4
Butyl acetate	A	B/4	B/4
N-Butylamine	B/4	B/4	B/4
Calcium chloride	A	C/1	A
Calcium hydroxide	B/1	C/1	B/1
Carbon disulfide	A	B/4	B/4
Carbon tetrachloride	A	B/4	B/4
Chlorine dioxide	B/2/4	B/4	B/2/4
Chlorobenzene	A	B/4	B/4
Chlorobutane	A	B/4	B/4
Chloroethanol	A	B/4	B/4
Chloroform	B/4	B/4	B/4
Chlorosulfuric acid 100%	B/3	B/3/4	B/3/4
Chromic acid 100%	B/3	B/3/4	B/3/4
Citric acid	A	B/1	A
Cyanoacrylate	C/1	C/1	C/1
Cyclohexane	A	B/4	B/4
Cyclohexanone	A	B/4	B/4
1,4-Dioxane (Diethylene dioxide)	A	B/4	B/4
Dichlorobenzene	A	A	A
Dichloroethane (DCE)	B/4	A	A
Diesel oil (Heating oil)	A	A	A
Diethylene glycol	A	A	A
Diethylether	A	B/4	B/4
Dimethyl sulfoxide (DMSO)	A	B/1/4	B/4
Dimethylformamide (DMF)	B/4	B/4	B/4
Ethanol	A	A	A
Ether	B/4	B/4	B/4
Ethyl acetate	A	B/4	B/4
Ethylenediamine	A	A	A
Ethylene glycol	A	A	A
Formaldehyde (Formalin)	A	A	A
Formamide	A	A	A
Formic acid	A	A	A
Gamma-butyrolactone	A	A	A
Gasoline	A	B/4	B/4
Glycerin <40%	A	A	A
Heptane	A	A	A
Hexane	A	A	A
Hydrochloric acid 20%	A	A	A
Hydrochloric acid 37% (HCl)	A	B/3	B/3
Hydrofluoric acid (HF)	C/5	C/5	C/5
Hydrogen peroxide	A	A	B/2

Chemicals I – Z	Calibrex™		
	520	525	530
Iodine	A	C/1	B/1
Iodine bromide / chloride	C/2/4	C/4	C/2/4
Isooctane	A	A	A
Isopropanol	A	A	A
Isopropylamine	A	B/4	B/4
Lactic acid	A	C/1	A
2-Methoxyethanol	A	A	A
Methanol	A	A	A
Methyl chloride (Chloromethane)	A	B/4	B/4
Methyl methacrylate (MMA)	A	B/4	B/4
Methyl propyl ketone (2-Pentanone)	B/4	A	A
Methylene chloride (Dichloromethane) (DCM)	B/2/4	B/4	B/2/4
Nitric acid 100%	B/3	C/3/4	C/2/3/4
Nitric acid dil. <30%	A	B/4	B/4
Nitro-hydrochloric acid (Aqua regia)	B/3	B/4	B/2/4
N-methyl-2-pyrrolidone (NMP)	A	A	A
Octane	A	A	A
Octanol	A	A	A
Oil, mineral (engine oil)	A	A	A
Oil, vegetable, animal	A	B/4	B/4
Oil of turpentine	A	B/4	B/4
Oxalic acid	A	C/1	A
Pentane	B/4	B/4	B/4
Perchloric acid 100%	B/3	B/4	B/4
Perchloric acid diluted	A	A	A
Petroleum	A	B/4	B/4
Petroleum ether / spirit	A	B/4	B/4
Phenol	A	A	A
Phenylhydrazine	A	B/1/4	B/4
Phosphoric acid 85%	A	A	A
Potassium chloride	A	C/1	A
Potassium dichromate	A	C/1	B/1
Potassium hydroxide	B/1	C/1	A
Potassium iodide	A	C/1	A
Potassium permanganate	A	C/1	B/1
Propionic acid (Propanoic acid)	A	A	A
Propylene glycol (Propane-1,2-diol)	A	A	A
Picric acid (Trinitrophenol)	A	B/4	B/4
Pyridine	B/4	B/4	B/4
Scintillation fluid	A	A	A
Silver nitrate	B/1	C/1	A
Sodium acetate	A	C/1	A
Sodium chloride (Kitchen salt)	A	C/1	A
Sodium hydroxide 30%	B/1	C/1	A
Sodium hypochlorite (Javel water)	A	C/1	B/4
Sodium thiosulfate	A	C/1	A
Sulfonitic acid 100%	B/2/3	B/3/4	B/2/3/4
Sulfuric acid 98%	B/2	B/4	B/2/4
Tetrachloroethylene	B/4	B/4	B/4
Tetrahydrofuran (THF)	B/2/4	B/4	B/2/4
Toluene	B/4	B/4	B/4
Trichlorethylene	B/4	B/4	B/4
Trichloroacetic acid	A	B/1/4	B/4
Trichloroethane	B/4	B/4	B/4
Trichloromethane (Chloroform)	B/4	B/4	B/4
Triethylene glycol	A	A	A
Trifluoroacetic acid (TFA)	B/3	B/4	B/4
Xylene	B/4	B/4	B/2/4

Compatibility statement

A = Good resistance – B = Acceptable with limitations – C = Not recommended

Technical risks

- 1 = Possible crystallization, valve or plunger blockage. Scratches on plunger coating if plunger/barrel dried and stick together.
- 2 = Swell of plunger coating, possible peeling.
- 3 = Release of acid vapours (risk increases with concentration). Do not leave dispenser on bottle.
- 4 = Damage, softening or discoloration of external parts through vapours. Do not leave dispenser on bottle.
- 5 = Chemical degradation of glass parts (plunger/barrel).

Model differentiation 520 / 525 / 530

Ordering information - Instruments

Calibrex™ universal 520

32 mm base thread. Supplied with 300 mm feed tubing, 80 mm delivery jet with stopper, 28, 40, 45 mm bottle neck adapters, QC certificate and operating instructions.

Volume mL	Division mL	Cat. No. Fixed volume*	Cat. No. Adj. Volume
0.25 - 2	0.05	520.F02	520.002
1 - 5	0.1	520.F05	520.005
1 - 10	0.25	520.F10	520.010

* Specify desired fixed volume within instrument range when ordering.

Calibrex™ organo 525

45 mm base thread. Supplied with feed tubing, delivery jet with stopper, bottle neck adapters, QC certificate and operating instructions.

1 and 2.5 mL: 300 mm feed tubing, 90 mm delivery jet, 25, 28 and 32 mm adapters.

5 and 10 mL: 300 mm feed tubing, 90 mm delivery jet, 28, 32 and 40 mm adapters.

As of 25 mL: 350 mm feed tubing, 120 mm delivery jet, 32, 38 and 40 mm adapters.

Volume mL	Division mL	Cat. No. Without stopcock	Cat. No. With stopcock
0.1 - 1	0.02	525.001	525.001FC
0.25 - 2.5	0.05	525.002.5	525.002.5FC
0.5 - 5	0.1	525.005	525.005FC
1 - 10	0.2	525.010	525.010FC
2.5 - 25	0.5	525.025	525.025FC
5 - 50	1.0	525.050	525.050FC
10 - 100	1.0	525.100	525.100FC

Calibrex™ solutae 530

45 mm base thread. Supplied with feed tubing, delivery jet with stopper, bottle neck adapters, QC certificate and operating instructions.

1 and 2.5 mL: 300 mm feed tubing, 90 mm delivery jet, 25, 28 and 32 mm adapters.

5 and 10 mL: 300 mm feed tubing, 90 mm delivery jet, 28, 32 and 40 mm adapters.

As of 25 mL: 350 mm feed tubing, 120 mm delivery jet, 32, 38 and 40 mm adapters.

Volume mL	Division mL	Cat. No. Without stopcock	Cat. No. With stopcock
0.1 - 1	0.02	530.001	530.001FC
0.25 - 2.5	0.05	530.002.5	530.002.5FC
0.5 - 5	0.1	530.005	530.005FC
1 - 10	0.2	530.010	530.010FC
2.5 - 25	0.5	530.025	530.025FC
5 - 50	1.0	530.050	530.050FC
10 - 100	1.0	530.100	525.100FC

► Bottle-neck adapters

Each Calibrex™ dispenser comes along with three additional PP material bottle neck adapters. To fit additional bottle neck diameter, order the appropriate adapter size or combine two adapters to reach the adequate solution.

Type	Ext. Ø of bottle neck	Cat. No.
Adapters for Calibrex™ 525/530 ①		
Threaded, PP	32 - 25 mm	1.525.GL25
Threaded, PP	32 - 28 mm	1.525.GL28
Threaded, PP	45 - 32 mm	1.525.GL32
Threaded, PP	45 - 38 mm	1.525.GL38
Threaded, PP	45 - 40 mm	1.525.GL40
Adapters for Calibrex™ universal 520 ②		
Threaded, PP	22, 25, 28, 30, 34, 36, 38, 40, 45 mm	GLP + Ø
Threaded, PTFE	36, 38, 40, 45 mm	GLT + Ø
Tapered, PP	18.8, 24, 29.2, 45 mm	NSP + Ø
Tapered, PTFE	18.8, 24, 29.2 mm	NST + Ø

bottle-top dispensers

Reagent bottles

Glass and polyethylene reagent bottles supplied with its PP screw cap. Suitable for all bottle top dispensers. Corresponding neck adapters supplied with the dispensers.

Shape	Volume	Neck, ext. Ø	Cat. No.
Amber glass¹⁾ ①			
Square	100 mL	32 mm	314.0100
Square	250 mL	32 mm	314.0250
Square	500 mL	32 mm	314.0500
Square	1000 mL	45 mm	314.1000
Round	2500 mL	45 mm	314.2500
Amber glass, with handle¹⁾ ②			
Round	2500 mL	45 mm	314.2500H
PE coated amber glass²⁾ ③			
Square	500 mL	32 mm	314.0500PE
Square	1000 mL	45 mm	314.1000PE
Round	2500 mL	45 mm	314.2500PE
Clear borosilicate glass¹⁾ ④			
Round	250 mL	45 mm	314.0250C
Round	500 mL	45 mm	314.0500C
Round	1000 mL	45 mm	314.1000C
Round	2000 mL	45 mm	314.2000C
Clear Pyrex glass, with connection neck¹⁾ ⑤			
Round	500 mL	32 mm	314.0500P
Polyethylene²⁾ ⑥			
Square	250 mL	25 mm	315.0250
Square	500 mL	25 mm	315.0500
Square	1000 mL	32 mm	315.1000
Square	2500 mL	45 mm	315.2500

¹⁾ autoclavable, ²⁾ not autoclavable

Your local distributor

Safety first

Refer to package inserts for safety precautions, operating instructions and warranty terms.

Mind risks involved in handling hazardous liquids with respect to personal, third party and environmental protection and safety.

QC and warranty

The Calibrex™ line is manufactured and tested to fully comply with current regulations. Each instrument bears its own serial number and passes strict performance control attested by an individual QC certificate. Products and specifications are subject to change without prior notice.

US Patent pending

www.socorex.com/patents-en.html

SOCOREX ISBA SA
 Champ-Colomb 7, P.O. Box
 1024 Ecublens/Lausanne
 Switzerland
socorex@socorex.com
www.socorex.com
 Tel. +41 (0)21 651 6000
 Fax +41 (0)21 651 6001

B.CLX.E – B316