


Dräger Fabius® MRI Anaesthesia Workstations

Increase the diagnostic capability of your MRI unit with the help of state-of-the-art ventilation in the Dräger Fabius® MRI anaesthesia system specially designed for use in MRI environments.



Benefits

Fully certified for MRI

Specifically designed for use in magnetic environments, the Dräger Fabius MRI is certified for field strengths of up to 40 mTesla or 400 Gauss. This provides reliable ventilation with 1.5 and even 3 Tesla MRI systems.

Electronic ventilator

The Dräger Fabius MRI incorporates a servo-controlled, piston ventilator that requires no drive gas. The piston is capable of delivering higher inspiratory flows than the traditional bellows ventilator, and provides a more accurate volume delivery.

Broad range of therapy modes

The Fabius MRI gives you access to a broad range of ventilation modes - including volume, pressure, pressure support and SIMV/PS - providing a maximum level of respiratory care regardless of the patient acuity level.

Compact design, intelligent ergonomics

The Fabius MRI has a compact breathing system, COSY 2.6 which can be mounted on the left or right hand side of the unit, depending on the room layout. Mounted on a short, 8" arm to provide positioning flexibility, it simplifies use in an often cramped MRI environment. The Fabius MRI also gives you the option of using the Dräger CLIC disposable soda lime absorber canister for an easy dust-free soda lime canister exchange.

High visibility alarms

The Fabius MRI is equipped with highly visible repeater LED optical alarms that alert caregivers in the central station to changing patient conditions.

Standard operating philosophy

Because it uses Dräger's uniform user interface, Fabius MRI is quick and easy to learn. All relevant information is displayed on a single, high resolution LCD color screen (6.5 inch diagonal).

System components



Dräger Vapor® 2000 and D-Vapor®

Dräger vaporisers have been the benchmark for quality for over 50 years. Quality trusted by doctors and nurses around the world: to date, over 400,000 Vapor units have been sold to hospitals around the world.

Accessories



VentStar® MRI 300

Disposable breathing circuit, consisting of 2 smoothbore hoses, Y-piece and LuerLock elbow.

Suitable for use in MRI environment.

Length: 300 cm (118 inch).

Latex free.



VentStar® MRI (N) 300

Disposable neonatal breathing circuit, consisting of 2 smoothbore hoses (Ø 10 mm), angled Y-piece with LuerLock.

Suitable for use in MRI environment.

Length: 300 cm (118 inch).

Latex free.

Accessories



VacuSmart® Gel

The VacuSmart Gel is a disposable cartridge with an integrated gelling agent for all Dräger Medical bronchial aspirators using 700 ml secretion jars. The gelling agent coagulates the bronchial secretion and prevents thereby the leakage e.g. in a case of disposal in a waste press.



CLIC Absorber 800+

Disposable CLIC Absorber 800+ is filled with Drägersorb 800+, 1.2 I (42.2 fl oz.).



WaterLock® 2

Perfect protection for precise gas measurement. Dräger WaterLock® 2 safely stops water from getting into the Multi-Gas Sensor. The membrane technology developed by Dräger for the WaterLock® 2 stops any bacteria or germs from getting into the gas measurement system. The WaterLock® 2 is also safe and simple to empty – with a further advantage in handling and hygiene.

Related Products



Dräger Fabius® plus XL

The Dräger Fabius plus XL combines proven German engineering you can count on with high performance ventilation therapy. Thanks to its scalable design concept, it allows you to choose the quality workstation you want now without losing sight of your future goals and needs.



Dräger Fabius® Plus

Combine quality ventilation, easy operation and maintenance with open architecture expandability. The Dräger Fabius® Plus combines quality ventilation with enhanced flexibility and integration capabilities. It was designed to accommodate a wide range of options and accessories, allowing you to customize your Dräger Fabius® Plus to suit your particular needs.



Dräger Fabius® GS premium

The Dräger Fabius® GS premium is an anaesthesia workstation that is simple to use, highly efficient and ready for the future. It features a solid design with modular architecture plus a wide range of ventilation capabilities. Customize your Fabius® GS premium exactly the way you need it.



Dräger Fabius® Tiro

Get the most out of even the smallest spaces with a compact yet fully featured anaesthesia solution designed for use in a variety of specialized environments.

Technical Data

BASE UNIT

Dimensions (W x H x D)	
Trolley Version (Cart) with COSY	approx. 39 x 55 x 35.5 in (99 x 140 x 90 cm)
Neight and load	
Fabius MRI Trolley (with COSY) without supplementary cylinders	365 lbs. (165,8 kg)
and vaporizers	
Power and battery backup	
Power Input	100 to 240 VAC, 50 / 60 Hz, 70 VA, including additional power outlets
Operation time with fully charged batteries	> 45 min
ANAESTHESIA GAS SUPPLY MODULE	
Range of fresh gas flow indicators	0.0 to 12.0 L/min
Total fresh gas flowmeter	0 to 10 L/min
O ₂ flush	at 87 psi (6 bar): max 75 L/min;
	at 41 psi (2.8 bar): min. 25 L/min
Vapor mount	Dräger Interlock or Selectatec for:
	Dräger Isoflurane Vapor 2000,
	Dräger Sevoflurane Vapor 2000,
	Dräger Halothane Vapor 2000
Hosts Auto Exclusion and Interlock Vapor.	
VENTILATOR OPERATING SPECIFICATIONS	Flactronically controlled electrically driven
VENTILATOR OPERATING SPECIFICATIONS Ventilator E-vent®	Electronically controlled, electrically driven
	Volume Controlled Ventilation Pressure Controlled Ventilation
/entilator E-vent®	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support
Ventilator E-vent®	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS
Ventilator E-vent®	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation
Ventilator E-vent®	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS
Ventilator E-vent® Operating modes Control input ranges	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing
Ventilator E-vent® Operating modes Control input ranges Breathing Frequency (rate)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa)
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa)
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiration pause (Tip:Ti)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 %
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiratory pause (Tip:Ti) SIMV Inspiratory time	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 % 0.3 - 4.0 sec
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiration pause (Tip:Ti) SIMV Inspiratory pressure (Pinsp)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 % 0.3 - 4.0 sec PEEP + 5 to 65 cm H ₂ O (hPa)
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiration pause (Tip:Ti) SIMV Inspiratory time	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 % 0.3 - 4.0 sec PEEP + 5 to 65 cm H ₂ O (hPa) 10 to 75 L/min in Volume and Pressure Control modes
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiratory pressure (Tip:Ti) SIMV Inspiratory time Inspiratory Flow (InspFlow) Inspiratory Flow (InspFlow)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 % 0.3 - 4.0 sec PEEP + 5 to 65 cm H ₂ O (hPa) 10 to 75 L/min in Volume and Pressure Control modes 10 to 85 L/min in Pressure Support and SIMV/PS modes
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory Flow (InspFlow) Inspiratory Flow (InspFlow) Inspiratory Support Level (PPS)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 % 0.3 - 4.0 sec PEEP + 5 to 65 cm H ₂ O (hPa) 10 to 75 L/min in Volume and Pressure Control modes 10 to 85 L/min in Pressure Support and SIMV/PS modes PEEP + 3 to 20 cmH ₂ O (hPa)
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory Flow (InspFlow) Pressure Support Level (PPS) Min. frequency for apnoe-ventilation (Freq. Min.)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 % 0.3 - 4.0 sec PEEP + 5 to 65 cm H ₂ O (hPa) 10 to 75 L/min in Volume and Pressure Control modes 10 to 85 L/min in Pressure Support and SIMV/PS modes PEEP + 3 to 20 cmH ₂ O (hPa) 3 to 20 bpm and "OFF"
Control input ranges Breathing Frequency (rate) Positive End Expiratory Pressure (PEEP) Inspiration/expiration ratio (Ti:Te) Pressure limiting (Pmax) Fidal Volume (Vt) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory pressure (Pinsp) Inspiratory Flow (InspFlow) Inspiratory Flow (InspFlow) Inspiratory Support Level (PPS)	Volume Controlled Ventilation Pressure Controlled Ventilation Pressure Support SIMV/PS Manual Ventilation Spontaneous Breathing 4 to 60 bpm 0 to 20 cmH ₂ O (hPa) 4:1 to 1:4 15 to 70 cm H ₂ O (hPa) 20 to 1400 mL in Volume Control 20 to 1100 mL in SIMV/PS 0 to 50 % 0.3 - 4.0 sec PEEP + 5 to 65 cm H ₂ O (hPa) 10 to 75 L/min in Volume and Pressure Control modes 10 to 85 L/min in Pressure Support and SIMV/PS modes PEEP + 3 to 20 cmH ₂ O (hPa)

Technical Data

	N ₂ O cut-off if O ₂ fresh gas valve is closed or if O ₂ flow is less than
	0.2 L/min.
	Audible and visual (flashing red LED) indication in case O_2 pressure drops below 1.38 bar (20 psi) \pm 0.27 bar (4 psi). In case of electricity and battery failure, manual ventilation, gas delivery and agent delivery are possible. Positive pressure relief valve opens at 75 \pm 5 cm H ₂ O. Negative pressure relief valve open at -7.5 to -9 cm H ₂ O.
Ventilator monitoring	
Monitoring	Continuous monitoring of inspiratory O ₂ concentration, breathing frequency, tidal volume (expiratory), minute volume (expiratory), peak airway pressure, PEEP, and selection of mean or plateau pressure. In addition, all fresh gas flow information is displayed as
	virtual flow tubes.
Expiratory Minute Volume range	0 to 99 L/min
Control Screen	6.5 in (16.5 cm) color screen
BREATHING SYSTEM	
Volume of entire compact breathing system	1.7 L + bag
	1.7 L + bag 1.5 L (standard) [option: Prefilled Dräger Sorb CLIC absorber with 1.2 Liter]
Volume of CO ₂ absorber	1.5 L (standard) [option: Prefilled Dräger Sorb CLIC absorber with
Volume of CO ₂ absorber GAS SUPPLY AND CONNECTION	1.5 L (standard) [option: Prefilled Dräger Sorb CLIC absorber with
Volume of CO ₂ absorber GAS SUPPLY AND CONNECTION Gas Supply	1.5 L (standard) [option: Prefilled Dräger Sorb CLIC absorber with 1.2 Liter]
Volume of CO ₂ absorber GAS SUPPLY AND CONNECTION Gas Supply Cylinder Yokes	1.5 L (standard) [option: Prefilled Dräger Sorb CLIC absorber with 1.2 Liter] O ₂ , N ₂ O & Air
Volume of entire compact breathing system Volume of CO ₂ absorber GAS SUPPLY AND CONNECTION Gas Supply Cylinder Yokes OTHER Writing surfaces	1.5 L (standard) [option: Prefilled Dräger Sorb CLIC absorber with 1.2 Liter] O ₂ , N ₂ O & Air

Notes

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA Moislinger Allee 53–55 23558 Lübeck, Germany www.draeger.com

Manufacturer:

Dräger Medical GmbH Moislinger Allee 53-55 23558 Lübeck, Germany

As of August 2015

Dräger Medical GmbH changes to Drägerwerk AG & Co. KGaA

Locate your Regional Sales Representative at: www.draeger.com/contact



REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH Moislinger Allee 53-55 23558 Lübeck, Germany Tel +49 451 882 0 Fax +49 451 882 2080 info@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial S. de R.L. Complejo Business Park, V tower, 10th floor Panama City Tel +507 377 9100 Fax +507 377 9130 contactcsa@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH Branch Office P.O. Box 505108 Dubai, United Arab Emirates Tel +971 4 4294 600 Fax +971 4 4294 699 contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Medical South East Asia Pte Ltd. 25 International Business Park #04-27/29 German Centre Singapore 609916, Singapore Tel +65 6572 4388 Fax +65 6572 4399 asia.pacific@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dImfr-contact@draeger.com

CANADA

Draeger Medical Canada Inc. 2425 Skymark Avenue, Unit 1 Mississauga, Ontario, L4W 4Y6 Tel +1 905 212 6600 Toll-free +1 866 343 2273 Fax +1 905 212 6601 Canada.support@draeger.com