Operation Manual

Auto Hematology Analyzer Reagent (M-52LH Lyse)

English (3 – 5)

Spanish (6 - 8)

French (9 – 11)

Portuguese (12 – 14)

Russian (15 – 18)

Turkish (19 – 21)

© 2012-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

Intellectual Property Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyright or patents and does not convey any license under the patent rights or copyright of Mindray, or of others.

Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

Release, amendment, reproduction, distribution, rental, adaptation, translation or any other derivative work of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

mindray

MINDRAY

are the trademarks, registered or otherwise, of Mindray in China and other countries. All other trademarks that appear in this manual are used only for informational or editorial purposes. They are the property of their respective owners.

Responsibility on the Manufacturer Party

Contents of this manual are subject to change without prior notice.

All information contained in this manual is believed to be correct. Mindray shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this manual.

Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- 1) Use of this product is conducted as instructed in this operation manual by Mindray authorized personnel or skilled/trained clinical professionals;
- 2) The product is stored and transported as instructed in the operation manual and the labeling;
- 3) The product is used only on the specified models.

Warranty

Free service will be provided to:

- Products that are covered by the warranty terms of Mindray.

Pay service will be provided to:

- 1) Expired products;
- 2) Products within the warranty period, but are damaged by improper use or human error, stored in improper environment, or damaged by force majeure or other reason not caused by the product itself.

Return Procedure

Follow the procedure below to return products to Mindray if needed.

Contact Mindray customer service department, inform them the product name and lot No. (the name and lot No. are labeled on the outer package of the product). If the name and lot No. are not legible, your return request will not be accepted. Please specify the product name, model, lot No. and the return reason.

M-52LH Lyse

1 Product

M-52LH Lyse

2 Model

M-52LH

3 Specification

100mL×4

4 Intended use

M-52LH Lyse dissolves RBCs, realizes WBC counting, basophil differential and hemoglobin measurement.

5 Principle

M-52LH lyse further processes the blood cells that have been processed by M-52DIFF lyse, it turns hemoglobin into methemoglobin, which then measured by the colorimetric method. It also processes WBCs to intensify the difference between basophils and other WBC sub-populations, and incorporates laser scatter and flow cytometry to realize basophil analysis.

6 Ingredients

Quaternary ammonium salts.....<50g/L

Isopropanol.....2-10g/L

7 Storage condition and shelf life

The product must be stored in the temperature of 2°C~30°C, and the relative humidity must not exceed 90%. The storage room must be well ventilated and without corrosive gas. The shelf life is 2 years.

The product must be used in the temperature of 10°C~30°C, and it stays valid for 60 days after being opened.

8 Applicable analyzer models

The product applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd..

9 SAMPLE REQUIREMENT

Fresh human whole blood anticoagulated. Do not use contaminated samples.

10 MATERIALS REQUIRED BUT NOT PROVIDED

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents.

11 Instructions for use

- 1) Restore the M-52LH lyse to usage temperature;
- 2) Open the external packing and insert the pickup tube into the reagent container based on the matching color of the reagent cap and the connector of analyzer cap assembly;
- 3) Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
- 4) Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis. See the operator's manual of the analyzer for details.

12 CUT-OFF VALUE/REFERENCE INTERVALS

Not applicable.

13 RESULT ELABORATION

Not applicable.

14 LIMITATIONS

Not applicable.

15 PRODUCT SPECIFICATIONS

- 1) Appearance: transparent liquid, no sediments, grains or flocculi.

2) Blank count results: the blank count results of M-52LH lyse tested on Mindray Analyzers meet the requirements in Table 1.

Table 1 Requirements of M-52LH lyse blank count

Test parameter	Blank count requirements
WBC	$\leq 0.2 \times 10^9/L$
HGB	$\leq 1g/L$

16 PRECAUTIONS

- 1) For professional use in vitro diagnosis only.
- 2) Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
- 3) Do not use the reagent if it is frozen.
- 4) If the blank count is abnormal after the reagent is transported, place it still for 24 hours at room temperature before using it.
- 5) If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
- 6) Dispose of the waste, residual and contaminated packing based on local regulations.
- 7) The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
- 8) Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.
- 9) Disposal of waste liquid and materials should be in accordance with local guidelines.
- 10) The Material Safety Data Sheet (SDS) is available upon request
- 11) All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
- 12) Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- 13) This product contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

 	
Warning	
H411	Toxic to aquatic life with long lasting effects.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
Prevention	
P273	Avoid release to the environment.
P280	Wear protective gloves, protective clothing, eye protection and face protection.
P264	Wash all exposed external body areas thoroughly after handling.
Response	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313	If eye irritation persists: Get medical advice/attention.
P391	Collect spillage.
P302+P352	IF ON SKIN: Wash with plenty of water.
P332+P313	If skin irritation occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation

17 GRAPHICAL SYMBOLS

		
Batch code	Use-by date	Temperature limit
		
HUMIDITY LIMITATION	Manufacturer	In vitro diagnostic medical device
		
CATALOGUE NUMBER	Authorized representative in the European Community	Consult instructions for use
		
European Conformity	Unique Device Identifier	

18 REFERENCES

Not applicable.

19 COMPANY CONTACT

Manufacturer:

Address:

E-mail Address:

Tel:

Fax:

Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Mindray Building, Keji 12th Road South, High-tech industrial park, Nanshan, Shenzhen 518057, P.R.China

service@mindray.com

+86 755 81888998

+86 755 26582680

EC-Representative:

Address:

Tel:

Fax:

Shanghai International Holding Corp. GmbH(Europe)

Eiffestraße 80, 20537 Hamburg, Germany

0049-40-2513175

0049-40-255726

20 Approval Date of the Operator's Manual

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