RealPCR* CSFV RNA Test

English Version

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REF 99-56022

 **Version**
06-56022-01

*This insert is not for a product available to U.S. or Canadian customers.
Please refer to the RealPCR insert webpage and select "United States
or Canada" to view the correct product insert.*

Test With Confidence™

IDEXX

MENU

PRINT

English version

RealPCR* CSFV RNA Test

For veterinary use only

Name and Intended Use

The RealPCR* CSFV RNA Test is used for the detection of classical swine fever virus (CSFV) RNA extracted from blood (EDTA), serum, plasma, oral fluids, swab samples or tissue (spleen, kidney, lymph node and tonsil) from swine, including wild boars. Tissue samples can be tested in pools of up to 10 samples and blood (EDTA), plasma and serum can be tested in pools of up to 20 samples. Oral fluids can be tested as a composite sample taken from pens of up to 30 pigs. Pools containing a single weak sample (for example, Ct > 32) may yield a negative result due to the dilution effect of pooling.

General Information

CSFV is the causative agent of classical swine fever (CSF), a highly contagious hemorrhagic disease affecting wild and domestic pigs. CSF is characterized by high morbidity and mortality rates, causing serious economic losses to the pig industry, and it is a World Organization for Animal Health (OIE)-listed disease.

Pigs infected with CSFV may shed a high amount of virus before showing clinical signs of the disease. If animals survive an acute or sub-acute infection, they can become chronically infected and excrete the virus intermittently or continuously until death. In pregnant sows, CSFV is able to cross the placenta and infect fetuses, causing abortions, fetal mummifications and stillborns. In mid-gestation infections (~50–70 days of pregnancy), weak or persistently viremic piglets can be born. These persistently infected piglets are able to shed high levels of virus for several months.

The IDEXX* RealPCR system is a modular format in which disease-specific target mixes are paired with a standardized DNA or RNA master mix and a single pooled positive control. Reagents are individually packaged and sold separately to allow for flexible reagent handling.

The RealPCR CSFV RNA Mix (CSFV RNA Mix) contains primers and probes for the detection of CSFV RNA when amplified with RealPCR RNA Master Mix (RNA MMx). The assay is a single-tube reverse transcriptase and polymerase reaction. The internal control for the test is based on the detection of an endogenous porcine RNA sequence present in the host sample and is referred to as the internal sample control (ISC) in this protocol. Detection of endogenous RNA in swine samples controls for sample addition, extraction and amplification. Primers and probe for detection of the internal sample control are included in the CSFV RNA Mix. An optional internal positive control, the RealPCR Internal Positive Control (IPC ≥ v1.1), is also available and should be used when endogenous host RNA is at low levels or unlikely to be present after extraction (such as environmental samples). The IPC contains a synthetic version of the swine ISC RNA target and is therefore compatible with the CSFV RNA Mix. Refer to the RealPCR Internal Positive Control (REF 99-56330) product insert for guidance.

Materials and Storage

Identification/ General Information	Cap color	Quantity	Storage		Freeze/Thaw cycles
		100 tests	At receipt	After reconstitution	
RealPCR[®] CSFV RNA Mix (CSFV RNA Mix), dried	Yellow	1 x 1.0 mL	-25 to 8°C	-25 to -15°C	≤6
<p>[REF] 99-56022</p> <p>Reconstitute to 1 mL in PCR Grade Water. Store the CSFV RNA Mix in the dark. The expiration date on the vial is valid for either the dry or reconstituted form. The CSFV RNA Mix tube label indicates the PC version that is compatible with that target mix. For example, PC ≥ v1.4 means the target mix can be used with PC versions greater or equal to 1.4.</p>					
RealPCR RNA Master Mix (RNA MMx)	Black	1 x 1.0 mL	-25 to -15°C (Long-term)	N/A	≤6
<p>[REF] 99-56280</p> <p>Concentrated master mix that includes reverse transcriptase and hot-start polymerase for use with RNA target mixes in the IDEXX RealPCR system. The RNA MMx is more viscous than most master mixes— see the Test Procedure section for handling recommendations. A reference dye (ROX[®]) has been added for normalizing volume inaccuracies. Protect the RNA MMx from light.</p>					
RealPCR Positive Control, dried (PC)	Blue	1 x 500 μL	-25 to 8°C	-25 to -15°C	≤6
<p>[REF] 99-56310</p> <p>Reconstitute to 500 μL in PCR Grade Water. The PC contains all IDEXX[®] RealPCR and ISC targets (including the target for CSFV) and is intended for use with all IDEXX RealPCR target mixes. The expiration date on the vial is valid for either the dry or reconstituted form. The PC has a version number (such as v1.3). When new target mixes are developed for the RealPCR product line, the target sequences are added to the PC and the PC version number increases (for example, v1.3 would be updated to v1.4).</p> <p>The PC includes the IDEXX Signature (a unique oligonucleotide sequence). Presence of the IDEXX Signature, in the work environment, indicates PC contamination. For laboratories that want to monitor for PC contamination, the IDEXX Signature can be detected using the RealPCR PC Tracker DNA Mix and RealPCR DNA MMx.</p>					
RealPCR PCR Grade Water	Clear	2 x 1.0 mL	-25 to 8°C		N/A
<p>[REF] 99-56350</p> <p>PCR Grade Water has been qualified for reverse transcription-PCR (RT-PCR) use. It is used for the reconstitution of RealPCR reagents. It is also used as the PCR negative control for each test run. Do not transfer PCR Grade Water vials between PCR work areas. Separate vials of water are needed for each area to avoid contamination risk.</p>					

Note: See table at the end of the insert for a description of symbols used on the insert and labels.

Materials Required but Not Provided

- Commercial RNA extraction kit
 - Optional—Centrifuge with a rotor and adapters for multi-well plates
 - Micro-centrifuge for 2 mL microtubes capable of reaching 1500 – 3000 x g
 - Appropriate personal protective equipment (e.g., gloves, lab coat)
 - Nuclease-free, aerosol-resistant pipette tips
 - Sterile microtubes for preparation of PCR Mix
 - Pipettes (5–1000 μ L); dedicated pipettes for preparation of PCR Mix
 - 96- or 384-well format PCR plates and optical adhesive film/plate covers
 -  Real-time PCR instrument (Applied Biosystems[®] 7500, Applied Biosystems[®] 7500 Fast System (Standard and Fast Mode), Applied Biosystems[®] ViiA[®] 7, Applied Biosystems QuantStudio 5, Applied Biosystems QuantStudio[®] Flex, Agilent Mx3000P[®], Agilent Mx3005P[®], Agilent AriaMx, Bio-Rad CFX96 Touch[®], Bio Molecular Systems Mic qPCR Cycler, QIAGEN[®] Rotor-Gene (72-Well Rotor only), Roche LightCycler[®] 480 or equivalent).
- Note-** the Roche LC480 instrument requires additional calibration and software settings. IDEXX[®] Technical Service can provide guidance for use of the instruments mentioned above with RealPCR[®] reagents.

Laboratory Practices and Warnings

- Do not use reagents past expiration date.
- The entire procedure must be performed under nuclease-free conditions.
- Wear powder-free gloves when working with the reagents and nucleic acids.
- To avoid cross-contamination, use nuclease-free, aerosol-resistant pipette tips for all pipetting, and physically separate the workplaces for nucleic acid extraction/handling, PCR setup and PCR.

Reconstitution of Dried Components

Reconstitute the CSFV RNA Mix and Positive Control by pipetting PCR Grade Water to the volume indicated on the component label. Allow to sit at 18 to 26°C for at least 10 minutes; mix and microcentrifuge briefly prior to use. Once the CSFV RNA Mix and the Positive Control are reconstituted, aliquot as appropriate and store the solutions frozen. When handling frozen components, thaw at 18 to 26°C for approximately 15 to 30 minutes, mix gently and then microcentrifuge briefly (~1,500 – 3,000 x g).

RNA Extraction

The CSFV RNA Mix has been validated using the commercial extraction methods listed below. Other extraction or lysis methods may also be used once validated by the laboratory.

-  RealPCR DNA/RNA Magnetic Bead Kit (IDEXX)
- RealPCR DNA/RNA Spin Column Kit (IDEXX)
- NucleoSpin RNA Kit (Machery Nagel)
- QIAamp Viral RNA Mini Kit (Qiagen)

Store the purified RNA at <-15°C if testing is not performed immediately after RNA extraction. It is recommended that a negative extraction control (mock sample) is included as a sample.

Test Procedure

- 1 Preparation of the PCR Mix.
 - Mix the thawed RNA MMx by inversion or gentle vortex.
 - The RNA MMx is a viscous solution; always pipette it slowly.
 - To prepare the PCR Mix add 10 μL CSFV RNA Mix and 10 μL RNA MMx for each reaction.
 - When preparing the PCR Mix, first pipette CSFV RNA Mix into the tube and then add the RNA MMx. Pipette up and down a few times to rinse the MMx pipette tip.
 - Gently vortex the solution to ensure the components are mixed well.
 - Pipette the PCR Mix slowly into the PCR plate.

Load the PCR plate within 20 minutes or store at 2 to 8°C for up to 4hrs. The PCR Mix can be stored at –25 to –15°C for up to 2 weeks. Protect from light.
- 2 Pipette 20 μL of the PCR Mix into the required wells of the multiwell plate.
- 3 Add 5 μL of sample RNA to each well. The final reaction volume is 25 μL .
- 4 Include the Positive Control (5 μL) and PCR Negative Control (5 μL PCR Grade Water) for each test run.
- 5 Cover the plate and briefly spin the plate, if necessary, to settle contents and remove air bubbles.
- 6 Set up thermal cycler with IDEXX® RealPCR® Standard DNA/RNA Cycling Program.

Settings for Reporter and Quencher

<u>Target</u>	<u>Reporter</u>	<u>Quencher</u>
CSFV	FAM*	BHQ* (none)
Internal Control (ISC/IPC)	HEX* (VIC)	BHQ (none)
Passive Reference	ROX*	N/A

RealPCR Standard DNA/RNA Cycling Program

	<u>Temperature</u>	<u>Time</u>	<u>Cycles</u>
Reverse transcription (RT)	50°C	15 min.	1
Denaturation	95°C	1 min.	1
Amplification [§]	95°C	15 sec.	45
	60°C	30 sec.	

[§]Ensure the instrument is set to record fluorescence following the 60°C amplification step.

7 Analyze data

When setting up the instrument software, assign a unique identifier to each target and internal control. For example, when target A and B are on the same plate, the A wells must be analyzed independently from the B wells. Refer to specific instrument's user manual for guidance on how to analyze data. To set the threshold, use the Auto Ct setting.

- Agilent Mx3000P* and Mx3005P* – ensure the background-based threshold fluorescence method is being used for analysis.
- QIAGEN* Rotor-Gene instrument – manually adjust the threshold line above background in the linear phase of exponential amplification. This is best done on log view plots and should be repeated for each reporter in the target mix.
- Applied Biosystems* instruments – In certain situations the auto threshold setting does not produce satisfactory results. In these situations, a manual adjustment of the threshold is necessary to determine Ct values. This is best done on log view plots and should be repeated for each reporter in the target mix.

Validity Criteria

	<u>FAM*</u> Ct value	<u>HEX* (VIC)</u> Ct value
Positive Control	<38	<38
PCR Negative Control	No Signal	No Signal

8 Interpretation of results

<u>Sample Result</u>	<u>FAM Signal</u>	<u>HEX (VIC) Signal</u>	<u>Other characteristics</u>
CSFV RNA detected	Yes	Yes/No	A positive Ct value and characteristic amplification curve in comparison to the PCR Negative Control. An internal control amplification curve in the HEX (VIC) channel is expected; some strong CSFV-positive RNA samples may result in a negative internal control result. [†]
CSFV RNA not detected	No	Yes	Amplification curve in the HEX (VIC) internal control channel.
Invalid [‡]	No	No	Absence of an amplification curve in the FAM and HEX (VIC) channels indicates an invalid result.

[†]The target mix is optimized for the detection of CSFV RNA; a strong positive RNA sample may out compete the detection of the internal control.

[‡]An invalid sample can be an indication of failed sample addition, extraction and/or PCR. It is recommended that the RNA be diluted five-fold into PCR grade water and retested; include the undiluted RNA as a sample. If the test is still not valid a new extraction is recommended.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contact/ptd

Patent information: idexx.com/patents

©2024 IDEXX Laboratories, Inc. All rights reserved.

 Dye compounds in this product are sold under license from Biosearch Technologies, Inc. and protected by U.S. and world-wide patents either issued or in application. This license covers veterinary and human applications, limited to research and development and diagnostic use.

*IDEXX, RealPCR and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories or its affiliates in the United States and/or other countries. All other product and company names and logos are trademarks or registered trademarks of their respective holders.

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IDEXX Montpellier SAS
IDEXX B.V.