Protectalon

Na Hyaluronate

%1.0

%1.2

%1.4

%1.6

%1.8

%2.0

(NL)	Instructies voor gebruik	1-2
(EN)	Instructions for use	3-4
(F)	Notice d'Utilisation	5-6
(DE)	Gebrauchsanweisung	7-8
(TR)	Kullanma Talimatı	9-10
(ES)	Instrucciones de uso	11-12
(IT)	Istruzioni per l'uso	13-14
(PT)	Instruções de Uso	15-16
(RU)	Инструкция по использованию	17-18
(AZ)	İstifadə üçün təlimat	19-20
(CN)	操作手册	21-22
(KK)	Пайдалану бойынша нұсқаулық	23-24
(HR)	Upustvo za uporabu	25-26
(NO)	Bruksanvisning	27-28
(THAI)	เอกสารกำกับเครื่องมือแพทย์	29-30

Protectalon

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Cohesive sodium hyaluronate for intraocular use

Product:

Protectalon is a viscoelastic sodium hyaluronate solution.

Each product consists of 1.1 ml of viscoelastic solution in a single-use glass syringe for intraocular use, together with one sterile 27G cannula.

Description:

Protectalon Viscosurgical Device is sterile, non-pyrogenic, clear, non-inflammatory, highly purified sodium hyaluronate of high molecular weight, dissolved in a buffered physiological saline solution. The highly-purified sodium hyaluronate is obtained from bacteria by fermentation.

Protectalon contains sodium hyaluronate, sodium chloride, disodium hydrogen phosphate, sodium dihydrogen phosphate and water for injection.

Protectalon is easy to instil. Sodium hyaluronate is viscous at rest, but its pseudo-plastic properties allow for great ease of flow when pressure is applied to expel the solution through the cannula. After instillation the viscous properties are restored.

	Protectalon 1.0 %	Protectalon 1.2 %	Protectalon 1.4 %	Protectalon 1.6 %	Protectalon 1.8 %	Protectalon 2.0 %
1.0 ml contains	10 mg sodium hyaluronate	12 mg sodium hyaluronate	14 mg sodium hyaluronate	16 mg sodium hyaluronate	18 mg sodium hyaluronate	20 mg sodium hyaluronate
pH	6.8-7.6	6.8-7.6	6.8-7.6	6.8-7.6	6.8-7.6	6.8-7.6
Osmolality	300-350 mOsm/kg	300-350 mOsm/kg	300-350 mOsm/kg	300-350 mOsm/kg	300-350 mOsm/kg	300-350 mOsm/kg
Zero shear viscosity (App.)	60.000 mPa.s.	130.000 mPa.s.	200.000 mPa.s.	400.000 mPa.s.	600.000 mPa.s.	900.000 mPa.s.
Molecular Weight (App.)	3.0M Da	3.0M Da	3.0M Da	3.0M Da	3.0M Da	3.0M Da

Indications:

- Intraocular surgery of the anterior segment, including cataract extraction and intraocular lens implantation.
- The viscoelastic properties of Protectalon cohesive allow lubrication, support and protection of ocular tissues during ophthalmic surgery.
- Maintains the depth of the anterior chamber during surgery, allows for efficient manipulation with reduced trauma to the corneal endothelium and other surrounding tissues.

Contraindications:

At present, there are no known contraindications to the use of Protectalon Ophthalmic Viscosurgical Device when used as recommended: Care should be taken in patients with hypersensitivity to any components in this solution.

Warnings and Precautions:

- Residual viscoelastic not removed at the end of the surgical procedure is flushed out naturally via the trabecular system and Schlemm's canal. There is, however, a risk of blockage of the drainage channels, which can lead to an increase in intraocular pressure. To prevent the risk of IOP, irrigation/aspiration is used at the end of the surgical procedure to remove residual viscoelastic. In case of IOP increase, treat with the preferred therapy. Patients who have been diagnosed with wide-angle glaucoma, severe myopia, diabetic retinopathy or uveitis are at particular risk of increased intraocular therapy.
- Do not reuse syringe and cannula. Any repeat usage of the syringe and cannula carries a risk of contamination and infection of the patient.
- Do not re-sterilize the pre-filled syringe and cannula.
- Do not use if package is damaged or opened.
- Do not use after the expiry date printed on the pack.
- Dispose of the syringe and cannula in accordance with accepted medical practice and applicable national, local and institutional requirements.
- There is no evidence concerning the safety of Protectalon in human pregnancy and lactation. Administration during pregnancy and lactation is at the discretion of the ophthalmic surgeon.

Interaction with other agents:

Sodium hyaluronate is incompatible with quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore, Protectalon pre-filled syringes must never come into contact with surgical instruments rinsed with these solutions, or with ophthalmic products containing quaternary ammonium compounds as preservatives.

Directions for Use:

FOR INTRAOCULAR USE. FOR SINGLE USE ONLY.

Protectalon should only be used by a specialist trained in intraocular surgery.

1. Twist and pull out

2. Insert and twist the cannula

WRONG!

(Backstop opening at front)

CORRECT!


(Backstop opening at back)

Press down on the plunger gently to expel a few drops of the product to prevent the introduction of air bubbles into the anterior chamber of the eye. The syringe is ready to use. The amount of the viscoelastic used is determined by the type of surgery. Administration is by injection into the anterior chamber at different stages of the surgical procedure. At the end of surgery, irrigation /aspiration is used at the end of the surgical procedure to remove residual viscoelastic.

Storage:

Store Protectalon between 2°C and 25°C. Protect from light and shocks. Do not freeze.

The graph below represents the Rheological Profile of Protectalon OVD at 25°C.

Symbols:

Manufacturer

Do not reuse

Do not resterilize

Sterilized using steam (Product)

Sterilized using ethylene oxide (Packaging and cannula)

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