oceanus 14 pro

PTA balloon dilatation catheter

To reach the **most challenging** lesions
with the **highest visibility**

Best flexibility
and trackability*

Smallest balloon
diameters

Enhanced
visibility*



Lowest
profiles

Proprietary hydrophilic
durable coating

oceanus 14 pro features

* Radiopaque markers	
	$\varnothing \leq 2.0\text{mm}$, Length $\leq 20\text{mm}$: 1 tungsten marker
	$\varnothing \leq 2.0\text{mm}$, Length $> 20\text{mm}$: 2 tungsten markers
	$\varnothing \geq 2.5\text{mm}$: 2 Pt/Ir markers

- > Over the wire catheter (OTW)
- > Usable Catheter length: 100 cm & 150 cm
- > Guidewire compatibility: 0.014"
- > Tip profile: 0.017"
- > Crossing profile: 0.021" to 0.033"
- > Nominal pressure: 7 atm
- > Rated Burst Pressure (RBP): 16 atm
- > Average Burst Pressure (ABP): 22 atm

Product with CE mark, certified by Notified Body 0318

Usable catheter length	Balloon Ø (mm)	Balloon Length (mm)								
		10	15	20	40	60	80	120	150	200
100 cm	1.25	BP PC14N 100 125 010	BP PC14N 100 125 015	BP PC14N 100 125 020	-	-	-	-	-	-
	1.5	BP PC14N 100 150 010	BP PC14N 100 150 015	BP PC14N 100 150 020	BP PC14N 100 150 040	BP PC14N 100 150 060	BP PC14N 100 150 080	-	-	-
	2.0	BP PC14N 100 200 010	BP PC14N 100 200 015	BP PC14N 100 200 020	BP PC14N 100 200 040	BP PC14N 100 200 060	BP PC14N 100 200 080	BP PC14N 100 200 120	BP PC14N 100 200 150	BP PC14N 100 200 200
	2.5	-	-	-	BP PC14N 100 250 040	BP PC14N 100 250 060	BP PC14N 100 250 080	BP PC14N 100 250 120	BP PC14N 100 250 150	BP PC14N 100 250 200
	3.0	-	-	-	BP PC14N 100 300 040	BP PC14N 100 300 060	BP PC14N 100 300 080	BP PC14N 100 300 120	BP PC14N 100 300 150	BP PC14N 100 300 200
	3.5	-	-	-	BP PC14N 100 350 040	BP PC14N 100 350 060	BP PC14N 100 350 080	BP PC14N 100 350 120	BP PC14N 100 350 150	BP PC14N 100 350 200
	4.0	-	-	-	BP PC14N 100 400 040	BP PC14N 100 400 060	BP PC14N 100 400 080	BP PC14N 100 400 120	-	-
150 cm	1.25	BP PC14N 150 125 010	BP PC14N 150 125 015	BP PC14N 150 125 020	-	-	-	-	-	-
	1.5	BP PC14N 150 150 010	BP PC14N 150 150 015	BP PC14N 150 150 020	BP PC14N 150 150 040	BP PC14N 150 150 060	BP PC14N 150 150 080	-	-	-
	2.0	BP PC14N 150 200 010	BP PC14N 150 200 015	BP PC14N 150 200 020	BP PC14N 150 200 040	BP PC14N 150 200 060	BP PC14N 150 200 080	BP PC14N 150 200 120	BP PC14N 150 200 150	BP PC14N 150 200 200
	2.5	-	-	-	BP PC14N 150 250 040	BP PC14N 150 250 060	BP PC14N 150 250 080	BP PC14N 150 250 120	BP PC14N 150 250 150	BP PC14N 150 250 200
	3.0	-	-	-	BP PC14N 150 300 040	BP PC14N 150 300 060	BP PC14N 150 300 080	BP PC14N 150 300 120	BP PC14N 150 300 150	BP PC14N 150 300 200
	3.5	-	-	-	BP PC14N 150 350 040	BP PC14N 150 350 060	BP PC14N 150 350 080	BP PC14N 150 350 120	BP PC14N 150 350 150	BP PC14N 150 350 200
	4.0	-	-	-	BP PC14N 150 400 040	BP PC14N 150 400 060	BP PC14N 150 400 080	BP PC14N 150 400 120	-	-

● Introducer compatibility: 4F


Distributed by:

Manufactured by:

Life Vascular Devices Biotech, S.L.
www.ivascular.global
info@ivascular.global

iVascular® oceanus 14 pro

EN	INSTRUCTIONS FOR USE. PTA BALLOON DILATATION CATHETER	02
ES	INSTRUCCIONES DE USO. CATÉTER BALÓN DE DILATACIÓN PARA ATP	03
PT	INSTRUÇÕES DE USO. CATETER-BALÃO DE DILATAÇÃO PARA ATP	04
IT	ISTRUZIONI PER L'USO. CATETERE A PALLONCINO PER PTA DILATAZIONE	05
FR	MODE D'EMPLOI. CATHÉTE DE DILATATION À BALLONNET POUR ATP	06
DE	GEBRAUCHSHINWEIS. BALLONKATHETER DILATATION FÜR PTA	07
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ. БАЛЛОННЫЙ ДИЛАТАЦИИ КАТЕТЕР ДЛЯ ПТА	08
CS	NÁVOD K POUŽITÍ. BALÓNKOVÝ DILATAČNÍ KATETR PRO PTA	09
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ . ΠΤΑ ΚΑΘΗΤΗΡΑΣ ΔΙΑΛΥΣΗ ΜΠΑΛΟΝΩΝ	10
PL	INSTRUKCJA UŻYTKOWANIA. BALONOWY CEWNIK ROZSZERZAJACY DO PTA	11
NL	GEBRUIKSAANWIJZING. BALLONKATHETER DILATIE VOOR PTA	12
SV	ANVÄNDNINGSIKTRUKTIONER. BALLONKATETER DILATION FÖR PTA	14
TR	KULLANIM İÇİN TALİMATLAR. PTA DİLAYON BALLOON KATETER	15
UK	ІНСТРУКЦІЇ З ВИКОРИСТАННЯ. БАЛООННИЙ КАТЕТЕР ДЛЯ ПТА	15

ENGLISH

1. Product description

The balloon catheter is a coaxial catheter (dual lumen from the connector to the tip, also known as OTW), designed for percutaneous transluminal angioplasty of small peripheral arteries.

The catheter has a coaxial dual-lumen body from the connector to the tip.

On the distal part of the catheter, before the tip, is located the balloon or inflatable segment, to different pressures. The balloon dilates the artery when is inflated by means of infusion of contrast medium inside it. The connector is Y-shaped, with two ports:

- The side port allows contrast medium passage to dilate the balloon.
- The straight port is intended for guidewire insertion. The guidewire is inserted through the tip of the catheter and comes out by the connector port to lead the catheter to the lesion.

Moreover, the inner lumen of the catheter is the one through which the guidewire runs, while the outer lumen is intended for contrast medium passage. The maximum diameter of the guidewire must not exceed 0.36mm=0.014".

The balloon is intended to reach different diameters at different pressures (according to the compliance curve included on the primary packaging).

To inflate the balloon, the connector inflation port must be connected to an inflation device. The balloon is inflated to predictable diameters with the pressure shown by the manometer.

One radiopaque marker is located at each end of the balloon in order to mark its length and help the user to see the catheter while inside the patient.

In the case of small diameter and short length balloons, a single marker is placed in the middle of the balloon to help know the position of the balloon when it is inside the patient.

At the distal end is the tip, which is rounded and atraumatic in shape in order to avoid damaging the arteries while it is being advanced. The distal tube of the catheter is coated with a durable hydrophilic coating to help the catheter go through tortuous arteries smoothly.

The useful length of the catheter is 100cm or 150cm depending on the item reference.

The balloon catheter is supplied in the following balloon diameters and lengths (pg. 19).

Contents

- One OTW catheter covered by protection sheath for the balloon and a stylet to protect the guidewire lumen. This is inserted into a dispenser in a sterile bag.
- One card with the compliance curve showing the nominal inflation pressure and the recommended maximum pressure.
- One leaflet with instructions for use.

2. Indications

The device is indicated to dilate stenosis in femoral, popliteal, infra-popliteal and renal arteries and for treating obstructive lesions of native or synthetic arteriovenous dialysis fistulae, with a reference diameter balloon from 1.25 and 4.0 mm and lengths from 10 up to 200 mm.

3. Contraindications

Inability to cross lesion with guidewire.

4. Warnings

- The device should only be used by experienced doctors familiar in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.

- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not dry with gauze.
- Do not expose the product to organic solvents or to contrast media not suitable for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Advance the product on the guidewire using fluoroscopy. Do not allow the product to advance without the guidewire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.

5. Precautions

- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting intravascular catheters.
- Use with standard 0.014" guidewires and standard introducers. If the introducer length goes above 25cm or the braided type is used, it may be necessary to increase its size to avoid catheter friction.
- Before use, check the size and suitability of the catheter for its intended use.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast medium (preferably 50/50, though the proportion of contrast medium can be reduced in large balloons, like 4mm diameter or longer ones).
- If any resistance is noticed during advance, stop and determine its cause before continuing.
- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guidewire and introducer be removed as a single unit.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.

6. Possible adverse effects/complications

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Ictus / embolism / thrombosis
- Haemodynamic deterioration
- Acute elastic recoil
- Arrhythmia
- Aneurysm or pseudo-aneurysm
- Infections
- Total occlusion of the artery
- Reocclusion of the treated area: restenosis
- Spasm
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Local bleeding with haematoma at the access site
- Hypo/hypertension
- Allergic reactions to materials.

7. How to use

7.1 Equipment Required

- Heparinized normal saline solution.
- Contrast medium (use contrast media suitable for intravascular use).
- 0.014" guidewire; do not use any other size of guidewire.
- Introducer with a haemostatic valve of the size specified on the label. Do

not use smaller sizes as these might damage the catheter.

- Three-way stopcock.
- Line extensions.
- Inflation device (manual pump with built-in manometer).
- Several standard 10-20cc syringes with saline solution for washing the system.

7.2 Preparing the catheter

- Remove the catheter from the protective dispenser. Check that it is the right size. Gently remove the protective sheath and stylet from the balloon.
- Connect a 10ml syringe containing sterile saline solution to the port for the guidewire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the whole air system:

1. Attach a three-way stopcock to the contrast liquid port (the side part of the connector). An extension line could be optionally placed between the balloon connector and the three-way stopcock.
2. Close the air passage through the balloon.
3. Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
4. Open the passage in the three-way stopcock between the syringe and the catheter.
5. With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
6. When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

DANGER: If bubbles do not stop coming from the catheter into the syringe after 3 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

- In order to activate the hydrophilic coating, it is recommended to wet the catheter with a sterile saline solution immediately prior to its insertion through the introducer sheath.

7.3 Technique for Insertion / Treatment

- Position the 0.014" guidewire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.
- Insert the proximal end of the guidewire into the distal tip of the balloon catheter.
- Carefully advance the balloon catheter on the guidewire, through the introducer and with the aid of fluoroscopy, until it reaches the area to be treated. In the event that the balloon catheter selected cannot be placed inside the lesion, take a smaller size to predilate the lesion.
- Connect the inflation device (1/2-1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Apply pressure for 30-60 seconds depending on the lesion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guidewire in position, carefully retract the catheter until it can be removed through the introducer.
- Inject contrast liquid through the introducer to check the treated lesion before removing the guidewire.

8. Warranty

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

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© iVascular - MP45021 Ed. 0 Dec 2017