Synthesis Preloaded

PRELOADED ABERRATION CUTTING EDGE
FREE ASPHERIC LENS
Model SIP and SIPY (yellow)

EN INTRAOCULAR LENS

DESCRIPTION:
The Synthesis Intraocular lens (IOL) is a single-piece ultraviolet (UV) absorbing posterior chamber IOL developed to replace the natural crystalline lens in pseudophakic adult patients. The Synthesis IOL has aspheric optics designed to be free of spherical aberration. The preloaded intraocular Synthesis lens is supplied with an Accusject Pro delivery System to provide a sterile, controlled and touch-free method of delivery of the lens to the ocular bag. The preloaded Synthesis IOL is supplied in 2 different models: Model-SIP or SIPY. See label on the cardboard box for type of the lens.

PHYSICAL CHARACTERISTICS:	
Lens / Haptic Material	Model SIP: Foldable hydrophilic acrylic copolymer containing a UV absorber. Model SIPY: Foldable hydrophilic acrylic copolymer containing a UV absorber and a blue light filtering chromophore.
Material Characteristics	Index of Refraction when wet at 20° C : 1.460 Index of Refraction when in the eye at 35° C : 1.459
Spectral Transmittance	Model SIP: 10% transmittance at 365 nanometers for +20.0 diopter IOL. Model SIPY: 10% transmittance at 374 nanometers for +20.0 diopter IOL. See FIG.7 with chart of X-Value = Wavelength (nm) and Y-Value = % Transmittance
Optic Type	Aspheric aberration free
Powers	0.0 to +30.0 Diopters (0.0 to +10.0 in 1.0 Diopter increments, +10.0 to +30.0 in 0.5 Diopter increments)
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 10.5 mm in 11 mm across the diopter range; Haptic Angle: 0°
Lens Orientation	The anterior side of the lens is facing up towards the anterior side of the eye. As illustrated in FIG.6, when the orientation features are top right (A) and bottom left (B), you are facing the anterior side of the lens.

INDICATIONS:
The Synthesis posterior chamber lens is indicated for primary implantation for the visual correction of aphakia in capsular bag patients. The lens is designed for implantation in the capsular bag.

CONTRA-INDICATIONS AND WARNINGS:
Physicians considering lens implantation in any of the following circumstances (but not limited to) must weigh the potential risk/benefit ratio, recurrent severe anterior or posterior segment inflammation or uveitis; patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases; surgical difficulties increasing the potential for complications (such as but not limited to: persistent bleeding, significant iris trauma, uncontrolled positive pressure, or persistent vitreous prolapse or loss); distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible; circumstances that would result in damage to the endothelium during implantation; patients in whom neither the posterior capsule nor zonule are intact enough to provide support; asymptomatic children under the age of 2 years are not suitable candidates for IOL.
As with any surgical procedure, there is risk involved. Potential complications of cataract surgery may include, but are not limited to the following: inflammation (iridocyclitis), pupillary membrane, vitreous inflammation (CME, TASS...), corneal endothelial damage, endophthalmitis, retinal detachment, corneal edema, pupillary block, iris prolapse, hypotony, glaucoma, wound leak, hyphema, iris atrophy, vitreous prolapse, IOL decentration or tilt, IOL silicone adhesions, IOL optic decentering or opacities, secondary surgical intervention (such as but not limited to wound leak repair, vitreous aspiration, iridectomy, IOL repositioning or replacement and retinal detachment repair).
The IOL must be implanted in accordance to the following instructions for use. Improper use may impose risk to the patient's health.
Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage.
Do not attempt to modify or alter the device or any of its components.
Do not use the lens after expiration date.
Do not reuse the lens and the preloaded system. Reuse may impose serious risk to the patient's health.
Do not resterilize as this can produce undesirable side effects.

PRECAUTIONS:

- Store at room temperature. Avoid temperatures below 5° C and above 45° C. Do not freeze.
- The IOL model, its power and expiration date should be verified before opening the protective packaging and before opening the individual sterile pouch.
- Use correct injector compatible.
- Sterility of the IOL and of the injector parts is guaranteed only if the individual sterile pouches have not been opened or damaged.
- IOLs must be carefully rinsed prior to implantation. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution
- Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens after implantation in the capsular bag.

PACKAGING / STERILIZATION:
The Preloaded Synthesis IOL is individually packaged in a pouch and a vital containing balanced salt solution that should be opened under sterile conditions. The Synthesis IOL is placed in a loading chamber, fixed by a lens holder, all immersed in BSS. A patient card and self-adhesive labels are supplied in the carton to provide traceability of the IOL. The Preloaded Synthesis IOL is steam sterilized.

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DIRECTIONS FOR USE:

- The physician should determine preoperatively the power of the lens to be implanted, using calculation methods. The recommended A-constants listed on the lens carton are guidelines only.
- Prior to implanting, examine the lens package for model, power, and proper configuration.
- Open the carton and remove the sterilized pouch containing the lens vial. Gradually peel the pouch apart to release the vial into the sterile field. Inspect the vial. Make sure it is not damaged and the seal is not broken. Hold the vial in one hand. Grasp the pull-tab of the foil lid and peel the foil lid away from you while holding the vial in a horizontal position to expose the loading chamber. Do not touch the lens at this time. Please refer to the enclosed figures to place the loading chamber into the compatible injector
- Remove the loading chamber containing the lens from the lens vial by pulling up on the upper loop of the lens holder. Do not remove the lens holder at this time.
- Prior to inserting the loading chamber into the injector ensure that the plunger is in withdrawn position. Insert the loading chamber containing the lens into the opening of the compatible injector body until it locks into position (FIG. 1).
- Rinse the lens with sterile balanced salt solution.
- Taking care not to touch the IOL, apply viscoelastic through the filling hole of the lens holder in the direction of the cartridge tip (FIG. 2) and at the end of the loading chamber (FIG. 3). Ensure that the viscoelastic fills the injector tip and the distal half way under the lens.
- Push the plunger forward until the silicone cushion approaches the lens holder.
- Remove the lens holder from the loading chamber by grasping the loop of the holder with two fingers and by flipping up the lower hinge of the lens holder with the thumb and lifting straight up (FIG. 4).
- Check that the lens is positioned centrally in the loading chamber and close the wings of the loading chamber (FIG. 5).
- When the "click-lock" mechanism is engaged, the lens is securely in place. The system is ready for injection. For detailed descriptions, refer to the instructions for use supplied with the compatible injector.
- Press the injector plunger forward to push the lens. Implant the lens with the bevel of the injector tip facing down to the patient's eye. The anterior side of the lens is facing up towards the cornea when the orientation features are top right (A) and bottom left (B) (FIG. 6)
- Only viscoelastics and delivery systems that have been validated and approved for use with this lens should be used (NOTE: Please refer to the Directions for Use with the viscoelastic and compatible injector system for additional information).
- There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient.

SYMBOLS USED ON LABELING:

Symbol	Description
D	Diopter
ØB	Body Diameter
ØT	Overall Diameter

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