 IMMULITE®
2000

HBsAg

**For the Qualitative Detection
of Hepatitis B Surface Antigen
in Human Serum**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 HBsAg

English

Intended Use For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the qualitative detection of *hepatitis B surface antigen* (HBsAg) in human serum or plasma (heparinized, sodium citrate or EDTA), as an aid in the determination of acute and chronic hepatitis B virus infections.

Catalog Number: **L2KHB2** (200 tests)

Test Code: **HBS** Color: **Orange**

Summary and Explanation

Hepatitis B virus (HBV) is the sole human pathogen in the family of hepatitis-associated DNA viruses, and is found world-wide. Distribution of HBV infection will vary among geographical areas and population groups. Transmission of the virus is due to parenteral contact, through the exchange of blood or blood products, sexual contact, and perinatal spread from mother to newborn.^{1,2} Clinical manifestations range from mild asymptomatic infections to severe fulminant hepatitis.^{1,2,3} Over 90% of infected adults will have an acute self-limiting infection,¹ with jaundice and abnormal liver function. Recovery occurs without any chronic sequelae.^{1,2}

Chronic liver disease, a condition in which infection persists for more than six months, a known sequela of a hepatitis B infection, is usually progressive.^{1,2} The risk of developing the chronic carrier state is more likely to follow infection acquired in childhood than as an adult.^{4,5} In chronic HBV carriers, there is no evidence of continued hepatic damage,^{1,2} however, the infection persists and the carrier maintains the ability to transmit the virus.²

Availability of recombinant HBV vaccines, and the recommendation of universal immunization for infants and other high-risk persons has aided in the prevention of HBV infections. In addition, treatment with alpha-interferon to relieve symptoms is available. Results have shown positive response to treatment in

40–50% of selected individuals with chronic active hepatitis B.^{4,5}

Classification of a hepatitis B infection requires the identification of several serological markers expressed during three phases (incubation, acute and convalescent) of the infection. The first marker to appear in serum is hepatitis B surface antigen (HBsAg). Presence of this antigen indicates an ongoing infection with HBV, and is detectable in the acutely ill and in chronic carriers.^{1,2,4}

Principle of the Procedure

IMMULITE 2000 HBsAg is a solid-phase, two-step chemiluminescent enzyme immunoassay. The solid phase, a polystyrene bead enclosed within a Reaction Tube, is coated with an antibody directed against the hepatitis B surface antigen (anti-HBs).

The patient specimen is added to the Reaction Tube containing a coated bead. An alkaline phosphatase-labeled anti-HBs antibody is also added to the Reaction Tube. After the wash and incubation steps, chemiluminescent substrate undergoes hydrolysis in the presence of alkaline phosphatase. The photon output, as measured by the luminometer, is related to the presence of HBsAg in the sample.

Incubation Cycles: 2 × 30 minutes

Time to First Result: 65 minutes

Specimen Collection

Lipemia may interfere with the assay. An ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving

anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 HBsAg has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 100 μ L serum or plasma (heparinized, sodium citrate or EDTA).

Storage: 2 days at room temperature (15–28°C).⁸

3 days at 2–8°C.^{6,8}

For longer storage: at –20°C.⁷

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.⁸⁻¹⁰

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

The HBsAg Adjustor, HBsAg Low Positive and HBsAg Positive Controls contain HBsAg which has been inactivated by proven, documented methods. However, always handle all controls as if capable of transmitting infectious agents.

The HBsAg results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 HBsAg EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

HBsAg Bead Pack (L2HB12)

With barcode. 200 beads coated with murine monoclonal anti-HBs. Stable at 2–8°C until expiration date.

L2KHB2: 1 pack

HBsAg Reagent Wedge (L2HBA2)

With barcodes. 11.5 mL of a protein-based buffer, with preservative. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-HBs in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KHB2: 1 wedge

Before use, tear off the top of the label at the perforations without damaging the barcode on the main label. Remove the foil seal from the top of the Reagent Wedge, and snap the sliding cover down into the ramps on the reagent lid.

HBsAg Adjustor (LHBR)

4 mL containing human serum with inactivated HBsAg, with preservative. The Adjustor serves as the assay's Cutoff. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHB2: 1 vial

HBsAg Controls (LHBC1, LHBC2, LHBC3)

Three vials (Negative, Low Positive and Positive) containing 8 mL each. **LHBC1 (Negative Control):** human serum without HBsAg, with preservative. **LHBC2, LHBC3 (Low Positive Control, Positive Control):** human serum with inactivated HBsAg, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHB2: 1 set

For the current control ratio ranges, please refer to the Control insert.

Aliquot Labels with barcodes are supplied with the kit, for use with the Adjustor and Controls. Before use, place the appropriate Aliquot Labels on standard test tubes, so the barcodes can be read by the barcode reader on the IMMULITE 2000.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Washing

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Distilled or deionized water; test tubes

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Adjustment Interval: 4 weeks

Quality Control Samples: The control(s) supplied with the kit should be used as quality control material to monitor the performance of the assay.

For the current control ratio ranges, please refer to the Control insert.

Calculation of Cutoff and S/CO Ratio:

The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second (mean cps) of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the "Low Adjustor CPS" and "Curve Parameter 1" fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry: Kit Entry.)

Calculation of a signal/cutoff (s/co) ratio is done by using the following formula:

$$\text{S/CO Ratio} = \frac{\text{Sample or Control cps}}{\text{Mean Adjustor cps} \times \text{P1}}$$

Calculation and reporting of qualitative (reactive/nonreactive) and s/co ratio results are handled automatically by the IMMULITE 2000.

The result is "Reactive" if the sample's counts are above the cutoff, and "Nonreactive" if below the cutoff.

Interpretation of Results

- A result of "**Reactive**" (ratio of ≥ 1.0) indicates that HBsAg is present and was detected in the patient sample. Specimens found to be initially reactive for HBsAg should be reassayed in duplicate to verify that the initially reactive result is repeatable.

If at least two of three result ratios are ≥ 100 for an individual sample, the patient sample may be reported as reactive and does not require testing with the IMMULITE 2000 HBsAg Confirmatory Assay (L2KCH).

If at least two of three result ratios are reactive, with ratios < 100 , the patient sample must be tested using the IMMULITE 2000 HBsAg Confirmatory Assay (L2KCH), a procedure to confirm the presence of HBsAg based on the neutralization of HBsAg-reactive samples with anti-HBs. Those samples with initial and repeated ratios ≥ 1.0 and < 100 ,

in which the HBsAg is neutralized by the confirmatory test procedure, are considered confirmed reactive for HBsAg. It is essential that results be confirmed before patients are informed of the test results.

If neither of the duplicates of the reassayed sample is reactive, the patient sample should be considered nonreactive for HBsAg.

- A result of “**Nonreactive**” (ratio of < 1.0) indicates that HBsAg is not present and was not detected in the patient sample.
- When HBsAg is used as a stand-alone assay for screening pregnant women to identify neonates who are at risk for acquiring HBV during the perinatal period, supplemental testing such as the IMMULITE 2000 HBsAg Confirmatory Kit must be used to confirm the result regardless of whether the reactive ratio results are above or below 100.
- Since all highly sensitive immunoassays may produce false reactive results in rare cases, it is essential that results be confirmed before patients are informed of the test results.

The results reported by the laboratory to the physician should include: “The following results were obtained with the IMMULITE 2000 HBsAg EIA. Results obtained from other manufacturers’ assay methods may not be used interchangeably.”

Expected Values

Individuals acutely infected with the hepatitis B virus will exhibit a rapid rise and decline of detectable levels of HBsAg between one and four months after exposure; usually through the course of the clinical illness. Of acutely infected adults, 5–10% will become chronic carriers of the virus.^{1,4} In chronic carriers, no significant decrease in HBsAg level is observed.

Limitations

The results of the test must be taken within the context of the patient’s clinical history, symptomology and other laboratory findings.

Lipemic, hemolyzed or grossly contaminated samples may give erroneous results.

A nonreactive result does not indicate that the patient was not infected with HBV. The patient sample should be tested for the presence of other serological markers.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay’s performance. Results are expressed as a signal-to-cutoff ratio. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See “Precision” table.)

Bilirubin: Presence of bilirubin in concentrations up to 400 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Hemolysis may cause a shift in values. (See “Hemolysis” table.)

Lipemia: Lipemia may cause a shift in values. (See “Lipemia” table.)

Alternate Sample Type: Samples (n = 60) were collected into plain, heparinized, sodium citrate and EDTA vacutainer tubes. 41 samples were spiked with HBsAg and all samples were then assayed by the IMMULITE 2000 HBsAg

procedure. Results are expressed as a signal-to-cutoff ratio.

(Heparin) = 0.99 (Serum) + 1.93
r = 1.00

(NaCitrate) = 0.96 (Serum) + 0.83
r = 1.00

(EDTA) = 0.98 (Serum) + 1.52
r = 1.00

Means:

219 (Serum)
218 (Heparin)
211 (NaCitrate)
216 (EDTA)

Method Comparison 1: The IMMULITE 2000 HBsAg procedure was compared to the Abbott AUSZYME® MONOCLONAL, a commercially available enzyme immunoassay for HBsAg, on 231 patient samples.

IML 2000				
Abbott	React.	Non-react.	Relative Sensitivity	Relative Specificity
Reactive	16	0		
Non-reactive	3	212	100%	99%

Agreement: 99%

Method Comparison 2: The IMMULITE 2000 HBsAg procedure was compared to the Elecsys HBsAg immunoassay on 2438 samples, consisting of seroconversion panels, potentially cross-reactive samples, samples from patients with acute and chronic HBV infection, from blood donors, hospitalized patients, and routine samples. The table below presents the results of this study.

	No. of samples					
	HBsAg reactive	HBsAg non-reactive	Sensitivity	Specificity	PPV	NPV
IMMULITE 2000 HBsAg						
Reactive	238	0				
Non-reactive	16	2184	93.7%	100%	100%	99.3%
Elecsys HBsAg						
Reactive	254	0				
Non-reactive	0	2184	100%	100%	100%	100%

Method Comparison 3: The IMMULITE 2000 HBsAg procedure was compared to the AxSYM HBsAg immunoassay on 438 samples, consisting of seroconversion panels and samples from patients with chronic HBV infection. The table below presents the results of this study.

IMMULITE 2000	React.	React.	Non-Non-		Agreement
			react.	react.	
AxSYM	React.	react.	React.	react.	
Seroconversion (n = 340)	148	3	7	182	97.1%
Chronic HBV infection (n = 98)	98	0	0	0	100%
Total (n = 438)	246	3	7	182	97.7%

References

- 1) Locarnini SA, Gust ID. Hepadnaviridae: hepatitis B virus and the delta virus. In: Balows A, et al, editors. Laboratory diagnosis of infectious diseases: principles and practices. New York: Springer-Verlag, 1988: 750-96.
- 2) Follett EAC. Diagnosis of hepatitis B infection. In: Young H, McMillan A, editors. Immunological diagnosis of sexually transmitted diseases. New York: Marcel Dekker, 1988: 433-49.
- 3) Hollinger FB, Dienstag JL. Hepatitis B and D viruses. In: Lennette EH, et al, editors. Manual of clinical microbiology. 6th ed. Washington, D.C.: American Society for Microbiology, 1995: 1033-49.
- 4) Nowicki MJ, Balistreri WF. Hepatitis A to E: building up the alphabet. Contemporary Peds 1992: 118-28.
- 5) Zuckerman AJ, et al. Hepatitis B virus and hepatitis D virus. In: Principles and practice of clinical virology. 2nd ed. New York: John Wiley & Sons, 1992: 153-72.
- 6) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:354-60.
- 7) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders, 1995:322-4.
- 8) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.
- 9) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- 10) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only.
For technical assistance, contact your
National Distributor.

siemens.com/healthcare

The Quality System of Siemens Healthcare
Diagnostics Products Ltd. is certified to ISO
13485.

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2019 Siemens Healthcare Diagnostics. All rights reserved.