
VERTEBROPLASTIC

EN	Instructions for use	p.2
IT	Istruzioni per l'uso	p.3
PT	Instruções de uso	p.4
FR	Mode d'emploi	p.5

ES	Instrucciones de uso	p.6
DE	Gebrauchshinweise	p.8
RO	Instrucțiuni de utilizare	p.9

INSTRUCTIONS FOR USE

- Before using the needle check that the package is not damaged. Check that the needle is not damaged and if necessary inform the manufacturer or your local representative. Remove the needle from its package, being careful to operate under aseptic conditions.
- Position the patient. The whole procedure must be followed in real time under the guide of combined CT and fluoroscopy. Even a digital angiography may be useful. For the skin and peripheral tissues a local anesthetic injection could be useful;
- Perform the necessary procedures of disinfection of the skin.
- Insert the needle in the selected point, holding the handle on the palm of your hand; advance the entire needle with a slight twisting motion until reaching the vertebral surface.
- If the needle Chibell for deep anesthesia is present in the kit: remove the stylet and insert the needle Chibell in order to perform local anesthesia in the periosteum, connecting an appropriate syringe in the luer lock connection. Remove the needle for anesthesia and reinsert the stylet into the cannula.
- Advance the needle with a slight twisting motion through the cortex. Enter with the needle into the vertebral body and remove the stylet.
- Before injecting the cement (polymethylmethacrylate, PMMA), there must be a venography to exclude needle placement directly in the venous vertebral complex. It is also necessary to verify the continuity of the posterior vertebral wall by the injection of a contrast medium.
- If spatula and bowl are available in the kit, prepare the cement according to the manufacturer instructions of use, put it into an appropriate syringe or aspirate it using the suction nozzle (which may be present in the kit) firmly connected to the syringe.
- If a mixing system called 'CEMIX' is available in the kit, prepare the cement according to its instructions, and fill 'CEMIX' using a disposable syringe, or using the device 'Easynject' of Biopsybell.
- Inject the necessary amount of cement using a syringe and exerting adequate pressure (or following the Instructions for use of the device Easynject). If necessary you can use a suitable extension to be connected to the syringe.
- Ensure uniform filling of the vertebrae. You should check that there is no leakage of cement outside the predetermined space for ex. in the veins, discs or epidural spaces.
- Complete the extrusion of cement using the dedicated extractor present in the kit.
- After injecting the PMMA reinsert the stylet into the cannula, make sure that the two connections (cannula and stylet) are fitted to each other and finally extract the needle. During the introduction of the stylet do not use a hammer only on the handle of such a stylet, in order not to compromise its integrity.

INDICATIONS

The 'VERTEBROPLASTIC' is intended for vertebroplasty. Use only under medical surveillance.

PRECAUTIONS

The system is reserved exclusively to physicians. After use, the needle may be a risk of contamination and / or infection. So handle with care and dispose of the product according to applicable regulations and medical practice. In order to minimize risks or other complications associated with the use of cement for vertebroplasty is recommended that :

- In all patients having acute osteoporotic fractures of the vertebral body should be considered a period of traditional therapy;
- Only qualified and properly trained physicians carrying out of vertebroplasty should practice this procedure;
- Use only cements indicated for vertebroplasty and read and follow carefully the instructions for the use of cement;
- Monitor the process with high quality imaging systems to enable the recognition of spills or leaks of PMMA
- Keep closely monitored the blood pressure of the patient during and immediately after the procedure, a treatment level multiple may increase the risk of sudden drop in blood pressure due to release of the monomer of PMMA in the circulation. In a single session should not be made treatment more than 3 vertebrae together;
- A careful diagnosis and special precautions should be taken whenever the procedure is performed in the treatment of patients with spinal tumors that have eroded the posterior wall of the vertebral body.
- It's also recommended not to use the device if are necessary particularly risky maneuvers when the device should be subjected to high bending or twisting forces that could cause the breaking of the needle.

PACKAGE

The system is packed in rigid transparent blister and medical paper and it is sterilized with ethylene oxide (EO) in accordance with the requirements of ISO 11135. The device is sold in packs of 5 or 10 pieces. The device is not reusable. In case of reuse there may be a risk of deterioration in device performance and / or the possibility of cross contamination. Do not use if packaging has been opened or damaged or polluted.

Biopsybell disclaims any liability for any damage caused by improper use or different from that reported in this instruction sheet.

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