

SPHERA DUO

HYDROCEPHALUS SHUNT SYSTEM



SPHERA DUO

HYDROCEPHALUS SHUNT SYSTEM



HYDROCEPHALUS SHUNT SYSTEM SPHERA DUO

Sphera Duo valves were developed for a precise control of cerebral intraventricular pressure. The Adult, Infantile and Neonatal Plus sizes design present low profile of implant and flexible body with anatomical shape to contour the curvature of the cranium.

The valves are made of transparent medical grade silicone with internal structure in polysulfone. They have pumpable central reservoir with needle guard that protects against perforations on puncturing or sampling procedures.


The Neonatal size has no pumping chamber, as the design prioritizes low profile and minimum implant volume. The system can be combined to a reservoir connected to the ventricular catheter to enable puncturing and CSF sampling.

SPHERA MECHANISM: PRECISE MANAGEMENT

The pressure control system is composed of ruby ball and conic seat and stainless steel spring. The perfect match between ball and seat can safely set the pressures of opening and closing of the system, providing precise control of intracranial pressure.

Adult, Infantile and Neonatal Plus sizes have a dual pressure control encased in the input and output occluders. The Neonatal model encases a single mechanism for pressure control.


All sizes are supplied in four pressure ranges: high, medium, low and extra low, to meet the individual requirements of patients. The indications of flow and pressure printed on the valve body are radiopaque and allow visualization imaging after implantation.



CATHETER: FLEXIBILITY AND RADIOPACITY

The shunt system presents the valve accompanied by cerebral ventricular catheter and peritoneal catheter. These are made of soft transparent medical grade silicone with radiopaque stripe, which ensures catheter visualization in imaging exams.


The hardness of the silicone used in the manufacture of catheters was tailored to allow adequate flexibility and at the same time, prevent the occurrence of unwanted kink in the subcutaneous route, which can cause obstruction or decrease the flow of drainage. For all valve sizes, the system can be supplied in different configurations of peritoneal and ventricular catheter with or without reservoir.



The graphic represents average rates. Consider a range $\pm 1,5$ cmH₂O

MARKS		
Pressure	Marks	Values (Flow of 21ml/h)
Extra low	○ ○ ○ ▶	1 to 3 cmH ₂ O
Low	● ○ ○ ▶	3 to 7 cmH ₂ O
Medium	● ● ○ ▶	7 to 11 cmH ₂ O
High	● ● ● ▶	11 to 14 cmH ₂ O

	Model	Pressure	Content
ADULT	ADE10215D ADB10215D ADM10215D ADA10215D	Extralow Low Medium High	01 Hydrocephalus Valve Adult size 01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide 01 Peritoneal Catheter 102 cm 01 90° Angle Former, 01 irrigation tubing
	ADE10223D ADB10223D ADM10223D ADA10223D	Extralow Low Medium High	01 Hydrocephalus Valve Adult size 01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide 01 Peritoneal Catheter 102 cm 01 90° Angle Former, 01 irrigation tubing
INFANTILE	INE9015D INB9015D INM9015D INA9015D	Extralow Low Medium High	01 Hydrocephalus Valve Infantile size 01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide 01 Peritoneal Catheter 90 cm 01 90° Angle Former, 01 irrigation tubing
	INE9023D INB9023D INM9023D INA9023D	Extralow Low Medium High	01 Hydrocephalus Valve Infantile size 01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide 01 Peritoneal Catheter 90 cm 01 90° Angle Former, 01 irrigation tubing
NEONATAL PLUS	NPE9015D NPB9015D NPM9015D NPA9015D	Extralow Low Medium High	01 Hydrocephalus Valve Neonatal Plus size 01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide 01 Peritoneal Catheter 90 cm 01 Angle former, 01 irrigation tubing
NEONATAL	NNE9015D NNB9015D NNM9015D NNA9015D	Extralow Low Medium High	01 Hydrocephalus Valve Neonatal size 01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide 01 Peritoneal Catheter 90 cm 01 90° Angle Former, 01 irrigation tubing
	NNE9015RD NNB9015RD NNM9015RD NNA9015RD	Extralow Low Medium High	01 Hydrocephalus Valve Neonatal size with reservoir 01 Straight Ventricular Cerebral Catheter 15 cm and stainless steel guide 01 Peritoneal Catheter 90 cm; 01 90° Angle Former, 01 irrigation tubing



The technical information for these products is not limited to the characteristics presented in this catalog. For complete information, request the Instructions for Use at info@hpbio.com.br



HP BIOPROTESES LTDA
 Rua Maria José Rangel, 83 • Vila São Paulo • CEP: 04650-180
 São Paulo - SP - Brasil • Tel.: 55 11 3853 7625 • Fax: 55 11 3853 7628
info@hpbio.com.br • www.hpbio.com.br