CEA

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 CEA

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of carcinoembryonic antigen (CEA) in serum, as an aid in the management of cancer patients and in the assessment of prognosis.

Catalog Number: **L2KCE2** (200 tests), **L2KCE6** (600 tests)

Test Code: **CEA** Color: **Light Green**

Summary and Explanation

Carcinoembryonic antigen (CEA) comprises a heterogeneous family of glycoproteins ranging in weight from 175,000 to 200,000 daltons due to varying carbohydrate and amino acid content. Its biological function is not well defined, but it may play a role in intercellular recognition, regulation of the immune response, and metastasis of colorectal cancer.⁶ The name derives from the previous belief that CEA occurred only in gastrointestinal carcinomas and the fetal digestive tract. Elevated levels have since been detected in a number of malignant and nonmalignant conditions of the gastrointestinal tract and other sites. These conditions include various hepatic diseases; inflammatory lesions, particularly of the gastrointestinal tract; infections; trauma; infarction; collagen vascular disease; renal impairment; and smoking (current and past).² CEA also occurs at low levels in normal colon and other tissues. Serum values in healthy adults are generally less than 5 ng/mL.⁷ Though serum values exceeding five times the normal range usually indicate malignancy, values seen in malignant and nonmalignant conditions can overlap considerably, ruling out the use of CEA as a screen for malignancy. The value of CEA measurement lies, instead, in patient prognosis, status assessment, and monitoring.

Until 1980, the only assay method for CEA used in the United States was a radioimmunoassay requiring extraction. Since then, several other nonextraction RIAs and sandwich-type EIAs have been developed for use on serum and plasma samples.

CEA levels at diagnosis of colorectal carcinoma (CRC) correlate with prognosis.⁵ Elevated preoperative serum CEA predicts increased risk of recurrence⁶ or hepatic metastases.⁴ The site of the primary cancer may influence interpretation of elevated CEA: in colon carcinomas, high CEA levels indicate a poor prognosis, but not necessarily so in rectal carcinomas.⁹ In addition, monitoring levels during chemotherapy before surgery can be informative, and failure of CEA to fall during preoperative radio-therapy usually indicates the presence of a tumor outside the field of radiation and a poorer prognosis.³ Preoperative elevated CEA has also been observed to correlate inversely with estimated mean time to recurrence,⁴ and to correlate directly with degree of tumor stage, thickness, and differentiation in CRC patients at diagnosis.^{3,4,9} CEA is not always elevated, however, in sera of those with poorly differentiated CRCs, which produce less CEA but which are more aggressive than moderately or well-differentiated carcinomas.^{3,6} Rising CEA is deemed the most accurate indicator of recurrence in CRC.^{5,8} Levels decrease to normal in nearly all patients after complete resection of CRCs, usually within 4 to 6 weeks after surgery.³ Failure to decline after surgery may suggest incomplete resection,⁵ and increasing levels can precede recurrent CRCs an average of 4 to 6 months before they are clinically evident in about two-thirds of cases.^{3,5,8} Rising CEA levels have also been used to assist in selecting candidates for second-look surgery, in conjunction with antibody imaging using radiolabeled antibodies against CEA to provide information on tumor location.³ Several studies have suggested a correlation between the rate of rise in CEA and the presence or probability of hepatic metastases: more rapid rises are seen with hepatic metastases than with

localized recurrences and resectable localized tumors. In some cases, it has been CEA that signaled the asymptomatic recurrences.⁴

Principle of the Procedure

IMMULITE 2000 CEA is a solid-phase, two-site sequential chemiluminescent immunometric assay.

Incubation Cycles: 2 × 30 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 CEA has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 15 µL serum

Storage: 7 days at 2–8°C or for longer periods at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹⁵⁻¹⁷

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H318,
P280, P264,
P305 +
P351 +
P310**

Danger! Causes serious eye damage.

Wear protective gloves/protective clothing/eye protection/face protection. Wash hands thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Immediately call a POISON CENTER or doctor/physician. **Contains:** tetrasodium ethylene diamine tetraacetate; CEA Reagent Wedge A

Harmful to aquatic life with long lasting effects.

**H412
P273, P501**

Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. **Contains:** sodium azide; CEA Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

CEA Bead Pack (L2CE12)

With barcode. 200 beads, coated with monoclonal murine anti-CEA antibody. Stable at 2–8°C until expiration date.

L2KCE2: 1 pack **L2KCE6:** 3 packs

CEA Reagent Wedge (L2CEA2)

With barcode. 11.5 mL buffer/murine serum matrix, with preservative; 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal rabbit anti-CEA antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KCE2: 1 wedge **L2KCE6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

CEA Adjustors (LCEL, LCEH)

Two vials (Low and High) of lyophilized CEA in a CEA-free human serum, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Let stand for 30 minutes. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 7 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KCE2: 1 set **L2KCE6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of CEA.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Based on its relationship to IMMULITE CEA (see Method Comparison), the IMMULITE 2000 CEA procedure can be expected to have essentially the same reference ranges.

A reference range study for IMMULITE CEA was performed using serum samples from adult volunteers, including both men and nonpregnant women ranging from approximately 20 to 70 years of age (central 95%: 22–64 years, median: 40 years). The subjects were in apparent good health, based on a questionnaire.

Blood samples were collected in France, Germany, the Netherlands and Portugal. Results were generated by an independent laboratory in the Netherlands using IMMULITE kits. (The samples were collected in plain glass tubes without anticoagulants, gel barriers, or clot-promoting additives, and assayed in singlicate.)

Medians and 95th percentiles for relevant subgroups are tabulated below. Percentiles were determined nonparametrically.

Group	Median	95%ile	Units	n
Male Smokers	2.1	6.2	ng/mL	153
Male Nonsmokers	1.1	3.4	ng/mL	226
Female Smokers	1.3	4.9	ng/mL	81
Female Nonsmokers	0.8	2.5	ng/mL	262

Results from another study with the IMMULITE CEA assay are shown in the two tables below. The study was conducted at three sites in the United States. Included are individual results from 153 subjects in apparent good health (both smokers and nonsmokers); individual results from 243 patients with various nonmalignant diseases; and a total of 1382 samples (some collected during serial monitoring) from 657 patients with various types of cancer.

Patient Category	n	Median	95 %ile
Normal subjects			
Nonsmokers	86	1.2	4.1
Smokers	67	2.1	9.8
Nonmalignant Diseases			
Pulmonary Diseases	32	2.8	11.2
Renal Diseases	19	1.5	6.3
Hepatitis	47	1.9	8.5
Thyroid Diseases	65	1.1	4.9
Other Nonmalignant Diseases	80	1.9	13.6
Malignant Diseases			
Bladder Cancer	27	1.5	83.1
Breast Cancer	46	4.1	2230
Colorectal Cancer	944	5.1	849
Esophageal Cancer	44	4.5	323
Lung Cancer	52	4.8	355
Ovarian Cancer	50	1.3	24.6
Renal Cancer	39	1	3.5
Pancreatic Cancer	19	19.4	495
Stomach Cancer	25	56.5	277
Prostate Cancer	29	1.2	13.2
Rectal Cancer	21	1.4	23.8
Other Cancers	86	2	233

Patient Category	n	Percent of Samples in Range (ng/mL)			
Group		0–3.0	3.1–5.0	5.1–10	>10
Normal Subjects					
Nonsmokers	86	90.7%	7.0%	2.3%	0.0%
Smokers	67	71.6%	14.9%	9.0%	4.5%
Nonmalignant Diseases					
Pulmonary Diseases	32	62.5%	12.5%	18.8%	6.3%
Renal Diseases	19	73.7%	15.8%	10.5%	0.0%
Hepatitis	47	74.5%	10.6%	12.8%	2.1%
Thyroid Diseases	65	92.3%	3.1%	4.6%	0.0%
Other Nonmalignant Diseases	80	70.0%	12.5%	10.0%	7.5%
Malignant Diseases					
Bladder Cancer	27	70.4%	18.5%	3.7%	7.4%
Breast Cancer	46	43.5%	13.0%	10.9%	32.6%
Colorectal Cancer	944	36.1%	13.1%	15.8%	35.0%
Esophageal Cancer	44	34.1%	22.7%	9.1%	34.1%
Lung Cancer	52	38.5%	13.5%	7.7%	40.4%
Ovarian Cancer	50	86.0%	2.0%	6.0%	6.0%
Renal Cancer	39	89.7%	7.7%	0.0%	2.6%
Pancreatic Cancer	19	31.6%	5.3%	0.0%	63.2%
Stomach Cancer	25	36.0%	0.0%	8.0%	56.0%
Prostate Cancer	29	69.0%	24.1%	3.4%	3.4%
Rectal Cancer	21	71.4%	14.3%	0.0%	14.3%
Other Cancers	86	58.1%	12.8%	8.1%	20.9%

Consider these limits as *guidelines only*. Each laboratory should establish its own reference ranges.

Limitations

Patients with confirmed carcinoma often have pretreatment CEA levels in the same range as healthy persons. Elevated CEA levels are frequently observed in smokers, in cancer patients and in patients with a variety of nonmalignant diseases and inflammatory conditions. Therefore, a serum CEA level, regardless of its value, should not be interpreted as absolute evidence for the presence or absence of malignant disease. The CEA value should

be used in conjunction with information available from clinical evaluation and other diagnostic procedures. The IMMULITE 2000 CEA assay should not be used as a cancer screening test.

Some individuals have antibodies to mouse protein which can cause interference in immunoassays that employ antibodies derived from mice. In particular, specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). These specimens may show erroneous results in such assays.¹²⁻¹⁴ Therefore, results should be interpreted with caution for such patients.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1998;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: Up to 550 ng/mL

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 0.15 ng/mL

High-dose Hook Effect:

None up to 250,000 ng/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three CEA solutions (2520, 4160 and 5580 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibodies used in the IMMULITE 2000 CEA procedure are highly specific for CEA. (See "Specificity" table.)

Bilirubin: Presence of conjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin:

Biotin Test Level (ng/mL)	CEA Concentration (ng/mL)	
	3.6	210
	Bias (%)	
2.0	2	3
5.0	12	0
9.0	7	4
19.0	17	0
32.0	19	16
100	23	22
1500	25	23

Specimens that contain biotin at a concentration of 2 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely elevated results for patient samples.

The recommended adult daily dietary intake for biotin is 30 µg/day. Over the counter dietary supplements promoted for use in hair, skin and nail health may contain 5–100 mg of biotin, with recommendations to take multiple pills per day. Pharmacokinetic studies in healthy adults have shown that, in subjects ingesting 5 mg, 10 mg, and 20 mg of biotin, serum concentrations of biotin can reach up to 73 ng/mL, 141 ng/mL, and 355 ng/mL, respectively.¹⁸ Subjects who take up to 300 mg of biotin per day may

have plasma biotin levels as high as 1160 ng/mL.¹⁹

Hemolysis: Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 10 volunteers into plain, heparinized, EDTA and Becton Dickinson SST® vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of carcinoembryonic antigen, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 CEA procedure.

(Heparin) = 0.96 (Serum) – 0.45 ng/mL
r = 0.992

(EDTA) = 1.0 (Serum) – 0.91 ng/mL
r = 0.995

(SST) = 0.99 (Plain Tubes) + 0.97 ng/mL
r = 0.990

Means:

40.9 ng/mL (Serum)
38.9 ng/mL (Heparin)
40.0 ng/mL (EDTA)
41.4 ng/mL (SST)

Method Comparison: The assay was compared to IMMULITE CEA assay on 48 samples. (Concentration range: approximately 2.1 to 550 ng/mL. See graph.) By linear regression:

(IML 2000) = 1.06 (IML) + 5.0 ng/mL
r = 0.989

Means:

204 ng/mL (IMMULITE 2000)
187 ng/mL (IMMULITE)

References

- 1) Beard DB, Haskell CM. Carcinoembryonic antigen in breast cancer. Am J Med 1986 Feb;80:241-5.
- 2) Begent R. The value of carcinoembryonic antigen measurement in clinical practice. Ann Clin Biochem 1984;21: 231-8.
- 3) Begent R, Rustin GJS. Tumour markers: from carcinoembryonic antigen to products of hybridoma technology. Cancer Surv 1989;8(1):107-21.
- 4) Cooper EH, de Mello Jr. JP, Giles GR. Biochemical markers in gastrointestinal malignancies. Arq Gastroenterol 1989;26(4):131-40.
- 5) Fantini GA, DeCosse JJ.

Surveillance strategies after resection of carcinoma of the colon and rectum. Surg Gynecol Obstet 1990;171:267-73. 6) Jessup JM, Thomas P. Carcinoembryonic antigen: function in metastasis by human colorectal carcinoma. Cancer Metastasis Rev 1989;8:263-80. 7) Nisselbaum JS, Smith CA, Schwartz D, Schwartz MK. Comparison of Roche RIA, Roche EIA, Hybritech EIA, and Abbott EIA methods for measuring carcinoembryonic antigen. Clin Chem 1988;34(4):761-4. 8) Sardi A, Agnone CM, Nieroda CA, Mojzisek C, Hinkle G, Ferrara P, et al. Radioimmunoguided surgery in recurrent colorectal cancer: the role of carcinoembryonic antigen, computerized tomography, and physical examination. South Med J 1989;82:1235-9. 9) Sener SF, Imperato JP, Chmiel J, Fremgen A, Sylvester J. The use of cancer registry data to study preoperative carcinoembryonic antigen level as an indicator of survival in colorectal cancer. CA 1989;39(1):50-7. 10) Shinkai T, Saijo N, Tominaga K, Eguchi K, Shimizu E, Sasaki Y, et al. Serial plasma carcinoembryonic antigen measurement for monitoring patients with advanced lung cancer during chemotherapy. Cancer 1986;57:1318-23. 11) National Committee for Clinical Laboratory Standards, Procedures for the collection of diagnostic blood specimens by venipuncture, 3rd ed, NCCLS Document H3-A3, 1991. 12) Primus FJ et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine antibody for diagnosis and therapy. Clin Chem 1988;34:261-4. 13) Hansen HJ et al. Solving the problem of antibody interference in commercial "sandwich"-type immunoassay of carcinoembryonic antigen. Clin Chem 1989;35:146-51. 14) Schroff RJ et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-85. 15) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 16) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 17) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030. 18) Grimsey P, et al. 2017. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. Int. J. Pharmacokinet. 2(4), 247-256. 19) Piketty ML, et al. 2017. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. Clin Chem Lab Med. 55(6):817-825.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs


Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	6.0	0.18	3.0%	0.58	9.8%
2	15.8	0.58	3.7%	1.2	7.6%
3	54	1.48	2.7%	2.45	4.5%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	12	—	—
	8 in 16	5.8	6.0	97%
	4 in 16	2.8	3.0	93%
	2 in 16	1.4	1.5	93%
	1 in 16	0.7	0.8	88%
2	16 in 16	26	—	—
	8 in 16	13	13	100%
	4 in 16	6.3	6.5	97%
	2 in 16	3.1	3.3	94%
	1 in 16	1.5	1.6	94%
3	16 in 16	44	—	—
	8 in 16	22	22	100%
	4 in 16	11	11	100%
	2 in 16	5.2	5.5	95%
	1 in 16	2.4	2.8	85%
4	16 in 16	145	—	—
	8 in 16	75	73	103%
	4 in 16	38	36	106%
	2 in 16	19	18	106%
	1 in 16	9.4	9.1	103%
5	16 in 16	243	—	—
	8 in 16	130	122	107%
	4 in 16	62	61	102%
	2 in 16	30	30	100%
	1 in 16	14	15	93%

Method Comparison

(IML 2000) = 1.06 (IML) + 5.0 ng/mL
r = 0.989

Recovery (ng/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	23	—	—
	A	141	148	95%
	B	220	230	96%
	C	283	301	94%
2	—	73	—	—
	A	189	195	97%
	B	274	277	99%
	C	333	348	96%
3	—	85	—	—
	A	200	207	97%
	B	276	289	96%
	C	328	360	91%
4	—	117	—	—
	A	238	237	101%
	B	329	319	103%
	C	388	390	99%

Specificity

Compound ¹	ng/mL Added ²	% Cross reactivity ³
Chondroitin sulfate	270,000	ND
Dextran sulfate	50,000	ND
Hyaluronic acid	500,000	ND
Heparin sodium	6,000,000	ND
AFP	10,000	ND
Ferritin	10,000	ND
HCG	10,000	ND
PAP	1000	ND
PSA	1000	ND
Prolactin	500	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2017-2020 Siemens Healthcare Diagnostics. All rights reserved.
CONTROL AB