Coagulation Control Plasmas

helena
Biosciences Europe

REF 5186
REF 5187
REF 5183
REF 5482

Routine Control N
Routine Control A
Routine Control SA
Routine Coagulation Control Set

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-0482P 2016/01 (16)

Coagulation Control Plasmas

en

INTENDED PURPOSE

The Coagulation Control Plasmas kit is intended for use as a quality control material.
Routine Control N and Routine Control A and Routine Control SA are for use as normal, moderately prolonged and markedly prolonged controls for PT and aPTT assays. They are also assayed for Fibrinogen, TCT and ATIII, and are prepared from normal human plasma.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.
Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg), Hepatitis C Antigen (HcAg), HIV-1 antibody, HIV-2 antibody.
However they should be handled with the same precautions as a human patient sample.

COMPOSITION

| REF | Component | Content | Description |
|------|----------------------------------|-----------|--|
| 5186 | Routine Control - N | 10 x 1 mL | Prepared from pooled normal plasma. |
| 5187 | Routine Control - A | 10 x 1 mL | Prepared from asclerotic human plasma. |
| 5183 | Routine Control - SA | 10 x 1 mL | Prepared from asclerotic human plasma. |
| 5482 | Routine Coagulation Control Set: | | |
| | Routine Control - N | 4 x 1 mL | |
| | Routine Control - A | 3 x 1 mL | |
| | Routine Control - SA | 3 x 1 mL | |

Each kit contains instructions for use.

Each kit contains 1 mL of buffered, lyophilised human plasma.
Reconstitute each vial of the appropriate control with 1 mL of distilled or deionised water. Swirl gently. Allow to stand for 10 minutes for complete dissolution and mix well before use.

ITEMS REQUIRED BUT NOT PROVIDED

Coagulation Control Plasmas may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with suitable commercial reagents.

STORAGE SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. The reconstituted controls are stable for 8 hours when kept at 2 – 8°C or 4 weeks at -20°C when flash frozen. Keep covered.

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Not applicable.

PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS

Routine Control N should give values within the laboratory normal range for PT, aPTT and fibrinogen assays. Routine Control A and Routine Control SA have been standardised to give prolonged and markedly prolonged PT and aPTT times respectively. Lot and routine specific expected values are provided with each pack of controls.

LIMITATIONS

The results obtained with Coagulation Control Plasmas depend on several factors strongly associated with instrumentation. Types of reagent, reagent substrate and laboratory to laboratory variations. Each laboratory should establish an expected range for the particular instrument being used.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility

| Sample | n | Intra-assay precision aPTT CV (%) | PT CV (%) |
|--------------------|---|--------------------------------------|-----------|
| Routine Control N | 5 | 2/83 | 1.01 |
| Routine Control A | 5 | 2/78 | 1.71 |
| Routine Control SA | 5 | 1/72 | 1.03 |

BIBLIOGRAPHY

1. Kirkwood TBL, et al. (1977) Identification of Sources of Variation in Factor VIII Assay. *British Journal of Haematology*, 37:555-568.
2. Goldfarb MD (1971) Reproducibility in Coagulation Assays. *AJCP* 55:561-564.
3. Pallitt HA and Longbery JR (1979) A Precision Study of Coagulation Factor Assay Techniques. *AJCP* 59:231-235.

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BIBLIOGRAFIE

1. Marion E.Reid & Christine Lomas-Francis, Blood Group Antigens & Antibodies, SBB Books, New York 2007; pagina 181.
2. Issitt PD. Applied Blood Group Serology, 3rd Edition. Montgomery Scientific, Miami 1985; capitolul 6.
3. Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationary Office.
4. AABB Technical Manual, 16th edition, AABB 2008.
5. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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
Advena Ltd. Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013, Malta

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SUMMARY

In 1900, Landsteiner discovered the serum of some people would agglutinate the red cells of others. Four common phenotypes are now recognised: O, A, B and AB. Subgroups of A and B have since been identified.

| Forward Group | | | Reverse Group | | | ABO Phenotype | Caucasians % ¹ | |
|---------------|---|-----|----------------|----------------|---|---------------|---------------------------|----|
| A | B | A,B | A ₁ | A ₂ | B | O | | |
| + | 0 | + | 0 | 0 | + | 0 | A | 43 |
| 0 | + | + | + | + | 0 | 0 | B | 9 |
| 0 | 0 | 0 | + | + | + | 0 | O | 44 |
| + | + | + | 0 | 0 | 0 | 0 | AB | 4 |

INTENDED PURPOSE

The ABO reagents are blood grouping reagents intended to be used to qualitatively determine the presence or absence of the A and/or B antigens on the red cells of blood donors or patients requiring a blood transfusion when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagents contain antibodies against the appropriate A and/or B antigen on human red cells and will cause direct agglutination (clumping) of red cells that carry the corresponding ABO antigen. No agglutination generally indicates the absence of the corresponding ABO antigen on human red cells (see **Limitations**).

REAGENT

Lorne Monoclonal IgM ABO blood grouping reagents contain mouse monoclonal antibodies diluted in a phosphate buffer containing sodium chloride, EDTA and bovine albumin. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. Each reagent is supplied at optimal dilution for use with all the recommended techniques stated below without the need for further dilution or addition. For lot reference number and expiry date see **Vial Label**.

| Product | Cell Line/Clone | Colour | Dye Used |
|----------|-------------------------|------------|-------------|
| Anti-A | 9113D10 | Blue | Patent Blue |
| Anti-B | 9621A8 | Yellow | Iartrazine |
| Anti-A,B | 152D12 + 9113D10 + ES15 | Colourless | None |

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PRECAUTIONS

- The reagents are intended for in vitro diagnostic use only.
- If a reagent vial is cracked or leaking, discard the contents immediately.
- Do not use the reagents past the expiration date (see **Vial Label**).
- Do not use the reagents if a precipitate is present.
- Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
- The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but are not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
- The reagents contain < 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
- No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

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For information on disposal of the reagent and decontamination of a spillage site see **Material Safety Data Sheets**, available on request.

CONTROLS AND ADVICE

- It is recommended a positive control and a negative control be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.

- Since these reagents do not contain macromolecular potentiators, it is very unlikely that false positive reactions are caused with IgG coated cells.
- Blood specimens of weak A or B subgroups (e.g Ax) may give rise to false negative or weak reactions when tested using slides, microtitre plates or gel cards. It is advisable to re-test weak subgroups using tube technique.
- Individuals older than six months should have their ABO blood-grouping results confirmed by testing their serum or plasma against known group A, and B cells before their ABO blood group can be confirmed.
- Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
- In the Recommended Techniques one volume is approximately 50µl when using the vial dropper provided.
- The use of the reagents and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the reagents are in use.
- The user must determine the suitability of the reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED

- Applicator sticks.
- Automatic plate reader.
- Bio-Rad ID-Cards (NaCl, enzyme test and cold agglutinins).
- Bio-Rad ID-Centrifuge.
- Bio-Rad ID-CellStab or ID-Diluent 2.
- Glass microscope slides or white card tiles.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- Microplate centrifuge.
- Ortho BioVue System Cassettes (Neutral).
- Ortho BioVue System Centrifuge.
- Ortho 0.8% Red Cell Diluent.
- Plate shaker.
- PBS solution (pH 6.8–7.2) or Isotonic saline solution (pH 6.5–7.5).
- Positive and negative control red cells:
- Anti-A: group A (positive control) and group O (negative control).
- Anti-B: group B (positive control) and group O (negative control).
- Anti-A,B: group A and group B (positive controls) and group O (negative control).
- Test tube centrifuge.
- Validated "U" well microplates.
- Volumetric pipettes.

RECOMMENDED TECHNIQUES

A. Tube Technique

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in a labelled test tube: 1 volume of Lorne Anti-ABO reagent and 1 volume of red cell suspension.
- Mix thoroughly and incubate at room temperature for 1 minute.
- Centrifuge all tubes for 10 seconds at 1000 rcf or for a suitable alternative time and force.
- Gently resuspend red cell button and read macroscopically for agglutination
- Any tubes, which show a negative or questionable result, should be incubated for 15 minutes at room temperature.
- Following incubation, repeat steps 4 and 5.

B. Bio-Rad-ID Technique (NaCl, enzyme test and cold agglutinins cards)

- Prepare a 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
- Remove aluminium foil from as many microtubes as needed.
- Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne Anti-ABO reagent.
- Centrifuge ID-Card(s) in the Bio-Rad gel card centrifuge.
- Read macroscopically for agglutination.

C. Ortho BioVue Technique (Neutral cassettes)

- Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
- Remove aluminium foil from as many reaction chambers as needed.
- Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne Anti-ABO reagent.
- Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
- Read macroscopically for agglutination.

D. Microplate Technique, using "U" wells

- Prepare a 2-3% suspension of red cells in PBS or Isotonic saline.
- Place in the appropriate well: 1 volume Lorne Anti-ABO reagent and 1 volume red cell suspension.
- Mix thoroughly, preferably using a microplate shaker, taking care to avoid cross-well contamination.
- Incubate at room temperature for 15 minutes (time dependant on user).
- Centrifuge the microplate for 1 minute at 140 rcf or for a suitable alternative time and force.
- Resuspend the cell buttons using carefully controlled agitation on a microplate shaker

7. Read macroscopically or with a validated automatic reader.
8. Any weak reactions should be repeated by the tube technique.

E. Slide Technique

1. Prepare a 35-45% suspension of red cells in serum, plasma or PBS or Isotonic saline or use anti-coagulated whole blood (in its own plasma).
2. Place on a labelled glass slide or card tile: 1 volume of Lorne Anti-ABO reagent and 1 volume of red cell suspension.
3. Using a clean applicator stick, mix reagent and cells over an area of about 20 x 40 mm.
4. Slowly tilt the slide back and forth for 30 seconds, with occasional further mixing during the 1-minute period, maintaining slide at room temperature.
5. Read macroscopically after 1 minute over a diffuse light and do not mistake fibrin strands as agglutination.
6. Any weak reactions should be repeated by the tube technique.

INTERPRETATION OF TEST RESULTS

1. Positive: Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the appropriate ABO antigen on the red cells.
2. Negative: No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the appropriate ABO antigen on the red cells.
3. Discrepancies: If the results obtained with reverse group don't correlate with forward group, further investigation is required.

STABILITY OF THE REACTIONS

1. Read all tube and microplate tests immediately after centrifugation.
2. Slide tests should be interpreted after a maximum of one minute to ensure specificity and to avoid the possibility a negative result may be incorrectly interpreted as positive due to drying of the reagent.
3. Caution should be exercised in the interpretation of results of tests performed at temperatures other than those recommended.

LIMITATIONS

1. ABO antigens are not fully developed at birth and so weaker reactions may therefore occur with cord or neonatal specimens.
2. When using Monoclonal Anti-A,B, blood specimens of weak A or B subgroups (e.g Ax) may give rise to false negative or weak reactions when tested using slides, microtitre plates or gel cards. It is advisable to re-test weak subgroups using the tube technique.
3. Lorne monoclonal Anti-A and monoclonal Anti-B are not validated to detect Ax and A3 or Bx and B3 antigens resp and we therefore do not claim reactivity of the monoclonal Anti-A or Anti-B reagent against these weak A and B sub-groups.
4. Stored blood may give weaker reactions than fresh blood.
5. False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques
 - Cord samples contaminated with Wharton's jelly

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1. Prior to release, each lot of Lorne ABO monoclonal reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the 'Guidelines for the Blood Transfusion Services in the United Kingdom' and the 'Common Technical Specifications'.
2. Specificity of source monoclonal antibodies is demonstrated using a panel of antigen-negative cells.
3. The potency of the reagents has been tested against the following minimum potency reference standards obtained from National Institute of Biological Standards and Controls (NIBSC): Anti-A reference standard 03/188 And / Or Anti-B reference standard 03/164
4. Lorne Anti-B does not react with "Acquired-B" red cells.
5. Lorne Monoclonal ABO reagents do not detect crypt antigens such as T, Tn or Cad.
6. The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

1. The user is responsible for the performance of the reagents by any method other than those mentioned in the Recommended Techniques.
2. Any deviations from the Recommended Techniques should be validated prior to use⁸.

BIBLIOGRAPHY

1. Marion E. Reid and Christine Lomas-Francis, Blood Group Antigens and Antibodies, SBB Books, New York 2007; Page 181.
2. Issitt PD. Applied Blood Group Serology, 3rd Edition. Montgomery Scientific, Miami 1985; Chapter 6
3. Guidelines for the Blood Transfusion Service in the United Kingdom 6th Edition 2002. The Stationery Office.

4. AABB Technical Manual, 16th Edition, AABB 2008.
5. British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

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

TABEL SIMBOLURI

| | | | |
|---|---------------------|---|----------------------------|
|  | Batch Number |  | <i>In-vitro</i> Diagnostic |
|  | Catalogue Reference |  | Store At |
|  | Expiry Date |  | Manufacturer |
|  | Read Pack Insert | | |

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Lorne Laboratories Limited

Unit 1 Cutbush Park Industrial Estate

Danehill

Lower Earley, Reading,

Berkshire, RG6 4UT

United Kingdom

Tel: +44 (0) 118 921 2264

Fax: +44 (0) 118 986 4518

E-mail: info@lornelabs.com

DIRECTIONS FOR USE

Anti-D Clone 1 and Clone 2 Monoclonal:

For Tube, Bio-Rad-ID, Ortho BioVue, Microplate and Slide Techniques.

SUMMARY

The Rh blood group system was discovered in 1940. The D antigen is the most clinically significant non-ABO red blood cell antigen and has been implicated in causing Haemolytic Transfusion Reactions and Haemolytic Disease of the Newborn.

| Anti-D | Phenotype | Caucasians % ¹ | Afro-Americans % ² |
|--------|-----------|---------------------------|-------------------------------|
| + | Rh D +ve | 83 | 92 |
| 0 | Rh D -ve | 17 | 8 |

INTENDED PURPOSE

The Anti-D reagents are blood grouping reagents intended to be used to qualitatively determine the presence or absence of the Rh D antigen on the red cells of blood donors or patients requiring a blood transfusion when tested in accordance with the recommended techniques stated in this IFU.

PRINCIPLE

The reagents contain antibodies against the D antigen on human red cells and will cause direct agglutination (clumping) of human red cells that carry the D antigen. No agglutination (no clumping) generally indicates the absence of the D antigen on human red cells (see **Limitations**).

REAGENT

Lorne Monoclonal IgM Anti-D Clone 1 and Clone 2 blood grouping reagents are low protein reagents containing a human monoclonal IgM antibody diluted with sodium chloride (0.9 g%), bovine albumin (2.0 g%) and macromolecular potentiators (1.5 g%). When typing patient samples, each reagent will directly agglutinate Rh D positive cells, including majority of variants (but not D^u) and a high proportion of weak D (D^w) phenotypes when using the recommended techniques. The reagents do not contain or consist of CMR substances, or endocrine disrupting substances or that could result in sensitisation or an allergic reaction by the user. Each reagent is supplied at optimal dilution for use on patient samples with all recommended techniques stated below without need for further dilution or addition. For lot reference number and expiry date see Vial Label.

| Product | Cell Line / Clone |
|----------------|-------------------|
| Anti-D Clone 1 | RUM-1 |
| Anti-D Clone 2 | MS-201 |

WEAKENED EXPRESSION OF THE RhD ANTIGEN

The collective term D^w is widely used to describe red cells which have a weaker expression of the D antigen than normal. The term weak D denotes individuals with a reduced number of complete D antigen sites per red cell. The term partial D denotes individuals with missing D antigen epitopes. D^w cells is a partial D category which misses most D epitopes. Both Clone 1 and Clone 2 reagents will detect most examples of partial and weak D red cells by direct agglutination, but will not detect D^w cells.

STORAGE

Reagent vials should be stored at 2 - 8°C on receipt. Prolonged storage at temperatures outside this range may result in accelerated loss of reagent reactivity. This reagent has undergone transportation stability studies at 37°C and -25°C as described in document BS EN ISO 23640:2015.

SAMPLE COLLECTION AND PREPARATION

Blood samples can be collected into EDTA, citrate, CPDA anticoagulants or as a clotted sample. The samples should be tested as soon as possible following collection. If a delay in testing should occur, store the samples at 2-8°C. Samples displaying gross haemolysis or microbial contamination should not be used for testing. Blood samples showing evidence of lysis may give unreliable results. It is preferable (but not essential) to wash all blood samples with PBS or isotonic saline before being tested.

PRECAUTIONS

1. The reagents are intended for in vitro diagnostic use only.
2. If a reagent vial is cracked or leaking, discard the contents immediately.
3. Do not use the reagents past the expiration date (see Vial Label).
4. Do not use the reagents if a precipitate is present.
5. Protective clothing should be worn when handling the reagents, such as disposable gloves and a laboratory coat.
6. The reagents have been filtered through a 0.2 µm capsule to reduce the bio-burden, but are not supplied sterile. Once a vial has been opened the contents should remain viable up until the expiry date as long as there is no marked turbidity, which can indicate reagent deterioration or contamination.
7. The reagents contain < 0.1% sodium azide. Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive metal azides. On disposal flush away with large volumes of water.
8. Materials used to produce the products were tested at source and found to be negative for HIV 1+2 and HCV antibodies and HBsAg using approved microbiological tests.

9. No known tests can guarantee that products derived from human or animal sources are free from infectious agents. Care must be taken in the use and disposal of each vial and its contents.

DISPOSAL OF REAGENT AND DEALING WITH SPILLAGES

For information on disposal of the reagent and decontamination of a spillage site see Material Safety Data Sheets, available on request.

CONTROLS AND ADVICE

1. It is recommended a positive control (ideally R,r cells), and a negative control (ideally rr cells) be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show expected results.
2. When typing red cells from a patient who is diagnosed with a disease that causes the red cells to become coated with antibody or other proteins (such as HDN, AIHA), it is important to test the patient's red cells using Lorne's Monoclonal D Negative Control (catalogue # 650010). Tests must be considered invalid if red cells are agglutinated using Lorne's Monoclonal D Negative Control (catalogue # 650010).
3. Weak and partial D antigen variants are poorly detected by the gel card, microtitre plate and slide technique. It is recommended that weak and partial D variants are tested using the tube test technique.
4. Before use, let the reagent warm up to room temperature. As soon as the reagent has been used, put the reagent back in storage at 2-8°C.
5. In the Recommended Techniques one volume is approximately 50µl when using the vial dropper provided.
6. The use of the reagents and the interpretation of results must be carried out by properly trained and qualified personnel in accordance with the requirements of the country where the reagents are in use.
7. The user must determine suitability of the reagents for use in other techniques.

REAGENTS AND MATERIALS REQUIRED

- Applicator sticks.
- Automatic plate reader.
- Bio-Rad ID-Cards (NaCl, enzyme test and cold agglutinins).
- Bio-Rad ID-Centrifuge.
- Bio-Rad ID-CellStab or ID-Diluent 2.
- Glass microscope slides or white card tiles.
- Glass test tubes (10 x 75 mm or 12 x 75 mm).
- Microplate centrifuge.
- Ortho BioVue System Cassettes (Neutral).
- Ortho BioVue System Centrifuge.
- Ortho 0.8% Red Cell Diluent.
- Plate shaker.
- PBS solution (pH 6.8-7.2) or isotonic saline solution (pH 6.5-7.5).
- Positive (ideally R,r) and negative (rr) control red cells.
- Test tube centrifuge.
- Validated "U" well microplates.
- Volumetric pipettes.

RECOMMENDED TECHNIQUES

- A. Tube Technique**
1. Prepare a 2-3% suspension of red cells in PBS or isotonic saline.
 2. Place in a labelled test tube: 1 volume of Lorne Anti-D reagent and 1 volume of red cell suspension.
 3. Mix thoroughly and centrifuge all tubes for 20 seconds at 1000 rcf or for a suitable alternative time and force.
 4. Gently resuspend red cell button and read macroscopically for agglutination
 5. Any tubes, which show negative or questionable result (as can happen with weak D samples), should be incubated for 15 minutes at room temperature.
 6. Following incubation, repeat steps 3 and 4.
- B. Bio-Rad-ID Technique (NaCl, enzyme test and cold agglutinins cards)**
1. Prepare a 0.8% suspension of red cells in ID-CellStab or ID-Diluent 2.
 2. Remove aluminium foil from as many microtubes as needed.
 3. Place in appropriate microtube: 50µl of red cell suspension and 25µl of Lorne Anti-D reagent.
 4. Centrifuge ID-Card(s) in a Bio-Rad gel card centrifuge.
 5. Read macroscopically for agglutination.
- C. Ortho BioVue Technique (Neutral cards)**
1. Prepare a 0.8% suspension of red cells in 0.8% Ortho Red Cell Diluent.
 2. Remove aluminium foil from as many reaction chambers as needed.
 3. Place in appropriate reaction chamber: 50µl of red cell suspension and 40µl of Lorne Anti-D reagent.
 4. Centrifuge cassette(s) in an Ortho BioVue System Centrifuge.
 5. Read macroscopically for agglutination.
- D. Microplate Technique, using "U" wells**
1. Prepare a 2-3% suspension of red cells in PBS or isotonic saline.
 2. Place in the appropriate well: 1 volume Lorne Anti-D reagent and 1 volume

- red cell suspension.
- Mix thoroughly, preferably using a microplate shaker, taking care to avoid cross-well contamination.
 - Incubate at room temperature for 15 minutes (time dependant on user).
 - Centrifuge the microplate for 1 minute at 140 rcf or for a suitable alternative time and force.
 - Resuspend the cell buttons using carefully controlled agitation on a microplate shaker
 - Read macroscopically or with a validated automatic reader.
 - Any weak reactions should be repeated by the tube technique.

E. Slide Technique

- Prepare a 35-45% suspension of red cells in serum, plasma or PBS or Isotonic saline or use anti-coagulated whole blood (in its own plasma).
- Place on a labelled glass slide or card tile: 1 volume of Lorne Anti-D reagent and 1 volume of red cell suspension.
- Using a clean applicator stick, mix reagent and cells over an area of about 20 x 40 mm.
- Slowly tilt the slide back and forth for 30 seconds, with occasional further mixing during the 1-minute period, maintaining slide at room temperature.
- Read macroscopically after 1 minute over a diffuse light and do not mistake fibrin strands as agglutination.
- Any weak reactions should be repeated by the tube technique.

INTERPRETATION OF TEST RESULTS

- Positive: Agglutination of the red cells constitutes a positive test result and within accepted limitations of test procedure, indicates the presence of the D antigen on the red cells.
- Negative: No agglutination of the red cells constitutes a negative result and within the accepted limitations of the test procedure, indicates the absence of the D antigen on the red cells.
- Test results of cells that are agglutinated using the reagent negative control shall be excluded, as the agglutination is most probably caused by the effect of the macromolecular potentiators in the reagent on sensitised cells.

STABILITY OF THE REACTIONS

- Read all tube and microplate tests immediately after centrifugation.
- Slide tests should be interpreted after a maximum of one minute to ensure specificity and to avoid the possibility a negative result may be incorrectly interpreted as positive due to drying of the reagent.
- Caution should be exercised in the interpretation of results of tests performed at temperatures other than those recommended.

LIMITATIONS

- Lorne Anti-D is not suitable for use with enzyme treated cells, cells suspended in LISS or for use in indirect antiglobulin (IAT) techniques.
- Stored blood may give weaker reactions than fresh blood.
- False positive agglutination may be seen due to the presence of macromolecular potentiators in the reagent when testing IgG sensitised cells, e.g. ALHA, HDN.
- False positive or false negative results may also occur due to:
 - Contamination of test materials
 - Improper storage, cell concentration, incubation time or temperature
 - Improper or excessive centrifugation
 - Deviation from the recommended techniques

SPECIFIC PERFORMANCE CHARACTERISTICS

- Prior to release, each lot of Lorne Anti-D monoclonal reagent was tested using the recommended test methods listed in this IFU. The tests complied with the test requirements as stated in the current version/issue of the 'Guidelines for the Blood Transfusion Services in the United Kingdom' and the 'Common Technical Specifications'.
- Anti-D grouping reagents for D grouping of patients should not react with D⁰ cells using the method(s) recommended for use.
- Specificity of source monoclonal antibodies is demonstrated using a panel of antigen-negative cells.
- The potency of the reagents has been tested against the following minimum potency reference standard obtained from National Institute of Biological Standards and Controls (NIBSC):
 - Anti-D reference 99/836.
- The Quality Control of the reagents was performed using red cells with phenotypes that were verified by a UK blood transfusion centre and had been washed with PBS or Isotonic saline prior to use.

DISCLAIMER

- The user is responsible for the performance of the reagents by any method other than those mentioned in the Recommended Techniques.
- Any deviations from the Recommended Techniques should be validated prior to use⁹.

BIBLIOGRAPHY

- Issitt PD. Applied Blood Group Serology, 3rd Edition, Montgomery Scientific, Miami, 1985, Chapter 10.
- AABB Technical Manual, 16th Edition, AABB 2008.

- Marion E. Reid and Christine Lomas-Francis, Blood Group Antigens and Antibodies, SBB Books, New York 2007; Page 192.
- Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine 1995. 5, 171-184
- Guidelines for the Blood Transfusion Service in the United Kingdom, 6th Edition 2002. The Stationery Office.
- British Committee for Standards in Haematology, Blood Transfusion Task Force. Recommendations for evaluation, validation and implementation of new techniques for blood grouping, antibody screening and cross matching. Transfusion Medicine, 1995, 5, 145-150.

AVAILABLE REAGENT SIZES

| | Vial Size | Catalogue Number | Test per vial |
|------------------------------|-----------|------------------|---------------|
| Anti-D Clone 1 Monoclonal | 10 ml | 730010 | 200 |
| | 1000 ml | 730000* | 20,000 |
| Anti-D Clone 2 Monoclonal | 5000 ml | 730000X5* | 100,000 |
| | 10 ml | 710010 | 200 |
| | 1000 ml | 710000* | 20,000 |
| | 5000 ml | 710000X5* | 100,000 |

*This size is For Further Manufacturing Use (FFMU) only and is therefore not CE marked.

Advena Ltd, Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013, Malta

Lorne Laboratories Limited

Unit 1 Cutbush Park Industrial Estate, Danehill, Lower Earley, Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264 Fax: +44 (0) 118 986 4518 Email: info@lornelabs.com www.lornelabs.com

Blood Reagents and Diagnostic Kits

Quality blood reagents and diagnostic kits delivered worldwide


Lorne Laboratories provides high quality blood-grouping reagents and diagnostic kits. We do so at affordable prices and back this up with excellent customer service.

For over forty four years, Lorne has seen that philosophy transform the company from a small UK operation to one that spans the globe. The Lorne name is known and respected in over 110 countries around the world and it has come to mean 'quality' to the many blood transfusion professionals who use the products, both in the UK and across the five continents in which we operate.

Lorne has achieved ISO 13485 and MDSAP accreditation and our reagents and diagnostic kits are all CE marked. A significant number of our blood-grouping reagents are also registered with Health Canada.

Contents

BLOOD GROUPING REAGENTS Page

| | |
|--|---|
| ABO system – Monoclonal Reagents..... | 4 |
| Rhesus System – Monoclonal Reagents..... | 4 |
| Grouping Reagents for M, N and S Blood Group Systems | 5 |
| Grouping Reagents for Kell Blood Group System..... | 5 |
| Grouping Reagents for Rare Blood Groups..... | 6 |
| Control Reagents..... | 7 |
| Anti-Human IgG Reagent..... | 7 |
| Anti-Human Globulin Reagent | 7 |
| Monoclonal Anti-C3d Reagent..... | 7 |

ENZYMES AND POTENTIATORS Page

| | |
|-------------------------------|---|
| Buffered Saline Tablets | 8 |
| Papenzyme-plus..... | 8 |
| Bromelite | 8 |
| LISS Concentrate..... | 8 |
| LISS Ready for use..... | 8 |
| LISS-ADD | 8 |
| PEG-ADD..... | 8 |
| Serological Albumin | 8 |
| Accessories | 9 |

RED CELLS Page

| | |
|--|----|
| Red Cells – Reverse Grouping Cells | 10 |
| Red Cells – Antibody Screening Cells..... | 10 |
| Red Cells – Identicells..... | 10 |
| Red Cells – Coombs Control Cells..... | 10 |
| Red Cells Preservative – Preservacell..... | 11 |
| Red Cells Preservative – ABO Preservacell..... | 11 |
| Alsevers Solution..... | 11 |
| Red Cell Elute | 11 |

DIAGNOSTIC KITS Page

| | |
|------------------------|----|
| Syphilis Kits | 12 |
| Latex Kits..... | 13 |
| Febrile Antigens | 14 |
| Rose Waaler..... | 14 |

NON CE MARKED BLOOD GROUPING REAGENTS Page

| | |
|--|----|
| ABO System - Monoclonal Reagents..... | 15 |
| Rhesus System - Monoclonal Reagents..... | 15 |

ABO SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM ABO blood grouping reagents contain mouse monoclonal antibodies diluted in a phosphate buffer containing sodium chloride, EDTA and bovine albumin. Each reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|---------------------|--------------|------|--------------------|
| Anti-A Monoclonal | 600010 | 10ml | 36 Months |
| Anti-B Monoclonal | 610010 | 10ml | 36 Months |
| Anti-A,B Monoclonal | 620010 | 10ml | 36 Months |

Lorne Anti-A₁ Lectin blood grouping reagent is prepared from an extract of *Dolichos biflorus* seeds, diluted with a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube technique.

Lorne Anti-H Lectin blood grouping reagent is prepared from an extract of *Ulex europaeus* seeds, diluted with a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube technique.

| | | | |
|----------------------------|--------|-----|-----------|
| Anti-A ₁ Lectin | 116005 | 5ml | 24 Months |
| Anti-H Lectin | 115002 | 2ml | 24 Months |

RHESUS SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM Anti-D Clone 1 and Clone 2 blood grouping reagents are low protein reagents containing a human monoclonal IgM antibody diluted with sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, each reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^v) and a high proportion of weak D (D^w) phenotypes when using the slide, tube, gel card and microplate techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|---------------------------|--------------|------|--------------------|
| Anti-D Clone 1 Monoclonal | 730010 | 10ml | 30 Months |
| Anti-D Clone 2 Monoclonal | 710010 | 10ml | 30 Months |

Lorne Monoclonal Anti-D Duoclone blood grouping reagent is a low protein, blended reagent containing human monoclonal IgM and IgG Anti-D, diluted in a phosphate buffer containing sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, this reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^v) and a high proportion of weak D (D^w) phenotypes when using the slide, tube, gel card and microplate techniques. It will agglutinate D^v cells in the IAT phase of testing.

| | | | |
|----------------------------|--------|------|-----------|
| Anti-D Duoclone Monoclonal | 740010 | 10ml | 30 Months |
|----------------------------|--------|------|-----------|

Lorne Monoclonal IgM Anti-Rh blood grouping reagents are low protein reagents containing human monoclonal antibodies diluted with sodium chloride, bovine albumin and macromolecular potentiators. Each reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

| | | | |
|--------------------------------|--------|------|-----------|
| Anti-C Monoclonal | 690005 | 5ml | 24 Months |
| Anti-C ^w Monoclonal | 750002 | 2ml | 24 Months |
| Anti-E Monoclonal | 691005 | 5ml | 24 Months |
| Anti-c Monoclonal | 692005 | 5ml | 24 Months |
| Anti-e Monoclonal | 693005 | 5ml | 24 Months |
| Anti-C+D+E Monoclonal | 700010 | 10ml | 24 Months |

All pack inserts are available on www.lornelabs.com

GROUPING REAGENTS FOR M, N AND S BLOOD GROUP SYSTEMS

Lorne Human Anti-M blood grouping reagent is prepared from human serum diluted in a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube and gel card techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|-------------------|--------------|------|--------------------|
| Anti-M Polyclonal | 311002 | 2ml | 24 Months |

Lorne Anti-N Lectin blood grouping reagent is prepared from an extract of *Vicia*, diluted with a sodium chloride solution containing bovine albumin. The reagent is supplied at optimal dilution for use by tube and gel card techniques.

| | | | |
|---------------|--------|-----|-----------|
| Anti-N Lectin | 312002 | 2ml | 24 Months |
|---------------|--------|-----|-----------|

Lorne Monoclonal Anti-S and Anti-s blood grouping reagents contain human monoclonal antibodies diluted in a phosphate buffer containing sodium chloride and bovine albumin. They are supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

| | | | |
|-------------------|--------|-----|-----------|
| Anti-S Monoclonal | 770002 | 2ml | 24 Months |
| Anti-s Monoclonal | 771002 | 2ml | 24 Months |

GROUPING REAGENTS FOR KELL BLOOD GROUP SYSTEM

Lorne Monoclonal Anti-K blood grouping reagent is a low protein reagent containing the monoclonal IgM antibody, diluted in a phosphate buffer containing sodium chloride, bovine albumin and macromolecular potentiators. The reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|-------------------|--------------|------|--------------------|
| Anti-K Monoclonal | 760010 | 10ml | 24 Months |

Lorne Monoclonal Anti-k (Cellano) blood grouping reagent contains monoclonal IgG antibody diluted in sodium chloride containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube and gel card techniques without the need for further dilution or addition.

| | | | |
|-----------------------------|--------|-----|-----------|
| Anti-k (Cellano) Monoclonal | 325002 | 2ml | 24 Months |
|-----------------------------|--------|-----|-----------|

Lorne Human Anti-Kp^a and Anti-Kp^b blood grouping reagents are prepared from human serum diluted in sodium chloride containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube and gel card techniques without the need for further dilution or addition.

| | | | |
|---------------------------------|--------|-----|-----------|
| Anti-Kp ^a Polyclonal | 321002 | 2ml | 24 Months |
| Anti-Kp ^b Polyclonal | 322002 | 2ml | 24 Months |

All products available in bulk quantities

GROUPING REAGENTS FOR RARE BLOOD GROUPS

Lorne Monoclonal Fy^a blood grouping reagent contains human monoclonal antibodies diluted in a phosphate buffer containing sodium chloride and bovine albumin.

Lorne Human Anti-Fy^b blood grouping reagent is prepared from human serum diluted in a sodium chloride solution containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|---------------------------------|--------------|------|--------------------|
| Anti-Fy ^a Monoclonal | 774002 | 2ml | 24 Months |
| Anti-Fy ^b Polyclonal | 317002 | 2ml | 24 Months |

Lorne Anti-Jk^a and Anti-Jk^b Monoclonal reagents contain human monoclonal IgM antibodies, sodium chloride, phosphate buffer and bovine albumin. When typing red cells in accordance with the instructions for use, the reagents will directly agglutinate either Jk^a positive cells or Jk^b positive cells. Both reagents are supplied at optimal dilution for use by the test tube technique.

| | | | |
|---------------------------------|--------|-----|-----------|
| Anti-Jk ^a Monoclonal | 775002 | 2ml | 24 Months |
| Anti-Jk ^b Monoclonal | 776002 | 2ml | 24 Months |

Lorne Human Anti-Lu^a and Anti-Lu^b blood grouping reagents are prepared from human serum diluted in a sodium chloride solution containing macromolecular potentiators and bovine albumin. Each reagent is supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

| | | | |
|---------------------------------|--------|-----|-----------|
| Anti-Lu ^a Polyclonal | 330002 | 2ml | 24 Months |
| Anti-Lu ^b Polyclonal | 331002 | 2ml | 24 Months |

Lorne Monoclonal IgM Anti-P₁ blood grouping reagent contains mouse monoclonal IgM antibodies, diluted in a solution containing sodium chloride and bovine albumin. The reagent is supplied at optimal dilution for use by tube and gel card techniques.

| | | | |
|--------------------------------|--------|-----|-----------|
| Anti-P ₁ Monoclonal | 315002 | 2ml | 24 Months |
|--------------------------------|--------|-----|-----------|

Lorne Monoclonal Anti-Le^a contains human monoclonal IgM and Anti-Le^b contains mouse monoclonal IgM antibodies, diluted in a phosphate buffer containing sodium chloride, EDTA, bovine albumin and macromolecular potentiators. Each reagent is supplied at optimal dilution for use by the tube technique. The Anti-Le^a reagent can also be used for gel card techniques.

| | | | |
|---------------------------------|--------|-----|-----------|
| Anti-Le ^a Monoclonal | 632002 | 2ml | 24 Months |
| Anti-Le ^b Monoclonal | 631002 | 2ml | 24 Months |

Lorne Human Anti-Di^a blood grouping reagent is prepared from human serum diluted in a sodium chloride solution containing macromolecular potentiators and bovine albumin. The reagent is supplied at optimal dilution for use by indirect Coombs tube test and gel card techniques.

| | | | |
|---------------------------------|--------|-----|-----------|
| Anti-Di ^a Polyclonal | 328002 | 2ml | 24 Months |
|---------------------------------|--------|-----|-----------|

All pack inserts are available on www.lornelabs.com

CONTROL REAGENTS

Lorne Monoclonal Rh Control is formulated with the same levels of potentiators and protein as Lorne Monoclonal Rh Grouping Reagents with the blood group antibodies omitted. The reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|-----------------------|--------------|------|--------------------|
| Monoclonal Rh Control | 640010 | 10ml | 24 Months |

Lorne Monoclonal D Negative Control is for the control of Monoclonal Anti-D reagents and is formulated with the same concentrations of phosphate buffer, sodium chloride, bovine albumin and macromolecular potentiators as Lorne Monoclonal Anti-D reagents with the antibodies omitted. The reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

| | | | |
|-------------------------------|--------|------|-----------|
| Monoclonal D Negative Control | 650010 | 10ml | 30 Months |
|-------------------------------|--------|------|-----------|

Lorne Inert AB Serum is prepared from pooled human serum. No potentiators or any other chemicals have been added to the reagent. This reagent is supplied at optimal dilution for use by the techniques recommended for the reagent to be controlled.

| | | | |
|----------------|--------|------|-----------|
| Inert AB Serum | 110010 | 10ml | 24 Months |
|----------------|--------|------|-----------|

Lorne Precise Weak Anti-D Control Reagent is prepared from pools of human serum containing low activity Anti-D. The pool is diluted in inert serum to give a final concentration of 0.09 IU/ml Anti-D. ABO antibodies are not absorbed. This polyclonal reagent is supplied at the optimal dilution, for use by the techniques recommended for the reagent to be controlled.

| | | | |
|---------------------|--------|-----|-----------|
| Precise Weak Anti-D | 209005 | 5ml | 24 Months |
|---------------------|--------|-----|-----------|

ANTI-HUMAN IgG REAGENT

Lorne Monospecific Anti-Human IgG Clear and Anti-Human IgG Green reagents contain anti-human IgG derived from rabbits. All non-specific activity is removed by absorption. These monospecific reagents are supplied at optimal dilution, for use by spin tube technique.

| Item | Product Code | Size | Maximum Shelf Life |
|------------------------|--------------|------|--------------------|
| Anti-Human IgG (Clear) | 401010 | 10ml | 24 Months |
| Anti-Human IgG (Green) | 402010 | 10ml | 24 Months |

ANTI-HUMAN GLOBULIN REAGENT

Lorne Polyspecific Anti-Human Globulin Elite Clear and Anti-Human Globulin Elite Green reagents contain anti-human IgG derived from rabbits with non-specific activity removed by absorption and mouse monoclonal IgM Anti-C3d, Clone BRIC-8. The antibodies are diluted in a buffered solution containing bovine albumin. These polyspecific reagents are supplied at optimal dilution, for use by spin tube technique.

| Item | Product Code | Size | Maximum Shelf Life |
|---------------------|--------------|------|--------------------|
| A.H.G Elite (Clear) | 415010 | 10ml | 24 Months |
| A.H.G Elite (Green) | 435010 | 10ml | 24 Months |

MONOCLONAL ANTI-C3d REAGENT

Lorne Monoclonal IgM Anti-C3d blood grouping reagent contains mouse monoclonal Anti-C3d, Clone BRIC-8. The reagent is supplied at optimal dilution, for use by direct tube technique.

| Item | Product Code | Size | Maximum Shelf Life |
|---------------------|--------------|------|--------------------|
| Anti-C3d Monoclonal | 427002 | 2ml | 24 Months |

All products available in bulk quantities

ENZYMES AND POTENTIATORS

Lorne Phosphate Buffered Saline Tablets provide a safe, standardised saline solution for transfusion serology. Each tablet makes 1 litre of solution.

| Item | Product Code | Size | Maximum Shelf Life |
|-------------------------|--------------|------|--------------------|
| Buffered Saline Tablets | 490025 | 25 | 24 Months |
| | 490250 | 250 | 24 Months |

Lorne Papanzyme-plus reagent is a ready to use liquid preparation of stabilised papain. The reagent is standardised by serological methods for use in blood group antibody investigations. The reagent is supplied at optimal dilution for use by tube technique.

| | | | |
|----------------|--------|------|-----------|
| Papanzyme-plus | 441010 | 10ml | 12 Months |
|----------------|--------|------|-----------|

Lorne Bromelite reagent is a ready to use liquid preparation of stabilised bromelin. The reagent is standardised by serological methods for use in blood group antibody investigations. The reagent is supplied at optimal dilution for use by tube technique.

| | | | |
|-----------|--------|------|-----------|
| Bromelite | 443010 | 10ml | 12 Months |
|-----------|--------|------|-----------|

Lorne LISS Concentrate is a solution of glycine, phosphate buffer and 0.3 M sodium chloride. The solution is supplied at a stronger concentration than needed for serological use. It must be diluted 10 times in de-ionised water before being used by all recommended techniques in the pack insert.

| | | | |
|------------------|--------|--------|-----------|
| LISS Concentrate | 460500 | 500ml | 24 Months |
| | 460025 | 2500ml | 24 Months |

Lorne LISS ready for use is a low ionic strength solution containing glycine, sodium chloride and phosphate buffer. The reagent is supplied at the optimal dilution ready for use by all recommended techniques in the pack insert.

| | | | |
|--------------------|--------|------------|-----------|
| LISS Ready for use | 470020 | 20 x 250ml | 12 Months |
| | 470250 | 4 x 250ml | 12 Months |
| | 470025 | 2500ml | 12 Months |

Lorne LISS-ADD is a low ionic strength solution containing glycine, sodium chloride, phosphate buffer and bovine albumin. The reagent is supplied at the optimal dilution, for use by all the recommended techniques in the pack insert.

| | | | |
|----------|--------|------|-----------|
| LISS-ADD | 480010 | 10ml | 24 Months |
|----------|--------|------|-----------|

Lorne PEG-ADD is a low ionic strength solution contains glycine, a phosphate buffer and polyethylene glycol. The reagent is supplied at optimal dilution for use by all the recommended techniques stated in the pack insert without the need for further dilution or addition.

| | | | |
|---------|--------|------|-----------|
| PEG-ADD | 485010 | 10ml | 24 Months |
|---------|--------|------|-----------|

Lorne 22% and 30% Serological Albumin is prepared from a mixture of bovine serum albumin and buffered saline. No artificial avidity enhancers or high molecular weight agglutination potentiators are added to any BSA preparation. None of the BSA reagents contain sodium caprylate. Each BSA reagent is supplied at optimal dilution for use by all recommended techniques in the pack insert.

| | | | |
|-------------------------|--------|------|-----------|
| Serological Albumin 22% | 451010 | 10ml | 24 Months |
| Serological Albumin 30% | 452010 | 10ml | 24 Months |

All products available in bulk quantities

ACCESSORIES

| Item | Product Code | Size |
|---|--------------|-----------|
|  <p>The Essex Blockfile is a sturdy box for the storage of up to 200 cassette type histology blocks in each box. Delivered flat in packs of 25, they are easily assembled when needed and provide convenient and economical long-term storage for tissue samples. Dimensions – (L) 390mm x (D) 200mm x (H) 45mm</p> | 882025 | 25 Boxes |
|  <p>The Brighton Slidestak is a sturdy box for the storage of over 1000 microscope slides upright in each box. Delivered flat in packs of 25, they are easily assembled and provide convenient and economical long-term storage. Dimensions – (L) 310mm x (D) 170mm x (H) 80mm</p> | 881025 | 25 Boxes |
|  <p>The Cardtiles are suitable for a wide range of haemagglutination tests. Each card has twenty 30mm square white test sites.</p> | 880100 | 100 Cards |
|  <p>The Mini-Cardtiles are specially laminated cards suitable for VDRL and similar agglutination tests requiring a white background. Each card has ten 20mm diameter white test sites that are slightly compressed to form shallow wells.</p> | 880120 | 25 Cards |
|  <p>The Latex Cardtiles are suitable for all latex agglutination tests. Each card has six 30mm diameter black reaction sites.</p> | 880130 | 25 Cards |
|  <p>The mini pipettes are 11.5cm in length with one sealed end. The diameter of the tube is 4mm.</p> | 044000 | 500 |
|  <p>We use these products in our own production. The 5ml vials are glass and the 10ml vials are plastic. Each vial is supplied with its own dropper. The 10ml vials are supplied in a sterile format in multiples of 1224 while all the other products are supplied individually and require autoclaving.</p> | LAB00002 | 5ml |
| Vials and Droppers | LAB00036 | 10ml |
|  <p>We use these products in our own production. They are not branded. We supply flatpacks that hold 5 or 10 vials</p> | LAB00027 | 5 Vial |
| Flatpacks | LAB00007 | 10 Vial |

All pack inserts are available on www.lornelabs.com

RED CELLS - REVERSE GROUPING CELLS

These Grouping Cells are made from red cells that have been washed to remove all traces of blood group antibodies and then resuspended in a preservative solution of buffered saline containing adenine, inosine, chloramphenicol and neomycin sulphate. The preservative solution does not interfere with complement-mediated haemolysis. Each group A₁, A₂ and B vial contains a 2.5-3.0% red cell suspension derived from the blood of a single donor, whereas the group O vial contains a 2.5-3.0% suspension of pooled group O red cells from two donors in equal proportions. Antigens for which the O cells have been typed are noted on the antigen profile accompanying each lot. NB: Some antigens are present on only 50% of the cells in each lot.

| Item | Product Code | Size | Minimum Shelf Life |
|--|--------------|----------|--------------------|
| Reverse Grouping A ₁ + B Cells | 910020 | 2 x 10ml | 30 Days |
| Reverse Grouping A ₁ + A ₂ + B Cells | 910030 | 3 x 10ml | 30 Days |
| Reverse Grouping A ₁ + A ₂ + B + O Cells | 910040 | 4 X 10ml | 30 Days |
| Reverse Grouping A ₂ Cells | 920002 | 1 x 10ml | 30 Days |

RED CELLS - ANTIBODY SCREENING CELLS

Maxi-Screen 3 red cells are for antibody screening. Each reagent vial contains a 2.5-3.0% suspension of red cells derived from the blood of a single group O donor. The donor red cells have been washed to remove blood group antibodies and then resuspended in a preservative solution containing adenine and inosine to help preserve carbohydrate metabolism and chloramphenicol and neomycin sulphate as preservatives.

| Item | Product Code | Size | Minimum Shelf Life |
|---------------|--------------|----------|--------------------|
| Maxi-Screen 3 | 950030 | 3 x 10ml | 30 Days |

RED CELLS - IDENTICELLS

Identicells are made up of 10 vials which each contain a 2.5-3.0% suspension of red cells derived from the blood of a single group O donor. The donor cells have been washed to remove blood group antibodies and then resuspended in a preservative solution containing adenine, inosine, chloramphenicol and neomycin sulphate as preservatives.

| Item | Product Code | Size | Minimum Shelf Life |
|-------------|--------------|----------|--------------------|
| Identicells | 960050 | 10 x 5ml | 30 Days |

RED CELLS - COOMBS CONTROL CELLS

Coombs Control Cells are made up of a 3.8-4.2% suspension of single donor group O red cells washed to remove all blood group antibodies and then resuspended in a preservative solution. The preservative solution contains neomycin sulphate and chloramphenicol as preservatives. The cells are then sensitised with IgG.

| Item | Product Code | Size | Minimum Shelf Life |
|----------------------|--------------|------|--------------------|
| Coombs Control Cells | 970010 | 10ml | 30 Days |

Deliveries take place every 28 days.

As well as the standard Red Cell products above, we can also supply other products for special standing orders only. These include 0.8% panels and screening cells for use with column agglutination systems and Papainised red cell panels. Contact us with full details of your requirements for a quotation.

All the above Red Cell products are CE marked and comply with the Red Book.

All pack inserts are available on www.lornelabs.com

They are not Lorne labelled products.

RED CELLS PRESERVATIVE – PRESERVACELL

Lorne Preservacell is a phosphate buffered solution containing glucose, calcium chloride and purine bases, with chloramphenicol, gentamycin sulphate and neomycin sulphate as antibiotics. The reagent is supplied at optimal dilution.

| Item | Product Code | Size | Maximum Shelf Life |
|--------------|--------------|-------|--------------------|
| Preservacell | 980500 | 500ml | 12 Months |

RED CELLS PRESERVATIVE – ABO PRESERVACELL

Lorne ABO Preservacell is a phosphate buffered solution containing glucose, calcium chloride and purine bases, with chloramphenicol, gentamycin sulphate and neomycin sulphate as antibiotics. EDTA is added to stop complement binding so that potent ABO haemolysins act as simple agglutinins. The reagent is supplied at optimal dilution.

| Item | Product Code | Size | Maximum Shelf Life |
|------------------|--------------|-------|--------------------|
| ABO Preservacell | 981500 | 500ml | 12 Months |

ALSEVERS SOLUTION

Alsevers Solution is an isotonic, balanced salt solution that is routinely used as an anticoagulant/blood preservative, which permits the storage of whole blood at refrigerator temperatures for approximately 10 weeks. The solution contains antibiotics and other chemicals that ensure maximum viability of stored red cells.

| Item | Product Code | Size | Maximum Shelf Life |
|-------------------|--------------|--------|--------------------|
| Alsevers Solution | 983000 | 1000ml | 12 Months |

RED CELL ELUTE

Lorne Red Cell Elute is an acid elution kit. The kit consists of Concentrated Wash Solution, which is used to minimise antibody dissociation during washing, Acid Eluting Solution, which is a low pH glycine buffer containing a colouring agent and a Base Buffering Solution, Tris solution containing bovine albumin. The Concentrated Wash Solution requires dilution and all the other solutions are supplied at optimal dilution for use by all recommended techniques in the pack insert.

| Item | Product Code | Size | Maximum Shelf Life |
|----------------|--------------|------|--------------------|
| Red Cell Elute | 930110 | Kit | 30 Months |

All products available in bulk quantities

SYPHILIS KITS

Lorne TPHA Kit detects antibodies to *T. pallidum*. Test Cells are preserved avian erythrocytes coated with antigenic components of pathogenic *T. pallidum* (Nichol's strain). Any non-specific reactions are detected using the Control Cells; avian erythrocytes not coated with *T. pallidum* antigens. Non-specific reactions can also be absorbed out using Control Cells. Antibodies to non-pathogenic treponemes are absorbed by an extract of Reiter's treponemes in the cell suspension. Reagents are supplied at optimal dilution for use by the recommended techniques in the pack insert.

| Item | Product Code | Size | Maximum Shelf Life |
|---------------------------|--------------|-----------|--------------------|
| TPHA Microtitre Plate Kit | 043100A | 100 Tests | 18 Months |

The VDRL test is a non-treponemal slide agglutination test for the qualitative and semi-quantitative detection of plasma reagins. The antigen suspension, a lipid complex, is agglutinated when mixed with samples containing reagins of patients affected by syphilis. The reagent is supplied ready to use.

| | | | |
|-----------------------------|---------|-----------|-----------|
| VDRL Stabilised Reagent Kit | 046511A | 250 Tests | 30 Months |
|-----------------------------|---------|-----------|-----------|

RPR Carbon Antigen contains micro particulate carbon, which aids in the microscopic reading of results. Lorne provides reagents, controls and kits. All the reagents are supplied at optimum dilution for use by all recommended techniques without the need for further dilution or addition.

| | | | |
|-----------------------------|---------|-----------|-----------|
| RPR Carbon Antigen | 045005A | 100 Tests | 30 Months |
| RPR Carbon Positive Control | 047001A | 1ml | 30 Months |
| RPR Carbon Kit | 044150A | 150 Tests | 30 Months |
| | 044500A | 500 Tests | 30 Months |

RPR Carbon Kit

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Strep Test Kit

All pack inserts are available on www.lornelabs.com

LATEX KITS

Lorne ASO Latex Kit is a serologic test for the detection of ASO antibodies. All the reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

| Item | Product Code | Size | Maximum Shelf Life |
|---------------|--------------|-----------|--------------------|
| ASO Latex Kit | 031100A | 100 Tests | 30 Months |

Lorne RF Latex Kit is for the detection of Rheumatoid Factor. The latex reagent is a suspension of latex particles coated with human gamma globulins, which agglutinate in the presence of Rheumatoid Factor. All latex reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

| | | | |
|--------------|---------|-----------|-----------|
| RF Latex Kit | 830100A | 100 Tests | 30 Months |
|--------------|---------|-----------|-----------|

Lorne CRP Latex Test Kit is for the detection of CRP. The test reagent consists of latex particles coated with either goat or rabbit Anti-CRP (IgG). All the latex reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

| | | | |
|---------------|---------|-----------|-----------|
| CRP Latex Kit | 850100A | 100 Tests | 30 Months |
|---------------|---------|-----------|-----------|

Lorne IM Latex Test Kit is for detection of the heterophile antibody associated with Infectious Mononucleosis. The test reagent consists of latex particles coated with partially purified glycoprotein from bovine red cells. All reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

| | | | |
|--------------|---------|----------|-----------|
| IM Latex Kit | 041050A | 50 Tests | 30 Months |
|--------------|---------|----------|-----------|

Lorne LE Latex Kit is for the detection of nuclear proteins in Systemic Lupus Erythematosus (SLE or LE). The test reagent consists of DNP coated latex particles. All the reagents are supplied at optimal dilution for use by recommended techniques in the pack insert.

| | | | |
|-------------------|--------|----------|-----------|
| LE Latex Test Kit | 840050 | 50 Tests | 24 Months |
|-------------------|--------|----------|-----------|

Lorne Strep Kit is a Latex agglutination grouping kit for the identification of Streptococci of Lancefield groups A, B, C, D, F and G by agglutination of specific antibody coated latex particles in the presence of enzymatically-extracted antigen. All the reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

| | | | |
|----------------|--------|--------------|-----------|
| Strep Test Kit | 860050 | 6 x 50 Tests | 18 Months |
|----------------|--------|--------------|-----------|

Lorne Staph Kit is a Latex agglutination kit for the identification of *Staph. aureus*. Includes latex reagent, control reagent and agglutination slides. All the reagents are supplied at optimal dilution for use by all recommended techniques in the pack insert.

| | | | |
|----------------|--------|-----------|-----------|
| Staph Test Kit | 870050 | 50 Tests | 18 Months |
| | 870100 | 100 Tests | 18 Months |

All products available in bulk quantities

FEBRILE ANTIGENS

Lorne Stained Febrile Antigens are for the detection of certain Salmonellae, Rickettsiae and Brucellae pathogens. The antigens are suspensions of killed bacteria, stained to enhance the reading of agglutination tests. The blue stained antigens are specific to the somatic "O" antigens and the red stained antigens are specific to the flagellar "H" antigens. Suspensions of Proteus OX2, OX19 and OXK are used to detect rickettsial antibodies.

Rose Bengal is a slide agglutination test for the qualitative and semi-quantitative detection of anti-Brucella antibodies in human and animal serum. The stained bacterial suspension agglutinates when mixed with samples containing specific IgG or IgM antibodies present in the patient sample.

| Item | Product Code | Size | Maximum Shelf Life |
|--------------------------------|--------------|-----------------------|--------------------|
| Salmonella Typhi H | 502005A | 100 Tests | 30 Months |
| Salmonella Paratyphi AH | 504005A | 100 Tests | 30 Months |
| Salmonella Paratyphi BH | 506005A | 100 Tests | 30 Months |
| Salmonella Paratyphi CH | 508005A | 100 Tests | 30 Months |
| Salmonella Typhi O | 510005A | 100 Tests | 30 Months |
| Salmonella Paratyphi AO | 512005A | 100 Tests | 30 Months |
| Salmonella Paratyphi BO | 514005A | 100 Tests | 30 Months |
| Salmonella Paratyphi CO | 516005A | 100 Tests | 30 Months |
| Brucella Abortus | 518005A | 100 Tests | 30 Months |
| Brucella Melitensis | 520005A | 100 Tests | 30 Months |
| Proteus OX2 | 522005A | 100 Tests | 30 Months |
| Proteus OX19 | 524005A | 100 Tests | 30 Months |
| Proteus OXK | 526005A | 100 Tests | 30 Months |
| Febrile Antigen Kit + Controls | 532042A | 8x100 Tests and 2x1ml | 30 Months |
| Febrile Positive Control | 536001A | 1ml | 30 Months |
| Febrile Negative Control | 537001A | 1ml | 30 Months |
| Rose Bengal | 155050A | 50 Tests | 30 Months |

ROSE WAALER

Rose Waaler is a technique of passive haemagglutination for qualitative and semi-quantitative detection of human serum Rheumatoid Factor. Sheep red blood cells are coated with a concentration of antiserum to sheep red blood cells that is too low to cause agglutination. The addition of serum from a patient with Rheumatoid Factor will cause agglutination.

| Item | Product Code | Size | Maximum Shelf Life |
|-------------|--------------|----------|--------------------|
| Rose Waaler | 156050A | 50 Tests | 30 Months |

Febrile Antigen Kit + Controls

All pack inserts are available on www.lornelabs.com

ABO SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM ABO Standard Grade blood grouping reagents contain mouse monoclonal antibodies diluted in a phosphate buffer containing sodium chloride, EDTA and bovine albumin. Each reagent is supplied at optimal dilution for use by slide, tube, gel card and microplate techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|--------------------------------------|--------------|------|--------------------|
| Anti-A Monoclonal - Standard Grade | 600010E | 10ml | 36 Months |
| Anti-B Monoclonal - Standard Grade | 610010E | 10ml | 36 Months |
| Anti-A,B Monoclonal - Standard Grade | 620010E | 10ml | 36 Months |

RHESUS SYSTEM – MONOCLONAL REAGENTS

Lorne Monoclonal IgM Anti-D Clone 1 Standard Grade blood grouping reagent is a low protein reagent containing a human monoclonal IgM antibody diluted with sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, the reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^{vi}) and a high proportion of weak D (D^u) phenotypes when using the slide, tube, gel card and microplate techniques.

| Item | Product Code | Size | Maximum Shelf Life |
|--|--------------|------|--------------------|
| Anti-D IgM Monoclonal - Standard Grade | 730010E | 10ml | 30 Months |

Lorne Monoclonal Anti-D Duoclone Standard Grade blood grouping reagent is a low protein, blended reagent containing human monoclonal IgM and IgG Anti-D, diluted in a phosphate buffer containing sodium chloride, bovine albumin and macromolecular potentiators. When typing patient samples, this reagent will directly agglutinate Rh D positive cells, including the majority of variants (but not D^{vi}) and a high proportion of weak D (D^u) phenotypes when using the slide, tube, gel card and microplate techniques. It will agglutinate D^{vi} cells in the IAT phase of testing.

| | | | |
|---|---------|------|-----------|
| Anti-D Duoclone Monoclonal - Standard Grade | 740010E | 10ml | 30 Months |
|---|---------|------|-----------|

All products available in bulk quantities

Lorne Laboratories Limited
Unit 1 Cutbush Park Industrial Estate
Danehill, Lower Earley
Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264

Email: info@lornelabs.com www.lornelabs.com

Published October 2018

In Vitro Diagnostic Medical Device
For professional use only

Hypochlorite solution 0.5%

| REF | Name | Packaging size |
|------|----------------------------|----------------|
| 3917 | Hypochlorite solution 0.5% | 1000ml |

Intended purpose

Hypochlorite solution 0.5% is a non – sterile reagent dedicated for intensive cleaning, rinsing and washing of hematology analyzers (capillaries, tubing, chambers).

Principle

The use of reagents based on sodium hypochlorite is recommended by the manufactures of hematology analyzers to maintain equipment in proper working condition.

Hypochlorite solution 0.5% is cleaning agent “emergency cleaner”, used to remove residual specimen and the remaining reagents contained in the elements of the measuring system of hematology analyzers. It removes any residual cellular, proteinaceous deposits (serum albumins) and the remaining reagents that may appear during the operation of the analyzer.

Specimens (collection and preparation)

Not applicable.

Reagent preparation

This reagent is ready to use and can be applied straight from the bottle, no special reagent preparation is needed unless Operators manual for used instrument include it.

Procedure (instruction for use)*

Due to differences in the rules for implementing the procedures for cleaning in different types of analyzers, Hypochlorite solution 0.5% should be used according to instrument manufacturer`s instructions for use and should be connected as listed in the Operators manual for instrument.

Recommended models of instruments:

| Hypochlorite 0.5% Emergency Cleaner | Model of instrument* |
|--|---|
| | Abbott Cell-Dyn 1800, 1700, 1600, 1300 |
| | ABX Pentra80, 60, 60C, 60C+, 60MS, ABX Micros 60, ES60, 45 |
| | Benesphera™ H32, Benesphera™ H32 VET |
| | Beckman Coulter AcT 5™, Beckman Coulter AcT Diff™, AcT Diff 2™, Beckman CoulterAcT 8™, AcT 10™ |
| | Diatron Abacus, Diatron Abacus Junior VET |
| | Drew Excell 18 (BT2100) |
| | Erma PCE-210 |
| | Hospitex Hemascreen 18 |
| | HTI Micros CC18 |
| | Medonic CA620-20, CA620-16, CA530-16 |
| | Melet-Schloesing MS9, MS4, MS8, MS8 VET |
| | Mindray BC-3200, BC-3000 Plus, BC-2800, BC-2300, BC-2000 |
| | Nihon Kohden Celltac™ F MEK-8222K, Celltac ES MEK-7300, Celltac E MEK-7222K, °Celltac™ α MEK-8118K + QA-810V |
| | Seac H20 Genius, SEAC H12 |
| | Sysmex K4500, Sysmex K1000, Sysmex KX21, KX21-N |

Composition (in water)

| Component | Concentration |
|---------------------|---------------|
| Sodium hypochlorite | < 0,6 % |
| Sodium hydroxide | < 0,5% |
| Fragrant | < 0.5 % |

Storage and shelf life

Store in temperature 2-30°C.

The shelf life of Hypochlorite solution 0.5% is 18 months from manufacturing date, if stored at the prescribed temperature range.

Do not use reagent beyond the expiration date printed on label.

Warnings and precautions

For in vitro diagnostic use
For professional use only

Hypochlorite solution 0.5% meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Avantor Performance Materials Poland S.A.
ul. Sowińskiego 11, 44-101 Gliwice, Poland
Tel: +48 32 2392 000

| | |
|---|---|
|  | H315: Causes skin irritation |
| | H319: Causes serious eye irritation |
| | H412: Harmful to aquatic life with long lasting effects. |
| PREVENTION | P264: Wash thoroughly after handling |
| | P280: Wear protective gloves/protective clothing/eye protection/face protection |
| | P273: Avoid release to the environment |
| RESPONSE | P332+P313: If skin irritation occurs: Get medical advice/attention |
| | P362+P364: Take off contaminated clothing and wash it before reuse |
| | P337+P313: If eye irritation persists: Get medical advice/attention |

For further information please refer to Master Safety Data Sheet.

Limitations of use

Do not use reagents with visible physical or chemical changes (color, turbidity) or in case of direct packaging damage.

Please refer to Operators manual for instrument for information about any additional limitation of use.

*The information contained herein has not been approved by analyzers manufacturers, it is recommendation for use only. Always refer to the user manual provided with the equipment at issue.

Disposal information

Dispose of contents/container to an appropriate treatment and disposal facility in accordance with applicable laws and regulations, and product characteristic at time of disposal.

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ELITECH CLINICAL SYSTEMS SAS
Zone Industrielle
61500 SEES FRANCE

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Design, production, control and sales of clinical chemistry products intended to be used for in vitro diagnostics. Validation of the combination reagents and analyzers. Distribution of clinical chemistry analyzers and products for in vitro diagnostics.

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ELITech Clinical Systems SAS
Zone industrielle - 61500 SEES - FRA

est conforme aux exigences des normes internationales
complies with the requirements of the international standards

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
GMED N° 10462-8

Ce certificat est délivré selon les règles de certification GMED / This certificate is issued according to the rules of GMED certification

Renouvelle le certificat 10462-7

On behalf of the President
Marjorie PERRIMON
Certification Director

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DECLARATION OF EC CONFORMITY

We, ELITech Clinical Systems SAS, Zone Industrielle 61500 SEES France, hereby certify, under our own responsibility, that the reagents such as listed attached (6 pages), conform to the essential requirements of appendices I and III of European Directive 98/79/EC, relating to in vitro diagnostic medical devices and to the public health code.

These devices are classified in the "other device" category since they do not belong neither to list A or list B of annex II nor to self-testing class.

This declaration is based on the contents of each technical file and is supported by the certification of our quality system according to the standard NF EN ISO 13485 : 2016 (Certification valid until July 27th, 2026).

We, ELITech Clinical Systems SAS, Zone Industrielle 61500 SEES France, hereby certify electrodes; conform to Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, including Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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Intended Use

Spectrum Diagnostics bilirubin reagent is intended for the in-vitro quantitative, diagnostic determination of bilirubin in human serum on both **automated and manual systems**.

Background

The average level of the bilirubin produced in humans from different sources ranges between 250 to 300 mg/day, of which 85% is derived from the heme moiety of the haemoglobin released from senescent erythrocytes that are destroyed in the reticuloendothelial system. The remaining 15 % is produced from erythrocytes destroyed in the bone marrow and from catabolism of other heme containing proteins such as cytochromes and myoglobin.

After it is produced in the peripheral tissues, bilirubin is transported to the liver in association with albumin. In the liver, bilirubin is conjugated with glucuronic acid for solubilization and subsequent transport through the bile duct and elimination via the digestive tract. Disease or conditions which, through hemolytic processes, produce bilirubin faster than the liver can metabolize it, cause the levels of unconjugated (indirect) bilirubin to increase in the circulation. Bile duct obstruction or damage to hepatocellular structure causes increases in the levels of both conjugated (direct) and unconjugated (indirect) bilirubin in the circulation.

Method

Colorimetric Diazo method.

Assay Principle

The total bilirubin concentration is determined in presence of caffeine by the reaction with diazotized sulphanilic acid to produce an intensely colored diazo dye (560-600 nm). The intensity of color of this dye formed is proportional to the concentration of total bilirubin.

Direct bilirubin is determined in absence of caffeine by the direct reaction with diazotized sulphanilic acid to form red-colored azobilirubin, the color intensity of which measured at 546 nm is proportional to the concentration of the direct bilirubin in the sample.

Reagents

Reagent 1 (R1)
Sulfanilic acid 31.0 mmol/l
HCL 0.20 N

Reagent 2 (R2)
Sodium nitrite 28.0 mmol/l

Reagent 3 (R3)
Caffeine 0.28 mol/l
Sodium benzoate 0.55 mol/l

Reagent 4 (R4)
Tartarate 0.99 mol/l
Sodium hydroxide 2.0 N

Reagent 4 contains caustic material.

Corrosive (C)

R35 Causes severe burns.

R41 Risk of serious damage to eyes.

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 After contact with skin, wash immediately with plenty of soap and water.

For further information, refer to the Bilirubin reagent material safety data sheet.

SYMBOLS IN PRODUCT LABELLING

| | | | |
|---|------------------------------|---|---------------------------------------|
|  | Authorised Representative |  | Use by/Expiration Date |
|  | For in-vitro diagnostic use |  | CAUTION. Consult instructions for use |
|  | Batch Code/Lot number |  | Manufactured by |
|  | Catalogue Number |  | (Xi) - Irritant |
|  | Consult instructions for use |  | Temperature Limitation |

Precautions and Warnings

Do not ingest or inhale. In case of contact with eyes or skin; rinse immediately with plenty of soap and water. In case of severe injuries; seek medical advice immediately.

Reagent Preparation, Storage and Stability

Spectrum bilirubin reagents are supplied ready-to-use and stable up to the expiry date labeled on the bottles when stored at room temperature. The opened vials are stable for 6 months at the specified temperature if contamination is avoided.

Deterioration

Do not use the Spectrum bilirubin reagents if precipitate forms. Failure to recover control values within the assigned range may be an indication of reagent deterioration.

Specimen Collection and Preservation

Avoid exposure of the specimen to light. If plasma is used, only heparin and oxalate plasma are suitable. Other anticoagulants should not be used. The average half-life of total bilirubin and direct bilirubin in serum is 17 days and few hours respectively.

Stability:

| | -20 °C | 4 – 8 °C | 20 – 25 °C |
|--------|----------|----------|------------|
| Total | 6 months | 7 days | 1 day |
| Direct | 6 months | 7 days | 2 days |

Procedure

Total Bilirubin

| | Sample blank | Sample |
|-----------|--------------|--------|
| Reagent 1 | 200 µl | 200 µl |
| Reagent 2 | ---- | 1 drop |
| Reagent 3 | 1.0 ml | 1.0 ml |
| Sample | 200 µl | 200 µl |

Mix and incubate for 10 minutes at 20 – 25 °C. then add;

| | | |
|-----------|--------|--------|
| Reagent 4 | 1.0 ml | 1.0 ml |
|-----------|--------|--------|

Mix and incubate for 5 minutes at 20 – 25 °C. Measure absorbance of sample (A_{sample}) against sample blank at 578 nm(560 - 600 nm) The color intensity is stable for 30 minutes.

Direct Bilirubin

| | Sample blank | Sample |
|------------------|--------------|--------|
| Reagent 1 | 200 µl | 200 µl |
| Reagent 2 | ---- | 1 drop |
| Saline 0.9% NaCl | 2.0 ml | 2.0 ml |
| Sample | 200 µl | 200 µl |

Mix and incubate for exactly 5 minutes at 20 – 25 °C. Measure absorbance of sample (A_{sample}) against sample blank at 546 nm (530 - 560 nm).

Calculation

Total bilirubin (mg/dl) = $A_{\text{Sample}} \times 10.8$

Direct bilirubin (mg/dl) = $A_{\text{Sample}} \times 14.4$

Quality Control

Normal and abnormal commercial control serum of known concentrations should be analyzed with each run.

Performance Characteristics

Precision

Within run (Repeatability)

| | Total | | Direct | |
|--------------|---------|---------|---------|---------|
| | Level 1 | Level 2 | Level 1 | Level 2 |
| n | 20 | 20 | 20 | 20 |
| Mean (mg/dL) | 0.79 | 4.37 | 0.299 | 0.77 |
| SD | 0.016 | 0.18 | 0.016 | 0.057 |
| CV% | 2.03 | 4.12 | 5.35 | 7.4 |

Run to run (Reproducibility)

| | Total | | Direct | |
|--------------|---------|---------|---------|---------|
| | Level 1 | Level 2 | Level 1 | Level 2 |
| n | 20 | 20 | 20 | 20 |
| Mean (mg/dL) | 0.82 | 4.52 | 0.32 | 0.82 |
| SD | 0.02 | 0.17 | 0.023 | 0.062 |
| CV% | 2.44 | 3.76 | 7.19 | 7.56 |

Methods Comparison

A comparison between Spectrum Diagnostics Bilirubin and a commercial reagent of the same methodology was performed on 20 human sera. A correlation of 0.975 was obtained.

Sensitivity

When run as recommended, the sensitivity of this assay is 0.1 mg/dL (1.7 µmol/L) for both total and direct bilirubin.

Linearity

The reaction is linear up to a total bilirubin concentration of 30 mg/dL (513 µmol/L) and a direct bilirubin concentration of 10 mg/dL (171 µmol/L). Specimens showing higher concentration should be diluted 1+4 with physiological saline and repeat the assay (result×5).

Interfering substances

Haemolysis

Avoid haemolysis since it interferes with the test.

Lipemia

Lipemic specimens interfere with the test.

Drugs

Theophylline and propranolol may cause artificially low total bilirubin levels.

Expected Values

Total Bilirubin

Adults and infants >1 month < 0.2-1.0 mg/dL (3.4 -17 µmol/l)
Newborns premature (3-5 d) 10-14 mg/dL (171-239 µmol/l)

Newborns:

(3-5 d) 4.0 - 8.0 mg/dL (68-137 µmol/l)
(<48 h) 6.0 - 10.0 mg/dL (103-171 µmol/l)
(<24 h) 2.0 - 6.0 mg/dL (34-103 µmol/l)

Direct Bilirubin 0 - 0.3 mg/dL (0 - 51 µmol/L)

Spectrum Diagnostics does not interpret the results of a clinical laboratory procedure; interpretation of the results is considered the responsibility of qualified medical personnel. All indications of clinical significance are supported by literature references.

Analytical Range

Total bilirubin : 0.1 – 30 mg/dL (1.7 – 513 µmol/L)
Direct bilirubin : 0.1 – 10 mg/dL (1.7 – 171 µmol/L)

Waste Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

S56: dispose of this material and its container at hazardous or special waste collection point.

S57: use appropriate container to avoid environmental contamination.

S61: avoid release in environment. refer to special instructions/safety data sheets.

References

- Balistreri WF, Shaw LM. Liver function. In: Tietz NW, ed. Fundamentals of clinical chemistry. 3rd ed. Philadelphia: WB Saunders; 1987:729-761.
- Malloy HT, Evelyn KA. The determination of bilirubin with the photoelectric colorimetric method. J Biol Chem. 1937;119:481-490.
- Tietz NW, ed. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders; 1995:268-273.

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Egyptian Company for Biotechnology (S.A.E)
MDSS GmbH