OCLUDERE DE DUCT AMPLATZER™

INDICAȚII

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AMPLATZER™ DUCT OCCLUDERS

INDICATIONS

The Amplatzer™ Duct Occluder and Amplatzer™ Duct Occluder II are percutaneous, transcatheter occlusion devices intended for the nonsurgical closure of a patent ductus arteriosus (PDA) in patients with a weight of 6 kg or larger.


The Amplatzer Piccolo™ Occluder is a percutaneous, transcatheter occlusion device intended for the nonsurgical closure of a patent ductus arteriosus (PDA) in patients with a weight 700 grams and up at time of the procedure.

MRI SAFETY INFORMATION

Through nonclinical testing, Amplatzer™ devices have been shown to be MR Conditional. Refer to the appropriate Instructions for Use to obtain more detailed MRI scanning information.

MR image quality may be compromised if the area of interest is in the same area as or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

(10, 11, 12)

DEVICE SPECIFICATIONS AND DELIVERY SYSTEM COMPATIBILITY


Amplatzer™ Duct Occluder

Amplatzer™ TorqVue™ 180° Delivery System

Model/ Reorder Number	Device Diameter at Descending Aorta (mm) A	Device Diameter at Pulmonary Artery (mm) B	Retention Skirt Diameter (mm) C	Device Length (mm) D	Model/Reorder Number	Minimum Recommended Sheath Size	Sheath Inner Diameter (mm [inch])	Sheath Outer Diameter (mm [inch])
9-PDA-003	5	4	9	5	9-ITV05F180/60	5 F	1.83 [0.072]	2.51 [0.099]
9-PDA-004	6	4	10	7				
9-PDA-005	8	6	12	7				
9-PDA-006	10	8	16	8				
9-PDA-007	12	10	18	8	9-ITV07F180/80	7 F	2.44 [0.096]	3.18 [0.125]
9-PDA-008	14	12	20	8				
9-PDA-009	16	14	22	8				

(13)

(2)

(10, 11, 12)

DEVICE SPECIFICATIONS AND DELIVERY SYSTEM COMPATIBILITY

Amplatzer™ Duct Occluder II

Amplatzer™ TorqVue™ LP Delivery System

Model/Reorder Number	Waist Diameter (mm) A	Device Length (mm) B	Disc Diameter (mm) C	Model/Reorder Number	Minimum Recommended Catheter Size	Catheter Inner Diameter (mm [inch])	Catheter Outer Diameter (mm [inch])
9-PDA2-03-04	3	4	9				
9-PDA2-03-06	3	6	9	9-TVLP4F90/060 or 9-TVLP4F90/080	4 F	1.17 [0.046]	1.40 [0.055]
9-PDA2-04-04	4	4	10				
9-PDA2-04-06	4	6	10				
9-PDA2-05-04	5	4	11				
9-PDA2-05-06	5	6	11	9-TVLP5F90/060 or 9-TVLP5F90/080	5 F	1.50 [0.059]	1.73 [0.068]
9-PDA2-06-04	6	4	12				
9-PDA2-06-06	6	6	12				

(13)

DEVICE SPECIFICATIONS AND DELIVERY SYSTEM COMPATIBILITY

Amplatzer™ Piccolo™ Occluder

Amplatzer™ TorqVue™ LP Catheter

Model/Reorder Number	Waist Diameter (mm) A	Length Between Retention Discs (mm) B	Disc Diameter (mm) C	Model/Reorder Number	Minimum Recommended Catheter Size	Catheter Inner Diameter (mm [inch])	Catheter Outer Diameter (mm [inch])
9-PDAP-03-02-L	4.00(0.157)	3.00(0.118)	2.00(0.079)				
9-PDAP-03-04-L	4.00(0.157)	3.00(0.118)	4.00(0.157)				
9-PDAP-03-06-L	4.00(0.157)	3.00(0.118)	6.00(0.236)				
9-PDAP-04-02-L	5.25(0.207)	4.00(0.157)	2.00(0.079)				
9-PDAP-04-04-L	5.25(0.207)	4.00(0.157)	4.00(0.157)	9-TVLPC4F90/080	4 F	1.17 [0.046]	1.4 [0.055]
9-PDAP-04-06-L	5.25(0.207)	4.00(0.157)	6.00(0.236)				
9-PDAP-05-02-L	6.50(0.256)	5.00(0.197)	2.00(0.079)				
9-PDAP-05-04-L	6.50(0.256)	5.00(0.197)	4.00(0.157)				
9-PDAP-05-06-L	6.50(0.256)	5.00(0.197)	6.00(0.236)				

ANCILLARY PRODUCT SPECIFICATIONS

Amplatzer™ Guidewire

Model/Reorder Number	Diameter (inch)	Body	Tip Description	Useable Length (cm)
9-GW-001	0.035	Super Stiff	7.5 mm, Modified J-tip	260

(14)

LATEX-FREE INFORMATION

These Amplatzer™ products do not contain latex.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Always check the regulatory status of the device in your region.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Abbott

3200 Lakeside Dr, Santa Clara, CA 95054, USA Tel: 1.800.227.9902
www.cardiovascular.abbott

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

© 2020 Abbott. All Rights Reserved.

37971 MAT-2002324 v2.0 | Item approved for Global OUS only.

THE INTAGLIO™ WIRE TREATMENT

Amplatzer™ Structural Interventions

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.
Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

Abbott Vascular International BVBA

Park Lake, Culliganlaan 2b, 1831 Diegem, Belgium
www.cardiovascular.abbott

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

© 2019 Abbott. All Rights Reserved. 34513-SJM-AMPLP-0519-0121c | Item approved for Chile only.
SL-16047-01 09/2019

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Clinical Insights

AMPLATZER™ STRUCTURAL INTERVENTIONS

THE INTAGLIO™ WIRE TREATMENT

(4)

KEY MESSAGES

- Intaglio is a chemical etch treatment of Nitinol wire used in the manufacturing of Amplatzer Devices.
- Although Amplatzer products have always met the standards for nickel leaching, they were improved through the use of Intaglio wire treatment.
- The Intaglio wire treatment reduces the amount of nickel that may be leached from the device by more than 95% depending on the model. See next page for Nickel Reduction from Intaglio in Table 1 and PFO example in Graph 1.
- Intaglio is part of Abbott's commitment to increasing quality of all products.

INTRODUCTION


The Amplatzer Occlusion Devices have historically been manufactured using Black Oxide Nitinol wire. In 2014, the Intaglio chemically etched finish was introduced to Amplatzer Occlusion devices.

SIGNIFICANCE OF INTAGLIO

Amplatzer products have been characterized utilizing a robust in-vitro chemical leaching study to assess the materials-related toxicity risks associated with nickel release after implantation. The conversion to the Intaglio chemical etching process has further reduced these toxicological risk factors with respect to nickel leaching. The mechanical performance of the device remains unchanged with the conversion to chemically etched wire.

INTAGLIO TREATMENT

The Intaglio treatment is a secondary step to the Nitinol material previously utilized in the production of Amplatzer products with a Black Oxide finish. The metal wire undergoes a chemical etch as a final process in the manufacturing of the product resulting in a thinner, more consistent oxide layer on the device. The mechanical properties of the wire and finished device are not affected by this process so the proven clinical performance remains unaltered. Devices treated with Intaglio have a slight blue or purple hue compared to previous devices with a Black Oxide finish.

CONCLUSIONS

Amplatzer Occluder Devices and the introduction of Intaglio chemical etching are part of Abbott's commitment to continuous product improvement. Although Amplatzer products have an excellent safety profile and have always met the standards for nickel leaching, they were improved through the use of Intaglio wire treatment. Testing of Amplatzer Occluder devices with Intaglio chemical etching has demonstrated in multiple devices more than 95% (varying depending on device model) further reduction in nickel leaching compared to devices manufactured with a Black Oxide finish. Lower levels of nickel is intuitively better for patients and may reduce the potential for an adverse reaction to nickel.


TABLE 1: NICKEL REDUCTION FROM INTAGLIO IN AMPLATZER OCCLUSION PRODUCTS*

Product	Black Oxide, Max Ni Leach 1st Day [µg]	Black Oxide, Max Ni Leach 10th Day [µg]	Intaglio, Max Ni Leach 1st Day [µg]	Intaglio, Max Ni Leach 10th Day [µg]	1st Day Reduction in Nickel Leach with Intaglio (%)	10th Day Reduction in Nickel Leach with Intaglio (%)
PFO	77.1	18.4	1.79	0.17	97.7%	99.1%
ADOII	45.2	7.0	0.74	0.08	98.4%	98.9%
ASD	78.1	21.8	3.36	1.05	95.7%	95.2%

- The Intaglio wire treatment reduces the amount of nickel that may be leached from the device by more than 95% depending on the model.
- Lower levels of nickel increase safety and may reduce the potential for an adverse reaction to nickel

*Chemical Characterization Testing data on file at Abbott. For product specific Nickel precautions, please check the Instructions for Use.

GRAPH 1: NICKEL LEACH REDUCTION WITH INTAGLIO - PFO

AMPLATZER™ TORQVUE™ LP DELIVERY SYSTEM

Ordering Information

(16)

Contents: 1 delivery cable, delivery catheter, hemostasis valve, loader and plastic vise

Model/ Reorder No	Cath Sz	Tip Angle	Usable Lgth (cm)	Inner Diam of Catheter (mm/in)	Outer Diam of Catheter (mm/in)	Delivery Cable Lgth (cm)
9-TVLP4F90/060	4 F	90°	60	1.17/0.046	1.40/0.055	160
9-TVLP4F90/080	4 F	90°	80	1.17/0.046	1.40/0.055	195
9-TVLP5F90/060	5 F	90°	60	1.50/0.059	1.73/0.068	160
9-TVLP5F90/080	5 F	90°	80	1.50/0.059	1.73/0.068	195

AMPLATZER™ GUIDEWIRES

Ordering Information

Contents: 1 guidewire

Model/ Reorder No	Diam (in)	Body	Floppy Tip Lgth (cms)	Tip Description	Usable Lgth (cm)
9-GW-001	0.035	Super Stiff	6	7.5 mm, Modified J-tip	260
9-GW-002	0.035	Super Stiff	5	1.5 mm, Modified J-tip	260
9-GW-003	0.035	Super Stiff	20	6 mm, J-tip	300
9-GW-004*	0.035	Soft Tip, Fixed Core	NA	6 mm, J-tip	300

*9-GW-004, also referred to as “Noodlewire,” is a soft tip, flexible guidewire recommended for establishing an arterial-venous loop, facilitating closure of ventricular septal defects.

STRUCTURAL HEART THERAPY

AMPLATZER™ FAMILY OF DUCT OCCLUDERS

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**.
Always check the regulatory status of the device in your region.

MULTIPLE DEVICE DESIGNS

FOR CONFIDENT AND COMPLETE PATENT DUCTUS ARTERIOSUS OCCLUSION

Patent Ductus Arteriosus (PDA) is one of the most common of all congenital anomalies. The Amplatzer™ Family of Duct Occluders offers an extensive line of devices and a simple treatment option for PDA occlusion. Each duct occluder has been designed to occlude PDAs of various shapes and sizes. With an easy-to-perform deployment procedure, the implanting physician can be confident in placing the device with precision.

AMPLATZER™ DUCT OCCLUDER

TREATMENT RANGE

Ductus Diameter: 2.0-12.0 mm

Ductus Length: Greater than 3.0 mm

CONFIDENCE IN LARGE PDA CLOSURE

The Amplatzer™ Duct Occluder is designed to accommodate large PDAs with a single device-which may reduce the complexity of the procedure.

DESIGNED TO REDUCE RISK

Retention Skirt: Built to provide secure positioning within the ductus and prevent embolization³

Cone Shape: Engineered to conform to tapering ducts and promote occlusion³

DESIGNED TO PROVIDE COMPLETE OCCLUSION

Complete Closure of PDA: With over 98% of patients showing complete occlusion at six months¹

Polyester Material: Promotes occlusion and tissue in-growth

AMPLATZER™ DUCT OCCLUDER II

TREATMENT RANGE

Max Ductus Diameter: 12.0 mm

Max Ductus Length: 5.5 mm

VERSATILITY IN PDA CLOSURE

The Amplatzer™ Duct Occluder II conforms to a variety of ducts while achieving complete closure from a pulmonary or aortic approach.²

DESIGNED FOR FLEXIBILITY

Compliant, Multi-layered Mesh and Dual Articulating Discs: Provide high conformability to treat most classifications of patent ductus arteriosus³

Symmetrical Design: Offers procedural flexibility to choose an aortic or pulmonary artery approach³

PROMOTES RAPID OCCLUSION

Six Planes of Occlusion: Create full crosssectional coverage designed for rapid², complete closure²

Fabric-free Technology: Allows for delivery through a low-profile catheter while maintaining a high rate of occlusion²

AMPLATZER PICCOLO™ OCCLUDER

TREATMENT RANGE

Ductus Diameter: Less than 4.0 mm

Ductus Length: Greater than 3.0 mm

PROVEN PDA CLOSURE FOR PATIENTS 700G AND UP

The Amplatzer Piccolo™ Occluder is the only PDA closure solution indicated for premature infants ($\geq 700\text{g}$ + ≥ 3 days old) and proven to deliver safe and effective closure.

RELIABILITY BY DESIGN

Tightly woven single-layer mesh: Designed to minimize residual shunt immediately after placement

Engineered for structure and strength: Occludes ducts while minimizing protrusion into surrounding vasculature

Intaglio wire treatment: Designed to reduce nickel leaching

SMOOTH DELIVERY IN EVEN THE MOST CHALLENGING MORPHOLOGIES

Predictable placement: Disc size and shape designed for predictable positioning in the duct

Deliverable in 4 F through systems: 4 F catheter facilitates delivery small vasculatures

Procedural flexibility: Symmetrical design offers procedural flexibility to choose an antegrade (venous) or retrograde (arterial) approach.

For infants $\leq 2\text{kg}$ a venous approach is recommended

ORDERING INFORMATION

AMPLATZER™ Duct Occluder

Model / Reorder Number	Device Diameter at Descending Aorta (mm)	Device Diameter at Pulmonary Artery (mm)	Retention Skirt (mm)	Length (mm)
9-PDA-003	5	4	9	5
9-PDA-004	6	4	10	7
9-PDA-005	8	6	12	7
9-PDA-006	10	8	16	8
9-PDA-007	12	10	18	8
9-PDA-008	14	12	20	8
9-PDA-009	16	14	22	8

AMPLATZER™ Duct Occluder II

Model / Reorder Number	Waist Diameter (mm)	Device / Nominal Length (mm)	Disc Diameter (mm)
9-PDA2-03-04	3	4	9
9-PDA2-03-06	3	6	9
9-PDA2-04-04	4	4	10
9-PDA2-04-06	4	6	10
9-PDA2-05-04	5	4	11
9-PDA2-05-06	5	6	11
9-PDA2-06-04	6	4	12
9-PDA2-06-06	6	6	12

AMPLATZER PICCOLO™

Model / Reorder Number	Waist Diameter (mm)	Length between discs (mm)	Disc Diameter (mm)
9-PDAP-03-02-L	3	2	4
9-PDAP-03-04-L	3	4	4
9-PDAP-03-06-L	3	6	4
9-PDAP-04-02-L	4	2	5.25
9-PDAP-04-04-L	4	4	5.25
9-PDAP-04-06-L	4	6	5.25
9-PDAP-05-02-L	5	2	6.50
9-PDAP-05-04-L	5	4	6.50
9-PDAP-05-06-L	5	6	6.50

REFERENCES:

1. AMPLATZER Duct Occluder Instructions for Use.
2. Saliba Z, El-Rassi I, Abi-Warde MT, et al. The Amplatzer Duct Occluder II: A New Device for Percutaneous Ductus Arteriosus Closure. J Interv Cardiol. 2009;22(6):496-502.
3. Test(s) performed by and data on file at Abbott.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**.
Always check the regulatory status of the device in your region.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

Abbott

3200 Lakeside Dr, Santa Clara, CA 95054, USA

www.structuralheartsolutions.com

™ Indicates a trademark of the Abbott group of companies.

© 2021 Abbott. All Rights Reserved. MAT-2004321 v1.0 | Item approved for Global OUS use only.