# ClinChem Multi Control (level 2)

## Order Information

Level	Cat. No.	Package size
level 2	105-009120-00	6×5 mL

## Intended Purpose

ClinChem Multi Control is used for quality control of quantitative determination of routine chemistry analytes by monitoring accuracy and precision of Mindray BS series chemistry analyzers and test ability of clinical laboratory.

## Summary

ClinChem Multi Control contains the following analytes:

Albumin (ALB)	Alkaline Phosphatase (ALP)	Alanine Aminotransferase (ALT)
$\alpha$ -Amylase( $\alpha$ -AMY)	Aspartate Aminotransferase (AST)	Direct Bilirubin (D-Bil)
Total Bilirubin (T-Bil)	Calcium (Ca)	Total Cholesterol (TC)
Creatine Kinase (CK)	Creatinine (CREA)	Glucose (Glu)
Gamma-Glutamyltransferase ( $\gamma$ -GT)	$\alpha$ -Hydroxybutyrate Dehydrogenase ( $\alpha$ -HBDH)	Lactate Dehydrogenase (LDH)
Magnesium (Mg)	Phosphorus (P)	Total Protein (TP)
Triglycerides (TG)	Urea (UREA)	Uric Acid (UA)
Iron (Fe)	Cholinesterase (CHE)	Lipase (LIP)
Sodium (Na <sup>+</sup> )	Potassium(K <sup>+</sup> )	Chlorine (Cl-)
Immunoglobulin A (IgA)	Immunoglobulin G (IgG)	Complement C3(C3)
Immunoglobulin M (IgM)	C-Reactive Protein (CRP)	Complement C4(C4)
HDL-Cholesterol (HDL-C)	LDL-Cholesterol (LDL-C)	Apolipoprotein A1 (ApoA1)
Apolipoprotein B (ApoB)	Prealbumin (PA)	Creatine Kinase-MB (CK-MB)
Antistreptolysin O (ASO)	Transferrin (TRF)	Ferritin (FER)
Unsaturated Iron Binding Capacity (UIBC)		

Mindray BS measurement system is composed of Mindray BS series chemistry analyzers, Mindray reagent kits, calibrators and control materials.

## Components

ClinChem Multi Control is a lyophilized control based on human serum with less than 0.095% sodium azide as preservative, which contains 2 levels. The

concentrations or activities of the control components are lot-specific with the matched models listed in the value sheet and almost at normal or/and pathological levels for the methods used.

### Storage and stability

The control is stable up to expiration date indicated on the label when stored in unopened vial at 2-8°C and protected from light.

Once reconstituted, please store control tightly capped and protected from light when not in use.

Once thawed to room temperature, the control should not be re-frozen. The stability information is as below.

Temperatures of storage	Stability		
	Universal components	Bil-T / Bil-D / UIBC	ALT
15-25°C	12 hours	8 hours	12 hours
2-8°C	5 days	1 day	5 days
(-25)-(-15)°C	4 weeks	2 weeks	2 weeks

\*The supplementary cautions for storage and stability:

1. Alkaline Phosphatase activity in the reconstituted control will rise over the stability period. It is recommended that the reconstituted control is allowed to stand for 1 hour at 15-25°C before measurement.
2. Alanine Aminotransferase activity in the reconstituted control is affected by storage temperature and storage time. It is recommended to use control as soon as possible after reconstitution and shorten the storage time.
3. Total bilirubin, direct bilirubin in the control is light sensitive and it is recommended that the control serum is stored in the dark.

### Preparation

1. It is recommended to take the vial out of the freezer and allow reaching room temperature.
2. Tap the vertically positioned vial gently and ensure that the lyophilized material is at the bottom of the vial.
3. Remove the screw cap and rubber stopper carefully, avoiding any loss of lyophilizate.
4. Reconstitute by adding 5.0 mL distilled/deionized water to the side of the vial slowly and exactly.
5. Carefully replace the rubber stopper and leave to stand for 30 minutes at room temperature.
6. During the standing time, mix contents by inverting several times and swirling gently to ensure that all components of lyophilizate are dissolved. Avoid the formation of foam.

7. According to the quality control procedure described in the operation manual for the analyzers.
8. Dispense the required volume into a sample cup, and then run QC test you have requested or have it analyzed in the same way as patient samples.
9. Quality control is recommended to run daily in parallel with the patient samples, after new calibration or reagent cartridge change, after specific maintenance or trouble shooting procedures.
10. Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

### **Materials required but not provided**

1. General laboratory materials: NaCl solution 9 g/L (saline), distilled/deionized water.
2. Mindray reagent kit and Mindray calibrator. The catalogue number of reagent and calibrator used with ClinChem Multi Control can be obtained from the corresponding Mindray reagent kit and calibrator instruction for use.
3. Mindray BS series chemistry analyzers, and General laboratory equipments.

### **Performance Characteristics**

#### **■ Assay Values and Accuracy**

The control assay values (the target value and range) determined by Mindray standard transfer procedure and routine method are listed in the value sheet. The target value was obtained from Mindray measurement system and the range was calculated as the target values  $\pm 3$  standard deviations. The control values are specific for different lot numbers and different model of chemistry analyzers, so please check the lot number and model of analyzer before use.

Mindray systems (BS series chemistry analyzers/Mindray reagents) are calibrated by Mindray calibrators, then test the control on Mindray systems, the assayed values should be within the defined range.

#### **Result interpretation**

If the control result falls outside the range, the measurement system should be checked. For instance, position of reagent or sample placed on the analyzer; expiration date or storage condition of the calibrator, reagent and control; the parameter setting in software; performance of the analyzer; or the calibration process. Each laboratory should establish its own internal quality control scheme and procedure for corrective action.

#### **Warnings and precaution**

1. For in vitro diagnostic use only. For laboratory professional use.

2. Please Take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the controls with damaged packages. The controls avoid direct exposure to sunlight. The results can't be assured when stored at inappropriate condition.
4. If unintentionally opened before used, please reconstitute to use as soon as possible.
5. Do not mix controls with different lots and bottles.  
Do not use the control beyond the expiration date and the in-use date.  
Do not mix new opened control with in-use control.  
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
8. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the control accidentally enters the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. The control was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen<sup>1</sup>.
13. All human material should be considered potentially infectious.
14. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## References

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington:2007.

## Graphical symbols

In Vitro Diagnostic  
medical device

Unique device  
identifier

European  
Conformity

Consult Instructions  
For use

Use-by  
date

Authorized representative in  
the European Community

Batch Code

Temperature  
limit

Manufacturer

Catalogue  
number

Control

Biological Risks

Keep away from sunlight  
Indicates a medical device that needs  
protection from light sources

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-81888998; **Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH(Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés
**Fabricant** : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
**Représentant en Europe** : Shanghai International Holding Corp. GmbH (Europe)