

GE Healthcare

# Xeleris 4.0

Product Data Sheet

SECRET DE AFACERI

CONFIDENTIAL

## Overview

Xeleris™ 4.0 functional imaging workstation is a Nuclear Medicine, PET, NM/CT, and PET/CT processing, analysis, and review system. Designed with productivity in mind, it can accelerate workflow and provides a powerful clinical diagnostic tool to the medical imaging community.

Combining streamlined workflow with a comprehensive clinical library and extensive networking capabilities on a functional imaging workstation, Xeleris 4.0 is at the nucleus of productivity in the clinical imaging department. Utilizing the GE Healthcare-wide graphical user interface, Xeleris 4.0 is the processing and review platform of the Discovery™, Optima™ and Brivo™ NM and NM/CT series, Infinia™ Hawkeye™ 4, Ventri, Discovery PET/CT 600 series, and all other molecular imaging cameras in GE Healthcare's current offering. Xeleris 4.0 provides the automated processing and connectivity necessary in today's demanding environment.

Xeleris 4.0 operates in a variety of client-server configurations. It features optional Xeleris Floating License (XFL) technology and Xeleris Suite for PACS. XFL allows patient studies and the Xeleris applications needed for their processing and review to be accessed across a LAN or a WAN network on a facility provided Windows®-based PC (minimum hardware and software requirements apply), thus providing an untethered solution for simultaneous use of Xeleris applications across the enterprise.

Xeleris Suite extends this power to GE Healthcare Centricity® PACS systems and other PACS systems supporting Open API integration.

Xeleris 4.0 is available in Full, Cardiac and Molecular Breast Imaging (MBI) configurations. These configurations of Xeleris can simplify workflow and may provide greater economy as only applications needed for specific imaging systems are provided. The Cardiac and Molecular Breast Imaging (MBI) configurations are fully upgradeable to full Xeleris 4.0 processing and review.

## Indications for Use

The system is intended for use by Nuclear Medicine (NM) or Radiology practitioners and referring physicians for display, processing, archiving, printing, reporting and networking of NM data, including planar scans (Static, Whole Body, Dynamic, Multi-Gated) and tomographic scans (SPECT, Gated SPECT, dedicated PET or Camera-based PET) acquired by gamma cameras or PET scanners.

The system can run on dedicated workstation or in a server-client configuration.

The NM or PET data can be coupled with registered and/or fused CT or MR scans, and with physiological signals in order to depict, localize, and/or quantify the distribution of radionuclide

tracers and anatomical structures in scanned body tissue for clinical diagnostic purposes.

DaTQUANT optional application enables visual evaluation and quantification of 123I-ioflupane (DaTscan™) images. The DaTQUANT Normal Database option enables quantification relative to normal population databases of 123I-ioflupane (DaTscan™) images.

These applications may assist in detection of loss of functional dopaminergic neuron terminals in the striatum, which is correlated with Parkinson disease.

Q.Lung application may aid physicians in:

- Diagnosis of Pulmonary Embolism (PE), Chronic Obstructive Pulmonary Disease (COPD), Emphysema and other lung deficiencies.
- Assess the fraction of total lung function provided by a lobe or whole lung for Lung cancer resection requiring removal of an entire lobe, bilobectomy or pneumonectomy.

Q.Brain\* allows the user to visualize and quantify relative changes in the brain's metabolic function or blood flow activity between a subject's images and controls, when used with radiopharmaceuticals approved by the regulatory authority in the country of use, which may be resulting from brain function alterations in:

- Epileptic seizures
- Dementia. Such as Alzheimer's disease, Lewy body dementia, Parkinson's disease with dementia, vascular dementia, and frontotemporal dementia.
- Inflammation
- Brain death
- Cerebrovascular disease such as Acute stroke, Chronic and acute ischemia
- Traumatic Brain Injury (TBI)

When integrated with the patient's clinical and diagnostic information, Q.Brain application may aid the physician in the interpretation of cognitive complaints, neuro-degenerative disease processes and brain injuries.

The Alcyone CFR application allows for the quantification of coronary vascular function by deriving Myocardial Blood Flow (MBF) and then calculating Coronary Flow Reserve (CFR) indices on data acquired on PET scanners and on stationary SPECT scanners with the capacity for dynamic SPECT imaging. These indices may add information to physicians using Myocardial Perfusion Imaging for the diagnosis of Coronary Artery Disease (CAD).

**SECRET DE AFACERI**  
**CONFIDENTIAL**

## Primary Benefits

### Streamlined workflow

#### Ignite

Xeleris 4.0 can enhance nuclear imaging productivity through Ignite operational flexibility and automated workflow. Xeleris 4.0 combines the speed of automated intuitive processing with the freedom to modify processing parameters (if necessary), helping to optimize study results without losing the benefits of automation.

Ignite, combined with GE nuclear imaging systems, can make most clinical scans as easy as 1, 2, 3:

1. Select the patient from the worklist on the camera
2. Position the patient and press Start (Ignite the process)
3. Review the results that will be automatically displayed without the need for further interaction.

Departmental efficiency is enhanced via workflow automation in the form of:

1. Enhanced Filtering - enables detailed and personalized data filtering and search functionality which enables faster access to patients' exams using smart filters and search options compared to manual search and filter. It allows user to find or filter patients' exams based on multiple criteria of modality, time, referring physician, data type, therefore helping improve productivity.
2. Auto Snap Shot- Screen capture (printing). This functionality enables automatic generation of DICOM Screen Capture (SC) and DICOM Multi-Frames SC (MFSC) from predefined application screens.
3. Auto Transfer and Archive - enables automatic transfer of patient data and Xeleris results series to other pre-defined multiple Xeleris Database (DB), DICOM or PACS destination, with no operator intervention.

These features improve workflow by reducing the required number of user interventions for performing data transfer or data printing. Manual data transfer or archiving or printing requires Multiple button clicks. Using the auto transfer/auto archive features enables transferring the data with no user intervention, thus workflow is improved. The user can still use the manual data transfer as well.

## Clinical library

Xeleris 4.0 features a comprehensive clinical library of user-friendly processing and review tools and protocols, covering all nuclear imaging needs, and providing the flexibility to customize protocols per user's requirements.

## Access and Integration

Xeleris 4.0 brings together the entire molecular imaging department, including GE Healthcare and most non-GE Healthcare nuclear imaging systems, providing effective solutions for current and future information technology needs. Xeleris 4.0 supports processing, archiving and review of data received from DICOM 3 compatible NM, PET, CT and MR data.

Xeleris 4.0 interfaces with nuclear imaging and PET systems, scanners, and workstations, including the Infinia, Optima, Brivo and Discovery line of hybrid imagers, Millennium and DS-series.

### DICOM Conformance

Xeleris 4.0 provides excellent connectivity tools allowing you to **share Nuclear Medicine/PET and Hybrid images with PACS and DICOM-conformant systems**, including full screen dynamic displays in Multiframe Secondary Capture format. See the Xeleris DICOM conformance statement for details.

### IHE Workflow

Xeleris 4.0 implements the IHE scheduled workflow and nuclear medicine image integration profiles (as image display and evidence creator actors) and the access to radiology information integration profile (as image display actor). See the Xeleris IHE integration statement for details.

### Remote Viewing

Xeleris can capture most images in PC compliant screen capture and movie formats (for easy drop in presentations, email or on the network).

### Common User Interface

Xeleris utilizes a highly advanced graphical user interface, common to both the processing and review workstation as well as the acquisition console. This common user interface can help shorten the operators' learning curve and helps promote departmental productivity and efficiency by helping to enhance the operators' performance through ease-of-use.

SECRET DEAFACER1

# Xeleris 4.0 Main Features

## Clinical Library

### Volumetrix™ MI

A consolidated tomographic data viewing and processing package for SPECT and PET data processing with or without anatomical data (CT or MR), capable of Multi-SPECT sessions pasting, follow-up on multiple studies and input of anatomical oblique slices. Its Integrated SPECT reconstruction and correction algorithms include interactive selection of optimal image reconstruction parameters and Hybrid QC (ACQC) tools to help improve alignment of SPECT and CT data sets for accurate attenuation correction and localization.

- Inputs anatomical oblique slices and features on-the-fly oblique creation/manipulation
- Facilitates dual isotope and study follow-up reading and accepts multiple SPECT, PET, CT and MR datasets
- Supports standard uptake value (SUV) calculation and display mode for PET (LBM, BSA, & BW methods)
- Flexible and user customized layouts
- Auto Condense mode to automatically condense slices to fit on one screen for printing or filming
- Multiple fusion modes, including weighted fusion, which allows the user to shift between the contribution of the PET and CT image to the fused results
- CT and PET images retain resolution
- Fast CT scroll

**Volumetrix MI includes the following packages:**

### MDC - Motion Detection and Correction

- Automated cardiac and general purpose SPECT motion correction
- Sinogram and linogram images for QC analysis
- Tools for manual adjustment and correction

### SPECT Options

- Adjustment of any of the reconstructions settings
- On-the-fly image correction setup including resolution recovery (option), attenuation and scatter.
- Filter selection with filter gallery review to determine the most effective filter and filter parameters. 2D Hanning, Metz, Butterworth, Hamming, and Wiener filters are available for FBP reconstruction. 3D Butterworth, Hanning and Gauss filters are available for OSEM reconstruction
- Chang uniform Attenuation Correction for brain SPECT
- Dual scatter energy window scatter correction for SPECT

## QC for hybrid NM/CT studies

Enabling effective alignment of hybrid NM and CT data sets in order to reduce potential inter-study misalignments. It can help reduce the need for study retakes and features:

- Automatic Hybrid QC (user overridden) of the NM FOV relative to Hybrid CT FOV
- Confirm or reject each FOV separately
- Correct or generate attenuation maps

## Multi-FOV Pasting

- Automatic pasting (user overridden) up to seven FOV's (covering the whole body)
- Preview pasted volume on SPECT/CT (NM & fused available)

## Registration

- Manual adjustment via shift/rotate
- Single landmark registration
- Regional registration
- Fusion & spots collection are supported on registered images

## Single CT Attenuation Correction

- Use of a single CT scan as opposed to multiple CT scans in some protocols (e.g. Lung Ventilation/Perfusion or consecutive imaging periods)
- Enhanced image quality, compared to non-attenuation corrected images.
- Patient dose reduction by eliminating the need for additional CT scans for purpose of Attenuation Correction (AC) of each SPECT Study.

Auto accept or user intervention are allowed with any one of the above packages

## Volumetrix Suite (option)

Optional software for Volumetrix MI that expands the review capabilities of GE Hawkeye 4 SPECT/CT and offers flexibility to the users of Discovery and Optima SPECT/CT and PET/CT systems, including:

SECRET DE AFACERI

CONFIDENTIAL

## Volumetrix 3D

Integrates 3D Fusion into nuclear medicine workflow, bringing together the benefits of advanced 3D visualization with productivity of traditional 2D image analysis.

- 3D VR of fused SPECT-CT, PET-CT or MR images
- NM segmentation tools to include or exclude portions of 3D rendered images, including removal of the table from the CT image and threshold based semi-automatic segmentation
- Clip & Cut Planes to integrate traditional Axial, Sagittal and Coronal slices simultaneously in to the 3D rendered objects

## Volumetrix IR

Offers the clinician a choice of virtually any CT to be registered to NM or PET data in the NM workflow:

- Automatic rigid registration of CT to other CT or MR
- Manual adjustment via shift/rotate
- Full or Regional registration (fusion & collect supported)
- Auto accept or user intervention are allowed

## Cardiology

### Myovation

- Provides side-by-side reconstruction and auto reformat of cardiac SPECT, gated SPECT, and PET data including Sestamibi, Thallium, Tetrofosmin, Dual Isotope, FDG, and Rb-82
- Side-by-side review of up to six datasets including beating heart slices
- Comparative review of primary and secondary reconstructions
- Supports reformat of iterative reconstructed data
- Integrated Hawkeye, TAC and ACuscan processing
- Calculation of TID and lung heart ratio
- Polar Plots and Reversibility Polar Plots for Technetium
- Volume based ejection fraction for Tc99m MPI
- 3D Transaxial Post-Filter for both SPECT and PET
- Masking of transaxial and oblique slices
- Enhanced Quality control tools include:
  - Automatic and/or manual QC and Review screens snapshots collection for quick and easy review
  - Gated QC information for Stress and Rest studies including average HR, accepted\Rejected beats and R-R distribution
- Filtered back projection, iterative reconstruction w/scatter correction and/or attenuation correction (where relevant data is provided)

- Three methods of SPECT motion correction are available:
  - Stasis, Hopkins -Cardiac Oriented and MDC-SPECT Motion Detection & Correction
- Well-suited cardiac review screen, including 3D, Polar Maps, four sets, five slices review, ED/ES review, beating slices, EF, side-by-side perfusion and beating gated review
- Synchronized gated slice beating
- Direct linkage with optional 3rd party packages
- Automatic gated and non-gated SPECT motion correction
- Clip normalization

### Alcyone Dynamic Analysis Tool (option)

The Dynamic Analysis Tool option assists with flow measurement of dynamic SPECT and planar images and allows display, fit, and export time activity curves (TAC) based on two dimensional region of interest (ROI) of reconstructed data.

The tool includes:

- ROI Drawing tools for Dynamic Images
- Interactive count calculation and display
- Online gamma fit calculation and display
- Data Export for post-processing calculations

Applicable for Discovery NM 530c and Discovery NM/CT 570c data only.

*Note: Results derived from this application are for technology exploratory use only and should not be used for patient diagnosis.*

### Alcyone Reconstruction Tools (option)

The Alcyone Reconstruction tools option includes five applications that improve productivity for image reconstruction:

- Batch Recon – process one or multiple studies using a variety of reconstruction and filter parameters to optimize image quality. Are able to set a different set of reconstruction parameters for Stress and Rest, Gated and Summed, and for attenuation corrected and non-corrected images
- Filt3D Batch – performs batch filtering of multiple transaxial images from Alcyone data
- Make SA – reformat transaxial Alcyone slice data to short access (SA) oblique slices according to reformat limits used in SA images that were created in Myovation
- Make Scatter – generates an image based on weighted sum of up to 5 projections of planar Alcyone data

SECRET DE AFACERI

CONFIDENTIAL

- System Matrix Generator - allows users to generate additional system matrices with different voxel and FOV sizes, thus providing additional control over reconstructed images

Applicable for Discovery NM 530c and Discovery NM/CT 570c data only.

*Note: Results derived from this application are for technology exploratory use only and should not be used for patient diagnosis.*

### Alcyone Lister Tool (option)

The Alcyone Lister Tool is an optional application for reframing data acquired in list mode on Alcyone technology cameras, using scan parameters that are different than those used to generate projection data during the acquisition.

The Lister tool has two modes of operation:

- Interactive mode: to modify reframing parameters, display the distributions of the basic reframing parameters, and view results.

There are two main usages:

- Interactive Processing of Specific Dataset: reframe a specific patient study (with embedded list mode data) for research purposes.
- Setting Up a Clone for Batch Processing: create an automatic Lister clone; useful for processing multiple sets of list mode data.
- Automatic Batch Processing mode: used for batch processing of multiple datasets with identical Lister settings, using automatic Lister clones

With Alcyone Lister tool, it is possible to modify the following reframing parameters: Time, Energy, and Number of Gates

Applicable for Discovery NM 530c and Discovery NM/CT 570c data only.

*Note: Results derived from this application are for technology exploratory use only and should not be used for patient diagnosis.*

### Alcyone Motion Management (option)

Discovery NM 530c and Discovery NM/CT 570c cameras incorporate a novel

multi-pinhole geometry which is not compatible with traditional motion correction techniques.

- Alcyone Motion Management (Alcyone MDC) application allows for motion detection and correction of studies acquired on Alcyone gamma cameras - Discovery NM 530c & Discovery NM/CT 570c
- Alcyone MDC application is applicable in low dose and high dose 99mTc and 201Tl scans

- Alcyone MDC is designed to enhance visual clarity of Alcyone cardiac studies.

- Simultaneous all view acquisition in Alcyone cameras enables high temporal resolution which may allow for detection and correction of respiratory motion by the Alcyone MDC application when adequate count rate is available.

- Alcyone MDC application does not require any H/W accessory devices for the detection and correction of patient and/or respiratory motion

*Note: MDC for Alcyone can generate only non-attenuation corrected SPECT images.*

### Alcyone to Multigated Views (option)

The Alcyone to Multigated Views tool performs parallel forward projection to the gated transaxial images to generate four standard multi-gated images that correspond to given/selected projection directions: anterior (ANT), left anterior oblique (LAO), left lateral oblique (LLO), left posterior oblique (LPO) viewing directions. These dynamic views can then be loaded to EF Analysis for processing and review.

### Alcyone External CTAC (option)

Alcyone External CTAC is an optional application where external GE CT scanner data can be used for the Attenuation Correction (AC) of the Emission data. The attenuation maps are built from the nominal CT slices by the Attenuation Correction QC (ACQC) application that is integrated in the Myovation application.

- Attenuation Correction is implemented as a part of the MLEM reconstruction of the NM planar data.
- The data is transferred/imported into Xeleris 4.0 via DICOM import or transfer

### Single CTAC (option)

Optional feature for Myovation that enables the use of single CTAC scan of same patient acquired by Hybrid GE camera or single stand alone CT for accurate attenuation correction of multiple SPECT scans, assuming patient position (including all movable organs) and body size are similar, yielding:

1. Enhanced image quality, compared to non-corrected image
2. Eliminating the need for acquisition of separate CT scans for AC of each SPECT study leading to simplification and potentially higher productivity
3. Potentially improved patient care (due to potentially reduced exposure to CT radiation) by eliminating the need for acquisition of separate CT scans for AC of each SPECT study

*Note: feature was validated using data from GE Healthcare 600's series, Infinia, Infinia Hawkeye 4 and Ventri.*

SECRET DE AFACERI  
CONFIDENTIAL

## AdreView™ Planar (Option)

AdreView Planar is an optional Xeleris application for analysis of cardiac (chest) 123I mIBG (Iobenguane 123I Injection) images that provides a guided workflow and interactive tools to generate accurate and repeatable quantitative Heart to Mediastinum ratio reports (as compared to manual processing).

### AdreView Planar provides:

- Semi-automated placement and sizing of standardized Regions-Of-Interest (ROIs) on 123I mIBG chest planar images,
- Automatic computation of Heart-to-Mediastinum (H/M) ratio,
- Saving of results series, including image, ROIs and H/M ratio.
- AdreView Planar H/M application is enabled for use in all countries where Iobenguane 123I Injection radiopharmaceutical is allowed to use in cardiac gamma imaging.

## First Pass and EF Analysis

- EF analysis with volume curve
- Left/right ventricle selection
- Auto/manual ROIs for systole and diastole
- Interactive beat selection
- Phase/Amplitude analysis
- Quality control review screen including transit times
- Supports list mode data from DS-Series cameras

## Peak Filling Rate

- The Peak Filling Rate protocol uses the EF Analysis results as input, which includes the left ventricle (LV) emptying index and determines several timing, ejection fraction and rate parameters, which provide additional information of the heart.

## L-R Shunt

- Calculates the pulmonary to system flow ratio. Used to determine the existence and size of the inter-cardiac shunts
- Supports list mode data from DS-Series cameras

## Cardiac Morphing (option)

- Delivers enhanced perfusion image quality for Gated SPECT data
- Overcomes non-uniform blurring of the myocardium due to cardiac motion, enhancing the visual clarity of the images
- Based on elastic summation of tomograms representing all bins in a cardiac cycle – rescaled to match end-diastolic LV size

- Cardiac Morphing shall be performed on gated Rest- Stress data acquired at the standard full time or full dose
- Processing of dual gated Rest-Stress study is likely to complete in approximately 5 minutes if processed using 64x64 image matrix or 20 minutes if processed using 128x128 image matrix.

## Evolution for Cardiac (option)<sup>1</sup>

An OSEM resolution recovery reconstruction algorithm which models the collimator-detector response developed at John Hopkins University & UNC Chapel Hill. It improves cardiac SPECT resolution, signal to noise ratios, and noise variability. The package provides non-inferior image quality with up to 50% reduction in count density of Cardiac SPECT studies, enabling:

1. Imaging at ½ acquisition time (6 – 7.5 min) compared to standard MPI protocols, or
2. Imaging at ½ the dose compared to standard MPI protocols, or
3. Various combinations of 1 and 2.

## MDC - Motion Detection and Correction

Refer to “General SPECT and Hybrid SPECT/CT”

## XT-ERNA (option)

An application used to analyze cardiac gated data to assess cardiac functionality. The application uses edge detection algorithms to generate the Left and Right Ventricle End Diastole and End Systole ROIs, from which it calculates Ejection Fraction, Peak-filling rate, Peak-emptying rate and other essential parameters.

## Emory Cardiac Toolbox™ (option)

A comprehensive set of nuclear cardiology protocols for advanced cardiac analysis, including a variety of databases of normal patients for comparison covering various acquisition protocols, including 1 and 2 day Sestamibi, Dual Isotope, Stress/Rest Tetrofosmin, Thallium, Myoview Pharmacologic Stress, and a 1 day Rest/Stress Tetrofosmin for GE Healthcare Discovery NM 530c and NM/CT 570c cameras.

Options available with the Emory Cardiac Toolbox software:

### PERFEX™ (option)

Optional module that provides an automatically generated text interpretation of Cequal™ perfusion results.

**SECRET DE AFACERI**  
**CONFIDENTIAL**

#### **Emory SyncTool™ (option)**

Optional module, which enables assessment of LV asynchrony by phase analysis of Gated Single Photon Emission Computed Tomography (GSPECT) Myocardial Perfusion Imaging (MPI) studies. The LV asynchrony quantification is a feasible predictor of response to CRT (Cardiac Resynchronization Therapy) and its applicability in clinical settings.

#### **AdreView™ Tools (option)**

Optional module which enables assessment of SPECT images acquired using AdreView (123I mIBG), a radiopharmaceutical for imaging myocardial innervation. The following tools are provided:

- Determination of uptake of mIBG in the heart relative to uptake of some SPECT perfusion agent such as tetrofosmin, sestamibi or thallium.
- Calculation of mIBG heart-to-mediastinum ratio from transaxial images

#### **PET database for ECToolbox™ (option)**

- Optional module that provides Rb-82 and NH3 databases, RB/ FDG and SPECT/FDG mismatch tools.

#### **Cedars Sinai Cardiac Packages (option)**

A comprehensive set of nuclear cardiology protocols for advanced cardiac analysis, including:

##### **Cedars Sinai Quantitative Perfusion SPECT™ (option)**

Automatic 3-Dimensional software approach to quantitative Perfusion SPECT.

##### **Cedars Sinai Quantitative Gated SPECT™ (option)**

An application calculating the ejection fraction of the left ventricle and a 3D surface display is generated.

##### **Cedars Sinai Quantitative Blood Pool GSPECT™ (option)**

Application for the quantitative analysis of gated cardiac blood pool datasets computing volumes and ejection fractions. The application supports Alcyone data as well as conventional angio cameras.

##### **Cedars Sinai Automatic Report Generator™ (ARG) (option)**

Optional module available for the QGS/QPS applications for generating reports.

##### **Cedars Sinai Companion™ (option)**

Optional module for QGS and QPS applications features

- 17 segment scores and templates in QPS
- Diastolic filling parameters in QGS
- Eccentricity ratio in QGS

#### **Pluspack lite Package (option):**

- Prone/Supine
- Stress/Rest registration and serial change
- Shape Index

#### **Invia Corridor 4DM™ (option)**

A comprehensive set of nuclear cardiology protocols for advanced cardiac analysis developed by Invia and the University of Michigan Medical Center in Ann Arbor, Michigan including 4DM-SPECT, 4DM-PET and 4DM-CT options.

- Automated Workflows
- Multi-Monitor support
- Fusion SUV
- Calcium Scoring Database
- Audit Logging
- HL7 Export

#### **4DM-SPECT (option)**

A comprehensive cardiac SPECT display and quantification program for gated and ungated SPECT perfusion studies which includes its own report generation package.

#### **4DM-CT Fusion (option)**

Optional module available for the 4DM SPECT application allows CT and SPECT fusion viewing, CT reporting options and coronary calcium analyses. Coronary arteries display is not available.

#### **4DM-PET (option)**

A comprehensive cardiac display and quantification program for gated and ungated PET perfusion studies which includes its own report generation package.

#### **Alcyone CFR (option)**

The optional Alcyone CFR application capitalizes the unique Dynamic SPECT capabilities of the Alcyone design. It allows for the quantification of coronary vascular function by deriving Myocardial Blood Flow (MBF) and then calculating Coronary Flow Reserve (CFR) indices on data acquired on PET scanners and on stationary SPECT scanners with the capacity for dynamic SPECT imaging. These indices may add information to physicians using Myocardial Perfusion Imaging for the diagnosis of Coronary Artery Disease (CAD).

**SECRET DE AFACERI**  
**CONFIDENTIAL**

Quantitation and Display of Alcyone CFR is powered by INVIA's Corridor 4DM application.

Alcyone CFR Application capabilities include:

- SPECT (Alcyone Technology) and PET support
- Automated CFR quantification
- Global and regional absolute blood flow results
- QA tools to track LV uptake and territory distributions
- Configurable Time Activity Curve (TAC) displays
- Interactive Polar Maps

## Nephrology

### Renal analysis

Tools for renal perfusion and function analysis and review which automatically identify the kidney ROIs, yet still allow adjustment via threshold, manual draw, and editing tools.

- Generate Renogram for diuretic and captopril acquisitions
- DTPA, MAG3 and LASIX renography
- DMSA Renal
- Pediatric kidney depth calculation
- Analysis methods include Gates GFR, QuantEM (Option),
- Modified Gates, Modified Schlegel
- Single Sample Clearance include Dubovsky (for Ortho Iodo hippuran and Tc99mMAG3) and Bubeck (Tc99m MAG3)
- Perfusion methods include Hilson, Peter, and Kirchner
- Relative uptake methods include Slope and Integral methods
- Excretion index (20 min, 30 min, residual)
- Manchester with Rutland slope for relative uptake support

### Renogram DMSA

- Relative and absolute function (%dose)
- Three methods of entry of dose (camera, well counter, and dose calibrator)
- Five methods of kidney depth estimation (manual, Raynaud, Taylor, Tonnensen, and conjugate view)

## Bone

### Volumetrix MI Evolution for Bone (option)<sup>1,2</sup>

Optional iterative reconstruction with resolution recovery well-suited for Bone SPECT studies. The Evolution for Bone algorithm, which was developed at John Hopkins University & UNC Chapel Hill, models the collimator-detector response, improves Bone SPECT resolution, signal to noise ratios and reduces noise variability

Use of the option enables:

1. Improved resolution of Bone SPECT studies acquired over standard acquisition time; or
2. Non-inferior image quality with up to 50% Reduction in count density, enabling:
  - a. Imaging at ½ acquisition time (6 – 7.5 min) compared to standard Bone SPECT protocols; or
  - b. Imaging at ½ the dose compared to standard Bone SPECT imaging protocols; or
  - c. Various combinations of (a) and (b).

The Evolution for Bone reconstruction is an additional module for the Volumetrix MI application.

### Whole Body and Bone Spots Review

- Review whole body and spot images
- Analysis of the perfusion and blood pool phases of three-phase or four-phase bone studies
- Dual zoom mode for spot generation from whole body images
- Predefined automatic layouts of whole body images/spots
- Regional masking
- ROI comparison based analysis
- Statistic and count display
- Sacroiliac ratio with manual or auto edge detection

### Evolution for Planar Bone – ASM/NL (option)<sup>1,2</sup>

An optional add-on feature to Whole Body and Spots Bone Review application. Evolution for Planar Bone includes a noise-reduction algorithm that preserves the fine structures in the image. Evolution for Planar Bone uses well-suited pixel size and optimal energy window settings. (Adaptive Structure Matching Non-Local Filter).

SECRET DE AFACERI  
CONFIDENTIAL

It enables non-inferior image quality with up to 50% reduction in count density compared to standard protocol, delivering any of the following:

- Improved WB planar image quality for the same scan time
- Shorter WB planar scan time while preserving image quality
- Reduced injected dose of the same scan time while preserving image quality.

## Oncology & Infection

### Evolution Toolkit (option)<sup>2,3</sup>

A package enabling improved resolution and reduced noise for SPECT studies of <sup>99m</sup>Tc, <sup>123</sup>I, <sup>111</sup>In, <sup>131</sup>I, <sup>67</sup>Ga by the use of the Evolution reconstruction technique with resolution-recovery.

Compared to standard FBP or iterative reconstruction, Evolution Toolkit can enable improved visual clarity. Included in the Volumetrix MI for Oncology Evolution and Statistical Resampling packages:

Volumetrix MI for Oncology Evolution is well-suited for the processing and review of SPECT and SPECT-CT data of <sup>67</sup>Ga, <sup>111</sup>In, <sup>131</sup>I and <sup>123</sup>I using iterative reconstruction technique with resolution-recovery.

Including Poisson and Angular re-sampling tools, which enable imaging simulation of various levels of count densities to test the impact of time or dose reduction on image quality.

### Q.Metrix (option)

A package enabling quantitative SPECT results in the form of MBq/ml and SPECT SUV (Standard Uptake Value). The

Q.Metrix application utilizes advanced Evolution reconstruction with compensation for Attenuation, Resolution and Scatter. Patient demographics and dose information are incorporated to provide accurate quantitative results.

Quantitative SPECT results are further enhanced with advance segmentation tools providing 2D and 3D organ and lesion characterization. The application provides for patient follow-up and longitudinal comparison. Q.Metrix supports data from GE Healthcare 600's Hybrid products using the following isotopes and collimators:

- <sup>99m</sup>Tc, <sup>201</sup>Tl, <sup>111</sup>In, <sup>123</sup>I, <sup>131</sup>I, and <sup>67</sup>Ga
- LEHR, MEGP, HEGP, ELEGP, LEHS

### Dosimetry Toolkit (option)

Dosimetry Toolkit uses consecutive patient WB SPECT/CT and/or WB planar datasets for quantifying changes in radiopharmaceutical uptake over time and calculating residence time per organ for Radio-Isotope Therapy (RIT) treatment planning purposes.

The application results can be used as input for RIT planning applications such as OLINDA or similar.

The Dosimetry toolkit includes tools that can help the user to perform in a quick and convenient way the following tasks:

- Reconstruction of all SPECT/CT data, including accurate SPECT/CT registration and quality control, along with patient motion detection and correction, attenuation, scatter and collimator blurring corrections. Image reconstruction is performed using iterative algorithm (OSEM).
- Registrations of all scans to one common reference with
- semi-automatic/manual tools that enable the user to perform local organ registration
- Segmentation of the different organs with semi-automatic/manual tools to differentiate between overlapping organs
- Create time activity curves for each of the organs
- Curve fitting
- Calculate imaging agent residence time in each organ
- Export results and all numerical values in format suitable for input to OLINDA or similar RIT planning applications

*Note: Dosimetry Toolkit is only available with GE Healthcare NM 600's Hybrid products data*

## Gastroenterology

### Gall Bladder EF

Determines ejection fraction (EF) from a dynamic dataset. The EF is computed from the gallbladder and background ROIs on individual images, and can help reduce possible errors due to patient motion during the study.

### Gastric Emptying

Determines the Gastric Emptying Curve as well as other parameters from up to 12 static images or from dynamic images in both posterior and anterior positions.

SECRET DE AFACERI

CONFIDENTIAL

## Esophageal Motility Analysis

- Provides analysis of esophageal dynamic studies via transit curves and compressed images
- Check the pass of a bolus along the esophagus in order to detect esophageal dysmotility and obstruction
- Help visualize refl of the stomach contents into the esophagus

## Pulmonary System

### Lung analysis

- Analysis of lung ventilation and perfusion spots lung studies
- Ventilation index and quantitative perfusion analysis including ratio
- Includes templates for lung segments visualization

### Q.Lung (option)

The optional Q.Lung application is used for processing and reading lung volumetric data, including NM SPECT and hybrid SPECT-CT, external CT (CT from a separate non-hybrid scan). The Q. Lung application provides processing, quantification, and multidimensional review for pulmonary scintigraphy for display and quantification of global and regional ventilation (V) and perfusion (P) on SPECT and SPECT/CT studies.

The combination of ventilation and perfusion studies enables the user to determine whether defects in blood flow are associated with ventilation defects.

Q.Lung application may aid physicians in:

- Diagnosis of Pulmonary Embolism (PE), Chronic Obstructive Pulmonary Disease (COPD), Emphysema and other lung deficiencies.  
Assess the fraction of total lung function provided by a lobe or whole lung for Lung cancer resection requiring removal of an entire lobe, bilobectomy or pneumonectomy.

Q.Lung capabilities include:

- Registration of SPECT Ventilation, SPECT Perfusion and CT studies
- Image Segmentation of lung structures based on CT scans, namely the lungs, airways and lobes.
- Fusion and display of SPECT Ventilation / SPECT Perfusion with CT / CTPA studies
- Display of Ventilation, Perfusion and V/P ratio (Quotient) images
- Computation and display of lobe statistics, including volumes and tracer uptake in Ventilation and Perfusion scans

*Note: Q.Lung is only available with GE Healthcare 600's series hybrid scanners data.*

## Neurology

### Brain SPECT

- Brain SPECT protocol includes automatic reconstruction and reformatting, 3D volume rendering, filter selection and filter gallery
- Motion correction algorithm, MDC - Motion Detection and Correction.
- Automated well-suited parameter selections including attenuation correction, filter type/values and reformat configuration for a variety of organs. Activated based on relevant clinical category: renal, bone, liver, oncology (gallium, indium), lungs and general purpose
- Comprehensive review of all slice planes, 3-plane spot collection and triangulation to show image interrelationships
- QC screen with sinogram and linogram analysis
- 3D masking definable at projection view as well as on reconstructed slices
- Filter selection with filter gallery review to determine the most effective filter and filter parameters. 2D Hanning, Metz, Butterworth, Hamming, and Wiener filters are available for FBP reconstruction. 3D Butterworth, Hanning and Gauss filters are available for OSEM reconstruction
- CBF Segmentation Protocol, including:
  - Trans-axial slice subdivided into equiangular sectors over 360 degrees
  - Brain oblique reformatting
  - OM line definition with external marker
  - Analysis of mean count, total count, maximum pixel value and standard deviation
  - Early/Late and Right/Left Mean ratios
  - Quality control of reconstructed data
  - Image slice comparisons
- Chang Attenuation Correction

An approximate Attenuation correction based on simple calculation of the average attenuation for photons traveling from each point in the organ at different angles.

Applicable to Brain SPECT scans under assumption of uniform attenuation across the entire brain tissue.

SECRET DE AFACERI

CONFIDENTIAL

## Q.Brain (Option)\*

The optional Q. Brain application is intended to aid clinicians in the assessment and quantification of pathologies derived from SPECT and SPECT-CT, PET and PET-CT brain scans.

Q. Brain application features automated analysis through quantification of tracer uptake and comparison with the corresponding tracer uptake in control subjects. The resulting quantification is presented using volume of interests, voxel-based and 3D stereotactic surface projection maps of the brain.

Q. Brain image analysis standardizes individual brain shapes into a standard atlas shape while preserving the functional information measured by SPECT and PET imaging.

SPECT/PET co-registration to MR and fusion display capabilities allows functional findings to be related to anatomy, and offers visualization of structural abnormalities.

The package allows the user to visualize and quantify relative changes in the brain's metabolic function or blood flow activity between a subject's images and controls, when used with radiopharmaceuticals approved by the regulatory authority in the country of use, which may be resulting from brain function alterations in:

- Epileptic seizures
- Dementia. Such as Alzheimer's disease, Lewy body dementia, Parkinson's disease with dementia, vascular dementia, and frontotemporal dementia.
- Inflammation
- Brain death
- Cerebrovascular disease such as Acute stroke, Chronic and acute ischaemia
- Traumatic Brain Injury (TBI)

When integrated with the patient's clinical and diagnostic information, Q.Brain application may aid the physician in the interpretation of cognitive complaints, neuro-degenerative disease processes and brain injuries.

Q.Brain Application provides a full suite of processing, review and reporting capabilities.

Processing capabilities include:

- Motion Detection and Correction (SPECT)
- FBP and OSEM Reconstruction (SPECT) including
  - CT-based or Chang attenuation and scatter corrections
  - OSEM reconstruction with Resolution Recovery (licensed option)
- Manual Re-orientation to AC/PC line.
- Hybrid Quality Control (SPECT/CT)

- Co-Registration (Licensed Option): Automatic and/or user-adjusted co-registration for
  - SPECT to SPECT, PET to PET (e.g. Early-Late studies).
  - SPECT to MR, PET to MR
  - SPECT/CT to MR, PET/CT to MR
- Proprietary, adaptive template registration technique with anatomic normalization to a standardized space.

Review Capabilities include:

- Comparison with normals databases for the following tracers: 18F-FDG, HMPAO, ECD and IMP
- Creation of user defined normal databases
- 3D SSP models for uptake ratio and Z-score images
- PET/SPECT only based quantitation. Not dependent on MR to deliver quantitative analysis.
- Regional quantitative results
- Quantitative comparison and display of longitudinal studies
- Same standardized space is used by the DaTQUANT and Q.Brain applications; this allows for fused display and review of registered Perfusion and 123I-ioflupane images

Reporting Capabilities include:

- Referring physician and patient based exam summary reports for streamlining communication of results
- Customizable interface and easy exporting capability of results

*Note: Q.Brain is only available with GE Healthcare 600's products data and PET products.*

## DaTQUANT™ (option)

DaTQUANT optional application enables visual evaluation and quantification of 123I-ioflupane images. DaTQUANT Normal Database enables quantification relative to normal population databases of 123I-ioflupane images.

This application may assist in detection of loss of functional dopaminergic neuron terminals in the striatum, which is correlated with Parkinson disease.

- DaTQUANT analysis gives objective method for the assessment of the extent and intensity of the striatal signal. Relative comparison of uptake ratios may be used for monitoring the progression of concentration of dopamine transporters in the synapses of striatal dopaminergic neurons. Therefore, DaTQUANT analysis may assist to the visual assessment of loss of presynaptic dopaminergic neurons, related to Parkinsonian syndromes, and may provide adjunct tool for Parkinson disease diagnosis.
- DaTQUANT provides accurate analysis of 123I-ioflupane images using pre-defined 123I-ioflupane VOI template, which is more consistent, objective and repeatable compared to manual ROI analysis.

- DaTQUANT provides tools to analyze binding ratios differences between multiple scans.  
DaTQUANT enables comparison with user chosen suitable database of age matched reference values, and may provide adjunct information to aid in diagnosis of Parkinson disease. DaTQUANT supports creation of user defined databases, which enable site specific population comparisons for potentially improved diagnostic accuracy.
- DaTQUANT PDF report provides an easy and convenient reporting tool intended to share results with referring physicians

*Note: DaTQUANT is available for use only in countries where <sup>123</sup>I-ioflupane pharmaceutical is approved for use.*

## Endocrinology

### Thyroid uptake index

Computes thyroid uptake index from Tc99m and I123 static studies for either the entire thyroid per thyroid lobe, or by region (cold or hot spot)

### Parathyroid imaging analysis

Analysis of parathyroid from dual isotope (Tc99m & Tl201) and dual phase MIBI studies.

## General Processing Tools

### Workspace

General application enabling manual processing and review of patient data, including:

- Customized review templates for screen layout and various display settings
- Creating and saving templates appropriate for various types of patient data. Once saved as Review Templates, these layouts can be applied to patient data of the same type for standardization and comparison purposes
- A full set of display tools is provided to manipulate the screen layouts
- An extensive set of image analysis tools
- Processing results can be saved, creating new datasets within the same series or within a new series
- Macro Record mode for the recording of operations and saving them as user applications

### Dynamic Motion Correction (DMC)

Provides tools to manually and automatically correct motion introduced during dynamic acquisition

## Display and analysis

- Single and multiple display of images, curves and ROIs
- Circular, elliptical, rectangular, irregular, or pixel intensity regions of interest are available
- Predefined display templates for easy image repositioning
- Operator control to display multi-frame data sets in user's choice of format
- Ability to use up to eight separate patient cards at one time for processing and display, dependant upon system resources available and utilized by each card
- Image manipulation (i.e., pan, zoom, flip, rotate) using icons or mouse control
- Image filtering (i.e., spatial convolution (9 or 25 point kernel, Gaussian), Fourier (Hanning, Butterworth, Metz or Winner), temporal convolution)
- 2D region-of-interest statistical analysis (total & average counts, ROI area, maximum & minimum counts, standard deviation)
- Cinematic display of dynamic / multi-frame dataset studies
- Synched cine
- Normalization of images to image set or image to itself
- Threshold and windowing control for current image or all displayed images including 2 methods of gamma control and CT windowing (using center and window)
- Saving of modified images for later display and analysis
- Text annotation: predefined annotation with interactive font size and color options
- Numerous color and monochrome map selections
- Scroll facility for multiple image displays
- Generate JPEG screen capture images and AVI cine for export
- Integrated web browser

## Curves, Graphs, Image processing

- Time activity curve, histogram analysis and statistical analysis from ROIs. Curve parameters include max. and min., time to peak, T1/2, slope, FWHM, centroid
- Storage of curves/graphs for later recall and comparison
- Curve math functions for curve processing:
  - Add, subtract, multiply, divide: curve or constant
  - Normalized curve subtract
  - Curve filter: convolution, Hanning, Butterworth
  - Curve fit: linear, gamma variant, exponential
  - List curve points

- Image processing and math functions including:
  - Image fusion
  - Statistics in the image and within a defined ROI
  - Image filter: spatial and temporal convolution, Fourier
  - Copy, combine
  - Add, subtract, multiply, divide: image or constant
  - Reframe to a composite or multiframe dataset
  - Magnify/minify image
  - Normalize to maximum count value
  - Threshold above, below, or between a value
  - Mask inside/outside ROI
  - Vertical & horizontal activity profiles
  - Interactive control of location and size

### User Commands

The extended functionality allows users to record a sequence of operations (Macro command like) using the intuitive Xeleris 4.0 user interface tools to easily customize applications to their needs.

The added user command functionalities to Xeleris include:

- Arithmetic and statistics calculations for images and ROIs
- Graphs, plots and curves manipulations
- Annotations and time stamps
- All Workspace application operations are now available for a range of 2D applications including image manipulation tools such as masking and reframing.

### Discovery NM/CT 670 CZT Lister Tool (option)

True list mode acquisition and processing to enable wide range of research exploration including:

- Energy window optimization – identify the best energy window for the digital detector which maximize signal to scatter ratio
- Dual Isotope – simultaneous dual isotope acquisitions and processing to harness the power of multi isotope imaging, Tc99m and Tl201, or either of these and I123
- Reconstruction optimization – optimize reconstruction parameters for the various isotopes

Applicable for Discovery NM/CT 670 CZT data only.

## Xeleris 4.0 for Molecular Breast Imaging configuration

Xeleris 4.0 for MBI (Discovery 750b) configuration shares the same basic functionality of the Xeleris 4.0 platform optimized to Molecular Breast Imaging (MBI). The Xeleris 4.0 for MBI (Discovery 750b) supports the following applications to provide simplified processing and review of molecular breast imaging exams.

- MBI Viewer
- Workspace
- Aladdin

The following applications are available as optional enhancements to the Xeleris 4.0 MBI configuration

- Xeleris DICOM Viewer
- Dual Monitor

### MBI Viewer application

Xeleris 4.0 for MBI includes the MBI Viewer application enabling:

1. Display of MBI images in customizable formats
2. Dynamic display of MBI Images
3. Reformat of Dynamic MBI Images
4. Registered top and bottom image display

In addition to the MBI viewer, general application enabling manual processing and review of MBI patient data is supported, including:

1. Customized review templates for screen layout and various display settings
2. Creating and saving templates appropriate for various types of patient data. Once saved as Review Templates, these layouts can be applied to patient data of the same type for standardization and comparison purposes
3. A full set of display tools is provided to manipulate the screen layouts
4. An extensive set of image analysis tools
5. Processing results can be saved, creating new datasets within the same series or within a new series
6. Customized Aladdin programming - Visual Basic® based user protocol development environment including protocol development tools, navigation editor, menu editor and results load/save

**SECRET DE AFACERI**  
**CONFIDENTIAL**

## Xeleris 4.0 Cardiac Configuration

This configuration of Xeleris can simplify workflow and may provide greater economy as only applications needed for specific cardiac workflow are provided, without altering the infrastructure functionalities of the Workstation. The Xeleris 4.0 Cardiac Configuration supports the following applications to provide simplified processing and review of Nuclear

Cardiac exams;

- Myovation
- First Pass
- EF Analysis
- Peak Filling Rate
- L-R Shunt
- MDC Motion Detection and Correction

The following applications are available as optional enhancements to the Xeleris 4.0 Cardiac configuration

- Evolution for Cardiac
- Cardiac Morphing
- Emory Cardiac Toolbox
- Cedars Cardiac
- Invia Corridor 4DM
- Xeleris DICOM Viewer
- Dual Monitor
- Access tools

## Customization

### Applications customization

- Create a user defined list of favorite applications quickly accessible through a single tab
- Allows different customization sets for the same application
- Save defaults from within most applications
- Standardized annotation templates
- User-customizable review templates may be created for each study type
- User Commands may be added to 2D planar applications specific for user needs creating personalized applications.

### Color map customization

- Create user-defined color map
- Allow the definition of the default color maps

### Customized Security tools

Tools to assist privacy and security measures compliance.

- Individual user accounts
- Account strength and policies
- Audit logging
- Virus protection software
- Protected configuration and security settings
- Node-to-node authentications
- Anonymized Export

The Anonymized Export application performs anonymization of selected patient data, in compliance with HIPAA regulations. For example, removing patient ID, patient name and institution name. The purpose of this application is to enable data sharing without providing patient identification (for example, for research purposes.)

### Launch Two

Invoke two Xeleris applications simultaneously:

- Same patient data with two different applications.
- Two instances of the same application with two different studies of the same patient for follow-up purposes.

The applications are invoked as if the user manually invoked them in a sequence.

### Customized Aladdin™ programming

Visual Basic® based user protocol development environment including protocol development tools, navigation editor, menu editor and results load/save.

Aladdin user-created protocols have not been validated by GE. Those who wish to create and use such protocols on Xeleris should verify their accuracy (including both display and quantitative results) before doing so. Key features include:

- Creation of protocols by recording actions on the processing card (record mode)
- Protocol text editor and syntax checker
- Protocol debugger including: syntax check, breakpoints, single step, parameter watch windows
- Full set of nuclear function calls added to a Visual Basic® like foundation
- Image and curve processing using standard Xeleris tools ROI, display, database access functions
- Compatible with most eTEGRA / GENIE P&R user Aladdin software
- Protocols have same philosophy as rest of Xeleris user interface providing the ability to incorporate all standard Xeleris tools, including reconstruction and reformatting
- On-line help text with examples for most functions

SECRET DE AFACERI  
CONFIDENTIAL

- Interactive development of programs: running, editing, testing, debugging
- Custom definition of user interface: menus, tabs and buttons
- Review template editor (RTE) to define display screen layouts
- Shortcut toolbar
- Split screen: program code and display screen
- Interface to Xeleris patient database
- Perform customized processing of data files (compression of avi files, reframing data to different frame size, etc.)
- Aladdin advanced programming training available

## Flexibility

- Integrated user interface for data archive/retrieval, image data deletion, filming/hard copy and data transfer
- Simultaneous operation of display, analysis, archive, network filming and reconstruction
- Flexible user configuration and customization
- Print templates allow the user to setup commonly used printers or file formats and print header layouts
- Multiple print headers can be defined, controlling the location and content of textual information including hospital name, hospital logo, and various data attributes
- Concurrent access to data allows the same patient study to be accessed from multiple Xeleris units simultaneously
- Merge series from different studies into a single one in case of split studies
- Easily correct errors entered when typing patient name, study name, and series name during acquisition of Nuclear Medicine, CT, MR, and PET data (not including results series)
- Easily change NM datasets names (excluding results series)

## Operate Xeleris with Native Language

Xeleris user interface (UI) is available translated to native languages other than English for improved customer convenience.

# Integration and Access

## Connectivity & networking

### Xeleris DICOM Viewer (option)

The Xeleris DICOM Viewer enhances access to molecular imaging data by providing a DICOM viewing application for general and hybrid Nuclear Medicine studies in a common Xeleris user interface. The application may be distributed on CD, DVD or network destination via Xeleris interface. Operates on Windows XP, Windows Vista, Windows 7, Windows 8

### Essential Viewer including

- Planar, Dynamic, DICOM SC, DICOM MFSC support
- Basic Workspace review tools including
  - Scroll, Pan and Zoom
  - Color map and window level
  - Measurement

### Volumetrix Viewer including

- 3D MIP display and triangulation to 2D images
- Synchronized triangulation between SPECT-CT and PET-CT images for optimum localization
- SPECT-CT and PET-CT fusion with blending tools
- Predefined layouts for optimum SPECT-CT and PET-CT display
- Oblique reformatting

## Xeleris Floating License (XFL) (option)

### XFL Remote - WAN (option)

Xeleris allows the user to load Xeleris client software across the network (over the LAN and WAN) in a Server-client configuration where the client can be any Windows-based PC that meets minimum hardware requirements.\*\* Included with XFL remote option is XFL ToGo. XFL ToGo provides the ability to work while disconnected from Local or Wide area networks for a period of 72 hours potentially improving workflow and access to clinical images. XFL enables the use of most Xeleris applications on the client remote PC (based on floating license availability on the Xeleris 4.0), see Supported Application Table below.

#### Notes:

1. Some data may require more than 2GB RAM to be used with XFL. For PET-CT and SPECT-CT 2GB RAM is the minimum and 4GB RAM is recommended

2. Performance may vary based on PC and network characteristics, data loaded, and other SW running concurrently

3. The XFL Software may not function properly in cases of conflicts with other software installed on the PC

4. Studies should be read on monitors that comply with the SMPTE and BWH standards

5. XFL software can only be installed on Windows XP SP3, Windows 7 SP1, Windows 8, and Windows 10.

\*\*Intel Based CPU computer, 2GB RAM, Intel multi-Core or Xeon® hyper threaded CPU (equivalent or better), minimum 500MB free disk space for XFL LAN client connected via at least 100BaseT LAN to a Xeleris 4.0 PC and 10GB free disk space for XFL WAN client, IE6.0.

OpenGL version 2.0 or higher is required for Volumetrix 3D. Acrobat Reader® is required for viewing user guide documentation. Any candidate client must be tested with XFLCheckSystem tool for compatibility of other restrictions, such as font size and screen resolution.

### **Xeleris Suite Package for PACS (option)**

Integrating the diagnostic power of Xeleris Nuclear Medicine clinical applications with the PACS functionality of Centricity RA1000 workstation, or any PACS or similar 3rd party Viewing system (Image Review) supporting Open Interface Plug-in architecture.

- Provides a single Functional imaging and radiology workstation environment by running a single departmental-shared worklist for Radiologists and Nuclear Medicine practitioners
- Enables seamless launch of Xeleris applications by up to 5 concurrent users simply by opening the patient study, while making other relevant radiology exams immediately available for inclusive multi-modality patient review
- Eliminates the need for multiple workstations
- Xeleris Suite clients must meet XFL clients minimum specifications
- Xeleris Suite display capabilities are extended to support a variety of screen resolutions and configurations typically used by PACS systems.

Validated screen resolutions:

- 1280 x 1024 (19")
- 1680 x 1050 (22")
- 1920 x 1080 (23")
- 1200 x 1600 (20.1")

### **NM and PET DICOM 3.0 Support**

- DICOM 3.0 support for Nuclear Medicine and PET image and secondary capture image storage and query/retrieve as an SCU and SCP  
Provides the ability to query/retrieve or receive DICOM 3.0 Nuclear Medicine Image data or secondary capture image data from GE Nuclear Medicine and PET-CT scanners
- Xeleris clinical applications are validated with data from GE NM, Discovery PET/CT 690, Discovery PET/CT 600, Discovery LS and Discovery ST scanners

### **CT/MR DICOM support**

Provides support for DICOM 3.0 CT and MR Image data storage as a SCP and SCU.

- Receive and send DICOM 3.0 CT/MR data from compatible CT, MR or PACS systems
- CT or MR data converted into Xeleris database format
- Send DICOM 3.0 Nuclear Medicine Image data or secondary capture image data to PACS or other valid DICOM compatible systems.
- Refer to Xeleris conformance statement

### **IHE compatibility**

Implements the components necessary for compliance with

the Scheduled Workflow profile, IHE nuclear medicine image profile, Portable Data for Imaging profile and Radiation Exposure Monitoring profile, including:

- Scheduled workflow as evidence creator and image display actors
- DICOM storage commit (as SCU) provides ability to automatically mark as archived data stored on PACS (requires PACS support of Storage Commit as SCP)
- Compliance with the nuclear medicine image profile as evidence creator including the result screen export option and image display including the MPR option, general nuclear medicine, and Cardiac nuclear medicine options
- DICOM Multiframe secondary capture (BW and Color)
- Use of DICOM acquisition context tags to identify cardiac data sets
- Synched Cine support
- Compliance with Portable Data for Imaging profile as PDI Exporter and Importer actors
- Compliance with Radiation Exposure Monitoring profile as dose consumer actor, able to receive and send dose reports.

### **DICOM part 10, import and export**

- Exchange of DICOM data via various electronic formats (CD, FTP, etc.)
- Run Xeleris applications on DICOM P10 files by enabling:
  1. Direct data access from the archived media without the need to import it to the local database first
  2. Sharing multiple data sources, providing virtually unlimited archiving space for processing and review

Most Xeleris applications can operate directly on DICOM P10 files, see Supported Application Table list of applications for XFL-Remote (WAN).

### **Remote access for iLinq™ support, InSite™ and data transfer (option)**

GE's remote serviceability features supported via Broadband connection.

### **TVA – TiP™ Virtual Assist (option)**

A secure and flexible way to receive training and applications support via iLinq Broadband connections. Allows the Xeleris user to share the display and control of their workstation with a remote GE Applications Specialist.

SECRET DE AFACERI

CONFIDENTIAL

### Legacy system STP connectivity protocol

Sopha Transfer Protocol release 3.0 allows either the DS Camera with XT 3.0 acquisition software or NXT systems to send data to the Xeleris. The Xeleris can also query and fetch data from either the DS or NXT systems.

### Remote file transfer

Provides a simple user interface to import or export Interfile and/or DICOM Part 10 files via the local network. Modem is no longer supported.

### ICANL/ACR compliance

Supports special layouts to meet ICANL/ACR requirements for accreditation (excluding patient orientation in Planar and Dynamic scans).

### Image file translation with network connectivity from Pegasys workstation (option)

Software option, which allows an existing Pegasys workstation to transfer images to Xeleris for processing and review.

- One way file transfer from Pegasys to Xeleris
- Unprocessed Nuclear Medicine studies only
- Further restrictions may apply to the types of data and hardware compatible with this option

### Film/hardcopy

- Queued digital filming and digital hard copy
- Support for DICOM 3.0 basic grayscale and color print management as an SCU (see conformance statement)
- Standard postscript level 2 printer support via TCP/IP Ethernet interface
- DASM LCAM printers are no longer supported

### Archiving

Local DB size is 100GB of clinical data available online for processing & review

### Re-writable combo CD-ROM drive with DVD reader

Allows user to archive nuclear medicine, PET, CT and MR on CD in DICOM Part 10 format

## Minimum Hardware Requirements

The Xeleris 4.0 workstation is available as "Xeleris", "Xeleris MBI" or "Xeleris Cardiac" and consists of:

- High-performance Intel based computer with high resolution color monitor, keyboard, mouse and CD-RW/ DVD-RW multi drive
- Windows™ 7 operating system
- 4GB RAM memory size
- 500 GB internal disk capacity for on-line image/data storage
- Ethernet network connection (100/1000 Base T)
- 19" LCD color monitor or 22" or 23" wide-screen LCD color monitor
- Keyboard and operator manual are available in various languages

## Options

### DVD-RAM Archiving (option)

Archiving of patient data on DVD-RAM without change in archived data characteristics, functionality or intended use. The DVD archiving is performed through standard Windows tools, increasing data storage versatility

- Read and write support of DVD-RAM format
- 4.0 GB per side; double-sided capacity for DVD-RAM

### Extended database (option)

Extend the Xeleris database with additional 700GB of data space enabling:

- 8 times more than the standard database size
- Logical partition to up to 5 different local databases
- Connectivity to other stations is standard (some limitations may apply with previous Xeleris versions)
- Does not require additional HW

*Note: the extended database does not replace other archive solutions and it is recommended to frequently backup any data.*

**SECRET DE AFACERI**  
**CONFIDENTIAL**

### Dual Monitor (option)

Optional dual monitor configuration available for desktop and XFL configurations, enabling:

- Viewing of two full cards side by side for enhanced compare mode
- Viewing of patient browser on the left screen while running applications on the right
- Viewing more images in dual screen well-suited layouts in Myovation, Volumetrix MI, and Workspace that utilize the 2 monitors
- Optional laser filming and/or color hardcopy

## Siting requirements

- Storage temperature: -40° to 140°F (-40° to 60° C)
- Storage humidity: 8% to 90% RH, non condensing
- Storage altitude: 30,000 ft (0 to 9,144m)
- 23" Wide Screen LCD Display monitor dimensions (with height adjustable stand): (WxHxD) 21.5 x 13.3 - 19.2 x 8.98" (54.3 x 33.7 - 48.7 x 22.8 cm)
- Keyboard Dimensions: (WxDxH) 18.3 x 6.3 x 1.4 inches (46.48 x 16 x 3.56cm)
- Minimum table/desk top space: 28 x 27" (71.0 x 68.6 cm)

Consult your GE representative for complete site planning information.

## Power Requirements

- Input voltage: 100 to 127V VAC and 200 to 240VAC and 118VAC
- Over 50Hz to 60Hz
- Rated input current: 10A@100-127 VAC, 6A@200- 240 VAC, 10A@118 VAC
- Maximum rated power: 475 watt custom power supply - (Wide Ranging, Active PFC)

## Environment Siting Requirements

- Operating temperature: 40° to 95° F (5° to 35° C)
- Operating humidity: 8% to 85% RH, non condensing
- Operating altitude: 10,000 ft (0 to 3,048m)

## Warranty

Please contact your sales representative for warranty information.

<sup>1</sup>In clinical practice, Evolution options (Evolution for Bone, Evolution for Cardiac, Evolution for Bone Planar) and Evolution Toolkit are recommended for use following consultation of a NM physician, physicist and/or application specialist to determine the appropriate dose or scan time reduction to obtain diagnostic image quality for a particular clinical task, depending on the protocol adopted by the clinical site.

<sup>2</sup>Evolution claims are supported by simulation of count statistics using default factory protocols and imaging of 99mTc based radiotracers with LEHR collimator on anthropomorphic phantom or realistic NCAT - SIMSET phantom followed by quantitative and qualitative images comparison.

<sup>3</sup>Evolution Toolkit claims are supported by simulation of full count statistics using lesion simulation phantom images based on various radiotracers and collimators and by showing that SPECT image quality reconstructed with Evolution Toolkit provide equivalent clinical information but have better signal-to-noise, contrast, and lesion resolution compare to the images reconstructed with FBP / OSEM.

<sup>4</sup>In the US, 18F-FDG is not indicated for the evaluation for Alzheimer's disease.

## Xeleris™ 4.0 Configurations

|                                                                                                                                               | Xeleris Full | XFL & XFL Suite | XFL Remote | Xeleris Cardiology | Xeleris MBI |
|-----------------------------------------------------------------------------------------------------------------------------------------------|--------------|-----------------|------------|--------------------|-------------|
| <b>General Tomography</b>                                                                                                                     |              |                 |            |                    |             |
| Single CTAC                                                                                                                                   | X            | X               | X          |                    |             |
| Volumetrix 3D & IR                                                                                                                            | X*           | X*              | X*         |                    |             |
| Workspace                                                                                                                                     | X            | X               | X          | X                  | X           |
| <b>Cardiac</b>                                                                                                                                |              |                 |            |                    |             |
| Myovation                                                                                                                                     | X            | X               | X          | X                  |             |
| Myovation Evolution (Evolution for Cardiac)                                                                                                   | X*           | X**             | X**        | X*                 |             |
| Cardiac Morphing                                                                                                                              | X*           | X**             |            | X*                 |             |
| Emory Cardiac Toolbox™                                                                                                                        | X*           | X*              | X*         | X*                 |             |
| ECTB – PerfEx™                                                                                                                                | X*           |                 |            | X*                 |             |
| ECTB – SyncTool™                                                                                                                              | X*           | X*              | X*         | X*                 |             |
| ECTB – NRP                                                                                                                                    | X*           | X*              | X*         | X*                 |             |
| ECTB – Heart Fusion™                                                                                                                          | X*           | X*              | X*         | X*                 |             |
| ECTB - Adreview™ Tools                                                                                                                        | X*           |                 |            | X*                 |             |
| Corridor-4DM™ packages                                                                                                                        | X*           | X*/*****        |            | X*                 |             |
| Cedars QGS™, QPS™ and Companion™                                                                                                              | X*           | X*              | X*         | X                  |             |
| Blood Pool Gated SPECT                                                                                                                        | X*           |                 |            | X*                 |             |
| Cedars ARG™                                                                                                                                   | X*           |                 |            | X*                 |             |
| EF Analysis                                                                                                                                   | X            | X               |            | X                  |             |
| First Pass                                                                                                                                    | X            | X               |            | X                  |             |
| L-R Shunt                                                                                                                                     | X            | X               |            | X                  |             |
| Peak Filling Rate                                                                                                                             | X            | X               |            | X                  |             |
| XT-ERNA                                                                                                                                       | X*           |                 |            | X*                 |             |
| AdreView™ Planar                                                                                                                              | X*           | X**             | X**        | X*                 |             |
| Alcyone Lister Tool                                                                                                                           | X*           |                 |            | X*                 |             |
| Alcyone Access Tools (includes Alcyone Dynamic Analysis Tool, Alcyone Reconstruction Tool, Alcyone Motion Management and Alcyone External AC) | X*           | X***            | X***       | X*                 |             |
| <b>Renal</b>                                                                                                                                  |              |                 |            |                    |             |
| Renal Analysis                                                                                                                                | X            | X               | X          |                    |             |
| QuantEM                                                                                                                                       | X*           |                 |            |                    |             |
| DMSA                                                                                                                                          | X            | X               | X          |                    |             |
| <b>Bone</b>                                                                                                                                   |              |                 |            |                    |             |
| WB & Spots Bone Review                                                                                                                        | X            | X               | X          |                    |             |
| Evolution for Planar Bone                                                                                                                     | X*           |                 |            |                    |             |
| Volumetrix MI Evolution for Bone                                                                                                              | X*           | X**             | X**        |                    |             |

|                                        | Xeleris Full | XFL & XFL Suite | XFL Remote | Xeleris Cardiology | Xeleris MBI |
|----------------------------------------|--------------|-----------------|------------|--------------------|-------------|
| <b>Oncology &amp; Infection</b>        |              |                 |            |                    |             |
| WB & Spots Oncology Review             | X            | X               | X          |                    |             |
| Evolution Toolkit                      | X*           | X**             |            |                    |             |
| Dosimetry Toolkit                      | X*           |                 |            |                    |             |
| Q.Metrix                               | X*           | X**             |            |                    |             |
| Volumetrix MI Evolution for Oncology   | X*           | X**             | X**        |                    |             |
| MBI Viewer                             | X*           | X**             | X**        |                    | X           |
| <b>Gastrointestinal</b>                |              |                 |            |                    |             |
| Gallbladder EF                         | X            | X               | X          |                    |             |
| Gastric Emptying                       | X            | X               | X          |                    |             |
| Esophageal Motility Analysis           | X            | X               |            |                    |             |
| <b>Pulmonary</b>                       |              |                 |            |                    |             |
| Lung Analysis                          | X            | X               | X          |                    |             |
| Q.Lung                                 | X*           | X*              | X*         |                    |             |
| <b>Brain</b>                           |              |                 |            |                    |             |
| Brain SPECT                            | X            | X               | X          |                    |             |
| Q.Brain                                | X*           | X*/****         | X*/****    |                    |             |
| DaTQUANT™                              | X*           | X**             | X**        |                    |             |
| <b>Miscellaneous</b>                   |              |                 |            |                    |             |
| Motion Detection and Correction (MDC)  | X            | X               | X          | X                  |             |
| Dynamic Motion Correction              | X            | X               | X          |                    |             |
| Thyroid Uptake Index                   | X            | X               | X          |                    |             |
| Parathyroid Imaging                    | X            | X               | X          |                    |             |
| Xeleris DICOM Viewer                   | X*           |                 |            | X*                 | X*          |
| Poisson Resampling (Evolution Toolkit) | X*           | X**             |            |                    |             |
| Angular Resampling (Evolution Toolkit) | X*           | X**             |            |                    |             |
| ROI comparative Analysis               | X            | X               |            |                    |             |
| Launch Two                             | X            | X               |            | X                  |             |
| 3D Filtering                           | X            | X               |            | X                  |             |
| Extract Tool                           | X            | X               |            |                    |             |
| WMV Encoder                            | X            | X               |            |                    |             |
| Color Map Creator                      | X            | X               |            | X                  |             |
| Aladdin Editor                         | X            | X               |            | X                  | X           |
| Flood Uniformity                       | X            |                 |            | X                  |             |
| Hawkeye Registration QC                | X            |                 |            | X                  |             |
| COR Test                               | X            |                 |            | X                  |             |
| SPECT-CT Registration QC               | X            |                 |            | X                  |             |
| Make 8 Gates                           | X            |                 |            |                    |             |

|                         | Xeleris Full | XFL & XFL Suite | XFL Remote | Xeleris Cardiology | Xeleris MBI |
|-------------------------|--------------|-----------------|------------|--------------------|-------------|
| Filt3DBatch             | X            |                 | X          | X                  |             |
| Make Scatter            | X            |                 | X          | X                  |             |
| Anonymized Export       | X            | X               | X          | X                  | X           |
| Interfile Import/Export | X            |                 |            | X                  | X           |
| Review Templates        | X            | X               | X          | X                  |             |
| User Commands           | X            | X               | X          | X                  |             |
| D600 CZT Lister         | X*           |                 |            |                    |             |
| <b>Other</b>            |              |                 |            |                    |             |
| Dual Monitor            | X*           | X**             | X**        | X*                 | X*          |
| Pegasys DirectConnect   | X*           |                 |            |                    |             |

\*Denotes optional software

\*\*Denotes software available on XFL if available on the Xeleris server

\*\*\*Denotes that only some of the Alcyone Access Tools are available on XFL due to Lister NOT being available on XFL. Available: Alcyone Count Metric, Alcyone Multigated Views, Alcyone motion detection.

\*\*\*\*Denotes software available on 64-bit operating systems only.

\*\*\*\*\* Denotes software available on Windows 7 or later operating systems only.

Note: The QuantEm is the software license only. It does not include the syringe holder.

© 2016 General Electric Company-All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue any products described at any time without notice or obligation. Please contact your GE representative for the most current information.

GE, GE Monogram and imagination at work, Xeleris, Discovery, Optima, Brivo, Infinia, Hawkeye, Genie, Entegra, Volumetrix, AdreView, DaTQUANT, Aladdin, iLinq, InSite, and TiP trademarks of General Electric Company.

All third party trademarks are the property of their respective owners.

GE Healthcare, a division of General Electric Company.


GE Healthcare  
3000 N Grandview Blvd  
Waukesha, WI 53188  
USA  
www.gehealthcare.com

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