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#### CYTO-BRUSH

*Not sterile, disposable brush for endocervical cells collection, in cytology and microbiology. It allows the drawing and laying on slides of a bigger quality of cells, preventing to damage their structure.*

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 *Medical Devices - manufactured in compliance with 93/42/EC Directive (Italian D.lgs 46 del 24/02/1997) + 2007/47/EC - Class I not sterile*

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*This device is intended as “MEDICAL DEVICE” Class I not sterile – CYTO-BRUSH is intended for gynaecological use for the collection of cells in obstetrics and gynaecology for subsequent diagnostic investigations. It is used for cell collection in the endocervical tract. **For professional use only.***

**National classification of medical devices (CND - For Italian law) > U089002 (Gynaecological cytology devices)**

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*For handling, before the use, wash your hands with disinfectant soap and use appropriate gloves. There are no other special handling conditions*

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*First make the sampling in the vaginal fornix then gently insert the brush side of CYTO-BRUSH in the endocervix leaving the bristles visible. Rotate one full turn and then remove. Caution: do not use CYTO-BRUSH for endometrial sampling. Be gentle with pregnant patients. For professional use only. Disposable - do not reuse: risk of infection and loss of functionality. Store in a cool and dry place, protect from dampness. Dispose after use according to the regulations in force.*

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*Do not use the units in case of damages or contamination. Do not use after expiry date.*

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


1. Proceed to drawing by delicately introducing the device into the endocervices, for the length of brushes;
2. Turn brushes at 360° by one single turn and then extract the device;
3. Lay drawn sample on a slide by turning brushes and slightly pressing them on it. Make a slight pressure to maintain the integrity of the cells and prevent degenerative processes. The material crawled must be thin because the presence of superimposed layers prevents an adequate reading of the preparation.
4. Fixing: fix the prepared immediately to prevent spontaneous cell degeneration consequent drying: even a short delay may cause damage to the cells.
5. If you are using a fixative spray form, the device must be kept at a distance of 15/20 cm from the slide: in fact, if carried out in the shortest distance, the power delivery of the cells moves to the periphery and induce the formation of air bubbles in the centre of the slide.
6. Allow the slide to dry in a horizontal position.
7. Dispose the CYTO-BRUSH according to current regulations.

**WARNING:** the use of the device in case of pregnancy must be extremely delicate. Do not use for endometrial drawing.

**For professional use only.**