# Protégé™ EverFlex™

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## Instructions for use

### Device description

The Protégé™ EverFlex™ self-expanding stent system is a self-expanding nitinol stent system intended for permanent implantation. The stent is made of a nickel titanium alloy (nitinol) and comes premounted on a 6 Fr (0.035 in / 0.89 mm) over-the-wire (OTW) delivery system. The stent is cut from a nitinol tube in an open lattice design and has tantalum radiopaque markers at the proximal and distal ends of the stent. After deployment, the stent achieves its predetermined diameter and exerts a constant, gentle outward force to establish patency.

The delivery system, as shown in **Figure 1**, is comprised of an inner shaft (1) and outer sheath (2), which are secured with a safety lock (3). The nylon inner shaft terminates distally in a flexible catheter tip (4) and originates proximally at the hub (5). The distal portion of the delivery system for the 20 to 150 mm stents, as shown in **Figure 2**, is comprised of two radiopaque markers, one marker distal (6) and one retainer marker proximal (7) to the constrained stent, on the inner shaft.

The distal portion of the delivery system for the 200 mm stents, as shown in **Figure 3**, is comprised of the same components as those in **Figure 2** except for the following radiopaque markers: one retainer marker distal (13) and one holder marker proximal (14) to the constrained stent, on the inner shaft.

As shown in **Figure 1**, the outer sheath connects proximally to the manifold (8). The self-expanding stent is constrained within the space between the inner shaft and the outer sheath. This space is flushed before the procedure through the stopcock (9). The outer sheath has a radiopaque marker at its distal end (10), as shown in **Figures 2 and 3**.

Stent positioning at the target lesion (peripheral indication) or the stricture (biliary indication) is achieved before deployment by using the two radiopaque markers on the inner shaft, which mark the location of the constrained stent. Before stent deployment, turn the safety lock counterclockwise to unlock the outer sheath.

Retract the outer sheath by pulling the distal grip (11) toward the proximal grip (12). The stent is completely deployed when the radiopaque marker on the outer sheath passes the proximal radiopaque marker on the inner shaft.


Figure 1. Stent on delivery system


Figure 2. Distal portion of 20 to 150 mm delivery system


Figure 3. Distal portion of 200 mm delivery system

## Indications for use

### Peripheral

The stent is indicated for use in the following types of occlusions or lesions:

- Occlusions or lesions at high risk for abrupt closure or threatened closure following percutaneous transluminal angioplasty (PTA)
- Lesions that appear to be at high risk for restenosis following PTA in the common iliac, external iliac, superficial femoral, proximal popliteal, or subclavian arteries

Stenting is intended to improve and maintain the artery luminal diameter.

### Biliary

The stent is indicated as a palliative treatment of malignant neoplasms in the biliary tree.

### Contraindications

The Protégé EverFlex self-expanding peripheral stent system is contraindicated under the following conditions:

#### Peripheral

- Patients in whom anticoagulant and/or antiplatelet therapy is contraindicated
- Patients who exhibit persistent acute intraluminal thrombus of the proposed lesion site, after thrombolytic therapy
- Perforation at the angioplasty site evidenced by extravasation of contrast medium
- Aneurysm of the artery to be treated
- Patients with known hypersensitivity to nickel-titanium
- Patients with uncorrected bleeding disorders
- All the customary contraindications for PTA

#### Biliary

- None known

### Warnings

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the stent system or the vessel lumen. If resistance is felt, carefully withdraw the stent system without deploying the stent.

- If resistance is felt when pulling back the distal grip, do not force deployment. Carefully withdraw the stent system without deploying the stent.
- If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the system as one unit.

### Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
- Do not use the stent if it is partially deployed.
- Support from a sheath is necessary to minimize lengthening or shortening of the stent during stent deployment.
- Always use a sheath during the implant procedure to protect the vessel (peripheral indication) or the biliary tract (biliary indication), and the puncture site.
- If the lesion is not predilated, it may be difficult to properly position or remove the stent system.
- The stent is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition (peripheral indication), or duct apposition (biliary indication).
- Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.
- The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.
- The effects of overlapping stents have not been evaluated.
- Use caution when crossing a deployed stent with any adjunct device.
- Do not expand the stent beyond its nominal diameter.

### Potential adverse events

Potential adverse events from the stent implantation procedure include the usual complications associated with any percutaneous procedure, PTA, or stent placement.

#### Table 1. Potential adverse events

##### Risks associated with any percutaneous procedure

- Arteriovenous fistula
- Allergic reaction to device materials or procedure medications
- Aneurysm
- Artery perforation or rupture
- Bleeding requiring transfusion
- Contrast medium reaction or renal failure
- Death
- Duct rupture resulting from overstretching the duct (biliary only)
- Embolism
- Fever
- Hematoma
- Hypertension or hypotension
- Infection (for biliary placement, secondary to contamination of the stent may lead to cholangitis, hemobilia, peritonitis, or abscess)
- Inflammation
- Intraluminal thrombus
- Myocardial infarction
- Pseudoaneurysm
- Sepsis
- Shock
- Stroke
- Transient ischemic attack
- Venous thromboembolism

##### Risks associated with PTA and stent placement

- Abrupt or subacute closure
- Amputation
- Persons with allergic reaction to nickel titanium (nitinol) may suffer an allergic response to this implant
- Device breakage
- Dissection or intimal flap
- Failure to deploy stent
- Gastrointestinal bleeding due to anticoagulation
- Partial stent deployment
- Renal insufficiency (new or worsening)
- Renal failure requiring dialysis
- Restenosis
- Stent collapse or fracture
- Stent migration
- Stent misplacement
- Surgical or endovascular intervention
- Thrombosis or occlusion of the stent
- Vessel spasm

### Storage

Store the stent system away from sunlight at room temperature in a dry place.

### Directions for use

**Warning:** This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.

### Selection and preparation

#### Required items

- 5 to 10 cc syringe filled with heparinized saline
- 0.035 in (0.89 mm) exchange guidewire
- Hemostatic sheath
- PTA balloon (peripheral only)

#### Select stent size for peripheral interventions

Measure the diameter of the reference vessel (proximal and distal to lesion). Refer to **Table 2** for stent diameter sizing. Measure the length of the target lesion. Choose a stent length that extends proximal and distal to the target lesion.

#### Preparation for biliary interventions

1. Inject contrast medium.
2. Perform a percutaneous cholangiogram using standard technique.
3. Evaluate and mark the stricture by using fluoroscopy, and observe the most distal level of the biliary stricture.

#### Select stent size for biliary interventions

Measure the diameter of the reference bile duct (proximal and distal to stricture). Refer to **Table 2** for stent diameter sizing. Measure the length of the target stricture. Choose a stent length that extends proximal and distal to the tumor to protect against obstruction from further tumor growth.

Table 2. Stent diameter sizing

Stent diameter (mm)	Reference vessel diameter (mm)	Bile duct diameter (mm)
5	3.5–4.5	3.5–4.5
6	4.5–5.5	4.5–5.5
7	5.5–6.5	5.5–6.5
8	6.5–7.5	6.5–7.5

### Preparation of stent delivery system

1. Open the shelf box to reveal the pouch containing the stent and delivery catheter.
2. After careful inspection of the pouch, looking for damage to the sterile barrier, carefully peel open the pouch and extract the tray with contents.
3. Set the tray on a flat surface. Carefully pull off the tray lid and remove the stent and delivery system. If sterility has been compromised or the device is damaged, do not use the device.  
**Caution:** Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
4. Verify that the device is locked by tightening the safety lock clockwise.  
**Caution:** Do not exceed 300 psi (2068 kPa) / 20 atm while flushing the delivery system.
5. Attach a 5 to 10 cc syringe filled with heparinized saline to the stopcock on the manifold. Open the stopcock and inject saline into the annular space between the inner shaft and the outer sheath until it comes out the outer sheath.
6. Attach a 5 to 10 cc syringe filled with heparinized saline to the proximal luer lock injection hub. Inject the saline solution through the guidewire lumen until it comes out the catheter tip.
7. Ensure that the stent is contained within the outer sheath by examining the distal end of the catheter. Do not use the stent system if the stent is partially deployed. If a gap exists between the catheter tip and the outer sheath, open the safety lock and gently pull the inner shaft in a proximal direction until the gap is closed. Secure the safety lock after making the adjustment by turning the lock clockwise.  
**Caution:** Do not use the stent if it is partially deployed.

### Stent deployment procedure

#### Insertion of sheath and guidewire

1. To gain access to the appropriate site, use a sheath with a hemostatic valve that is compatible with a 6 Fr (2.0 mm) delivery system. If the iliac artery is treated, the access site is the common femoral artery. If the subclavian artery is treated, the access site is the brachial or axillary artery.  
**Caution:** Support from a sheath is necessary to minimize lengthening or shortening of the stent during stent deployment.
2. Insert an 0.035 in (0.89 mm) exchange guidewire of appropriate length across the target lesion (peripheral indication) or the stricture (biliary indication) via the sheath.  
**Caution:** Always use a sheath during the implant procedure to protect the vessel (peripheral indication) or the biliary tract (biliary indication), and the puncture site.

#### Dilation of lesion

##### Peripheral

Predilate the lesion using standard PTA techniques. Remove the PTA balloon from the patient while maintaining lesion access with the guidewire.

**Caution:** If the lesion is not predilated, it may be difficult to properly position or remove the stent system.

##### Biliary

Predilation is not typically performed on malignant strictures. However, if predilation is necessary, use standard balloon dilation techniques. Remove the balloon catheter from the patient while maintaining stricture access with the guidewire.

#### Introduction of stent delivery system

Advance the device over the guidewire through the hemostatic valve and sheath.

**Warning:** If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may cause damage to the stent system or the vessel lumen. If resistance is felt, carefully withdraw the stent system without deploying the stent.

#### Stent deployment

1. Advance the delivery system until the distal (leading) radiopaque inner shaft marker is distal to the target lesion (peripheral indication) or the stricture (biliary indication).  
**Note:** When deploying the 200 mm stent delivery system, the stent moves back approximately 5 mm from the distal retainer after the initial release.
2. Pull back on the delivery system until there is no slack in the delivery system and the radiopaque inner shaft markers extend distal and proximal to the target lesion (peripheral indication) or the stricture (biliary indication).
3. Open the safety lock by turning it counterclockwise.
4. Initiate stent deployment by pinning down the inner shaft (proximal grip) in a fixed position and pulling the outer sheath (distal grip) toward the proximal grip as shown in **Figure 4**.
5. Once initial deployment is visible and before achieving vessel apposition (peripheral indication) or duct apposition (biliary indication), reposition the stent as needed using the radiopaque markers.  
**Note:** To ensure that there is no relative movement between the grips during repositioning, secure the safety lock.  
**Caution:** The stent is not designed for recapturing at any time, and is not designed for repositioning after establishing vessel apposition (peripheral indication), or duct apposition (biliary indication).
6. During release of the stent, keep the entire length of the flexible deployment system as straight as possible. To ensure that no slack is introduced into the delivery system, hold the proximal grip stationary and fixed. Deployment is complete when the outer sheath marker passes the proximal inner shaft stent marker and the stent is released.  
**Warning:** If resistance is felt when pulling back the distal grip, do not force deployment. Carefully withdraw the stent system without deploying the stent.


Figure 4. Stent deployment

**Caution:** Failure to hold the proximal grip in a fixed position during deployment of the stent may result in partial deployment, foreshortening, lengthening, or increased deployment force.

**Caution:** The stent is not designed to be lengthened or shortened relative to its nominal length. Excessive stent lengthening or shortening may increase the risk of stent fracture.

**Note:** If a second stent is needed, place the more distal stent first. If overlap of sequential stents is necessary, keep the amount of overlap to a minimum.

**Caution:** The effects of overlapping stents have not been evaluated.

#### After stent deployment

1. While using fluoroscopy, withdraw the entire delivery system as one unit, over the guidewire, into the access sheath and out of the body. Remove the delivery system from the guidewire.  
**Warning:** If resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the system as one unit.
2. Using fluoroscopy, visualize the stent to verify full deployment.
3. If the stent is not completely expanded at any point along the lesion (peripheral indication) or stricture (biliary indication), balloon dilation (standard PTA technique) can be performed after deploying the stent. To dilate the stent, select an appropriately sized PTA balloon catheter and dilate the stent by using conventional technique. Select a PTA balloon with an inflation diameter that approximates the diameter of the reference vessel. Remove the PTA balloon from the patient.  
**Caution:** Use caution when crossing a deployed stent with any adjunct device.  
**Caution:** Do not expand the stent beyond its nominal diameter.
4. Remove the guidewire and the sheath from the body.
5. Close the entry wound as appropriate.
6. Discard the delivery system, guidewire, and sheath.  
**Note:** Physician experience and discretion determines the appropriate drug regimen for each patient.

#### MRI Conditional

#### Stent migration and image artifact

Comparison to published test results indicate that the Protégé EverFlex stent will not exhibit substantial magnetic field interactions or migrate in an MR environment of 3.0 Tesla or less. MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent.

#### RF heating

Radio frequency (RF) heating in an MR environment was evaluated for the Protégé EverFlex stent. Testing was conducted on two stent/MR scanner scenarios: a 200 mm length stent with a 1.5 Tesla scanner and a 120 mm stent with a 3.0 Tesla scanner. These combinations were chosen since they represent the theoretical worst-case conditions for RF heating.

Protégé EverFlex Stent, single configuration, 8 x 200 mm length, SFA position (or overlapped stents totaling approximately 200 mm in length) was tested for RF heating at 1.5 Tesla. 3.0 Tesla RF heating testing was not performed on the 8 x 200 mm.

#### Recommendations:

If performing an MR procedure with a patient with a 200 mm Protégé EverFlex stent (or overlapping stents approximately 200 mm in length), exclude the stent area from RF exposure.

The local body SAR should be limited to  $\leq 0.9$  W/kg if an MR body coil must be used (the recommended SAR level is based on theoretical extrapolation to maintain localized temperature increase  $\leq 2.0$  °C based on in vitro test results, valid for test configuration only). Note that the whole body averaged SAR displayed by the MR system is inappropriate to scale exact local temperature increases.

Gradient magnetic fields: a stimulation level of 63.80% was used (1.5 T Siemens) in respect to heating issues. No tests have been performed on possible nerve or other tissue stimulation.

Protégé EverFlex stent, single configuration, 8 x 120 mm length, SFA position, was tested for RF heating at 3.0 Tesla. 1.5 Tesla RF heating testing was not performed on the 8 x 120 mm.

#### Recommendations:

Local-body SAR  $\leq 3.1$  W/kg with MR body coil (the recommended SAR level is based on theoretical extrapolation to maintain localized temperature increase  $\leq 2.0$  °C based on in vitro test results, valid for test configuration only). Note that the whole-body-averaged SAR displayed by the MR system is inappropriate to scale exact local temperature increases. Gradient magnetic fields: stimulation level of 74.19% was used (3T Siemens) in respect to heating issues. No tests have been performed on possible nerve or other tissue stimulation.

#### Warranty

This product is subject to standard Medtronic warranty terms.

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© 2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™\* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company
**ev3 Inc.**
Medtronic Australasia Pty Ltd