## INSTRUCTIONS

### EXTRACORPOREAL BLOOD CIRCUIT FOR BLOOD PURIFICATION EQUIPMENT

**EN** Extracorporeal Blood Circuit for Blood Purification Equipment

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
### SCHEMATIC DIAGRAM OF EXTRACORPOREAL BLOOD CIRCUIT FOR PURIFICATION EQUIPMENT

# **EN** Extracorporeal Blood Circuit for Blood Purification

## **Equipment**

### **PRODUCT NAME**

Extracorporeal Blood Circuit for Blood Purification Equipment

### **MODELS**

JRHLL-040, JRHLL-041, JRHLL-010, JRHLL-011, JRHLL-020, JRHLL-021, JRHLL-022, JRHLL-023, JRHLL-024, JRHLL-025, JRHLS-001

### **STRUCTURE AND PERFORMANCE OF PRODUCT**

Tubings: arterial line, venous line, heparin line, pressure sensor line, sampling line, priming line, fluid infusion line, fluid discharge line, and secondary circulation line.

Accessories: two-way luer connector, rotating connector, clamp, tubing, sampling connector, trifurcated adapter, luer-lock connector, male luer connector, roller clamp, drip chamber, spike, contamination shield, monitoring pillow, bifurcated adapter, pump tubing, luer-lock connector cap, dropping bottle, dialysis connector, waste bag, filter screen and sensor protective connector.

The above tubing and accessories should be combined to form a complete set of extracorporeal blood circuit as per different clinical treatment modes.

All PVC raw materials are medical grade.

### **STORAGE CONDITIONS**

This product should be stored under temperature between 0°C-35°C, and in cool, dry, well-ventilated and clean place at 80% RH maximum without corrosive gas. Avoid pressing heavily and hitting and getting wet in the rain or snow.

### **STERILIZATION METHOD**

The product is sterile and non-pyrogen after ethylene oxide sterilization.

### **INTENDED USE**

This product is designed for single-use extracorporeal blood circuit for blood purification equipment, to provide blood access channel for blood purification.

### **CONTRAINDICATIONS**

There are no special contraindications for the extracorporeal blood circuit for blood purification equipment during regular dialysis.

### **PRECAUTIONS AND WARNINGS**

1. Read the instruction manual carefully before use.
2. The product is for single use only. Use promptly after the packaged is opened and discard after use.
3. Do not use the product if internal package is damaged, contaminated or damp.
4. Do not use the product if blood circuit is leaked or tubing and components damaged during priming and treatment.
5. The maximum positive pressure is 66.7 kPa. The maximum negative pressure is 62.2 kPa. The maximum blood flow is 400 mL/min.
6. The product is intended to be used only with the dialysis machine and the dialyzer.
7. There may be 10% difference between the actual blood flow and the value displayed on the machine.
8. Do not kink or knot the tubing during haemodialysis. Doing so may apply excessive local pressure to blood circuit.
9. Air may enter blood circuit during priming or operation. Evacuate all air promptly. It is recommended to keep the blood tubing and components at the air monitoring position close to the air detector of the dialysis machine to improve the monitoring accuracy of the detector.
10. Where the external or internal parts of the connector contact the lubricant, the adapter may be separated.
11. If the air is introduced via syringe at the sampling site located at the back end of the air detector, the air detector may fail to induct the air.
12. Biohazard waste should be treated seriously according to national MOH regulation.

### **INDICATIVE EXPLANATION**

1. In case sensor protective connector is wet by normal saline or blood during dialysis, fasten the clamp on the side tube and remove sensor protective connector from dialysis machine. Remove sensor protective connector from wing-type luer-lock connector and replace with a new one.
2. The extracorporeal blood circuit is intended to be used with hollow fiber dialyzer, plasma separator, perfusion apparatus, internal fistula puncture needle, infusion apparatus, and waste bag.
3. During blood drawing or components replacement, sterilize the surface of the puncture site or replaced components with medical alcohol or iodine.
4. The maximum outer diameter of the needle used for blood sampling should be 0.8mm.
5. The shelf-life is 3 years for sterility.
6. The product contains DEHP.
7. For high-risk group, such as the newborns, preadolescent males, pregnant and

breast-feeding females, the physician should use the alternatives as practical as possible.

8. The product is not applicable for infusion of fat-soluble solution such as the fat emulsion.

9. Clinical medical staff should note that the efficacy may be changed due to the interaction between the PVC tubing and the medicine infused.

10. Do not use the product to infuse medicines that are incompatible with PVC.

### **INSTRUCTIONS FOR INSTALLATION AND USE**

1. Check the expiry date prior to use. Inspect if the package is damaged and contents completed.
2. Take out the tubing from the package, and connect the arterial line (in red) and the venous line (in blue) to the two corresponding ends of the dialyzer.
3. Use 1L heparin normal saline to prime blood circuit and the dialyzer (the dose of heparin should be determined by the physician-in-charge).
  - Use heparin normal saline to prime the arterial line of blood circuit, the dialyzer and the venous line of the blood circuit successively. Note: install the pump tubing as per the instructions for use. Improper operation may result in damaged blood circuit.
  - Keep the priming volume inside the arterial and venous chambers no more than 2/3 of the total volume.
  - Prime the blood circuit sufficiently and evacuate all air from blood circuit and dialyzer.
4. After priming blood circuit, connect the arterial line of blood circuit to the arterial channel of the patient. Pump the blood to flow through the dialyzer, the venous line, the venous channel and finally back to the vein of the patient. Haemodialysis is initiated.
  - Check whether the lines and dialyzer are correctly installed, whether the lines are damaged or falsely installed, whether all ports are tightly connected, whether the tubing is kinked, whether the switch clamp is in normal condition, and whether the sensor protective connector is firmly installed before haemodialysis is initiated.
  - Note that the connection port should be tightened and locked as it connects to patient's fistula needle.
5. After completion of dialysis, connect the arterial line of blood circuit to normal saline. Pump the blood inside the blood circuit back to the vein of the patient with normal saline. Stop the blood pump. Clamp the blood inlet and outlet port and the fluid infusion port. Disconnect venous line of blood circuit from patient's venous cannula.
6. It is recommended that the product should be used with sphygmomanometer during haemodialysis.

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-  **CE** 0123 CE Mark
-  Do not re-use
-  Caution
-  **STERILE** EO Sterilization using ETO
-  Date of manufacture
-  Use-by date
-  **LOT** Batch code
-  Do not use if package is damaged
-  Consult instructions for use
-  Temperature limit
-  **EC REP** Authorized representative in the European Community
-  Manufacturer

B00E00/21FI

**EC REP MedNet EC-REP GmbH**

Borkstrasse 10, 48163 Muenster, GERMANY.

**Shandong Weigao Blood Purification Products Co., Ltd.**

No.20 Xingshan Road, Weihai Torch Hi-tech Science Park, 264210 Weihai,

Shandong Province, PEOPLE'S REPUBLIC OF CHINA.

Tel: +86-631-5660598 Fax: +86-631-5660598

**WEGO**

THE BLOOD PURIFICATION EXPERT