

BACTICID WIPES

ALCOHOL BASED RAPID DISINFECTION WIPES FOR MEDICAL DEVICES

- ready-to-use disinfection wipes pre-soaked with Bactacid
- aldehyde, colourant and fragrance free
- rapid disinfection with broad biocidal activity
- residue free



Bactacid wipes are ethanol-based rapid disinfection wipes for non-invasive medical devices and alcohol resistant surfaces. Disposable wipes are pre-soaked with Bactacid solution and come in a convenient and ready to use packaging. Perfect product in conditions with higher hygiene requirements, where rapid and regular surface disinfection is significant. It does not contain aldehydes, fragrance and colorants and is suitable for all alcohol-resistant surfaces.

Bactacid wipes have a broad spectrum of bactericidal (incl. MRSA, VRE), mycobactericidal (incl. Mycobacterium avium and Mycobacterium terrae), fully virucidal (incl. all enveloped viruses, Rota, Adeno), fungicidal and yeasticidal activity.

APPLICATION AREAS:

Alcohol-based rapid disinfection wipes for medical devices and other alcohol resistant surfaces, where short contact times are required, for example in the following areas:

- Hospitals and other medical institutions
- Operating theatres and procedure rooms
- Intensive care unit
- In accident and emergency departments
- Dialysis departments
- Family doctor offices, first-aid stations
- Ambulances
- Laboratories and diagnostics.

- Dentist's offices
- For home care of patient
- Food industry and food processing areas
- Pharmaceutical and cosmetics industry

USAGE:

BACTICID wipes are ready for use. Open the cover of the bucket and pull out the impregnated wipe up to the perforation. To tear off the wipe, pull it suddenly and at a small angle. After tearing off the wipe, close the cover immediately to avoid drying of the wipes.

Thoroughly wipe the surface with the impregnated wipe, ensuring complete coverage. The surfaces must remain visibly moistened for 30 seconds to take effect.

Rinsing of the surfaces is not required. Mark the date of the first use on the product packaging.

Bactacid wipes is a perfect product for rapid disinfection of surfaces like patient nearby surfaces, operating tables, examination couches, stretchers, patient chairs, wheelchairs, medical equipment, trays, trolleys, hospital and maternity beds.

If the surface tolerance is unknown, try it out on an inconspicuous area first. This product is not suitable for acrylic glass, plexiglas and alcohol soluble varnish.

MICROBIOLOGY:

Efficacy	Concentration	Contact time
Escherichia coli	Ready to use	30 sec
Enterococcus hirae	Ready to use	30 sec
Staphylococcus aureus	Ready to use	30 sec
Pseudomonas aeruginosa	Ready to use	30 sec
Bacteria (incl. MRSA)	Ready to use	30 sec
Mycobacteria: (M.Avium, M.Terrae)	Ready to use	30 sec
Candida albicans	Ready to use	30 sec
Aspergillus Niger	Ready to use	5 min
Aspergillus brasiliensis	Ready to use	60 sec
All enveloped viruses (incl. HIV, HBV, HCV, Vaccinia)	Ready to use	30 sec
Rotavirus, Adeno	Ready to use	30 sec
Noro virus	Ready to use	15 sec

BACTICID wipes are tested according to the following European standards: EN 1040, EN 1275, EN 1276, EN 1650, EN 13624, EN 13697, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563, EN16615.




INGREDIENTS:

100 g of product contains active ingredients: ethanol 72 g, Alkyl (C12-18) dimethylbenzyl ammonium chloride (ADBAC (C12-18)) 0,0125 g, Alkyl (C12-C14) ethylbenzylammonium chloride (ADEBAC (C12-C14)) 0,0125 g.


PRECAUTIONS:

Flam. Liq. 2: H225 - Highly flammable liquid and vapour P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking P233: Keep container tightly closed P403+P235: Store in a well-ventilated place. Keep cool P501 - Dispose of contents/container in accordance with local/ regional/ national/ international regulation.

PACKAGING:



100 pcs bucket
16x20cm wipes



150 pcs bucket
20x20cm wipes



200 pcs bucket
20x20cm wipes



600 pcs bucket
18x20cm wipes



AS Chemi-Pharm
Tänassilma tee 11, Harju maakond, 76406, Eesti
www.chemi-pharm.com
chemi-pharm@chemi-pharm.com