Perifix®









GB	Set for continuous epidural anaesthesia	4
DE	Set zur kontinuierlichen Epiduralanästhesie	5
BG	Набор за непрекъсната епидурална анестезия	6
CN	持续硬膜外麻醉包	7
CZ	Souprava pro kontinuální epidurální anestézii	8
DK	Sæt til kontinuerlig epiduralanæstesi	9
EE	Püsiva epiduraalanesteesia komplekt	10
ES	Set para una anestesia epidural continuada	11
FI	Jatkuvan epiduraalipuudutuksen sarja	12
FR	Set pour anesthésie épidurale continue	13
GR	Σετ για συνεχή επισκληρίδια αναισθησία	14
HR	Set za kontinuiranu epiduralnu anesteziju	15
HU	Készlet folytonos epiduralis érzéstelenítéshez	16
ID	Set untuk anestesi epidural secara terus-menerus	17
IT	Set per anestesia epidurale continua	18
KR	지속형 경막외 마취 세트	19
LT	Nenutrūkstamos epiduralinės anestezijos priemonių rinkinys	20
LV	Ilgstošās epidermālās anestēzijas komplekts	21
NL	Set voor continue epidurale anesthesie	22
NO	Sett for kontinuerlig epiduralanestesi	23
PL	Zestaw do ciągłego znieczulania zewnątrzoponowego	24
PT	Conjunto para anestesia epidural, contínua	25
RO	Set pentru anestezie epidurală continuă	26
RS	Set za kontinuiranu epiduralnu anesteziju	27
RU	Набор для продленной эпидуральной анестезии	28
SE	Set för kontinuerlig epiduralanestesi	29
SI	Set za kontinuirno epiduralno anestezijo	30
SK	Súprava na kontinuálnu epidurálnu anestéziu	31
TH	ชุดสำหรับระงับความรู้สึกผ่านไขสันหลัง	32
TR	Kontinü epidural anestezi seti	33
UA	Набір для продовженої епидуральної анестезії	34
VN	Bộ dành cho việc gây mê ngoài màng cứng liên tục	35

B | BRAUN

GB	Line drawings of the individual components
DE	Strichzeichnungen der Komponenten
BG	Скици на отделните съставни части
CN	单独组件的图示
CZ	Perokresby jednotlivých částí
DK	Stregtegninger af de enkelte dele
EE	Üksikute komponentide skeemid
ES	Esquemas de los componentes individuales
FI	Osien merkit
FR	Croquis des différents constituants
GR	Γραμμικά σχέδια των ξεχωριστών στοιχείων

HR	Prikaz pojedinih komponenti
HU	Az egyes alkotórészek vonalas rajzai
ID	Ilustrasi setiap komponen
IT	Schemi dei singoli componenti
KR	각 구성요소의 라인 드로잉
LT	Kiekviena detalė pažymėta brūkšniu kodu
LV	Atsevišķo materiāla sastāvdaļu svītrkodi
NL	Tekeningen van de verschillende onderdelen
NO	Strektegninger av de enkelte bestanddelene
PL	Rysunki pojedynczych części składowych
PT	Traços da marca de cada componente







RO	Desene liniare ale componentelor
RS	Skice pojedinačnih komponenti
RU	Рисунки отдельных компонентов
SE	Bilder på de enskilda komponenterna
SI	Shematski prikazi posameznih sestavnih delov
SK	Perokresbay jednotlivých částí
TH	ภาพส่วนประกอบต่าง ๆ
TR	Komponentlerin ayrı çizimleri
UA	Графічне зображення окремих компонентів
VN	Bản vẽ nét của từng thành phần riêng lẻ

GB Directions
DE Anwendungshinweise
BG Указания за употреба
CN 说明
CZ Pokyny pro použití
DK Betjeningsvejledning
EE Kasutusjuhendid
ES Aplicaciones
FI Käyttöohje
FR Conseils d'utilisation
GR Υποδείξεις εφαρμογής

HR Upute
HU Használati útmutatás
ID Petunjuk
IT Modalità d'uso
KR 방향
LT Įspėjimai
LV Norādījumi par lietošanu
NL Gebruiksaanwijzingen
NO Bruk
PL Przepisy stosowania
PT Instruções de uso

RO Instrucțiuni de folosire
RS Uputstva
RU Указания по применению
SE Användningsanvisningar
SI Napotki za uporabo
SK Aplikačné pokyny
TH วิธีการใช้
TR Kullanım talimatları
UA Вказівки
VN Chỉ dẫn

GB Instructions for use

Contents

(depending on set)

1. Periflex® Epidural Needle, Tuohy bevel (Fig. A)
2. Periflex® Catheter, 1000 mm (Fig. D)
3. Periflex® Catheter Connector (Fig. C)

Optional:

1. Periflex® Filter 0.2 µm (Fig. F)
2. Periflex® PinPad (Fig. J)
3. Periflex® L.O.R. Syringe, latex-free (Fig. B)
4. Single-use syringes (Fig. G, H, I)
5. Single-use needles (Fig. G, H, I)
6. Catheter label (Fig. K)

Materials used

PE, PC, PS, PA, ABS, PP, PEBA, TPE, SAN, stainless steel

Indications

Epidural anaesthesia is a central neuraxial block procedure with use in surgery, obstetrics and pain control.

Epidural anaesthesia (lumbar, thoracic, cervical) for fast and reliable anaesthesia during surgery. If necessary, anaesthesia can be extended and prolonged during surgery. Analgesia is possible at any time after surgery.

Contraindications and restrictions on use

The product should only be used by physicians that have been adequately trained in this technique. Absolute contraindications include patient refusal, severe uncorrected hypovolemia, increased intracranial pressure, infection at the site of injection and known hypersensitivity to local anaesthetics. Relative contraindications include coagulation disorders, anticoagulant therapy, anatomical abnormalities of the vertebral column, chronic low back pain, neurological diseases (e.g. nerve lesions, polyneuropathy, multiple sclerosis), cardiovascular diseases (inability to increase cardiac output) and sepsis. Additional contraindications for cervical epidural anaesthesia are sympathetic blockade in patients with congestive heart failure or severe aortic or mitral stenosis or congestive lung disease.

CAUTION:

Do not perform cervical epidural anaesthesia for paediatric anaesthesia!

Please consult medical textbooks for standard contraindications to performing continuous epidural anaesthesia procedures.

Risks

Known risks of epidural anaesthesia techniques include, in particular, hypotension, dyspnoea / apnoea, cardiac arrest, bradycardia, hypothermia, urinary retention, toxic effects of the local anaesthetic / adjuvant (please refer to the Pharma-SPC's), inadequate block, temporary or permanent neurological complications and local or systemic infections (postdural puncture headache (PDPH), paraesthesia, meningitis / cauda equine syndrome, transverse myelitis or anterior spinal artery syndrome with permanent paralysis, epidural haematoma/epidural space infections, abscesses, skin and subcutaneous tissue infections), inadvertent puncture of blood vessels.

Additional risks for cervical epidural anaesthesia include spinal cord damage, T1 – T4 sympathetomy, phrenic nerve paresis or paralysis, decrease in heart rate by blocking sympathetic heart innervations and modulation of the vagal nerve and blood flow reduction in the vertebralbasilar artery system.

Trained staff and sufficient emergency equipment, in particular for artificial ventilation and haemodynamic support, must be nearby. During and after cervical epidural anaesthesia, cardiovascular and respiratory monitoring is important to exclude cardio-pulmonary insufficiency.

Warnings

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Should not be used in patients with known hypersensitivity to any of the materials employed.

Do not use if packaging is damaged.

Do not resterilise.

Do not wet Periflex® Filter 0.2 µm with disinfectants containing alcohol as this may result in cracks in the filter housing.

Do not use syringes smaller than 10 ml as they may cause the filter housing to be damaged as a result of excessive pressure from the syringe.

Periflex L.O.R. Syringe:

Special syringe for the Loss of Resistance technique for regional anaesthesia.

Syringe, only suitable to a limited extent for aspiration of liquids. Do not aspirate via the mechanical stop. For hygiene reasons, do not fill the syringe until shortly before use.

Duration of application

Inspect the insertion site of the catheter daily. Remove the catheter in the case of any signs of local or systemic inflammation of unknown origin.

Directions

CAUTION:

– Use aseptic techniques.

– Never pull the catheter through the needle as it may otherwise shear off.

1. To make the puncture use the Loss of Resistance technique (technique of choice), or use the hanging drop technique or use the Macintosh balloon.
2. After making the puncture and identifying the epidural space, use the threading assist guide (E) (Fig. 1) and insert the Periflex® epidural catheter (D) through the epidural needle into the epidural space to the required position. (The catheter is marked on the tip and at distances of 1 and 5 cm. The place where the catheter tip leaves the epidural needle has a broad mark.)
3. Pull the epidural needle back over the catheter, insert the catheter as far as possible into the Periflex® Catheter Connector (C) and close (Fig. 2–5). The Catheter-Connector works according to the Click & Ready principle.
 - a) Open the catheter connector lid as far as possible (Fig. 2).
 - b) Insert the catheter into the funnel-shaped opening on the connector tip (Fig. 3).
 - c) Push the catheter well into the connector until it clearly reaches a stop (Fig. 3).
 - d) Close the connector lid by pressing down firmly on the tip of the lid until the lid engages with a clear click (Fig. 4).
 - e) After removal of the protective cap attach the connector to the Periflex® Filter (if provided) with a twist action (Luer Lock) (Fig. 5).

If necessary: Re-opening the connector: Hold the connector firmly between two fingers and press the connector tip on the upper side of the lid firmly downwards. The lid will spring open with an audible click (Fig. 6).

4. Remove air from the Periflex® Filter 0.2 µm (F) using physiological saline solution, connect the filter to the catheter and rinse the catheter with 1–2 ml physiological saline solution.

CAUTION:

The catheter must be inserted into the connector as far as possible as otherwise injection is impossible and it can become disconnected.

Do not use if the packaging is damaged.

Do not resterilise.

The catheter connector and epidural filter should not be wet with disinfectants containing alcohol as this may result in cracks in the housing.

Test dose

Before any epidural injection, it is recommended that an adequate test dose is administered. This test dose may contain a suitable cardiovascular active substance in order to check whether the catheter is inadvertently in a blood vessel. Monitor heart rate continuously in order to detect complications, i.e. tachycardia, immediately. The catheter allows the continuous or intermittent administration of compounds that have been approved for epidural anaesthesia to the epidural space over a period of hours or days.

CAUTION:

Do not use an adhesive spray which contains organic solvents to fix the catheter.

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Verwendete Materialien

PE, PC, PS, PA, ABS, PP, PEBA, TPE, SAN, rostfreier Stahl

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[AU] For sterilization type see primary
packaging
B. Braun Australia Pty. Ltd.
Level 5, 7 – 9 Irvine Place
Bella Vista NSW 2153
Australia

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[CA] Distributed by: B. Braun Medical Inc.
Bethlehem, PA USA 18018-3524

[CO] Importado por B. Braun Medical S.A.
Calle 44 N°8-08/31 Bogotá

[ID] Imported by:
PT. B. Braun Medical Indonesia
Jakarta, Indonesia.

[IN] Imported by:
B. Braun Medical India Pvt. Ltd.,
Bldg: B, Gala No:1-10, Prerna Comp.,
Anjur Phata, Dapoda Rd, Val Vill: B
Taluka-Bhiwandi (Thane-Zone 5),
Maharashtra, Pin 421302-India
Email ID: cs.bind@bbraun.com,
Customer Care No.: 022-6668 2232

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[MY] Authorised representative:
B. Braun Medical Industries Sdn. Bhd.
Bayan Lepas Free Industrial Zone,
11900 Penang, Malaysia.
www.bbraun.com

[RS] Nosilac upisa med. sredstva u Registar
med. sredstava i distributer:
B. Braun Adria RSRB d.o.o.
Milutina Milankovića 11g, Beograd
Br. rešenja o upisu u Registar:
515-02-01391-18-001 od 24.07.2018.

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[US] Rx only • Pictorial is for reference only
Distributed by:
B. Braun Medical Inc. Bethlehem,
PA 18018-3524
Phone 1-800-227-2862

B. Braun Medical AG
Seesatz 17
6204 Sempach
Switzerland

Manufacturing site/Fabricado por:
B. Braun Melsungen AG,
Carl-Braun-Str. 1,
34212 Melsungen, Germany.

CE 0123

1221

12410888

B. Braun Melsungen AG
34209 Melsungen
Germany
www.bbraun.com