



*Nasal CPAP generator
for use with ventilators*

Miniflow



Distributed by:
SLE

When the smallest thing matters

SLE Miniflow

The Miniflow

The Miniflow is a patient interface for the CPAP mode on ventilators and is connected to the inspiratory and expiratory limbs of a dual limb patient circuit, allowing a seamless transition from ventilation to nasal CPAP.

Functions and advantages

- The ventilator provides the necessary gas flow and controls the CPAP Pressure.
- An included T-piece allows easy connection if the ventilator needs an external pressure measurement.
- A wide range of masks and prongs allow easy fitting to babies of any size giving a comfortable, leak-free seal.


Special features include

- Minimal dead space
- Adjustable hinge (45°-60°) that adapts to the baby's anatomy
- Separation of inspiratory and expiratory gases
- Flexible prong adapter for optimal positioning
- Quiet
- Patented non-kinking tubes
- Simple to set-up and use


Tips for use

- Always use the largest prong size possible to avoid leakage.
- Regularly switch between prongs and mask to avoid irritation to the patient's nose.
- Do not insert prongs more than half way into the nose.
- Recommended weight range for Miniflow is 450 g to 4 kg.
- **IMPORTANT:** This is a single use, single patient product.

Incorporating a Miniflow into a patient breathing circuit:



Parts of the Miniflow:



- (A) Inspiratory limb (B) Expiratory limb
(C) Separation wall (D) Adjustable prong adaptor

Miniflow nCPAP Generator	Specifications
Material	Terlux TPE-S Makrolon
Weight	10 g
Connectors	flow (expiratory and inspiratory)
CPAP level setting	0 to > 20 cmH ₂ O
Noise Level (@ 5 cmH ₂ O)	44 dBa
Nasal Connector	Changeable 45° and 60°

Miniflow Parts

Miniflow interface & Starter Kit			
Part No.	Description	Size	Qty
BM4000/10SLE	Miniflow nCPAP starter kit	n/a	1
BM4000SLE	Miniflow nCPAP interface	n/a	20

Miniflow Prongs (100% Silicone)						
Part No.	Description	Size	Qty	Weight (g)	Prong Diameter (mm)	Prong Distance (mm)
BM1026	Nasal prongs	XS	10	0.8	3	5
BM1200/21	Nasal prongs	S	10	0.9	3.5	6
BM12000/02	Nasal prongs	M	10	1.0	4.1	7.5
BM12000/22	Nasal prongs	L	10	1.1	4.75	7.25
BM12000/03	Nasal prongs (Suggested baby weight: 2200-3000g)	XL	10	1.2	5.5	9.5
BM12000/32	Nasal prongs (Wide)	WM	10	1.6	3.7	8.2
BM12000/33	Nasal prongs (Wide)	WL	10	1.5	5	9.7

Miniflow Masks (100% Silicone)					
Part No.	Description	Size	Qty	Weight (g)	Mask Width (mm)
BM1200/08	Nasal mask	XS	10	1.9	14
BM1200/04	Nasal mask	S	10	1.5	18.5
BM1200/05	Nasal mask	M	10	2	21
BM1200/06	Nasal mask	L	10	2.6	23
BM1200/07	Nasal mask (Suggested baby weight: 3000-4000g)	XL	5	4	26

Miniflow Bonnets			
Part No.	Description	Size	Qty
BM1213/10	Bonnet	XXS	10
BM1214/10	Bonnet	XS	10
BM1215/10	Bonnet	S	10
BM1216/10	Bonnet	M	10
BM1217/10	Bonnet	L	10
BM1218/10	Bonnet	XL	10
BM1219/10	Bonnet	XXL	10
BM1220/10	Bonnet	XXXL	10
BM1212/13	Bonnet fixation strip	13 cm	10
BM1212/15	Bonnet fixation strip	15 cm	10
BM1212/18	Bonnet fixation strip	18 cm	10
BM1212/20	Bonnet fixation strip	20 cm	10

Miniflow Accessories		
Part No.	Description	Qty
BM2150	Measuring tape	1 block, 20 strips
BM4010	Pressure line + T-piece	10


Please consult your current price list for all part numbers and pricing, or contact your local SLE representative for more info.

Starter Kit includes:

CPAP Interface with Foam x10



Measuring Tape
Block of 10



Prongs
3 x S 3 x L
3 x M 3 x XL



Masks
1 x S 1 x L
1 x M




Single Use Bonnet
2 x S 3 x L
3 x M 1 x XL



Bonnets are supplied in a range of colours with each accompanied by two fixation strips of appropriate length.

 +44 (0)20 8681 1414

 +44 (0)20 8649 8570

 sales@sle.co.uk

 www.sle.co.uk

SLE Limited
Twin Bridges Business Park
232 Selsdon Road
South Croydon
Surrey
CR2 6PL
UK



When the smallest thing matters