EKVITESTLAB LLC

Velyka Vasylkivska St. 114
03150 Kyiv, Ukraine
Tel. 0-800-31-89-87
e-mail: info@equitest.com.ua
www.equitest.com.ua

STATEMENT

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
Borrelia burgdorferi IgG

**ELISA kit for the qualitative and semiquantitative
detection of IgG antibodies to *Borrelia burgdorferi***

Instructions for use

IVD

REF
EI-801

Σ 96
tests

EQUI *Borrelia burgdorferi* IgG

ELISA kit for the qualitative and semiquantitative detection of IgG antibodies to *Borrelia burgdorferi*

1. INTENDED USE

The «EQUI *Borrelia burgdorferi* IgG» ELISA is a kit intended to qualitatively and semi-quantitatively detect anti-*Borrelia burgdorferi* sensu lato IgG in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) in order to diagnose Lyme disease. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated open immunoassay analysers.

Target group: patients with non-specific infectious symptoms, people who spend lots of time outdoors in wooded, grassy areas or tick bites, inhabitants of endemic areas.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Lyme disease, or Lyme borreliosis, is the most common bacterial infection in the northern hemisphere that is transmitted through the bite of Ixodes mites. Spirochetes of the genus *Borrelia* cause this disease. Lyme disease causes polysystemic affections and without proper treatment in the chronic form causes a number of complications.

Lyme borreliosis is caused by several species and serotypes of bacteria in North America, Europe and Asia, which are grouped into the group *Borrelia burgdorferi* sensu lato. The most common are *Borrelia burgdorferi* sensu stricto, *Borrelia afzelii*, *Borrelia garinii*. These are mobile gram-negative spiral spirochetes 10-25 µm long and about 0.2 µm thick.

After infection with *Borrelia*, an inflammatory reaction occurs in response to surface antigens of *Borrelia burgdorferi*. The development of primary migratory erythema is associated with the interaction of mast cells with lipoproteins of the pathogen. During the first stage of the disease, 1-3 weeks after infection, specific IgM antibodies to *Borrelia burgdorferi* antigens begin to be produced. They are found in the vast majority of people with a history of primary migratory erythema and tick bites. Specific IgG antibodies begin to be synthesized about a month after infection and are detected in high titers in the second stage of Lyme disease, as well as in chronic borreliosis in the absence of proper treatment. After successful therapy, the titer of IgG antibodies decreases after a few months, but they may continue to be low for a long time. The presence of specific antibodies to *Borrelia burgdorferi* antigens does not provide immunity to re-infection.

Asymptomatic course of the disease in the early stages, as well as nonspecific manifestations of tick-borne borreliosis complicate its timely detection. The detection of specific antibodies plays a special role to *Borrelia* antigens in serum or cerebrospinal fluid by ELISA and immunoblotting. However, false-positive results of serodiagnosis can be observed in the presence of other spirochert infections

(syphilis, leptospirosis, etc.). Bacteriological methods of detection of the pathogen and its genetic material by polymerase chain reaction are used to verify the diagnosis. In addition, there is a seronegative course of the disease. Therefore, when diagnosing Lyme borreliosis, both the data of the epidemiological anamnesis and the clinical picture of the infection are taken into account.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Borrelia burgdorferi* specific IgG in «EQUI *Borrelia burgdorferi* IgG» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Recombinant antigens of *Borrelia burgdorferi* sensu lato are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Borrelia burgdorferi* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species IgG monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

		Microplate
		Each plate well is coated with anti- <i>Borrelia burgdorferi</i> sensu lato: <i>B.burgdorferi</i> sensu stricto, <i>B.afzelii</i> , <i>B.garinii</i> . The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months
STRIPS	1 x 96 wells	
		Positive control
		Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C
CONTROL +	1 x 0,25 ml	
		Negative control
		Negative human serum with a preservative (yellow). Store at 2-8 °C
CONTROL -	1 x 0,6 ml	
		Serum dilution solution
		Buffer solution with a milk extract, a detergent and a preservative (violet). Store at 2-8 °C
DIL SAMPLE	1 x 13 ml	
		Conjugate solution (ready to use)
		Buffer solution of monoclonal antibodies to human IgG, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
SOLN CONJ	1 x 13 ml	

SOLN TMB	1 x 13 ml	TMB solution (ready to use) TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
TWEEN WASH 20x	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
SOLN STOP	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use reagents from different lots of the ELISA kit during the assay; do not mix reagents from different lots of the ELISA kit; do not mix reagents from ELISA kits of different nosology; do not use reagents of other manufacturers along with «EQUI Borrelia burgdorferi IgG» kits;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA ki;
- SOLN|TMB solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of SOLN|TMB with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for SOLN|CONJ and SOLN|TMB;
- do not evaluate the test results visually (without a reader);

- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls of ELISA kit «EQUI Borrelia burgdorferi IgG» were tested and found to be negative for anti-HIV1/2, anti-HCV, anti-*T.pallidum* antibodies and HBsAg. Nevertheless, all controls and patient samples should be regarded and handled as potentially infectious;
- some components of the ELISA kit contain low concentrations of harmful agents and could cause irritation of the skin and the mucosa. In the case of contact with [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with skin or mucosa the place of contact should be immediately rinsed with large amounts of water;
- in case of spilling solutions that do not contain acid, e.g. sera, disinfect the surface thoroughly, then dry it with absorbent paper. If the spilling fluid is an acid, it must be initially neutralized with sodium bicarbonate and then use the mechanism described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.

9.2. Fill in the sample application plan.

9.3. Prepare the detergent as per para. 8.2.

9.4. Dispense 90 µl of **DIL|SAMPLE** into each well.

9.5. Add 10 µl of controls and test samples into the wells:

CONTROL|+ – into well A1,

CONTROL|- – into wells B1, C1, D1,

and test samples into the remaining wells.

At the time of adding, colour of the sample diluent changes from violet to blue. Pipette mixture in the wells gently, avoiding foaming.

9.6. Cover strips with an adhesive film and incubate for 30 min at 37°C.

9.7. Remove and discard the adhesive film and wash all wells 6 times with automatic washer or 8-channel pipette as follows:

- aspirate the content of all wells into a liquid waste container;
- add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
- aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
- repeat the washing step 5 more times;
- after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.

9.8. Add 100 µl of [SOLN|CONJ] into all wells. Cover strips with a new adhesive film and incubate for 30 min at 37°C.

9.9. Remove and discard the adhesive film and wash all wells 5 times as described above in 9.7.

9.10. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.

9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.

9.12. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.

9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD of the negative control (\bar{Nc}), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL + OD $\geq 1,5$

CONTROL - OD $\leq 0,150$

CONTROL - $\bar{Nc} \times 0,5 \leq Ncn \leq \times 2,0$ where Ncn is the OD for each Nc run

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and \bar{Nc} is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE*
$IP_{\text{sample}} < 0,9$	NEGATIVE

* Borderline samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

Use of the positivity index allows for a semi-quantitative comparative analysis of the level of specific antibodies in the dynamics of paired samples obtained from patients with an interval of 2-4 weeks. It should be borne in mind that IP_{sample} within 1,1 – 7,0 proportional to the content of specific antibodies. If IP_{sample} is higher 7,0, to correctly assess the content of specific antibodies, it is recommended to re-analyze the sample previously diluted 10 times. When determining the final result in this case, multiply the value obtained IP_{sample} on the degree of dilution (x10).

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 32 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
2	0,449	1,3	7,2
13	0,675	2,0	7,1

Inter assay reproducibility

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
2	0,414	1,2	7,8
13	0,714	2,1	6,9

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 µmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

To determine the sensitivity and specificity of the «EQUI *Borrelia burgdorferi* IgG» ELISA kit in comparison with a similar commercial ELISA kit, 122 samples of donor blood sera (random sampling) and 212 samples of children's blood sera (a total of 334 samples), as well as samples of a commercial panel of sera, were used «Lyme Disease (Anti-*Borrelia burgdorferi*) Mixed Titer Performance Panel PTL202» (SeraCare Life Sciences, USA). The results of determination of antibodies of the IgG class to *Borrelia burgdorferi*, obtained in the «EQUI *Borrelia burgdorferi* IgG» ELISA kit, were compared with the results obtained in similar commercial tests. The relative sensitivity of the «EQUI *Borrelia burgdorferi* IgG» ELISA kit was 99.5%, the relative specificity was 98.8%.

12. LIMITATIONS OF ASSAY

A positive result in the «EQUI *Borrelia burgdorferi* IgG» ELISA kit is evidence of the presence of IgG class antibodies in the patient, specific to *Borrelia burgdorferi* sensu lato, which are produced by the body when a person is infected with the causative agent of Lyme borreliosis.

It should be noted that in the case of early infection, the ELISA result may be negative due to the absence of antibodies at the initial stage of the disease. In the presence of clinical manifestations of the disease, it is recommended to conduct repeated testing after two to four weeks, as well as examine the patient's sample for specific antibodies of the IgM class (for example, in the «EQUI *Borrelia burgdorferi* IgM» ELISA kit).

It is not possible to completely exclude false-positive results, which may be caused by the presence of specific antibodies in the blood in diseases caused by spirochetes (syphilis, typhus, leptospirosis, etc.).

To confirm a positive ELISA result, it is recommended to conduct an additional study of the sample for the presence of antibodies to individual *Borrelia burgdorferi* proteins by immunoblotting. However, the final diagnosis cannot be established only on the basis of serological test results. For the correct diagnosis of tick-borne borreliosis, research should also be conducted to identify the causative agent, for example, using PCR or cultural methods. When establishing a diagnosis, the results of a complex of laboratory and instrumental studies, as well as clinical manifestations of the disease, should be taken into account.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

REFERENCES

1. Bernard Q., Wang Z. et al. Interaction of primary mast cells with *Borrelia burgdorferi* (sensu stricto): role in transmission and dissemination in C57BL/6 mice // *Parasites & Vectors*. - 2017. - 10:313.
2. CDC. Lyme Disease // <https://www.cdc.gov/lyme/index.html>.
3. Cerar T., Strle F., Stupica D. et al. Differences in Genotype, Clinical Features, and Inflammatory Potential of *Borrelia burgdorferi* sensu stricto Strains from Europe and the United States // *Emerging Infectious Diseases*. - 2016. - Vol. 22(5). - P. 818-827.
4. Embers M. E., Hasenkampf N. R. et al. Dynamic Longitudinal Antibody Responses during *Borrelia burgdorferi* Infection and Antibiotic Treatment of Rhesus Macaques // *Clinical & Vaccine Immunology*. – 2012. – Vol. 19, No. 8. – P. 1218-1226.
5. Lantos P. M., Auwaerter P. G., Wormser G. P. A Systematic Review of *Borrelia burgdorferi* Morphologic Variants Does Not Support a Role in Chronic Lyme Disease // *Clinical Infectious Diseases*. - 2014. - Vol. 58(5). - P. 663–671.
6. Lyme Borreliosis (Lyme disease). In: *International travel and health*. Geneva: World Health Organization; 2014 (<http://www.who.int/ith/diseases/lyme/en/>).
7. Tilly K., Rosa P. A. and Stewart P. E. Biology of Infection with *Borrelia burgdorferi* // *Infectious Disease Clinics of North America*. - 2008. - Vol. 22(2). - P. 217–234.
8. Shapiro E. D. *Borrelia burgdorferi* (Lyme Disease) // *Pediatrics in Review*. - 2014. - Vol. 35(12). - P. 500–509.
9. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU.
10. Закон України «Про відходи» // *Відомості Верховної Ради України*. - 1998. - №36-37.
11. Наказ МОЗ України №325 від 08.06.2015 «Про затвердження Державних санітарно-протиепідемічних правил і норм щодо поводження з медичними відходами».
12. Постанова КМУ від 02 жовтня 2013р. №754 «Про затвердження технічного регламенту щодо медичних виробів для діагностики in vitro».
13. Hanna Tolonen, Kari Kuulasmaa, Tiina Laatikainen, Hermann Wolf and the European Health Risk Monitoring Project. Recommendation for indicators, international collaboration, protocol and manual of operations for chronic disease risk factor surveys Part 4.Storage and transfer of serum/plasma samples// Finnish National Public Health Institute 2002// https://thl.fi/publications/ehrm/product2/part_iii4.htm.
14. Surveillance Guidelines for Measles, Rubella and Congenital Rubella Syndrome in the WHO European Region. Annex 3.Collection, storage and shipment of specimens for laboratory diagnosis and interpretation of results// Geneva: World Health Organization; 2012 Dec.

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
	Date of manufacture
	Use by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Keep away from sunlight
	Mark of compliance with technical regulations

Edition 7, 13.12.2022

For questions and suggestions regarding the ELISA kit contact:

Ekvittestlab LLC
Velyka Vasylkivska St. 114, Kyiv, Ukraine, 03150

Tel: 0(800)31-89-87, +38 (044)334-89-87,
e-mail: info@equitest.com.ua, www.equitest.com.ua

ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL|SAMPLE] into the wells (violet)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],
other wells – examined samples
(change of colour from violet to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|CONJ] into all wells (green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN|STOP] into all wells (change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at 450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,3;$$

$$IP_{\text{sample}} = OD_{\text{sample}} / CO$$

\bar{Nc} - the average value of OD 3-x [CONTROL|-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE
$IP_{\text{sample}} < 0,9$	NEGATIVE

Borrelia burgdorferi IgM

ELISA kit for the qualitative detection of IgM
antibodies to *Borrelia burgdorferi*

Instructions for use

IVD

REF
EI-802

Σ 96
tests

EQUI *Borrelia burgdorferi* IgM

ELISA kit for the qualitative detection of IgM antibodies to *Borrelia burgdorferi*

1. INTENDED USE

The «EQUI *Borrelia burgdorferi* IgM» ELISA kit is intended for the qualitative detection of IgM antibodies to *Borrelia burgdorferi* sensu lato in human serum or blood plasma by enzyme-linked immunosorbent assay (ELISA) for the purpose of diagnosing acute Lyme borreliosis. The analysis procedure is designed both for manual setup with automatic pipettes and standard equipment, and for an automatic enzyme immunoassay of the «open» type.

Target group: patients with non-specific infectious symptoms, visits to the forest or a history of tick bites, summer residents, residents of endemic areas.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Lyme disease, or Lyme borreliosis, is the most common bacterial infection in the northern hemisphere that is transmitted through the bite of Ixodes mites. Spirochetes of the genus *Borrelia* cause this disease. Lyme disease causes polysystemic affections and without proper treatment in the chronic form causes a number of complications.

Lyme borreliosis is caused by several species and serotypes of bacteria in North America, Europe and Asia, which are grouped into the group *Borrelia burgdorferi* sensu lato. The most common are *Borrelia burgdorferi* sensu stricto, *Borrelia afzelii*, *Borrelia garinii*. These are mobile gram-negative spiral spirochetes 10-25 µm long and about 0.2 µm thick.

After infection with *Borrelia*, an inflammatory reaction occurs in response to surface antigens of *Borrelia burgdorferi*. The development of primary migratory erythema is associated with the interaction of mast cells with lipoproteins of the pathogen. During the first stage of the disease, 1-3 weeks after infection, specific IgM antibodies to *Borrelia burgdorferi* antigens begin to be produced. They are found in the vast majority of people with a history of primary migratory erythema and tick bites. Specific IgG antibodies begin to be synthesized about a month after infection and are detected in high titers in the second stage of Lyme disease, as well as in chronic borreliosis in the absence of proper treatment. After successful therapy, the titer of IgG antibodies decreases after a few months, but they may continue to be low for a long time. The presence of specific antibodies to *Borrelia burgdorferi* antigens does not provide immunity to re-infection.

Asymptomatic course of the disease in the early stages, as well as nonspecific manifestations of tick-borne borreliosis complicate its timely detection. The detection of specific antibodies plays a special role to *Borrelia* antigens in serum or cerebrospinal fluid by ELISA and immunoblotting. However, false-positive results of serodiagnosis can be observed in the presence of other spirochete infections (syphilis, leptospirosis, etc.). Bacteriological methods of detection of the pathogen

and its genetic material by polymerase chain reaction are used to verify the diagnosis. In addition, there is a seronegative course of the disease. Therefore, when diagnosing Lyme borreliosis, both the data of the epidemiological anamnesis and the clinical picture of the infection are taken into account.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Borrelia burgdorferi* sensu lato specific IgM in «EQUI *Borrelia burgdorferi* IgM» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Recombinant antigens of *Borrelia burgdorferi* are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Borrelia burgdorferi* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species IgM monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS	1 x 96 wells	<p>Microplate Each plate well is coated with anti-<i>Borrelia burgdorferi</i> sensu lato: <i>B.burgdorferi</i> sensu stricto, <i>B.afzelii</i>, <i>B.garinii</i>. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months</p>
CONTROL +	1 x 0,25 ml	<p>Positive control Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C</p>
CONTROL -	1 x 0,6 ml	<p>Negative control Negative human serum with a preservative (yellow). Store at 2-8 °C</p>
DIL SAMPLE	1 x 13 ml	<p>Serum dilution solution Buffer solution with a milk extract, a detergent and a preservative (brown). Store at 2-8 °C</p>
SOLN CONJ	1 x 13 ml	<p>Conjugate solution (ready to use) Buffer solution of monoclonal antibodies to human IgM, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C</p>

[SOLN TMB]	1 x 13 ml	TMB solution (ready to use) TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
[TWEEN WASH 20x]	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
[SOLN STOP]	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use reagents from different lots of the ELISA kit during the assay; do not mix reagents from different lots of the ELISA kit; do not mix reagents from ELISA kits of different nosology; do not use reagents of other manufacturers along with «EQUI Borrelia burgdorferi IgM» kits;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA ki;
- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for [SOLN|CONJ] and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);

- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls of ELISA kit «EQUI Borrelia burgdorferi IgM» were tested and found to be negative for anti-HIV1/2, anti-HCV, anti-*T.pallidum* antibodies and HBsAg. Nevertheless, all controls and patient samples should be regarded and handled as potentially infectious;
- some components of the ELISA kit contain low concentrations of harmful agents and could cause irritation of the skin and the mucosa. In the case of contact with [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with skin or mucosa the place of contact should be immediately rinsed with large amounts of water;
- in case of spilling solutions that do not contain acid, e.g. sera, disinfect the surface thoroughly, then dry it with absorbent paper. If the spilling fluid is an acid, it must be initially neutralized with sodium bicarbonate and then use the mechanism described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with

patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.

9.2. Fill in the sample application plan.

9.3. Prepare the detergent as per para. 8.2.

9.4. Dispense 90 µl of **DIL|SAMPLE** into each well.

9.5. Add 10 µl of controls and test samples into the wells:

CONTROL|+ – into well A1,

CONTROL|- – into wells B1, C1, D1,

and test samples into the remaining wells.

At the time of adding, colour of the sample diluent changes from brown to blue. Pipette mixture in the wells gently, avoiding foaming.

- 9.6. Cover strips with an adhesive film and incubate for 30 min at 37°C.
- 9.7. Remove and discard the adhesive film and wash all wells 6 times with automatic washer or 8-channel pipette as follows:
- aspirate the content of all wells into a liquid waste container;
 - add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
 - aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
 - repeat the washing step 5 more times;
 - after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100µl of [SOLN|CONJ] into all wells. Cover strips with a new adhesive film and incubate for 30 min at 37°C.
- 9.9. Remove and discard the adhesive film and wash all wells 5 times as described above in 9.7.
- 9.10. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for 30 minutes in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD of the negative control (\bar{Nc}), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,25$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL + $OD \geq 1,5$

CONTROL - $OD \leq 0,150$

CONTROL - $\bar{N}_C \times 0,5 \leq N_{cn} \leq \bar{N}_C \times 2,0$

where N_{cn} is the OD
for each N_c run

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and \bar{N}_C is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE*
$IP_{\text{sample}} < 0,9$	NEGATIVE

* Borderline samples are recommended to be tested again in two wells of the ELISA kit. If the results are again indeterminate, a new sample should be collected and analyzed after 7-14 days. In case of repeated reception of indeterminate results, such samples should be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 44 replicates on one series of ELISA kits.

Sample No.	OD_{av}	IP_{av}	CV, %
269	2,186	7,7	5,6
809	0,387	1,4	8,6

Inter assay reproducibility

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD_{av}	IP_{av}	CV, %
269	2,073	7,2	8,1
809	0,364	1,3	7,5

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 $\mu\text{mol/L}$), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

To determine the sensitivity of the «EQUI *Borrelia burgdorferi* IgM» ELISA kit, we used serum samples that were characterized as positive by the Western blot method (59 samples), as well as samples of the commercial serum panel «Lyme Disease (Anti-*Borrelia burgdorferi*) Mixed Titer Performance Panel PTL202 (SeraCare Life Sciences, USA). According to the results of the analysis, the sensitivity of the ELISA kit was 96.9%. To assess the specificity, samples of donor blood sera (122 samples) and 212 samples of children's blood sera (334 samples in total) were used. The results of the determination of IgM antibodies to *Borrelia burgdorferi* obtained in the «EQUI *Borrelia burgdorferi* IgM» ELISA kit were compared with the results obtained in similar commercial tests. The relative specificity of the «EQUI *Borrelia burgdorferi* IgM» ELISA kit was 100.0%.

12. LIMITATIONS OF ASSAY

A positive result in the «EQUI *Borrelia burgdorferi* IgM» ELISA kit is evidence of the presence of IgM class antibodies in the patient, specific to *Borrelia burgdorferi* sensu lato, which are produced by the body when a person is infected with the causative agent of Lyme borreliosis.

It should be noted that in the case of early infection, the ELISA result may be negative due to the absence of antibodies at the initial stage of the disease. In the presence of clinical manifestations of the disease, it is recommended to conduct repeated testing after at least two weeks. A two- or three-fold increase in the level of antibodies indicates the activity of the infectious process.

The results obtained in immunosuppressed individuals should be interpreted with caution.

It is not possible to completely exclude false-positive results, which may be caused by the presence of specific antibodies in the blood in diseases caused by spirochetes (syphilis, typhus, leptospirosis, etc.).

To confirm a positive ELISA result, it is recommended to conduct an additional study of the sample for the presence of antibodies to individual *Borrelia burgdorferi* proteins by immunoblotting. However, the final diagnosis cannot be established only on the basis of serological test results. For the correct diagnosis of tick-borne borreliosis, research should also be conducted to identify the causative agent, for example, using PCR or cultural methods. When establishing a diagnosis, the results of a complex of laboratory and instrumental studies, as well as clinical manifestations of the disease, should be taken into account.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use


Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

REFERENCES

1. Bernard Q., Wang Z. et al. Interaction of primary mast cells with *Borrelia burgdorferi* (sensu stricto): role in transmission and dissemination in C57BL/6 mice // *Parasites & Vectors*. - 2017. - 10:313.
2. CDC. Lyme Disease // <https://www.cdc.gov/lyme/index.html>.
3. Cerar T., Strle F., Stupica D. et al. Differences in Genotype, Clinical Features, and Inflammatory Potential of *Borrelia burgdorferi* sensu stricto Strains from Europe and the United States // *Emerging Infectious Diseases*. - 2016. - Vol. 22(5). - P. 818-827.
4. Embers M. E., Hasenkampf N. R. et al. Dynamic Longitudinal Antibody Responses during *Borrelia burgdorferi* Infection and Antibiotic Treatment of Rhesus Macaques // *Clinical & Vaccine Immunology*. – 2012. – Vol. 19, No. 8. – P. 1218-1226.
5. Lantos P. M., Auwaerter P. G., Wormser G. P. A Systematic Review of *Borrelia burgdorferi* Morphologic Variants Does Not Support a Role in Chronic Lyme Disease // *Clinical Infectious Diseases*. - 2014. - Vol. 58(5). - P. 663–671.
6. Lyme Borreliosis (Lyme disease). In: *International travel and health*. Geneva: World Health Organization; 2014 (<http://www.who.int/ith/diseases/lyme/en/>).
7. Tilly K., Rosa P. A. and Stewart P. E. Biology of Infection with *Borrelia burgdorferi* // *Infectious Disease Clinics of North America*. - 2008. - Vol. 22(2). - P. 217–234.
8. Shapiro E. D. *Borrelia burgdorferi* (Lyme Disease) // *Pediatrics in Review*. - 2014. - Vol. 35(12). - P. 500–509.
9. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU.
10. Закон України «Про відходи» // *Відомості Верховної Ради України*. - 1998. - №36-37.
11. Наказ МОЗ України №325 від 08.06.2015 «Про затвердження Державних санітарно-протиепідемічних правил і норм щодо поводження з медичними відходами».
12. Постанова КМУ від 02 жовтня 2013р. №754 «Про затвердження технічного регламенту щодо медичних виробів для діагностики in vitro».
13. Hanna Tolonen, Kari Kuulasmaa, Tiina Laatikainen, Hermann Wolf and the European Health Risk Monitoring Project. Recommendation for indicators, international collaboration, protocol and manual of operations for chronic disease risk factor surveys Part 4.Storage and transfer of serum/plasma samples// Finnish National Public Health Institute 2002// https://thl.fi/publications/ehrm/product2/part_iii4.htm.
14. Surveillance Guidelines for Measles, Rubella and Congenital Rubella Syndrome in the WHO European Region. Annex 3.Collection, storage and shipment of specimens for laboratory diagnosis and interpretation of results//Geneva: World Health Organization; 2012 Dec.

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
	Date of manufacture
	Use by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Keep away from sunlight
	Mark of compliance with technical regulations

Edition 7, 14.11.2022

For questions and suggestions regarding the ELISA kit contact:

Ekvittestlab LLC
Velyka Vasylkivska St. 114, Kyiv, Ukraine, 03150

Tel: 0(800)31-89-87, +38 (044)334-89-87,
e-mail: info@equitest.com.ua, www.equitest.com.ua

ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL|SAMPLE] into the wells (brown)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL|+], B1, C1, D1 – [CONTROL|-],
other wells – examined samples
(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|CONJ] into all wells (green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN|TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN|STOP] into all wells
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at 450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,25;$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO$$

Nc - the average value of OD 3-x [CONTROL|-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE
$IP_{\text{sample}} < 0,9$	NEGATIVE


Date: 23.03.2021

STATEMENT

We, Vitrotest Bioreagent LLC, having a registered office at _М. Boychuka street 18b/56, Kyiv 01103 Ukraine, assign SRL SANMEDICO having a registered office at A. Corobceanu street 7A, apt. 9, Chişinău MD-2012, Moldova, as authorized representative.

We declare that the company mentioned above is authorized to register, notify, renew or modify the registration of medical devices on the territory of the Republic of Moldova.

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- Baughman A.L., Bisgard K.M., Edwards K.M. et al. Establishment of diagnostic cutoff points for levels of serum antibodies to pertussis toxin, filamentous hemagglutinin, and fimbriae in adolescents and adults in the United States // Clinical and Diagnostic Laboratory Immunology. - 2004. - Vol.11, No.6. - P. 1045-1053.
- European Centre for Disease Prevention and Control. Guidance and Protocol for the serological diagnosis of human infection with *Bordetella pertussis*. Stockholm: ECDC; 2012.
- European Centre for Disease Prevention and Control. Laboratory diagnosis and molecular surveillance of *Bordetella pertussis* - Stockholm: ECDC; 2022.
- Guiso N., Berbers G., Fry N.K. He Q., Riffelmann M., Wirsing von König C.H.; EU Pertstrain group. What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories // European Journal of Clinical Microbiology & Infectious Diseases. - 2011. - Vol.30, No.3. - P. 307-312.

5. Paradowska-Stankiewicz I., Rumik A., Bogusz J. et al. Duration of protection against Bordetella pertussis infection elicited by whole-cell and acellular vaccine priming in Polish children and adolescents // *Vaccine*. – 2021. – Vol.39, No.41. – P. 6067-6073.
6. Riffelmann M., Thiel K., Schmetz J. et al. Performance of commercial enzyme-linked immunosorbent assays for detection of antibodies to Bordetella pertussis // *Journal of Clinical Microbiology*. – 2010. – Vol.48, No.12. – P. 4459-4463.
7. World Health Organization. Laboratory Manual for the diagnosis of whooping cough caused by bordetella pertussis/bordetella parapertussis : update 2014.
8. Xing D., Wirsing von König C.H., Newland P. et al. Characterization of reference materials for human antiserum to pertussis antigens by an international collaborative study // *Clinical and Vaccine Immunology*.- 2009. – Vol.16, No.3. – P. 303–311.

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1. Baughman A.L., Bisgard K.M., Edwards K.M. et al. Establishment of diagnostic cutoff points for levels of serum antibodies to pertussis toxin, filamentous hemagglutinin, and fimbriae in adolescents and adults in the United States // *Clinical and Diagnostic Laboratory Immunology*. - 2004. – Vol.11, No.6. -P. 1045-1053.
2. European Centre for Disease Prevention and Control. Guidance and Protocol for the serological diagnosis of human infection with *Bordetella pertussis*. Stockholm: ECDC; 2012.
3. European Centre for Disease Prevention and Control. Laboratory diagnosis and molecular surveillance of *Bordetella pertussis* - Stockholm: ECDC; 2022.
4. Guiso N., Berbers G., Fry N.K. He Q., Riffelmann M., Wirsing von König C.H.; EU Pertstrain group. What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories // *European Journal of Clinical Microbiology & Infectious Diseases*. – 2011. – Vol.30, No.3. - P. 307-312.
5. Paradowska-Stankiewicz I., Rumik A., Bogusz J. et al. Duration of protection against *Bordetella pertussis* infection elicited by whole-cell and acellular vaccine priming in Polish children and adolescents // *Vaccine*. – 2021. – Vol.39, No.41. - P. 6067-6073.
6. Riffelmann M., Thiel K., Schmetz J. et al. Performance of commercial enzyme-linked immunosorbent assays for detection of antibodies to *Bordetella pertussis* // *Journal of Clinical Microbiology*. – 2010. – Vol.48, No.12. – P. 4459-4463.
7. World Health Organization. Laboratory Manual for the diagnosis of whooping cough caused by *bordetella pertussis/bordetella parapertussis* : update 2014.
8. Xing D., Wirsing von König C.H., Newland P. et al. Characterization of reference materials for human antiserum to pertussis antigens by an international collaborative study // *Clinical and Vaccine Immunology*.- 2009. – Vol.16, No.3. – P. 303–311.

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