Flow sensor
Flow-Sensor
Capteur de débit
Sensor de flujo
Sensor fluxo
Sensore di flusso
流量传感器
フローセンサ
Датчик потока
Akış sensörü

Accessories

Instructions for use
Gebrauchsanweisung
Instructions d'utilisation
Instrucciones de uso
Instruções de uso
Istruzioni per l'uso
使用说明
使用説明書
Инструкции по эксплуатации
Kullanım talimatları

REF 950185

PN 624610/04 | 2017-09-30

CE 0197

HAMILTON
MEDICAL
Intelligent Ventilation since 1983

Table of Contents

English 3

Deutsch 6

Français 9

Español 12

Português 15

Italiano 18

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Türkçe 30

© 2017 Hamilton Medical AG. All rights reserved. Printed in Switzerland. No part of this publication may be reproduced or stored in a database or retrieval system, nor transmitted, in any form or by any means, electronic, mechanical, by photocopying, recording, or otherwise, without the prior written permission of Hamilton Medical.

This manual may be revised or replaced by Hamilton Medical at any time and without notice. Ensure that you have the most current applicable version of this manual; if in doubt, contact Hamilton Medical AG Marketing Department. While the information set forth is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this manual shall limit or restrict in any way Hamilton Medical's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Hamilton Medical has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated and serviced by trained professionals only.

 **WARNING**

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

 **CAUTION**

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

Table of Contents

Intended use

4

Connecting the flow sensor to the ventilator

4

Connecting the flow sensor to the patient

4

Attaching the flow sensor tubes to the breathing circuit tubes

4

Performing a tightness test and visual inspection

5

Calibrating the flow sensor

5

Cleaning and disinfection

5

Disposal

5

Environmental and storage conditions

5

Follow the instructions for use. Consult the appropriate ventilator operator's manual for more information.

Do not use if package is damaged

Quantity

MR unsafe

Autoclavable

Used for pediatric and adult patient groups

⚠ WARNING

- **Hamilton Medical cannot be held responsible for the correct functioning of flow sensors that are not cleaned and sterilized according to Hamilton Medical's reprocessing instructions for the autoclavable flow sensor.**
- **Replace the flow sensor in accordance with hospital infection control procedures, or depending on patient's secretions and nebulization of medication, if required.**

- **Follow hospital infection control procedures for all tests and use.**
- **Visually inspect the flow sensor body, tubing, and internal flap prior to use.**
- **Discard the flow sensor if there is any sign of damage.**

⚠ CAUTION

- *Calibrate the flow sensor prior to patient use. If calibration fails, it may be repeated once. The flow sensor must be discarded if the calibration fails the second time.*
- *Note that the use of rinse aids will reduce the lifespan of the flow sensor. Rinse aids, such as Neodisher Mediklar, lead to early failure and cracks in the plastic flow sensor body.*
- *If reprocessing flow sensors manually, take care that no damage is caused by the use of hard brushes, scouring agents, or by the exertion of too much force. Do not insert brushes into the flow sensor, as this will damage the flap.*
- *(USA only) Federal law restricts this device to sale by or on the order of a physician.*

Intended use

The Hamilton Medical flow sensor is a device to measure patient air flow with

Hamilton Medical ventilators. The autoclavable flow sensor is not to be used in an environment with magnetic resonance imaging (MRI) equipment.

Connecting the flow sensor to the ventilator

Connect the small tubes of the flow sensor to the corresponding color-coded connectors on the ventilator panel.


Connecting the flow sensor to the patient

Insert the flow sensor between the Y-piece of the patient circuit and the patient connection. The end of the flow sensor that should face the patient is marked by "Patient" and an arrow indicating flow direction. It is recommended that you use a short section of flexible tubing between the flow sensor and the actual patient connection.

Attaching the flow sensor tubes to the breathing circuit tubes

The flow sensor tubes must be attached to the breathing circuit tubes to avoid entangling or hooking the flow sensor tubes in other equipment.

To attach the flow sensor tubes to the breathing circuit

- 1 Place the tube attachment strap around both the flow sensor tubes and the breathing circuit tubes as shown. Choose a length that suits your needs, but make sure you do not cut off the breathing current.
- 2 Insert the pointed end with the balls into the opening on the other end of the tube attachment strap.

Performing a tightness test and visual inspection

Carry out a visual inspection and a tightness test prior to operation as described in the *ventilator Operator's Manual*. If the test fails, replace leaking parts and repeat the tightness test.

Calibrating the flow sensor

The flow sensor must be calibrated prior to patient use. To calibrate, follow the instruc-

tions in the appropriate ventilator *Operator's Manual*.

Connect the supplied calibration connector for the calibration. After successful completion of the procedure, the flow sensor is ready for patient use. If calibration fails, it may be repeated once. The flow sensor must be discarded if the flow calibration fails the second time.

Cleaning and disinfection

The flow sensor is delivered clean and ready for clinical use. After use, the flow sensor has to be cleaned, disinfected, and sterilized immediately. After processing, the flow sensor must be visually inspected for damage and optical changes. Worn or damaged flow sensors must be removed from use, as these no longer fulfill their function. Find more information on cleaning and disinfection in the *reprocessing instructions for the autoclavable flow sensor*.

Disposal

The flow sensor must be handled as a contaminated item. Follow all local, state, and federal regulations with respect to environmental protection when disposing of used flow sensors.

Environmental and storage conditions

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