COD 12500 10 x 50 mL

Only for in vitro use in the clinical laboratory

INTENDED USE

Reagents for the measurement of protein (total) concentration in human serum or plasma. The obtained values are useful as an aid in the evaluation of the states of dehydration and water intoxication and salt retention syndromes, and to control the evolution of multiple myeloma, and to control Wadenström macroglobulinemia.

This reagent is for use in the BioSystems A25 and A15 analyzers.

CLINICAL SIGNIFICANCE

Most of the plasma proteins are synthesized by the liver. The major exception to this is the immunoglobulins which are produced by plasma cells found in the spleen, lymph nodes and bone marrow.

The two general causes of alterations of serum total protein are a change in the volume of plasma water and a change in the concentration of one or more of the serum proteins.

Hyperproteinemia can be caused by dehydration (inadequate water intake, severe vomiting, diarrhea, Addison's disease, diabetic acidosis) or as a result of an increase in the concentration of specific proteins (immunoglobulins in chronic infections, multiple myeloma)^{1,2}

Hypoproteinemia may be caused by hemodilution (salt retention syndromes, massive intravenous infusions), by an impaired synthesis (severe malnutrition, chronic liver disease, intestinal malabsorptive disease), or by an excessive protein loss due to a chronic kidney disease or severe burns^{1,2}

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Protein in the sample reacts with copper (II) ion in alkaline medium forming a coloured complex that can be measured by spectrophotometry³.

CONTENTS AND COMPOSITION

- A. Reagent. 10 x 50 mL. Copper (II) acetate 6 mmol/L, potassium iodide 12 mmol/L, sodium hydroxide 1.15 mol/L, detergent.
 - DANGER: H314: Causes severe skin burns and eye damage. P280: Wear protective gloves/protective clothing/eye protection/face protection. P303+P361+P353: IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

For further warnings and precautions, see the product safety data sheet (SDS).

STORAGE AND STABILITY

Store at 2-30°C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 2 months

Indications of deterioration: Absorbance of the blank over the limit indicated in "Test Parameters"

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044)

REAGENT PREPARATION

Reagent is provided ready to use

SAMPLES

Serum or heparinized plasma collected by standard procedures. Stable for 4 weeks at 4-8°C4. Anticoagulants other than heparin should not be used.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 60 days, after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, 18009 and 18042) and II (cod. 18007, 18010 and 18043) to verify the accuracy of the measurement procedure

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits

REFERENCE VALUES

Serum, adults1:



06/2020 8435287110TPROT00000JK

Trueness: Results obtained with this reagent did not show systematic differences when

compared with reference reagents. Details of the comparison experiments are available on request

Concentrations are lower in child. Plasma total protein concentration is 2 to 4 g/L higher due to

These ranges are given for orientation only; each laboratory should establish its own reference

The metrological characteristics described below have been obtained using an A25 analyzer.

Linearity limit: 150 g/L. For samples with higher values, dilute manually or refer to the Test

Parameterization for Automatic dilution (note that all these samples will be diluted with the

Repeatability (CV)

1.0 %

1.2 %

the presence of fibrinogen as well as some other trace proteins¹

LIMITATIONS OF THE PROCEDURE

METROLOGICAL CHARACTERISTICS

Mean concentration

51.8 g/L

82.1 a/L

Results are similar with A15.

Detection limit: 1.6 a/L.

same dilution ratio).

Precision

Interferences: bilirubin (up to 30 mg/dL), hemolysis (hemoglobin up to 500 mg/dL) and lipemia (triglycerides up to 975 mg/dL) do not interfere. Other drugs and substances may interfere⁵

BIBLIOGRAPHY

- 1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
- 2 Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
- 3. Gornall AG, Bardawill CS, David MM. Determination of serum proteins by means of the Biuret reaction. J Biol Chem 1949; 177: 751-766.
- 4. Word Health Organization (WHO). Use of anticoagulants in diagnostic investigations Document WHO/DIL/LAB/99.1, Rev.2; laboratory 2002. http://www.who.int/iris/handle/10665/65957.
- 5. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

TEST PARAMETERS

R1: use Reagent A

R2: use R

GENER

PROCE

eagent B.	

PROTEIN TOTAL	PROTEIN TOTAL	
SER	SER	
endpoint mon.	endpoint mon.	
g/L	g/L	
no	no	
0	0	
increasing	increasing	
bichrom.	bichrom.	
535	535	
670	670	
4	4	
300	300	
-	-	
1.2	1.2	
21	14	
-	-	
-	-	
-	-	
	PROTEIN TOTAL SER endpoint mon. g/L no 0 increasing bichrom. 535 670 4 300 - 1.2 21 - 21 -	PROTEIN TOTAL SER PROTEIN TOTAL SER endpoint mon. g/L g/L g/L no no 0 0 increasing increasing bichrom. bichrom. 535 535 670 670 4 4 300 300 - - 1.2 1.2 21 14 - - - - - -

A25

Predilution factor	-	-
Predilution reduced factor	2	2
CALIBRATION AND BLANK		
Calibration type	multiple	multiple
Number of calibrators	-	-
Calibration curve	-	-
OPTIONS		
Blank absorbance limit	0.150	0.150
Kinetic blank limit	-	-
Linearity limit	150	150
Substrate depletion	-	-





ranges

_

PROTEIN (TOTAL)

Within-laboratory (CV)

1.1 %

1.5 %

A15

BIURFT