## **40814, 40818**

Anticoagulant Citrate Dextrose Solution Ph Eur (ACD) Solution A

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## Description

---

The solution is presented as a single-use container with a 500 mL or 750 mL volume in a clear polyolefin bag covered by a laminated polypropylene overwrap.

Each 100 mL contains:

- Glucosum Monohydricum Ph Eur - 2.45 g
- Natrii Citras Dihydricum Ph Eur - 2.2 g
- Acidum Citricum Monohydricum Ph Eur - 0.8 g
- Aqua ad Iniectabilia Ph Eur

## Storage Recommendations

---

Store up to 25 °C. Avoid excessive heat and protect from freezing.

The expiration date of the solution is printed on the bag and the product should not be used after this date.

## Indications

---

The Anticoagulant Citrate Dextrose Solution A is intended for use ONLY with automated apheresis devices for the collection of human blood components.

Please see the apheresis device operator's manual for a complete listing of warnings, cautions and undesirable effects.

## Warnings

---

1. Patients or donors with impaired or abnormal citrate and/or calcium metabolism (e.g. liver and renal diseases) can present an increased risk of citrate sensitivity. The attending physician should assess the appropriateness of such patients or donors for apheresis and prescribe how they should be monitored during the apheresis procedure.
2. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
  - Fluid leaks
  - Parts that are warped or deformed
  - Plastics that are brittle and discolored
  - Filters that have reduced filtration capabilities
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

## Caution

---

1. This product is sterilized with steam.
2. When you break frangible connectors, bend them in both directions to ensure that you break them completely. Failure to do so may result in restricted flow.

## Procedures for Use

---

See the apheresis device operator's manual for additional information.

- Not for direct intravenous infusion.
- Do not use unless the solution is clear and the container is intact.
- Discard unused portion.

## Symbols and Certifications

---

| Symbol                                                                              | Definition                                                                                                                  |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
|    | Indicates the temperature limitations for product storage when accompanied by temperature(s).                               |
|    | Indicates that the product is sterilized using steam or dry heat.                                                           |
|    | Indicates that the product should not be reused.                                                                            |
|    | Indicates the product was manufactured in accordance with Annex II of the European Council Directive 93/42/EEC, as amended. |
|   | Indicates the lot of the product when accompanied by a number.                                                              |
|  | Indicates the expiration date of the product when accompanied by a specific date.                                           |

## Return of Used Product

---

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department.

**IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.**

**Please contact your local representative for information regarding returned goods and product complaints.**

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©2019 Terumo BCT, Inc.