Medtronic


C315H20, C315H40, C315HIS, C315J, C315S4, C315S5, C315S10

Delivery catheter
Cathéter d'introduction
Einführungskatheter
Catéter introductor
Plaatsingskatheter
Catetere di posizionamento
Innføringskateter
Införingskateter
Καθετήρας τοποθέτησης
Fremføringskateter
Cateter introdutor
Vezetőkatéter
Cewnik wprowadzający
Zaváděcí katétr
Доставочный катетер
Απλικαčný katéter
Taşıma kateteri

Instructions for use ▪ Mode d'emploi ▪ Gebrauchsanweisung
▪ Instrucciones de uso ▪ Gebruiksaanwijzing ▪ Istruzioni per
l'uso ▪ Bruksanvisning ▪ Bruksanvisning ▪ Οδηγίες χρήσης
▪ Brugsanvisning ▪ Instruções de utilização ▪ Használati
útmutató ▪ Instrukcja użytkowania ▪ Návod k použití
▪ Инструкция по эксплуатации ▪ Pokyny na používanie
▪ Kullanım talimatları

! USA **Caution:** Federal Law (USA) restricts this
device to sale by or on the order of a
physician.

CE
0050
2010

The following are trademarks of Medtronic:
 Medtronic

Explanation of symbols on package labeling

Refer to the package label to see which symbols apply to this product

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive AIMD 90/385/EEC.

Date of manufacture

Manufacturer

Use by

Lot number

Sterilized using ethylene oxide

Do not reuse

Maximum storage temperature

Open here/Peel here

Catalogue number

Do not use if package is damaged

For US audiences only

Consult instructions for use

Outer diameter

Manufactured in


EC authorized representative

Inner diameter

Package contents

Guide-catheter dilator

Product documentation

Contents of Package

This package contains product documentation, 1 vessel dilator, and 1 of the following C315 catheters:

C315 H20 & H40

C315 HIS

C315 J

C315 S4



C315 S5

C315 S10

1 Device description

The C315 is a single-use percutaneous catheter intended to introduce various types of pacing or defibrillator leads and catheters. The C315 is packaged with a vessel dilator for introduction into the vasculature. Proximally, the C315 is equipped with a hemostatic valve, and the distal tip is radiopaque to facilitate imaging under fluoroscopy. The C315 is designed to be slittable, thereby allowing its removal after device placement. A variety of curves are available to accommodate various anatomies and different lead locations.

1.1 Contents of package

See page 3 for package contents.

2 Indications for use

The C315 is indicated for the introduction of various types of pacing or defibrillator leads and catheters.

3 Contraindications

Use of the C315 percutaneous catheter is contraindicated in patients with obstructed or inadequate vasculature.

4 Warnings and precautions

Necessary hospital equipment – Keep external defibrillation equipment nearby for immediate use during insertion, placement, acute lead system testing, or whenever arrhythmias are possible or intentionally induced. Backup pacing should be readily available during implant. Use of a delivery system and/or leads may cause heart block. For further information about this and other potential adverse events or complications, refer to the technical manual packaged with the appropriate product.

Vessel and tissue damage – Use care when passing a catheter through vessels and tissue.

- Avoid damage, such as perforations and dissections, to vessels and cardiac tissue, during catheter passage and positioning.
- Do not push, pull, or rotate the catheter against resistance. If resistance is met, discontinue movement, determine the reason for resistance, and take appropriate action before continuing. Damage to the catheter may prevent it from performing with accurate torque response and control, and may cause vessel damage.

Device compatibility – Use the catheter only with compatible transvenous devices. No test data is available to demonstrate compatibility of any non-Medtronic device with the catheter.

Consequences of using the catheter with incompatible devices may include the inability to deliver the transvenous device or damage to the transvenous device or catheter during delivery.

Slitting evaluations on this product were conducted using Medtronic Universal II or Medtronic Adjustable Slitters. These are the only slitters recommended for use with the C315 catheter.

Inspecting the sterile packages – Inspect the packages prior to opening.

- Contact your local Medtronic representative and provide the lot number from the outer package label if the seal or package is damaged or if there are any concerns regarding product integrity.
- Do not use the product after its expiration date.
- The recommended storage temperature for this product is at or below 40 °C (104 °F).
- The catheter set has been sterilized with ethylene oxide prior to shipment.

For single use only – This device is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize this device for purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Handling the catheter – Handle the catheter with care at all times.

- Do not kink, stretch, or severely bend the catheter.
- Do not use surgical instruments to grasp the catheter.
- Do not use excessive force when inserting a catheter into a vessel.
- Ensure the catheter is thoroughly flushed and free of air prior to use.
- Avoid contact with liquids other than blood, saline, or contrast solution.
- Use in conjunction with fluoroscopic guidance and proper anti-coagulation agents

Flushing the catheter – Use the side port on the hub to flush the catheter. The catheter must be thoroughly flushed and free of air prior to use. When using the side port, consider sealing the proximal opening of the hemostatic valve with a thumb or forefinger.

5 Potential adverse events

- air embolism
- allergic reaction to contrast media
- arteriovenous fistula formation
- bleeding at the insertion site
- brachial plexus injury
- cardiac tamponade
- dislodgement
- dissection
- endocarditis
- heart block
- hematoma formation
- hemothorax
- infection
- irregular heart beat
- mediastinal widening
- perforation
- pneumothorax
- subclavian artery puncture
- thrombophlebitis
- thrombosis
- valve damage
- vascular occlusion
- vessel damage

6 Directions for use

6.1 Device preparation

1. Flush the C315 with saline.
2. Place the dilator in the C315.
3. Flush the dilator with saline.

6.2 Recommended procedure

1. Obtain venous access
2. Advance the C315 and dilator together over a guidewire into the right atrium using fluoroscopy.
3. Remove the dilator.
4. Position the C315 catheter in the desired location.
5. Remove the guidewire.
6. Advance a lead or catheter through the C315 to the desired location.
7. When a lead has been positioned in the desired location, remove the C315 by slitting.
8. For instructions on removing the C315 by slitting, refer to the technical manual packaged with the slitter.

7 Detailed device description

7.1 Specifications (nominal)

Parameter	Models C315
Guide-catheter dilator	Material: Polyethylene
	Inner diameter: 0.96 mm (2.8 Fr, 0.038 in) minimum
	Outer diameter: 1.75 mm (5.2 Fr, 0.068 in) maximum
Guide catheter	Material: Polyether block amide,
	Inner diameter: 1.8 mm (5.4 Fr, 0.072 in) minimum
	Outer diameter: 2.4 mm (7.0 Fr, 0.092 in) maximum
Compatible Leads	1.4 mm (4.2 Fr, 0.055 in) or smaller
Compatible Inner Transvenous Devices	1.75 mm (5.2 Fr, 0.068 in) or smaller


Model number	H20	S4, S5, J	S10, H40	HIS
Guide catheter	20 cm	30 cm	40 cm	43 cm
usable length	(7.9 in)	(11.8 in)	(15.7 in)	(16.9 in)
Guide-catheter dilator	28.5 cm	38.5 cm	51.5 cm	51.5 cm
usable length	(11.2 in)	(15.2 in)	(20.3 in)	(20.3 in)

8 Medtronic disclaimer of warranty

For complete disclaimer of warranty information, see the accompanying disclaimer of warranty document.

9 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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© Medtronic, Inc. 2010
Medtronic B.V.
Medtronic Australasia Pty Ltd