INgezim® ADV gE PLUS

Prod Ref: 11.GEP.K3

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KIT COMPOSITION

REAGENT	2 Plates		5 Plates		25 Plates	
	Uni.	Vol.	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5		25	
Vials with Positive Control	1	1,5mL	2	1,5mL	2	1,5mL
Vials with Negative Control	1	1,5mL	2	1,5mL	2	1,5mL
Vials with Peroxidase Conjugate, ready to use	1	25mL	1	60mL	3	100mL
Bottles with Washing Solution 25x concentrated	1	100mL	1	100mL	2	250mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL	5	60mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL	2	250mL
Bottles with Diluent (DE31)	1	60mL	1	60mL	1	250mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (450 nm filter)

TECHNICAL BASIS

The gE test uses a monoclonal antibody (Mab) specific for gE glycoprotein of ADV and microtiter plates coated with ADV antigen (gE). After sample addition, if antibodies are absent in the test serum, the gE antigen will remain free and the Mab peroxidase-conjugated specific for gE, will bind. After substrate addition a colour reaction is produced. When gE antibodies are present in the test sample they will bind to the gE antigen of the well and will block the binding of the conjugated Mab. In this case no colour reaction will be produced.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.
Control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them in aliquots at -20°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

The sera sample needs to be at 1/2 dilution to be tested. The dilution of sera samples could be done directly on the wells by adding 50 µL of serum and 50 µL of sera diluent to each well.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (40 mL of concentrate washing solution into 960 mL of distilled water). This solution remains stable when is stored between +2°C and +8°C.

Positive and Negative Controls:

Controls need to be at 1/2 dilution to be tested. The dilution of control samples could be made directly on the wells by adding 50 µL of serum and 50 µL of sera diluent to each well.

TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add 50 µL of diluent and 50 µL serum (controls and samples) to each well. Shake carefully to homogenate the mixture, seal the plate and incubate for **2 hours at room temperature (20°C to 25°C)** (preferably with agitation) or overnight between **+2°C and +8°C**.
3. Wash 3 times following the described procedure.
4. Add 100 µL of conjugate to each well. Seal the plate and incubate for **30 minutes at room temperature (20°C to 25 °C)**.
5. Wash 4 times following the procedure.
6. Add 100 µL of substrate, to each well. Keep the plate for **10-15 min at room temperature in a dark place**.
7. Add 100 µL of stop solution to each well.
8. Read the OD at 450 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

TEST VALIDATION:

$$\frac{P}{N} < 0.3$$

- P= OD of positive control
- N= OD of negative control

Negative control's OD must be higher than 0.800

CUT OFF VALUES CALCULATION:

- Positive Cut off = $N - [(N - P) \times 0.4]$
- Negative Cut off = $N - [(N - P) \times 0.35]$
- N = OD of Negative control (average)
- P = OD of positive control (average).

RESULTS INTERPRETATION:

- If you are running samples in duplicate wells, Optical Density (OD) value must be calculated as the arithmetic average of the two wells.
- Sample is positive when its OD value is lower than or equal to the Positive cut off value.
- Sample is negative when its OD value is higher than or equal to the negative cut off value.
- Samples with OD values between both cut off values must be scored as doubtful. In these cases, we recommend to run a second sample of the same animal collected 2 weeks later.

INHIBITION PERCENTAJE CALCULATION (IPM)

$$IPM = \frac{OD\ CN - OD\ M}{OD\ CN - OD\ CP} \times 100$$

- IP ≥ 40: Sample positive to gE antibodies.
- IP ≤ 35: Sample negative to gE antibodies.
- 40 > IP > 35: Doubtful sample to gE antibodies.

PRACTICAL EXAMPLE:

Obtained OD average values:

- NC = 1.160 (N); PC = 0.110 (P)
- Sample 1 = 0.373; Sample 2 = 1.395

Assay validation:

- P/N ratio = $0.110 / 1.160 < 0.3$
- $1.160 > 0.800$

CUT OFF Calculation:

- Cut Off (+) = $1.160 - \{(1.160 - 0.110) \times 0.4\} = 0.740$
- Cut Off (-) = $1.160 - \{(1.160 - 0.110) \times 0.35\} = 0.792$

Results:

- Sample 1 = Positive for gE
- Sample 2 = Negative for gE

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

CERTIFIED

MANAGEMENT SYSTEM
QUALITY - ISO 9001
ENVIRONMENTAL - ISO 14001

9191.INGE

9175.INGE

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Distributed in

por

by