
ACCUVAC Pro suction device

Experience at its best – for Professional Demands

ACCUVAC Pro

Experience at its best – for professional demands

The first life-saving step involves clearing the patient's airways from obstructions so that effective ventilation can be administered. In an emergency situation the patient's survival depends on keeping the airways clear. ACCUVAC Pro ensures free and clear airways by suctioning the mouth and throat, endotracheal or bronchial areas. The pump's high suction capacity and its easy handling assist medical personnel in carrying out this effective and proven patient procedure in professional emergency medical services.

Broad application range, flexible operation

With the choice of four pre-defined suction levels, the user can adjust the suction capacity to every situation. ACCUVAC Pro can be used on babies, children and adults. In addition to suctioning airways, ACCUVAC Pro also deflates vacuum splints and mattresses.

Your benefits at a glance

- Variable uses
- Safe and hygienic
- Simple operation, robust
- EN 1789-compliant, air worthy as per RTCA DO-160G

Convenient operation, broad application spectrum

- High suction capacity of about 34 liters/minute at -0.8 bar (at device inlet)
- Rechargeable battery: user can easily change battery without tools
- Modern lithium-ion rechargeable battery with greater run time of more than 60 minutes
- Four pre-defined vacuum levels can be selected by user
- Automatic function check with rapid visual and acoustic feedback
- Safe secretion collection in auto-clavable reusable secretions canister with disposable bacteria filter and overflow protection or disposable Serres® suction bag with integrated bacteria filter
- Robust housing made of impact-resistant materials
- Non-tip housing with low center of gravity

Ideal operation in everyday medical services

- Optimum accommodation of suction tube in holder on side of device
- Single-handed release from the wall mounting with the press of a button
- Compatible with existing ACCUVAC wall mounting
- Optionally available: large accessories bag, protective bag and/or shoulder strap

For every patient and every emergency

- Suctioning of secretions and food particles
- Endotracheal and bronchial suction
- Particularly suitable for pre-hospital care
- Deflate vacuum splints and mattresses
- Suitable for use on babies, children and adults

Safety and reliability day after day

In the event of a malfunction, you can read out the service files from your ACCUVAC Pro and send them to WEINMANN Emergency. At best, the data are enough for our specialists to


 **Service Hotline**
+49 40 88 18 96 122

provide a remote diagnosis and work with you to resolve the problem. If not, we will take a closer look at your device and, as needed, we will give you a replacement device to use while yours is out of service.

Service data: ACCUVAC Pro

Manufacturer's warranty	2 years
Safety check interval	no check required
Maintenance schedule	no maintenance required
Function check	automatic
Time required for function check	about 15 seconds
Battery change	simple, from outside device, no tools required
Rechargeable battery	no maintenance, no calibration required
Deep discharge protection for battery	✓
Protection from electrical surge and reverse polarity	✓
Conversion of disposable system/reusable system	simple, without tools
Software update can be made by operator/user	✓
Remote diagnosis of device malfunction	✓

Intuitive Operation Ensures Safety

1. Holder for suction tube
for ideal storage

2. Four pre-defined suction levels
quick push-button selection

3. Backlit
to make suction levels visible in poor lighting conditions

4. Release for wall mounting
single-handed release with press of button

5. Battery compartment cover
can be opened without tools for fast and easy access to battery

6. Rechargeable battery
modern lithium-ion battery with longer service life and rapid charging

7. Device base with tube guide

8. Automatic function check
with quick visual and acoustic feedback

9. Canister holder
impact-resistant materials increase safety for the canister systems

10. Secretions canister Serres® 1000 ml
with disposable suction bag, integrated bacteria filter and solidifying agents. When the bag is full, ACCUVAC Pro stops automatically

Technical Data

ACCUVAC Pro

Product class as per Directive 93/42/EEC	IIa	
RoHS-conformity as per Directive 2011/65/EU (RoHS II)	Yes	
Classification as per EN ISO 10079-1	High vacuum/high flow	
Standards used	EN 60601-1, EN 60601-1-11, EN 60601-1-12, EN ISO 10079-1, EN 1789, RTCA DO-160 G	
Dimensions (W x H x D)	• with canister system • with canister system + accessories bag	370 x 277 x 146 mm 370 x 277 x 152 mm
Weight	• Device with battery, no canister system and holder • Reusable canister system with holder • Disposable canister system with holder	3.65 kg 1.00 kg 0.65 kg
Temperature range: Operation Transport/storage	-5 °C to +50 °C -40 °C to +70 °C	
Maximum current consumption	3.8 A	
Rated voltage	12 V DC nominal (min. 10 V, max. 15 V)	
Suction capacity at device inlet (without canister system) at -0.8 bar, fully charged battery and 21 °C/1013 hPa	34 liters/min	
Suction capacity at device inlet to reusable canister system at -0.8 bar, fully charged battery and 21 °C/1013 hPa	30 liters/min	
Vacuum settings	Pre-defined levels: -0.1 bar, -0.2 bar, -0.5 bar and -0.8 bar, electronically regulated	
Type of operation	S2 60 min	
Degree of protection	IP34D	
Type of rechargeable battery	L-ion	
Rechargeable battery	• Charging time • Service life • Battery operating time	Battery status 80%: 2 hrs, 40 min (at 20 °C while not in use) Battery status 100 %: about 4 hrs about 500 charging cycles 60 min (during continuous operation with fully charged/new battery)
Volume of secretions canister	1000 ml	
Reusable suction tube	10 mm ID, 1300 mm long	
Disposable suction tube	7 mm ID, 1800 mm long	

CE 0124

Certified Quality Management System meeting EC directive 93/42/EEC, Annex II (EN ISO 9001:2008/EN ISO 13485)

All rights to design and specification modifications reserved.

Designed in Germany

Select your preferred ACCUVAC Pro combination – designed to meet your needs:

ACCUVAC Pro with disposable canister system, WM 11605

ACCUVAC Pro with reusable canister system* WM 11600

ACCUVAC Pro with disposable canister system with accessories bag, WM 11645

ACCUVAC Pro with reusable canister system* with accessories bag, WM 11640

Combine your ACCUVAC Pro with additional accessories from page 6

*Please note that the reusable canister system is equipped with a disposable bacteria filter.

Accessories and Replacement Parts




Accessories

- | | |
|--|----------|
| 1. Protective bag
(cannot be combined with accessories bag) | WM 11692 |
| 2. Shoulder strap | WM 11693 |
| 3. 12-V connection cable | WM 10650 |
| 4. Power supply unit und charger
for 100V to 240V alternating current to ACCUVAC
device plug, including power cable in compliance
with EN 50075 (Europlug) | WM 2620 |
| 5. Wall mounting for ACCUVAC
including installation set | WM 15208 |
| 6. Holding plate for standard hospital rails | WM 15845 |
| 7. Attachment set with
two adapters for standard hospital rails | WM 15805 |
| 8. Wall mounting for power supply unit
and charger | WM 15844 |
| 9. Installation set pole bracket
dia. 19 mm - 40 mm | WM 15806 |
| 10. Conversion set for reusable canister system,
consisting of:
- Set, reusable canister system (WM 17821)
- Holder for reusable canister system (WM 11654) | WM 17820 |
| 11. Conversion set for disposable canister system,
consisting of:
- Holder for disposable canister system (WM 11754)
- Vacuum tube for Serres® secretions
canister (WM 11761)
- Serres® secretions canister, complete (WM 10790) | WM 17825 |
| 12. Conversion kit for accessories bag,
consisting of:
- Accessories bag (WM 11691)
- Battery compartment cover for accessories bag (WM 11614) | WM 17829 |

Consumables

- | | |
|--|----------|
| 13. Serres® suction bag 1000 ml, with
hydrophobic filter and solidifying agent
(Packaging unit: 32 items in carton) | WM 17800 |
| 14. Disposable suction tube with fingertip
control for use with disposable system,
180 cm long, 7 mm ID, also available in
sets of 10, 20, 32 and 50 tubes | WM 10778 |
| 15. Reusable suction tube 130 cm long,
10 mm ID, also available in sets of
10, 20, and 50 tubes | WM 10662 |
| 16. Fingertip control for reusable suction tube
10 mm ID, also available in sets of
10, 20 and 50 pieces | WM 10666 |
| 17. Set of 10 disposable bacteria filters
for reusable secretions canister | WM 17830 |

Replacement Parts for Reusable System

- | | |
|--|----------|
| 1. Upper part of secretions canister cover | WM 11657 |
| 2. Filter holder | WM 11661 |
| 3. O-Ring for filter holder | WM 11663 |
| 4. Disposable bacteria filter | |
| 5. Lower part of secretions canister cover | WM 11658 |
| 6. Ball float valve | WM 11662 |
| 7. Reusable secretions canister 1000 ml | WM 11653 |
| 8. Reusable suction tube 10 mm ID | WM 10662 |
| 9. Fingertip control for reusable suction tube 10 mm ID | WM 10666 |
| 10. Holder for reusable canister system | WM 11654 |
| • Secretions canister cover, complete, consisting of replacement parts 1 – 6 | WM 17822 |
| • Set, reusable canister system consisting of replacement parts 1 – 9 | WM 17821 |
| • Conversion set for reusable canister system consisting of replacement parts 1 – 10 | WM 17820 |

Replacement Parts for Serres® Disposable System

- Serres® secretions canister, complete, consisting of: WM 10790
 - Secretions canister, 1000 ml (WM 10775)
 - Suction bag
 - Disposable suction tube with fingertip control (WM 10778)
- Serres® secretions canister, 1000 ml WM 10775
- Vacuum tube for Serres® secretions canister WM 11761
- Holder for disposable canister system WM 11754
- Set disposable canister system, consisting of: WM 17826
 - Vacuum tube for Serres® secretions canister (WM 11761)
 - Serres® secretions canister, complete (WM 10790)
- Angled connector for Serres® canister, reusable WM 10799

Replacement Parts for ACCUVAC Pro

- Set, catch release WM 17827
- Lithium-ion rechargeable battery WM 11603
- Battery compartment cover (without eyelets for accessories bag) WM 11604
- Holder for suction tube WM 11664
- Device base WM 11677

Simply Professional

WEINMANN Emergency is a family-owned, internationally active medical technology company. With our mobile system solutions for emergency, transport and disaster medicine, we set standards for saving human lives. In close collaboration with professional users in emergency medical services, hospitals and military medical corps, we develop innovative medical products for ventilation and defibrillation. For more than 100 years we have offered our customers a high degree of reliability, extensive experience and quality made in Germany.

Headquarter

WEINMANN Emergency
Medical Technology GmbH + Co. KG
Frohbösestraße 12
22525 Hamburg
Germany

T: +49 40 88 18 96-0
F: +49 40 88 18 96-480
T: +49 40 88 18 96-120 Customer Service
T: +49 40 88 18 96-122 After-Sales Service
E: info@weinmann-emt.de

Center for Production, Logistics, Service

WEINMANN Emergency
Medical Technology GmbH + Co. KG
Siebenstücken 14
24558 Henstedt-Ulzburg
Germany

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Conformity assessment procedure: Directive 93/42/EEC Annex II on medical products, passed by the commission on 14th June 1993, last amended on 5th September 2007

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Lenzkirch, 01.04.2020
Place and date of issue

i.V. Andreas Heer
Head of Quality Management

i.V. Steffi Focke
Quality Management

EC Certificate
Directive 93/42/EEC Annex II, excluding Section 4
Full Quality Assurance System
Medical Devices

Registration No.: HD 60148646 0001

Report No.: 21201240 016

Manufacturer: WEINMANN Emergency Medical
Technology GmbH + Co. KG
Frohbösestr. 12
22525 Hamburg
Deutschland

Products: Medical devices for emergency and transport medicine
(see attachment for products and sites included)
Replaces Certificate, Registration No.: HD 60129559 0001

Expiry Date: 2024-05-26

The Notified Body hereby declares that the requirements of Annex II, excluding section 4 of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex II, section 5 of the aforementioned directive. For placing on the market of class III devices covered by this certificate an EC design-examination certificate according to Annex II, section 4 is required.

Effective Date: 2020-04-16

Date: 2020-04-16

Notified Body

Dipl.-Ing. I. Munkler

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

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
Products included:

Emergency and transport medicine:

- Suction pumps
- Ventilators
- Module systems
- Resuscitators
- Ventilation masks
- Patient hose systems
- Pressure reducers
- Click dial flowmeter
- Defibrillators-/Monitoring systems
- SpO2 sensors

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
including the locations according to annex

Scope: Design and development, production, distribution
and servicing of medical devices and digital services
for emergency and transport medicine

Proof has been furnished by means of an audit that the
requirements of ISO 9001:2015 are met.

Validity: The certificate is valid from 2022-07-01 until 2025-06-30.

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WEINMANN Emergency Medical Technology GmbH + Co. KG
TÜV Rheinland Cert GmbH