Instruction for use

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«Mycoplasma IgG-ИФА»

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUALITATIVE DETERMINATION
OF IgG ANTIBODIES TO MYCOPLASMA
IN HUMAN SERUM OR PLASMA**

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
Authorized Representative in EU:
Polmed.de
Steinacker 20, D-73773
Aichwald, Germany
e-mail: info@polmed.de

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CONTENT

1. INTENDED USE	9
2. SUMMARY AND EXPLANATION	9
3. PRINCIPLE OF THE TEST	9
4. WARNINGS AND PRECAUTIONS	10
5. KIT COMPONENTS	11
6. SPECIMEN COLLECTION AND STORAGE	12
7. TEST PROCEDURE	12
8. QUALITY CONTROL	14
9. CALCULATION OF RESULTS	14
10. EXPECTED VALUES	14
11. LITERATURE	15

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IgG ANTIBODIES TO MYCOPLASMA IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the qualitative determination of IgG antibodies to Mycoplasma in blood serum or plasma.

This kit is designed for measurement of IgG antibodies to Mycoplasma in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 45 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Mycoplasma represent a separate class of microorganisms. Unique metabolic properties of Mycoplasmae determine their poor growth on standard microbiological media and require the application of serological methods in diagnostics.

Among large variety of species, *M.hominis*, *M.genitalium*, *M.pneumoniae* and closely immunologically related *Ureaplasma urealyticum*, play the most considerable role in medical practice. All these microorganism share common antigenic epitopes.

M.pneumoniae causes pneumonia, bronchitis and bullous meningitis; other mycoplasmae can cause acute or chronic pelvic inflammations and may contribute to male and female infertility.

3. PRINCIPLE OF THE TEST

This test is based on indirect enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific Ig, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0 % H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5. KIT COMPONENTS 5.1. Contents of the Kit

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	Mycoplasma IgG EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL	Calibrator set, 1 ml	1	pcs	blue	2 months
3 CONTROL-CONTROL+	Control sera (1 ml and 1 ml, resp.)	2	pcs	colourless and red	2 months
4 CONJ HRP	Conjugate, 14 ml	1	pcs	red	until exp.date
5 DIL	EIA buffer 14 ml	1	pcs	blue	until exp.date
6 SUBS TMB	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7 BUF WASH 26X	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
8 STOP	Stop solution, 14 ml	1	pcs	colourless	until exp.date
9 N003	Plate sealing tape	2	pcs		N/A
10 K106I	Instruction Mycoplasma IgG EIA	1	pcs		N/A
11 K106Q	QC data sheet Mycoplasma IgG EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 80–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 20–250 µl;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 6 wells for the calibrators CAL and control samples CONTROL - , CONTROL + and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 80 µl of buffer DIL into the wells (EIA buffer).
3	Pipet 20 µl of calibrators, control samples CONTROL - , CONTROL + and 20 µl of unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit). After add of samples (without heparin or EDTA), buffer change colour from blue to red. Do not dilute control sample and calibrators.
4	Incubate 30 minutes at +18...+25 °C.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +18...+25 °C.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10–20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.
14	Apply coefficient method for data reduction. 1. Calculate a cutoff value by multiplying mean OD of Calibrator by the Q value (see Lot QC insert). Cutoff = OD(Calibrator) * Q value 2. Divide mean OD of each sample by Cutoff. K = OD(sample) / Cutoff 3. If the K value is greater than 1.1, the result is POSITIVE. If the K value is less than 0.9, the result is NEGATIVE. If the K value is between 0.9 and 1.1, the result is EQUIVOCAL.

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS

Some laboratories, based on their population studies, set up a second cutoff, which stands between anamnestic ('normal') IgG antibody level and 'high' IgG antibody level characteristic of reactivation or late period of primary infection. Recommended values for this second cut-off for two age groups (8 months – 3 year, > 3 years) are presented in the table below.

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Mycoplasma IgG. Based on data obtained by XEMA, the following normal range is recommended (see below).

Sex, age	Units, K	
	Lower limit	Upper limit
Seronegative	<0.1	0.9
Seropositive > 3 years	1.1	3.9
newborn*	<0.1	1.3
under 8 months*	<0.1	2.6
8 months – 3 years	<0.1	2.1

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