<sup>(BIC)</sup> **Sphera Duo Hydrocephalus Valve<sup>1</sup>**  
**Sphera Duo Hydrocephalus Shunt System<sup>2</sup>**

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---

# Sphera Duo Hydrocephalus Shunt System

## Sphera Duo Hydrocephalus Valve


### DESCRIPTION

The Sphera Duo Valve is made of transparent medical grade silicone and transparent polysulfone input and output connectors encasing the pressure control mechanism, comprised by a synthetic ruby ball and conic seat and 316L stainless steel spring (doesn't interfere in CT scan or Magnetic Resonance imaging). The valve is available in sizes Adult, Infantile and Neonatal. All sizes are supplied in four pressure ranges: High, Medium, Low and Extra-Low.


The adult and infantile sizes design presents low profile of implant and flexible body with anatomical shape to contour the curvature of the cranium. They have pumpable central reservoir with needle guard that protects against perforations on puncturing or sampling procedures.

The Neonatal size has no pumping chamber, as the design prioritizes low profile and minimum implant volume. The system can be combined to a reservoir connected to the ventricular catheter to enable puncturing and CSF sampling.

The valves are individually multitestated for flow rate and pressure, assuring its quality and functionality. The shunt system is sterile and non-pyrogenic if the packing integrity is maintained, and during the validity term of the sterilization. It is packed in double blister and sterilized by ethylene oxide. The external blister, after removal of the sealing foil allows access to the internal blister and the sterile surgical field. The later, after removal of the sealing foil allows access to the components.


## NEONATAL


## SIZES


## ADULT

## INFANTILE

## NEONATAL

The arrow printed on the valve body indicates the flow direction, and the circular markings symbolize the pressure according to the legend below. The indication markings of flow and pressure are radiopaque and detectable by X-ray image.

### LEGEND


| Pressure  | Marks   | Values (Flow of 20 ml/h)        |
|-----------|---------|---------------------------------|
| Extra Low | ○ ○ ○ ➤ | 1 to 3 cm of H <sub>2</sub> O   |
| Low       | ● ○ ○ ➤ | 3 to 7 cm of H <sub>2</sub> O   |
| Medium    | ● ● ○ ➤ | 7 to 11 cm of H <sub>2</sub> O  |
| High      | ● ● ● ➤ | 11 to 14 cm of H <sub>2</sub> O |

The nominal pressure showed on the package is determined by “in vitro” standard test with continuous flow of 20ml/hour and 4cm-water pulse with 60-cycle/minute amplitude and frequency. Each valve is individually tested to assure conformity to the pressure/flow rate features defined for each rate.

The bands in the graphics represent the range where the pressure curves operate when submitted to a progressive increasing flow from 5 to 50 ml/h. They reflect the “in vitro” results obtained and are referential values for the surgeon’s option to attend to the clinical conditions of the patient.

The pressure/flow performance is described on Graphic:

### PRESSURE / FLOW

The graphic represents average rates. Consider a range of  $\pm 1,5$  cm H<sub>2</sub>O.

Neither the valve nor any item of the shunt system has any component which can cause interferences in CT scan or Magnetic Resonance imaging, and the shunt flow characteristics are not affected by these exams.

The valve should be implanted with complete shunt system, including ventricular catheter and peritoneal catheter. The models SDVAE, SDVAB, SDVAM, SDVAA, SDVINE, SDVINB, SDVINM, SDVINA, SDVNNNE, SDVNNB, SDVNNM, SDVNNNA, SDVNNER, SDVNNBR, SVNNMR and SDVNNAR don't present these components, and they should be provided separately.

In addition to the valve the Shunt System contains: ventricular catheter with stainless steel guide and angle former, peritoneal catheter and tubing for irrigation with adapter.

The ventricular cerebral catheter is made with medical grade silicone with radiopaque fillet (allowing visualization through x-ray) always accompanied by the stainless steel stylet for introduction. The catheters have marks at every 1 cm to direct introduction.

The peritoneal catheter is made with medical grade silicone with radiopaque fillet. The catheter shows marks to direct introduction at every 5 cm. The catheter have opened distal tip. Other drainage components such as connectors and reservoirs may be separately purchased and connected to the system. All connections are universal.

### INDICATION

This product is indicated for intracranial pressure control in post-surgical drainage and Hydrocephalus with ventricle-peritoneal or ventricle-atrial derivation.

*NOTE: The atrial catheter is not provided and should be acquired separated.*

*NOTE: In the models SDVAE, SDVAB, SDVAM, SDVAA, SDVINE, SDVINB, SDVINM, SDVINA, SDVNNE, SDVNNB, SDVNNM, SDVNNA, SDVNNER, SDVNNBR, SVNNMR and SDVNNAR the proximal (ventricular) and distal (peritoneal or atrial) catheters are not provided and should be acquired separated. Catheters from different manufacturers can be used if they have an inner diameter from 1.2 to 1.6 mm for use with Adult and Infantile valve models, and inner diameter from 1.0 to 1.4 to be used with Neonatal valve.*

### CONTRAINDICATION

The Shunt System implant is contraindicated in patients with knowledge or suspicion of infection on the implant region of anyone of the components (meningitis, ventriculitis, skin infections, bacteremia, septicemia, peritonitis, etc) or infection anywhere in the body. The implant of the ventricle-atrial derivation system is contraindicated in cases of pulmonary embolism or pressure increase on the right atrium, and congenital cardiopathies.

### WARNINGS AND PRECAUTIONS

The implant of a shunt in cases with hemorrhagic liquor or high protein content (above 200 mg/dL) must be carefully assessed due to increased risk of obstruction. Avoid excessive pumping of the valve after the implant due to the risk of rapidly draining the ventricles, which may result in ventricular collapse or subdural hematoma. Avoid pumping and unnecessary tests. Frequent pumping may cause excessive ventricle drainage, especially small-volume ones.

We recommend that the Neonatal model for low weight patients and/or with chronic malnutrition, subject to difficulty of healing and scalp erosion at the site of the valve implant be carefully assessed, studying also the possibility of implanting in another place in the derivation system path.

The post-surgical of patients implanted with derivation system must be closely observed for detection and evaluation of obstruction symptoms or overdrainage.

The overdrainage may be discerned by the signs and symptoms of reduced intracranial pressure and to the development of a subdural hematoma, subdural hygroma, or collapse of the lateral ventricular wall.

The obstruction may lead to the signs and symptoms of increased intracranial pressure.

In children the signs may be increased tension of the fontanelle, congestion of the scalp veins, apathy, drowsiness, rigidity of the nape of the neck, irritability and vomiting.

In adult patients, the signs may be: vomiting, headaches, loss of conscience, rigidity of the nape of the neck, blurred vision. If some of these occur, an immediate surgical evaluation is needed.

The surgeon is the best person to select the procedure to be used.

The product implant is restricted to surgeon doctors with specific training and qualification in this procedure.

Do not use this device for other purposes rather than the intended and stipulated by the manufacturer.

This product must be handled by qualified personnel only.

Patients or the persons in charge of them must be informed of possible complications and warned about the precautions involving the implant of this product.

Handle the material under strict aseptic technique.

The silicone elastomer is very electrostatic and susceptible of contamination from air or surface particles. Avoid contact with surfaces releasing threads or particles. Handle the material with sterile and clean gloves, avoiding contact with frayed fabrics, talc or other contaminants. The safety of the connections can be reinforced with ligatures. Care must be taken to make the connections and binding in such a way that the catheters do not suffer cuts or occlusion. Care must be taken also not to kink catheters at any point, because there is risk of obstruction of the system.

Avoid contact with piercing or cutting materials when opening the package and handling the components of the Shunt System.

When using hypodermic needles for injection or puncturing of the reservoir, calibers greater than 0.5mm should not be used.

Perform the puncturing so as not to exceed a maximum of five procedures over the life cycle of the valve to prevent damage to the chamber.

Assure the absence of air and the correct introduction depth of the ventricular catheter through observation of liquor output on the proximal tip.

Excessive pumping of the valve must be avoided.

The maximum positive pressure that can be applied to the system without impairing its performance is 2m H<sub>2</sub>O. Negative pressures applied to the system can't impair its performance.

Do not introduce stylets through the valve connectors.

In the early and immediate postoperative period, the proper functioning of the valve can be checked by pumping and observing the valve response, which will indicate if there are proximal or distal obstruction. In the case of distal obstruction, the pumping chamber will offer great resistance to being pressed. In the case of proximal obstruction, the camera will not return to spherical shape after being pumped, staying down.

As an after implantation care, the patient should be instructed not to lay the head down on the site of valve implant.

## **TRANSPORTATION AND STORAGE**

Always keep the product in cool and dry place, free of contamination. It is advisable to inspect the products storage at regular intervals to assure the usage before the expiration of sterilization validity date.

## **REESTERILIZATION**

This is a single-use product. Do not resterilize this device.

## **DISPOSAL**

The product should be disposed as hospital waste, or according regional regulation.

If the sterilization validity of the product has expired, or in case of damaged product packaging, please return to the manufacturer or discard as convenient.

## **IMPLANT PROCEDURES**

The surgical procedure for component implants, including positioning of the ventricular catheter and valve over the cranium can be directed by different implant techniques, at the discretion of the surgeon.

The components assembly in the shunt system is made through the connection between proximal (ventricular) catheter, valve and distal (peritoneal) catheter.

The Unishunt models do not need connection of valve and peritoneal catheter, because they are supplied pre-assembled.

The Unitized models do not need connection of valve and catheters ventricular and peritoneal, because they are supplied pre-assembled.

## INSTRUCTIONS OF USE

### PRE-IMPLANT PROCEDURES

ENG

POR

SPA

**ATTENTION:** Pumping and / or irrigation tests are intended only to remove air bubbles and evaluate the permeability of the valve for implant safety, and doesn't indicate its functional pressure range. There is no suitable method to evaluate the functional pressure in pre-implantation procedure or during the implantation. Hpbio does not recommend testing to verify the functional range of valve, because this demands an extensive procedure, specific equipment and extreme caution to avoid contamination of the valve, hindering its execution in the surgical theater.

### VALVE

#### ADULT AND INFANTILE SIZES

It is advisable to perform the pumping or irrigation procedure to assure the air bubbles withdraw from the components and assure the absence of obstruction.

#### NEONATAL SIZE

The removal of air bubbles is performed by the Irrigation Procedure.

#### PUMPING PROCEDURE

1. Immerge the valve in sterile isotonic solution.
2. Pump slightly (do not apply excessive pressure/frequency), observing the flow output up to the elimination of air bubbles.

#### IRRIGATION PROCEDURE

1. Fill a 5 to 10ml syringe with sterile isotonic solution.
2. Adapt the syringe to the provided transparent test tube. Attach the other tip of the tube on the valve input (observe the flow direction trough the printed arrow).
3. Hold the valve in vertical position with the arrow pointing up, and slightly irrigate the valve with the solution, until it completely fills the valve, observing the outflow at the output connector.

#### PERITONEAL CATHETER

1. Fill a 5 to 10mL syringe with sterile isotonic solution.
2. Connect the syringe to the adapter of the irrigation tubing, and the tubing to one end of peritoneal catheter.
3. Inject the solution in the catheter until totally full, observing the outflow at the other end.

#### UNISHUNT SYSTEM (VALVE AND PERITONEAL CATHETER PRE-ASSEMBLED)

##### IRRIGATION PROCEDURE

1. Fill a 5 to 10mL syringe with sterile isotonic solution.
2. Attach the syringe to the irrigation tubing (provided) and the tubing to valve's input connector (observe the flow direction trough the printed arrow).
3. Hold the valve on vertical position, with the arrow pointing down and press the syringe plunger gently and slowly, observing the line water to fill the valve and eliminate air bubbles until the output of flow at the distal end tip of peritoneal catheter.

#### UNITIZED SYSTEM (VALVE PRE-ASSEMBLED TO VENTRICULAR AND PERITONEAL CATHETER)

##### IRRIGATION PROCEDURE

1. Attach the provided connector to a 5 to 10mL syringe filled with sterile isotonic solution.
2. Press the syringe plunger gently and slowly with the connector inserted through one of the ventricular catheter's orifices, until completely fill the system and eliminate air bubbles.

## COMPLICATIONS

The use of this device may cause complications similar to those associated with other surgical procedures, or to any foreign object implanted in the body, depending on the degree of reaction and sensitivity of the patient.

The most common complications when using shunt systems and their components are infection, obstruction, system disconnection, mechanical failure of components, overdrainage and insufficient drainage.

The use of the ventricular catheter may cause infection, subcutaneous infiltration of liquor, neurological sequela or intraventricular obstruction by particles (including hemorrhagic liquor or with high protein content).

The ventricular catheter can be obstructed by the adhesion to choroid plexus, if inserted in the brain or by the coaptation of the ventricular walls.

The ventricle-peritoneal derivation may cause inflammation of the abdominal cavity, skin erosion, intestinal perforation, cysts or pseudocysts, umbilical fistula and ascitis.

The Peritoneal Catheter is subject to obstruction by blood clots, debris, bacterial colonization, aggregated cells of tumor, and may suffer fracture.

The ventricle-atrial derivation may cause infection and predisposition to infection dissemination to other areas or organs, skin erosion, pulmonary embolism, thrombosis, cardiac perforation and cardiac tamponade, cardiac arrhythmia, pulmonary hypertension.

The presence of implanted devices can result in lesions originated by degradation of the tissue or epidermis around the implant area and may lead to development of infections. Infections can also be originated by contaminants presented in surgical incision site or by pathogenic agents presented in the bloodstream.

The squeezing of the catheters may cause system obstruction or drainage decrease.

Growth of the patient may cause displacement of the Peritoneal Catheter from the peritoneal cavity, deriving liquor to areas of difficult absorption.

The disconnection of the catheters or fracture of a component can cause migration of the components into the peritoneal cavity, lateral ventricle of the brain or adjacent areas.

The use of a hydrocephalus shunt system has been related to the incidence of epilepsy, and incidence of apoplexy in cases of constant revision of the catheter.

In children, the extreme reduction of the intracranial pressure will cause marked depression of the anterior fontanelle, overriding of the cranial bones and can transform a communicating into a non-communicating hydrocephalus.

Sudden postural changes can cause feeling of unease, discomfort and / or dizziness.

## IMPLANT TRACEABILITY

This derivation system is provided with a PATIENT IDENTIFICATION CARD. We recommend that it be filled out and that the patient or person in charge keep it. This card is useful for clinical follow-up and on emergencies.

Extra labels of product identification are also provided; they must be attached to the card, to the patient records, and to the surgical description card.

## COMMUNICATION FORM

Please fill the communication form provided with the packing, if any comment on the product is necessary before or after the implant, and send it directly to HpBio or through the supplier or representative. If possible, attach the part or component that generated the comment with serial number and lot.

| MODEL                           | CATALOGUE N°                                     | PRESSURE                          | CONTENT   |
|---------------------------------|--|-----------------------------------|---|
| <b>SPHERA<br/>DUO<br/>ADULT</b> | ADE10210D<br>ADB10210D<br>ADM10210D<br>ADA10210D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size<br>01 Right Angled Ventricular Cerebral Catheter 10 cm with stainless steel guide<br>01 Peritoneal Catheter 102 cm<br>01 irrigation tubing                    |
|                                 | ADE10215D<br>ADB10215D<br>ADM10215D<br>ADA10215D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 Peritoneal Catheter 102 cm<br>01 90° Angle Former, 01 irrigation tubing   |
|                                 | ADE10223D<br>ADB10223D<br>ADM10223D<br>ADA10223D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size<br>01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide<br>01 Peritoneal Catheter 102 cm<br>01 90° Angle Former, 01 irrigation tubing |
|                                 | ADE12015D<br>ADB12015D<br>ADM12015D<br>ADA12015D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 Peritoneal Catheter 120 cm<br>01 90° Angle Former, 01 irrigation tubing   |
|                                 | ADE12023D<br>ADB12023D<br>ADM12023D<br>ADA12023D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size<br>01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide<br>01 Peritoneal Catheter 120 cm<br>01 90° Angle Former, 01 irrigation tubing |

| MODEL  | CATALOGUE N°   | PRESSURE                          | CONTENT   |
|--|--|-----------------------------------|---|
| <b>SPHERA<br/>DUO<br/>ADULT<br/>UNISHUNT</b> | UADE10210D<br>UADB10210D<br>UADM10210D<br>UADA10210D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Right Angled Ventricular Cerebral Catheter 10 cm with stainless steel guide<br>01 irrigation tubing                    |
|  | UADE10215D<br>UADB10215D<br>UADM10215D<br>UADA10215D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing   |
|  | UADE10223D<br>UADB10223D<br>UADM10223D<br>UADA10223D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing |
|  | UADE12015D<br>UADB12015D<br>UADM12015D<br>UADA12015D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size pre-assembled to<br>01 Peritoneal Catheter 120 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing   |
|  | UADE12023D<br>UADB12023D<br>UADM12023D<br>UADA12023D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Adult size pre-assembled to<br>01 Peritoneal Catheter 120 cm;<br>01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing |

| MODEL                            | CATALOGUE N° | PRESSURE | CONTENT  |
|----------------------------------|--------------|----------|--|
| <b>SPHERA DUO ADULT UNITIZED</b> | UNADE10213D  | Extralow | 01 Hydrocephalus Valve Adult size pre-assembled to                           |
|                                  | UNADB10213D  | Low      | 01 Straight Ventricular Cerebral Catheter 13,5 cm with                       |
|                                  | UNADM10213D  | Medium   | 01 Peritoneal Catheter 102 cm; 01 90° Angle Former; 01 stainless steel guide |
|                                  | UNADA10213D  | High     | 01 irrigation tubing; 01 connector for irrigation                            |

| MODEL                       | CATALOGUE N° | PRESSURE   | CONTENT  |
|-----------------------------|--------------|--|--|
| <b>SPHERA DUO INFANTILE</b> | INE907D      | Extralow   | 01 Hydrocephalus Valve Infantile size  |
|                             | INB907D      | Low  | 01 Right Angled Ventricular Cerebral Catheter 07 cm with stainless steel guide |
|                             | INM907D      | Medium   | 01 Peritoneal Catheter 90 cm   |
|                             | INA907D      | High   | 01 irrigation tubing   |
|                             | INE9010D     | Extralow   | 01 Hydrocephalus Valve Infantile size  |
|                             | INB9010D     | Low  | 01 Right Angled Ventricular Cerebral Catheter 10 cm with stainless steel guide |
|                             | INM9010D     | Medium   | 01 Peritoneal Catheter 90 cm   |
|                             | INA9010D     | High   | 01 irrigation tubing   |
|                             | INE9015D     | Extralow   | 01 Hydrocephalus Valve Infantile size  |
|                             | INB9015D     | Low  | 01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide     |
|                             | INM9015D     | Medium   | 01 Peritoneal Catheter 90 cm   |
|                             | INA9015D     | High   | 01 90° Angle Former, 01 irrigation tubing                                      |
|                             | INE9023D     | Extralow   | 01 Hydrocephalus Valve Infantile size  |
|                             | INB9023D     | Low  | 01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide   |
|                             | INM9023D     | Medium   | 01 Peritoneal Catheter 90 cm   |
|                             | INA9023D     | High   | 01 90° Angle Former, 01 irrigation tubing                                      |
|                             | INE10215D*   | Extralow   | 01 Hydrocephalus Valve Infantile size  |
|                             | INB10215D*   | Low  | 01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide     |
|                             | INM10215D*   | Medium   | 01 Peritoneal Catheter 102 cm  |
|                             | INA10215D*   | High   | 01 90° Angle Former, 01 irrigation tubing                                      |
|                             | INE10223D*   | Extralow   | 01 Hydrocephalus Valve Infantile size  |
|                             | INB10223D*   | Low  | 01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide   |
|                             | INM10223D*   | Medium   | 01 Peritoneal Catheter 102 cm  |
|                             | INA10223D*   | High   | 01 90° Angle Former, 01 irrigation tubing                                      |
|                             | INE12015D*   | Extralow   | 01 Hydrocephalus Valve Infantile size  |
|                             | INB12015D*   | Low  | 01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide     |
|                             | INM12015D*   | Medium   | 01 Peritoneal Catheter 120 cm  |
|                             | INA12015D*   | High   | 01 90° Angle Former, 01 irrigation tubing                                      |
| INA12023D*                  | Extralow     | 01 Hydrocephalus Valve Infantile size  |  |
| INE12023D*                  | Low          | 01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide |  |
| INB12023D*                  | Medium       | 01 Peritoneal Catheter 120 cm  |  |
| INM12023D*                  | High         | 01 90° Angle Former, 01 irrigation tubing                                    |  |

\* These configurations are not included in the CE Certificate

| MODEL  | CATALOGUE N°   | PRESSURE                          | CONTENT   |
|--|--|-----------------------------------|---|
| <b>SPHERA<br/>DUO<br/>INFANTILE<br/>UNISHUNT</b> | UINE907D<br>UINB907D<br>UINM907D<br>UINA907D             | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Right Angled Ventricular Cerebral Catheter 07 cm with stainless steel guide<br>01 irrigation tubing                     |
|  | UINE9010D<br>UINB9010D<br>UINM9010D<br>UINA9010D         | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Right Angled Ventricular Cerebral Catheter 10 cm with stainless steel guide<br>01 irrigation tubing                     |
|  | UINE9015D<br>UINB9015D<br>UINM9015D<br>UINA9015D         | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing    |
|  | UINE9023D<br>UINB9023D<br>UINM9023D<br>UINA9023D         | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Straight Ventricular Cerebral Catheter 23,5 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing  |
|  | UINE10215D*<br>UINB10215D*<br>UINM10215D*<br>UINA10215D* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing   |
|  | UINE10223D*<br>UINB10223D*<br>UINM10223D*<br>UINA10223D* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing |
|  | UINE12015D*<br>UINB12015D*<br>UINM12015D*<br>UINA12015D* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 120 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing   |
|  | UINE12023D*<br>UINB12023D*<br>UINM12023D*<br>UINA12023D* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Peritoneal Catheter 120 cm;<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 90° Angle Former, 01 irrigation tubing |

| MODEL  | CATALOGUE N°   | PRESSURE                          | CONTENT  |
|--|--|-----------------------------------|--|
| <b>SPHERA<br/>DUO<br/>INFANTILE<br/>UNITIZED</b> | UNINE9011D<br>UNINB9011D<br>UNINM9011D<br>UNINA9011D | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Infantile size pre-assembled to<br>01 Straight Ventricular Cerebral Catheter 11,5 cm with<br>01 Peritoneal Catheter 90 cm; 01 90° Angle Former; 01 stainless steel guide<br>01 irrigation tubing; 01 connector for irrigation |

\* These configurations are not included in the CE Certificate

| MODEL                      | CATALOGUE N°   | PRESSURE                          | CONTENT   |
|----------------------------|--|-----------------------------------|---|
| <b>SPHERA DUO NEONATAL</b> | NNE905D<br>NNB905D<br>NNM905D<br>NNA905D                 | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Right Angled Ventricular Cerebral Catheter 05 cm with stainless steel guide<br>01 Peritoneal Catheter 90 cm<br>01 irrigation tubing                                    |
|                            | NNE905RD<br>NNB905RD<br>NNM905RD<br>NNA905RD             | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Right Angled Ventricular Cerebral Catheter 05 cm, with stainless steel guide<br>01 Peritoneal Catheter 90 cm;<br>01 irrigation tubing                   |
|                            | NNE907D<br>NNB907D<br>NNM907D<br>NNA907D                 | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Right Angled Ventricular Cerebral Catheter 07 cm with stainless steel guide<br>01 Peritoneal Catheter 90 cm<br>01 irrigation tubing                                    |
|                            | NNE907RD<br>NNB907RD<br>NNM907RD<br>NNA907RD             | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Right Angled Ventricular Cerebral Catheter 07 cm, with stainless steel guide<br>01 Peritoneal Catheter 90 cm; 01 irrigation tubing                      |
|                            | NNE9015D<br>NNB9015D<br>NNM9015D<br>NNA9015D             | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 Peritoneal Catheter 90 cm<br>01 90° Angle Former, 01 irrigation tubing                   |
|                            | NNE9015RD<br>NNB9015RD<br>NNM9015RD<br>NNA9015RD         | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Straight Ventricular Cerebral Catheter 15 cm and stainless steel guide<br>01 Peritoneal Catheter 90 cm;<br>01 90° Angle Former, 01 irrigation tubing    |
|                            | NNE9023D*<br>NNB9023D*<br>NNM9023D*<br>NNA9023D*         | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 Peritoneal Catheter 90 cm<br>01 90° Angle Former, 01 irrigation tubing                 |
|                            | NNE9023RD*<br>NNB9023RD*<br>NNM9023RD*<br>NNA9023RD*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Straight Ventricular Cerebral Catheter 23.5 cm and stainless steel guide<br>01 Peritoneal Catheter 90 cm;<br>01 90° Angle Former, 01 irrigation tubing  |
|                            | NNE10215D*<br>NNB10215D*<br>NNM10215D*<br>NNA10215D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 Peritoneal Catheter 102 cm<br>01 90° Angle Former, 01 irrigation tubing                  |
|                            | NNE10215RD*<br>NNB10215RD*<br>NNM10215RD*<br>NNA10215RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Straight Ventricular Cerebral Catheter 15 cm and stainless steel guide<br>01 Peritoneal Catheter 102 cm;<br>01 90° Angle Former, 01 irrigation tubing   |
|                            | NNE10223D*<br>NNB10223D*<br>NNM10223D*<br>NNA10223D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 Peritoneal Catheter 102 cm<br>01 90° Angle Former, 01 irrigation tubing                |
|                            | NNE10223RD*<br>NNB10223RD*<br>NNM10223RD*<br>NNA10223RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Straight Ventricular Cerebral Catheter 23.5 cm and stainless steel guide<br>01 Peritoneal Catheter 102 cm;<br>01 90° Angle Former, 01 irrigation tubing |
|                            | NNE12015D*<br>NNB12015D*<br>NNM12015D*<br>NNA12015D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 Peritoneal Catheter 120 cm<br>01 90° Angle Former, 01 irrigation tubing                  |
|                            | NNE12015RD*<br>NNB12015RD*<br>NNM12015RD*<br>NNA12015RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Straight Ventricular Cerebral Catheter 15 cm and stainless steel guide<br>01 Peritoneal Catheter 120 cm;<br>01 90° Angle Former, 01 irrigation tubing   |
|                            | NNE12023D*<br>NNB12023D*<br>NNM12023D*<br>NNA12023D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 Peritoneal Catheter 120 cm<br>01 90° Angle Former, 01 irrigation tubing                |
|                            | NNE12023RD*<br>NNB12023RD*<br>NNM12023RD*<br>NNA12023RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 Straight Ventricular Cerebral Catheter 15 cm and stainless steel guide<br>01 Peritoneal Catheter 120 cm;<br>01 90° Angle Former, 01 irrigation tubing   |

\* These configurations are not included in the CE Certificate

| MODEL   | CATALOGUE N°   | PRESSURE                          | CONTENT   |
|---|--|-----------------------------------|---|
| <b>SPHERA<br/>DUO<br/>NEONATAL<br/>UNISHUNT</b> | UNNE905D<br>UNNB905D<br>UNNM905D<br>UNNA905D                 | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Angled 90° Ventricular Cerebral Catheter 05 cm with stainless steel guide<br>01 Irrigation tubing                                      |
|   | UNNE905RD<br>UNNB905RD<br>UNNM905RD<br>UNNA905RD             | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Right Angled Ventricular Cerebral Catheter 05 cm<br>with stainless steel guide; 01 Irrigation tubing                    |
|   | UNNE907D<br>UNNB907D<br>UNNM907D<br>UNNA907D                 | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Right Angled Ventricular Cerebral Catheter 07 cm with stainless steel guide<br>01 Irrigation tubing                                    |
|   | UNNE907RD<br>UNNB907RD<br>UNNM907RD<br>UNNA907RD             | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Right Angled Ventricular Cerebral Catheter 07 cm<br>with stainless steel guide; 01 Irrigation tubing                    |
|   | UNNE9015D<br>UNNB9015D<br>UNNM9015D<br>UNNA9015D             | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 Irrigation tubing                   |
|   | UNNE9015RD<br>UNNB9015RD<br>UNNM9015RD<br>UNNA9015RD         | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm<br>and stainless steel guide; 01 90° Angle Former, 01 Irrigation tubing    |
|   | UNNE9023D*<br>UNNB9023D*<br>UNNM9023D*<br>UNNA9023D*         | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 90° Angle Former, 01 Irrigation tubing                 |
|   | UNNE9023RD*<br>UNNB9023RD*<br>UNNM9023RD*<br>UNNA9023RD*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 90 cm;<br>01 Straight Ventricular Cerebral Catheter 23.5 cm<br>and stainless steel guide; 01 90° Angle Former, 01 Irrigation tubing  |
|   | UNNE10215D*<br>UNNB10215D*<br>UNNM10215D*<br>UNNA10215D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 Irrigation tubing                  |
|   | UNNE10215RD*<br>UNNB10215RD*<br>UNNM10215RD*<br>UNNA10215RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm<br>and stainless steel guide; 01 90° Angle Former, 01 Irrigation tubing   |
|   | UNNE10223D*<br>UNNB10223D*<br>UNNM10223D*<br>UNNA10223D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 102 cm;<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 90° Angle Former, 01 Irrigation tubing                |
|   | UNNE10223RD*<br>UNNB10223RD*<br>UNNM10223RD*<br>UNNA10223RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 102 cm<br>01 Straight Ventricular Cerebral Catheter 23.5 cm<br>and stainless steel guide; 01 90° Angle Former, 01 Irrigation tubing  |
|   | UNNE12015D*<br>UNNB12015D*<br>UNNM12015D*<br>UNNA12015D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 120 cm<br>01 Straight Ventricular Cerebral Catheter 15 cm with stainless steel guide<br>01 90° Angle Former, 01 Irrigation tubing                   |
|   | UNNE12015RD*<br>UNNB12015RD*<br>UNNM12015RD*<br>UNNA12015RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 120 cm;<br>01 Straight Ventricular Cerebral Catheter 15 cm<br>and stainless steel guide; 01 90° Angle Former, 01 Irrigation tubing   |
|   | UNNE12023D*<br>UNNB12023D*<br>UNNM12023D*<br>UNNA12023D*     | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size pre-assembled to<br>01 Peritoneal Catheter 120 cm;<br>01 Straight Ventricular Cerebral Catheter 23.5 cm with stainless steel guide<br>01 90° Angle Former, 01 Irrigation tubing                |
|   | UNNE12023RD*<br>UNNB12023RD*<br>UNNM12023RD*<br>UNNA12023RD* | Extralow<br>Low<br>Medium<br>High | 01 Hydrocephalus Valve Neonatal size with reservoir pre-assembled to<br>01 Peritoneal Catheter 120 cm;<br>01 Straight Ventricular Cerebral Catheter 23.5 cm<br>and stainless steel guide, 01 90° Angle Former, 01 Irrigation tubing |

\* These configurations are not included in the CE Certificate

## SINGLE COMPONENTS

| MODEL                   |          | PRESSURE | CONTENT   |
|-------------------------|----------|----------|---|
| ADULT                   | Extralow | SDVAE    | 01 Hydrocephalus Valve Adult size<br>01 irrigation tubing                   |
|                         | Low      | SDVAB    |   |
|                         | Medium   | SDVAM    |   |
|                         | High     | SDVAA    |   |
| INFANTILE               | Extralow | SDVINE   | 01 Hydrocephalus Valve Infantile size<br>01 irrigation tubing               |
|                         | Low      | SDVINB   |   |
|                         | Medium   | SDVINM   |   |
|                         | High     | SDVINA   |   |
| NEONATAL                | Extralow | SDVNNE   | 01 Hydrocephalus Valve Neonatal size<br>01 irrigation tubing                |
|                         | Low      | SDVNNB   |   |
|                         | Medium   | SDVNNM   |   |
|                         | High     | SDVNNA   |   |
| NEONATAL WITH RESERVOIR | Extralow | SDVNNER  | 01 Hydrocephalus Valve Neonatal size with reservoir<br>01 irrigation tubing |
|                         | Low      | SDVNNBR  |   |
|                         | Medium   | SDVNNMR  |   |
|                         | High     | SDVNNAR  |   |

| PRODUCT                                | FEATURES    |                 |                 |             |
|--|-------------|-----------------|-----------------|-------------|
|  | Length (cm) | ø Internal (mm) | ø External (mm) | Catalog no. |
| Cerebral Catheter<br>Neonatal Straight | 7           | 1,0             | 2,0             | VN7         |
|  | 10          | 1,0             | 2,0             | VN10        |
|  | 15          | 1,0             | 2,0             | VN15        |
|  | 22          | 1,0             | 2,0             | VN22        |
| Infantile Straight                     | 15          | 1,2             | 2,5             | VI15        |
|  | 20          | 1,2             | 2,5             | VI20        |
|  | 23,5        | 1,2             | 2,5             | VI23.5      |
| Adult Straight                         | 15          | 1,4             | 2,7             | VA15        |
|  | 20          | 1,4             | 2,7             | VA20        |
|  | 23,5        | 1,4             | 2,7             | VA23.5      |
| Neonatal Right-Angled                  | 5           | 1,0             | 2,0             | VNG6        |
|  | 6           | 1,0             | 2,0             | VNG5        |
|  | 7           | 1,0             | 2,0             | VNG7        |
| Neonatal Right-Angled with Reservoir   | 5           | 1,0             | 2,0             | VNG5R       |
|  | 6           | 1,0             | 2,0             | VNG6R       |
|  | 7           | 1,0             | 2,0             | VNG7R       |
| Infantile Right-Angled                 | 5           | 1,2             | 2,5             | VIG5        |
|  | 6           | 1,2             | 2,5             | VIG6        |
|  | 7           | 1,2             | 2,5             | VIG7        |
| Infantile Right-Angled with Reservoir  | 5           | 1,2             | 2,5             | VIG5R       |
|  | 6           | 1,2             | 2,5             | VIG6R       |
|  | 7           | 1,2             | 2,5             | VIG7R       |
| Adult Right-Angled                     | 7 x 3       | 1,4             | 2,7             | VAG7        |
|  | 7 x 2       | 1,4             | 2,7             | VAG7.2      |
|  | 9 x 3       | 1,4             | 2,7             | VAG9        |
|  | 9 x 2       | 1,4             | 2,7             | VAG9.2      |
|  | 10 x 3      | 1,4             | 2,7             | VAG10       |
| Adult Right-Angled with Reservoir      | 7 x 3       | 1,4             | 2,7             | VAG7R       |
|  | 10 x 3      | 1,4             | 2,7             | VAG10R      |

| PRODUCT                       | FEATURES    |                 |                 |             |
|-------------------------------|-------------|-----------------|-----------------|-------------|
| Peritoneal Catheter           | Length (cm) | ø Internal (mm) | ø External (mm) | Catalog no. |
| Infantile Peritoneal Catheter | 52          | 1,0             | 2,3             | PI 52       |
|                               | 90          | 1,0             | 2,3             | PI 90       |
|                               | 102         | 1,0             | 2,3             | PI 102      |
|                               | 120         | 1,0             | 2,3             | PI 120      |
| Adult Peritoneal Catheter     | 90          | 1,2             | 2,6             | PA 90       |
|                               | 102         | 1,2             | 2,6             | PA 102      |
|                               | 120         | 1,2             | 2,6             | PA 120      |

| PRODUCT                                  | FEATURES    |                 |                 |             |
|--|-------------|-----------------|-----------------|-------------|
| Cerebral Reservoir                       | Length (cm) | ø Internal (mm) | ø External (mm) | Catalog no. |
| Infantile Ventricular Drainage Reservoir | 13          | 4               | 6               | RHM-1       |
| Adult Ventricular Drainage Reservoir     | 17          | 5,5             | 9,5             | RHM-3       |

| PRODUCT                  | FEATURES                         |                                   |             |
|--------------------------|----------------------------------|-----------------------------------|-------------|
| Atrial Catheter (venous) | Length (mm)<br>Tube of 2,5 x 1,2 | Length (mm)<br>Tube of 1,5 x 0,75 | Catalog no. |
| Cateter Atrial Tipo A    | 450                              | -                                 | AA          |
| Cateter Atrial Tipo B    | 220                              | 150                               | AB          |
| Cateter Atrial Tipo C    | 380                              | 30                                | AC          |
| Cateter Atrial Tipo E    | 220                              | 220                               | AE          |
| Cateter Atrial Tipo H    | 30                               | 380                               | AH          |

| PRODUCT          | FEATURES    |                 |                 |             |
|------------------|-------------|-----------------|-----------------|-------------|
| Conector         | Length (mm) | ø Internal (mm) | ø External (mm) | Catalog no. |
| Straight A Short | 11          | 1,0             | 1,8             | CAC         |
| Straight A Long  | 17          | 1,0             | 1,8             | CAL         |
| Straight B       | 08          | 0,8             | 1,4 / 1,8       | CAB         |
| 90°              | 08          | 1,0             | 1,8             | C 90        |
| Y Short          | 09          | 1,0             | 1,8             | CYC         |
| Y Long           | 12          | 1,0             | 1,8             | CYL         |

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