Operation Manual

Auto Hematology Analyzer Reagent (M-52 series)

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Kullanma Kılavuzu (35~40)

© 2014 - 2017 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.
For this Operation Manual, the issue date is 2017.

Intellectual Property Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyright or patents and does not convey any license under the patent rights or copyright of Mindray, or of others.

Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

Release, amendment, reproduction, distribution, rental, adaptation, translation or any other derivative work of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

mindray

MINDRAY

are the trademarks, registered or otherwise, of Mindray in China and other countries. All other trademarks that appear in this manual are used only for informational or editorial purposes. They are the property of their respective owners.

Responsibility on the Manufacturer Party

Contents of this manual are subject to change without prior notice.

All information contained in this manual is believed to be correct. Mindray shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this manual.

Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- 1) Use of this product is conducted as instructed in this operation manual by Mindray authorized personnel or skilled/trained clinical professionals.;
- 2) The product is stored and transported as instructed in the operation manual and the labeling;
- 3) The product is used only on the specified models.

Warranty

Free service will be provided to:

- Products that are covered by the warranty terms of Mindray.

Pay service will be provided to:

- 1) Expired products;
- 2) Products within the warranty period, but are damaged by improper use or human error, stored in improper environment, or damaged by force majeure or other reason not caused by the product itself.

Return Procedure

Follow the procedure below to return products to Mindray if needed.

Contact Mindray customer service department, inform them the product name and lot No. (the name and lot No. are labeled on the outer package of the product). If the name and lot No. are not legible, your return request will not be accepted. Please specify the product name, model, lot No. and the return reason.

Company Contact

Manufacturer:	Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address:	Mindray Building, Keji 12th Road South, Hi-tech industrial park, Nanshan, Shenzhen 518057, P.R. China
Website:	www.mindray.com
E-mail Address:	service@mindray.com
Tel:	+86 755 81888998
Fax:	+86 755 26582680

EC-Representative	Shanghai International Holding Corp. GmbH(Europe)
Address	Eiffestraße 80, Hamburg 20537, Germany
Tel	0049-40-2513175
Fax	0049-40-255726

Warning

If the reagent accidentally gets into your mouth, seek medical assistance immediately; if the reagent spills on your skin, wash it off with plenty of water immediately; if the reagent spills into your eyes, wash it off with plenty of water and seek medical assistance immediately.

Note:

This product must only be used by Mindray authorized personnel or skilled/trained clinical professionals.

Diluent(M-52D Diluent)

1 Product name

Product name: diluent

Commodity name: M-52D diluent

2 Model

M-52D

3 Package Specification

20L×1

4 Principle

M-52D diluent provides stable solution environment for blood cells, dilute blood sample, generates sheath fluid and provide conductance in blood cell analysis process.

5 Intended use

M-52D diluent applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Mindray, it is used in sample dilution, blood cell counting, volume measurement and hemoglobin measurement.

6 Ingredients

Borate Buffer..... $\leq 0.5\%$

Sodium Chloride..... $\leq 0.1\%$

7 Performance characteristics

- 1) Appearance: transparent liquid, no sediments, grains or flocculi.
- 2) Blank count results: the blank count results of M-52D diluent tested on Mindray Analyzers meet the requirements in Table 1.

Table 1 Requirements of M-52D diluent blank count

Test parameter	Blank count requirements
WBC	$\leq 0.20 \times 10^9/L$
RBC	$\leq 0.02 \times 10^{12}/L$
HGB	$\leq 1g/L$
PLT	$\leq 10 \times 10^9/L$

8 Storage condition and shelf life

The product must be stored in the temperature of 2°C~30°C, and the relative humidity must not exceed 90%. The storage room must be well ventilated and without corrosive gas. The shelf life is two years.

The product must be used in the temperature of 10°C~30°C, and it stays valid for 60 days after being opened.

9 Applicable analyzer models

The product applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd..

10 Instructions for use

- 1) Restore the M-52D diluent to usage temperature;
- 2) Open the external packing and insert the pickup tube into the reagent container based on the matching color of the reagent cap and the connector of analyzer cap assembly;
- 3) Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
- 4) Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis. See the operator's manual of the analyzer for details.

11 Warning

- If the reagent accidentally gets into your mouth, seek medical assistance immediately.
- If the reagent spills on your skin, wash it off with plenty of water immediately.
- If the reagent spills into your eyes, wash it off with plenty of water immediately.

12 Note

- 1) This product can only be used for Mindray Auto Hematology Analyzers. Read the operation manual carefully before use, and do not use expired products.
- 2) Do not use frozen reagents.
- 3) If background results of the reagent are abnormal due to transportation, keep the reagent still in room temperature for 24 hours before use.
- 4) If reagents are contaminated or become abnormal due to other reasons, stop using and then replace them at once.
- 5) Be sure to dispose of reagents, waste and contaminated packing materials according to government regulations.
- 6) The following conditions may produce interference to test results: expired reagent; reagent contaminated by dusts; unclean sample dilution; reagent mixed with those of other manufacturers; residue reagent mixed with newly opened reagent; reagent not used as required.

Lyse(M-52DIFF Lyse)

1 Product name

Product name: Lyse

Commodity name: M-52DIFF Lyse

2 Model

M-52DIFF

3 Package Specification

500mLx4, 500mLx1

4 Principle

M-52DIFF lyse dissolves RBCs in the blood cell analysis process, it processes WBCs to intensify the difference of the WBC sub-populations, and incorporates laser scatter and flow cytometry to realize DIFF channel WBC analysis.

5 Intended use

M-52DIFF applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Mindray, it dissolves RBCs and differentiates WBCs.

6 Ingredients

Surfactant.....<35g/L

7 Performance characteristics

- 1) Appearance: transparent liquid, no sediments, grains or flocculi.
- 2) Blank count results: the blank count results of M-52DIFF lyse tested on Mindray Analyzers meet the requirements in Table 2.

Table 2 Requirements of M-52DIFF lyse blank count

Test parameter	Blank count requirements
WBC	$\leq 0.2 \times 10^9/L$
HGB	$\leq 1g/L$

8 Storage condition and shelf life

The product must be stored in the temperature of 2℃~30℃, and the relative humidity must not exceed 90%. The storage room must be well ventilated and without corrosive gas. The shelf life is two years.

The product must be used in the temperature of 10℃~30℃, and it stays valid for 60 days after being opened.

9 Applicable analyzer models

The product applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd..

10 Instructions for use

- 1) Restore the M-52DIFF lyse to usage temperature;
- 2) Open the external packing and insert the pickup tube into the reagent container based on the matching color of the reagent cap and the connector of analyzer cap assembly;
- 3) Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
- 4) Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis. See the operator's manual of the analyzer for details.

11 Warning

- If the reagent accidentally gets into your mouth, seek medical assistance immediately.
- If the reagent spills on your skin, wash it off with plenty of water immediately.
- If the reagent spills into your eyes, wash it off with plenty of water immediately.

12 Note

- 1) This product can only be used for Mindray Auto Hematology Analyzers. Read the operation manual carefully before use, and do not use expired products.
- 2) Do not use frozen reagents.
- 3) If background results of the reagent are abnormal due to transportation, keep the reagent still in room temperature for 24 hours before use.
- 4) If reagents are contaminated or become abnormal due to other reasons, stop using and then replace them at once.
- 5) Be sure to dispose of reagents, waste and contaminated packing materials according to government regulations.
- 6) The following conditions may produce interference to test results: expired reagent; reagent contaminated by dusts; unclean sample dilution; reagent mixed with those of other manufacturers; residue reagent mixed with newly opened reagent; reagent not used as required.

Lyse(M-52LH Lyse)

1 Product name

Product name: Lyse

Commodity name: M-52LH Lyse

2 Model

M-52LH

3 Package Specification

100mLx4, 100mLx1

4 Principle

M-52LH lyse further processes the blood cells that have been processed by M-52DIFF lyse, it turns hemoglobin into methemoglobin, which then measured by the colorimetric method. It also processes WBCs to intensify the difference between basophils and other WBC sub-populations, and incorporates laser scatter and flow cytometry to realize basophil analysis.

5 Intended use

M-52LH Lyse applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Mindray, it dissolves RBCs, realizes WBC counting, basophil differential and hemoglobin measurement.

6 Ingredients

Quaternary ammonium salts.....<50g/L

Isopropanol.....2-10g/L

7 Performance characteristics

- 1) Appearance: transparent liquid, no sediments, grains or flocculi.
- 2) Blank count results: the blank count results of M-52LH lyse tested on Mindray Analyzers meet the requirements in Table 3.

Table 3 Requirements of M-52LH lyse blank count

Test parameter	Blank count requirements
WBC	$\leq 0.2 \times 10^9/L$
HGB	$\leq 1g/L$

8 Storage condition and shelf life

The product must be stored in the temperature of 2°C~30°C, and the relative humidity must not exceed 90%. The storage room must be well ventilated and without corrosive gas. The shelf life is 2 years.

The product must be used in the temperature of 10°C~30°C, and it stays valid for 60 days after being opened.

9 Applicable analyzer models

The product applies to BC-5000, BC-5150, BC-5120, BC-5130 and BC-5140 Auto Hematology Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd..

10 Instructions for use

- 1) Restore the M-52LH lyse to usage temperature;
- 2) Open the external packing and insert the pickup tube into the reagent container based on the matching color of the reagent cap and the connector of analyzer cap assembly;
- 3) Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
- 4) Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis. See the operator's manual of the analyzer for details.

11 Warning

- If the reagent accidentally gets into your mouth, seek medical assistance immediately.
- If the reagent spills on your skin, wash it off with plenty of water immediately.
- If the reagent spills into your eyes, wash it off with plenty of water immediately.

12 Note

- 1) This product can only be used for Mindray Auto Hematology Analyzers. Read the operation manual carefully before use, and do not use expired products.
- 2) Do not use frozen reagents.
- 3) If background results of the reagent are abnormal due to transportation, keep the reagent still in room temperature for 24 hours before use.
- 4) If reagents are contaminated or become abnormal due to other reasons, stop using and then replace them at once.
- 5) Be sure to dispose of reagents, waste and contaminated packing materials according to government regulations.
- 6) The following conditions may produce interference to test results: expired reagent; reagent contaminated by dusts; unclean sample dilution; reagent mixed with those of other manufacturers; residue reagent mixed with newly opened reagent; reagent not used as required.

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