Prolactin

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Prolactin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of prolactin in serum, as an aid in the diagnosis and treatment of pituitary disorders.

Catalog Number: **L2KPR2** (200 tests),
L2KPR6 (600 tests)

Test Code: **PRL** Color: **Dark Blue**

Summary and Explanation

Human prolactin is a polypeptide hormone of the anterior pituitary with a molecular mass of about 22,800. It plays an essential role in the secretion of milk and has the ability to suppress gonadal function. The very existence of prolactin as a substance distinct from human growth hormone was established only in 1970. Since that time, determination of prolactin has become an important tool in the investigation of amenorrhea, galactorrhea and hypothalamic-pituitary disorders.

As a reference range for circulating prolactin, the literature suggests concentrations up to approximately 20 ng/mL. Values are distinctly elevated at birth but decline to adult levels in less than three months.

Women are reported to have slightly higher mean levels than men, with a slight rise at puberty — apparently estrogen related — and a corresponding fall at menopause. During pregnancy, the prolactin level climbs steadily to ten or twenty times its former value, then drops back down to normal after delivery — within three weeks in nonnursing mothers. In those who breast-feed, the decline to normal is more gradual because of the prompt and dramatic surges in prolactin release induced by suckling. Women taking oral contraceptives or under estrogen treatment may have prolactin levels higher than normal.

In assessing the significance of moderate elevations, it is important to keep in mind

that prolactin is a stress hormone. Not only surgery, but events no more distressing than venipuncture or a clinical interview have been reported to occasion a transient rise. Moreover, the release of prolactin is inherently episodic, and day-to-day fluctuations with CVs as high as 30% have been encountered. Finally, there is a sleep-related diurnal variation: prolactin levels increase during sleep and reach their lowest a few hours after waking. The advice sometimes given to draw samples "between nine and noon" is based on the assumption that subjects observe reasonably normal waking hours.

Principle of the Procedure

IMMULITE 2000 Prolactin is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes

Time to first result: 35 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Prolactin has not been tested with all possible variations of tube types.

Volume Required: 25 µL serum

Storage: 7 days at 2–8°C, or
3 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.
	Contains: 2-methyl-2H-isothiazol-3-one; Prolactin Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Prolactin Bead Pack (L2PR12)

With barcode. 200 beads, coated with monoclonal murine anti-prolactin. Stable at 2–8°C until expiration date.

L2KPR2: 1 pack **L2KPR6:** 3 packs

Prolactin Reagent Wedge (L2PRA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-prolactin, with preservative. Stable at 2–8°C until expiration date.

L2KPR2: 1 wedge **L2KPR6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Prolactin Adjustors (LPRL, LPRH)

Two vials (Low and High), of lyophilized human prolactin in a serum/buffer matrix. Reconstitute each vial with **2.0 mL** distilled or deionized water, and mix by gentle inversion. Stable after reconstitution for 60 days (aliquotted) at –20°C.

L2KPR2: 1 set **L2KPR6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Prolactin Sample Diluent (L2PRZ)

For on-board dilution of high samples. 25 mL of concentrated (ready-to-use) processed, prolactin-free serum/buffer matrix. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2PRZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
4 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of prolactin.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Based on its relationship to IMMULITE Prolactin (see Method Comparison), the assay can be expected to have essentially the same reference range.

Group	n	Median	96% Range
Adult Males	19	6.2 ng/mL	2.5 – 17
		131 mIU/L	53 – 360

A study¹⁷ performed with IMMULITE Prolactin yielded the following results.

Group	n	Median	95% Range
Adult Females	115	9.4 ng/mL	1.9 – 25
		199 mIU/L	40 – 530

As summarized in technical report ZB157,¹⁷ in a study which followed normal ovulating women on a daily basis throughout one complete cycle, somewhat higher values were obtained — an outcome consistent with the increasingly stressful sample collection process entailed by the study design.

A cross-sectional study of pediatric fertility values performed with IMMULITE Prolactin at a "wellness" clinic in the southwestern United States yielded the following results.

Group	Age (yr)	n	Prolactin, ng/mL	
			Median	Central 95%
Females	Cord	28	380	200 – 675
	0.1 – 0.5	28	15	1 – 140
	0.6 – 9	55	11	2 – 43
Males	Cord	27	295	150 – 565
	0.1 – 0.5	36	19	4 – 65
	0.6 – 9	55	8	0.6 – 29
Combined	Cord	55	340	160 – 665
	0.1 – 0.5	64	117	2 – 125
	0.6 – 9	110	9	1 – 40

Group	Age (yr)	n	Prolactin, mIU/L	
			Median	Central 95%
Females	Cord	28	8056	4240 – 14,310
	0.1 – 0.5	28	318	21 – 2968
	0.6 – 9	55	233	42 – 912
Males	Cord	27	6254	3180 – 11,978
	0.1 – 0.5	36	403	85 – 1378
	0.6 – 9	55	170	13 – 615
Combined	Cord	55	7208	3392 – 14,098
	0.1 – 0.5	64	2480	42 – 2650
	0.6 – 9	110	191	21 – 848

Pediatric: Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Prolactin assay in accordance with CLSI guideline EP28-A3C.¹⁸ Samples were collected prospectively from apparently healthy pediatric subjects, using predefined inclusion criteria. Reference values were generated for subpopulations based on age and Tanner stage subgroups based on physiological development. The study was designed to establish reference values across genders, and to include approximately equal numbers of males and females within each age or Tanner stage subgroup. The subject's Tanner stage was assessed based on pubic hair and genitalia/breast development. The scale proposed by Neinstein and Kaufman was used for the determination of the Tanner stages.¹⁹

The reference intervals and Tanner values are based on the central 90% (5th and 95th percentiles). Where sample sizes were insufficient to calculate the 5th or

95th percentile, the minimum or maximum observed values are presented in the Reference Intervals tables.

IMMULITE 2000/2000 XPi Prolactin Pediatric Reference Intervals

Male

Age (Years)	n	Median	Range
		ng/mL	
2 – 3	10	7.8	< 3.5* – > 21.9 [†]
4 – 9	57	6.7	3.1 – 17.3
10 – 16	203	5.8	2.7 – 11.8
17 – 21	37	7.1	4.6 – 15.2

Age (Years)	n	Median	Range
		mIU/L	
2 – 3	10	164.3	< 74.2* – > 464.3 [†]
4 – 9	57	142.0	67.2 – 366.1
10 – 16	203	123.0	57.7 – 249.3
17 – 21	37	150.5	96.9 – 322.2

Female

Age (Years)	n	Median	Range
		ng/mL	
2 – 3	17	8.0	< 2.8* – > 13.9 [†]
4 – 9	47	6.3	2.7 – 14.3
10 – 12	93	7.0	3.2 – 17.6
13 – 21	127	8.1	3.8 – 22.1

Age (Years)	n	Median	Range
		mIU/L	
2 – 3	17	169.6	< 59.4* – > 294.7 [†]
4 – 9	47	133.6	58.1 – 302.7
10 – 12	93	148.4	68.9 – 372.9
13 – 21	127	171.7	80.6 – 467.7

*Value presented is the minimum reportable value observed; insufficient sample size to calculate a 5th percentile limit.

[†] Value presented is the maximum value observed; insufficient sample size to calculate a 95th percentile limit.

IMMULITE 2000/2000 XPi Prolactin Pediatric Reference Intervals by Tanner Stage

Male

Tanner Stage	n	Median	Range
		ng/mL	
1	73	6.0	3.4 – 17.8
2	64	5.3	2.3 – 11.4
3	63	5.2	3.2 – 11.3
4	59	6.6	2.8 – 15.2
5	48	6.7	4.3 – 13.2

Tanner Stage	n	Median	Range
		mIU/L	
1	73	127.2	70.8 – 377.6
2	64	112.4	48.8 – 242.2
3	63	110.2	67.8 – 238.7
4	59	139.9	59.4 – 322.2
5	48	142.0	91.2 – 280.4

Female

Tanner Stage	n	Median	Range
		ng/mL	
1	73	6.5	2.8 – 17.1
2	47	6.8	3.1 – 19.9
3	65	8.0	3.3 – 17.7
4	47	8.2	3.4 – 22.1
5	52	7.3	3.8 – 21.9

Tanner Stage	n	Median	Range
		mIU/L	
1	73	137.8	59.4 – 361.9
2	47	144.2	67.0 – 422.3
3	65	169.6	69.7 – 375.5
4	47	173.8	72.9 – 467.7
5	52	155.8	81.1 – 464.4

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely

exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

ng/mL $\times$ 21.2 $\rightarrow$ mIU/L [3rd IS 84/500]

Reportable Range: 0.5–150 ng/mL (10.6–3180 mIU/L) in terms of the 3rd IS 84/500

Analytical Sensitivity: 0.5 ng/mL (10.6 mIU/L)

High-dose Hook Effect: None up to 20,500 ng/mL (434,600 mIU/L)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with four prolactin solutions (190, 356, 928 and 1740 ng/mL), were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for prolactin. (See "Specificity" table.)

Alternate Sample Type: A *limited* study of matched samples from laboratory volunteers yielded comparable results for serum and heparinized plasma.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 375 mg/dL has no

effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to IMMULITE Prolactin on 148 samples. (Concentration range: approximately 1 to 140 ng/mL. See graph.) By linear regression:

(IML 2000) = 1.00 (IML) – 0.96 ng/mL
 $r = 0.988$

Means:

44 ng/mL (IMMULITE 2000)

45 ng/mL (IMMULITE)

References

- 1) Cowden E, et al. Laboratory assessment of prolactin status. *Ann Clin Biochem* 1979;16:113–21.
- 2) Cowden E, et al. Tests of prolactin secretion. *Lancet* 1979;1155–8.
- 3) Djursing H. Short- and long-term fluctuations in plasma prolactin concentration in normal subjects. *Acta Endocrinol* 1981; 97:1–6.
- 4) Greer ME, et al. Prevalence of hyperprolactinemia in anovulatory women. *Obstet Gynecol* 1980;56:65–9.
- 5) Harrington RA, et al. Metoclopramide: an updated review. *Drugs* 1983;25:451–94.
- 6) Healy DL, et al. Pituitary autonomy in hyperprolactinemia secondary amenorrhea: results of hypothalamic-pituitary testing. *J Clin Endocrinol Metab* 1977;44:809–18.
- 7) Kleinberg DL, Noel GL, Frantz AG. Galactorrhea: a study of 235 cases. *N Engl J Med* 1977;296:589–600.
- 8) Knuth UA, Friesen HG. Prolactin and pregnancy. *Curr Top Exp Endocrinol* 1983;4:69–96.
- 9) Maxson WS, Hammond CB. Hyperprolactinemia. *Contemp Ob/Gyn* 1982 Jan;19(1):49–63, 19:67–83.
- 10) Pepperell RJ. Prolactin and reproduction. *Fertil Steril* 1981;35:267–74.
- 11) Babson, AL. The IMMULITE Automated Immunoassay System. *J Clin Immunoassay* 1991;14:83–8.
- 12) Thorner MO. Prolactin: clinical physiology and management of hyperprolactinemia. In: Martini L, Besser GM, et al, editors. *Clinical endocrinology*. New York: Academic Press, 1977;319–61.
- 13) Tolis G, Franks S. Physiology and pathology of prolactin secretion. In: Tolis G, et al, editors. *Clinical endocrinology: a pathophysiological approach*. New York: Raven Press, 1979;291–317.
- 14) Tolis G, et al. Aspects of prolactin pathophysiology. *Periodicum Biologorum* 1983;85 Suppl 1:29–36.
- 15) Yen SSC. Neuroendocrine regulation of gonadotropin and prolactin secretion in women: disorders in reproduction. In: Vaitukaitis J, editor. *Clinical reproductive endocrinology*. New York: Elsevier Biomedical, 1982;137–76.
- 16) Zacur HA. Use of the human prolactin immunoassay. *J Clin Immunoassay* 1983

Spring;6(1):63–67. 17) Vankrieken L. IMMULITE reproductive hormone assays: multicenter reference range data. Los Angeles: Diagnostic Products Corporation, 2000. Document No. ZB157-D. 18) Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C. 19) Neinstein LS and Kaufman FR, Chapter 1: Normal Physical Growth and Development in Neinstein L.S. *Adolescent Health Care: A Practical Guide*, 4th ed.

Technical Assistance

Contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	5.5	0.18	3.3%	0.27	4.9%
2	11.9	0.39	3.3%	0.57	4.8%
3	22.3	0.60	2.7%	0.90	4.0%
4	64.5	1.8	2.8%	3.3	5.1%
5	121	3.4	2.8%	6.4	5.3%
6	127	4.3	3.4%	5.8	4.6%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	34	—	—
	4 in 8	17	17	100%
	2 in 8	8.2	8.5	96%
	1 in 8	4.0	4.3	93%
2	8 in 8	73	—	—
	4 in 8	36	37	97%
	2 in 8	17	18	94%
	1 in 8	8.1	9.1	89%
3	8 in 8	84	—	—
	4 in 8	41	42	98%
	2 in 8	20	21	95%
	1 in 8	9.6	11	87%
4	8 in 8	120	—	—
	4 in 8	59	60	98%
	2 in 8	28	30	93%
	1 in 8	13	15	87%

Specificity


Compound ¹	Amount Added ² (ng/mL)	% Cross- reactivity ³
FSH	1000	ND
	200	ND
HCG	1000	ND
	100	ND
hPL	200,000	ND
LH	1000	ND
	500	ND
TSH	1000	0.14%

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Method Comparison: Prolactin: Prolactine.

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determinabile. **Method Comparison.** Prolactin:

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.

© 2015 Siemens Healthcare Diagnostics.
All rights reserved.

Made in: UK

Siemens Healthcare
Diagnostics Products Ltd.
Glyn Rhonwy, Llanberis,
Gwynedd LL55 4EL
United Kingdom

2017-09-27

PIL2KPR – 23

Changes in this Edition:

cc#EU23186: 1) In Expected Values, added Pediatric information. 2) Added Pediatric references. 3) Updated Technical Assistance.

Understanding the Symbols

Understanding the Symbols	En English
Erklärung der Symbole	De Deutsch
Descripción de los símbolos	Es Español
Explication des symboles	Fr Français
Definizione dei simboli	It Italiano
Descrição dos símbolos	Pt Português

The following symbols may appear on the product labeling: / Die folgenden Symbole können auf dem Produktetikett erscheinen: / Los siguientes símbolos pueden aparecer en la etiqueta del producto: / Les symboles suivants

peuvent apparaître sur les étiquettes des produits : / Sull'etichetta del prodotto possono essere presenti i seguenti simboli: / Os seguintes símbolos podem aparecer no rótulo dos produtos:

Symbol Definition

En: *In vitro* diagnostic medical device

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Fr: Fabricant

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Pt: Fabricante

En: Authorized Representative in the European Community

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Pt: Marca CE

Symbol Definition

En: CE Mark with identification number of notified body

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