



**Medtronic**

THERAPIES

# CD HORIZON<sup>®</sup> SOLERA<sup>™</sup> Fenestrated Screw Spinal System

Surgical Technique

4.75mm and 5.5/6.0mm

CD HORIZON<sup>®</sup> SOLERA<sup>™</sup> Fenestrated Screws provide immediate enhanced fixation in poor bone quality through the application of Fenestrated Screw Cement

SERVICES & SUPPORT

ENABLING TECHNOLOGIES

Surgical **Synergy**<sup>™</sup>



## Table of Contents

Implant Features and Instrument Set	4
The Locking Mechanism of the Lock Sleeve	7
Assembling the Adapter Driver	8
Preparation of the Spine	10
Attaching the Fenestrated Screw to the Adapter Driver	10
Screw Insertion	12
Mixing the Cement	13
Option A. Using the Bone Filler Device for Delivery	
Step 1: Loading the Bone Filler Device	16
Step 2: Checking the Viscosity of the Cement	16
Step 3: Injecting the Cement	17
Option B. Using the CDS Cement Delivery System	
Step 1: Preparing the CDS Gun	20
Step 2: Checking the Viscosity of the Cement	23
Step 3: Priming the Bone Filler Device	23
Step 4: Injecting the Cement	25
Removal of the Adapter Driver	27
Completing the Procedure	28
Disassembling the Adapter Driver	28
Screw Explanation	29
Product Ordering Information	30
Important Information on the CD HORIZON® SOLERA™ Fenestrated Screw Spinal System	32
Explanation of Symbols	34
Important Information on Resin for Fenestrated Screw (FS) Cement High Viscosity, Radiopaque Bone Cement	35




## Implant Features and Instrument Set

CD HORIZON® SOLERA™ Fenestrated Screw Spinal System offers the same design and capabilities of the CD HORIZON® SOLERA® Cannulated Multi-Axial Screw, with the addition of six fenestrations near the tip. These features allow cement to pass directly through the screw shaft and fenestrations, providing immediate enhanced fixation in poor bone quality.

### IMPLANT CD HORIZON® SOLERA™ FENESTRATED SCREW

Available for CD HORIZON® SOLERA®  
4.75mm and 5.5/6.0mm Spinal Systems



### INSTRUMENTS ADAPTER DRIVER

The CD HORIZON® SOLERA™ Fenestrated Screw Adapter Driver comes in three parts.

Lock Sleeve 4.75mm or 5.5/6.0mm  
(reusable item)



Adapter Driver QC Shaft  
(reusable item)

Adapter Driver T25 Tip  
(disposable item)

## Implant Features and Instrument Set

### CEMENT MIXING

Kyphon® Mixer



Luer Lock



Cartridge



Syringe



Paddle



Funnel



Plunger



Fenestrated Screw (FS) Cement  
powder packet and liquid bottle



## Implant Features and Instrument Set

### CEMENT INJECTION



## The Locking Mechanism of the Lock Sleeve

### Locked position

Latch is closed (Figure 1.1)

### Engaged position

Latch is opened until the hard stop (Figure 1.2)

### Released position

Push the latch down and continue rotating until fully opened (Figure 1.3)



Figure 1.1 Locked position

Adapter Driver and screw are fully secured and cannot be disengaged



Figure 1.2 Engaged position

Adapter Driver is secure, screw can be attached and disconnected




Figure 1.3 Released position

Adapter Driver can be assembled/disassembled

## Assembling the Adapter Driver

» Place the Lock Sleeve in the released position (**Figure 1.3**)

» Insert the Adapter Driver T25 Tip into the top of the Lock Sleeve ensuring the T25 portion is inserted first (**Figure 2.1**)

The T25 Tip will extend out of the bottom of the Lock Sleeve (**Figure 2.2**)




Figure 2.1 Insert the Tip with T25 portion leading



Figure 2.2

» Insert the Adapter Driver QC Shaft into the top of the Lock Sleeve (**Figure 3.1**)



Figure 3.1

» Ensure that the QC Section extends above the Lock Sleeve (**Figure 3.2**)




Figure 3.2

## Assembling the Adapter Driver

» Rotate the latch until it clicks into the engaged position (**Figure 4**)



Figure 4

» Press the tip of the Adapter Driver against a hard surface to move the inner shafts upwards (**Figure 5.1**) and rotate the latch to the locked position (**Figure 5.2**)



Figure 5.1



Figure 5.2

## Preparation of the Spine

Prepare the spine according to the appropriate CD HORIZON® Spinal System Surgical Technique.

## Attaching the Fenestrated Screw to the Adapter Driver

- » Ensure the assembled Adapter Driver is in the locked position (refer back to **Figure 5.2**)
- » Thread the Adapter Driver to the selected Fenestrated Screw from the screw caddy (**Figure 6**)
- » Ensure the Fenestrated Screw is securely tightened onto the Adapter Driver to prevent cement leakage and screw toggling

### ✓ Note

If the Adapter Driver is not engaging with the screw, ensure the inner drive was properly stacked-up before locking the latch (**Figure 5.1**)

### ! Important

After tightening the screw, ensure the latch is fully closed (**Figure 5.2**)



Figure 6

## Attaching the Fenestrated Screw to the Adapter Driver

- » Break off the VERIFY!® Implant Tracking Tag as shown (**Figure 7**)



Figure 7

## Screw Insertion

- » Attach the Quick Connect Handle to the Adapter Driver (Figure 8)




Figure 8

- » Using the assembled Fenestrated Screw Adapter Driver and Quick Connect Handle, insert screws into levels to be instrumented. **The Adapter Driver must remain attached to the screw to allow for cement injection (Figure 9)**

**! Important**

To reduce the risk of cement leakage, it is essential to verify with intraoperative imaging that the tip of each screw is within the confines of the vertebral body, just beyond the midpoint




Figure 9 Screw insertion

## Mixing the Cement

- » To prepare the FS Cement, insert the funnel into the Kyphon® Mixer. Make sure the purple latches are facing up so they are in the **closed position**. If using the CDS, thread the Luer Lock adapter onto the spout at the bottom of the mixer body (Figure 10)



Figure 10

- » Pour entire contents of the FS Cement powder packet into the Kyphon® Mixer (Figure 11)



Figure 11

## Mixing the Cement

- » Take the syringe and blunt needle and thread together. Snap off the cap of the FS Cement liquid bottle and withdraw the contents of the FS Cement liquid with the syringe and needle assembly (Figure 12.1, 12.2)



Figure 12.1



Figure 12.2

- » Inject the liquid into the Kyphon® Mixer (Figure 13)



Figure 13

## Mixing the Cement

- » Remove the funnel and insert the paddle (Figure 14.1)



2 min.

- » Mix the cement for approximately two minutes by agitating the paddle (Figure 14.2)

### Note

The mixing and curing characteristics of the cement will vary according to temperature and other factors. For example, in cooler environments it may be desired to mix for slightly longer to achieve a more viscous consistency before filling the Bone Filler Devices. Refer to the FS Cement instruction for use



Figure 14.1



Figure 14.2

- » Remove the paddle and insert the plunger (Figure 15)



Figure 15



## Option A. Using the Bone Filler Device for Delivery

### STEP 1: LOADING THE BONE FILLER DEVICE

- » Attach the Bone Filler Device to the Kyphon® Mixer by threading it onto the Luer Lock at the base of the mixer. Lower both valve levers and fill the Bone Filler Devices by depressing the plunger (Figure 16)

**Note**

Each Bone Filler Device distributes 1.5cc of cement



Figure 16

### STEP 2: CHECKING THE VISCOSITY OF THE CEMENT



8 min.

Begin checking the viscosity of the cement after approximately eight minutes. To check the viscosity, extrude a small amount from the tip of a filled Bone Filler Device.

**If the cement is sticky, it is not ready and should be allowed to stand further. The cement should retain its shape when compressed into a ball (Figure 17).**

**Note**

Do not roll between finger tips; use a pinching motion

**Note**

Time should be measured beginning with initial mixing of the cement, not from the loading of the Bone Filler Device. Check to make sure that the cement has lost its sheen, a sign that it has become doughy and is ready for application




Figure 17

## Option A. Using the Bone Filler Device for Delivery

### STEP 3: INJECTING THE CEMENT

- » Remove the Quick Connect Handle from the Adapter Driver (Figure 18)



Figure 18

- » Insert the Bone Filler Device (Figure 19)



Figure 19

## Option A. Using the Bone Filler Device for Delivery

### STEP 3: INJECTING THE CEMENT

**Note**

As an option, the Stabilisation Clip can be used to hold the Bone Filler Device in place during cement injection (Figure 20)



Figure 20

- » Insert the plunger portion of the Bone Filler Device into the outer sheath (Figure 21) and inject the cement by depressing the plunger until the desired volume of cement is delivered (Figure 22)

**Important**

Cement injection should only be performed under fluoroscopic control

**Note**

Each Bone Filler Device can distribute 1.5cc of cement. The Fenestrated Screw and Adapter Driver cannula will fill with approximately 0.4cc of cement before any is ejected into the vertebral body. To achieve adequate fixation, it is recommended that 0.8cc of cement be implanted in the vertebral body for each screw in the thoracic spine (except for T11 and T12) and 1.8cc of cement be used for each screw in the lumbar spine along with T11 and T12

Based on internal cadaveric testing (Source: Medtronic TR08-426)  
Testing is not indicative of human clinical outcome



Figure 21

## Option A. Using the Bone Filler Device for Delivery

### STEP 3: INJECTING THE CEMENT

- » Repeat injection for each Fenestrated Screw

**Important**

- » Remove the Stabilisation Clip and Bone Filler Device from the Adapter Driver before cement is cured
- » Leave the Adapter Driver in place until the cement is fully cured to avoid pulling cement into the saddle of the screw

Once the cement is set, see page 27 for removal of the Adapter Driver and completing the procedure.



Figure 22

## Option B. Using the CDS Cement Delivery System

### STEP 1: PREPARING THE CDS GUN

- » Ensure the gun is working properly by removing the protective rubber stopper from the metal part at the end of the tubing (Figure 23)



Figure 23

- » Squeeze the handle several times to ensure the saline is flowing out of the rubber tubing (Figure 24.1, 24.2)

At this point you can mix the cement (see page 13 for instructions)



Figure 24.1

**Note**

The cartridge should be loaded immediately once the cement is mixed



Figure 24.2

## Option B. Using the CDS Cement Delivery System

### STEP 1: PREPARING THE CDS GUN

- » Attach the cartridge to the Kyphon® Mixer by threading it onto the Luer Lock at the base of the mixer (Figure 25)

Lower both valve levers and fill the cartridge by depressing the plunger



Figure 25

- » Return the valve levers to the locked position and unthread the cartridge from the mixer (Figure 26). Wipe away any excess cement from the end of the cartridge
- » Attach the second cartridge and fill as previous



Figure 26

## Option B. Using the CDS Cement Delivery System

### STEP 1: PREPARING THE CDS GUN

- » Firmly thread the connector at the end of the tubing to the cartridge (**Figure 27**)



Figure 27

### STEP 1: PREPARING THE CDS GUN

- » A tight connection is necessary to maintain the internal pressure. Prime the gun several times until cement flows through the cartridge (**Figure 28**)

- » When cement starts flowing out the end of the cartridge you can begin to test the viscosity of the cement

Ensure that you keep the tip of the cartridge clean from excess cement



Figure 28

## Option B. Using the CDS Cement Delivery System

### STEP 2: CHECKING THE VISCOSITY OF THE CEMENT



8 min.

Begin checking the viscosity of the cement after approximately eight minutes. To check the viscosity, extrude a small amount from the tip of the cement cartridge.

**If the cement is sticky, it is not ready and should be allowed to stand further. The cement should retain its shape when rolled into a ball (Figure 29).**

✓ **Note**

Do not roll between fingertips; use a pinching motion

✓ **Note**

Time should be measured beginning with initial mixing of the cement, not from the filling of the cartridge. Check to make sure that the cement has lost its sheen, a sign that it has become doughy and is ready for application



Figure 29

### STEP 3: PRIMING THE BONE FILLER DEVICE

- » Thread the Bone Filler Device onto the end of the cartridge (**Figure 30**)



Figure 30

## Option B. Using the CDS Cement Delivery System

### STEP 3: PRIMING THE BONE FILLER DEVICE

- » Pull the trigger several times until cement is flowing from the end of the Bone Filler Device (**Figure 31**)



Figure 31

- » To stop cement flow, depress the purple button for 1-2 seconds (**Figure 32**)



Figure 32

## Option B. Using the CDS Cement Delivery System

### STEP 4: INJECTING THE CEMENT

- » Remove the Quick Connect Handle from the Adapter Driver and insert the Bone Filler Device (**Figure 33**)



Figure 33

- » Take the Stabilisation Clip and attach this to the Adapter Driver and Bone Filler Device (**Figure 34.1, 34.2**)
- The Bone Filler Device is now locked onto the Adapter Driver and will not eject upwards during cement injection



Figure 34.1

Figure 34.2

## Option B. Using the CDS Cement Delivery System

### STEP 4: INJECTING THE CEMENT

- » Pull the trigger until the desired volume of cement is delivered (**Figure 35**). To stop cement flow, depress the purple button for 1-2 seconds (refer back to **Figure 32** if required)

#### ! Important

Cement injection should only be performed under fluoroscopic control

#### ✓ Note

Each squeeze of the CDS trigger will deliver 0.2cc of cement. The Fenestrated Screw and Adapter Driver cannula will fill with approximately 0.4cc of cement before any is ejected into the vertebral body. To achieve adequate fixation, it is recommended that 0.8cc of cement be implanted in the vertebral body for each screw in the thoracic spine (except for T11 and T12) and 1.8cc of cement be used for each screw in the lumbar spine along with T11 and T12

Based on internal cadaveric testing (Source: Medtronic TR08-426)  
Testing is not indicative of human clinical outcome



Figure 35

- » Disconnect the Stabilisation Clip from the Adapter Driver and move the Bone Filler Device and cartridge to the next screw. Reconnect the Stabilisation Clip and repeat injection

If you require the second cartridge, unscrew the first cartridge and secure firmly the new cartridge to the tubing and the connector (see **Figure 27**).

#### ! Important

- » Remove the Stabilisation Clip and Bone Filler Device from the Adapter Driver before cement is cured
- » Leave the Adapter Driver in place until the cement is fully cured to avoid pulling cement into the saddle of the screw

## Removal of the Adapter Driver

### DISENGAGING THE ADAPTER DRIVER FROM THE FENESTRATED SCREW

- » Once the cement is fully cured, place the Adapter Driver in the engaged position (**Figure 36**)

#### ! Important

Firmly hold the outer sleeve to avoid rotation while opening the latch (**Figure 36 – zoom**)



Figure 36

- » With your thumb in-between the latch and the driver, unthread the driver in an anticlockwise direction (**Figure 37**). Your thumb will prevent the latch from re-locking

The Adapter Driver can then be removed and disassembled on the back table

#### ! Important

Do not rotate the QC attachment at this stage. This will unthread the screw from the bone



Figure 37

## Completing the Procedure

See appropriate CD HORIZON® Spinal System Surgical Technique for guidance in rod, set screw and transverse link placement, final tightening and bone grafting.

## Disassembling the Adapter Driver

- » Place the Adapter Driver in the released position (**Figure 38**)



Figure 38 Adapter Driver in released position

- » Remove the Adapter Driver QC Shaft and the Adapter Driver T25 Tip from the Lock Sleeve by tipping the assembly upside down (**Figure 39.1, 39.2**)

The Adapter Driver QC Shaft and the Lock Sleeve can be reprocessed

### ! Important

The Adapter Driver T25 Tip is discarded




Figure 39.1 Removal of the Adapter Driver QC Shaft

Figure 39.2 Removal of the Adapter Driver T25 Tip

## Screw Explanation

If removal of a Fenestrated Screw is necessary, attach a standard Medtronic Quick Connect Handle to the T25 screwdriver.

Next, fully engage the T25 end of the driver into the screw head; then thread the instrument sleeve into the screw head. Turn anticlockwise until the Fenestrated Screw has been removed.

### ✓ Note

The cured cement will break away from the screw under this torsional force to allow for screw removal. Refer to the appropriate CD HORIZON® System Surgical Technique for proper explanation of the other CD HORIZON® implants

## Product Ordering Information

### COBALT CHROME/TITANIUM MULTI-AXIAL SCREWS

4.75mm	5.5/6.0mm	Description	4.75mm	5.5/6.0mm	Description
54840024530*	55840024530*	4.5mm x 30mm	54840027570*	55840027570*	7.5mm x 70mm
54840024535*	55840024535*	4.5mm x 35mm	54840027575*	55840027575*	7.5mm x 75mm
54840024540*	55840024540*	4.5mm x 40mm	54840027580*	55840027580*	7.5mm x 80mm
54840024545*	55840024545*	4.5mm x 45mm	54840028530*	55840028530*	8.5mm x 30mm
54840024550*	55840024550*	4.5mm x 50mm	54840028535*	55840028535*	8.5mm x 35mm
54840024555*	55840024555*	4.5mm x 55mm	54840028540*	55840028540*	8.5mm x 40mm
54840024560*	55840024560*	4.5mm x 60mm	54840028545	55840028545	8.5mm x 45mm
54840024565*	55840024565*	4.5mm x 65mm	54840028550	55840028550	8.5mm x 50mm
54840024570*	55840024570*	4.5mm x 70mm	54840028555	55840028555	8.5mm x 55mm
54840024575*	55840024575*	4.5mm x 75mm	54840028560*	55840028560*	8.5mm x 60mm
54840024580*	55840024580*	4.5mm x 80mm	54840028565*	55840028565*	8.5mm x 65mm
54840025530*	55840025530*	5.5mm x 30mm	54840028570*	55840028570*	8.5mm x 70mm
54840025535	55840025535	5.5mm x 35mm	54840028575*	55840028575*	8.5mm x 75mm
54840025540	55840025540	5.5mm x 40mm	54840028580*	55840028580*	8.5mm x 80mm
54840025545	55840025545	5.5mm x 45mm	–	55840029530*	9.5mm x 30mm
54840025550	55840025550	5.5mm x 50mm	–	55840029535*	9.5mm x 35mm
54840025555	55840025555	5.5mm x 55mm	–	55840029540*	9.5mm x 40mm
54840025560*	55840025560*	5.5mm x 60mm	–	55840029545*	9.5mm x 45mm
54840025565*	55840025565*	5.5mm x 65mm	–	55840029550*	9.5mm x 50mm
54840025570*	55840025570*	5.5mm x 70mm	–	55840029555*	9.5mm x 55mm
54840025575*	55840025575*	5.5mm x 75mm	–	55840029560*	9.5mm x 60mm
54840025580*	55840025580*	5.5mm x 80mm	–	55840029565*	9.5mm x 65mm
54840026530*	55840026530*	6.5mm x 30mm	–	55840029570*	9.5mm x 70mm
54840026535*	55840026535*	6.5mm x 35mm	–	55840029575*	9.5mm x 75mm
54840026540	55840026540	6.5mm x 40mm	–	55840029580*	9.5mm x 80mm
54840026545	55840026545	6.5mm x 45mm	–	55840021540*	10.5mm x 40mm
54840026550	55840026550	6.5mm x 50mm	–	55840021545*	10.5mm x 45mm
54840026555	55840026555	6.5mm x 55mm	–	55840021550*	10.5mm x 50mm
54840026560*	55840026560*	6.5mm x 60mm	–	55840021555*	10.5mm x 55mm
54840026565*	55840026565*	6.5mm x 65mm	–	55840021560*	10.5mm x 60mm
54840026570*	55840026570*	6.5mm x 70mm	–	55840021565*	10.5mm x 65mm
54840026575*	55840026575*	6.5mm x 75mm	–	55840021570*	10.5mm x 70mm
54840026580*	55840026580*	6.5mm x 80mm	–	55840021575*	10.5mm x 75mm
54840027530*	55840027530*	7.5mm x 30mm	–	55840021580*	10.5mm x 80mm
54840027535*	55840027535*	7.5mm x 35mm			
54840027540*	55840027540*	7.5mm x 40mm			
54840027545	55840027545	7.5mm x 45mm			
54840027550	55840027550	7.5mm x 50mm			
54840027555	55840027555	7.5mm x 55mm			
54840027560*	55840027560*	7.5mm x 60mm			
54840027565*	55840027565*	7.5mm x 65mm			

\* Available as Add-Ons

## Product Ordering Information

### INSTRUMENTS (NON-STERILE, MULTI-USE)

Reference	Description
7480741	Adapter Driver QC Shaft
7480745	Lock Sleeve 4.75
7480755	Lock Sleeve 5.5/6.0
7480760	Stabilisation Clip

### INSTRUMENTS (STERILE, DISPOSABLE)

Reference	Description
F04B	Size 3 KyphX® Bone Filler Device
A07A	Kyphon® Mixer
CDS3A	CDS Cement Delivery System
CC02A	2x Cement Cartridges

### INSTRUMENTS (NON-STERILE, DISPOSABLE)

Reference	Description
7480740	Adapter Driver T25 Tip














### CEMENT (STERILE)

Reference	Description
7480724	Fenestrated Screw Cement





## Explanation of Symbols

	Authorised representative in the European Community		Non-sterile
	<b>CAUTION:</b> Federal law (USA) restricts these devices to sale by or on the order of a physician		For US Audiences only
	Consult instructions for use		The device complies with European Directive MDD 93/42/EEC
	Do not reuse		Sterilised by irradiation
	Batch/Lot code		Sterilised by ethylene oxide
	Manufacturer		Use-by date
	Catalogue number		

©2013 Medtronic Sofamor Danek USA, Inc. All rights reserved.

## Important Information on Resin for Fenestrated Screw (FS) Cement High Viscosity, Radiopaque Bone Cement

	<b>Distributed by</b>
<b>ELMDOWN SRL</b> Unipersonale Via Manfredo Camperio 14 20123 Milano ITALIA Tel: +39 02 87166971 Fax: +39 02 87152991	<b>Medtronic Sofamor Danek USA, Inc.</b> 1800 Pyramid Place Memphis, TN 38132 Tel: +1 800 933 2635 (In USA) +1 901 396 3133 (Outside of USA) Fax: +1 901 396 0356

**INSTRUCTIONS FOR USE**  
For the attention of the surgeon. Carefully read all instructions prior to use.

**INDICATION FOR USE**  
FS Cement is a radiopaque acrylic resin to be used in screw augmentation surgical procedures to increase the fixation strength of bone screws. In particular this cement is used to augment pedicle screw fixation in patients diagnosed with osteoporosis.

**STERILE**  
FS Cement is supplied sterile. The powder and the package are sterilised with gamma irradiation. The liquid is sterilised using filtration and is contained in a glass vial. The outside of the glass vial is sterilised with ethylene oxide gas. This device is intended for single use only. Do not re-sterilise. Do not use if package is open or damaged.

**DESCRIPTION**  
FS Cement is a polymethylmethacrylate (PMMA) that contains approximately 30% barium sulphate. It is designed for delivery in a highly viscous state.

**COMPOSITION OF FS CEMENT**

POWDER – 20 g of sterile powder	20 g packet of sterile powder
Methylmethacrylate – Styrene copolymer	68.0% w/w
Barium sulphate	30.0% w/w
Benzoyl peroxide	2.0% w/w
LIQUID – 9,0 g phial of sterile liquid	Gravity displacement
Methylmethacrylate	99.1% w/w
NN-dimethyl-p-toluidine	0.9% w/w
Hydroquinone	75ppm

**CONTRAINDICATIONS**

- Active infections
- Haemorrhagic diathesis
- Fractured or partially resected vertebral bodies
- Severe Osteopenia
- Fever

The use of the device must be carefully considered in presence of hypersensitivity to monomer or to any of the other components of the resin

**WARNINGS**

**NOTE: The use of this product is not recommended in patients that do not exhibit a pathological condition, such as primary or secondary osteoporosis or a tumour, which would impair the ability of the patient to heal using conservative treatment methods**

- The sterile FS Cement is designed and intended for one time use only. DO NOT reuse, reprocess, or re-sterilise this product. This device is designed for single patient use only. Reuse, reprocessing, or re-sterilisation may compromise the properties of the cement and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. Avoid the partition of the product in two or more portions to be used in different moments. This would also be classified as a reuse which could lead to a ratio error between powder and liquid components and loss of sterility with subsequent risk of infection to the patient. Moreover there would be the risk of inadequate cement volume to complete the surgical operation thus potentially lengthening the surgery time. The residual material must be considered surgical waste and therefore it must be eliminated at the end of the surgical procedure. Sterility is assured only if the unit container is not damaged
- Maintain aseptic surgical technique to prevent possible infection including treatment site infection
- Always use live imaging when injecting the material
- Monitor patients carefully for any change in blood pressure during and immediately following the application of bone cement. Adverse patient reactions affecting the cardiovascular system, including Bone Cement Implantation Syndrome (BCIS), have been associated with the use of bone cements. Hypotensive reactions have occurred between 10 and 165 seconds following application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement, especially those potentially at increased risk for peri-operative death, including elderly patients, patients with underlying cardiac or pulmonary compromise, and patients being treated for multiple vertebral body fractures in one procedure
- Caution should be exercised during the mixing of the two components to prevent excessive exposure to the concentrated vapours of the monomer, which may produce irritation of the respiratory tract, eyes, and possibly the liver
- Polymerisation of the bone cement is an exothermic reaction, which occurs while the cement is hardening in situ. The released heat may damage bone or other tissues surrounding the implant. The heat produced in situ has a long term effect, and positioning of the patient should be maintained securely to obtain proper fixation
- The liquid component is a powerful lipid solvent. It should not be allowed to come into contact with rubber or latex gloves. Should contact occur, the gloves may dissolve and tissue damage may occur.

Wearing a second pair of gloves and adherence to mixing instructions may diminish the possibility of hypersensitivity reactions. The mixed bone cement should not make contact with the gloved hand until the bone cement is being tested for the consistency of dough

- Personnel wearing permeable contact lenses should NOT be near or involved in mixing the bone cement; soft contact lenses are permeable
- The handling characteristics of bone cements are affected by operating room conditions, including the room temperature, temperature of the cement components prior to mixing, humidity, the geometry of the mixing apparatus, time spent mixing, and the geometry of the delivery device. Any change in one or more of these conditions can alter the handling characteristics of the bone cement, including the time it takes for bone cement to reach the doughy state (the handling period), the time the bone cement remains in the doughy state (the working time), and the time the bone cement hardens (the setting time). The user must be aware of these factors and adjust technique to account for variability in operating room conditions
- Avoid over-pressurisation of the bone cement and do not insert the bone cement into the cavity of the vertebral body until the bone cement has reached the doughy state. Maintain patient positioning until the end of the polymerisation process. Failure to do so may lead to extravasation of the bone cement beyond the site of its intended application and may or could damage the surrounding tissues
- The safety of the bone cement in pregnant women, breastfeeding women or in children has not been established. Bone cement may adversely affect bone growth and foetal health
- Difficulty in swallowing and blistering of the throat, alleged to be an allergic reaction to the use of bone cement, have also been reported
- DO NOT mix more than one vial of liquid and one packet of powder together at any one time. Never modify the ratios between the liquid and solid components. Doing so could affect bone cement properties, including handling characteristics
- Long-term follow-up is advised for all patients on a regularly scheduled basis. Inadequate filling or unanticipated postoperative events may affect the stability of the bone cement and/or the bone cement-bone interface. A fibrous tissue layer may develop between the cement and the bone, and loosening of the bone cement may occur leading to failure
- FS Cement is provided in finished form with all the necessary components for use. The addition of radiopacifier, e.g. barium; antibiotics; or other drugs or materials to FS Cement is not recommended. The safety and effectiveness of adding such drugs or materials has not been evaluated and may cause patient harm
- Modifying the polymerisation time by either warming or cooling the bone cement and/or associated delivery devices has not been tested and could affect bone cement properties, including handling characteristics
- If bone cement is seen outside of the pedicle or in the circulatory system during the procedure, immediately stop the injection

- PRECAUTIONS**
- Store product at room temperature below 25°C and away from light. Do not use the product if the powder has a yellowish or brownish color or if the liquid is syrupy. These two conditions indicate that the product has not been stored correctly. Do not use monomer if the liquid component shows any sign of thickening or premature polymerisation. Keep the product at a temperature of 23 ± 1°C for a period of 24 hours prior to use
  - Never deliver bone cement into the vertebral body without the use of high quality fluoroscopic guidance capable of visualising movement of bone cement
  - The operator should have specific training and be familiar with the properties, handling characteristics, and application of the bone cement and adhere to the instructions for use
  - Do not use after the expiration date printed on the package. The device may not be safe or effective beyond its expiration date
  - Follow the mixing and handling instructions to avoid contact dermatitis. Strict adherence to the instructions for mixing the powder and liquid components may reduce the incidence of this complication
  - Adequately ventilate the operating room to eliminate as much monomer vapour as possible. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported
  - Dispose of the polymer component in an authorised waste facility. The liquid component can be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal

- ADVERSE EVENTS**
- Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements include:
- Cardiac arrest
  - Cerebrovascular accident
  - Myocardial infarction
  - Pulmonary embolism
- Adhere strictly to good surgical principles and techniques. Deep wound infection is a serious postoperative complication and may require total removal of the embedded cement. Deep wound infection may be latent and not manifest itself even for several years postoperatively.
- Other reported adverse events relevant to the anatomy being treated with acrylic bone cements include:
- Fistula
  - Haematoma
  - Hemorrhage
  - Heterotopic new bone formation
- Extravasation of bone cement potentially resulting in but not limited to:
- Compression or irritation of nerve structures, such as the spinal cord or nerve roots, causing radiculopathy, paresthesia, paraplegia or paralysis and/or;
  - Introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae
  - Pyrexia due to allergy to bone cement
  - Short-term conduction irregularities
  - Thrombophlebitis
  - Transitory fall in blood pressure
  - Migration of hardened bone cement bolus
  - Fracture of the pedicle





THERAPIES

SERVICES & SUPPORT

ENABLING TECHNOLOGIES

#### NOT FOR DISTRIBUTION IN THE USA OR ITS TERRITORIES

The surgical technique shown is for illustrative purposes only. The technique(s) actually employed in each case will always depend upon the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient.

Please see the package insert for the complete list of indications, warnings, precautions, and other important medical information.



**Medtronic Sofamor Danek USA, Inc.**  
1800 Pyramid Place  
Memphis, TN 38132  
Tel: +1 800 933 2635 (In USA)  
+1 901 396 3133 (Outside of USA)  
Fax: +1 901 396 0356

[www.medtronic.eu](http://www.medtronic.eu)

**Medtronic International Trading Sàrl**  
Route du Molliau 31  
Case Postale  
CH-1131 Tolochenaz  
Tel: +41 (0)21 802 70 00  
Fax: +41 (0)21 802 79 00

