

# Thermo Scientific DreamTaq Hot Start DNA Polymerase

Pub. No. MAN0015972

Rev. Date 29 July 2016 (Rev. A.00)

Lot: \_\_\_\_ Expiry Date: \_\_\_\_

## Ordering Information

Catalog No.	DreamTaq Hot Start DNA Polymerase, 5 U/μL	10X DreamTaq Buffer*
EP1701	200 U	1.25 mL
EP1702	500 U	2 × 1.25 mL
EP1703	2500 U	10 × 1.25 mL
EP1704	4 × 2500 U	40 × 1.25 mL

\* includes 20 mM MgCl<sub>2</sub>

Store at **-20°C**

www.thermofisher.com

For Research Use Only. Not for use in diagnostic procedures.

## DESCRIPTION

Thermo Scientific™ DreamTaq™ Hot Start DNA Polymerase is an enhanced hot start *Taq* DNA polymerase optimized for most PCR applications. It ensures higher sensitivity, specificity, and yields compared to conventional hot start *Taq* DNA polymerase. It is capable of amplifying long amplicons such as 6 kb genomic DNA and 20 kb λ DNA.

DreamTaq Hot Start DNA Polymerase combines *Taq* DNA polymerase and a specific antibody that inhibits the DNA polymerase activity at ambient temperatures, thus preventing the amplification of non-specific products. At polymerization temperatures, the antibody molecule is released, rendering the polymerase fully active.

DreamTaq Hot Start DNA Polymerase uses the same reaction set-up and cycling conditions as conventional *Taq* DNA polymerases, but the antibody-based hot start allows the reactions to be set up at room temperature. Because the enzyme is supplied with the optimized DreamTaq buffer, which includes 20 mM MgCl<sub>2</sub>, extensive optimization of reaction conditions is not required.

DreamTaq Hot Start DNA Polymerase generates PCR products with 3'-dA overhangs. The enzyme tolerates dUTP and can incorporate modified nucleotides.

## FEATURES

- High specificity due to antibody based hot start.
- Robust amplification with minimal optimization.
- High yields of PCR products.
- Higher sensitivity compared to conventional hot start *Taq* DNA polymerases.
- Amplification of long targets up to 6 kb from genomic DNA and up to 20 kb from viral DNA.
- Generates 3'-dA overhangs.
- Incorporates dUTP and modified nucleotides.

## APPLICATIONS

- Routine PCR amplification of DNA fragments up to 6 kb from genomic DNA and up to 20 kb from viral DNA.
- RT-PCR.
- Genotyping.
- Generation of PCR products for TA cloning.

## CONCENTRATION

5 U/μL

## DEFINITION OF ACTIVITY UNIT

One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 minutes at 74°C.

## 10X DREAMTAQ BUFFER

DreamTaq Buffer is a proprietary formulation, which contains KCl and (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> at a ratio optimized for robust performance of DreamTaq Hot Start DNA Polymerase in PCR. DreamTaq Buffer also includes MgCl<sub>2</sub> at a concentration of 20 mM.

## INHIBITION AND INACTIVATION

- Inhibitors: Ionic detergents (deoxycholate, sarkosyl, and SDS) at concentrations higher than 0.06, 0.02, and 0.01%, respectively.
- Inactivated by phenol/chloroform extraction.

## PROTOCOL

To set up parallel reactions and to minimize the possibility of pipetting errors, prepare a PCR master mix by mixing water, buffer, dNTPs, primers, and DreamTaq Hot Start DNA Polymerase. Prepare sufficient master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes, then add template DNA.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. For each 50-μL reaction, add the following components into a thin-walled PCR tube:

<b>10X DreamTaq Buffer*</b>	5 μL
<b>dNTP Mix, 2 mM each (#R0241)</b>	5 μL (0.2 mM of each)
<b>Forward primer</b>	0.1–1.0 μM
<b>Reverse primer</b>	0.1–1.0 μM
<b>Template DNA</b>	10 pg–1 μg
<b>DreamTaq Hot Start DNA Polymerase</b>	1.25 U
<b>Water, nuclease-free (#R0581)</b>	to 50 μL
<b>Total volume</b>	50 μL

\*10X DreamTaq Buffer contains 20 mM MgCl<sub>2</sub>, which is optimal for most applications. If further optimization is required, additional MgCl<sub>2</sub> can be added to the master mix. The volume of water should be reduced accordingly.

Volumes of 25 mM MgCl<sub>2</sub> (#R0971), required for specific final MgCl<sub>2</sub> concentration:

Final concentration of MgCl <sub>2</sub>	2 mM	2.5 mM	3 mM	4 mM
Volume of 25 mM MgCl <sub>2</sub> to be added for 50-μL reaction	0 μL	1 μL	2 μL	4 μL

3. Gently vortex the samples and briefly centrifuge.
4. When using a thermal cycler that does not contain a heated lid, overlay the reaction mixture with 25 μL of mineral oil.

5. Place the reactions in a thermal cycler. Perform PCR using the recommended thermal cycling conditions outlined below:

Step	Temperature, °C	Time	Number of cycles
Initial denaturation	95	1–3 min	1
Denaturation	95	30 s	25–40
Annealing	T <sub>m</sub>	30 s	
Extension*	72	1 min	1
Final Extension	72	5–15 min	

\* The recommended extension step is 1 minute for PCR products up to 2 kb. For longer products, the extension time should be prolonged by 1 minute/kb.

## GUIDELINES FOR PREVENTING CONTAMINATION OF PCR REACTION

During PCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. Follow the general recommendations below to lower the risk of contamination.

- Prepare your DNA sample, set up the PCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up PCR mixtures in a laminar flow cabinet equipped with an UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive displacement pipettes, or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.
- Use PCR-certified reagents, including high quality water (e.g., Water, nuclease-free, #R0581).
- Always perform “no template control” (NTC) reactions to check for contamination.

DreamTaq Hot Start DNA Polymerase incorporates dUTP; therefore, you can control carry-over contamination using Uracil-DNA Glycosylase (#EN0361).

## GUIDELINES FOR PRIMER DESIGN

Use special design software or follow the general recommendations for PCR primer design as outlined below to design optimal primers:

- Use PCR primers that are 15–30 nucleotides long.
- Optimal GC content of the primer is 40–60%. Ideally, C and G nucleotides should be distributed uniformly along the primer.
- Avoid placing more than three G or C nucleotides at the 3'-end to lower the risk of non-specific priming.
- If possible, the primer should terminate with a G or C at the 3'-end.

- Avoid self-complementary primer regions, and complementarities between the primers and direct primer repeats to prevent hairpin formation and primer dimerization.
- Check for possible sites of undesired complementarity between primers and template DNA.
- When designing degenerate primers, place at least 3 conserved nucleotides at the 3'-end.
- Differences in melting temperatures (T<sub>m</sub>) between the two primers should not exceed 5°C.

## ESTIMATION OF PRIMER MELTING TEMPERATURE

For primers containing less than 25 nucleotides, the approximate melting temperature (T<sub>m</sub>) can be calculated using the following equation:

$$T_m = 4(G + C) + 2(A + T),$$

where G, C, A, T represent the number of respective nucleotides in the primer.

If the primer contains more than 25 nucleotides, we recommend using specialized computer programs to account for interactions of adjacent bases, effect of salt concentration, etc.

## COMPONENTS OF THE REACTION MIXTURE

### Template DNA

Optimal amount of template DNA for a 50- $\mu$ L reaction volume is 1 pg–1 ng for both plasmid and phage DNA, and 100 pg–1  $\mu$ g for genomic DNA. Higher amounts of template increase the risk of non-specific PCR products. Lower amounts of template reduce the accuracy of the amplification.

All routine DNA purification methods are suitable for template preparation; e.g., Thermo Scientific™ GeneJET™ Genomic DNA Purification Kit (#K0721) or GeneJET Plasmid Miniprep Kit (#K0502). Trace amounts of certain agents used for DNA purification, such as phenol, EDTA, and proteinase K, can inhibit DNA polymerases. Ethanol precipitation and repeated washes of the DNA pellet with 70% ethanol normally removes trace contaminants from DNA samples.

### MgCl<sub>2</sub> concentration

DreamTaq Hot Start DNA Polymerase is provided with an optimized 10X DreamTaq Buffer, which includes MgCl<sub>2</sub> at a concentration of 20 mM. A final MgCl<sub>2</sub> concentration of 2 mM is generally ideal for PCR. MgCl<sub>2</sub> concentration can be further increased up to 4 mM by the addition of 25 mM MgCl<sub>2</sub> (#R0971).

If the DNA samples contain EDTA or other metal chelators, Mg<sup>2+</sup> ion concentration in the PCR mixture should be increased accordingly (1 molecule of EDTA binds 1 Mg<sup>2+</sup>).

## dNTPs

The recommended final concentration of each dNTP is 0.2 mM. In certain PCR applications, higher dNTP concentrations may be necessary. It is essential to have equal concentrations of all four nucleotides (dATP, dCTP, dGTP and dTTP) in the reaction mixture.

To obtain a 0.2 mM concentration of each dNTP in the PCR mixture, refer to the table below.

Volume of PCR mixture	dNTP Mix, 2 mM each (#R0241)	dNTP Mix, 10 mM each (#R0191)	dNTP Mix, 25 mM each (#R1121)
50 $\mu$ L	5 $\mu$ L	1 $\mu$ L	0.4 $\mu$ L
25 $\mu$ L	2.5 $\mu$ L	0.5 $\mu$ L	0.2 $\mu$ L
20 $\mu$ L	2 $\mu$ L	0.4 $\mu$ L	0.16 $\mu$ L

Use 200  $\mu$ M of each dNTP. dUTP or dITP can be added up to 200  $\mu$ M. For longer amplicons, a lower dUTP concentration (20–100  $\mu$ M) may be required for high yields.

## Primers

The recommended concentration range of the PCR primers is 0.1–1  $\mu$ M. Excessive primer concentrations increase the probability of mispriming and generation of non-specific PCR products.

For degenerate primers and primers used for long PCR, we recommend higher primer concentrations in the range of 0.3–1  $\mu$ M.

## CYCLING PARAMETERS

### Initial DNA denaturation and enzyme activation

DreamTaq Hot Start DNA polymerase is inactive at room temperature during the reaction set up and is activated during the 1–3 minute initial denaturation/enzyme activation step.

It is essential to completely denature the template DNA at the beginning of the PCR run to ensure efficient utilization of the template during the first amplification cycle. If the GC content of the template is 60% or less, an initial 1–3 minute denaturation at 95°C is sufficient. For GC-rich templates this step can be prolonged.

### Denaturation

A DNA denaturation time of 30 seconds per cycle at 95°C is normally sufficient. For GC-rich DNA templates, this step can be prolonged to 3–4 minutes. DNA denaturation can also be enhanced by the addition of 5–10% glycerol, 5% DMSO, 1% formamide, or 1–1.5 M betaine. The melting temperature of the primer-template complex decreases significantly in the presence of these reagents. Therefore, the annealing temperature has to be adjusted accordingly.

Note that higher than 10% DMSO or 5% formamide in the reaction mix inhibit DNA polymerases. Therefore, it may be necessary to increase the amount of the enzyme in the reaction if these additives are used.

## Primer annealing

The annealing temperature should be equal to the melting temperature (T<sub>m</sub>) of the primers. Annealing for 30 seconds is normally sufficient. If non-specific PCR products appear, the annealing temperature should be optimized stepwise in 1–2°C increments. When additives that change the melting temperature of the primer-template complex are used (glycerol, DMSO, formamide and betaine), the annealing temperature must also be adjusted.

## Extension

The optimal extension temperature for DreamTaq Hot Start DNA Polymerase is 70–75°C. The recommended extension step is 1 minute at 72°C for PCR products up to 2 kb. For longer products, the extension time should be increased by 1 minute/kb. For amplification of templates >6 kb, we recommend reducing the extension temperature to 68°C.

## Number of cycles

The number of cycles may vary depending on the amount of template DNA in the PCR mixture and the expected PCR product yield.

If less than 10 copies of the template is present in the reaction, about 40 cycles are required. For higher template amounts, 25–35 cycles are sufficient.

## Final extension

After the last cycle, we recommend incubating the PCR mixture at 72°C for an additional 5–15 minutes to fill in any possible incomplete reaction products. If the PCR product will be cloned into TA vectors such as the Thermo Scientific™ InsTAclone™ PCR Cloning Kit (#K1213), the final extension step may be prolonged to 15 minutes to ensure the complete 3'-dA tailing of the PCR product. If the PCR product will be used for cloning using Thermo Scientific™ CloneJET™ PCR Cloning Kit (#K1231), the final extension step can be omitted.

## TROUBLESHOOTING

For troubleshooting, visit [www.thermofisher.com](http://www.thermofisher.com).

## CERTIFICATE OF ANALYSIS

### Endodeoxyribonuclease Assay

No detectable conversion of supercoiled plasmid DNA to a nicked form was observed.

### Residual Activity Assay


No detectable extension of labeled double stranded oligonucleotide with 5'-overhangs after incubation in the presence of dNTPs.

### E. coli DNA Assay

No detectable E.coli DNA was observed.

### Functional Assay

Performance in PCR is tested by the amplification of a 594 bp and 7.5 kb fragments of human genomic DNA.

Quality authorized by:  Jurgita Zilinskiene

### LIMITED USE LABEL LICENSE No. 593: Newcastle License for Modified DNA Polymerase

Notice to Purchaser: This product is licensed under patents owned by University of Newcastle upon Tyne

### LIMITED USE LABEL LICENSE No. 599: Internal Research and Development Use Only

Notice to Purchaser: The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackaging, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights.

For information on purchasing a license for uses other than internal research and development purposes, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

### PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

### For Research Use Only. Not for use in diagnostic procedures.

Refer to [www.thermofisher.com/support](http://www.thermofisher.com/support) for the Safety Data Sheets.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

## Safety Data Sheet

Safety Data Sheet according to Regulation (EC) No. 1907/2006 (REACH)  
Classification according to Regulation (EC) No. 1272/2008 [CLP]

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### Product identifier

<b>Product code</b>	C4221
<b>Product name</b>	10X DreamTaq Buffer
<b>Chemical Name</b>	Not Applicable
<b>REACH registration number</b>	No registration number is given yet for this substance / substances in this mixture since the annual import quantity is less than one tonnage per annum or the transition period for its registration according to Article 23 of REACH has not yet expired.

#### Relevant identified uses of the substance or mixture and uses advised against

<b>Relevant identified uses</b>	For research use only. Not for use in diagnostic procedures
<b>Use Description Code</b>	SU22 - Professional uses: Public domain (administration, education, entertainment, services, craftsmen), PROC15 - Use as laboratory reagent, PC21 - Laboratory chemicals, SU24 - Scientific research and development
<b>Uses advised against</b>	Not for consumer use.

#### Details of the supplier of the safety data sheet

##### Manufacturer / Supplier

Thermo Fisher Scientific Baltics UAB  
V.Graiciuno 8  
LT-02241 Vilnius  
Lithuania  
Tel.: +370 5 2602131  
Fax.: +370 5 2602142

LIFE TECHNOLOGIES EUROPE BV  
KWARTSWEG 2  
2665 NN BLEISWIJK  
NETHERLANDS  
31-(0)180 392 400  
Email: MSDS@lifetech.com

Life Technologies Limited  
3 Fountain Drive  
Inchinnan Business Park  
Paisley  
PA4 9RF, UK  
+44 (0)141 814 6100

**24 hour Emergency Response for Hazardous Materials [or Dangerous Goods] Incident. Spill, Leak, Fire, Exposure, or Accident. Call CHEMTREC** Within the USA + Canada: 1-800-424-9300 and 1-703-527-3887  
Outside the USA + Canada: 1-703-741-5970

#### Country Specific Emergency Number (if available):

CHEMTREC Ireland (Dublin)	+(353)-19014670 (Greeting Language: English and Irish)
CHEMTREC UK (London)	+(44)-870-8200418 (Greeting Language: English)

**Classification of the substance or mixture**

**Classification according to Regulation (EC) No. 1272/2008 [CLP]**

**Physical hazards**

Not Hazardous

**Health hazards**

Not Hazardous

**Environmental hazards**

Not Hazardous

**Additional information**

Not Applicable

**Label elements**

**Labelling according to Regulation (EC) No 1272/2008 [CLP]**

**Hazard pictograms**

None

**Signal Word**

None

**Hazard Statements**

Not Applicable

**EU Specific Hazard Statements**

Not Applicable

**Precautionary Statements**

**Prevention**

Not Applicable

**Response**

Not Applicable

**Storage**

Not Applicable

**Disposal**

Not Applicable

**Other hazards**

Contains a known or suspected endocrine disruptor

### SECTION 3: Composition/information on ingredients

We recommend handling all chemicals with caution.

Component	CAS No	EINECS-No.	Weight-%	REACH registration number	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Nonylphenol ethoxylates 9016-45-9 (<0.9)	9016-45-9	500-024-6	<0.9	-	Acute Tox. 4 - H302 Eye Dam. 1 - H318 Chronic Aquatic 2 - H411

### SECTION 4: First aid measures

#### Description of first aid measures

<b>Skin contact</b>	Rinse skin with water. Immediate medical attention is not required.
<b>Eye contact</b>	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
<b>Ingestion</b>	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice.
<b>Inhalation</b>	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. Consult a physician if necessary.
<b>Notes to Physician</b>	Treat symptomatically.

#### Most important symptoms and effects, both acute and delayed

Not Applicable

#### Indication of any immediate medical attention and special treatment needed

None.

### SECTION 5: Firefighting measures

#### Extinguishing media

**Suitable extinguishing media**  
**Unsuitable extinguishing media**

Water spray. Carbon dioxide (CO<sub>2</sub>). Foam. Dry chemical.  
No information available.

#### Special hazards arising from the substance or mixture

Not known

#### Protective equipment and precautions for firefighters

Standard procedure for chemical fires.

## SECTION 6: Accidental release measures

### **Personal precautions, protective equipment and emergency procedures**

Ensure adequate ventilation

Always wear recommended Personal Protective Equipment

Use personal protection equipment

See section 8 for more information

### **Environmental precautions**

Prevent product from entering drains. Do not allow material to contaminate ground water system.

### **Methods and material for containment and cleaning up**

Soak up with inert absorbent material.

### **Reference to other sections**

See section 8 for more information.

## SECTION 7: Handling and storage

### **Precautions for safe handling**

Use personal protective equipment as required. No special handling advices are necessary.

### **Conditions for safe storage, including any incompatibilities**

Keep in a dry, cool and well-ventilated place. Keep in properly labelled containers.

### **Specific end use(s)**

For research use only. Not for use in diagnostic procedures.

## SECTION 8: Exposure controls/personal protection

### Control parameters

Chemical Name	EU OEL (TWA)	EU OEL (STEL)	EU Skin Notation
Nonylphenol ethoxylates 9016-45-9	None	None	None

Chemical Name	Austria	Belgium (TWA)	Czech Republic
Nonylphenol ethoxylates 9016-45-9	None	None	None

Chemical Name	Denmark (TWA)	Finland OEL (TWA)	France OEL (VME)
Nonylphenol ethoxylates 9016-45-9	None	None	None

Chemical Name	Germany OEL (TWA)	Ireland (TWA)	Italy OEL (TWA)
Nonylphenol ethoxylates 9016-45-9	None	None	None

Chemical Name	Lithuania OEL (TWA)	Netherlands OEL (MAC)	Norway
Nonylphenol ethoxylates 9016-45-9	None	None	None

Chemical Name	Poland	Portugal	Spain OEL (TWA)
Nonylphenol ethoxylates 9016-45-9	None	None	None

Chemical Name	Sweden - Occupational Exposure Limits - TLVs (LLVs)	Switzerland	United Kingdom
Nonylphenol ethoxylates 9016-45-9	None	None	None

**Engineering Measures**      Ensure adequate ventilation, especially in confined areas.

### Exposure controls

#### **Personal protection equipment**

**Respiratory protection**      In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

**Hand protection**              Wear suitable gloves Glove material: Compatible chemical-resistant gloves.

**Eye protection**                Tight sealing safety goggles.

**Skin and Body Protection**    Wear suitable protective clothing.

**Hygiene Measures**            Handle in accordance with good industrial hygiene and safety practice.

#### **Environmental exposure controls**

Prevent product from entering drains. Do not allow material to contaminate ground water system.

## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

<b>Appearance</b>	liquid	
<b>Odour</b>	No data available	
<b>Odour Threshold</b>	No data	
<b>Molecular Weight</b>	No data	
<b>pH</b>	No data available	
<b>Melting point / melting range</b>	°C No data	°F No data
<b>Boiling point / boiling range</b>	°C No data	°F No data
<b>Flash point</b>	°C No data	°F No data
<b>Autoignition Temperature</b>	°C No data	°F No data
<b>Decomposition temperature</b>	°C No data	°F No data
<b>Evaporation rate</b>	No data	
<b>Flammability (solid, gas)</b>	No data available	
<b>Upper explosion limit</b>	No data	
<b>Lower explosion limit</b>	No data	
<b>Vapour Pressure</b>	No data	
<b>Vapour density</b>	No data	
<b>Relative density</b>	No data	
<b>Specific gravity</b>	No data	
<b>Solubility</b>	No data available	
<b>Partition coefficient: n-octanol/water</b>	No data available	
<b>Viscosity</b>	No data	
<b>Explosive properties</b>	No data	
<b>Oxidising properties</b>	No data	

### Other information

No data available.

## SECTION 10: Stability and reactivity

<b>Reactivity</b>	None known.
<b>Chemical stability</b>	Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	Hazardous reaction has not been reported.
<b>Conditions to avoid</b>	No information available.
<b>Incompatible materials</b>	No dangerous reaction known under conditions of normal use.
<b>Hazardous decomposition products</b>	No data available.



## SECTION 11: Toxicological information

### Information on toxicological effects

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Nonylphenol ethoxylates	= 1310 mg/kg (Rat) = 2590 mg/kg (Rat) = 1300 mg/kg (Rat) = 1410 µL/kg (Rat)	No data available	No data available

### Principal Routes of Exposure

**Skin corrosion/irritation** Data are conclusive but insufficient for classification

**Serious eye damage/irritation** Data are conclusive but insufficient for classification

**Respiratory or skin sensitisation** Data are conclusive but insufficient for classification

**Specific target organ toxicity (STOT) – single exposure** Data are conclusive but insufficient for classification

**Specific target organ toxicity (STOT) – repeated exposure** Data are conclusive but insufficient for classification

**Carcinogenicity** Data are conclusive but insufficient for classification

**Germ cell mutagenicity** Data are conclusive but insufficient for classification

**Reproductive Toxicity** Data are conclusive but insufficient for classification

**Aspiration Hazard** Data are conclusive but insufficient for classification

## SECTION 12: Ecological information

### Ecotoxicity

Hazardous to the Aquatic Environment.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Nonylphenol ethoxylates	No data available	No data available	No data available	No data available	No data available

**Persistence and degradability** No information available.

**Bioaccumulative potential** No information available.

### Results of PBT and vPvB assessment

No information available.

### Other adverse effects

Contains a known or suspected endocrine disruptor.

Chemical Name	EU - Endocrine Disruptors Candidate List
Nonylphenol ethoxylates	Group III Chemical; Group III Chemical (9<EO<19); Group III Chemical (with EO>19); Group III Chemical (with EO<9) Group III Chemical

## SECTION 13: Disposal considerations

### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

Do not empty into drains

Do not dispose of waste into sewer.

## SECTION 14: Transport information

### IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

<b>UN number</b>	Not Applicable
<b>UN proper shipping name</b>	Not Applicable
<b>Transport hazard class(es)</b>	Not Applicable
<b>Packing group</b>	Not Applicable

### Environmental hazards

Not Applicable

### Special precautions for user

Not Applicable

### Transport in bulk according to Annex II of MARPOL and the IBC Code

Not Applicable.

## SECTION 15: Regulatory information

### Safety, health and environmental regulations/legislation specific for the substance or mixture

None.

### Substances of Very High Concern

Chemical Name	Weight-%	EU - REACH (1907/2006) - Article 59(1) - Candidate List of Substances for Eventual Inclusion in Annex XIV
Nonylphenol ethoxylates	<0.9	Reason for inclusion Endocrine disrupting properties, Article 57f - environment

### Substance subject to authorisation per REACH Annex XIV

None

Chemical Name	Weight-%	Substance subject to authorisation per REACH Annex XIV
Nonylphenol ethoxylates	<0.9	Intrinsic properties: Endocrine disrupting properties (Article 57(f) - environment) Application date: July 4, 2019 Sunset date: January 4, 2021 Exempted uses: None

### Restricted substances under EC 1907/2006, Annex XVII

Revision date 15-Mar-2019  
Product code C4221

Page 8 / 10  
Product name 10X DreamTaq Buffer

None.

Chemical Name	Weight-%	EU - REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances
Nonylphenol ethoxylates	<0.9	Use restricted. See item 46[b]. Use restricted. See item 46a.

#### Substances listed under Annex I of Regulation (EC) No 689/2008

None.

#### Restricted substances under Annex V of Regulation (EC) No 689/2008

None.

#### Substances under Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC

None.

#### German Water hazard classes (Wassergefährdungsklassen)

Chemical Name	Weight-%	Water hazard class (WGK)
Nonylphenol ethoxylates	<0.9	hazard class 3 - highly hazardous to water hazard class 2 - obviously hazardous to water

#### Other International Inventories

Chemical Name	EINECS (European Union)	ELINCS (European List of Notified Chemical Substances)	ENCS (Japan)	PICCS (Philippines)
Nonylphenol ethoxylates	-	-	Listed	Listed

Chemical Name	AICS (Australia)	South Korea (KECL)	Canada (DSL)	NDSL
Nonylphenol ethoxylates	Listed	Listed	Listed	-

#### Chemical safety assessment

No Chemical safety assessment has been carried out.

## SECTION 16: Other information

<b>Reason for revision</b>	Update according to Commission Regulation (EU) No 830/2015
<b>Revision number</b>	2
<b>Revision date</b>	15-Mar-2019

### References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC)

#### 1272/2008 [CLP]:

Not classified

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

## CERTIFICATE OF ANALYSIS

**EN0581 Exonuclease I (ExoI)**

**Packaging Lot: 01016122**

**Expiry Date: 31.10.2024 (DD.MM.YYYY)**

**Storage: at -20±5°C**

### Filling lots for components in package:

Lot	Quantity	Description
00988399	4 ku	Exonuclease I (ExoI)
00977842	1 mL	10X React Buffer for Exo I

### QUALITY CONTROL

Parameter	Method	Requirement	Result
Unit definition	One unit of the enzyme catalyzes the release of 10 nmol of acid soluble nucleotides in 30 min at 37 °C.	20 U/μL ± 10%	Conforms
Endodeoxyribonucleases (nicking activity)	Incubation of supercoiled plasmid DNA with enzyme.	Not detectable	Conforms
Single-stranded Endodeoxyribonucleases	Incubation of circular single-stranded DNA with enzyme.	Not detectable	Conforms
Double-stranded Exodeoxyribonucleases	Incubation of double-stranded DNA fragments with enzyme.	Not detectable	Conforms
Ribonucleases	Incubation of [3H]-RNA with enzyme.	Not detectable	Conforms

#### ISO CERTIFICATION

Manufactured by Thermo Fisher Scientific Baltics UAB, in compliance with ISO 9001 and ISO 13485 certified quality management system.

Quality authorized by QC: **J. Žilinskienė**



# Exonuclease I (Exo I)

Catalog Number EN0581, EN0582

Pub. No. MAN0012007 Rev. C.00



**WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Sheets (SDSs) are available from [thermofisher.com/support](http://thermofisher.com/support).

## Contents and storage

Cat. No.	Contents	Amount	Storage
EN0581	Exonuclease I (Exo I)	4000 U, 20 U/ $\mu$ L	-25 °C to -15 °C
	10X Reaction Buffer	1 mL	
EN0582	Exonuclease I (Exo I)	20000 U, 20 U/ $\mu$ L	
	10X Reaction Buffer	5 x 1 mL	

## Description

Exonuclease I (Exo I) degrades single-stranded DNA in a 3'→5' direction, releasing deoxyribonucleoside 5'-monophosphates in a stepwise manner and leaving 5'-terminal dinucleotides intact. It does not cleave DNA strands with terminal 3'-OH groups blocked by phosphoryl or acetyl groups (1).

## Applications

- Primer removal from PCR mixtures:
  - prior to PCR product sequencing (2),
  - for one-tube “megaprimer” PCR mutagenesis (3).
- Removal of single-stranded DNA containing a 3'-hydroxyl terminus from nucleic acid mixtures.
- Assay for the presence of single-stranded DNA with a 3'-hydroxyl terminus (4).

## Source

*E. coli* cells with a cloned *E. coli sbcB* gene.

## Definition of Activity Unit

One unit of the enzyme catalyzes the release of 10 nmol of acid soluble nucleotides in 30 min at 37 °C.

Enzyme activity is assayed in the following mixture: 67 mM glycine-KOH (pH 9.5), 6.7 mM MgCl<sub>2</sub>, 1 mM DTT and 0.17 mg/mL single-stranded [<sup>3</sup>H]-DNA.

## Storage Buffer

The enzyme is supplied in: 20 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM DTT and 50 % (v/v) glycerol.

## 10X Reaction Buffer

670 mM glycine-KOH (pH 9.5 at 25 °C), 67 mM MgCl<sub>2</sub>, 10 mM DTT.

## Inhibition and Inactivation

- Inhibitors: 20 % (w/v) PEG 8000 (5).
- Inactivated by heating at 80 °C for 15 min.

## Note

The enzyme is not suitable for removing 3'-overhangs of dsDNA.

## Protocol for PCR product clean-up prior to sequencing

The clean-up reaction removes unincorporated primers and degrades unincorporated nucleotides. The resulting PCR product is ready to use for sequencing without additional purification, e.g., using column purification kits.

1. Prepare the following reaction mixture:

Components	Volume
PCR mixture (directly after completion of PCR)	5 $\mu$ L
Exonuclease I	0.5 $\mu$ L (10 U)
Thermo Scientific™ FastAP™ Thermosensitive Alkaline Phosphatase (#EF0651)	1 $\mu$ L (1 U)

- Mix well and incubate at 37 °C for 15 min.
- Stop the reaction by heating the mixture at 85 °C for 15 min.

### Note

- Up to 5  $\mu$ L of purified PCR products can be used directly for DNA sequencing without further purification.
- For reliable sequencing results there should not be nonspecific PCR products.
- The protocol may be applied for clean-up of PCR products, generated by any thermophilic DNA polymerase or polymerase mix.
- The procedure is not recommended for downstream cloning applications.

### Reference

- Lehman, I.R., Nussbaum A.L., The deoxyribonucleases of *Escherichia coli*. V. On the specificity of exonuclease I (phosphodiesterase), J. Biol. Chem., 239, 2628-2636, 1964.
- Werle, E., et al., Convenient single-step, one tube purification of PCR products for direct sequencing, Nucleic Acids Res., 22, 4354-4355, 1994.
- Nabavi S., Nazar R.N., Simplified one tube "megaprimer" polymerase chain reaction mutagenesis, Anal Biochem., 2, 346-348, 2005.
- Rosamond, J., et al., Modulation of the action of the recBC enzyme of *Escherichia coli* K-12 by Ca<sup>2+</sup>, J. Biol. Chem., 254, 8646-8652, 1979.
- Sasaki, Y., Miyoshi, D. and Sugimoto, N., Regulation of DN nucleases by molecular crowding., Nucleic Acids Res., 35, 4086-4093, 2007.

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania  
For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](http://thermofisher.com/symbols-definition).

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

[thermofisher.com/support](http://thermofisher.com/support) | [thermofisher.com/askaquestion](http://thermofisher.com/askaquestion)  
[thermofisher.com](http://thermofisher.com)

**ThermoFisher**  
SCIENTIFIC

## CERTIFICATE OF ANALYSIS

**EF0651      FastAP Thermosensitive AP**

**Packaging Lot: 01017842**

Expiry Date: 31.08.2023 (DD.MM.YYYY)

Storage: at -20±5°C

### Filling lots for components in package:

Lot	Quantity	Description
00991580	1 ku	FastAP Thermosensitive AP
01003773	2 x 1.5 mL	10X Buffer for fastAP

### QUALITY CONTROL

Parameter	Method	Requirement	Result
Unit definition	One unit is the amount of the enzyme required to dephosphorylate 5'-termini of 1 µg of linearized pUC57 DNA in 10 min at 37 °C in FastAP buffer.	1.0 U/µL ± 15%	Conforms
Endodeoxyribonucleases (nicking activity)	Incubation of supercoiled plasmid DNA with enzyme.	Not detectable	Conforms
Ribonucleases	Incubation of [3H]-RNA with enzyme.	Not detectable	Conforms
Endo- and Exodeoxyribonucleases	Incubation of single stranded and double stranded radiolabeled oligonucleotides with enzyme.	Not detectable	Conforms

### ISO CERTIFICATION

Manufactured by Thermo Fisher Scientific Baltics UAB, in compliance with ISO 9001 and ISO 13485 certified quality management system.

Quality authorized by QC: **J. Žilinskienė**





# Dephosphorylation of Proteins

This protocol is for the Dephosphorylation of Proteins

## Reaction Mixture:

1X Thermo Scientific™ FastAP™ reaction buffer,  
0.1 - 0.2 mg/mL of phosphoprotein,  
10 u of FastAP Thermosensitive Alkaline Phosphatase.  
Incubate at 37 °C for 1 hour.

**For example:** if you are doing a 20 µL reaction setup you need 2 µL 10x FastAP buffer, 2-4 µg of protein (to be in the range of 0.1 - 0.2 mg/mL) and 10 U of FastAP Thermosensitive Alkaline Phosphatase (1 u/µL).

## Note:

- The reaction can be stopped by addition of a final concentration of 50 mM EDTA (#R1021) or by addition of a final concentration of 10 mM sodium orthovanadate (Na<sub>3</sub>VO<sub>4</sub>).
- The optimal incubation time and the enzyme concentration must be determined experimentally for each substrate

[www.thermoscientific.com/onebio](http://www.thermoscientific.com/onebio)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

### Europe

Customer Service  
cs.molbio.eu@thermofisher.com

Technical Support  
ts.molbio.eu@thermofisher.com

Tel 00800 222 00 888  
Fax 00800 222 00 889

### United States

Customer Service  
cs.molbio@thermofisher.com

Technical Support  
ts.molbio@thermofisher.com

Tel 800 235 9880  
Fax 800 292 6088

### Canada

Customer Service  
cs.molbio@thermofisher.com

Technical Support  
ts.molbio@thermofisher.com

Tel 800 340 9026  
Fax 800 472 8322

**Thermo**  
S C I E N T I F I C

Part of Thermo Fisher Scientific



# SAFETY DATA SHEET

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### Identification of the substance or mixture

**Product code** 4482777  
**Product name** TAQMAN QSY PROBE 6,000 PMOLES

### Company/undertaking identification

Life Technologies  
5791 Van Allen Way  
PO Box 6482  
Carlsbad, CA 92008  
+1 760 603 7200

LIFE TECHNOLOGIES LIMITED  
3 FOUNTAIN DRIVE  
INCHINNAN BUSINESS PARK  
PAISLEY, PA4 9RF  
SCOTLAND  
44-141 814-6100

### 24 hour Emergency Response:

866-536-0631  
301-431-8585  
Outside of the U.S. ++1-301-431-8585

**For research use only. Not intended for human or animal diagnostic or therapeutic uses.**

## SECTION 2: Hazards identification

### In accordance with local and national regulations

In accordance with local and national regulations Regulation (EC) No 1272/2008 Classification according to Directive 67/548/EEC or 1999/45/EC

### GHS - Classification

#### Signal word

Not Hazardous

#### Health hazards

Not Hazardous

#### Physical hazards

Not Hazardous

**Revision date** 23-Jan-2015  
**Product code** 4482777

**Product name** TAQMAN QSY PROBE 6,000 PMOLES **Page** 1 / 5

## European Union

### EU Specific Hazard Statements

#### **R-phrase(s)**

None

#### **S phrases**

None

### Principle Routes of Exposure/

#### Potential Health effects

**eyes**

May cause eye irritation with susceptible persons.

**Skin**

May cause skin irritation in susceptible persons.

**Inhalation**

May be harmful by inhalation.

**INGESTION**

May be harmful if swallowed.

#### Specific effects

**Carcinogenic effects**

None

**Mutagenic effects**

None

**Reproductive toxicity**

None

**Sensitisation**

None

**Target Organ Effects**

No known effects under normal use conditions

### SECTION 3: Composition/information on ingredients

The product contains no substances which at their given concentration, are considered to be hazardous to health. We recommend handling all chemicals with caution.

### SECTION 4: First aid measures

**Skin contact**

Rinse cautiously with water for several minutes. If symptoms occur, obtain medical advice.

**Eye contact**

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. If symptoms persist, call a doctor.

**INGESTION**

Never give anything by mouth to an unconscious person. If symptoms persist, call a doctor. Do not induce vomiting without medical advice.

**Inhalation**

Remove to fresh air. If symptoms persist, call a doctor. If not breathing, give artificial respiration.

**Notes to Physician**

Treat symptomatically.

### SECTION 5: Firefighting measures

**Suitable extinguishing media**

Water spray. Carbon dioxide (CO<sub>2</sub>). Foam. Dry chemical.

**Special protective equipment for firefighters**

Wear self-contained breathing apparatus and protective suit.

### SECTION 6: Accidental release measures

Revision date 23-Jan-2015  
Product code 4482777

Page 2 / 5  
Product name TAQMAN QSY PROBE 6,000 PMOLES

**Personal precautions** Use personal protection equipment.  
**Methods for cleaning up** Soak up with inert absorbent material.

### Environmental precautions

Prevent further leakage or spillage if safe to do so.

See Section 12 for more information.

## SECTION 7: Handling and storage

**Handling** Always wear recommended Personal Protective Equipment. Wear personal protective equipment.  
**Storage** Keep in a dry, cool and well-ventilated place.

## SECTION 8: Exposure controls/personal protection

### Exposure Limits

We are not aware of any national exposure limit.

**Engineering measures** Ensure adequate ventilation, especially in confined areas.

### Personal protective equipment

Personal Protective Equipment requirements are dependent on the user institution's risk assessment and are specific to the risk assessment for each laboratory where this material may be used.

**Respiratory protection** In case of insufficient ventilation, wear suitable respiratory equipment.

**Hand protection** Impervious gloves.  
**Eye protection** Safety glasses with side-shields.  
**Skin and body protection** Lightweight protective clothing.  
**Hygiene measures** Handle in accordance with good industrial hygiene and safety practice.

**Environmental exposure controls** Prevent product from entering drains.

## SECTION 9: Physical and chemical properties

### General information

**Form** liquid  
**Appearance** No information available  
**Odour** No information available  
**Boiling point / boiling range** °C No data available °F No data available

Revision date 23-Jan-2015  
Product code 4482777

Page 3 / 5  
Product name TAQMAN QSY PROBE 6,000 PMOLES

<b>Melting point / melting range</b>	°C No data available	°F No data available
<b>flash point</b>	°C No data available	°F No data available
<b>Autoignition temperature</b>	°C No data available	°F No data available
<b>Oxidising properties</b>	No information available	
<b>Water solubility</b>	soluble	

#### SECTION 10: Stability and reactivity

<b>Stability</b>	Stable under normal conditions.
<b>Materials to avoid</b>	No dangerous reaction known under conditions of normal use.
<b>Hazardous decomposition products</b>	None under normal use conditions.
<b>polymerisation</b>	Hazardous polymerisation does not occur.

#### SECTION 11: Toxicological information

##### Acute toxicity

Not Hazardous.

##### Principle Routes of Exposure/ Potential Health effects

<b>eyes</b>	May cause eye irritation with susceptible persons
<b>Skin</b>	May cause skin irritation in susceptible persons
<b>Inhalation</b>	May be harmful by inhalation
<b>INGESTION</b>	May be harmful if swallowed

<b>Carcinogenic effects</b>	None
<b>Mutagenic effects</b>	None
<b>Reproductive toxicity</b>	None.
<b>Sensitisation</b>	None

**Target Organ Effects** No known effects under normal use conditions

#### SECTION 12: Ecological information

<b>Ecotoxicity</b>	No information available
<b>Mobility</b>	No information available
<b>Biodegradation</b>	Inherently biodegradable.
<b>Bioaccumulation</b>	Material does not bioaccumulate

#### SECTION 13: Disposal considerations

Dispose of contents/containers in accordance with local regulations.

#### SECTION 14: Transport information

##### IATA

Revision date 23-Jan-2015  
Product code 4482777

Page 4 / 5  
Product name TAQMAN QSY PROBE 6,000 PMOLES

<b>Proper shipping name</b>	Not classified as dangerous in the meaning of transport regulations.
<b>Hazard Class</b>	None
<b>Subsidiary class</b>	None
<b>Packing group</b>	None
<b>UN-No</b>	None

SECTION 15: Regulatory information

**International Inventories**

Complies

SECTION 16: Other information

**Reason for revision**                      SDS sections updated

For research use only. Not intended for human or animal diagnostic or therapeutic uses.

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRENTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PUPOSE"

**End of Safety Data Sheet**


# TaqMan multiplex real-time PCR

## Get more data out of your sample

- A complete multiplex real-time PCR (qPCR) solution for gene expression and genotyping applications
- Applied Biosystems™ ABY™ and JUN™ dyes, QSY™ quencher, and a multiplex master mix for optimal amplification performance
- Up to 4-plex reactions—as sensitive as singleplex reactions, decreases the starting material required, and minimizes optimization processes

Obtaining the maximum amount of genetic information from an important but small amount of sample can be challenging. This is particularly true with formalin-fixed, paraffin-embedded (FFPE) samples or tumor biopsies that are used for translational research studies. Singleplex qPCR is frequently used for these clinical research samples, but this typically has a higher cost per sample than running in multiplex format. The additional time and materials required to set up multiple single-assay reactions could also significantly increase the cost of a complex project.

Multiplex qPCR, a strategy where more than one target in a sample is amplified and quantified in a single tube, can decrease the quantity of sample material and reagents required. A complete solution for multiplex qPCR is presented here,



**Figure 1. Fluorescence emission spectra of FAM, VIC, ABY, and JUN dyes used for multiplex real-time PCR.** Grey zones represent the filters available on Applied Biosystems™ real-time PCR systems: 1 through 6 for the QuantStudio™ 7 or 12K Flex Real-Time PCR Systems; 1 through 5 for the QuantStudio™ 6 Flex Real-Time PCR System, ViiA™ 7 Real-Time PCR System, and 7500 or 7500 Fast Real-Time PCR System. MP = Mustang Purple™ dye.

with components designed to work together for better data quality and less time for optimization. The solution consists of the following:

- Applied Biosystems™ TaqMan® probes using QSY quencher, providing maximal PCR efficiency in a multiplex format. These probes can be ordered with Applied Biosystems™ FAM™ and VIC™ dyes and also with the ABY and JUN dyes, allowing amplification of up to 4 targets in a single reaction. These reporter dyes are optimized to work together with minimal spectral overlap for improved performance (Figure 1). In addition, the QSY quencher is fully compatible with probes that have minor-groove binder (MGB) quenchers.


- The Applied Biosystems™ TaqMan® Multiplex Master Mix was developed to allow amplification of 4 targets simultaneously, without competition between targets. This master mix contains the Applied Biosystems™ Mustang Purple™ dye, a passive reference used for normalization instead of the Applied Biosystems™ ROX™ dye, allowing for measurement of JUN dye in the channel previously used to measure ROX dye.

- Off-the-shelf, predesigned assays— an RNase P assay using an ABY-QSY probe and a GAPDH assay using a JUN-QSY probe. Both assays are available in limited and nonlimited primer concentrations.
- Calibration plates for ABY, JUN, and Mustang Purple dyes, available in 96-well, 96-well Fast, and 384-well formats.
- Additional services provided through our custom services program— save time and let our Applied Biosystems™ TaqMan® Assay experts design your multiplex assays.

This multiplex solution is compatible with the Applied Biosystems™ QuantStudio™ 6, 7, and 12K Flex Real-Time PCR Systems, as well as the Applied Biosystems™ ViiA™ 7 Real-Time PCR System and the Applied Biosystems™ 7500 and 7500 Fast Real-Time PCR Systems.

### Multiplexing without compromise

The multiplex format enables cost savings and preservation of limited sample, but it's important to obtain the same sensitivity as in the singleplex format. Figure 2 demonstrates comparable results between reactions performed in individual tubes or in 4-plex reactions for a gene quantification experiment.



**Figure 2. Comparison of singleplex and multiplex gene expression assays.** (A) EGFR assay, FAM dye; (B) BRCA1 assay, VIC dye; (C) ESR1 assay, JUN dye; (D) RNase P assay, ABY dye. Amplification was performed on the QuantStudio 7 Real-Time PCR System using TaqMan Multiplex Master Mix. The figure shows amplification plots (left) and linear curves (right) for 4 assays amplified in singleplex (blue) and 4-plex reactions (red) in a dilution series from 20,000 pg to 2 pg of reference colon cDNA per 10  $\mu$ L reaction. Average  $C_q$  value (lines) and average standard deviation (bars) for the dilution series are represented in their respective graphs and show the concordance between singleplex and 4-plex reactions. PCR efficiencies are: 96.09% for EGFR singleplex and 96.39% for EGFR 4-plex; 93.56% for BRCA1 singleplex and 94.93% for BRCA1 4-plex; 97.13% for ESR1 singleplex and 95.81% for ESR1 4-plex; 96.91% for RNase P singleplex and 98.1% for RNase P 4-plex.




### Improved probe performance


Introduction of ABY and JUN reporter dyes and Mustang Purple passive reference dye allows for optimal 4-color multiplex assays when used with our FAM and VIC reporter dyes. Please note that ABY and JUN reporter dyes are available only with QSY quencher, while FAM and VIC dyes are available with either MGB or QSY quencher. A comparison with a set of dyes from another supplier shows that our combination of dyes provides an earlier  $C_t$  for the majority of assays (Figure 3).

### Optimized multiplex master mix

In multiplex PCR, it's important to have a robust master mix that allows for amplification of each target in a highly competitive environment. Our new master mix composition was developed to provide optimal multiplex performance for each target in the reaction. A comparison of our master mix and a master mix from another supplier in a 4-plex reaction shows an earlier  $C_t$  for 3 of the targets amplified with our new master mix and a lower standard deviation for most of the dilution points, demonstrating the excellent performance of our solution (Figure 4).



**Figure 3. Comparison of our new dye combination with a dye combination from another supplier.** Probes for (A) B2M, (B) HPRT, (C) GAPDH, and (D) RNase P gene expression assays were synthesized with FAM, VIC, ABY, and JUN dyes with QSY quencher (green bars and diamonds) and with another commercially available dye combination (blue bars and triangles). All possible gene-dye combinations were tested. Reactions were prepared with TaqMan Multiplex Master Mix using 900 nM of primer, 250 nM of probe, and 10 ng of cDNA. Amplification was performed on the QuantStudio 7 Real-Time PCR System using TaqMan Multiplex Master Mix. Bars represent average standard deviation. Triangles and diamonds represent average  $C_t$  values.



**Figure 4. Comparison of TaqMan Multiplex Master Mix with another commercially available master mix.** (A) B2M assay, FAM dye; (B) RNase P assay, VIC dye; (C) GAPDH assay, ABY dye; (D) HPRT assay, JUN dye. All assays used QSY quencher. The graph shows average standard deviation (bars) and average  $C_t$  values (cross and triangle) for 4-plex reactions using a dilution series from 100 ng to 10 pg of cDNA per 10  $\mu$ L reaction. All amplifications were performed on the ViiA 7 Real-Time PCR System using the cycling conditions recommended for each master mix. Green represents TaqMan Multiplex Master Mix, and blue represents 4-plex reactions with another commercially available master mix.

## Optimized to minimize time-to-results

Developing a multiplex PCR assay requires time to correctly design the assay and optimize the reaction. Using our complete solution, for which all components were developed to work together, helps increase your chances of success and limits your

development time. A new multiplex PCR user guide was developed to guide you through the development and optimization process [1], and our custom services will allow you to delegate assay design to our experienced team to minimize your efforts.

## References

1. Multiplex PCR User Guide. Available at [thermofisher.com/multiplexqpcr](http://thermofisher.com/multiplexqpcr)
2. TaqMan multiplex qPCR: Accurate, sensitive, and as efficient as traditional singleplex qPCR. Application note available at [lifetechnologies.com/multiplexqpcr](http://lifetechnologies.com/multiplexqpcr)

## Ordering information

Product	Cat. No.
<b>TaqMan QSY probes</b>	
TaqMan QSY Probe, 6,000 pmol	4482777
TaqMan QSY Probe, 20,000 pmol	4482778
TaqMan QSY Probe, 50,000 pmol	4482779
<b>Control kits</b>	
TaqMan GAPDH Assay, JUN-QSY 20X	4485712
TaqMan GAPDH Assay, JUN-QSY PL 20X	4485713
TaqMan RNaseP Assay, ABY-QSY 20X	4485714
TaqMan RNaseP Assay, ABY-QSY PL 20X	4485715
<b>Multiplex master mixes</b>	
TaqMan Multiplex Master Mix, 1 mL	4461881
TaqMan Multiplex Master Mix, 5 mL	4461882
TaqMan Multiplex Master Mix, 50 mL	4486295

Other formats are available at [lifetechnologies.com/multiplexqpcr](http://lifetechnologies.com/multiplexqpcr)

## Calibration plates

96-Well Calibration Plate, Mustang Purple dye	4461599
96-Well Calibration Plate, JUN dye	A24737
96-Well Calibration Plate, ABY dye	A24738

Calibration plates are also available for 96-well Fast and 384-well plate formats.

Visit [thermofisher.com/multiplexqpcr](http://thermofisher.com/multiplexqpcr) for more information.

Find out more at [thermofisher.com/multiplexqpcr](http://thermofisher.com/multiplexqpcr)

**ThermoFisher**  
SCIENTIFIC

# Custom Primers and TaqMan® Probes shipped at ambient temperature reduce environmental impact and retain their quality and stability

## Abstract

To minimize the adverse environmental impact of packaging and shipping products on gel or dry ice, Thermo Fisher Scientific investigated the feasibility of shipping its Custom Primers and TaqMan® Probes at ambient temperatures. This report describes stability testing of dye-labeled primers and MGB, TAMRA™, and QSY® probes after subjecting them to simulated summer shipping conditions. Analytical and stability testing demonstrated that Custom Primers and TaqMan® Probes that underwent simulated summer ambient-temperature shipping conditions maintained the same integrity and functionality as primers and probes that were kept at the recommended storage condition. By shipping at ambient temperatures, the need for expanded polystyrene (EPS) coolers and added refrigerant is eliminated and the fuel consumption and greenhouse gas emissions from transporting the product are significantly reduced.

## Introduction

The adverse environmental impact of shipping refrigerated or frozen products is tremendous. The annual carbon footprint to manufacture EPS and convert it into coolers for our Custom Primers and TaqMan® Probes is approximately 6 tons CO<sub>2</sub>-equivalents (CO<sub>2</sub>e) [1]. Factoring in the number of shipments, the average distance traveled per package, and the fact that most packages are shipped via air, the annual total carbon footprint for transporting Custom Primers and TaqMan® Probes is 32 tons CO<sub>2</sub>e [2].

There are other facts to consider beyond the greenhouse gas emissions. When a cooler arrives at the laboratory, the researcher is often put in the untenable position of deciding whether to burn additional fossil fuels to transport the empty cooler across country for reuse/recycling or to dispose of the cooler in a landfill. The best way to address the total environmental impact of “cold-chain” transport is to follow the hierarchy of “reduce, reuse, recycle”: 1) Design the product for stability to ensure it can withstand the rigors of ambient shipping conditions without added refrigerant or insulation; 2) Design the packaging to be reusable, without increasing source material consumption; and 3) Recycle locally. We have opted to reduce whenever possible, reuse when it is an environmentally preferable option, and to encourage our customers to recycle locally.

Thermo Fisher Scientific has been systematically evaluating novel ways to minimize the impact of shipping Life Technologies™ products on gel or dry ice, and the CO<sub>2</sub> footprint left by these products during distribution. Here we demonstrate that selected Custom Primers and TaqMan® Probes are stable at ambient temperatures during shipping. By avoiding the cooler and refrigerant, the product can be shipped in a smaller, corrugated cardboard box, which improves the carrier's freight density (less fuel and emissions per box) and reduces the amount of packaging materials requiring disposal or recycling. By eliminating the cooler and gel or dry ice for these products, Thermo Fisher Scientific is helping to divert an annual total of nearly 1,826 kg (5,062 cubic feet)

of EPS from landfills and incinerators by replacing it with recyclable corrugated paper packaging, and to reduce the annual total carbon footprint by 38 tons CO<sub>2</sub>e [1,2].

In 2009, we investigated the stability of five TaqMan<sup>®</sup> Assays: TaqMan<sup>®</sup> Gene Expression, Custom TaqMan<sup>®</sup> Gene Expression, TaqMan<sup>®</sup> MicroRNA, TaqMan<sup>®</sup> Drug Metabolism Genotyping, and TaqMan<sup>®</sup> SNP Genotyping Assays [3]. These assays comprise a preformulated set of unlabeled locus-specific oligonucleotide primers and minor groove binder–nonfluorescent quencher (MGB–NFQ) probes labeled with a fluorescent dye (VIC<sup>®</sup> or FAM<sup>™</sup> dye), and are supplied in liquid form. A total of 42 different TaqMan<sup>®</sup> Assays were selected to represent the widest range of performance as well as chemical, sequence, and structural motifs. Assays were subjected to simulated summer ambient shipping conditions and subsequently analyzed for physical integrity and functional performance. Stressed samples were compared to controls in analytical HPLC and functional real-time PCR assays. In all cases, simulated ambient shipping of the assays had no effect on their quality, integrity, or functional performance. This study provided ample evidence for the stability of a wide range of structural motifs and oligonucleotide sequences under ambient shipping conditions and also demonstrated the stability of the VIC<sup>®</sup> and FAM<sup>™</sup> dyes and the MGB moiety at the concentrations found in the assays.

For many years, Custom Primers and TaqMan<sup>®</sup> Probes have been shipped refrigerated on gel ice (with storage after shipping at +4°C or –20°C, depending on the product). Building on our 2009 study, this paper describes results from stability testing carried out after the Custom Primers and TaqMan<sup>®</sup> Probes were exposed to established summer shipping temperature profiles. These experiments demonstrate that by shipping selected Custom Primers and TaqMan<sup>®</sup> Probes under ambient conditions, not only can we supply researchers with the same superior-quality product they are used to receiving, but we can also reduce our environmental footprint in the process. This is a win for our customers (eliminating packaging waste and extra costs associated with refrigerated shipments), a win for our planet (reducing resource consumption and total carbon footprint), and a win for our company (eliminating the need to manage cold-chain transport).

## Materials and methods

**Products tested.** Custom Primers are 5'-labeled oligos that come with a choice of six dyes: 6-FAM<sup>™</sup>, TET<sup>™</sup>, VIC<sup>®</sup>, HEX<sup>™</sup>, NED<sup>™</sup>, or PET<sup>™</sup> dye. The Custom Primer Pairs also come with an unlabeled oligo in a separate tube. Primers and Primer Pairs may be HPLC purified and can be ordered in two or three different quantities, with the largest having the highest concentration. For this study, four different labeled primers at the highest concentration were selected to represent the variety of primer types and dyes available (Table 1). The FAM<sup>™</sup> and VIC<sup>®</sup> dyes were not tested with the primers because the 2009 study demonstrated the stability of these dyes in the TaqMan<sup>®</sup> Assays under ambient shipping conditions. Additionally, the Sequence Detection Primers, which are unlabeled, were not tested because the 2009 study established that unlabeled oligos are not affected by simulated ambient shipping conditions. All primers tested were formulated in Tris-EDTA (TE) buffer and were not HPLC purified. HPLC purification has no impact on the stability of the oligo or dye. Formulations in water were not evaluated because the pH of Tris buffers is known to vary inversely with temperature [4,5], something that does not occur in water, making TE a higher-risk formulation for ambient shipping.

TaqMan<sup>®</sup> MGB Probes incorporate a 5' reporter dye (FAM<sup>™</sup>, VIC<sup>®</sup>, TET<sup>™</sup>, or NED<sup>™</sup> dye) and a 3' nonfluorescent quencher, with the MGB moiety attached to the quencher molecule. The TAMRA<sup>™</sup> probes incorporate a 5' reporter dye (FAM<sup>™</sup>, TET<sup>™</sup>, or VIC<sup>®</sup> dye) and a 3' TAMRA<sup>™</sup> quencher dye. The TaqMan<sup>®</sup> QSY<sup>®</sup> Probes can be ordered with a 5' reporter dye (FAM<sup>™</sup>, VIC<sup>®</sup>, ABY<sup>®</sup>, or JUN<sup>®</sup> dye) and the QSY<sup>®</sup> quencher. All TaqMan<sup>®</sup> Probes are HPLC purified and supplied at a single concentration in TE buffer. The MGB, TAMRA<sup>™</sup>, and QSY<sup>®</sup> probes were each tested with their respective dyes, with the exception of the MGB probe. This probe was not tested with VIC<sup>®</sup> dye because our own unpublished studies have shown that FAM<sup>™</sup> is more labile than VIC<sup>®</sup> at elevated temperatures; therefore, FAM<sup>™</sup> was used to represent a “worst-case” scenario. Because the 2009 study showed that variation in sequence and length did not affect oligo stability, a single sequence was chosen for all primers and probes:

5' - TGGACAGCCACCGACGAGAGCCTGG - 3'

**Table 1. Custom Primers and TaqMan® Probes represented in this study.**

Product Description	Reporter dye	Cat. No.
<b>Custom Primers</b>		
Sequence Detection Primers, 10,000 pmol, 80,000 pmol, 130,000 pmol	None	4304970, 4304971, 4304972
<b>Custom 5'-Labeled Primer Pair Di-Repeats</b> , 10,000 pmol, 80,000 pmol, <b>300,000 pmol</b>	HEX™, NED™, PET®, 6-FAM™, VIC®, TET™	4304976, 4304977, <b>4304978</b>
<b>Custom 5'-Labeled Primer</b> , 10,000 pmol, 80,000 pmol, <b>300,000 pmol</b>	HEX™, NED™, PET®, 6-FAM™, VIC®, TET™	450007, 450006, <b>450017</b>
<b>Custom 5'-Labeled Primer Pair</b> , 10,000 pmol, 80,000 pmol, <b>300,000 pmol</b>	HEX™, NED™, PET®, 6-FAM™, VIC®, TET™	450056, 450059, <b>450062</b>
<b>Custom 5'-Labeled Primer Pair Di-Repeat + Tail</b> , 10,000 pmol, 80,000 pmol, <b>300,000 pmol</b>	HEX™, NED™, PET®, 6-FAM™, VIC®, TET™	4304979, 4304981, <b>4304982</b>
<b>TaqMan® Custom Probes</b>		
<b>TaqMan® MGB Probe</b> , 6,000 pmol, 20,000 pmol, <b>50,000 pmol</b>	6-FAM™, VIC®, NED™, TET™	4316034, 4316033, <b>4316032</b>
<b>TaqMan® TAMRA™ Probe</b> , 6,000 pmol, 20,000 pmol, <b>50,000 pmol</b>	VIC®, 6-FAM™, TET™	450025, 450024, <b>450003</b>
<b>TaqMan® QSY® Probe</b> , 6,000 pmol, 20,000 pmol, <b>50,000 pmol</b>	6-FAM™, VIC®, ABY®, JUN®	4482777, 4482778, <b>4482779</b>


Products tested are in **bold**

**Creating replicates.** To help eliminate manufacturing lot variability when creating the replicates, individual tubes of the primers and probes were manufactured, pooled, and aliquoted into the same packaging tube at the same fill volume as specified for the manufactured product. A total of 10 lots for each primer and probe were used to create five replicate stress tubes and five replicate control tubes. The control tubes were kept at -20°C for the duration of the study.

**Simulated shipping conditions.** To simulate temperatures experienced during shipping, samples were placed in a cycling environmental chamber (Thermotron® S-16) programmed to reproduce a “worst-case” 288-hour (12-day) summer temperature profile (Figure 1). This profile is adopted from one developed and

validated by Amgen to simulate global ambient shipping conditions and mimics product temperature extremes encountered during transit of over 2,500 shipments during summer months between the latitudes of 59.9° N and 37.8° S [6]. Testing of winter ambient conditions was not considered, due to the low risk of exposing the Custom Primers and TaqMan® Probes to cold conditions.

**Stability and integrity testing.** Structural integrity changes in stressed samples compared to controls were measured by reverse-phase HPLC (RP-HPLC) and MALDI mass spectrometry. RP-HPLC samples were analyzed using an Agilent® 1200 HPLC. The HPLC column used was a Phenomenex® Clarity® 3 µm Oligo-RP, 2.0 mm ID x 50 mm. Mobile phases used were 0.1 M TEAA (triethylamine acetate) in water and 0.1 M TEAA in



**Figure 1. 288 hr summer temperature profile used to simulate shipping conditions.** The summer temperature profile was used to mimic average high temperature extremes between the latitudes of 59.9° N and 37.8° S.

50% water/40% acetonitrile/10% methanol for the primers, MGB, and TAMRA™ probes, and 0.1 M TEAA in water and 2.0 M TEAA in 5% water/95% methanol for the QSY® probes. Absorbance was monitored at 260 nm for the oligonucleotide and at the maximum absorbance wavelength of the dye. Samples for MALDI mass spectrometry were analyzed on an AB Sciex® 4800 Plus MALDI TOF/TOF™ Analyzer.

## Results

**RP-HPLC.** RP-HPLC was used to create peak profiles of the dye-labeled primers using UV/Vis absorbance detection. Matched test and control tubes from each assay were analyzed. An example of the data is shown in Figure 2. Test and control peak profiles were compared, and the purity (peak areas) were calculated (data not shown). For all samples analyzed, test samples were judged as identical to matched controls (no degradation), confirming that the simulated shipping stress did not affect product integrity.




Figure 2. Simulated summer ambient shipping does not affect oligonucleotide stability—representative HPLC data. The effect of simulated summer ambient shipping on oligonucleotide integrity was measured by comparing RP-HPLC profiles of matched test and control samples. The HPLC chromatogram profiles of the test samples are comparable to the profiles of the control samples. There was no indication of probe or primer degradation in the simulated ambient-shipped 5'-Labeled Primer Pair Di-Repeats with the NED™ dye (red) compared to the matched control (blue).




Figure 3. Simulated summer ambient shipping does not affect oligonucleotide stability—representative MALDI mass spectrometry data. The effect of simulated summer ambient shipping on oligonucleotide integrity was measured by comparing mass spectrum profiles of matched test and control samples. The profiles of the test samples are comparable to the profiles of the control samples. There was no indication of probe or primer degradation in the simulated ambient-shipped TaqMan® TAMRA™ Probe with the VIC™ dye (bottom) compared to the matched control (top).

**MALDI mass spectrometry.** MALDI mass spectrometry was used to generate mass profiles of the dye-labeled primers and probes. Again, matched test and control assays were analyzed and compared to each other. An example mass spectrum is shown in Figure 3. Test and control samples showed the same mass profiles, indicating that no degradation of the oligo, dye, or quencher occurred during the shipping simulation, further confirming that the simulated shipping stress did not affect product integrity.

## Conclusions

The data described in this paper demonstrate that ambient shipping conditions have no effect on the quality and stability of Custom Primers and TaqMan® Probes. For each dye-labeled primer and probe tested, we were able to clearly demonstrate that ambient-temperature shipping conditions do not affect the product quality or integrity.

These results substantiate the change to ambient shipping conditions, and provide the researcher with confidence that when shipped under ambient conditions, their Custom Primers and TaqMan® Probes will exhibit no difference in function or stability compared to dry or gel ice-shipped products. In addition to ensuring our customers will continue to receive the highest quality possible, this study enables us to reduce the impact of transport of these products by 32 tons CO<sub>2</sub>e. Our customers will see a reduction of 1,826 kg of EPS waste. Our planet will collectively see CO<sub>2</sub> emissions reduced by 38 tons every year.

## References

1. Data derived from Bousted, I. 2006. Eco-profiles of the European Plastics Industry POLYSTYRENE (Expandable) (EPS). *Plastics Europe*, June.
2. Data derived from US Environmental Protection Agency. 2008. Climate leaders, greenhouse gas inventory protocol core module guidance: Optional emissions from commuting, business travel and product transport. US EPA document EPA-430-R-08-006. Available at [http://www.epa.gov/climateleadership/documents/resources/commute\\_travel\\_product.pdf](http://www.epa.gov/climateleadership/documents/resources/commute_travel_product.pdf).
3. Chang H-R, Crespo J, Lee K, McGall K, MacWhorter S, Russell P, Schatz J, Seid D, Walworth C, Yu J. 2009. TaqMan® Assays shipped at ambient temperature reduce environmental impact and retain their quality and stability. Life Technologies publication O-090071 0410. Available at [http://tools.invitrogen.com/content/sfs/brochures/cms\\_081489.pdf](http://tools.invitrogen.com/content/sfs/brochures/cms_081489.pdf).
4. Sambrook J, Russell DW. 2001. *Molecular Cloning: A Laboratory Manual*. 3rd Ed.; CSHL Press, New York, pp A1.3.
5. AppliChem GmbH. Biological Buffers. 2011. AppliChem GmbH, Darmstadt, Germany. Available at <https://www.applichem.com/fileadmin/Broschueren/BioBuffer.pdf>.
6. Cowland R. 2007. Developing ISTA Cold Chain Environmental Standards. Paper presented at the Dimensions '07 Conference, Orlando, Florida. Available at [www.ista.org/forms/COWLAND\\_RAY\\_Dimensions07.pdf](http://www.ista.org/forms/COWLAND_RAY_Dimensions07.pdf)

Find out more at [lifetechnologies.com/ecotaqman](http://lifetechnologies.com/ecotaqman)

 life  
technologies

A Thermo Fisher Scientific Brand



## Safety Data Sheet

Safety Data Sheet according to Regulation (EC) No. 1907/2006 (REACH)  
Classification according to Regulation (EC) No. 1272/2008 [CLP]

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### Product identifier

**Product code** 4351379  
**Product name** TAQMAN SNP ASSAY MTO HUMAN SM  
**Chemical Name** Not Applicable  
**REACH registration number** Formamide : 01-2119496064-35-XXXX

#### Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** For research use only  
**Use Description Code** SU22 - Professional uses: Public domain (administration, education, entertainment, services, craftsmen), PROC15 - Use as laboratory reagent, PC21 - Laboratory chemicals, SU24 - Scientific research and development  
**Uses advised against** Not for consumer use.

#### Details of the supplier of the safety data sheet

##### Manufacturer / Supplier

LIFE TECHNOLOGIES EUROPE BV  
KWARTSWEG 2  
2665 NN BLEISWIJK  
NETHERLANDS  
31-(0)180 392 400  
Email: MSDS@lifetech.com

Life Technologies Limited  
3 Fountain Drive  
Inchinnan Business Park  
Paisley  
PA4 9RF, UK  
+44 (0)141 814 6100

**24 hour Emergency Response for Hazardous Materials** Within the USA + Canada: 1-800-424-9300 and  
**[or Dangerous Goods] Incident. Spill, Leak, Fire,** 1-703-527-3887  
**Exposure, or Accident. Call CHEMTREC** Outside the USA + Canada: 1-703-741-5970

#### Country Specific Emergency Number (if available):

CHEMTREC Ireland (Dublin) +(353)-19014670 (Greeting Language: English and Irish)  
CHEMTREC UK (London) +(44)-870-8200418 (Greeting Language: English)

### SECTION 2: Hazards identification

**Revision date** 17-Apr-2020  
**Product code** 4351379

**Page** 1 / 10  
**Product name** TAQMAN SNP ASSAY MTO HUMAN SM



## Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

### Physical hazards

Not Hazardous

### Health hazards

Carcinogenicity	Category 2
Reproductive Toxicity	Category 1B
Specific target organ toxicity - Repeated exposure	Category 2

### Environmental hazards

Not Hazardous

### Additional information

No information available

## Label elements

Labelling according to Regulation (EC) No 1272/2008 [CLP]

### Hazard pictograms



### Signal Word

Danger

### Hazard Statements

H360 - May damage fertility or the unborn child if swallowed

H351 - Suspected of causing cancer if swallowed

H373 - May cause damage to organs through prolonged or repeated exposure

### Precautionary Statements

#### Prevention

P201 - Obtain special instructions before use

P202 - Do not handle until all safety precautions have been read and understood

P260 - Do not breathe dust/fume/gas/mist/vapours/spray

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P281 - Use personal protective equipment as required

#### Response

P308 + P313 - IF exposed or concerned: Get medical advice/attention

P314 - Get medical advice/attention if you feel unwell

#### Storage

Not Applicable

#### Disposal

Revision date 17-Apr-2020  
Product code 4351379

Page 2 / 10  
Product name TAQMAN SNP ASSAY MTO HUMAN SM

P501 - Dispose of contents/ container to an approved waste disposal plant

**Other hazards**

Not Applicable

SECTION 3: Composition/information on ingredients

Component	CAS No	EINECS-No.	Weight-%	REACH registration number	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Formamide 75-12-7 ( 1-5% )	75-12-7	200-842-0	1-5%	01-2119496064-35-X XXX	Repr. 1B - H360 Carc. 2 - H351 STOT RE 2 - H373

SECTION 4: First aid measures

**Description of first aid measures**

<b>Skin contact</b>	Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Immediate medical attention is required.
<b>Eye contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.
<b>Ingestion</b>	Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice. If swallowed, rinse mouth with water (only if the person is conscious). Risk of serious damage to the lungs (by aspiration). Get medical attention if symptoms occur.
<b>Inhalation</b>	Remove to fresh air. If not breathing, give artificial respiration. If symptoms persist, call a doctor.
<b>Notes to Physician</b>	Treat symptomatically.

**Most important symptoms and effects, both acute and delayed**

H360 - May damage fertility or the unborn child if swallowed H351 - Suspected of causing cancer if swallowed H373 - May cause damage to organs through prolonged or repeated exposure

**Indication of any immediate medical attention and special treatment needed**

IF exposed or concerned: Get medical advice/attention. Get medical advice/attention if you feel unwell.

SECTION 5: Firefighting measures

**Extinguishing media**

<b>Suitable extinguishing media</b>	Foam. Dry powder. Dry chemical. Carbon dioxide (CO <sub>2</sub> ). Water spray.
<b>Unsuitable extinguishing media</b>	No information available.

**Special hazards arising from the substance or mixture**

None known

**Protective equipment and precautions for firefighters**

Wear self-contained breathing apparatus and protective suit.

## SECTION 6: Accidental release measures

### **Personal precautions, protective equipment and emergency procedures**

Avoid exposure to vapour  
Avoid breathing vapours or mists  
Ensure adequate ventilation  
Avoid contact with skin, eyes or clothing  
Use personal protection equipment  
See section 8 for more information

### **Environmental precautions**

Should not be released into the environment. Prevent product from entering drains.

### **Methods and material for containment and cleaning up**

Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder, sawdust). Sweep up and shovel into suitable containers for disposal. Clean contaminated surface thoroughly.

### **Reference to other sections**

See section 8 for more information.

## SECTION 7: Handling and storage

### **Precautions for safe handling**

Always wear recommended Personal Protective Equipment. Wash hands before breaks and immediately after handling the product. Do not get in eyes, on skin, or on clothing. Avoid breathing vapours or mists. If during normal use the material presents a respiratory hazard, use adequate ventilation and/or wear appropriate respirator. See section 8 for more information.

### **Conditions for safe storage, including any incompatibilities**

Keep in properly labelled containers. Keep in a dry, cool and well-ventilated place. Store in accordance with local regulations. Keep away from combustible material.

### **Specific end use(s)**

For research use only.

SECTION 8: Exposure controls/personal protection

**Control parameters**

Chemical Name	EU OEL (TWA)	EU OEL (STEL)	EU Skin Notation
Formamide 75-12-7	None	None	None

Chemical Name	Austria	Belgium (TWA)	Czech Republic
Formamide 75-12-7	9 ppm 16 mg/m <sup>3</sup>	10 ppm 18 mg/m <sup>3</sup>	None

Chemical Name	Denmark (TWA)	Finland OEL (TWA)	France OEL (VME)
Formamide 75-12-7	10 ppm 18 mg/m <sup>3</sup>	10 ppm 19 mg/m <sup>3</sup>	20 ppm 30 mg/m <sup>3</sup>

Chemical Name	Germany OEL (TWA)	Ireland (TWA)	Italy OEL (TWA)
Formamide 75-12-7	None	10 ppm 18 mg/m <sup>3</sup>	None

Chemical Name	Lithuania OEL (TWA)	Netherlands OEL (MAC)	Norway
Formamide 75-12-7	10 ppm 20 mg/m <sup>3</sup>	None	10 ppm TWA 18 mg/m <sup>3</sup> TWA 15 ppm STEL 27 mg/m <sup>3</sup> STEL

Chemical Name	Poland	Portugal	Spain OEL (TWA)
Formamide 75-12-7	23 mg/m <sup>3</sup> TWA Skin Notation	10 ppm TWA skin - potential for cutaneous exposure	10 ppm 19 mg/m <sup>3</sup>

Chemical Name	Sweden - Occupational Exposure Limits - TLVs (LLVs)	Switzerland	United Kingdom
Formamide 75-12-7	10 ppm TLV NGV; 20 mg/m <sup>3</sup> TLV NGV	10 ppm TWA 18 mg/m <sup>3</sup> TWA	20 ppm TWA; 37 mg/m <sup>3</sup> TWA

**Engineering Measures**      Ensure adequate ventilation, especially in confined areas.

**Exposure controls**

**Personal protection equipment**

**Respiratory protection**      In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

**Hand protection**              Glove material: Nitrile rubber with thickness (mm) :5 Break through time (hours) :>1  
Recommended glove type has not been tested for use with this product.  
Information is based on professional knowledge

**Eye protection**                 Tight sealing safety goggles.

**Skin and Body Protection**      Wear laboratory coat for body protection.

**Hygiene Measures**              Handle in accordance with good industrial hygiene and safety practice.

**Environmental exposure controls**  
No special environmental precautions required.

## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

<b>Appearance</b>	liquid	
<b>Odour</b>	No data	
<b>Odour Threshold</b>	No data	
<b>Molecular Weight</b>	No data	
<b>pH</b>	8	
<b>Melting point / melting range</b>	°C 0-2.5	°F 32-36.5
<b>Boiling point / boiling range</b>	°C 100-210	°F 212-410
<b>Flash point</b>	°C >120	°F >248
<b>Autoignition Temperature</b>	°C >500	°F >932
<b>Decomposition temperature</b>	°C No data	°F No data
<b>Evaporation rate</b>	No data	
<b>Flammability (solid, gas)</b>	No data	
<b>Upper explosion limit</b>	No data	
<b>Lower explosion limit</b>	No data	
<b>Vapour Pressure</b>	No data	
<b>Vapour density</b>	No data	
<b>Relative density</b>	No data	
<b>Specific gravity</b>	No data	
<b>Solubility</b>	Soluble in water	
<b>Partition coefficient: n-octanol/water</b>	No data	
<b>Viscosity</b>	No data	
<b>Explosive properties</b>	No data	
<b>Oxidising properties</b>	No data	

### Other information

No data.

## SECTION 10: Stability and reactivity

<b>Reactivity</b>	None known.
<b>Chemical stability</b>	Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	Hazardous reaction has not been reported.
<b>Conditions to avoid</b>	High temperature. Thermal decomposition of the finish can take place above (>140 °C) >284 °C.
<b>Incompatible materials</b>	Oxidising agent. Acids. Bases. Sulphur trioxide. Iodine.
<b>Hazardous decomposition products</b>	Carbon monoxide. Hydrogen cyanide (hydrocyanic acid). Nitrogen oxides (NOx).

## SECTION 11: Toxicological information

### Information on toxicological effects

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Formamide	3200 mg/kg	13500 mg/kg	3900 ppm/6H

### Principal Routes of Exposure

<b>Skin corrosion/irritation</b>	Data are conclusive but insufficient for classification
<b>Serious eye damage/irritation</b>	Data are conclusive but insufficient for classification
<b>Respiratory or skin sensitisation</b>	Data are conclusive but insufficient for classification
<b>Specific target organ toxicity (STOT) – single exposure</b>	Data are conclusive but insufficient for classification
<b>Specific target organ toxicity (STOT) – repeated exposure</b>	Target organ(s) : Cardiovascular System Hematopoietic System
<b>Carcinogenicity</b>	Contains a known or suspected carcinogen
<b>Germ cell mutagenicity</b>	Data are conclusive but insufficient for classification
<b>Reproductive Toxicity</b>	May cause adverse reproductive effects - such as birth defect, miscarriages, or infertility
<b>Aspiration Hazard</b>	Data are conclusive but insufficient for classification

## SECTION 12: Ecological information

### Ecotoxicity

The environmental impact of this product has not been fully investigated.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Formamide	Desmodesmus subspicatus EC50>500 mg/L (72 h) Desmodesmus subspicatus EC50>500 mg/L (96 h)	Daphnia magna EC50>500 mg/L (48 h)	No data available	No data available	logPow-0.82

**Persistence and degradability** Readily biodegradable.

**Bioaccumulative potential** Material does not bioaccumulate.

### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

### Other adverse effects

No information available.

## SECTION 13: Disposal considerations

### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

## SECTION 14: Transport information

### IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

<b>UN number</b>	Not Applicable
<b>UN proper shipping name</b>	Not Applicable
<b>Transport hazard class(es)</b>	Not Applicable
<b>Packing group</b>	Not Applicable

### Environmental hazards

Not Applicable

### Special precautions for user

Not Applicable

### Transport in bulk according to Annex II of MARPOL and the IBC Code

Not Applicable.

## SECTION 15: Regulatory information

### Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Substances of Very High Concern

Chemical Name	Weight-%	EU - REACH (1907/2006) - Article 59(1) - Candidate List of Substances for Eventual Inclusion in Annex XIV
Formamide	1-5%	Reason for inclusion Toxic for reproduction, Article 57c

### Substance subject to authorisation per REACH Annex XIV

None

### Restricted substances under EC 1907/2006, Annex XVII

Chemical Name	Weight-%	EU - REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances
Formamide	1-5%	Use restricted. See item 30.

### Substances listed under Annex I of Regulation (EC) No 689/2008

None.

Revision date 17-Apr-2020  
Product code 4351379

Page 8 / 10  
Product name TAQMAN SNP ASSAY MTO HUMAN SM

**Restricted substances under Annex V of Regulation (EC) No 689/2008**

None.

**Substances under Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC**

None.

**German Water hazard classes (Wassergefährdungsklassen)**

Chemical Name	Weight-%	Water hazard class (WGK)
Formamide	1-5%	hazard class 1 - slightly hazardous to water

**Other International Inventories**

Chemical Name	EINECS (European Union)	ELINCS (European List of Notified Chemical Substances)	ENCS (Japan)	PICCS (Philippines)
Formamide	Listed	-	Listed	Listed

Chemical Name	AICS (Australia)	South Korea (KECL)	Canada (DSL)	NDSL
Formamide	Listed	Listed	Listed	-

**Chemical safety assessment**

No Chemical safety assessment has been carried out.



## SECTION 16: Other information

**Reason for revision** Update according to Commission Regulation (EU) No 830/2015  
**Revision number** 8  
**Revision date** 17-Apr-2020

### References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

Carcinogenicity	Category 2	Calculation method
Reproductive Toxicity	Category 1B	Calculation method
Specific target organ toxicity - Repeated exposure	Category 2	Calculation method

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"


# TaqMan SNP Genotyping Assays

Applied Biosystems™ TaqMan™ SNP Genotyping Assays from Thermo Fisher Scientific provide a highly flexible technology for detection of polymorphisms within any genome. Applied Biosystems™ TaqMan™ Assays have the simplest workflow available and are the quickest way to generate genotyping data. Based on powerful Applied Biosystems™ TaqMan™ probe and primer chemistry and designs, and coupled to dependable Applied Biosystems™ instruments and software, these made-to-order assays produce high-confidence results. TaqMan Assays are ideal for genotyping applications, including screening, association, candidate region, candidate gene, and fine-mapping studies.

Content-rich marker selection tools simplify study design and help you select from a library of human and mouse assays. This library includes over 7 million genome-wide human assays (of which 3.5 million are HapMap SNP-based assays, 160,000 are validated assays, and over 950,000 are coding region assays) and 10,000 mouse assays. We also offer 2,700 inventoried drug metabolism genotyping assays. Additionally, with Applied Biosystems™ Custom TaqMan™ SNP Genotyping Assays you can confidentially submit target SNP sequences for any genome to create your own assays. Let TaqMan SNP Genotyping Assays accelerate the pace of your discovery by eliminating time-consuming experimental design and optimization.

## Powerful, proven chemistry

Whether your genotyping studies require targeted detection of essential SNPs, or the flexibility to choose SNPs for mapping, TaqMan SNP Genotyping Assays are the technology of choice. Proven TaqMan



**Figure 1. Allelic discrimination is achieved by the selective annealing of TaqMan MGB probes.**

probes, which incorporate minor groove binder (MGB) technology at the 3' end, deliver superior allelic discrimination. The MGB molecule binds to the minor groove of the DNA helix, improving hybridization-based assays by stabilizing the MGB probe–template complex. This increased binding stability permits the use of probes as short as 13 bases for improved mismatch discrimination and greater flexibility when designing assays for difficult or variable sequences. In addition to SNP detection, TaqMan probes can be designed to detect multiple nucleotide polymorphisms (MNPs) and insertion/deletions (indels).

Detection is achieved with proven 5' nuclease chemistry by means of exonuclease cleavage of an allele-specific 5' dye label, which generates the permanent assay signal (Figure 1). All MGB probes include a nonfluorescent

quencher (NFQ) that virtually eliminates the background fluorescence associated with traditional quenchers, and provides a greater signal-to-noise ratio for superior assay sensitivity.

## TaqMan SNP Genotyping Assays collection

TaqMan SNP Genotyping Assays are the world's largest collection of single-tube, ready-to-use SNP assays available. The TaqMan SNP Genotyping Assays library consists of two collections of human assays and one of mouse assays, and can be supplemented with assays designed using our Custom TaqMan SNP Genotyping Assays Service.

### Over 7 million human SNP genotyping assays

This assay group contains over 7 million genome-wide SNPs, providing unprecedented marker coverage. Included in this collection are 160,000 validated assays that have approximately 10 kb spacing across gene regions. These assays were subjected to an extensive minor allele frequency test in 2–4 ethnic populations (45 individual samples per ethnic group) and as a result, offer the highest success rate. Also included are over 600,000 assays for the detection of nonsynonymous SNPs in coding regions, including many putative functional SNPs. Visit [thermofisher.com/taqmansnp](http://thermofisher.com/taqmansnp) for more information.

### Over 10,000 mouse SNP genotyping assays

The Applied Biosystems™ Mouse TaqMan™ Predesigned SNP Genotyping Assays collection consists of over 10,000 assays, and can be supplemented with assays designed using our Custom TaqMan SNP Genotyping Assays Service.

### TaqMan Drug Metabolism Genotyping Assays

The collection of Applied Biosystems™ TaqMan™ Drug Metabolism Genotyping Assays includes 2,700 assays that target high-value polymorphisms in 221 drug metabolism genes. These assays have proven performance in four different ethnic population samples, consisting of 45 individuals each. To enable easy identification, these assays have been mapped to the common public allele nomenclature websites where possible. Visit [thermofisher.com/taqmandme](http://thermofisher.com/taqmandme) for more information.

All TaqMan SNP Genotyping Assays are generated using next-generation algorithms from the Thermo Fisher Scientific bioinformatics pipeline. For all predesigned assays, bioinformatics evaluation of target SNP sequences includes the masking of adjacent SNPs and ambiguous bases so that assay design and subsequent performance is not affected by the poor quality of the underlying sequence. Lastly,

the assay designs are aligned to the human genome using BLAST to ensure that each assay binds uniquely to the intended polymorphism. As the Custom TaqMan SNP Genotyping Assay Service is confidential and secure, you simply perform your own bioinformatics analysis prior to submitting your sequence for assay design.

### Custom assay service for any possible SNP

Custom TaqMan SNP Genotyping Assays can be developed for any SNP in any organism. This service can generate assays for the detection of SNPs, MNPs, indels of up to 6 bases, or QSY™-labeled probes for multiplexing SNP assays.

Custom TaqMan SNP Genotyping Assays provide you with a complete service that includes secure and confidential ordering, assay design and manufacturing, and quality-control testing for synthesis accuracy and formulation completeness. Additionally, custom human assays are subjected to a functional test on 20 unique DNA samples.

Use the free Applied Biosystems™ Custom TaqMan™ Assay Design Tool to input and submit your sequence for assay design. This easy-to-use online resource lets you quickly submit your sequence information and start the ordering process securely and confidentially. Access the Custom TaqMan Assay Design Tool at [thermofisher.com/snpcadt](http://thermofisher.com/snpcadt)


### Quality design and manufacturing

Probes and primers used in TaqMan SNP Genotyping Assays are designed using our rigorous bioinformatics pipeline. This proprietary group of algorithms has generated millions of TaqMan Assay designs by utilizing heuristic design rules deduced from both manufacturing and assay performance data. All assays are designed to perform under universal reaction conditions, as calculated probe and primer melting temperatures are consistent and include contributions from associated probe conjugates (i.e., dyes and MGB).

After manufacturing, assay components undergo extensive laboratory testing at our state-of-the-art manufacturing facility. Quality-control testing includes mass spectrometry for sequence verification and formulation assessments of probe and primer concentrations. Additionally, all human SNP genotyping assays are functionally tested to ensure allelic discrimination.

### Simple workflow for quick results

TaqMan SNP Genotyping Assays constitute the simplest SNP genotyping technology available. We deliver your ready-to-use SNP genotyping assay at ambient temperature in a convenient, single-tube format. The rest is easy. Just combine the assay with Applied Biosystems™ TaqMan™ Genotyping Master Mix or TaqMan™ Universal PCR Master Mix and your purified DNA sample (Figure 2). There is no need to optimize probe, primer, salt concentrations, or temperature because all assays use universal reagent concentrations and thermal cycling conditions.



**Figure 2. A simple workflow and reliable instruments combine to generate fast, high-confidence results.**

After generating an endpoint read using a thermal cycler or real-time PCR instrument, no transfers, washes, or additional reagents are required, and the plate remains sealed; just read the plate and analyze the genotypes. This reduces the chance of contamination, sample mix-up, and sample loss. The simplicity of the chemistry allows you to easily automate the reaction for massively parallel genotyping studies, readily increasing the number of assays, number of samples, or both. Additionally, the analysis software allows you to auto-call genotypes, minimizing manual intervention.

### Reliable real-time PCR platforms

A suite of superior Applied Biosystems instrument platforms is available for processing and analyzing TaqMan SNP Genotyping Assays (Table 1). These instruments, which meet all throughput needs and budgets, include the 7500, 7500 Fast, 7900HT Fast, ViiA™ 7, StepOne™, and StepOnePlus™ Real-Time PCR Systems, and the QuantStudio™ 3, 5, 6, 7, and 12K Flex Real-Time PCR Systems (Figure 3). Following PCR amplification, an endpoint read can be performed on any Applied Biosystems real-time PCR system. All of these dependable instruments offer the advanced multicolor detection capabilities required for highly accurate and reproducible allelic discrimination assays.

### Data analysis software


The sophisticated SDS software package provided with all Applied Biosystems real-time PCR systems facilitates experimental setup, data collection, and assay performance analysis. The SDS software uses an advanced multicomponent algorithm to calculate the distinct signal contribution of each allele of a marker from the fluorescence measurements of each sample well during the assay plate read. The multicomponent data collected from the plate read are stored as SDS files, ready for genotype determination by the SDS software or optional Applied Biosystems™ TaqMan™ Genotyper Software (Figure 4).

**Table 1. Applied Biosystems instrument capacities.**

Instrument	Capacity
7500/7500 Fast Real-Time PCR System	96-well block (standard or Fast)
7900HT Fast Real-Time PCR System	96- and 384-well blocks (standard or Fast)
ViiA 7 Real-Time PCR System	96-well (standard or Fast), 384-well, and TaqMan Array Card blocks
StepOne Real-Time PCR System	48-well block (standard or Fast)
StepOnePlus Real-Time PCR System	96-well block (standard or Fast)
QuantStudio 3 Real-Time PCR System	96-well block (standard or Fast)
QuantStudio 5 Real-Time PCR System	96-well (standard or Fast) and 384-well blocks
QuantStudio 6 Real-Time PCR System	96-well (standard or Fast) and 384-well blocks
QuantStudio 7 Real-Time PCR System	96-well (standard or Fast), 384-well, and TaqMan Array Card blocks
QuantStudio 12K Flex Real-Time PCR System	96-well (standard or Fast), 384-well, TaqMan Array Card, and OpenArray plate blocks



**Figure 3. The QuantStudio 5 Real-Time PCR System (left) and the QuantStudio 12K Flex Real-Time PCR System (right), which offers the highest throughput of all Applied Biosystems real-time PCR instruments.**



**Figure 4. TaqMan Genotyper Software automatically determines sample genotypes and displays data.**

TaqMan Genotyper Software is a great resource for fast and accurate genotype calling. It is a free SNP genotyping data analysis tool for use with TaqMan SNP Genotyping Assays performed in 48-, 96-, or 384-well microtiter plates or on Applied Biosystems™ TaqMan™ OpenArray™ Plates. It has a state-of-the-art genotype-calling algorithm, an intuitive user interface, and enhanced study-based analysis features. The software enables multi-plate data analysis for high-throughput workflows and improved accuracy in genotype calling; versatile export features and comprehensive quality-control features facilitate streamlining of the entire workflow. TaqMan Genotyper Software can be downloaded at [thermofisher.com/taqmangenotyper](http://thermofisher.com/taqmangenotyper)

### Simple ordering

Selecting and ordering TaqMan SNP Genotyping Assays is as simple as “point and click.” Use SNPbrowser Software to select the most informative SNPs for your genotyping studies. As you identify SNPs of interest, simply upload your selected TaqMan SNP Genotyping Assays to our online ordering tool.

The TaqMan Assay online ordering tool (Figure 5) enables you to search, select, and order from our catalog of over 7 million made-to-order predesigned TaqMan SNP Genotyping Assays. You can search for SNPs using any of several criteria: National Center for Biotechnology Information (NCBI) gene ID, NCBI SNP reference ID (rs#), or gene symbol. You can further refine your search by using SNP type (i.e., intragenic, 5′ or 3′ UTR, chromosome, etc.).

Our Custom TaqMan SNP Genotyping Assays supply you with SNPs that are not available from our predesigned assay collection, including those from any nonhuman organism. This service designs assays for all possible SNP, MNP, and indel targets but without the up-front bioinformatic preparation used for the predesigned made-to-order assays. Our complementary Custom TaqMan Assay Design Tool conveniently formats your target sequence for submission to our manufacturing facilities. To order custom assays, simply prepare your target sequence according to the Design and Ordering Guide, and upload your submission file at [thermofisher.com/snpcaadt](http://thermofisher.com/snpcaadt)

The screenshot displays the TaqMan Assay search and order tool interface. On the left, there are search filters for Species (Human), Gene (BRCA1), Assay Type (Functionally Tested), and SNP Type (Intragenic, Silent Mutation, Transition Substitution). The main area shows search results for SNP ID rs34844365, located on Chromosome 17 at position 41251811. The SNP Type is Silent Mutation, and the Assay Type is Functionally Tested. The interface includes options to view details, add to cart, and view the assay on a map. Below the main result, there are sections for Product Details, Gene Details, and a table of related SNPs.

Transcript Accession	SNP Location	SNP Type	Codon Change	Amino Acid Change	Protein ID
NM_007254.3	760	Silent Mutation	ACA.ACG	T176T	NP_009225.1
NM_007257.3	760	Silent Mutation	ACA.ACG	T129T	NP_009228.2
NM_007258.3	760	Silent Mutation	ACA.ACG	T176T	NP_009229.2
NM_007299.3	760	Silent Mutation	ACA.ACG	T176T	NP_009230.3
NM_007350.3	760	Silent Mutation	ACA.ACG	T176T	NP_009231.2

**Figure 5.** Our TaqMan Assay search and order tool makes online ordering easy. For convenient online ordering and multiple search options for all our genotyping assays, including keyword, batch, and location searches, visit [thermofisher.com/taqmansnp](http://thermofisher.com/taqmansnp)



## Ordering information

Size	Human assays (Cat. No.)	Nonhuman assays (Cat. No.)	Number of SNPs	No. of 5 $\mu$ L reactions (384-well plate)	No. of 25 $\mu$ L reactions (96-well plate)	Assay mix formulation	Assay type
<b>Predesigned TaqMan SNP Genotyping Assays for Human and Mouse</b>							
Small	4351379	4351384*	>4.5 million	1,500	300	40X	made-to-order
Medium	4351376	4351382*	>4.5 million	5,000	1,000	40X	made-to-order
Large	4351374	4351380*	>4.5 million	12,000	2,400	80X	made-to-order
<b>Custom TaqMan SNP Genotyping Assays</b>							
Small	4331349	4332077	$\infty$	1,500	300	40X	made-to-order
Medium	4332072	4332075	$\infty$	5,000	1,000	40X	made-to-order
Large	4332073	4332076	$\infty$	12,000	2,400	80X	made-to-order
<b>TaqMan Drug Metabolism Genotyping Assays</b>							
Small	4362691	NA	2,700	750	150	20X	inventoried

All assays are quality-control tested using a mass spectrometer to verify sequence and yield. All assays have a VIC™ dye-labeled probe, a FAM™ dye-labeled probe, and two target-specific primers. All assays, excluding Custom TaqMan SNP Genotyping Assays, undergo bioinformatics evaluation of target SNP sequences.

Functional testing against 20 unique genomic DNA samples is performed on all custom and predesigned made-to-order human TaqMan SNP Genotyping Assays. Validation testing against four populations (45 samples/population) was performed on all 160,000 validated TaqMan SNP Genotyping Assays, and all TaqMan Drug Metabolism Genotyping Assays.

\* Over 10,000 mouse assays available.

applied  
biosystems

Find out more at [thermofisher.com/taqmansnp](https://thermofisher.com/taqmansnp)

**ThermoFisher**  
S C I E N T I F I C

**For Research Use Only. Not for use in diagnostic procedures.** © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a trademark of Roche Molecular Systems, Inc., used under permission and license. **CO020271 0116**

# SAFETY DATA SHEET

(In accordance with COMMISSION REGULATION (EU) No 830/2015)

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### Product identifier

**Product code** 4349763  
**Product name** 7500 Spectral Calibration Plate with SYBR® Dye  
**Chemical Name** Not Applicable  
**REACH registration number** No registration number is given yet for this substance / substances in this mixture since the annual import quantity is less than one tonnage per annum or the transition period for its registration according to Article 23 of REACH has not yet expired.

### Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses** For research use only  
**Use Description Code** SU22 - Public domain (administration, education, entertainment, services, craftsmen), PROC15 - Use as a laboratory reagent, PC21 - Laboratory chemicals, SU24 - Scientific research and development  
**Uses advised against** Not for consumer use.

### Details of the supplier of the safety data sheet

**Manufacturer/Supplier**  
LIFE TECHNOLOGIES EUROPE BV  
KWARTSWEG 2  
2665 NN BLEISWIJK  
NETHERLANDS  
31-(0)180 392 400  
Email: MSDS@lifetech.com

**24 hour Emergency Response for Hazardous Materials [or Dangerous Goods] Incident. Spill, Leak, Fire, Exposure, or Accident. Call CHEMTREC** Within the USA + Canada: 1-800-424-9300 and +1 703-527-3887  
Outside the USA + Canada: +1 703-741-5970

### Country Specific Emergency Number (if available):

CHEMTREC Ireland (Dublin) +(353)-19014670 (Greeting Language: English and Irish)  
CHEMTREC UK (London) +(44)-870-8200418 (Greeting Language: English)



**Classification of the substance or mixture**

**Classification according to Regulation (EC) No 1272/2008 [CLP]**

**Physical hazards**

Not Hazardous

**Health hazards**

Not Hazardous

**Environmental Hazards**

Not Hazardous

**Additional information**

Not Applicable

**Label elements**

**Labelling according to Regulation (EC) No 1272/2008 [CLP]**

**Hazard pictograms**

No Pictogram

**Signal word**

None

**Hazard Statements**

Not Applicable

**Precautionary Statements**

**Prevention**

Not Applicable

**Response**

Not Applicable

**Storage**

Not Applicable

**Disposal**

Not Applicable

**Other hazards**

Not Applicable

### SECTION 3: Composition/information on ingredients

Component	CAS-No.	EINECS-No.	Weight percent	REACH registration number	Classification according to Regulation (EC) No 1272/2008 [CLP]
Glycerin 56-81-5 ( 7-13 )	56-81-5	200-289-5	7-13	01-2119471987-18-X XXX	Not Classified

### SECTION 4: First aid measures

#### Description of first aid measures

<b>Skin contact</b>	Rinse with plenty of water. Immediate medical attention is not required.
<b>Eye contact</b>	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
<b>INGESTION</b>	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice.
<b>Inhalation</b>	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. Consult a physician if necessary.
<b>Notes to Physician</b>	Treat symptomatically.

#### Most important symptoms and effects, both acute and delayed

Not Applicable

#### Indication of any immediate medical attention and special treatment needed

None.

### SECTION 5: Firefighting measures

#### Extinguishing media

**Suitable Extinguishing Media**  
**Unsuitable Extinguishing Media**

Water spray. Carbon dioxide (CO<sub>2</sub>). Foam. Dry chemical.  
No information available.

#### Special hazards arising from the substance or mixture

Not Known.

#### Advice for fire-fighters

Standard procedure for chemical fires.

### SECTION 6: Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.

#### Environmental precautions

No special environmental precautions required. Avoid discharge into drains and waterways whenever possible.

#### Methods and material for containment and cleaning up

Soak up with inert absorbent material.

#### Reference to other sections

See section 8 for more information.

## SECTION 7: Handling and storage

### **Precautions for safe handling**

Use personal protective equipment as required. No special handling advices are necessary.

### **Conditions for safe storage, including any incompatibilities**

Keep in a dry, cool and well-ventilated place. Keep in properly labelled containers.

### **Specific end use(s)**

For research use only.

## SECTION 8: Exposure controls/personal protection

### Control parameters

Chemical Name	EU OEL (TWA)	EU OEL (STEL)	EU Skin Notation	
Glycerin 56-81-5	None	None	None	
Chemical Name	Austria	Belgium (TWA)	Denmark (TWA)	Finland OEL (TWA)
Glycerin 56-81-5	None	10 mg/m <sup>3</sup>	None	None
Chemical Name	France OEL (VME)	Germany OEL (TWA)	Ireland (TWA)	Italy OEL (TWA)
Glycerin 56-81-5	10 mg/m <sup>3</sup>	200 mg/m <sup>3</sup> exposure factor 2	10 mg/m <sup>3</sup>	None
Chemical Name	Sweden - Occupational Exposure Limits - TLVs (LLVs)	Netherlands OEL (MAC)	Spain OEL (TWA)	United Kingdom
Glycerin 56-81-5	None	None	10 mg/m <sup>3</sup>	10 mg/m <sup>3</sup> TWA (mist)
Chemical Name	European Union	France OEL (VME)	Germany OEL (TWA)	
Glycerin 56-81-5	None	10 mg/m <sup>3</sup>	200 mg/m <sup>3</sup> exposure factor 2	
Chemical Name	Italy OEL (TWA)	Portugal	Netherlands OEL (MAC)	Finland OEL (TWA)
Glycerin 56-81-5	None	None	None	None
Chemical Name	Austria	Denmark	Poland	Switzerland
Glycerin 56-81-5	None	None	None	None
Chemical Name	Ireland	Norway	Lithuania OEL (TWA)	Spain OEL (TWA)
Glycerin 56-81-5	None	None	None	10 mg/m <sup>3</sup>

**Engineering measures**      Ensure adequate ventilation, especially in confined areas.

### Exposure controls

#### Personal protection equipment

**Respiratory protection**      In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

**Hand Protection**              Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.

**Eye protection**                 Tight sealing safety goggles.

**Skin and body protection**      Wear suitable protective clothing.

**Hygiene measures**              Handle in accordance with good industrial hygiene and safety practice.

#### Environmental exposure controls

No special environmental precautions required.

## SECTION 9: Physical and chemical properties

### Information on basic physical and chemical properties

<b>Appearance</b>	Liquid	
<b>Odour</b>	no data available	
<b>pH</b>	Mixture has not been tested	
<b>Melting point / melting range</b>	°C >0	°F >32
<b>Boiling point / boiling range</b>	°C >100	°F >212
<b>Flash point</b>	°C >90	°F >160
<b>Autoignition temperature</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Decomposition temperature</b>	°C Mixture has not been tested	°F Mixture has not been tested
<b>Evaporation rate</b>	No data available	
<b>Flammability (solid, gas)</b>	Not applicable	
<b>Upper explosion limit</b>	Mixture has not been tested	
<b>Lower explosion limit</b>	Mixture has not been tested	
<b>Vapour Pressure</b>	Mixture has not been tested	
<b>Relative density</b>	Mixture has not been tested	
<b>Specific gravity</b>	No data available	
<b>Solubility</b>	Soluble in water	
<b>Partition coefficient: n-octanol/water</b>	No data available	
<b>Explosive properties</b>	Mixture has not been tested	

### OTHER INFORMATION

No data available.

## SECTION 10: Stability and reactivity

<b>Reactivity</b>	None known.
<b>Chemical stability</b>	Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	Hazardous reaction has not been reported.
<b>Conditions to Avoid</b>	None under normal processing.
<b>Incompatible Materials</b>	Strong acids. oxidising agents. Acetic anhydride. Isocyanates. Ammonia. Bases.
<b>Hazardous decomposition products</b>	Carbon oxides.

## SECTION 11: Toxicological information

### Information on toxicological effects

Chemical Name	LD50 (oral, rat/mouse)	LD50 (dermal, rat/rabbit)	LC50 (inhalation, rat/mouse)
Glycerin	= 12600 mg/kg Oral	No data available	>570mg/m <sup>3</sup> (Rat)

### Principal Routes of Exposure, Potential health effects

<b>Irritation</b>	Conclusive but not sufficient for classification
<b>Corrosivity</b>	Conclusive but not sufficient for classification
<b>Sensitisation</b>	Conclusive but not sufficient for classification
<b>STOT - Single Exposure</b>	Conclusive but not sufficient for classification
<b>STOT - Repeated Exposure</b>	Conclusive but not sufficient for classification
<b>Carcinogenicity</b>	Conclusive but not sufficient for classification
<b>Mutagenicity</b>	Conclusive but not sufficient for classification
<b>Reproductive Toxicity</b>	Conclusive but not sufficient for classification
<b>Aspiration Hazard</b>	Conclusive but not sufficient for classification

## SECTION 12: Ecological information

### Toxicity

Contains no substances known to be hazardous to the environment or not degradable in waste water treatment plants.

Chemical Name	Freshwater Algae Data	Water Flea Data	Freshwater Fish Species Data	Microtox Data	log Pow
Glycerin	No data available	Daphnia magna EC50>500 mg/L (24 h)	No data available	No data available	logPow-1.76

**Persistence and degradability** Inherently biodegradable.

**Bioaccumulative potential** Material does not bioaccumulate.

### Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

**Other adverse effects** No information available.

## SECTION 13: Disposal considerations

### Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

SECTION 14: Transport information

**IATA / ADR / DOT-US / IMDG**

Not classified as dangerous in the meaning of transport regulations.

<b>UN Number</b>	Not Applicable
<b>UN proper shipping name</b>	Not Applicable
<b>Transport hazard class(es)</b>	Not Applicable
<b>Packing group</b>	Not Applicable

**Environmental Hazards**

Not Applicable

**Special precautions for user**

Not Applicable

**Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

Not Applicable.

SECTION 15: Regulatory information

**Safety, health and environmental regulations/legislation specific for the substance or mixture**

**Substances of Very High Concern**

None.

**Restricted substances under EC 1907/2006, Annex XVII**

None.

**Substances listed under Annex I of Regulation (EC) No 689/2008**

None.

**Restricted substances under Annex V of Regulation (EC) No 689/2008**

None.

**Substances under Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC**

None.

**German Water hazard classes (Wassergefährdungsklassen)**

Chemical Name	Weight percent	Germany - Water Classification (VwVwS) - Annex 1	Germany - Water Classification (VwVwS) - Annex 2 - Water Hazard Classes	Germany - Water Classification (VwVwS) - Annex 3
Glycerin	7-13		hazard class 1 - low hazard to waters	

**Other International Inventories**

Chemical Name	EINECS (European Union)	ELINCS (European List of Notified Chemical Substances)	ENCS (Japan)	PICCS (Philippines)
Glycerin	Listed	-	Listed	Listed

Chemical Name	AICS (Australia)	South Korea (KECL)	Canada (DSL)	NDSL
Glycerin	Listed	Listed	Listed	-

**Chemical Safety Assessment**

No Chemical safety assessment has been carried out.



## SECTION 16: Other information

**Reason for revision** Update according to Commission Regulation (EU) No 830/2015  
**Revision number** 5  
**Revision date** 11-Dec-2017

### References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC)

#### 1272/2008 [CLP]:


Not classified

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

# 7500 Real-Time PCR Systems Spectral Calibration Kit I

Catalog Number 4349180

Pub. No. 4350071 Rev. D

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](http://thermofisher.com/support).

## Contents and storage

Contents	Amount	Storage
Background Plate sealed with an optical cover	1	-25°C to -15°C
Spectral Calibration Plates sealed with optical covers	7	
Region of Interest (ROI) Calibration Plate sealed with an optical cover	1	

## Related Documentation

For detailed information on instrument setup and the calibration process, refer to the *Applied Biosystems™ 7300/7500/7500 Fast Real-Time PCR System Installation and Maintenance Guide* (Pub. no. 4347828).

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

For support visit [thermofisher.com/support](http://thermofisher.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

The information in this guide is subject to change without notice.

### DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

**For Research Use Only. Not for use in diagnostic procedures.**

# Certificate of Analysis

**ThermoFisher**  
SCIENTIFIC

## BigDye® XTerminator™ Kit 100 Reactions

Product No. **4376486**  
Lot No. **2010180**  
Date of Manufacture **24NOV2020**  
Expiration Date **20MAY2021**

TEST	SPECIFICATION	RESULT
<b>4376495, SAM™ Solution (9mL)</b>		
<b>PRODUCT APPEARANCE</b>		
*Note: Solution may precipitate upon storage. See product protocol for guidance if precipitation is observed.		
Appearance	Colorless, Particulate Free Solution*	Pass
<b>pH TEST</b>		
pH Range	5.5 – 7.5	7.087
<b>CONDUCTIVITY TEST</b>		
Conductivity	<50µS/cm	8.507 µS/cm
<b>4376492, XTerminator™ Solution (2mL)</b>		
<b>PRODUCT APPEARANCE</b>		
Appearance	Light Brown, Opaque Suspension	PASS
<b>USE TEST</b>		
A DNA sequencing reaction prepared with BigDye® Terminator v.3.1 chemistry and 100ng of pGEM template was purified with XTerminator™ Solution/SAM™ Solution. The purified reaction was analyzed using an AB 3130xl Genetic Analyzer with a 50 cm capillary and POP-7™ polymer.		
% of capillaries exhibiting zero blobs between 60 and 100 bp	>=85%	PASS


SAM™ Solution Lot No.: 2010091  
XTerminator™ Solution Lot No.: 2010112

**For Research Use Only. Not for use in diagnostic procedures.**

Manufactured in compliance with our ISO 13485 certified quality management system.

**Bedford Site:** 2 Preston Court Bedford, MA, USA

Life Technologies  
5781 Van Allen Way  
Carlsbad, CA, USA 92008  
[www.thermofisher.com](http://www.thermofisher.com)  
For inquiries, contact us at  
[cofarequests@thermofisher.com](mailto:cofarequests@thermofisher.com)

  
Quality Assurance  
Issued 30NOV2020

Doc: 4376735 Rev F

# BigDye® XTerminator™ Purification Kit

- Single-well plate—No liquid transfer
- Complete dye blob removal
- Stabilized samples
- Improves any sequencing workflow
- Manual or automated protocols
- Rapid, reliable, and reproducible

The BigDye® XTerminator™ Purification Kit is a fast, simple purification method for DNA sequencing reactions that improves the sequencing workflow (Figure 1) and removes unincorporated BigDye terminators. No more dye blobs! Cleanup is complete in under 40 minutes and requires less than 10 minutes of labor.

More effective than ethanol precipitation (no disappearing pellets!), and cheaper than column purification and other commercial kits, the BigDye XTerminator Purification Kit removes dye blobs by capturing unincorporated dye

terminators, salts and other charged molecules that may interfere with base calling and electrokinetic sample injection (Figures 2–4). The kit is compatible with all terminator sequencing chemistries, reaction volumes, and template types.

### Simple Purification Process


Traditional purification methods, such as ethanol precipitation, require the addition of multiple reagents along with decanting and centrifuging steps.

The BigDye XTerminator Purification Kit requires the addition of only two reagents, which can be added sequentially or premixed:


- XTerminator™ Solution—Scavenges unincorporated dye terminators and free salts from the post-sequencing reaction
- SAM™ Solution—Enhances the performance of the XTerminator™ Solution and stabilizes the post-purification reactions

### Simplified Sequencing Workflow


Perform cycle sequencing with BigDye® terminator chemistry



**Figure 1.** Simplified and faster sequencing workflow using the BigDye® XTerminator™ Purification Kit.



**Figure 2.** High Quality, Long Read Sequence Data Using the BigDye XTerminator Purification Kit. High quality dye blob-free data can be produced using BigDye XTerminator Purification Kit without sacrifice to sequence length of read. Q20 base read = 1044 bp. Data kindly provided by the DNA Technology Unit, Plant Biotechnology Institute (NRC-PBI), Saskatoon, Canada.



**Figure 3.** Purification using BigDye XTerminator Kit versus Ethanol Precipitation (A) BigDye XTerminator Purification Kit yields dye blob-free data. (B) Ethanol purification generally yields sequence data with dye blobs.

### Simplified Sequencing Workflow

**STEP 1.** Perform cycle sequencing with BigDye® terminators.

**STEP 2.** After cycle sequencing, centrifuge the reaction plate briefly, then pipette the SAM™ Solution into each well, per the table below.

Note: XTerminator™ and SAM™ solutions may be premixed and added in a single step. For details, please see Big Dye® XTerminator™ Purification Kit Protocol (P/N 4374408).

Sequencing Reaction Volume	SAM Solution (µL/well)
10 µL/well (96-well plate)	45
20 µL/well (96-well plate)	90

**STEP 3.** Pipette the correct volume of XTerminator Solution into each well.

- Vortex the XTerminator Solution bulk container briefly.
- Using a wide-bore pipette tip, aspirate the solution per the table below and add to each well.

Sequencing Reaction Volume	BigDye XTerminator Solution (µL/well)
10 µL/well (96-well plate)	10
20 µL/well (96-well plate)	20

### Product Description

BigDye XTerminator Purification Kit

### Package

- 2 mL (~100 x 20 µL purifications)
- 20 mL (~1,000 x 20 µL purifications)
- 50 mL (~2,500 x 20 µL purifications)

### Part Number


- 4376486
- 4376487
- 4376484

For Research Use Only. Not for use in diagnostic procedures.

This product is compatible with DNA sequencing or fragment analysis methods covered by patents owned or licensed by Applied Biosystems. No license under these patents to use the DNA sequencing or fragment analysis methods is conveyed to the purchaser expressly, by implication, or by estoppel.

© Copyright 2006, Applied Biosystems. All rights reserved. Applera, AB (Design), Applied Biosystems, and BigDye are registered trademarks and SAM and XTerminator are trademarks of Applera Corporation or its subsidiaries in the U.S. and/or certain other countries. All other trademarks are the property of their respective owners.

Printed in the USA, 12/2006 Publication 107PR04-01



**Figure 4.** (A) BigDye XTerminator Purification yields high quality data and superior short fragment recovery and removes dye blobs from unincorporated dye terminators. (B) Traditional methods yield sequence data with dye blobs. Data kindly provided by the DNA Technology Unit, Plant Biotechnology Institute (NRC-PBI), Saskatoon, Canada.

- Repeat steps “a” through “b” to pipette the complete volume of XTerminator Solution into each well.

**STEP 4.** Seal the plate using a heat seal or Clear Adhesive Film (P/N 4306311). Vortex for 30 minutes, then centrifuge the reaction plate briefly.

**STEP 5.** Place the reaction plate in the Applied Biosystems DNA analyzer. Select run module and run plate.

**Note:** This protocol can be automated on a Biomek® FX Laboratory Automation Workstation (Beckman Coulter). For details please see Big Dye XTerminator Purification Kit Protocol (P/N 4374408).

### Software

Applied Biosystems provides downloadable run modules for use with the Big Dye XTerminator Purification Kit and Data Collection Software. These modules are designed for the 3100, 3100–Avant, 3130/3130xI, and 3730/3730xI Analyzers.

The software modules for the BigDye XTerminator Purification Kit are compatible with the Windows® NT, Windows® 2000, and Windows® XP operating systems and are available for download at [www2.appliedbiosystems.com/support/software/](http://www2.appliedbiosystems.com/support/software/).



# ExoSAP-IT™ PCR Product Cleanup

Brief Protocol

Catalog Number 78200, 78201, 78202, 78205, and 78250

Doc. Part No. 78200b Pub. No. MAN0016836 Rev. A.0 (02/2017)

**WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](http://thermofisher.com/support).

## Product description


ExoSAP-IT™ reagent treats PCR products ranging in size from less than 100 bp to over 20 kb with absolutely no sample loss by removing unused primers and nucleotides. Add ExoSAP-IT™ reagent directly to the reaction products following PCR. ExoSAP-IT™ PCR Product Cleanup is active in commonly used PCR buffers, so no buffer exchange is required. After treatment, ExoSAP-IT™ reagent is inactivated by heating to 80°C for 15 minutes. The treated PCR products are now ready for subsequent analysis in applications that require DNA to be free of excess primers and nucleotides.

## PCR cleanup protocol

**Note:** Store ExoSAP-IT™ reagent at -20°C in a non-frost-free freezer.

1. Remove ExoSAP-IT™ reagent from -20°C freezer and keep on ice throughout this procedure.
2. Mix 5 µL of a post-PCR reaction product with 2 µL of ExoSAP-IT™ reagent for a combined 7 µL reaction volume.  
When treating PCR product volumes greater than 5 µL, simply increase the amount of ExoSAP-IT™ reagent proportionally.
3. Incubate at 37°C for 15 minutes to degrade remaining primers and nucleotides.
4. Incubate at 80°C for 15 minutes to inactivate ExoSAP-IT™ reagent.
5. The PCR product is now ready for use in DNA sequencing, SNP analyses, or other primer-extension applications. Treated PCR products may be stored at -20°C until required.

## ExoSAP-IT™ Cleanup product overview



## Customer and technical support

Visit [thermofisher.com/support](http://thermofisher.com/support) for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
  - Product FAQs
  - Software, patches, and updates
  - Training for many applications and instruments
- Order and web support
- Product documentation, including:
  - User guides, manuals, and protocols
  - Certificates of Analysis
  - Safety Data Sheets (SDSs; also known as MSDSs)

**Note:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

---

**Manufacturer's address:** Affymetrix Inc. | 3450 Central Expressway | Santa Clara, CA 95051 | USA

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

**Corporate entity:** Life Technologies Corporation | Carlsbad, CA 92008 USA | Toll Free in USA 1 800 955 6288

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

## CERTIFICATE OF ANALYSIS

**78201.1.ML ExoSAP-IT™**

**Packaging Lot: 01013900**

**Expiry Date: 30.11.2022 (DD.MM.YYYY)**

**Storage: at -20±5°C**

### Filling lots for components in package:

Lot	Quantity	Description
01010236	1 mL	ExoSAP-IT™

## QUALITY CONTROL

Parameter	Method	Requirement	Result
Functional Assay	Test and Control ExoSAP-IT™ samples are used to eliminate unincorporated dNTP and unused primers from the PCR product. Exonuclease I/SAP-treated PCR product is sequenced. Quality of the sequence is evaluated visually.	There are no extra or missing bands, absence of stops in Test sequences. Test and Control sequences have uniform band intensity.	Conforms

### ISO CERTIFICATION

Manufactured by Thermo Fisher Scientific Baltics UAB, in compliance with ISO 9001 and ISO 13485 certified quality management system.

Quality authorized by QC: **J. Žilinskienė**





## CERTIFICATE OF ANALYSIS

**K0702 GeneJET PCR Purification Kit**

**Packaging Lot: 01017618**

**Expiry Date: 31.12.2023 (DD.MM.YYYY)**

**Storage: at 5±3°C**

### Filling lots for components in package:

Lot	Quantity	Description
00996508	60 mL	Binding Buffer
00990020	45 mL	Wash Buffer (concentrated)
01015366	30 mL	Elution Buffer
01013510		PCR Purification Columns

### QUALITY CONTROL

Parameter	Method	Requirement	Result
Functional testing	The kit was tested in the purification of PCR products. The quality of the purified DNA was evaluated spectrophotometrically, by agarose gel electrophoresis, digestion with FastDigest restriction enzymes and automated fluorescent sequencing.	Conforms	Conforms

#### ISO CERTIFICATION

Manufactured by Thermo Fisher Scientific Baltics UAB, in compliance with ISO 9001 and ISO 13485 certified quality management system.

Quality authorized by QC: **J. Žilinskienė**





PRODUCT INFORMATION

**Thermo Scientific**  
**GeneJET PCR Purification Kit**  
**#K0701, #K0702**

[www.thermoscientific.com/onebio](http://www.thermoscientific.com/onebio)

**#K0701**

**Lot**

50 preps

**Expiry Date**

### **CERTIFICATE OF ANALYSIS**

The kit was tested in the purification of PCR products according to the protocol described in the manual. The quality of the purified DNA was evaluated spectrophotometrically, by agarose gel electrophoresis, digestion with Thermo Scientific FastDigest restriction enzymes and automated fluorescent sequencing.

**Quality authorized by:**



Jurgita Zilinskiene

<b>CONTENTS</b>	page
COMPONENTS OF THE KIT .....	2
STORAGE AND STABILITY .....	2
DESCRIPTION.....	2
PRINCIPLE .....	3
IMPORTANT NOTES.....	3
ADDITIONAL MATERIALS AND EQUIPMENT REQUIRED .....	3
PURIFICATION PROTOCOLS.....	4
Protocol A. DNA purification using centrifuge .....	4
Protocol B. DNA purification using vacuum manifolds.....	5
TROUBLESHOOTING .....	6
SAFETY INFORMATION .....	7

## COMPONENTS OF THE KIT

GeneJET PCR Purification Kit	50 preps #K0701	250 preps #K0702
Binding Buffer	12 mL	60 mL
Wash Buffer (concentrated)	9 mL	45 mL
Elution Buffer (10 mM Tris-HCl, pH 8.5)	15 mL	30 mL
GeneJET Purification Columns (preassembled with collection tubes)	50	250

## STORAGE AND STABILITY

Thermo Scientific GeneJET PCR Purification Kit should be stored at room temperature (15-25 °C). For columns we recommend 4 °C storage for periods greater than 1 year. Any precipitate that forms in the buffers during storage can be redissolved by incubating briefly at 37 °C, then cooling to room temperature before use.


**Note. Close the bag with GeneJET Purification Columns tightly after each use!**

## DESCRIPTION

The GeneJET™ PCR Purification Kit is designed for rapid and efficient purification of DNA from PCR and other enzymatic reaction mixtures.

The kit utilizes a proprietary silica-based membrane technology in the form of a convenient spin column, eliminating the need for tedious resin manipulations or toxic phenol-chloroform extractions.

The GeneJET PCR Purification Kit effectively removes primers, dNTPs, unincorporated labeled nucleotides, enzymes and salts from PCR and other reaction mixtures. The kit can be used for purification of DNA fragments from 25 bp to 20 kb. The recovery rates are 90-100% in a 100 bp – 10 kb DNA fragment size range (see Fig. 1). Each GeneJET purification column has a total binding capacity of up to 25 µg of DNA and the entire procedure takes just 5 min. The purified DNA can be used in common downstream applications such as sequencing, restriction digestion, labeling, ligation, cloning, *in vitro* transcription, blotting or *in situ* hybridization.



**Fig. 1.** Recovery dependence on DNA fragment size

## PRINCIPLE

A reaction mixture containing DNA is combined with the binding buffer and added to a purification column. A chaotropic agent in the binding buffer denatures proteins and promotes DNA binding to the silica membrane in the column. As an added convenience, the binding buffer contains a color indicator that allows for easy monitoring of the solution pH for optimal DNA binding. Impurities are removed with a simple wash step. Purified DNA is then eluted from the column with the elution buffer. The recovered DNA is ready for use in downstream applications.

## IMPORTANT NOTES

- Prior to the initial use of the kit, dilute the **Wash Buffer** (concentrated) with ethanol (96-100%):

	<b>50 preps #K0701</b>	<b>250 preps #K0702</b>
Wash Buffer (concentrated)	9 mL	45 mL
Ethanol	45 mL	225 mL
Total Volume	54 mL	270 mL

After the ethanol has been added, mark the check box on the bottle to indicate the completed step.

- Examine the **Binding Buffer** for precipitates before each use. Re-dissolve any precipitate by warming the solution to 37 °C and cooling to 25 °C.
- Wear gloves when handling the **Binding Buffer** as this solution contains irritants (see p.7 for SAFETY INFORMATION).

## ADDITIONAL MATERIALS AND EQUIPMENT REQUIRED

- Ethanol 96-100%.
- Isopropanol.
- 3 M sodium acetate, pH 5.2 (may be necessary).
- Microcentrifuge.
- 1.5 or 2 mL microcentrifuge tubes.
- Heating block or water bath (may be necessary).

## PURIFICATION PROTOCOLS

### Note

- Read IMPORTANT NOTES on p. 3 before starting.
- All purification steps should be carried out at **room temperature**.
- All centrifugations should be carried out in a table-top microcentrifuge at **>12000 × g** (10 000-14 000 rpm, depending on the rotor type).

### Protocol A. DNA purification using centrifuge

Step	Procedure
1	Add a <b>1:1 volume</b> of <b>Binding Buffer</b> to completed PCR mixture (e.g. for every 100 µL of reaction mixture, add 100 µL of Binding Buffer). Mix thoroughly. Check the color of the solution. A yellow color indicates an optimal pH for DNA binding. If the color of the solution is orange or violet, add 10 µL of 3 M sodium acetate, pH 5.2 solution and mix. The color of the mix will become yellow.
2 for DNA ≤500 bp	<i>Optional:</i> if the DNA fragment is ≤500 bp, add a 1:2 volume of 100% isopropanol (e.g., 100 µL of isopropanol should be added to 100 µL of PCR mixture combined with 100 µL of Binding Buffer). Mix thoroughly. <b>Note.</b> If PCR mixture contains primer-dimers, purification without isopropanol is recommended. However, the yield of the target DNA fragment will be lower.
3	Transfer up to 800 µL of the solution from step 1 (or optional step 2) to the GeneJET purification column. Centrifuge for 30-60 s. Discard the flow-through. <b>Notes.</b> If the total volume exceeds 800 µL, the solution can be added to the column in stages. After the addition of 800 µL of solution, centrifuge the column for 30-60 s and discard flow-through. Repeat until the entire solution has been added to the column membrane. <b>Close the bag with GeneJET Purification Columns tightly after each use!</b>
4	Add <b>700 µL</b> of <b>Wash Buffer</b> (diluted with the ethanol as described on p. 3) to the GeneJET purification column. Centrifuge for 30-60 s. Discard the flow-through and place the purification column back into the collection tube.
5	Centrifuge the empty GeneJET purification column for an additional 1 min to completely remove any residual wash buffer. <b>Note.</b> This step is essential as the presence of residual ethanol in the DNA sample may inhibit subsequent reactions.
6	Transfer the GeneJET purification column to a clean 1.5 mL microcentrifuge tube (not included). Add <b>50 µL</b> of <b>Elution Buffer</b> to the center of the GeneJET purification column membrane and centrifuge for 1 min. <b>Note</b> <ul style="list-style-type: none"> <li>• For low DNA amounts the elution volumes can be reduced to increase DNA concentration. An elution volume between 20-50 µL does not significantly reduce the DNA yield. However, elution volumes less than 10 µL are not recommended.</li> <li>• If DNA fragment is &gt;10 kb, prewarm Elution Buffer to 65 °C before applying to column.</li> <li>• If the elution volume is 10 µL and DNA amount is ≥5 µg, incubate column for 1 min at room temperature before centrifugation.</li> </ul>
7	Discard the GeneJET purification column and store the purified DNA at -20 °C.

## Protocol B. DNA purification using vacuum manifolds

Step	Procedure
1	Perform DNA binding stage according to steps 1 - 2 in Protocol A on page 4.
2	Prepare the vacuum manifold according to the supplier's instructions. Place the GeneJET purification column(s) onto the manifold. <b>Close the bag with GeneJET Purification Columns tightly after each use!</b>
3	Transfer up to 800 $\mu\text{L}$ of the solution (from step 1 or 2 as in protocol A) to the GeneJET purification column. <b>Note.</b> If the total volume exceeds 800 $\mu\text{L}$ , the solution can be added to the column in stages. After each application, apply the vacuum and discard the flow-through. Repeat until the entire volume has been applied to the column membrane.
4	Apply the vacuum to draw the sample through the column. Switch off the vacuum after the entire sample has passed through the column.
5	Add <b>700 <math>\mu\text{L}</math> of Wash Buffer</b> (diluted with the ethanol as described on p. 3) to the GeneJET purification column.
6	Apply the vacuum to draw the solution through the column. Switch off the vacuum after the solution has passed through the column. Place the purification column back into the collection tube.
7	Centrifuge the empty GeneJET purification column for an additional 1 min to completely remove any residual wash buffer. <b>Note.</b> This step is essential as the presence of residual ethanol in the DNA sample may inhibit subsequent reactions.
8	Transfer the GeneJET purification column to a clean 1.5 mL microcentrifuge tube (not included). Add <b>50 <math>\mu\text{L}</math> of Elution Buffer</b> to the center of the GeneJET purification column membrane and centrifuge for 1 min. <b>Note.</b> <ul style="list-style-type: none"><li>• For low DNA amounts the elution volumes can be reduced to increase DNA concentration. An elution volume between 20-50 <math>\mu\text{L}</math> does not significantly reduce the DNA yield. However, elution volumes less than 10 <math>\mu\text{L}</math> are not recommended.</li><li>• If DNA fragment is &gt;10 kb, prewarm Elution Buffer to 65 °C before applying to column.</li><li>• If the elution volume is 10 <math>\mu\text{L}</math> and DNA amount is <math>\geq 5 \mu\text{g}</math>, incubate column for 1 min at room temperature before centrifugation.</li></ul>
9	Discard the GeneJET purification column and store the purified DNA at -20 °C.



## TROUBLESHOOTING

Problem	Possible Cause and Solution
<b>Low DNA yield</b>	<p><b>Inefficient DNA binding</b> Verify that a 1:1 volume of Binding Buffer is added to the reaction mixture. Ensure the solutions are mixed well. Check the color of the solution. A yellow color indicates an optimal pH for DNA binding. If the solution color is orange or violet, add 10 <math>\mu\text{L}</math> of 3 M sodium acetate, pH 5.2 solution and mix. The color of the mix will become yellow.</p> <p><b>Inefficient membrane wash</b> Ensure that the recommended volume of ethanol has been added to the Wash Buffer (concentrated) prior first use (see p. 3).</p> <p><b>Inefficient DNA elution</b> Add Elution Buffer directly to the center of the membrane and not to the side of the GeneJET purification column. Use 20-50 <math>\mu\text{L}</math> of Elution Buffer and ensure that the volume completely covers the surface of the membrane. Increase the Elution Buffer volume twice or perform two elution cycles when purifying larger amounts of DNA. (e.g., &gt;15 <math>\mu\text{g}</math>). In step 5 of Protocol A (step 7 of Protocol B), ensure all residual wash buffer is removed from the membrane. Longer centrifugation time (extra minute) can aid in removal of wash buffer.</p> <p><b>PCR reaction mixture does not contain DNA</b> Check for the presence and yield of the PCR product by running an aliquot of the reaction on an agarose gel.</p>
<b>Downstream reactions are unsuccessful</b>	<p><b>Presence of residual ethanol</b> In step 5 of Protocol A (step 7 of Protocol B), ensure all residual wash buffer is removed from the membrane. Longer centrifugation time can aid in removal of wash buffer.</p> <p><b>Inefficient membrane wash</b> Ensure that the collection tube is not overfilled during the wash step and that any of the wash buffer has remained in the bottom of the GeneJET purification column. Always discard the flow-through after centrifugation.</p> <p><b>Eluate contains excess salt</b> Ensure that the wash step 4 of Protocol A is effective. Incubate the GeneJET purification column with the Wash Buffer for several minutes before proceeding to centrifugation.</p>
<b>DNA does not remain in an agarose gel well</b>	<p>In step 5 of Protocol A (step 7 of Protocol B), ensure all residual wash buffer is removed from the membrane. Longer centrifugation time can aid in removal of wash buffer.</p>

## SAFETY INFORMATION



### Binding Buffer

Hazard-determining component of labeling: **guanidinium thiocyanate**

**Xn Harmful**

#### Risk phrases

R20/21/22 Harmful by inhalation, in contact with skin and if swallowed.

R32 Contact with acids liberates very toxic gas.

R52/53 Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

#### Safety phrases

S9 Keep container in a well-ventilated place.

S23 Do not breathe gas/fumes/vapour/spray.

S36/37 Wear suitable protective clothing and gloves.

S60 This material and its container must be disposed of as hazardous waste.

S61 Avoid release to the environment. Refer to special instructions/safety data sheets.



### Binding Buffer

#### Warning

##### Hazard statements:

H302 Harmful if swallowed.

H412 Harmful to aquatic life with long lasting effects.

##### Precautionary statements:

P273 Avoid release to the environment.

P264 Wash thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.

P301+P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

P330 Rinse mouth.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

#### References

1. Vogelstein, B. and Gillespie, D., Preparative and analytical purification of DNA from agarose, Proc. Natl. Acad. Sci. USA, 76, 615-619, 1979.
2. Marko, M.A., Chipperfield, R. and Birnboim, H.C., A procedure for the large-scale isolation of highly purified plasmid DNA using alkaline extraction and binding to glass powder, Anal. Biochem., 121, 382-387, 1982.
3. Boom, R., Sol, C.J.A., et al., Rapid and simple method for purification of nucleic acids, J. Clin. Microbiol., Mar, 495-503, 1990.

#### **LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.**

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

#### **PRODUCT USE LIMITATION**

This product is developed, designed and sold exclusively for research purposes and *in vitro* use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to [www.thermoscientific.com/onebio](http://www.thermoscientific.com/onebio) for Material Safety Data Sheet of the product.

© 2015 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific, Inc. and its subsidiaries.

**CERTIFICATE OF ANALYSIS  
#A27828**

**MagMAX™ mirVana™ Total RNA Isolation Kit Box 2 of 2**

**Quantity:** 96 react.  
**Packaging Lot:** 00302599  
**Storage:** at Room Temperature (15 to 25 °C)

**Expiry Date:** 05.2016 (MM.YYYY)

**The kit is composed of:**

	<b>Lysis Buffer</b>
<b>Quantity:</b>	115 mL
<b>Filling Lot:</b>	00289359
	<b>Wash Solution 1 Concentrate</b>
<b>Quantity:</b>	10 mL
<b>Filling Lot:</b>	00289361
	<b>Wash Solution 2 Concentrate</b>
<b>Quantity:</b>	12 mL
<b>Filling Lot:</b>	00289362
	<b>Rebinding Buffer</b>
<b>Quantity:</b>	4.8 mL
<b>Filling Lot:</b>	00278705
	<b>PK Digestion Buffer</b>
<b>Quantity:</b>	4.4 mL
<b>Filling Lot:</b>	00289363
	<b>Elution Buffer</b>
<b>Quantity:</b>	9.6 mL
<b>Filling Lot:</b>	00289364
	<b>RNA Binding Beads</b>
<b>Quantity:</b>	2 mL
<b>Filling Lot:</b>	00299967
	<b>MagMAX™ TURBO DNase™ Buffer</b>
<b>Quantity:</b>	4.6 mL
<b>Filling Lot:</b>	00289365

## QUALITY CONTROL

Parameter	Method	Requirement	Result
Functional Testing	RNA is isolated from mouse tissue and mouse plasma following the applicable manual extraction procedure described in the MagMAX™ mirVana™ Total RNA Isolation Kit User Guide (publication# MAN0011131). RNA is checked for minimal yield, purity and integrity requirements and suitability for real-time RT-PCR.	Conforms	Conforms

### ISO CERTIFICATION

This product is produced in accordance to Thermo Fisher Scientific Baltics UAB, a wholly owned subsidiary of Thermo Fisher Scientific Inc. (hereinafter Thermo Fisher) certified Quality Management (ISO 9001) and Environmental Management (ISO 14001) systems. Thermo Fisher Scientific Baltics UAB is also certified according to the standards ISO 13485 (Quality management system for IVD medical devices).

### PRODUCT USE LIMITATION

This product has been developed for research purposes and *in vitro* use only. This product has not been tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

### DISCLAIMER OF WARRANTY

Thermo Fisher warrants solely to the Buyer that the Products materially conform to Thermo Fisher's published specifications for such Products at the time of order submission. The warranty period is from the time Thermo Fisher delivers the Product until the earlier of either the Product's expiry or "use by" date or its specified number of or 12 months from the date Thermo Fisher delivers the Product. THERMO FISHER, ITS AFFILIATES AND DISTRIBUTORS EXPRESSLY DISCLAIM ALL OTHER REPRESENTATIONS AND WARRANTIES, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, AND ALL OTHER WARRANTIES, EXPRESS OR IMPLIED INCLUDING, WITHOUT LIMITATION, SATISFACTORY QUALITY AND NON-INFRINGEMENT. THIS LIMITED WARRANTY PROVIDED HEREIN DOES NOT APPLY TO ANY DEFECT CAUSED BY FAILURE TO PROVIDE A SUITABLE STORAGE, USE OR OPERATING ENVIRONMENT; USE OF NON-RECOMMENDED REAGENTS OR BIOCHEMICALS; THE USE OF ANY PRODUCT FOR A PURPOSE OR IN A MANNER OTHER THAN THAT FOR WHICH THEY WERE DESIGNED; MODIFICATIONS DONE BY BUYER; OR ANY OTHER ABUSE, MISUSE OR NEGLIGENCE OF THE PRODUCTS BY BUYER. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. The warranty provided herein is valid only when used by properly trained individuals. All technical advice, documentation and information provided by Thermo Fisher, whether by phone, e-mail, website or any other channel is provided "AS IS" and without any warranty of any kind. It is the responsibility of the Buyer to determine if a Product is suitable for a specific purpose and to apply the necessary safety precautions. Buyer's exclusive and sole remedy under this warranty is: (i) replacement of the Product that failed to conform to the warranty above or (ii) at Thermo Fisher's option or, where applicable, at Thermo Fisher's distributor's option, a credit for that portion of the Product which is not conforming.

### LIMITATION OF LIABILITY

NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, THE LIABILITY OF THERMO FISHER (WHETHER BY REASON OF BREACH OF CONTRACT, TORT, INDEMNIFICATION, OR OTHERWISE, BUT EXCLUDING LIABILITY OF THERMO FISHER FOR BREACH OF WARRANTY (THE SOLE REMEDY FOR WHICH WILL BE AS PROVIDED UNDER SECTION ABOVE)) WILL NOT EXCEED AN AMOUNT EQUAL TO THE TOTAL PURCHASE PRICE PAID BY CUSTOMER TO THERMO FISHER WITH RESPECT TO THE PRODUCT(S) GIVING RISE TO SUCH LIABILITY. NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, IN NO EVENT WILL THERMO FISHER BE LIABLE FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE SALE, USE OR PERFORMANCE OF ANY PRODUCTS (INCLUDING WITHOUT LIMITATION DAMAGES FOR LOSS OF USE OF FACILITIES OR EQUIPMENT, LOSS OF REVENUE, LOSS OF DATA, LOSS OF PROFITS OR LOSS OF GOODWILL, ANTICIPATED OR OTHERWISE), REGARDLESS OF WHETHER THERMO (a) HAS BEEN INFORMED OF THE POSSIBILITY OF SUCH DAMAGES OR (b) IS NEGLIGENT. This limitation of liability does not limit Thermo Fisher's liability for death or personal injury caused by Thermo's negligence or any other liability that cannot be excluded by law.

Rev. 1

**QC responsible:** J. Žilinskienė

**Date:** 24 July 2015