GB Cold light laryngoscopes

Safety Information

- The instrument is delivered unsterilized and must therefore be cleaned, disinfected and/or sterilized before using for the first time as well as before any further use.
- The lamp may be hot so allow it to cool down.
- Observe precisely the rules and procedures for cleaning, sterilization and care of instruments.
- Only use the procedures and products approved by the instrument manufacturer for cleaning, sterilization and care.
- Observe precisely the instructions provided by the manufacturers of these products concerning proper dilution and immersion time.
- Never clean the batteries, insert with water or a solution containing water; protect contacts against moisture.
- If batteries become totally discharged they may leak and thus damage the handle sleeve. For this reason, if they are not used for any length of time remove the battery container or rechargeable battery insert from the handle sleeve.

Cleaning preparation

Begin preparing instruments for use as soon as possible because any incrustations can cause irreparable damage. Remove any coarse dirt using soap and water and a plastic brush. Always rinse using completely demineralized water and then dry afterwards.

Disassembly

If necessary, remove the intubating blade from handle. To clean and sterilize, unscrew the battery from the handle sleeve. In the case of take-apart intubation laryngoscopes remove fiber optic light carrier from the intubating blade.

Manual cleaning, disinfection

Rub down standard batteries using a cloth moistened with alcohol. Remove any traces of dirt from the intubating blade and handle sleeve using a soft brush. Completely immerse and soak the intubating blade, fiber optic light carrier, and handle sleeve in cleaning/disinfectant solution. Afterwards thoroughly rinse with completely demineralized water and dry with sterile compressed air.

Machine cleaning, disinfection

The choice of rinsing programs must be agreed with the manufacturer of the washer. Usually there are special programs for endoscopic instruments (possibly even for telescopes) which are more suitable for cleaning and protecting materials than general programs. Clean the incubating blade and handle sleeves in the washer and disinfectant at up to 93°C.

Ultrasound cleaning

The ultrasound bath is suitable for particularly thorough gentle cleaning of heavily soiled and sensitive parts. Never use cleaning agents which contain chlorine. After cleaning, thoroughly rinse with completely demineralized water and dry with sterile compressed air.

Sterilization

NOTE: If it is necessary to deactivate pathogens which cause Creutzfeldt-Jakob Disease, the steam sterilization method recommended by the WHO (WHO/CDS/PH/92.104) can be used (18 minutes at 134°C). The sterilization time can be extended up to an hour at 134°C in special cases.

WARNING: Remember that sterilization can only be carried out successfully if recognized and validated methods are used.

Method

Steam sterilization: 134°C, P ₀ = 2 bar	
prevacuum method	5 min
fractionated prevacuum	5 min
gravitation method	30 min
Gas sterilization (formaldehyde, ethylene oxide)	
Plasma sterilization (e.g. STERRAD®)	
Chemical sterilization (e.g. STERIS®)	

Assembly

Check whether the batteries or rechargeable battery insert are in good condition: replace batteries which are empty, leaking or damaged (observe polarity); if in doubt, replace. Check whether the lamp is firmly in its holder, and tighten up if necessary. Screw the battery insert into the handle sleeve. Fit the intubating blade and check operation.

Replacement of nonrepairable parts

When replacing nonrepairable parts such as the lamp always use genuine spare parts. Fitting parts of other makes can lead to malfunctions for which the manufacturer is not liable.

Lamp

Unscrew the battery insert from the handle sleeve. Remove faulty lamp. Check to make sure the new lamp has the correct voltage and the insulation collar on the threaded base of the new lamp is present. Insert new lamp. Fit battery into the handle sleeve. Fit the intubating blade and check operation.

Batteries

Unscrew the battery from the handle sleeve. Always replace both batteries. New batteries observing correct polarity. Insert the battery into the handle sleeve. Fit the intubating blade and check operation. Dispose of old batteries in an environmentally friendly way.

Operating/storage conditions

	Temperature	rel. humidity
storage	-20°C ... +60°C	10% ... 90%
operation	+10°C ... +40°C	30% ... 70%

Disposal

To dispose these medical products, no special measures are necessary. National laws and regulations must be observed. Dispose of old batteries which can no longer be used in an environmentally friendly way: pass on to an official collection center.

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