# Medtronic

## Sentrant™

Introducer Sheath with Hydrophilic Coating  
Интродьюсерно дезиле с хидрофилно покритие  
Zaváděcí pouzdro s hydrofilní vrstvou  
Indføringsheath med hydrofil belægning  
Schleuse mit hydrophiler Beschichtung  
Θηκάρι εισαγωγής με υδρόφιλη επίστρωση  
Introductor con revestimiento hidrófilo  
Hüdrofiilse kattedga sisestuskanüül  
Sisäänvientiholkki, jossa on hydrofiilinen pinnoite  
Gaine d'introducteur avec revêtement hydrophile  
Uvodnica s hidofilnim premazom  
Hidofil bevonattal ellátott bevezetőhüvely  
Selubung Pengantar dengan Pelapis Hidrofilik  
Introduttore con rivestimento idrofilo  
Гидрофильдик қабат бар интродьюсер қабықшасы  
친수 코팅된 인트로듀서 쉬스  
Įvedimo prietaiso mova su hidrofiline dangą  
Ievadslūžas ar hidrofilu pārklājumu  
Обвивка за воведување со хидрофилен слој  
Introducersheath met hydrofiele coating  
Innføringshylse med hydrofilt belegg  
Koszulka naczyniowa z powłoką hydrofilną  
Bainha introdutora com revestimento hidrofílico  
Teacă de introducere cu strat hidofil  
Интродьюсер с гидрофильным покрытием  
Zavádzacie puzdro s hydrofilnou vrstvou  
Uvajalni tulec s hidofilno prevleko  
Omotáč uvodníka sa hidofilnom oblogom  
Införingshylsa med hydrofil beläggning  
Hidrofilik Kaplamalı İntrodüser Kılıf  
Інтрод'юсер із гідрофільним покриттям  
Ống dẫn có lớp phủ ưa nước

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# Contents

<b>1 Device description</b>	<b>12</b>
<b>2 Intended use</b>	<b>12</b>
2.1 Intended users	12
2.2 Indications for use	12
2.3 Contraindications	13
2.4 Clinical benefits	13
2.5 Intended patient population	13
2.6 Device performance characteristics	13
<b>3 Warnings and precautions</b>	<b>13</b>
<b>4 Adverse events</b>	<b>13</b>
4.1 Potential adverse events	13
<b>5 How supplied</b>	<b>14</b>
5.1 Sterility	14
5.2 Contents	14
5.3 Storage	14
5.4 Device disposal	14
<b>6 Clinical use information</b>	<b>14</b>
6.1 Physician training requirements	14
6.2 Device inspection	14
6.3 Additional recommended equipment (not included with this device)	14
<b>7 Instructions for use</b>	<b>14</b>
<b>8 Disclaimer of warranty</b>	<b>16</b>


## 1 Device description

The Sentrant introducer sheath with hydrophilic coating (*Figure 1*) is a single-use, disposable, hydrophilic catheter that provides a flexible and hemostatic conduit for the insertion of endovascular devices and to minimize blood loss associated with vascular procedures. The system is comprised of 2 components: a dilator and an introducer sheath, which accommodate a 0.035 in (0.89 mm) guidewire.

The dilator is radiopaque and has a tapered, flexible tip that facilitates atraumatic tracking through the vasculature. A female Luer taper fitting is located on the proximal end of the dilator grip. The distal end of the dilator grip is threaded to allow the dilator to be secured to the sheath seal housing.

The introducer sheath is comprised of a hydrophilic, coil-reinforced catheter that is attached to a rigid seal housing containing the hemostatic valve assembly. A sideport extension with a 3-way valve is permanently attached to the seal housing. A radiopaque markerband is located at the distal tip of the sheath. The device also has a suture loop for attaching it to the patient and a strain relief to prevent kinking of the catheter where it joins to the seal housing.

**Figure 1.** Sentrant introducer sheath with hydrophilic coating

- |  |                                       |
|--|---------------------------------------|
| 1 Dilator  | 6 RO marker band                      |
| 2 Dilator grip   | 7 Sideport extension with 3-way valve |
| 3 Tapered tip  | 8 Suture loop                         |
| 4 Introducer sheath  | 9 Strain relief                       |
| 5 Seal housing with hemostatic valve assembly (sheath hub) | 10 3-way valve handle                 |

**Note:** Graphical representation not drawn to scale.

This device does not contain natural rubber latex; however, during the manufacturing process, it may have incidental contact with latex.

## 2 Intended use

The intended use of the Sentrant introducer sheath is to provide a conduit for the insertion of diagnostic or endovascular devices into the vasculature.

### 2.1 Intended users

The Sentrant introducer sheath should only be used by physicians and teams trained in vasculature interventional techniques.

This device should only be used in a clinical setting in a sterile environment.

### 2.2 Indications for use

The Sentrant introducer sheath with hydrophilic coating is intended to provide a conduit for the insertion of diagnostic or endovascular devices into the vasculature and to minimize blood loss associated with such insertions.

## 2.3 Contraindications

There are no known contraindications with this device. It is not intended for use except as indicated.

## 2.4 Clinical benefits

This device can minimize blood loss associated with the insertion of endovascular devices or diagnostics.

**Note:** After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 0763000B00007398C.

## 2.5 Intended patient population

The intended patient population for the Sentrant introducer sheath is any patient undergoing a vascular procedure in which an endovascular device or diagnostic will be inserted.

This device has not been tested on a pediatric population.

## 2.6 Device performance characteristics

The Sentrant introducer sheath is a flexible hemostatic conduit to a patient's vasculature. The product is introduced to the body using a dilator with a tapered tip which is passed over a 0.035 in (0.89 mm) guidewire. Once in the desired position, the dilator is removed through the Sentrant introducer sheath allowing for the introduction of endovascular stent graft delivery systems and complementary accessories.

## 3 Warnings and precautions

**Caution:** Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Do not alter this device. Alterations may impair device function.
- Do not cut or puncture the catheter. This could result in major blood loss and vessel trauma.
- Do not use after the Use-by date printed on the label.
- Do not use if the device is damaged.
- After use, dispose of the device in accordance with hospital, administrative, or government policies.
- The catheter is only compatible with a 0.89 mm (0.035 in) guidewire.
- Verify the size compatibility of the sheath, device, catheter, and accessory components prior to use. The maximum diameter and length of the instrument or catheter to be introduced should be determined to ensure its passage through the sheath.  
**Note:** The working length of the sheath does not include the strain relief feature on the catheter.
- If scar tissue is present, predilatation may be required.
- Adequate vessel access is required to introduce the sheath into the vasculature. Careful evaluation of vessel size, anatomy, tortuosity, and disease state (including calcification, plaque and thrombus) is required to ensure successful sheath introduction and subsequent withdrawal. If the vessel is not adequate for access or the vessel size is smaller than the introducer sheath's outer diameter, major bleeding, vessel damage, serious injury, or death may result.
- When inserting, manipulating, or withdrawing a device through the sheath, always maintain sheath position.
- When inflating a balloon, ensure the balloon is completely outside of the Sentrant device.
- When puncturing, suturing, or incising the tissue near the sheath, use caution to avoid damaging the sheath.
- Do not attempt to insert or withdraw the introducer sheath if resistance is felt. The cause of the resistance must be determined before proceeding.
- Do not attempt to simultaneously insert multiple catheters or interventional devices into the sheath.
- Advance and withdraw the sheath (with the dilator fully inserted and locked in) only under fluoroscopic guidance.
- Individual patient anatomy and physician technique may require procedural variation.

## 4 Adverse events

### 4.1 Potential adverse events

Adverse events or complications associated with use of the Sentrant introducer sheath include, but are not limited to:

- Allergic reaction
- Blood loss, bleeding, or hematoma

- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Infection
- Vascular trauma (eg, dissection, rupture, perforation, or tear)
- Death

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

## 5 How supplied

### 5.1 Sterility

The Sentrant introducer sheath is individually packaged. It is supplied sterile (via ethylene oxide or electron beam) for single use only.

- Do not reuse or attempt to resterilize.
- If the device is damaged or the integrity of the sterile barrier has been compromised, do not use the product and contact your Medtronic representative for return information.

### 5.2 Contents

- One Sentrant introducer sheath with dilator
- One web reference leaflet

### 5.3 Storage

Store the Sentrant introducer sheath at room temperature in a dark, dry place to avoid extended exposure to light and moisture.

### 5.4 Device disposal

Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

## 6 Clinical use information

### 6.1 Physician training requirements

**Caution:** The Sentrant introducer sheath should only be used by physicians and teams trained in vascular interventional techniques.

### 6.2 Device inspection

Carefully inspect the device and packaging for damage or defects prior to use. If the Use-by date has elapsed, the device is damaged, or the sterilization barrier has been compromised, do not use the device; contact a Medtronic representative for return or replacement.


### 6.3 Additional recommended equipment (not included with this device)

- 0.035 in (0.89 mm) exchange guidewire
- Heparinized saline solution

## 7 Instructions for use


1. Remove the device from the package and ensure that the inner diameter (ID) of the sheath is appropriate for the maximum diameter of the instrument or catheter to be introduced.
2. Verify that the vessel is of adequate diameter and the tortuosity is not too severe to accommodate the introducer sheath.
3. Prepare the sheath for flushing and ensure that the 3-way valve is open. There are two options to do this. See *Figure 2*.

**Figure 2.** 3-way valve in open positions

4. Remove the cap from the open port on the 3-way valve.
5. Connect a syringe of heparinized saline solution to the open port and flush the sheath (tapping the introducer sheath to aid in releasing air bubbles).
6. Close the 3-way valve. See *Figure 3*.

**Figure 3.** 3-way valve in closed position

7. Connect a syringe of heparinized saline solution to the Luer fitting on the dilator and flush the dilator.

8. Insert the dilator completely into the introducer sheath and secure it by rotating the dilator grip clockwise, approximately ¼ turn.
9. Follow accepted clinical practice for vessel access and guidewire insertion.
10. Activate the hydrophilic coating by wetting the outer surface of the sheath with heparinized saline.  
**Note:** For best results, ensure the device remains wet during placement.
11. Under fluoroscopic guidance, advance the sheath and dilator as a unit over the guidewire to the desired location. Stop advancement if there is resistance and investigate the cause of the resistance before proceeding.
12. While holding the sheath steady and maintaining guidewire position, unlock and withdraw the dilator from the sheath until the dilator is completely removed from the sheath and guidewire.
13. While maintaining guidewire position, advance the selected interventional device over the guidewire, through the valve, and into the sheath. Stop advancing the interventional device if there is resistance and investigate the cause of the resistance before proceeding.
14. Interventional devices can be interchanged as necessary.
15. Upon removal of the sheath, take care to prevent bleeding, vessel damage, or other serious injury. Advancing the dilator into the sheath prior to withdrawing the sheath may ease withdrawal from the vasculature.

## 8 Disclaimer of warranty

**Note: This disclaimer of warranty does not apply in Australia.**

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
**Medtronic, Inc.**
**Medtronic B.V.**
Medtronic Limited