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DiaClon Rh-Subgroups + K	ID-Card
English	B002124 03.13
C, c, E, e, K, ctl monoclonal Determination of Rh phenotypes and K	
Product-Identification: 50110	

INTRODUCTION
Besides the RhD (RH1) antigen, other important antigens of the Rh system (ISBT, number 004) are: C (RH2), E (RH3), c (RH4) and e (RH5).
According to Issitt [1] their frequencies in the Caucasian population are as follows:
C.....70% c.....80% E.....30% e.....98%

Appropriate antigen-positive red cells may stimulate antibody production in antigen C (RH2), c (RH4), E (RH3) and e (RH5) negative individuals. The determination of the Rh phenotypes can, therefore be important during pregnancy, for previously transfused patients and for patients with known irregular antibodies [2].

Approximately 9% of the Caucasian population are K (KEL1, K, K1) positive. The K antigen is strongly immunogenic. Anti-K has been reported as the cause of hemolytic transfusion reactions, both immediate and delayed, and hemolytic disease of the newborn.


The ID-Card "DiaClon Rh-Subgroups + K" offers a complete profiling of the Rh phenotype and includes K typing.

REAGENTS

IVD

ID-Card "DiaClon Rh-Subgroups + K" contains monoclonal antibodies anti-C (cell line MS-24), anti-c (cell line MS-33), anti-E (cell line MS-260), anti-e (cell lines MS-16, MS-21, MS-63) and anti-K (cell line MS-56) within each gel matrix. The microtube (ctl) is the negative control. Preservative: < 0.1% NaN₃.

Caution: All reagents should be treated as potentially infectious.

ADDITIONAL REAGENTS REQUIRED

ID-Diluent 2: modified LISS for red cell suspensions.

(see related package insert)

FURTHER MATERIALS REQUIRED

- ID-Dispenser
- ID-Pipetor
- ID-Tips (pipetor tips)
- Suspension Tubes
- ID-Working table
- ID-Centrifuge 6, 12 or 24

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.

PREPARATION OF BLOOD SAMPLE

Prepare a 5% red cell suspension in ID-Diluent 2 as follows:

Allow the diluent to reach room temperature before use.

- Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
- Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

The cell suspension may be used immediately.

CONTROLS

Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance.

TEST PROCEDURE

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.

- Identify the ID-Card with the unique patient or donor number/details as appropriate.
- Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
- Add 10 or 12.5 µl of the patients' red cell suspension to all microtubes of the ID-Card.
- Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
- Read and record the results.

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