WIP'ANIOS Excel

The choice of **safety** into a **suitable format**

Make your choice!

CLEANSING DISINFECTANT WIPES

Wip'Anios Excel Travel 25

- 1 24 packets 25 wipes each
Size : 130 x 210 mm
Texture : 100% viscose
Ref. 2617.466

Wip'Anios Excel 50 et 100

- 3 12 packets 50 wipes each
Size : 180 x 200 mm
Texture : 100% viscose
Ref. 2446.424

Wip'Anios Excel King Size

- 5 1 bucket containing 1 refill of 350 wipes
Size : 200 x 360 mm
Texture : 70% viscose 30% polyester
Ref. 2446.672

Wip'Anios Excel 60

- 2 12 packets 60 wipes each
Size : 240 x 300 mm
Texture : 40% viscose 60% polyester
Ref. 2446.492

Wip'Anios Excel 100

- 4 6 packets 100 wipes each
Size : 180 x 200 mm
Texture : 100% viscose
Ref. 2446.655

Wip'Anios Excel 350

- 6 2 refills of 350 wipes
Size : 200 x 360 mm
Texture : 70% viscose 30% polyester
Ref. 2446.673

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EN2016-08 - Non contractual pictures. Dangerous- respect the precautions of use. Medical Devices class I/II. Use biocides with caution. Before using, read the label and information relative to the product. Biocidal product.

WIP'ANIOS Excel

WIP'ANIOS EXCEL wipes combine cleaning and disinfection in a single step making it possible to ensure **macroscopic and microscopic hygiene of high surfaces and medical devices**, at any time, in any environment (operating room, dialysis department, environment near the patients, medical office, interiors of medical vehicles, etc.).

✓ Cleansing disinfectant wipes

- Impregnation solution: Surfa'safe Premium.
- Efficiency ensured until the expiry date.
- Complies with EN 16615 phase 2/step 2.

✓ Ready to use

- Uniform and controlled impregnation, no dilution in qualitative water which is not controlled.
- Assured traceability.

✓ Suitable for all types of surfaces

- In particular, haemodialysis generators, including those that may come into contact with foodstuffs (Group 1 - PT2 & PT4).
- Class IIb non-immersible medical devices.
- Wiping endoscope sheaths in preliminary treatment.
- Wide compatibility with a wide range of materials.

✓ Formats suited to different uses

- Available in packets and buckets.
- 4 different sizes.

SAFETY GUARANTEED

MULTI-USE

ENSURED VERSATILITY

For the safety of your patients, operators and equipment.

SAFETY GUARANTEED

FOR THE USER

The impregnation solution is formulated :

- Without CMR substances
- **Alcohol free**
- Fragrance free

for **optimal use tolerance and comfort.**

FOR EQUIPMENT

The wipes are **impregnated with a Surfa'safe Premium solution**: its neutral pH and lack of amphoteric surfactants ensure maximum compatibility with a wide range of materials (metals, polymers, polycarbonate, aluminium, etc.).

FOR THE ENVIRONMENT

The **100% biodegradable viscose** of WIP'ANIOS EXCEL (50 or 100) and WIP'ANIOS Excel Travel wipes makes it possible to limit the impact on the environment.

VALIDATED ANTIMICROBIAL EFFICACY

- > Cleansing disinfectant wipes validated according to EN 16615.
- > Bactericide, fungicide and virucide impregnation solution.

ACTIVE AGAINST	STANDARDS	CONTACT TIME
BACTERIAS	EN 13727	2 minutes
	EN16615	2 minutes
	EN 13727 on VRE, A.baumannii β LSE, E.coli OXA 48, E.cloacae OXA 48, K.pneumoniae OXA 48, MRSA	1 minute
	EN 13697 sur ERV, A.baumannii β LSE, E.coli OXA 48, MRSA	2 minutes
YEASTS	EN 13624	2 minutes
	EN16615	2 minutes
VIRUS	PRV, (surrogate of HBV), BVDV (surrogate of HBV), Vaccinia virus, HSV, Rotavirus	1 minute

STANDARD EN 16615

> The application standard phase 2/step 2 **EN 16615** (May 2015) makes it possible to claim activity combining the action of the wipe and the solution by measuring the **bactericidal and yeasticidal efficacy** of a wipe pre-impregnated with a disinfectant cleansing solution.

Test condition: Quantitative surface test with mechanical action using wipes (4 zone test) in 1 to 5 minutes (or 60 minutes), at a temperature between +4°C and +30°C, diluted in 30°F hard water for concentrated products and in the presence of interfering substances (CP: 0.3g/L; CS: 3g/L albumin + 3mL/L erythrocytes).

A significant reduction must be established to claim the "cleansing disinfectant wipe" status: 5 log reduction for bacterias - 4 log reduction for yeast.

MAKE YOUR CHOICE !

Wip'Anios Excel Travel 25

Size : 130 x 210 mm
Texture : 100% viscose

Wip'Anios Excel 50 ou 100

Size : 180 x 200 mm
Texture : 100% viscose

Wip'Anios Excel 60

Size : 240 x 300 mm
Texture : 40% viscose 60% polyester

Wip'Anios Excel King Size

Size : 200 x 360 mm
Texture : 70% viscose, 30% polyester

➤ Microbiological quality managed :
Radio-sterilized wipes for an optimum safety.

Particularly suitable for haemodialysis generators

PERFECT FOR ...

- **Small medical devices** (pressure sensors, blood sugar testers, tensiometres, reflex hammers, conventional ultrasound probes)
- **Wiping endoscopes** preliminary treatment before disinfection
- **Surfaces in the vicinity of patients** (adaptable, bed, stretcher)

PERFECT FOR ...

- **Medical equipment and devices**
- Surfaces requiring wipes in a larger format (external device structures, examination table, dialysis equipment, etc.)

READY-TO-USE PRODUCTS

WIP'ANIOS EXCEL

- » 100% Biodegradable viscose wipes.
- » Impregnating solution formulated without fragrance nor alcohol.
- » Broad virucidal spectrum.

INSTRUCTIONS FOR USE

SINGLE USE WIPES

Active from 2 min.

INDICATIONS

Surface treatment (e.g. external structures of equipment: respirators, dialysis generators, etc.), non-immersible and non-invasive* medical devices, and wiping of endoscopy sheaths in preliminary treatment, prior to immersion for cleaning and disinfection. *Examples: conventional ultrasound probes, pressure sensors, glucose testers, tensiometres reflex hammers, etc.

COMPOSITION

Wipes impregnated with a cleaning and disinfecting solution: **Didecylmethylammonium chloride**, non-ionic surfactants, sequestering agent, excipients.

PACKAGING - REF.

Wip'Anios Excel (180 x 200 mm)	
6 packs, 100 wipes each.....	2446.655
12 packs, 50 wipes each.....	2446.424
Wall holder.....	425.185

WALL HOLDER

THE IMPREGNATION SOLUTION IS:

ACTIVE AGAINST	ACCORDING TO STANDARDS
BACTERIA	EN 13727, EN 1276, EN 13697, EN 16615 Active against MRB and XRB according to EN 13727
YEASTS	EN 13624, EN 16615, EN 13697
MOULDS	EN 13624, EN 13697
VIRUSES	Adenovirus, Norovirus MNV according to EN 14476 HIV, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus, HSV, Rotavirus, VRS, Polyomavirus SV 40, Coronavirus

MORE DETAILS IN THE DEDICATED BROCHURE

SURFA'SAFE PREMIUM

- » Formulated without alcohol, nor perfume.
- » Fast antimicrobial efficacy: bacteria including XRB, Rotavirus...
- » High compatibility with materials including polymers and polycarbonates...
- » Dense foam to limit aerosol droplets during diffusion.

INSTRUCTIONS FOR USE

READY-TO-USE
Active from 2 min.

INDICATIONS

Disinfecting detergent foam for cleaning and disinfection of surfaces, external equipment structures (respirators, haemodialysis generators, etc.) and non-immersible and non-invasive medical devices (conventional ultrasound probes, pressure sensors, blood sugar testers, blood pressure monitors, reflex hammers ...).

COMPOSITION

Didecylmethylammonium chloride, excipients.

PACKAGING - REF.

12 bottles, 750ml each, with foaming dispensers.....	2419.544
16 bottles, 500ml each	2419.105

ACTIVE AGAINST	ACCORDING TO STANDARDS
BACTERIA	EN 13727, EN 13697, EN 16615 MRB and XRB according to EN 13727, EN 13697
YEASTS	EN 13697, EN 13624, EN 16615
MOULDS	EN 13624, EN 13697 Aspergillus fumigatus and Tricophyton mentagrophytes according to EN 13624, EN 13697
VIRUSES	EN 14476 : Adenovirus, Norovirus MNV HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus, HSV, Rotavirus, VRS, Polyomavirus SV 40, Coronavirus according to EN 14476

Our advice

For large surfaces, preferably use the 500ml bottle.

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