


# Quark C12x-T12x

Stress ECGs

Diagnostic Quality 12-lead Stress Testing ECGs available in telemetry or cable configuration

“When heart and lung meet together..”



**COSMED**  
The Metabolic Company

“  
*Electrocardiographic (ECG) recording system for continuous monitoring of heart rhythm and evaluation of ischemic ECG changes during exercise and recovery*<sup>(1)</sup>”

- **Resting and Exercise ECG interpretation**
- **High resolution ECG processing**
- **Full Disclosure and Scroll back during the test**
- **ST segment, ST slopes, Trends**
- **Retrospective ECG and arrhythmia analysis even during test**
- **Standard and user defined Exercise Protocols**
- **Integration with COSMED Gas exchange equipment**
- **Patient Database & Network Compatibility**



COSMED Stress ECGs are the ideal diagnostic quality PC-based ECG devices for patient care from prevention to rehabilitation. Ultimate quality ECG signal allows an efficient diagnosis of coronary artery disease, the evaluation of patients after artery bypass surgery or angioplasty, and allows to document response to therapies and prevent from myocardial infarction.

Combined with any COSMED gas exchange analysis equipment, COSMED Stress ECGs become the complementary tools that help clinician to quantify patient's functional capacity and evaluate cardiac response in athletes. COSMED PC-based stress testing systems offer the ultimate in stress testing flexibility, meeting the most advanced requirements for clinical exercise testing. Two versions are currently available:

- **Quark C12x** - 12-lead Patient Cable Stress Test ECG, with USB communication interface.
- **Quark T12x** - 12-lead Wireless Stress Test ECG (Telemetry). State of the art Digital Radio Frequency technology, immune to environmental interference. Optimized for working in noisy EM Environments (Bluetooth, Wi-Fi Networks, Cellphones etc.) . It includes an automatic Channel Shift Feature for selecting the optimal Transmission Channel.

## Design

- **12 leads ECG** 12-lead simultaneous recording, display and printout on plain paper. Detachable 10 lead wires conform to AAMI (AHA) or IEC color codes. Every single lead can be replaced independently
- **High resolution of signal processing** (up to 500 samples/sec) guarantees exceptionally clear tracing for visual view and on screen measurements
- **On-line filtering** Advanced digital filters and automated baseline correction for real-time and stored ECGs to create exemplary ECG traces
- **Integration with COSMED gas exchange equipment** stress ECG and ergospirometry parameters can be thus recorded synchronously, identifying cardiac, pulmonary and vascular function in just a single measurement



*Integrated Cardio Pulmonary Exercise Testing: Quark C12x with a Quark CPET*



*Quark C12x and T12x can also be used as standalone Stress Test ECG Measuring Stations*

(1) AHA Guidelines for Clinical Exercise Testing Laboratories. Circulation. 1995; 91: 912-921



## Technical Specifications

Product	Description	REF
Quark C12x	Diagnostic quality 12 lead stress testing ECG with patient cable	C09080-01-99
Quark T12x	Diagnostic quality wireless 12 lead stress testing ECG	C09081-01-99
<b>Standard Tests</b>		
ECG	Stress Testing ECG, Resting ECG (w/ Interpretation), Full Disclosure, Scroll Back, St Segment, St Slopes, Trends, Arrhythmia Analysis w/stress	
<b>Features</b>		
Sampling rate	Processing: 500 s/sec - Acquisition: 16000 s/sec (C12x), 4000 s/sec (T12x)	
Defibrillator protection	to 360J	
CMRR	>100 Db	
Frequency range	0.05 - 300 Hz	
Patient leakage current	<10 µA	
Input impedance	>100 Mohm	
Pacemaker detection	From 0.1 to 2ms to 2-700mv	
<b>Electrodes</b>		
Number of electrodes	10 AHA / IEC	
Acquisition	Simultaneous 12-leads (I,II,III, aVR, AVL, aVF, V1, V2, V3, V4, V5, V6)	
Speed of display	5, 12.5, 25, 50, 100 mm/sec	
<b>Communication</b>	<b>Quark C12x</b>	<b>Quark T12x</b>
Communication	USB cable	Digital Radio Frequency
Radio frequency range (Auto-select)	2400 - 2483.5 Mhz	
Radio frequency output power	0.4 mW	
Maximal operational distance	~ 10 m	
<b>Hardware</b>	<b>Quark C12x</b>	<b>Quark T12x</b>
Dimensions & Weight	17 x 9 x 3 cm / 300 gr	16 x 10 x 4.5 cm / 350 gr
Power supply	Via USB (5V) (Consumption: < 200 mA ± 10%)	Battery (2 AA Alkaline Or NiMH Rechargeable)
PC Interface	Cable (USB 2.0)	Telemetry (Wi-Fi)
<b>Software</b>	<b>ECG Suite</b>	
Available languages	English, French, Spanish, German, Dutch, Italian	
PC Configuration	i5 processor required, RAM: minimum 2 GB, recommended 4 GB. OS: Windows (from XPSP3 to Windows 10)	
User defined exercise protocols	Bruce, Balke, Ellestad, Kattus and USAFSAM	
Display mode	3x4 / 6x2 / 12 leads with or without other parameters	
Export formats	HL7, XML, JPEG, GDT, PDF, DICOM	
<b>Accessories &amp; Options</b>	<b>Description</b>	<b>REF</b>
Additional Software	Software Stress, Rest ECG with Interpretation.& Patient Database	A-661-200-058
<b>Safety &amp; Quality Standards</b>		
MDD (93/42 EEC); FDA 510(k); EN 60601-1 (safety) / EN 60601-1-2 (EMC)		



**COSMED Srl**

Via dei Piani di Monte Savello 37  
Albano Laziale - Rome 00041, Italy

+39 (06) 931-5492 Phone  
+39 (06) 931-4580 Fax

info@cosmed.com | cosmed.com

To know more:

