

Exadrop®

en	Instructions for use	4
de	Gebrauchsanweisung	5
bg	Инструкция за употреба	7
cs	Návod k použití	8
da	Brugsanvisning	10
el	Οδηγίες χρήσης	11
es	Instrucciones de uso	13
et	Kasutusjuhised	14
fi	Käyttöohje	16
fr	Mode d'emploi	17
hr	Upute za uporabu	19
hu	Használati utasítás	20
it	Istruzioni per l'uso	22
kk	Пайдалану нұсқаулығы	23
ko	사용 지침	25
lt	Naudojimo instrukcija	26
lv	Lietošanas instrukcija	27
nl	Gebruikersinformatie	29
no	Bruksanvisning	30
pl	Instrukcja użytkowania	31
pt	Instruções de utilização	33
ro	Instrucțiuni de utilizare	34
ru	Инструкция по применению	36
sk	Návod na použitie	38
sl	Navodila za uporabo	39
sr	Uputstvo za upotrebu	40
sv	Bruksanvisning	42
tr	Kullanım Kılavuzu	43
uk	Інструкції для застосування	44
vi	Hướng dẫn sử dụng	46
zh	使用说明	47

B | BRAUN

Exadrop®

en Instructions for use

Device description

Articles equipped with PVC tube:

ID: 3 mm
OD: 4.1 mm
Shore A: 74

Articles equipped with PVC-free tube (incl. UV-protected tubes):

ID: 3 mm
OD: 4.1 mm
Shore A: 87

Exadrop® fulfils the standard EN ISO 8536-13 / ISO 8536-13 for graduated flow regulators which enables the setting or adjustment of a certain flow rate.

Flow rate accuracy:

Setting	Wheel position	Minimum (ml/h)	Maximum (ml/h)
High	"250"	175	325
Medium	"70"	45,5	95,5
Low	"20"	12	28

IV administration sets can be used with devices which are in compliance with EN ISO 80369-7 / ISO 80369-7, EN ISO 80369-20 / ISO 80369-20.

IV administration sets are in accordance with EN ISO 8536-4 / ISO 8536-4.

Precision flow control device with/ for I.V. administration set. Only for gravity infusion.

Sterility

Sterilized EO. Please refer to the label of primary packaging

Do not use if package is damaged.

Visually inspect the product packaging to ensure that the sterile barrier system is intact.

Intended purpose

Single-use, sterile flow control device with/ for I.V. administration set of infusion/drugs by gravity – indicated by symbol "G" (Gravity) on primary packaging.
Do not re-use.

Indication

Infusion sets are used for infusion therapy e.g.: Fluid volume maintenance, fluid volume replacement, infusion of total parenteral nutrition, infusion of medications (e.g. antibiotics, chemotherapeutic agents) with bottles and bags according to the SPC (Summary of product characteristics) of the drugs/ solutions.

Patient population

The device can be used for all patients for which infusion therapy is prescribed. No gender or age-related limitations. Exadrop® can be used for adults, pediatric and neonates.

Intended user

To be used by country-specific authorized persons e.g.: Healthcare professionals licensed and/or certified physicians, healthcare practitioners, care givers after instruction of a healthcare professional (according to local regulations). Exadrop® should only be used by healthcare worker who have been adequately trained within the educational background in this technique. Upon medical assessment and adequate instruction, patients and/or caregivers may also be allowed to take over defined handling steps

according to national guidelines. Healthcare workers are responsible for ensuring that patients and/or caregivers are instructed in the correct use of Exadrop® and that the instruction of caregiver and/or patient is documented by the healthcare worker according to local documentation guidelines.

Contraindications

Drugs/ Solutions which are proven incompatible to one of the used materials.

No administration of blood products or blood components.

SPCs of the drugs/solutions for contraindications should be considered.

Do not use for pressure infusion > 0,5 bar (50 kPa).

Residual Risks/Side effects

General risk occurring during infusion therapy including air embolism, local- to systemic infection, particulate contamination. Leakage of drugs/fluids can be an additional hazard to patients, visitors or staff.

Precautions

Do not sterilize.

Do not use with pressure infusion devices (infusion pumps).

Visual check of the device for damages before use.

Do not use if protective caps are missing or loose.

Please also consider the corresponding SPCs of the drugs/ solutions regarding possible incompatibilities of drugs/ solutions.

Mind the general guidelines for parallel infusions.

Monitor all connections for tightness over the whole infusion procedure.

The figures on the scale (= app. ml/h) are only approximate ones. They only refer to the infusion of 0.9 % NaCl solution administered by means of an Vasofix® Braunüle® G 18 being approximately 100 cm static difference of level without CVP.

Therefore, the figures on the scale must always be controlled by verification. Otherwise, there is a risk of over- or underdosage for the patient with all subsequent consequences.

Use aseptic techniques in accordance to national guidelines and/or hospital protocols for variants with Safeflow.

After medication administration, flush per institutional protocol, attaching completely between flushes, for variants with Caresite®.

Warning

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/ or limited functionality of the device may lead to injury, illness or death of the patient. Administration of light sensitive drugs requires UV-protected IV administration sets. Always check, if incompatibilities with other drugs or with materials is stated in the SPC (Summary of product characteristics) of the drug.

Do not re-attach the protective cap of the patient connector. Do not re-sterilize. Do not use after expiration date.

Variants with Safeflow:

Only use connection devices in compliance with ISO 80369-7. Syringes or connectors which are not in compliance with ISO 80369-7/ EN ISO 80369-7 can damage the swabable valve. Syringes and male luer connectors have a large variety of configurations and can vary significantly in design and dimensions.

Do not use, in any situation, needles in the needle-free injection port.

Operating Instructions

1. Open the single package and take out the product.
2. Unfold the line & inspect the IV-Set or the graduated flow regulator in-line set and check whether the slide clamp and flow regulator is opened. Check whether the green air vent cap is closed (only for IV-Set with graduated flow regulator). (Fig. 1 & 2)
3. Connect to accessories (e.g. 3 -way stopcocks, manifolds, injection sites) or connect graduated flow regulator in-line set to IV-Set.
4. Remove protective cap from spike (only for IV-Set with graduated flow regulator)
5. Insert the spike into the container (only for IV-Set with graduated flow regulator)
Insert administration set with closed snap vent in fluid container. (Fig. 3)
6. Close Flow regulator or graduated flow controller make sure the slide clamp is open. Rotate wheel of graduated flow controller towards "OPEN" in order to release initial position. Set graduated flow controller to "OFF". (Fig. 4)
The loop on the Exadrop®-handle is used as tube holding device when the connection of the I.V. administration set to the I.V. cannula cannot be established immediately.
8. Hang the container on the stand
9. Set up fluid level in drip chamber (Fig. 5)
10. Open Flow regulator (if applicable) and graduated flow regulator.
Set graduated flow regulator to "OPEN". (Fig. 6)
11. Open venting cap (only applies to glass bottles and device variants with venting flaps) (Fig. 7 & 8)
12. Take patient connector from hook and prime the IV line (including accessories (e.g. stopcocks or injection sites) and graduated flow regulator) – air free (Fig. 9)
13. If applicable, follow special priming procedures for accessories e.g. for IV Filter and Flushing sets
14. Close the flow regulator and/or the graduated flow regulator AND venting cap (for glass bottle use only and for device variants with venting flaps). (Fig. 10)
15. Connect to prepared and primed accessories if applicable (e.g. extension line, flushing sets secondary line)
16. Remove patient connector / tubing from flow regulator (only for IV-Set with flow regulator)
17. Remove protective cap from patient connector
18. Lock patient connector to the venous access (Fig. 11)

19. Open Flow Regulator (if applicable) and/or graduated flow regulator AND venting cap (for glass bottle use only), and start infusion (Fig. 12)
20. Adjust flow rate by turning the handle and monitor infusion rate
Set desired flow rate on the markings of the handle by working back from "OPEN".
21. Control the drop rate during infusion
22. Important note:
When the prepared IV-system is laid down horizontally, e.g. on a tablet, or the spiked glass container is placed on a table, it is necessary to close the air vent cap beforehand! Only re-open the air vent cap when the infusion is running.
Always ensure to establish a 20 cm siphon (Fig. 12 Et 13).
23. Stop of therapy: interrupt infusion by closing flow regulator or by closing slide clamp or by closing graduated flow regulator
With the slide clamp above the control the infusion can be interrupted without setting back the control device to "OFF".
After having opened the clamp the flow rate remains the same if the pressure conditions in the infusion system have not changed in the meantime. (Fig. 14)

Duration of use

Duration of use is dependent on the intended therapy according to the SPC of the drug or solution. In general, change of administration sets should occur according to national guidelines (e.g. CDC guideline) and/or hospital protocols.

Disposal

Disposal acc. to local guidelines and/or clinical protocols.

If applicable use integrated spike protection at flow regulator housing for disposal.

Storage and handling conditions

Storage conditions:

Keep away from sunlight.

Keep dry.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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B. Braun Melsungen AG, Carl-Braun-Str. 1, 34212 Melsungen, Germany.
B. Braun Medical AG
B. Braun Melsungen AG