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
EN DISPOSABLE MULTIFUNCTION PLATES INSTRUCTIONS FOR USE

IMPORTANT

The product is intended for use in non-sterile environment by authorized personnel only. Before using the product, the user should deeply understand these instructions.

DESCRIPTION

The disposable multifunction electrodes LESSA DEFI-PAD are constituted by a pair of adhesive pads provided with gels and direct connection to cables and defibrillators which may be used in place of the manual reusable paddles [1].

They are made in both adult and paediatric models, radiopaque or radiolucent and pre-connectable. Some paediatric models to be used with automatic defibrillators (AEDs) provide a system for reducing power.

PACKAGING

Each pair of LESSA disposable multifunction electrodes is packaged in sealed envelopes of opaque material suitable to protect the gel from light and moisture. Envelopes are included in the carton sales package along with a copy of the operating instructions.

INDICATIONS

The disposable multifunction electrode pads LESSA are indicated for:

- Trans thoracic external defibrillation.
- Trans thoracic synchronized cardioversion.
- Trans thoracic ECG Monitoring.
- Temporary trans thoracic cardiac pacing (non-invasive).

LESSA disposable multifunction electrode pads allow the user to effectively operate in the treatment of rhythm disorders related to the above-mentioned applications, without the risk of accidental electrocution related to the use of normally available reusable paddles.

CONTRAINDICATIONS

- LESSA DEFI-PAD disposable multifunction electrodes for adults are generally contraindicated in patients younger than 8 years old [2] (weighing less than 25kg [3][4]), but can be used if the size of the chest allows it, taking care that one pad does not come in contact with the other one. Follow the operating instructions of the defibrillator for energy to be delivered.
- The use of LESSA DEFI-PAD paediatric disposable multifunction electrodes is generally contraindicated in patients older than 8 years old [2] (weighing more than 25kg [3][4]); in fact the power release above 100J may increase the risk of burns.
- Do not apply on skin that shows signs of irritation or injury.

MODE OF USE

External Defibrillation and synchronized cardioversion: the LESSA DEFI-PAD disposable multifunction electrodes are able to transfer to the patient the electrical energy supplied by the defibrillator up to a maximum value of 360 Joules in the adult version [1] and of 100 Joules in the paediatric version [5].

The depolarization of the critical mass of the myocardium, which is essential for the success of the therapy, is only possible if it is crossed by a current of appropriate intensity: the active surface of the electrodes is optimized for this purpose. It is therefore appropriate, in addition to a targeted selection of the positioning sites, to apply the adhesive pads in such a way that their contact surface with the skin is maximum. The choice of power to supply is at the discretion of the operator. In paediatric applications the Guidelines for cardiopulmonary resuscitation recommend a supply of energy of 2-4J / kg; the recommended starting level is of 2 J / kg and it is preferable not to exceed 100J in order to avoid burns [5].

ATTENTION LESSA DEFI-PAD disposable multifunction electrode pads can withstand up to 50 defibrillation shocks; however it is good practice to replace them after 24 hours.

ATTENTION Do not supply a shock with manual metal paddle above the disposable electrode pads or ECG electrode pads.

Non-invasive trans thoracic pacing: LESSA DEFI-PAD disposable multifunction electrode pads can be used for non-invasive trans thoracic pacing. To minimize the threshold of pacing it is appropriate

to apply the adhesive plates in the manner described above. It is also necessary to have a good understanding of the equipment you want to use and follow the manufacturer's instructions.

ATTENTION It is good practice to replace LESSA DEFI-PAD disposable multifunction electrode pads after 8 hours, checking, in case of prolonged pacing (greater than 30 minutes), the skin of the patient for signs of irritation.

ATTENTION Replace LESSA DEFI-PAD disposable multifunction electrode pads after 30 minutes if the supplied pulses are monophasic and longer than 20ms.

ECG monitoring: LESSA DEFI-PAD disposable multifunction electrode pads can also be used for ECG monitoring.

ATTENTION If the tracing is not sufficiently clear, use an ECG patient cable, if it is present, and a separate set of ECG electrodes.

MODE OF APPLICATION

The multifunction electrode pads can be applied to the patient even in the mere suspicion that a severe arrhythmic disease may develop.

The points where it is possible to apply the adhesive electrode pads are listed in "PLACEMENT AND POLARITY".

- Uncover the chest and prepare the skin. Remove excessive hair [1]. Slightly abrade the skin surface to reduce the contact impedance. Avoid applying the adhesive pad on the nipple or breast tissue[6].
- Remove any debris (dirty, greasy and debris), using non-flammable cleaners. Finally, make sure the application sites are clean and dry.
- Open the package and remove multifunction electrode pads
- Gently remove the protective liner, starting from the tab to expose the adhesive and conduction areas.
- In case of multifunctional electrode pads with clips remove the protective support
- Apply the adhesive pads one-to-one starting with one side and pressing progressively over the entire surface to avoid the formation of air bubbles and ensure complete adhesion to the skin. Keep the adhesive pads well separated one from the other and be careful not to overlap them with other objects (ECG electrode pads, cables, transdermal patches, clothing etc.) [7][8].
- Do not replace the adhesive pads once applied. If the position must be changed, remove and replace with new multifunction electrode pads.
- To remove the adhesive plate without irritating the patient's skin, lift an edge and gently pull back. Hold at the same time the skin with the other hand.
- For multifunctional electrode pads without clip: connect the electrodes to the defibrillator or the adapter cable by following the instructions for use of the defibrillator.
- For multifunctional electrode pads with clips: connect the clip to the cable of the defibrillator for the correct polarity, observing the instructions for use of the defibrillator.
- For on demand pacing, separately connect ECG monitoring electrode pads.

POSITIONING AND POLARITY

The international guidelines indicate various placements as equally effective for the treatment of atrial or ventricular arrhythmias[1][5].

The following figures show the application sites commonly used and recommended by most manufacturers of defibrillators. Choose the most appropriate points of application of the therapy according to manufacturer's instructions for use of the defibrillator to be used.

For ease of placement and for training purposes, the anterior-lateral side (Fig.1) is preferred for arrhythmias defibrillation and cardioversion; the anterior-posterior side (Fig.2) is more common in hemodynamics and in trans thoracic pacing.


Fig. 1
• Defibrillation
• Cardioversion
• Pacing
• Monitoring (it provides a track Lead II)

Fig. 2
• Pacing
• Monitoring
• Defibrillation
• Cardioversion

To maintain the proper signal polarity, apply the electrode pads in the indicated positions (the apex is identified by the symbol of the heart). However note that for the purpose of the therapy, it is not relevant which electrode pad (apex / sternum) is placed in one of the two positions [2].

SIDE EFFECTS

- Plate adhesive may cause light cutaneous irritation.
- The prolonged trans thoracic stimulation or the repeated administration of defibrillation shock may cause more or less noticeable skin reddening according to the supplied energy.
- A lack in adhesion and/or air presence under the electrode may cause burnings.

PRECAUTIONS AND WARNINGS

- Use the product only on defibrillators brands, as indicated on the labels.
- Check that the product is compatible with the specific model of defibrillator intended to be used.
- Read the instructions for use of the defibrillator, with particular attention on the placement of multifunction electrode pads, their polarity and the power to be supplied.
- In paediatrics and for some models of automatic defibrillators the use of specific reducing power devices or the adoption of special precautions may be required. Always pay a special attention to energy levels set on the defibrillator and that can be delivered to the paediatric patient (see section "MODE OF USE").

ATTENTION

Do not use paediatric multifunction electrode pads marked with the symbol shown beside with automatic defibrillators. Paediatric multifunction electrode pads marked with the symbol shown beside are indicated for use with automatic defibrillators.

- After an extended period of trans thoracic pacing the ability to detect the evoked ECG signal can be reduced. In this case it is necessary to provide for the collection of the evoked signal by a separate set of electrode pads for ECG.
- Check the expiration date on the package. Do not use after this date.
- Do not use multifunction electrode pads if removed from the envelope for more than 24 hours. The adhesive pads are to be applied within 30 minutes after removal of the protective coating.
- Check that the packaging is intact. Do not use the product otherwise.
- Do not use the multifunction electrode pads if the gel is removed from the support or if it is ripped, torn or dry. Any discoloration localized on gel or on conductive foil does not affect the functionality of the product.
- Do not use the multifunction electrode pads if during removal of the protective coating the product is damaged (eg. The insulating coating of the contact has detached or there are tears in the foam support and / or in the electrode).
- Do not bend, do not cut and do not squash the adhesive pads.

- Do not use the multifunction electrode pads if the connector, the cable or the clips appear to be damaged.
- Check on the operating instructions of the defibrillator at which safety distances the devices (surgeon's electric knife, RF ablaters, diathermy equipment, mobile phones, etc.) that emit strong electromagnetic interferences must be placed.
- To prevent accidental damage from electric shock, ensure that during discharge operators are not in contact with the electrode pads, with the patient, or with conductive parts close to the patient.
- When defibrillators are used near oxygen sources or other flammable gases, use extreme care to avoid risk of fire or explosion.
- The product neither is sterile nor can be sterilized.
- The product is disposable. For use on a single patient. Discard after use.

POTENTIAL COMPLICATIONS

There are no complications related to the use of multifunction electrode pads.

ATTENTION: The defibrillator discharge may cause irregularities in the operation of an implanted pacemaker / defibrillator [1]; apply the multifunction electrode pads at a distance of at least 8cm [9]. After defibrillator discharge check its operation.

ATTENTION: If the chosen energy level is insufficient the success of therapy may be jeopardizing. On the contrary, higher levels may modify the enzyme structure without actual evidence of myocardial damage.

PRODUCT LIFE AND STORAGE

30 months from date of production. Check the expiration date printed on the package. The product should be stored in its original packaging in rooms with environment conditions of temperature and relative humidity, specified on the label. The storage at extreme temperatures must be limited to short periods (24 hours at -30°C or +65°C). Prolonged storage at extreme temperatures can shorten the life of the product.

ATTENTION: The overlap of weights on the packaging could damage the product.

DISPOSAL

Refuses deriving from health structures must be disposed in according to the regulation in force.

WARRANTY AND LIMITATIONS


AB MEDICA GROUP S.A. guarantees that the product complies with Directive 93/42/EEC. No responsibility may be ascribed to the producer who shall not be held liable for medical costs, director indirect damage due to lacking function or malfunction of the above product, when used differently from the instruction for use. We recommend to report opportunely any malfunction or defect of the product to AB MEDICA GROUP S.A. Quality Assurance Service.

CE 0051	European Conformity.	Warning: read the enclosed documentation	Consult instructions for use
REF	Catalogue Number	Serial/ Batch number	Use before
MAX	Manufactured by	Do not reuse	Production date
MIN	Keep away from sun light	Temperature limits	Operating temperature limits
%	Humidity limits	Keep away from humidity	Pes. per box/pack

	Latex free		Compatible with Corpus3 (Defibrillator Unit with reusable pads)		Do not use with Corpus3 with Defib-Unit SLIM
	Non-sterile		1 Pouch/2 Pads		

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