# Rheumatoid Factor Kit (Immunoturbidimetric Method)

## Order Information

| Cat. No.      | Package size              |
|---------------|---------------------------|
| 105-004632-00 | R1: 1×40 mL + R2: 1×11 mL |
| 105-004633-00 | R1: 2×40 mL + R2: 2×11 mL |
| 105-004634-00 | R1: 2×40 mL + R2: 2×11 mL |

## Intended Purpose

In vitro test for the quantitative determination of Rheumatoid Factor (RF) concentration in human serum on Mindray BS series chemistry analyzers. It is intended to be used for aiding to diagnose and monitoring the therapeutic effect of rheumatic diseases.

## Summary<sup>1-3</sup>

Rheumatoid factors are a heterogeneous group of autoantibodies directed against the antigenic determinants on the Fc-region of IgG molecules. They are important in the diagnosis of rheumatoid arthritis, but can also be found in other inflammatory-rheumatic diseases and in various non-rheumatic diseases. They are also found in clinically healthy persons over 60 years of age. Despite these restrictions, the detection of rheumatoid factors is a diagnostic criterion of the American College of Rheumatology for classifying rheumatoid arthritis.

## Assay Principle

Immunoturbidimetric Method

Denatured IgG + RF  $\rightleftharpoons$  Immunocomplex (agglutination)

Determination of the concentration of RF through photometric measurement of immunocomplex between denatured IgG and RF present in the sample, the absorbency increase is directly proportional to the concentration of RF.

## Reagents Components

|            |                  |          |
|------------|------------------|----------|
| <b>R1:</b> | Phosphate buffer | 13.6 g/L |
|            | Sodium chloride  | 11.4 g/L |
|            | Sodium azide     | 0.09%    |
| <b>R2:</b> | Tris buffer      | 6.05 g/L |
|            | Human IgG        | >1 g/L   |
|            | Sodium azide     | 0.09%    |

**Storage and stability**

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

On board in use, the reagents are stable for 28 days when refrigerated on the analyzer.

Contamination must be avoided.

Do not freeze the reagent.

**Specimen collection and preparation****■ Specimen types**

Serum is suitable for samples.

**■ Preparation for Analysis**

1. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
2. Centrifuge samples containing precipitate before performing the assay.
3. Specimens should be tested as soon as possible after sample collection and pre-analytical treatment.

**■ Sample Stability**

3 days at 15-25°C

7 days at 2-8°C

3 months at (-25)-(-15) °C<sup>4</sup> (Only frozen once)

For longer storage periods, samples should be frozen at (-20°C)<sup>5</sup>. Sample stability claims were established by manufacturer and/or based upon references, each laboratory should establish its own sample stability criteria.

**Reagent Preparation**

R1 and R2 are ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

**Materials required but not provided**

1. General laboratory materials: NaCl solution 9 g/L (saline), distilled/deionized water.
2. Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.
3. Mindray BS series chemistry analyzers and General laboratory equipment.

**Assay procedure**

| Parameters Item                                                                                              | BS-800 chemistry analyzers |
|--------------------------------------------------------------------------------------------------------------|----------------------------|
| Assay type                                                                                                   | Endpoint                   |
| Wavelength                                                                                                   | 340 nm                     |
| Reaction direction                                                                                           | Increase                   |
| R1                                                                                                           | 180 $\mu$ L                |
| Sample or Calibrator                                                                                         | 6 $\mu$ L                  |
| Mix, incubate at 37°C for 3-5 min, read the absorbance A1, then add:                                         |                            |
| R2                                                                                                           | 36 $\mu$ L                 |
| Mix thoroughly, incubate at 37°C for 5 min, read the absorbance A2,<br>Then calculate $\Delta A = (A2 - A1)$ |                            |

Parameters may vary in different chemistry analyzers, may adjust in proportion if necessary. For Mindray BS series chemistry analyzers, Reagent Parameters is available on request. Please refer to the appropriate operation manual for the analyzers.

## Calibration

1. It is recommended to use the Mindray Calibrator (Rheumatoid Factor Calibrator: 105-004645-00 or other suitable calibrators) for multi-point calibration. Traceability of the Mindray Rheumatoid Factor Calibrator can refer to the calibrator instructions for use of Mindray Company.

### 2. Calibration frequency

Calibration is stable for approximately 7 days on BS-800 chemistry analyzers. The calibration stability may vary on different instruments, each laboratory should set a calibration frequency in the instrument parameters appropriate to their usage pattern.

Recalibration may be necessary when the following occur:

- As changed reagent lot.
- As required following quality control procedures or out of control.
- As executes specific maintenance or troubleshooting procedure of chemistry analyzers.

3. The calibrator values are lot-specific with the matched models listed in the value sheet.

## Quality control

1. It is recommended to use the Mindray Control (ASO/CRP/RF Triple Control: 105-004650-00 or other suitable controls) to verify the performance of the measurement procedure; other suitable control material can be used in addition.

2. Two levels of control material are recommended to analyze each batch of

samples. In addition, the control should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or trouble shooting procedures as detailed in the appropriate system manual.

- Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control doesn't recover within the acceptable tolerances.

### Calculation

$\Delta A = (A_2 - A_1)$  sample or calibrator.

$C \text{ sample} = (\Delta A \text{ sample} / \Delta A \text{ calibrator}) \times C \text{ calibrator}$ .

The BS series chemistry analyzer detects the change of absorbance ( $\Delta A$ ) and calculates the RF concentration of each sample automatically after calibration.

### Dilution

If the value of sample exceeds 500 IU/mL, the sample should be diluted with 9 g/L NaCl solution (saline) (e.g. 1+19) and the result should be multiplied by 20.

### Expected values

| Sample Type | Units                   |
|-------------|-------------------------|
| Serum       | $\leq 18 \text{ IU/mL}$ |

The expected value above was established by Mindray using 211 samples of people.

Each laboratory should establish its own reference intervals based upon its particular locale and population characteristics since expected values may vary with geography, race, sex and age.

### Performance Characteristics

#### ■ Analytical Sensitivity

The Rheumatoid Factor Kit has an analytical sensitivity of 10 IU/mL on BS-800. Analytical sensitivity is defined as the lowest concentration of analyte that can be distinguished from a sample that contains no analyte. It is calculated as the value lying 3 standard deviations above that of the mean from 20 replicates of an analyte-free sample.

#### ■ Measuring range

The Mindray BS series systems provide the following linearity range:

| Sample Type | Units        |
|-------------|--------------|
| Serum       | 10-500 IU/mL |

A high RF concentration sample (approximately 500 IU/mL) is mixed with a low concentration sample ( $<10 \text{ IU/mL}$ ) at different ratios, generating

a series of dilutions. The RF concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient  $r \geq 0.990$ . The reportable range is 10-2000 IU/mL.

### ■ High-dose hook effect

Hook effect has not occurred up to a RF concentration of 2000 IU/mL on BS-800\*.

*\*Representative data, results in different instruments, laboratories may vary*

### ■ Precision

Precision was determined by following CLSI Approved Guideline EP05-A3<sup>6</sup>, each sample was assayed 2 times per run, 2 runs per day, a total of 20 days. The precision data of controls on BS-800 are summarized below\*.

| Specimen Type<br>(N=80) | Mean<br>(IU/mL) | Repeatability |      | Within-Lab    |      |
|-------------------------|-----------------|---------------|------|---------------|------|
|                         |                 | SD<br>(IU/mL) | CV % | SD<br>(IU/mL) | CV % |
| Control L               | 38.62           | 1.76          | 4.55 | 3.48          | 9.01 |
| Control H               | 83.50           | 2.27          | 2.72 | 3.80          | 4.55 |

*\*Representative data, results in different instruments, laboratories may vary.*

### ■ Analytical Specificity

The samples with different concentration interfering substance were prepared by addition of interferent to human serum pools, and recovers are within  $\pm 10\%$  of the corresponding control value to be considered as no significant interference.

No significant interference was observed when the following substances were tested for interference with this methodology. The data of interference studies on BS-800 are summarized below.

| Interfering Substance | Interferent Concentration<br>(mg/dL) | Analyte Concentration<br>(IU/mL) | Relative Deviation<br>(%)* |
|-----------------------|--------------------------------------|----------------------------------|----------------------------|
| Ascorbic acid         | 30                                   | 87.73                            | -0.99                      |
| Hemoglobin            | 500                                  | 88.76                            | -5.18                      |
| Bilirubin             | 40                                   | 89.03                            | -5.95                      |
| Intralipid            | 500                                  | 87.49                            | -3.99                      |

*\*Representative data, results in different instruments, laboratories may vary.*

In very rare cases gammopathy, in particular type IgM, may cause unreliable results<sup>7</sup>.

**■ Method Comparison**

Correlation studies were performed using CLSI Approved Guideline EP09<sup>8</sup>. The Mindray System (Mindray BS-2000/Mindray RF II Reagent) (y) was compared with comparison system (Hitachi 7180/ Diasys RF Reagent) (x) using the same serum specimens. The statistical data obtained by linear regression are shown in the table below \*:

| Regression Fit     | Correlation Coefficient (r) | Sample (N) | Concentration Range (IU/mL) |
|--------------------|-----------------------------|------------|-----------------------------|
| $y=0.9952x+1.2385$ | 0.9988                      | 100        | 10.11-488.39                |

*\*Representative data, results in different instruments, laboratories may vary.*

**Result interpretation**

The results could be affected by drugs, disease, or endogenous substances<sup>7,9</sup>. When the reaction curve is abnormal, it is recommended to retest and check the result.

**Warnings and precaution**

1. For in vitro diagnostic use only. For laboratory professional use.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The reagents avoid direct exposure to sunlight and freezing. The results can't be assured when stored at inappropriate condition.
4. If unintentionally opened before used, store the reagents tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.
5. Do not mix reagents with different lots and bottles.  
Do not use the reagents beyond the expiration date and the in-use date.  
Do not mix fresh reagents with in-use reagents.  
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
8. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the reagents accidentally enter the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.

10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. The reagent was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen<sup>10</sup>.
13. All human material should be considered potentially infectious.
14. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
15. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## References

1. Allard-Chamard H, Boire G. Serologic Diagnosis of Rheumatoid Arthritis. Clin Lab Med, 2019, 39:525-537.
2. Westwood O, Nelson P, Hay F. Rheumatoid factors: what's new? Rheumatology, 2006;45:379-385.
3. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 810-813.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:41pp.
5. CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007;45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-706.
10. Biosafety in Microbiological and Biomedical Laboratories; U.S.

Department of Health and Human Services; US Government Printing Office; Washington: 2007.

### Graphical symbols

In Vitro Diagnostic  
medical device

Unique device  
identifier

European  
Conformity

Consult Instructions  
For use

Use-by  
date

Authorized representative in  
the European Community

Batch Code

Temperature  
limit

Manufacturer

Catalogue  
number

Biological Risks

Keep away from sunlight

Indicates a medical device that needs protection from light sources

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-81888998; **Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH(Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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1. Allard-Chamard H, Boire G. Serologic Diagnosis of Rheumatoid Arthritis. Clin Lab Med, 2019, 39:525-537.
2. Westwood O, Nelson P, Hay F. Rheumatoid factors: what's new? Rheumatology, 2006;45:379-385.
3. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 810-813.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:41pp.
5. CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007,45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-706.
10. Biosafety in Microbiological and Biomedical Laboratories; U.S.

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4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:41pp.
5. CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007;45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-706.
10. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

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**Fabricante:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Morada:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, R. P. China

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- 4.Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:41pp.
- 5.CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
- 6.CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- 7.Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007, 45(9):1240-1243.
- 8.CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- 9.Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-706.
- 10.Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

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1998; 810-813.

4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:41pp.
5. CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007,45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-706.
10. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

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**Indirizzo:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, Repubblica popolare cinese

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## Referanslar

1. Allard-Chamard H, Boire G. Serologic Diagnosis of Rheumatoid Arthritis. Clin Lab Med, 2019, 39:525-537.
2. Westwood O, Nelson P, Hay F. Rheumatoid factors: what's new? Rheumatology, 2006;45:379-385.
3. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 810-813.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:41pp.
5. CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007,45(9):1240-1243.

8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-706.
10. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

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**Üretici:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Adres:** Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. Çin

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**Références**

1. Allard-Chamard H, Boire G. Serologic Diagnosis of Rheumatoid Arthritis. Clin Lab Med, 2019, 39:525-537.
2. Westwood O, Nelson P, Hay F. Rheumatoid factors: what's new? Rheumatology, 2006;45:379-385.
3. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 810-813.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:41pp.
5. CLSI. Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
6. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
7. Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med, 2007;45(9):1240-1243.
8. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-706.
10. Biosafety in Microbiological and Biomedical Laboratories (Sécurité biologique en laboratoires de microbiologie et laboratoires médicaux) ; U.S. Department of Health and Human Services (Ministère américain de la santé et des Services sociaux) ; US Government Printing Office (Bureau des publications du gouvernement américain) ; Washington : 2007

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© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés
**Fabricant** : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
**Représentant en Europe** : Shanghai International Holding Corp. GmbH (Europe)