
INSTRUCTIONS

LOW FLUX SERIES

[EN] Hollow Fiber Dialyzer/[DE] Hohlfaser Dialysator/[FR] Dialyseur àFibres Creuses/[ES] Dializador de Fibra Hueca
[PT] Dialisador De Fibra Oca/[RU] Капиллярный диализатор/[AR] وجلاف ايدلأابل ايدلأاز ايجافء/[JP] 中空糸型透析器

Specification		F12	F13	F14	F15	F16	F17	F18	F19	F20
Q _B (mL/min)		200/300/400	200/300/400	200/300/400	200/300/400	200/300/400	200/300/400	200/300/400	200/300/400	200/300/400
Clearances Q _D =500mL/min Q _F =0mL/min	Urea	180/230/257	185/238/267	190/245/276	193/251/284	195/256/291	196/260/297	197/263/302	198/266/306	199/269/310
	Creatinine	167/213/231	171/219/239	174/224/246	177/228/252	179/232/258	181/235/263	182/238/267	183/240/270	184/242/273
	Phosphate	147/175/195	152/183/205	156/190/213	159/196/219	162/201/225	165/205/230	167/208/235	168/211/239	169/214/243
	Vitamin B ₁₂	97/102/108	101/107/114	105/112/119	109/116/124	113/120/129	117/124/134	121/128/138	125/132/142	128/136/146
UF coefficient TMP=100mmHg	mL/h·mmHg	16.0	17.0	18.0	19.0	20.0	22.0	23.0	24.0	25.0
Blood Compartment Volume	mL	≤70	≤73	≤75	≤78	≤80	≤90	≤100	≤105	≤110
Pressure Drop Blood Q _B =200mL/min	kPa	≤8	≤8	≤8	≤8	≤8	≤9	≤9	≤10	≤10
Pressure Drop Dialysate Q _D =500mL/min	kPa	≤5	≤5	≤5	≤5	≤5	≤6	≤6	≤6	≤6
Effective Surface Area	m ²	1.2	1.3	1.4	1.5	1.6	1.7	1.8	1.9	2.0
Wall Thickness/ Internal Diameter	μm	40/200								
Effective Fiber Length	mm	232								
Maximum Transmembrane Pressure (Max. TMP)		500mmHg (66.7kPa)								
Sterilization Method		Radiation Sterilization								

SCHEMATIC DIAGRAM OF CONNECTION
OF HOLLOW FIBER DIALYZER

EN Hollow Fiber Dialyzer

INDICATIONS

WEGO hollow fiber dialyzers are designed for single use in haemodialysis therapy for acute and chronic renal failure.

CONTRAINDICATIONS

Generally, there are no absolute contraindications for the application of dialyzer in haemodialysis. Do not use the dialyzer if the patient is allergic to membrane materials.

SIDE EFFECTS

In rare cases, hypersensitivity reactions may occur during haemodialysis treatment. In severe cases, dialysis must be discontinued and the appropriate medication initiated as directed by the physician. In these cases, never transfuse extracorporeal blood back to the patient.

ANTICOAGULATION

It is recommended to introduce an anticoagulant to the extracorporeal circulation. Type, dose and method of application of an anticoagulant must be prescribed by the physician. Qualified hospitals may employ coagulation index to direct the application of anticoagulant.

MATERIALS

Membrane: polysulfone membrane; Housing and end cover: polycarbonate; Potting material: polyurethane; Sealing ring: medical grade silicone; Protective cap: polyethylene.

PRECAUTIONS AND WARNINGS

1. Read the instruction manual carefully before use.
2. This product is designed for single use only. Do not reuse the dialyzer. The manufacturer assumes no liability for any adverse consequences caused by reuse.
3. Inspect the package to ensure that it has not already been opened or damaged before use. Check the expiry date. Do not use dialyzer if package is opened or damaged. Do not use damaged or expired dialyzer.
4. The series of dialyzers are radiation sterilized. Preprocess the dialyzer and the associated materials according to the methods stated in the instruction manual before use.
5. The recommended maximum transmembrane pressure (TMP) does not exceed 500 mmHg (66.7 kPa), the maximum blood flow rate does not exceed 350 mL/min, and the maximum dialysate flow rate does not exceed 800 mL/min.
6. This product is to be used with haemodialysis machine and extracorporeal blood circuit.
7. In the course of regular treatment, place the dialyzer in a vertical position. Placing in a horizontal position may lead to flow deviation and performance degradation.
8. Use the Hansen quick disconnect coupling to connect the dialyzate port. Utilisez le raccord de déconnexion rapide Hansen pour connecter le port de dialyse.
9. Do not use dialyzate supply system without deaeration unit during treatment.
10. The absence of ultrafiltration during dialysis may lead to the reverse ultrafiltration.

INSTRUCTIONS FOR USE

PREPARATION

1. Check the specifications and model. Inspect the package to ensure that it has not been damaged or exposed to moisture. Check the sterilization method and expiry date.
2. Read the operating instruction carefully before use and be aware of the operation requirements for different dialysis membranes and sterilization methods.

PRIMING

1. Open the external packaging of the dialyzer. Place the dialyzer in the dialyzer holder in the vertical position (the blood outlet port facing upwardly). Connect the dialyzer to extracorporeal blood circuit in order. Remain protective caps for dialyzate ports in place. Connect arterial blood line to 1,000 mL physiological saline.
2. Start the blood pump. Pump speed should not exceed 100 mL/min. Evacuate all air from arterial blood line, dialyzer and venous blood line successively. Discharge the priming fluid from the venous blood line to the waste bag. (Intramembrane Priming)
3. Priming volume should not be less than 800 mL. Additional priming volume is appropriately allowed for patients with allergies and hypercoagulability.
4. Connect the arterial and venous blood lines. Connect Hansen quick disconnect couplings to dialyzer to form a closed circulation system. The recommended pump speed is 250-300 mL/min and the minimum cycle is 5 minutes. For patients with allergies, the cycle time may be appropriately extended to 30 minutes or more. (Extramembrane Priming)
5. In closed circulation, set the ultrafiltration volume to 100-200 mL. (Transmembrane Priming)
6. Stop the blood pump. Turn off the fluid infusion line or roller clamp. Reverse the dialyzer to keep the blood inlet port facing upwardly. Set the machine to treatment mode and get ready to dialyze.
7. Note: Do not prime reversely. Ensure that closed circulation should always be started after finishing intramembrane priming with no less than 800 mL physiological saline. It is recommended to inject 5-10 mg heparin to closed circulation via arterial port.

HAEMODIALYSIS INITIATION

1. Inject the heparin via venous port according to the prescription of the physician.
2. Introduction of blood: Connect the arterial blood line to the patient's arterial cannula. Pump speed should not exceed 100 mL/min. Ensure the blood slowly flow into the dialyzer. Stop the blood pump and clamp blood line while the blood flows into the venous chamber. Connect blood line to patient's venous cannula. Ensure that the air is eliminated and the device fixed properly.
3. Adjust the therapy parameters. Gradually increase the Vpump speed if the patient experiences no discomfort.
4. Conduct various monitoring activities according to nursing protocol for haemodialysis and provide treatment accordingly.

HAEMODIALYSIS TERMINATION

1. Stop the blood pump. Clamp the arterial blood line. Disconnect patient's arterial cannula. Connect arterial blood line to physiological saline. Open the clamp. Turn on blood pump and the pump speed should not exceed 100 mL/min.
2. Turn over the dialyzer. Keep the blood outlet port of the dialyzer to face upwardly to facilitate blood discharge from the blood compartment. Gently rub the dialyzer to facilitate blood discharge. Stop the blood pump when the fluid in the venous blood line becomes light pink or almost colorless. Clamp venous blood line and disconnect patient's venous cannula.
3. Provide nursing after dialysis and sterilize the medical devices according to nursing protocol for haemodialysis.

STORAGE CONDITIONS

1. This product should be stored in cool, dry, well-ventilated and clean place without corrosive gas at 80% RH maximum.
2. This product can be stored for three years at the temperature of 0-35°C without any change of quality and performance.

WARRANTY

1. The manufacturer warrants all its dialyzers have been manufactured, sterilized, packaged and inspected in compliance with MDD 93/42/EEC.
2. In case of defect dialyzer (damaged package or dialyzer), return the dialyzer to the manufacturer for free replacement. The manufacturer shall not be liable for any damage incurred subsequent to the manufacturer's delivery of the dialyzer.
3. In case of the use of dialyzer resulting in injury to patient or user, the manufacturer shall not be liable for any injury that cannot be determined to be caused by defect product.
4. The manufacturer shall assume no liability for any injury or adverse effects caused by the reuse of the dialyzer.
5. The manufacturer shall assume no liability for any injury caused by the use of expired dialyzer.
6. The manufacturer shall assume no liability for any injury to the patient or user due to misuse, improper handling, non-compliance with cautionary notes and instructions for use.

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	CE Mark
	Do not re-use
	Caution
	Sterilization using irradiation
	Date of manufacture
	Use-by date
	Batch code
	Do not use if package is damaged
	Consult instructions for use
	Temperature limit
	Authorized representative in the European Community
	Manufacturer

Version 4.0

MedNet GmbH

Borkstrasse 10, 48163 Muenster, GERMANY.

Weihai Weigao Blood Purification Products Co., Ltd.

No.20 Xingshan Road, Weihai Torch Hi-tech Science Park, 264210 Weihai,
Shandong Province, PEOPLE'S REPUBLIC OF CHINA.

Tel: +86-631-5660598 Fax: +86-631-5660598

WEGO
THE BLOOD PURIFICATION EXPERT