# HP IgG

## Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Helicobacter pylori in human serum and plasma

- for "in vitro" diagnostic use only -

**DIA.PRO**  
Diagnostic Bioprobes Srl  
Via G. Carducci n° 27  
20099 Sesto San Giovanni  
(Milano) - Italy  
Phone +39 02 27007161  
Fax +39 02 26007726  
e-mail: [info@dia.pro.it](mailto:info@dia.pro.it)

REF HP.G.CE  
96 Tests

### HP IgG

**A. INTENDED USE**  
Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Helicobacter pylori in human plasma and sera. The product is intended for the follow-up of patients showing gastrointestinal pathologies potentially related to HP-infection.  
For in vitro diagnostic use only.

**B. INTRODUCTION**  
Helicobacter pylori (HP) is a Gram negative bacterium, firstly isolated in gastric mucosa by Marshall and Warren in 1983. HP has been recognized to be the agent responsible of most of cases of gastric mucosal damage and to play a role in the evolution of gastric diseases to carcinoma.  
HP causes an immunological response during infection and specific antibodies of the different classes of IgG, IgM and IgA are produced by the patient.  
ELISAs are currently used to screen patients affected by gastritis or peptic ulcers for acute active infection due to some Helicobacter pylori (HP) strains.  
In particular the presence of IgA and IgM antibodies is reported to be correlated to the acute phase of illness, while IgG antibodies remain present at different times shortly after primary infection and should be used for many years.  
Quantitative ELISA are also useful in monitoring IgG titer undergoing antibiotic therapy, useful in monitoring IgG titer variations during and after the pharmaceutical treatment.

**C. PRINCIPLE OF THE TEST**  
Microplates are coated with H.pylori immunodominant antigens derived from tissue culture of a virulent strain.  
In the 1<sup>st</sup> incubation, the solid phase is treated with diluted samples and anti-HP IgG are captured, if present, by the antigens. Washing out all the other components of the sample in the 2<sup>nd</sup> incubation, anti-HP IgG are detected by the addition of an anti-IgG antibody, labeled with peroxidase (HRP).  
The enzyme captured on the solid phase, acting on the substrate/substrate mixture, generates an optical signal that is proportional to the amount of anti-HP IgG antibodies present in the sample.  
IgG in the sample may therefore be quantified by means of a standard curve calibrated in arbitrary units per milliliter (arbitrary) as no international standard is available.

- D. COMPONENTS**  
Each kit contains sufficient reagents to perform 96 tests.
1. Microplate: **MICROPLATE**
  2. strips: x 9 microtiter wells coated with HP specific immunodominant antigens derived from tissue culture of a virulent strain. Plates are sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.
  3. Calibration Curve: **CAL N°...**  
Ready to use and color coded standard curve ranging:  
4 ml CAL 1 = 0 arbu/ml  
4 ml CAL 2 = 5 arbu/ml  
4 ml CAL 3 = 10 arbu/ml  
2 ml CAL 4 = 20 arbu/ml  
2 ml CAL 5 = 50 arbu/ml  
4 ml CAL 6 = 100 arbu/ml.

Standards are calibrated against an Internal Gold Standard or IGS as no international one is defined.  
Contains human serum proteins, 2% casein, 10 mM Na-dicrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. Standards are blue colored.

3. Control Serum: **CONTROL...m**  
1 val. Lyophilized. It contains fetal bovine serum proteins, human IgG antibodies to HP at about 20 arbu/ml +/-20%, 0.3 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.
  4. Wash buffer concentrate: **WASHBUF 20X**  
1x60ml/bottle. 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0 +/-0.2, 0.05% Tween 20 and 0.05% Kathon GC.
  5. Enzyme conjugate: **CONJ**  
1x16ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8 +/-0.1, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.
  6. Chromogen/substrate: **SUBS TMB**  
1x16ml/vial. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine (or TMB) and 0.02% hydrogen peroxide (or H<sub>2</sub>O<sub>2</sub>).  
*Note: To be stored protected from light as sensitive to strong illumination.*
  7. Sulphuric Acid: **H2SO4 0.3 M**  
1x15ml/vial, contains 0.3 M H<sub>2</sub>SO<sub>4</sub> solution.  
Attention: Irritant (H315, H319, P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).
  8. Specimen Diluent: **DILSPE**  
2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.
  9. Plate sealing foils n°2  
10. Package insert n°1
- E. MATERIALS REQUIRED BUT NOT PROVIDED**
1. Calibrated Micropipettes (1000, 100 and 10ul) and disposable plastic tips.
  2. EtA grade water (distilled or deionised, chlorox treated to remove oxidizing chemicals used as disinfectants).
  3. Timer with 60 minute range or higher.
  4. Absorbent paper/tissues.
  5. Calibrated ELISA microplate thermoslastic incubator (dry or wet) set at +37°C (+/-0.5°C tolerance).
  6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking, strongly recommended) filters.
  7. Calibrated ELISA microplate washer.
  8. Vortex or similar mixing tools.
- F. WARNINGS AND PRECAUTIONS**
1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
  2. All the personnel involved in performing the assay have to wear protective laboratory clothes, lab-coats, gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in bio-safety procedures, as recommended by the Center for

Disease Control, Atlanta, U.S., and reported in the National Institute of Health's publication, *Biosafety in Microbiological and Biomedical Laboratories*, ed. 1984.

- All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
- The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born, microbial agents, when opening kit vials and microplates and when performing the test steps of the Chromagen (TM) from strong light and avoid vortices of the specimen during these test steps undertaken.
- Upon the receipt of the kit at 2,00°C into a temperature controlled refrigerator or cold room.
- Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
- Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, replace the laboratory supervisor to initiate the necessary procedures for kit replacement.
- Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
- Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
- Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 3 months.
- Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication, *Biosafety in Microbiological and Biomedical Laboratories*, ed. 1984.
- The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
- Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
- The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
- Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

**G. SPECIMEN: PREPARATION AND WARNINGS**

1. Blood is drawn aseptically, by venopuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No inhibitors has been observed in the preparation of the sample with citrate, EDTA and heparin.

2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.

3. Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of foam or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

4. Sera and plasma can be stored at +2...8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.

5. If particles are present, centrifuge at 2,000 rpm for 20 min or filter using 0.2-0.5µm filters to clean up the sample for testing.

**H. PREPARATION OF COMPONENTS AND WARNINGS**

**Microplate:**  
Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not in the container. Do not discard a desiccant or conservation. Unused strips have to be supplied, firmly zipped and stored at +2°-8°C. When opened the first time, residual strips are stable till the indicator of humidity inside the desiccant bag turns from yellow to green.

**Calibration Curve**  
Ready to use component. Mix carefully on vortex before use.

**Control Serum**  
Add the volume of ELISA grade water, reported on the label, to the lyophilised powder: fully dissolve and then gently mix on vortex.

**Note:** The control after dissolution is not stable. Store frozen in aliquots at -20°C.

**Wash buffer concentrate:**

The whole content of the concentrated solution has to be diluted 20X with distilled water and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

**Note:** Once diluted, the wash solution is stable for 1 week at +2, 8° C.

**Enzyme conjugate:**

Ready to use. Mix well on vortex before use.  
Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.  
If this component has to be transferred use only plastic, possibly sterile disposable containers.

**Chromogen/Substrate:**

Ready to use. Mix well on vortex before use.  
Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.  
Do not expose to strong illumination, oxidizing agents and metallic surfaces.  
If this component has to be transferred use only plastic, possibly sterile disposable container

**Sample Diluent**

Ready to use component. Mix carefully on vortex before use.

**Sulphuric Acid:**

Ready to use. Mix well on vortex before use.  
Attention: (H315, H319; P280, P302+P552, P332+P313, P305+P351+P338, P337+P313, P362+P363).  
Legend:  
Warning H statements:  
H315 Causes skin irritation.  
H319 Causes serious eye irritation.

*Precautions/ P statements:*  
P280 – Wear protective gloves/protective clothing/eye protection/face protection.  
P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.  
P312 + P313 – If skin irritation occurs: Get medical advice/attention.  
P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.  
P312 + P313 – If eye irritation persists: Get medical advice/attention.  
P302 + P352 – Take off contaminated clothing and wash it before use.

**1. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT**

- Microplates have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +1,2%. Decontamination of spills or residues of kit components should also be carried out regularly.
- The ELISA incubator has to be set at +37°C (tolerance of +/- 0,5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The ELISA washer is extremely important to the overall performance of the assay. The washer must be carefully validated and correctly optimised using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section "O Thermal quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of) needles of the washer has to be carried out according to the instructions of the manufacturer.
- Individual times have a tolerance of ±5%.
- The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (630-650nm, strongly recommended) for the detection of the absorbance range from 0 to 2,00 (0) linearly to 2, 200 absorbance, ± 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in this section "Thermal quality Control". The assay protocol has to be installed in the operating system of the unit, and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation, and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.

- Dia Pro's customer services offers support to the user in the setting and checking of instruments used in combination with the kit. In order to assure compliance with the requirements described, support is also provided for the installation of new instruments to be used with the kit.

**1. PRE-ASSAY CONTROLS AND OPERATIONS**

- Check the expiration date of the kit printed on the external label (but the liquid components are not contaminated by visible particles or aggregates (TM) is colorless or pale blue by separating a small volume of it with a sterile plastic pipette).
- Check that no leakage occurred in transportation and no spillage of liquid is present inside the box (primary microplate is not punctured or damaged).
- Dissolve the content of the Control Serum as reported.
- Dilute all the content of the 20X concentrated Wash Solution as described above.
- Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the microplates are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

**M. ASSAY PROCEDURE**

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.  
The kit may be used for quantitative and qualitative determinations as well.

**M1. QUANTITATIVE DETERMINATION:**

- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex, and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave the A1 and B1 empty for the operation of blanking.
- Dispense 100 µl of Calibrators and 100 µl Control Serum in duplicate. Then dispense 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for 60 min at +37°C.

**Important note:** Strips have to be sealed with the adhesive sealing foil, supplied only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.

- Wash the microplate with an automatic as reported previously (section 1.3).
- Pipette 100 µl Enzyme Conjugate into each well, except A1-B1 blanking wells, and cover with the sealer. Check that this red coloured component has been dispersed in all the wells, except A1 and B1.

**Important note:** Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

- Inclube the microplate for 60 min at +37°C.
- Wash microwells as in step 5.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank wells A1 and B1 included. Then incubate the microplate at room temperature (18-24°C) for 20 minutes.

**Important note:** Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section 1.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1 or B1 or both.

**M2. QUALITATIVE DETERMINATION**

If only a qualitative determination is required, proceed as described below:

- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
- Dispense 100 µl of Calibrator 0, arbutin and Calibrator 5 arbutin in duplicate and Calibrator 100 arbutin in single. Then dispense 100 µl of diluted samples in each properly identified well.
- Inclube the microplate for 60 min at +37°C.

**Important note:** Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.

- Wash the microplate with an automatic washer as reported previously (section 1.3).
- Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispersed in all the wells, except A1.

**Important note:** Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

- Inclube the microplate for 60 min at +37°C.
- Wash microwells as in step 5.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at room temperature (18-24°C) for 20 minutes.

**Important note:** Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section 1.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

**General important notes:**

- If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.

**N. ASSAY SCHEME**

Method	Operations
Calibrators & Control (*)	100 µl
Samples diluted 1:101	100 µl
1 <sup>st</sup> incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme conjugate	100 µl
2 <sup>nd</sup> incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H2O2	100 µl
3 <sup>rd</sup> incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

**(\*) Important Notes:**  
 • The Control Serum (CS) does not affect the test's results calculation.  
 • The Control Serum (CS) used only if a laboratory internal quality control is required by the Management.

An example of dispensation scheme for Quantitative Analysis is reported below:

	Microplate											
	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	CAL4	S1									
B	BLK	CAL4	S2									
C	CAL1	CAL5	S3									
D	CAL1	CAL5	S4									
E	CAL2	CAL6	S5									
F	CAL2	CAL6	S6									
G	CAL3	CS1	S7									
H	CAL3	CS1	S8									

Legend: BLK = Blank  
 CS = Control Serum - Not mandatory S = Sample  
 CAL = Calibrator

An example of dispensation scheme in qualitative assays is reported below:

	Microplate											
	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S3	S11									
B	CAL1	S4	S12									
C	CAL1	S5	S13									
D	CAL2	S6	S14									
E	CAL2	S7	S15									
F	CAL6	S8	S16									
G	S1	S9	S17									
H	S2	S10	S18									

Legend: BLK = Blank  
 S = Sample  
 CAL = Calibrators

**O. INTERNAL QUALITY CONTROL**

A validation check is carried out on the controls and the calibrator any time the kit is used in order to verify whether the performances of the assay are as expected and required by the IVD directive 98/79/EC.

Control that the following data are matched:

Check	Requirements
Blank well	< 0.100 OD450nm value
CAL 1	< 0.150 mean OD450nm value after blanking
CAL 2	coefficient of variation ≤ 30%
CAL 5	OD450nm > OD450nm CAL1 + 0.100
CAL 6	OD450nm > 1.000

If the results of the test match the requirements stated above, proceed to the next section.  
 If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/substrate solution has not got contaminated during the assay
CAL 1 > 0.150 OD450nm	1. that the washing procedure and the washer settings are as validated in the pre qualification study/ 2. that the proper washing solution has been used before use/ 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator or of their wells has occurred due to a contamination of the negative one: 4. that no contamination of the enzyme wells of positive samples or the enzyme samples that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 5. that the washer needles are not blocked or partially obstructed.
coefficient of variation > 30%	1. that the Chromogen/substrate solution has not got contaminated during the assay 2. that the proper washing solution has been used before use/ 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator or of their wells has occurred due to a contamination of the negative one: 4. that no contamination of the enzyme wells of positive samples or the enzyme samples that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 5. that the washer needles are not blocked or partially obstructed.

CAL 2 5 arbutin OD450nm < 0.000 OD450nm CAL1 + 0.000	1. that the procedure has been correctly executed/ 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred
CAL 6 100 arbutin < 1.000 OD450nm	1. that the procedure has been correctly executed/ 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred

Should one of these problems have happened, after checking, report to the supervisor for further actions.

**\*\* Note:**  
 If Control Serum has used, verify the following data:

Check	Requirements
Control Serum	Mean OD450nm CAL4 +20%

If the results of the test doesn't match the requirements stated above, operate as follows:


Problem	Check
Different from expected value	1. that the procedure has been correctly executed/ 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

Anyway, if all other parameters (Blank, CAL1, CAL2, CAL 6), match the established requirements, the test may be considered valid.

**P. RESULTS**

**P.1 Quantitative method**

If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm. (4-Parameters interpolation is suggested).  
 Then on the calibration curve calculate the concentration of anti Hyvion IgG antibody in samples.  
 An example of Calibration curve is reported below.

**Important Note:**  
Do not use the calibration curve above to make calculations.  
P.2 Quantitative method  
In the quantitative method, calculate the mean OD450nm values for the Calibrators 0 and 5 arbu/ml and then check that the assay is valid.

Example of calculation:

**Note:** The following data must not be used instead of real figures obtained by the user:

- Calibrator 0 arbu/ml: 0,020 – 0,024 OD450nm
- Mean Value: 0,022 OD450nm
- Lower than 0,150 – Accepted
- Calibrator 5 arbu/ml: 0,250 – 0,270 OD450nm
- Mean Value: 0,260 OD450nm
- Higher than Cal 0 + 0,100 – Accepted
- Calibrator 100 arbu/ml: 2,045 OD450nm
- Higher than 1,000 – Accepted

The OD450nm of the Calibrator 5 arbu/ml is considered the cut-off (or Co) of the system.  
The ratio between the OD450nm value of the sample and the OD450nm of the Calibrator 5 arbu/ml (or S/Co) can provide a semi-quantitative estimation of the content of specific IgG in the sample.

**Q. INTERPRETATION OF RESULTS**

Samples with a concentration lower than 5 arbu/ml are considered negative for anti H pylori IgG antibody.  
Samples with a concentration higher than 5 arbu/ml are considered positive for anti H pylori IgG antibody.

**Important notes:**

1. H pylori IgG results alone are not enough to provide a clear diagnosis of Helicobacter pylori infection. Other tests for Helicobacter pylori (supplied by DiaPro Diagnostic Bioprobes s.r.l. at code n° HPA.G.CE, HPA.CE and HPA.M.CE), should be carried out.

2. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
3. When test results are transmitted from the laboratory to the referring facility, attention must be paid to avoid erroneous data transfer; diagnosis must be done and released to the patient by a suitably qualified medical doctor.

**R. PERFORMANCE CHARACTERISTICS**

Evaluation of Performance has been conducted on panels of positive and negative samples in an external clinical laboratory with reference to a FDA approved reference kit.

**1. Limit of detection**

No international standard for HP IgG Antibody detection has been defined so far by the European Community.  
In its absence, an Internal Gold Standard (or GS) derived from a patient with an history of past, non-mucosaless, infection, has been defined in order to provide the device with a constant and excellent sensitivity.

**2. Diagnostic Sensitivity and Specificity:**

The diagnostic performances were evaluated in a performance evaluation study, conducted in an external center, with excellent experience in the diagnosis of infectious diseases.  
The diagnostic sensitivity was studied on more than 50 samples, pre-tested positive with the reference kit of European origin in use at the laboratory. Positive samples were collected from patients with a clinical history of H.pylori infection.  
The diagnostic specificity was determined on panels of more than 100 negative samples from normal individuals and blood donors, classified negative with the reference kit, including potentially interfering specimens.  
Both plasma, derived with different standard techniques of preparation (citrate, EDTA, and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.  
Frozen specimens have also been tested to check whether samples freezing interferes with the performance of the test. No interference was observed on clean and particle free samples.  
The Performance Evaluation provided the following values :

Sensitivity	> 98 %
Specificity	> 98 %

**3. Reproducibility:**

A study conducted on three samples of different HP IgG reactivity, examined in 16 replicates in three separate runs has shown CV% values ranging 2-18% depending on the OD450nm readings.  
The variability shown in the tables did not result in sample misclassification.

**S. LIMITATIONS**

False positivity has been assessed as less than 2% of the normal population.  
Frozen samples containing film particles or aggregates may generate false positive results.

**REFERENCES**

1. Lazzaroni M, et al., Medicina (1989), 9, 9-18.
2. Vaina D, et al., Federazione Medica ALI (1988), 7, 548-555.
3. Oberaid G, Et al., The Lancet (1989), vol.5, 7, 358-360.
4. Corleaid H, et al., The Lancet (1989) vol.5, 10, 354-355
5. Vaina D, et al., British Medical Journal (1989), vol.9, 43, 374-375.
6. Oberaid G, et al., Gut (1989), vol. 30, 7, 912-916.
7. Vaina D, et al., Ital.J.Gastroenterol. (1988), 20, 299-304.
8. Vaina D, et al., Current Opinion in Gastroenterology (1989), 5, 617-623.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Produced by  
Dia.Pro. Diagnostic Bioprobes Srl.  
Via G. Carducci, n° 27 – Sesto San Giovanni (MI) – Italy

**DECLARATION OF CONFORMITY**

1) **Manufacturer** (Name, department): Monobind Inc.

Address: 100 North Pointe, LAKE FOREST, CA 92630, UNITED STATES and

2) **European authorized representative**: CEPATNER4U BV,

Address: ESDOORNLAAN 13, 3951DB MAARN, THE NETHERLANDS;  
 (on product labels printed as:  
 CEPATNER4U, ESDOORNLAAN 13, 3951DB MAARN, THE NETHERLANDS, www.cepatner4u.eu)

3) **Product(s)** (name, type or model/batch number, etc.):

Immunoassay products;  
 ELISA,  
 CLIA,  
 Control,  
 Instruments (see appendix)

4) **The product(s) described above is in conformity with:**

Title	Document No.
In vitro Diagnostic Medical Devices Directive	98/79/EC

5) **Additional information** (Conformity procedure, Notified Body, CE certificate, Registration nr., etc.):

Conformity assessment procedure for CE marking: In vitro Diagnostic Medical Device Directive, Annex III

Registration nr.: NL-CA002-22758 and NL-CA002-22762

Lake Forest, USA, 2013-09-16

*A Shatola*

(Place & date of issue (YYYY-MM-DD))  
 Tony Shatola, QA Director, Monobind Inc.

Maarn, NL, 2013-09-16

*[Signature]*  
 Olga Teirlnck, Consultant, CEPATNER4U BV

(Place & date of issue (YYYY-MM-DD))  
 (name, function and signature of authorized representative)

**Appendix**

Date: 2013-09-16

**List of devices:**

Device types	Item# Accelife® ELISA, Microvelis	Item# CLIA Microvelis	Item# Control	Item# erit	EDMS code	Risk Class	First date of CE-marking
<b>Thyroid</b>							
Total Triiodothyronine (T3) Test System	125-300	175-300			12.04.01.05.00	Low	2005-11-11
Free Triiodothyronine (T3) Test System	1325-300	1375-300			12.04.01.01.00	Low	2005-11-11
Thyroxine (T4) Test System	225-300	275-300			12.04.01.07.00	Low	2005-11-11
Free Thyroxine (T4) Test System	1225-300	1275-300			12.04.01.02.00	Low	2005-11-11
Thyrotropin (TSH) Test System	325-300	375-300			12.04.01.11.00	Low	2005-11-11
Rapid TSH Test System	6025-300	6075-300			12.04.01.11.00	Low	2010-06-29
T3-Uptake (T3U) Test System	525-300	575-300			12.04.01.08.00	Low	2005-11-11
Thyroxine-Binding Globulin (TBG) Test System	3525-300	3575-300			12.04.01.09.00	Low	2005-11-11
Thyroglobulin (Tg) Test System	2225-300	2275-300			12.04.01.08.00	Low	2005-11-11
Total Thyroxine (T4), Total Triiodothyronine (T3) & Thyroid Stimulating Hormone (TSH) Thyroid Panel (VAST) Test System	8025-300	8075-300			12.04.01.01.00	Low	2005-11-11
Total Triiodothyronine (T3 SRS) Test System	8125-300	8175-300			12.04.01.01.00	Low	2010-06-29
Total Thyroxine (T4 SRS) Test System	8225-300	8275-300			12.04.01.01.00	Low	2010-06-29
Free Thyroxine (fT4), Free Triiodothyronine (fT3) & Thyroid Stimulating Hormone Free Thyroid Panel (VAST) Test System	7025-300	7075-300			12.04.01.01.00	Low	2010-06-29
<b>Neonatal Thyroid &amp; Genetics</b>							
Neonatal TSH (N-TSH) Test System	3425-300	3475-300			12.04.01.90.00	Low	2005-11-11
Neonatal (N-T4) Thyroxine Test System	2625-300	2675-300			12.04.01.12.00	Low	2005-11-11
Neonatal 17OHP (N-17OHP) Test System	5525-300	5575-300			12.05.01.07	Low	2008-02-01
Neonatal TBG (N-TBG) Test System	8925-300	8975-300			12.04.01.09.00	Low	2013-09-16
<b>Auto Immune Thyroid</b>							
Anti-Thyroglobulin (Anti-Tg) Test System	1025-300	1075-300			12.10.03.04.00	Low	2005-11-11
Anti-Thyropoxidase (Anti-TPO) Test System	1125-300	1175-300			12.10.03.01.00	Low	2005-11-11
<b>Fertility &amp; Prenatal</b>							
Luteinizing Hormone (LH) Test System	625-300	675-300			12.05.01.05.00	Low	2005-11-11
Folicle Stimulating Hormone (FSH) Test System	425-300	475-300			12.05.01.04.00	Low	2005-11-11
Prolactin Hormone (PRL) Test System	725-300	775-300			12.05.01.08.00	Low	2005-11-11
Prolactin Hormone Sequential (PRLs) Test System	6025-300	6075-300			12.05.01.08.00	Low	2005-11-11
B-Human Chorionic Gonadotropin (hCG) Test System	825-300	875-300			12.05.02.05.00	Low	2005-11-11
B-Human Chorionic Gonadotropin Extended Range (Ext. Range hCG) Test System	8825-300	8875-300			12.05.02.05.00	Low	2013-09-16
Rapid B-Human Chorionic Gonadotropin (Rapid	3325-300				12.05.02.05.00	Low	2005-11-11

Item#	Item#	Item#	Item#	Item#	Item#	Risk	First date of
Acculite®	Acculite®	QSurge®	System	EDMS code	Class	CE-Marking	
ELISA	CLIA	Control	ent				
Microvels	Microvels						
Human Chorionic Gonadotropin (hCG), Human Progactin (hPRL), Human Luteinizing Hormone (LH), Follicle Stimulating Hormone (FSH) Fertility Panel (VAST) Test System	8325-300	8375-300		12.05.01.90.00	Low	2006-08-24	
Alpha-Fetoprotein (AFP), Human Chorionic Gonadotropin (hCG), Unconjugated Estriol (u-E3) Triple Screen (VAST) Test System	8525-300	8575-300		12.05.01.90.00	Low	2010-06-29	
Pregnancy Associated Plasma Protein - A (PAPP-A) Test System	7925-300	7975-300		12.05.02.10.00	Low	2013-09-16	
Steroid							
Cortisol Test System	3625-300	3675-300		12.06.02.04.00	Low	2005-11-11	
DHEAS Test System	5125-300	5175-300		12.05.01.02.00	Low	2010-06-29	
Dehydroepiandrosterone (DHEA) Test System	7425-300	7475-300		12.05.01.02.00	Low	2011-09-26	
Estradiol (E2) Test System	4925-300	4975-300		12.05.01.03.00	Low	2010-06-29	
Unconjugated Estriol (u-E3) Test System	5025-300	5075-300		12.05.02.02.00	Low	2010-06-29	
Progesterone Test System	4825-300	4875-300		12.05.01.06.00	Low	2010-06-29	
Sex Hormone Binding Globulin (SHBG) Test System	9125-300	9175-300		12.05.01.09.00	Low	2013-09-16	
Testosterone Test System	3725-300	3775-300		12.05.01.10.00	Low	2007-11-01	
Free Testosterone Test System	5325-300	5375-300		12.05.01.10.00	Low	2010-06-29	
17α-OH Progesterone Test System	6225-300	6275-300		12.05.01.07.00	Low	2010-06-29	
17α-OH Progesterone - SI Test System	9925-300	9975-300		12.05.01.07.00	Low	2010-10-18	
Growth & Bone Metabolism							
Growth Hormone (hGH) Test System	1725-300	1775-300		12.06.04.02.00	Low	2005-11-11	
Parathyroid Hormone (PTH) Test System	9225-300	9275-300		12.06.03.13.00	Low	2011-09-26	
25-Hydroxyvitamin D3 (Vitamin D3) Test System	7725-300	7775-300		12.06.03.10.00	Low	2011-09-26	
Diabetes							
Insulin Test System	2425-300	2475-300		12.06.01.03.00	Low	2005-11-11	
Rapid Insulin Test System	5625-300	5675-300		12.06.01.03.00	Low	2010-06-29	
C-Peptide Test System	2725-300	2775-300		12.06.01.01.00	Low	2005-11-11	
Insulin - C-Peptide (VAST)	7325-300	7375-300		12.06.01.03.00	Low	2005-11-11	
Cardiac Markers							
CK-MB Test System	2925-300	2975-300		12.13.01.02.00	Low	2005-11-11	
Troponin I (cTnI) Test System	3825-300	3875-300		12.13.01.07.00	Low	2005-11-11	
Digoxin (DIG) Test System	925-300	975-300		12.06.01.01.00	Low	2005-11-11	
High Sensitivity CRP (hs-CRP) Test System	3125-300	3175-300		12.13.01.90.00	Low	2005-11-11	
Myoglobin Test System	3225-300	3275-300		12.13.01.05.00	Low	2005-11-11	

Item#	Item#	Item#	Item#	Item#	Item#	Risk	First date of
Acculite®	Acculite®	QSurge®	System	EDMS code	Class	CE-Marking	
ELISA	CLIA	Control	ent				
Microvels	Microvels						
Infectious Diseases							
Anti-H. Pylori Igs Test System	1425-300	1475-300		15.01.04.03.00	Low	2005-11-11	
Anti-H. Pylori IGM Test System	1525-300	1575-300		15.01.04.03.00	Low	2005-11-11	
Anti-H. Pylori Iga Test System	1625-300	1675-300		15.01.04.03.00	Low	2005-11-11	
Cancer Markers							
Alpha-Fetoprotein (AFP) Test System	1925-300	1975-300		12.03.90.01.00	Low	2005-11-11	
CA-125 Test System	3025-300	3075-300		12.03.01.06.00	Low	2005-11-11	
CA 15-3 Test System	5625-300	5675-300		12.03.01.02.00	Low	2010-06-29	
CA -19-9 Test System	3925-300	3975-300		12.03.01.03.00	Low	2005-11-11	
Carcinoembryonic Antigen (CEA) Test System	1825-300	1875-300		12.03.01.31.00	Low	2005-11-11	
Next Generation Carcinoembryonic Antigen (CEA) Test System	4625-300	4675-300		12.03.01.31.00	Low	2010-06-29	
Free β-Subunit Human Chorionic Gonadotropin (hHCG) Test System	2025-300	2075-300		12.03.01.90.00	Low	2005-11-11	
Allergy & Anemia							
Ferritin Test System	2825-300	2875-300		12.07.01.02.00	Low	2005-11-11	
Folate Test System	7525-300	7575-300		12.07.01.03.00	Low	2010-06-29	
Immunoglobulin E (IgE) Test System	2525-300	2575-300		12.02.01.02.00	Low	2005-11-11	
Transferrin Soluble Receptor (sTFR) Test System	8625-300	8675-300		12.07.01.06.00	Low	2010-06-29	
Vitamin B12 (B12) Test System	7625-300	7675-300		12.07.02.04.00	Low	2011-09-26	
Folate, Vitamin B-12 (VAST) Test System	7825-300	7875-300		12.07.01.00.00	Low	2013-09-16	
Miscellaneous Controls							
Anti-Thyroglobulin (Anti-Tg), Anti-Thyroglobulin (Anti-TPO) Control - Positive & Negative				AIT-101	Low	2010-06-29	
High Level Fertility Control - Single Level - Progesterone, Estradiol, Human Chorionic Gonadotropin				FC-300	Low	2010-06-29	
Material Control - Tri Level - Human Chorionic Gonadotropin, Free Beta Human Chorionic Gonadotropin Subunit, Alpha Fetal Protein, Estriol				MC-300	Low	2010-06-29	
Thyroglobulin (Tg) Control - Tri Level				TG-300	Low	2010-06-29	
H. Pylori Igs Control - Positive & Negative				HPY-IgG-300	Low	2010-06-29	
H. Pylori IGM Control - Positive & Negative				HPY-IgM-300	Low	2013-09-16	
H. Pylori Iga Control - Positive & Negative				HPY-IgA-300	Low	2013-09-16	
Thyroid Binding Globulin (TBG) Control - Tri Level				TBSG-300	Low	2013-09-16	
Miscellaneous Instruments							
Autopex ELISA & CLIA Analyzer				IND06	Low	2010-06-29	
Autopex Generation 2 ELISA & CLIA Analyzer				IND06-2	Low	2013-09-16	
Lumex CLIA Analyzer				IND01	Low	2006-08-24	
Nee-Lumax CLIA Analyzer				IND10	Low	2011-09-26	

Device types	Item# Account# ELISA Microwells	Item# Account# CLIA Microwells	Item# OSure# Control	Item# Instrum Kit	EDMS code	Risk Class	First date of CE-Marking
Impulse 2 CLIA Analyzer				IN005	21.02.10.01	Low	2006-08-24
Impulse 3 CLIA Analyzer				IN007	21.02.10.01	Low	2010-06-29
Lumax96 CLIA Analyzer				IN004	21.02.10.01	Low	2007-03-01
Lumatic CLIA Analyzer				IN008	21.02.10.01	Low	2011-09-26
Eldex 3 8 ELISA Analyzer				IN003	21.02.10.01	Low	2007-09-10
Neo-Eldex ELISA Analyzer				IN009	21.02.10.01	Low	2011-09-26
PriSMatic ELISA Analyzer				IN013	21.02.10.01	Low	2013-09-16
Plate Washer Microplate Washer				IN002	21.02.10.01	Low	2010-06-29

**Monobind Inc.**  
Lake Forest, CA 92630, USA

**AccuBind**  
ELISA Microwells

## Total Prostate Specific Antigen (tPSA) Test System

**Product Code: 2125-300**

### 1.0 INTRODUCTION

**Intended Use:** The Quantitative Determination of Total Prostate Specific Antigen (tPSA) Concentration in Human Serum by a Microplate Immunoassay.

### 2.0 SUMMARY AND EXPLANATION OF THE TEST

Prostate Specific Antigen (PSA) is a serine protease with chymotrypsin-like activity (1,2). The protein is a single chain glycoprotein with a molecular weight of 28.4 kDa (3). PSA derives its name from the observation that it is a normal antigen of the prostate but is not found in any other normal or malignant tissue.

PSA is found in benign, malignant, and metastatic prostate cancer. Since prostate cancer is the second most prevalent form of male malignancy, the detection of elevated PSA levels plays an important role in the early diagnosis. Serum PSA levels have been found to be more useful than prostatic acid phosphatase (PAP) in the diagnosis and management of patients due to increased sensitivity (4).

In this method, PSA calibrator, patient specimen or control is first added to a streptavidin coated well. Biotinylated monoclonal and enzyme labeled antibodies (directed against distinct and different epitopes of PSA) are added and the reactants mixed. Reaction between the various PSA antibodies and native PSA forms a sandwich complex that binds with the streptavidin coated to the well.

After the completion of the required incubation period, the enzyme-PSA antibody bound conjugate is separated from the unbound enzyme-PSA conjugate by aspiration or decantation. The activity of the enzyme present on the surface of the well is quantified by reaction with a suitable substrate to produce color.


The employment of several serum references of known prostate specific antigen (PSA) levels permits the construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with PSA concentration.

### 3.0 PRINCIPLE

#### Immunoenzymometric assay (TYPE 3):

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitopes recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on

the well and exogenously added biotinylated monoclonal anti-PSA antibody.  
Upon mixing monoclonal biotinylated antibody, the enzyme-labeled antibody and a serum containing the native antigen, reaction results between the native antigen and the antibodies, without competition or steric hindrance, to form a soluble sandwich complex. The interaction is illustrated by the following equation:

$Ag_{(s)}$  = Native Antigen (Variable Quantity)

$Enz^{Ab_{(b)}}$  = Enzyme labeled Antibody (Excess Quantity)

$Enz^{Ab_{(b)}}-Ag_{(s)}-En^{Ab_{(m)}}$  = Antigen-Antibodies Complex

$k_a$  = Rate Constant of Association

$k_d$  = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

Immobilized complex = complex bound to the solid surface  
After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antigen-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

### 4.0 REAGENTS

#### Materials Provided:

**A. Prostate Specific antigen (PSA) 1ml/vial - Icons A-F**  
Six (6) vials of references PSA Antigen at levels of 0.0A, 5.0B, 10.0C, 25.0D, 50.0E and 100.0F ng/ml. Store at 2-8°C. A preservative has been added.

**Note:** The calibrators, human serum based, were calibrated using a reference preparation, which was assayed against the 1<sup>st</sup> IS 96/670.

#### B. PSA Enzyme Reagent - 13ml/vial - Icon

One (1) vial containing enzyme labeled antibody, biotinylated monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

#### C. Streptavidin Coated Plate - 96 wells - Icon

One 96-well microplate coated with streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

#### D. Wash Solution Concentrate - 20 ml - Icon

One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-30°C.

#### E. Substrate A - 7ml/vial - Icon S<sup>A</sup>

One (1) bottle containing tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

#### F. Substrate B - 7ml/vial - Icon S<sup>B</sup>

One (1) bottle containing hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) in buffer. Store at 2-8°C.

#### G. Stop Solution - 8ml/vial - Icon

One (1) bottle containing a strong acid (1N HCl). Store at 2-30°C.

#### H. Product Instructions.

**Note 1:** Do not use reagents beyond the kit expiration date.

**Note 2:** Avoid extended exposure to heat and light. Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.

**Note 3:** Above reagents are for a single 96-well microplate.

#### 4.1 Required But Not Provided:

1. Pipette(s) capable of delivering 25 & 50µl volumes with a precision of better than 1.5%.
2. Dispenser(s) for repetitive deliveries of 0.100ml and 0.350ml volumes with a precision of better than 1.5%.
3. Microplate washers or a squeeze bottle (optional).
4. Microplate Reader with 450nm and 620nm wavelength absorbance capability.
5. Absorbent Paper for blotting the microplate wells.
6. Plastic wrap or microplate cover for incubation steps.
7. Vacuum aspirator (optional) for wash steps.
8. Timer.
9. Quality control materials

### 5.0 PRECAUTIONS

#### For In Vitro Diagnostic Use

#### Not for Internal or External Use in Humans or Animals

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA licensed reagents. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

Safe disposal of kit components must be according to local regulatory and statutory requirements.

### 6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood, serum in type and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain red top venipuncture tube without additives or anti-coagulants. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells. Samples may be refrigerated at 2-8°C for a maximum period of five (5) days, if the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.050ml of the specimen is required.

### 7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and elevated range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh serum reference should be used to determine the reason for the variations.

### 8.0 REAGENT PREPARATION:

#### 1. Wash Buffer

Dilute contents of wash concentrate to 1000ml with distilled or deionized water in a suitable storage container. Store at room temperature 20-27°C for up to 60 days.

#### 2. Working Substrate Solution

Pour the contents of the amber-wal-labeled Solution-A into the clear vial labeled Solution-B. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2-8°C.  
**Alert:** Do not use the working substrate if it looks blue.

**Note 2:** Do not use reagents that are contaminated or have bacteria growth.

### 9.0 TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum references and controls to room temperature (20 - 27°C).

**\*\*Test Procedure should be performed by a skilled individual or trained professional\*\***

1. Format the microplates' wells for each serum reference, control and patient specimen to be assayed in duplicate. Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.
2. Pipette 0.025 ml (25µl) of the appropriate serum reference, control or specimen into the assigned well.
3. Add 0.100 ml (100µl) of the PSA Enzyme Reagent to each well. It is very important to dispense all reagents close to the bottom of the coated well.
4. Swirl the microplate gently for 20-30 seconds to mix and cover.
5. Incubate 30 minutes at room temperature.
6. Decant the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
7. Add 350µl of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.
8. Add 0.100 ml (100µl) of working substrate solution to all wells (see Reagent Preparation Section). Always add reagents in the same order to minimize reaction time differences between wells.

#### DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION

9. Incubate at room temperature for fifteen (15) minutes.
10. Add 0.050ml (50µl) of stop solution to each well and mix gently for 15-20 seconds. Always add reagents in the same order to minimize reaction time differences between wells.
11. Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. The results should be read within thirty (30) minutes of adding the stop solution.

### 10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of PSA in unknown specimens.

1. Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
2. Plot the absorbance for each duplicate serum reference versus the corresponding PSA concentration in ng/ml on linear graph paper (do not average the duplicates or the serum references before plotting).
3. Draw the best-fit curve through the plotted points.
4. To determine the concentration of PSA for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (1.142) intersects the dose response curve at (23.6 ng/ml) PSA concentration (See Figure 1).


Note: Computer data reduction software designed for ELISA assays may also be used for the data reduction, if such software is utilized, the validation of the software should be ascertained.

**EXAMPLE 1**

Sample ID	Well	Abs (A)	Abs (B)	Mean	Value (ng/ml)
Cal A	A1	0.019	0.019	0	
	B1	0.019			
Cal B	C1	0.279	0.276	5	
	D1	0.273			
Cal C	E1	0.567	0.563	10	
	F1	0.559			
Cal D	G1	1.248	1.213	25	
	H1	1.179			
Cal E	A2	2.051	1.999	50	
	B2	1.947			
Cal F	C2	2.892	2.833	100	
	D2	2.775			
Patient	E3	1.186	1.142	23.6	
	F3	1.099			

\*The data presented in Example 1 and Figure 1 is for illustration only and should not be used in lieu of a dose response curve prepared with each assay.

**Figure 1**

**11.0 Q.C. PARAMETERS**

In order for the assay results to be considered valid the following criteria should be met:

1. The absorbance (OD) of calibrator F should be  $\geq 1.3$ .
2. Four out of six quality control pools should be within the established ranges.

**12.0 RISK ANALYSIS**

The MSDS and Risk Analysis Form for this product is available on request from Monobind Inc.

**12.1 Assay Performance**

1. It is important that the time of reaction in each well is held constant to achieve reproducible results.
2. Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.

3. Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
4. If more than one (1) plate is used, it is recommended to repeat the dose response curve.
5. The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.

**12.2 Interpretation**

6. Plate readers measure vertically. Do not touch the bottom of the wells.
7. Failure to remove adhering solution adequately in the application or decantation wash step(s) may result in poor replication and spurious results.
8. Use components from the same lot. No intermixing of reagents from different batches.
9. Patient specimens with PSA concentrations above 100 ng/ml may be diluted (for example 1/10 or higher) with normal female serum (PSA = 0 ng/ml) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor (10).
10. Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind's IFU may yield inaccurate results.
11. All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
12. It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
13. Risk Analysis- as required by CE Mark IVD Directive ISO 14971:2009 - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

**13.0 EXPECTED RANGES OF VALUES**

Healthy males are expected to have values below 4 ng/ml (4).

**TABLE 1**

Expected Values for the PSA Elisa Test System	Healthy Males
	<4 ng/ml

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal" persons is dependent upon a multiplicity of factors- the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by

the analysis using the method with a population indigenous to the area in which the laboratory is located.

**14.0 PERFORMANCE CHARACTERISTICS**

**14.1 Precision**

The within and between assay precisions of the PSA Accubind™ ELISA test system were determined by analyses on three different levels of control sera. The number, mean value, standard deviation and coefficient of variation for each of these control sera are presented in Table 2 and Table 3.

**TABLE 2**  
Within Assay Precision (Values in ng/ml)

Sample	N	X	S.D.	C.V.
Level 1	20	0.7	0.05	7.1%
Level 2	20	4.5	0.20	4.4%
Level 3	20	28.3	1.07	3.7%

**TABLE 3**  
Between Assay Precision\* (Values in ng/ml)

Sample	N	X	S.D.	C.V.
Level 1	10	0.8	0.09	11.3%
Level 2	10	4.3	0.25	5.8%
Level 3	10	27.5	1.42	5.2%

**14.2 Sensitivity**

The PSA Accubind™ ELISA test system has a sensitivity of 0.012 ng. This is equivalent to a sample containing 0.5 ng/ml PSA concentration.

**14.3 Accuracy**

The PSA Accubind™ ELISA method was compared with a reference Elisa method. Biological specimens from low, normal, and elevated concentrations were assayed. The total number of such specimens was 241. The least square regression equation and the correlation coefficient were computed for the PSA Accubind™ ELISA test method in comparison with the reference method. The data obtained is displayed in Table 4.

**TABLE 4**

Method	Mean	Least Square Regression Analysis	Correlation Coefficient
This Method (X)	5.82	$y = -0.0598x + 0.98(x)$	0.987
Reference (Y)	5.57		

Only slight amounts of bias between the PSA Accubind™ ELISA method and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

**14.4 Specificity:**

No interference was detected with the performance of PSA Accubind™ ELISA test system upon addition of massive amounts of the following substances to a human serum pool.

Acetylsalicylic Acid	100 µg/ml
Ascorbic Acid	100 µg/ml
Caffeine	100 µg/ml
CEA	10 µg/ml
AFP	10 µg/ml
CA-125	10,000 U/ml
hCG	10,000 IU/ml
hLH	100 IU/ml
hTSH	100 mIU/ml
hPRL	100 µg/ml

**15.0 REFERENCES**

1. Christenson A, Laurell CB, Lilja H. *Eur J Biochem.* 194, 755-63 (1990).
2. Wart KW, et al. *Proc Natl Acad Sci USA.* 83, 3166-70 (1986).
3. Chen Z, Prestigiacomo A, Stamey T, Chin Chem, 41, 1273-82 (1995).
4. Wild D. *The Immunoassay Handbook.* Stockton Press, 452, (1994).

5. Junker R, Brandt B, Zechel C, Assmann G. *Chin Chem.* 43, 1588-94 (1997).
6. Prestigiacomo A, Stamey T. "Physiological variations of serum prostate antigen in the (4-10 ng/ml) range in male volunteers." *J Urol.* 155, 1977-80 (1996).
7. Stamey T, McNeal JE, Yamoto CM, Sigal BM, Johnstone IM, "Biological determinants of cancer progression in men with prostate cancer." *JAMA* 281, 1395-1400 (1999).
8. Chen Z, Prestigiacomo A, Stamey T. "Purification and characterization of Prostate Specific Antigen (PSA) Complexed to  $\alpha$ -Antihymotrypsin: Potential reference Material for International Standardization of PSA Immunoassays." *Chin Chem.* 41(9), 1273-1282 (1995).
9. Horton GL, Balmson RR, Dart WJ, Chan KK, Catalona WJ and Landerson JH. Differences in values obtained with two assays of Prostate Specific Antigen." *J Urol.* 139, 762-72 (1988).
10. Steman JH, Lemonen J, Aitthan H, Rannikko S, Tuukkanen K and Aitthan O. "A complex between prostate specific antigen and  $\alpha$ -1-antitrypsin is the major form of prostate specific antigen in serum of patients with prostate cancer: assay of complex improves clinical sensitivity for cancer." *Cancer Res.* 51, 222-26 (1991).

Revision: 3 Date: 08/15/11 DCO: 0524

Cat #: 2125-300

Site	96(A)	96(B)
A1	1ml set	1ml set
B1	1 (13ml)	2 (13ml)
C1	1 plate	2 plates
D1	1 (20ml)	1 (20ml)
E1	1 (7ml)	2 (7ml)
F1	1 (7ml)	2 (7ml)
G1	1 (6ml)	2 (6ml)

For Orders and Inquiries, please contact


**Monobind Inc.**  
100 North Pointe Drive  
Lake Forest, CA 92630 USA

Tel: 949-951-2655 Email: [info@monobind.com](mailto:info@monobind.com)  
Fax: 949-951-3539 Web: [www.monobind.com](http://www.monobind.com)


Please visit our website to learn more about our other interesting products and services.

IVD   

EC REP **CEPARTNER** 3951 DB, 13, NL  
Tel: +31 (0) 6-516-536-26

**NSAI**

**NSAI**

Annex to Certificate Number: MD19.4585

# Certificate of Registration of Quality Management System to I.S. EN ISO 13485:2012

The National Standards Authority of Ireland certifies that:  
**Monobind Inc.**  
100 North Pointe Drive  
Lake Forest, CA 92630  
USA

has been assessed and deemed to comply with the requirements of the above standard in respect of the scope of operations given below:

**The Design, Manufacture and Distribution of In-Vitro Diagnostic Medical Device Immunoassays, and Related Reagent Products and Equipment**

**Scope of Registration:**  
The Design, Manufacture and Distribution of In-Vitro Diagnostic Medical Device Immunoassays, and Related Reagent Products and Equipment

Activity	Location
Headquarters, Design, Manufacture	Monobind Inc. 100 North Pointe Drive Lake Forest, CA 92630 USA
Manufacture, Design	File No.: MD19.4585 Monobind Inc. 103 North Pointe Drive Lake Forest, CA 92630 USA
	File No.: MD19.4585/A

Additional sites covered under this multi-site certification are listed on the Annex (File No. MD19.4585)

Verified by:  
Operations Manager

Approved by:  
Geraldine Larkin  
Chief Executive Officer

Approved by:  
Susan Murphy  
European Medical Device  
Operations Manager

Registration Number: MD19.4585  
Certification Granted: May 18, 2010  
Effective Date: Oct 29, 2017  
Expiry Date: Oct 28, 2020

National Standards Authority of Ireland, 1 Swire Square, Northwood, Santry, Dublin 9, Ireland T +353 1 607 3800

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**TO WHOM IT MAY CONCERN**

**This declaration has been established for GBG-MLD Moldova**

We, MACHEREY-NAGEL GMBH & CO KG, Neumann-Neander-Str. 6-8, 52355Düren, GERMANY, are Original Manufacturers of

- *Filtration products*
- *Rapid Test products*
- *MEDI-TEST products*
- *Chromatography products*

Our Authorized / (Non-)Exclusive Distributor in Moldova for the above mentioned products under MACHEREY-NAGEL brand, is the company

**"GBG-MLD" S.R.L.**  
**Tighina str.65, office 607**  
**MD-2001,Chisinau,**  
**Republic of Moldova**


Explicitly, GBG-MLD is allowed to take part and submit bids in official tenders for the goods manufactured by us in GBG-MLD's own name, risk and on their own account. MACHEREY-NAGEL gives no warranty regarding the fulfilment of any signed contracts or conditions granted by GBG-MLD

This declaration will remain valid up to 31.12.2019 and will automatically end this date without any termination or expiration notice given. In no event shall this declaration be extended automatically.

Düren, 25.02.2019

**MACHEREY-NAGEL GMBH & CO KG**

*Christos Evangelakakis*  
 Dr. Christos Evangelakakis  
 International Sales Manager

TÜVRheinland

### EC Certificate

Directive 98/79/EC Annex IV, excluding Sections 4 and 6  
Full Quality Assurance System  
In Vitro Diagnostic Medical Devices

Registration No.: HL 60119814 0001

Report No.: 21265422 001

**Manufacturer:** Macherer-Nagel GmbH & Co. KG  
Neumann-Neander-Str. 6-8  
52355 Düren  
Deutschland

**Products:**

Products for self-testing  
(see attachment for products and sites included)  
Replaces Certificate, Registration No.: HL 60076697 0001

**Expiry Date:**

2022-05-28

The Notified Body hereby declares that the requirements of Annex IV, excluding section 4 and 6 of the directive 98/79/EC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex IV, section 5 of the aforementioned directive. For placing on the market of List A devices covered by this certificate an EC design-examination certificate according to Annex IV, section 4 and a verification of manufactured products according to section 6 is required.

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TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 98/79/EC concerning in vitro diagnostic medical devices with the identification number 0197.

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Attachment to Certificate  
Registration No.: HL 60119814 0001  
Report No.: 21265422 001

Manufacturer: Macherer-Nagel GmbH & Co. KG  
Neumann-Neander-Str. 6-8  
52355 Düren  
Deutschland

Products for self-testing:  
- Single and multi-parameter disposable test strips for urine analysis  
- Indicator test strips and papers for measurement of pH in urine

Additional site for warehousing and logistics:

Bahnstr. 120  
52355 Düren, Germany

Notified Body

Dipl.-Ing. Sven Hoffmann

Date: 2017-05-29

CISQ is a member of

*IQNet, the association of the world's first class certification bodies, is the largest provider of management System Certification in the world. IQNet is composed of more than 30 bodies and counts over 150 subsidiaries all over the globe.*

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*Design and production of test tubes with predetermined vacuum for collection of haematological samples, biological liquids and urine samples. Production of test tubes for micro-collection of haematological samples. Design and production of Holders for vacuum sampling. Design and production of diagnostic kits for blood and biological liquids analysis. Design and production of culture media for microbiology. Design and production of sterile needles and devices for collection of haematological samples. Trading of the products of the Group: diagnostic kits, culture media for microbiology, plastic disposable labware, test tubes with predetermined vacuum and sterile needles. Design and production of moulds for plastic labware. Injection moulding of thermoplastic materials for medical devices.*

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For timely and updated information about any changes in the certification status referred to in this certificate, please contact the number +39 02 725341 or email address info@icim.it.

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"GBG-MLD" SRL  
STR. TIGHINA 65  
2001 CHISINAU  
MOLDOVA

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### DISTRIBUTOR AGREEMENT

To whom it may concern, we hereby declare that:

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appoints "GBG-MLD" SRL – STR. TIGHINA 65. - 2001 CHISINAU –MOLDOVA

as authorized distributor of KIMA plastic labware products in the territory of MOLDOVA

GBG MLD has the right to import and distribute KIMA plastic labware products.

This Agreement is valid one (2) years from the present date.

The Distributor does not have any possibility to oblige the company KIMA sas with quantities or delivery time as well as prices without prior written authorization from KIMA sas.

KIMA sas keeps the right to modify the prices according to the market of the raw materials.

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**"VACUTEST KIMA" vacuum blood and biological liquids  
collection tubes in plastic.**

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**Hereby we declare**

under our sole responsibility that the above mentioned devices meet the applicable provisions of the Directive 98/79/EC as amended on "In Vitro Diagnostic Medical Devices".

All the supporting documents, as required by Annex III, in order to prove conformity to the Essential Requirements as listed in Annex I, are retained under the premises of the Manufacturer

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0035188

## Capillary blood collection system

### Capillary blood collection GK

#### Capillary blood collection GK

It consists of a prepared test vessel and a prepared plastic end-to-end collection capillary with stopper.

- Easy handling
- Plastic capillary\* with exact filling volume; complete inner surface prepared, unbreakable
- Collection vessel serves as centrifugal vessel; prepared for all common tests
- Light protected, tinted vessels for bilirubin analyses

\* with conformity certificate in accordance with the Weights and Measurement Regulations

Fill the capillary with capillary blood from a horizontal position.

After filling let the blood flow into the vessel from a vertical position (shake out remaining blood).

Remove the capillary, press on the attached stopper and mix or centrifuge.


077003 EDTA 100 µl

Packing unit:

100 pcs in bag,

4,000 pcs in box

Ø 11/8 x 39 mm

# Certificate

The Certification Body of  
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization  
**KABE LABORTECHNIK GmbH**  
Jägerhofstr. 17  
51588 Nümbrecht  
Deutschland

has established and applies a quality management system for medical devices  
for the following scope:

**Design and development, production and distribution of  
in vitro diagnostic devices and consumption materials  
for sample withdrawal, preparation and storage  
as well as single-use medical devices**

Proof has been furnished that the requirements specified in

**EN ISO 13485:2016**

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2018-10-16  
Certificate Registration No.: SX 60133221 0001  
An audit was performed. Report No.: 21234760 009  
This Certificate is valid until: 2021-10-15

Certification Body

Date 2018-10-12

Dipl.-Ing. I. Munkler

**TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg**  
Tel.: +49 221 806-1371 Fax: +49 221 808-3935 e-mail: cert-validity@de.tuv.com http://www.tuv.com/safety

TÜVRheinland

### EC Certificate

Directive 93/42/EEC Annex II, excluding Section 4

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#### Products:

- Cannulas for blood collection
  - MBU Capillaries
- (see attachment for details)

Replaces Approval, Registration No.: HD 60034211 0001

#### Expiry Date:

2020-10-15

The Notified Body hereby declares that the requirements of Annex II, excluding section 4 of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex II, section 5 of the aforementioned directive. For placing on the market of class III devices covered by this certificate an EC design-examination certificate according to Annex II, section 4 is required.

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TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

TÜVRheinland

Doc. 1/1, Rev. 0

TÜV Rheinland

LGA Products GmbH

Tillystraße 2, 90431 Nürnberg

Attachment to Certificate  
Registration No.: HD 60105393 0001  
Report No.: 21234760 001

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#### Products included:

- Cannulas for blood collection

For the following devices the scope covers only the aspects of the manufacture concerned with the securing and maintaining sterile conditions:


- MBU Capillaries

Notified Body

Dipl.-Ing. I. Munkler

Date: 2015-10-16

Certificado ES10/81672

The management system of

# DELTALAB, S.L.

Pol. Ind. La Llana, Plaza De La Verneda, 1  
08191 Rubí (Barcelona)

has been assessed and certified as meeting the requirements of

## ISO 9001:2015

For the following activities

**Design, manufacture and sale of laboratory material for the collection, transport and conservation of samples for microbiological, molecular biology, haematology, biochemistry, histology, microscopy and colorimetric analysis.**  
**Commercialization of equipment for the storage of prepared samples, cryogenic stored samples, general labware and industrial packages.**

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in/ from the following sites

**Pol. Ind. La Llana, Plaza De La Verneda 1 - 08191 Rubí (Barcelona)**

This certificate is valid from  
29 November 2017 until 11 October 2019.  
Issue 7. Certified since October 2010.

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Page 1 of 1

This document is issued by the Company subject to its General Conditions of Certification Services accessible at [www.sgs.com/terms\\_and\\_conditions.htm](http://www.sgs.com/terms_and_conditions.htm). Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at <http://www.sgs.com/en/certified-clients-and-products/certified-client-directory>. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Certificado ES10/81672

The management system of

# DELTALAB, S.L.

Pol. Ind. La Llana, Plaza De La Verneda, 1  
08191 Rubí (Barcelona)

has been assessed and certified as meeting the requirements of

## ISO 9001:2015

For the following activities

**Design, manufacture and sale of laboratory material for the collection, transport and conservation of samples for microbiological, molecular biology, haematology, biochemistry, histology, microscopy and colorimetric analysis.**  
**Commercialization of equipment for the storage of prepared samples, cryogenic stored samples, general labware and industrial packages.**

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**UNI CEI EN ISO 13485-2016 (ISO 13485-2016)**

per i seguenti Processi  
*concerning the following kinds of Processes*

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*Management of manufacturing and placing on the market of sterile tampons for sampling of biological specimens in natural orifice and in surgical field. Design and manufacturing of diagnostic medical devices for analysis laboratories.*

*Marketing of medical and diagnostic devices in vitro.*

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
*Management of the manufacturing and placing on the market of sterile tampons for sampling of biological specimens in natural orifice and in surgical field. Design and manufacturing of diagnostic medical devices for laboratories of analysis. Marketing of medical and diagnostic devices in vitro. Marketing of laboratory articles.*

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## DECLARATION OF CONFORMITY FOR MATERIALS AND OBJECTS INTENDED TO COME IN CONTACT WITH FOOD

Hereby we declare that

1. Nuova Aptaca Srl in Vitro Medical Diagnostic Devices (Directive 98/79/CE) and Medical Device (93/42/CE):

Are conform to the following CE community legislation:

- 1935/2004/CE Rule ( objects intended to come in contact with food) and following updates and changes
- 10/2011 Regulation (specific migration limits) and following updates and changes
- 1895/2005/CE Rule (substances use restriction for food contact) and following updates and changes
- 2002/72/CE Directive and following updating and changes (food contact) and following updates and changes
- 2005/84/CE Directive (substances use restriction – phthalates, sulphates) and following updates and changes
- 1272/2008 Regulation ( labeling and use of dangerous substances) and following updates and changes
- 2011/65/UE Directive (heavy metals, RoHS) and following updating and changes

Conform to the following Italian legislation:

- 21/03/1973 Ministerial Decree and following updating and changes
  - 777/82 DPR and following updating and changes
2. During Devices manufacturing no materials containing natural rubber, latex, synthetic rubber are used (except for Articles of latex)
  3. Devices are produced with materials that do not contain substances submitted to restrictions as per the above legislations and respect the global and specific migration limits in accordance with the following conditions:
    - Simulant A ( distilled water) -40°C for 10 days
    - Simulant B ( acetic acid solution 3% p/v) – 40°C for 10 days
    - Simulant C ( Ethyl alcohol solution 10% v/v) – 40°C for 10 days
    - Simulant D1 (ethyl alcohol solution at 50% v/v) - 40°C for 10 days
    - Simulant D2 (Vegetable oil - Try substitute made with 95% ethyl alcohol as indicated by the Italian Ministerial Decree 34 of 21.03.1973) - 40°C for 10 days

The global migration limit, together with all other specific restrictions which monomers and/or additives present in the material can be exposed to, are respected in the use conditions here above.

This is supported with analytical tests made in accordance with 1935/2004/CE, 2023/2006/CE, 2007/19/CE, 10/2011/CE and following updating and changes and Italian Ministerial Decree n° 34 21 March 1973, where applicable.

Notes and/or simulant used for migration tests, indicated on the previous point of the present declaration, allow to fix the food or the group of food, admitted to the contact with food.

This declaration is issued in conformity with the legislation mentioned here above.

The use of the material, indicated in the present declaration, in industrial or commercial premises does not exclude the control of its conformity to the regulations in force and the control of its technological suitability for the target which it is designed for.

This declaration has a validity starting from the date mentioned here below and will be substituted only if there will be important changes in the material production able to change some requirements essential for the conformity and only if the legislative references mentioned in this declaration will be changed and updated so that a new control for the conformity will be needed.

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## CE Registration Certificate

This is to certify that, in accordance with the In Vitro Diagnostic Medical Device Directive 98/79/EC, Emergo Europe agrees to perform all duties and responsibilities as the Authorized Representative for

**ERMA Inc.**  
**2-31-6 Yushima, Bunkyo-ku**  
**Tokyo, 113-0034**  
**Japan**

as stipulated and demanded by the aforementioned Directive. The Dutch Competent Authorities have received the In Vitro Diagnostic Medical Device Registrations on the following dates:

**21 September 2007**  
**See attached product listing**

**Emergo Europe Registration Number: NL/CA01/601529**

The Manufacturer has provided Emergo Europe with the appropriate Declaration(s) of Conformity confirming that the In Vitro Diagnostic Medical Devices fulfill the applicable requirements of Directive 98/79/EC.

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# MANAGEMENT SYSTEM CERTIFICATE

Certificate No:  
59878-2009-AQ-MCW-FINAS

Initial certification date:  
20 December 2000

Valid:  
21 June 2018 - 31 August 2021

This is to certify that the management system of

## **THERMO FISHER SCIENTIFIC, JCS**

Kubinskaya 73, liter A, build.1, Saint-Petersburg, Russian Federation, 196240

has been found to conform to the Quality Management System standard:  
**ISO 9001:2015**

This certificate is valid for the following scope:

**MANUFACTURING OF LIQUID HANDLING PRODUCTS AND SPECIAL  
DIAGNOSTIC PLASTICS.**

Place and date:  
**Moscow, 21 June 2018**

**FINAS**  
Finnish Accreditation Service  
S001 (EN ISO/IEC 17021)

For the issuing office:  
**DNV GL – Business Assurance**  
Trekhpudny per. 9 build. 2, office 406,  
Moscow, Russian Federation

*S. Groubine*

**Serguei Groubine**  
Management Representative

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