MENDEC[®] AQUA

MENDEC[®] AQUA CARTRIDGE

TECRES[®]
ADVANCING HIGH TECHNOLOGY

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
To the attention of the surgeon

DESCRIPTION

MENDEC AQUA is a sterile single-use device for extrusion into the bone of acrylic resins or bone void fillers (not supplied with the kit) intended for vertebral augmentation procedures (vertebroplasty, kyphoplasty, screw-augmentation).

MENDEC AQUA is composed of: a hydraulic gun equipped with a tube (1,2,3) with a final connector (4), a cartridge (6,7,8) for loading the acrylic resin/bone void filler and an aspirating needle (5) of the saline solution (not supplied with the kit).

If the acrylic resin/bone void filler loaded in a single cartridge is not sufficient to complete the surgical procedure, a second cartridge is available singly packaged (Mendec Aqua Cartridge).

Utilisable tube length = 155 cm, Utilisable cartridge volume = 9.9 cc

INTENDED USE

Delivery of acrylic resins/bone void fillers in vertebral augmentation procedures (vertebroplasty, kyphoplasty, screw-augmentation).

PRECAUTIONS FOR USE:

- Read the instructions for use carefully before using the MENDEC AQUA device.
- Strictly aseptic surgical techniques must be observed.
- Make sure that the package is intact and that the components have not been tampered with in any way.
- Do not use after the expiry date printed as the product might be unsuitable for use.
- Before using MENDEC AQUA make sure that in the 24 hours prior to use it has been stored at a temperature of $23^{\circ}\text{C} \pm 1^{\circ}\text{C}$.
- For correct use of the device it is necessary to use only sterile saline solution (not supplied with the kit) as a thrust means for extrusion of the acrylic resin/bone void filler.
- Before using MENDEC AQUA for extrusion of acrylic resins or bone void fillers, tests must be carried out to verify the *working time* of the resins or bone void fillers.
- Sterility is not guaranteed if the unit container is damaged or opened. In these cases do not use the device.
- The device is single-use and intended for use on a single patient.

- Do not re-sterilize and/or re-utilise any of the components. The components should not be re-utilised as this involves a risk of infection for the patient and the performance level of the components is not guaranteed. Furthermore, re-sterilization can alter the device's morphology and/or its chemical/physical/mechanical properties, leading to a malfunctioning of the device with serious risks for the patient's health. Any residual material must be considered as surgical waste and therefore must be eliminated on conclusion of the surgical procedure.

WARNINGS

- Use only with extrusion devices provided with a Luer-Lock fitting for spine augmentation procedures (vertebroplasty, kyphoplasty, screw-augmentation).
- Do not use with bone void fillers containing solvents of PP (Polypropylene) and/or PA (Polyamide).

METHODS FOR USE OF THE DEVICE

The correct use of the device comprises 4 steps:

STEP 1: Loading the MENDEC AQUA hydraulic gun with saline solution

1. Immerse the connector in a container with sterile saline solution. If it is not possible to immerse the connector into a container with sterile solution, tightly fix the needle, supplied with the kit, to the connector.
2. Then unblock the knob by squeezing the lever (3) of the grip and pull the knob (1) to create a vacuum inside the syringe (2) so as to aspirate the sterile saline solution until the syringe is completely filled.
3. Newly unblock the knob by squeezing the lever on the grip and with the syringe on the gun facing upwards, push the knob to bleed the syringe of air.
4. Repeat the preceding steps up to when the syringe is completely filled with sterile saline solution. Lastly, remove the needle if present and temporarily replace the gun on the instrument table.

STEP 2: Loading the cartridge/s with acrylic resin/bone void filler

The cartridge/s can be loaded in 2 ways depending on whether the resin/bone void filler has been supplied in a mixing and extrusion device or whether it is to be manually mixed.

Note: In case two cartridges are needed, it is recommended to transfer the acrylic resin/bone void filler in the two cartridges soon after the mixing following the instructions below.

A) **Resin supplied in a mixing and extrusion device:** tightly screw the cartridge cap (8) onto the cartridge (6) without moving the piston (7); after having mixed the acrylic resin/bone void filler, in accordance with the IFUs of the selected device, connect the device to the cartridge of the Mendec Aqua and extrude the acrylic resin/bone void filler into it up to when the piston has reached the bottom of the cartridge. Refer to the instructions for use of the selected device; reduce the pressure of the cement and disconnect the device.

Warning: The devices compatible with the cartridge are those provided with a female Luer-Lock. If the device is not provided with a female Luer connection, it will not be possible to connect it to the Mendec Aqua cartridge and it will be necessary to fill the cartridge as if the cement had been manually mixed (see below).

B) **Manually-mixed resin:** manually push the piston onto the bottom of the cartridge, taking care to use sterile gloves; mix the acrylic resin/bone void filler in accordance with its IFUs and pour it into the cartridge up to about 2 mm from the top. Lastly tightly screw the cap onto the cartridge.

Warning: there will be some spillage of acrylic resin/bone void filler when the cap is screwed onto the cartridge. This is normal, though it is recommended to remove the excess of material.

STEP 3: Extrusion of the resin/bone void filler


1. Tightly fit the gun connector (previously loaded with sterile saline solution) to the cartridge, taking care to fully rotate the connector in a clockwise direction.
2. Test the partly-assembled system (on the instrument table) by extruding a small amount of acrylic resin/

bone void filler from the cartridge by rotating the pistol knob in a clockwise direction. This test guarantees that the acrylic resin/bone void filler is exiting correctly and safely; it further verifies that the connections have been correctly connected. **WARNING:** If there is any leakage of sterile solution or acrylic resin/bone void filler from parts other than the Luer connection, make a further verification of the preceding steps.

3. The extrusion with Mendec Aqua of the acrylic resin/bone void filler must be in accordance with its viscosity curve. Rotate the knob to newly extrude a small quantity of acrylic resin/bone void filler so as to verify that it has reached the necessary degree of viscosity.
4. Connect the cartridge to the extrusion device suitable for the preselected procedure (e.g. vertebroplasty needle, bone filler, cannulated pedicle screw) already positioned in the patient.
5. Retreat from the operating area up to the distance allowed by the length of the tube, while maintaining a clear view of the monitor.
6. Extrude the acrylic resin/bone void filler by slowly rotating the knob in a clockwise direction, while monitoring using radioscopic control. Slowly and smoothly dispense the cement up until the necessary filling amount in the vertebra has been achieved. A quantity of about 0.6 cc is dispensed for each complete cycle of the knob.
7. Dispensing of the cement can be halted at any time by reducing the pressure by rotating the knob in an anti-clockwise direction. To resume dispensing, newly rotate the knob in the clockwise direction.

Warning: The device is provided with an internal pressure limiter. If the flow of acrylic resin/bone void filler unexpectedly stops, the user should reduce the force and rotation velocity of the knob rotation by sliding the gun knob. If this is not sufficient to resume dispensing of the cement; the user should determine the cause of the blockage and if possible correct it as indicated in the following:

Possible causes of blockage:

- A) If the cement exits freely from the cartridge after having disconnected it from the extrusion device, the blockage is internal of the extrusion device. Replace the extrusion device.
- B) If the cement does not freely exit the cartridge when it is disconnected from the extrusion device, the probable cause of the blockage is an erroneous mixture of the acrylic resin/bone void filler or trapped air in the cartridge. Replace the device with a new Mendec Aqua.
- C) If following these actions it is not possible to continue dispensing, this means that the cement has reached an excessive viscosity.

STEP 4 (optional when using Mendec Aqua Cartridge):

1. reduce the pressure of Mendec Aqua by rotating the knob in an anti-clockwise direction;
2. disconnect the gun connector from the first cartridge, taking care to fully rotate the connector in a counterclockwise direction;
3. remove the first cartridge from the extrusion device (e.g. vertebroplasty needle, bone filler, cannulated pedicle screw);
4. load Mendec Aqua gun with sterile saline solution, as described in STEP 1;
5. extrude the resin/bone void filler from the second cartridge following the instructions described in STEP 3.

STEP 5: Disconnecting Mendec Aqua from the extrusion device

When the procedure has been completed, reduce the pressure inside the Mendec Aqua connected to the extrusion device by rotating the knob in an anti-clockwise direction; disconnect the connector from the cartridge and remove the cartridge from the extrusion device.

DISPOSAL


Disposal of the device should be in accordance with local waste regulations.

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TECRES S.p.A