

**HBsAg EIA Test Kit
Package Insert**

English

For professional in vitro diagnostic use only.

100

100

... ..

serum or plasma

[illegible][illegible]

100

- 100

FORMATION

- 100

- 1.0% Sodium Iodide

- Do not store St

10/25/2011

- covering the tra

1001

5 Stop Solution

Package Insert

- Disposable gloves

100

Configuration may d

absorbent pap

5	<ul style="list-style-type: none"> Add 50 µL of Substrate A to each well. (Clear Reagent) Add 50 µL of Substrate B to each well. (Clear Reagent) Substrate A and Substrate B can also be mixed together as a Working Substrate Solution before use. Add a volume of Substrate A to an equal volume of Substrate B in a clean glass or plastic vessel, mix well. Add 100 µL Working Substrate Solution to each well. Note: The Working Substrate Solution should be used within 30 minutes. Then a blue color should develop in wells containing Positive specimens. 	<ul style="list-style-type: none"> Add 50 µL of Substrate A to each well Add 50 µL of Substrate B to each well Or add 100 µL Working Substrate Solution to each well.
6	<ul style="list-style-type: none"> Mix gently then cover microwell plate with Plate Sealer and incubate in water bath or incubator using one of the following procedures: Standard Procedure: Incubate at 37°C ± 2°C for 10 minutes ± 1 minute. Enhanced Procedure: Incubate at 37°C ± 2°C for 30 minutes ± 2 minutes. 	<ul style="list-style-type: none"> Mix then cover microwell plate with Plate Sealer and incubate using one of the following procedures: Standard Procedure: Incubate at 37°C for 10 min Enhanced Procedure: Incubate at 37°C for 30 min
7	<ul style="list-style-type: none"> Remove the Plate Sealer. Add 50 µL of Stop Solution to each well. (Clear Reagent) Then a yellow color should develop in wells containing Positive specimens. 	<ul style="list-style-type: none"> Remove the Plate Sealer Add 50 µL of Stop Solution to each well
8	<ul style="list-style-type: none"> Read at 450/630-700 nm within 30 minutes. Note: Microwell plate can also be read at 450 nm, but it is strongly recommended to read it at 450/630-700 nm for better results. 	<ul style="list-style-type: none"> Read at 450/630-700 nm within 30 min

AUTOMATED PROCESSING

Automatic EIA microplate processors may be used to perform the assay after validating the results to ensure they are equivalent to those obtained using the manual method for the same specimens. Incubation times may vary depending on the processors used but do not program less incubation times than the procedure listed above. When automatic EIA microplate processors are used, periodic validation is recommended to ensure proper results.

VALIDATION REQUIREMENTS AND QUALITY CONTROL

1. Calculate the Mean Absorbance of Negative Control and Positive Control by referring to the table below.

Example of Negative Control Calculation	
Item	Absorbance
Negative Control: Well B1	0.023
Negative Control: Well C1	0.021
Total Absorbance of Negative Control	$0.023 + 0.021 = 0.044$
Mean Absorbance of Negative Control	$0.044/2 = 0.022$
Blank Absorbance: Well A1	0.002
NCx: Mean Absorbance of Negative Control - Blank Absorbance	$0.022 - 0.002 = 0.020$

2. Check the validation requirements below to determine if the test results are valid.

Validation Requirements	
Item	Validation Requirements
Blank Well	Blank Absorbance should be < 0.050 if read at 450/630-700 nm Note: It should be < 0.100 if read at 450 nm
Negative Control	Mean Absorbance after subtraction of Blank Absorbance should be < 0.100
Positive Control	Mean Absorbance after subtraction of Blank Absorbance should be > 1.000

NOTE: The test results are considered invalid if the above validation requirements are not met. Repeat the test or contact your local distributor.

3. Calculate the Cut-Off Value using the following formula if the test results are valid.

Example of Cut-Off Value Calculation

Item	Absorbance
NCx	0.020
Cut-Off Value: NCx + 0.070	$0.020 + 0.070 = 0.090$

INTERPRETATION OF RESULTS

NON-REACTIVE: Specimens with absorbance less than the Cut-Off Value are non-reactive for HBSAg and may be considered negative.

REACTIVE: Specimens with absorbance greater than or equal to the Cut-Off Value are considered initially reactive for HBSAg. The specimen should be retested in duplicate before final interpretation. Specimens that are reactive in at least one of the re-test are presumed to be repeatedly reactive and should be confirmed using other HBV markers or confirmatory testing. Specimens that are non-reactive on both retests should be considered non-reactive.

NOTE: Specimens with values within ±10% of the Cut-Off Value should be retested in duplicates for final interpretation.

LIMITATIONS

- The HBSAg EIA Test Kit is used for the detection of HBSAg in human serum or plasma. Diagnosis of an infectious disease should not be established based on a single test result. Further testing, including confirmatory testing, should be performed before a specimen is considered positive. A non-reactive test result does not exclude the possibility of exposure. Specimens containing precipitate may give inconsistent test results. Mutated HBSAg may not be detected by the test.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- As with other sensitive immunoassays, there is the possibility that non-repeatable reactive results may occur due to inadequate washing. The results may be affected due to procedural or instrument error.
- False positive results may occur due to high titers of Heterophilic Anti Mouse Antibodies (HAMAs).
- Erroneous result may also be due to fibrin particles and microbial contamination.
- False negative results may occur if the quantity of HBSAg present in the specimen is lower than the analytical sensitivity of the test, or if HBSAg is not present during the stage of disease when the specimen was collected.
- The Positive Control in the test kit is not to be used to quantify assay sensitivity. The Positive Control is used to verify that the test kit components are capable of detecting a reactive specimen provided the procedure is followed as defined in the kit and the storage conditions have been strictly adhered to.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

The analytical sensitivity of the HBSAg EIA Test Kit has been determined using reference HBSAg standards, ad and sy subtypes. The analytical sensitivity is 0.2 IU/mL using the standard procedure and 0.1 IU/mL using the enhanced procedure, which were all confirmed using the WHO NISBC International Standard with code number 01476-011-WIL for HBSAg.

Clinical Sensitivity and Specificity

The HBSAg EIA Test Kit has correctly identified specimens of a seroconversion panel and has been compared with a leading commercial HBSAg EIA test using clinical specimens. The results show that the clinical sensitivity of the HBSAg EIA Test Kit is >99.9%, and the clinical specificity is 99.9%.

HBSAg EIA vs. Other EIA

Method	Results		Other EIA		Total Results
	Positive	Negative	Positive	Negative	
HBSAg EIA	562	3	562	3	565
	0	5,234	0	5,234	5,234
Total Results	562	5,237	562	5,237	5,799

Clinical Sensitivity: >99.9% (99.4-100.0%)*
Overall Agreement: 99.9% (99.9-100.0%)*

Clinical Specificity: 99.9% (99.8-100.0%)*
95% Confidence Interval

Reproducibility

Intra-Assay: Within-run precision has been determined by using 10 replicates of three specimens: a low positive, a medium positive, and a high positive.

Inter-Assay: Between-run precision has been determined by 3 independent assays on the same three specimens: a low positive, a medium positive, and a high positive. Three different lots of the HBSAg EIA Test Kit have been tested using these specimens over a 5-day period.

Specimen	Intra-Assay			Inter-Assay		
	Mean Absorbance/ Cut-Off	Standard Deviation	Coefficient of Variation (%)	Mean Absorbance/ Cut-Off	Standard Deviation	Coefficient of Variation (%)
1	1.467	0.008	6.061	1.367	0.011	8.943
2	12.622	0.060	5.282	13.378	0.091	7.568
3	26.722	0.096	3.992	25.467	0.225	9.817

BIBLIOGRAPHY

- Blumberg, B.S. The Discovery of Australian Antigen and its Relation to Viral Hepatitis. *Virology*, 1971;7:223.
- Krugman, S., Giles J.P. Viral Hepatitis, Type B (MS-2-Strain). *Further Observations on Natural History and Prevention*. New England Journal of Medicine, 288, 755.
- Krugman, S., Oestreich, L.R., et al. *Viral Hepatitis Type B Studies On Natural History and Prevention Re-examined*. New England Journal of Medicine, 300, 101.

Index of Symbols

	Consult instructions for use		Tests per kit		Manufacturer
	diagnostic use only		Lot Number		Catalog #
	Store between 2-8°C		Substrate A		Conjugate
	Wash Buffer (25%)		Plate Sealer		Stop Solution
	Negative Control		Package Insert		Positive Control



ACON Laboratories, Inc.
10125 Mesa Rim Road,
San Diego, CA 92121, USA



An enzyme immunoassay (EIA) for the qualitative detection of IgG antibodies to Hepatitis C Virus (HCV) in human serum or plasma.
For professional in vitro diagnostic use only.

INTENDED USE

The HCV Antibody ELIA Test Kit is a qualitative enzyme immunoassay for the detection of IgG antibodies to Hepatitis C Virus (HCV) in human serum or plasma. It is intended for screening and as an aid in the diagnosis of possible Hepatitis C infection.

SUMMARY

Hepatitis C Virus is a small, enveloped, positive-sense, single-stranded RNA virus. HCV is now known to be the major cause of parenterally transmitted non-A, non-B hepatitis. HCV infection causes a wide variety of chronic liver disease, cirrhosis and liver cancer. The main route of transmission of the virus is via transfusion of blood and blood products, organ transplantation, and sharing contaminated needles and syringes. Antibodies to HCV are found in over 80% of patients with well-documented non-A, non-B hepatitis. Cloning the viral genome has made it possible to develop serologic assays that use recombinant antigens.^{1,2} Compared to the first generation HCV ELIA tests using single recombinant antigen, new serologic tests incorporate recombinant protein and/or synthetic peptide antigens to avoid nonspecific cross-reactivity and to increase the sensitivity.³ The HCV Antibody ELIA Test Kit is a third generation immunoassay for the qualitative detection of the presence of IgG antibodies to HCV in serum or plasma specimen. The test utilizes recombinant HCV antigens encoded by the genes for both structural (nucleocapsid) and non-structural proteins to selectively detect antibodies to HCV in serum or plasma.

PRINCIPLE

The HCV Antibody ELIA Test Kit is a solid phase qualitative indirect simultaneous enzyme immunoassay for the detection of IgG antibodies to HCV in human serum or plasma. The microwell plate is coated with HCV recombinant antigens. During testing, the specimen diluent and the specimens are added to the antigen coated microwell plate and then incubated. If the specimens contain antibodies to HCV, it will bind to the antigens coated on the microwell plate to form immobilized antigen-HCV antibody complexes. If the specimens do not contain antibodies to HCV, the complexes will not be formed. After initial incubation, the microwell plate is washed to remove unbound materials. The enzyme-conjugated anti-human IgG antibodies are added to the microwell plate and incubated. The enzyme-conjugated anti-human IgG antibodies will bind to the microwell antigen-HCV antibody complexes present. After the second incubation, the microwell plate is washed to remove unbound materials. Substrate A and substrate B are added and then incubated to produce a blue color indicating the amount of HCV antibodies present in the specimen. Sulfuric acid solution is added to the microwell plate to stop the reaction producing a color change from blue to yellow. The color intensity, which corresponds to the amount of HCV antibodies present in the specimen, is measured with a microplate reader at 450/630-700 nm or 450 nm.

PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after expiration date.
- Do not mix reagents from other kits with different lot numbers.
- Avoid cross contamination between reagents to ensure valid test results.
- Follow the wash procedure to ensure optimum assay performance.
- Use Plate Sealer to cover microwell plate during incubation to minimize evaporation.
- Use a new pipet tip for each specimen assayed.
- Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate. Do not allow wells to dry out during the assay procedure.
- Do not touch the bottom of the wells with pipette tips. Do not touch the bottom of the microwell plate with fingertips.
- Do not allow sodium hypochlorite fumes from chlorine bleach or other sources to contact the microwell plate during the assay as the color reaction may be inhibited.
- All equipment should be used with care, calibrated regularly and maintained following the equipment manufacturer's instructions.

HEALTH AND SAFETY INFORMATION

- Some components of this kit contain human blood derivatives. No known test method can offer complete assurance that products derived from human blood will not transmit infectious agents. Therefore, all blood derivatives should be considered potentially infectious. It is recommended that these reagents and human specimens be handled using established good laboratory working practices. Wear disposable gloves and other protective clothing such as laboratory coats and eye protection while handling kit reagents and specimens. Wash hands thoroughly when finished.
- ProClin™ 300 is included as a preservative in the Conjugate, Concentrated Wash Buffer, Specimen Diluent, Substrate and Controls. Avoid any contact with skin or eyes.
- Do not eat, drink or smoke in the area where the specimens or kits are handled. Do not pipette by mouth.

- Avoid any contact of Substrate A, Substrate B, and Stop Solution with skin or mucosa. The Stop Solution contains 0.5M sulfuric acid which is a strong acid. If spills occur, wipe immediately with large amounts of water. If the acid contacts the skin or eyes, flush with large amounts of water and seek medical attention.
- Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for one hour at 121°C. Disposable should be autoclaved or incinerated. Do not autoclave materials containing sodium hypochlorite.
- Handle and dispose all specimens and materials used to perform the test as if they contained infectious agents. Observe established precautions against microbiological hazards throughout all the procedures and follow the standard procedures for proper disposal of specimens.
- Observe Good Laboratory Practices when handling chemicals and potentially infectious material. Discard all contaminated material, specimens and reagents of human origin after proper decontamination and by following local, state and federal regulations.
- Neutralized acids and other liquids should be decontaminated by adding sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to a 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.

STORAGE AND STABILITY

- Unopened test kits should be stored at 2-8°C upon receipt. All unopened reagents are stable through the expiration date printed on the box if stored between 2-8°C. Once opened, all reagents are stable for up to 3 months after the first opening date if stored between 2-8°C. Return reagents to 2-8°C immediately after use.
- Allow the sealed pouch to reach room temperature before opening the pouch and remove the required number of strips to prevent condensation of the microwell plate. The remaining unused strips should be stored in the original resealable pouch with desiccant supplied at 2-8°C and can be used within 3 months of the opening date. Return the remaining unused strips and supplied desiccant to the original resealable pouch. Firmly press the seal closure to seal the pouch completely and immediately store at 2-8°C.
- Concentrated Wash Buffer may be stored at room temperature to avoid crystallization. If crystals are present, warm up the solution at 37°C. Working Wash Buffer is stable for 2 weeks at room temperature.
- Do not expose reagents especially the Substrate to strong light or hypochlorite fumes during storage or incubation steps.
- Do not store Stop Solution in a shallow dish or return it to the original bottle after use.

SPECIMEN COLLECTION AND PREPARATION

- The HCV Antibody ELIA Test Kit can be performed using only human serum or plasma collected from venipuncture whole blood.
- EDTA, sodium heparin, and ACD collection tubes may be used to collect venipuncture whole blood and plasma specimens. The preservative sodium azide neutralizes horseradish peroxidase and may lead to erroneous results.
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Grossly hemolytic, lipidic or turbid samples should not be used. Specimens with extensive particulate should be clarified by centrifugation prior to use. Do not use specimens with thin particles or contaminated with microbial growth.
- Serum and plasma specimens may be stored at 2-8°C for up to 7 days prior to assaying. For long term storage, specimens should be kept frozen below -20°C.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

REAGENTS AND COMPONENTS

Materials Provided

No.	Reagent	Component Description	Quantity
1	HCV Conjugate	Anti-human IgG antibody bound to peroxidase.	1 x 12 ml
2	Concentrated Wash Buffer (25x)	Tris-HCl buffer containing 0.1% Tween 20 Preservative, 0.1% ProClin™ 300	1 x 50 ml
2A	Specimen Diluent	Tris buffer Preservative, 0.1% ProClin™ 300	1 x 12 ml
3	Substrate A	Citra-phosphate buffer containing hydrogen peroxide. Preservative, 0.1% ProClin™ 300	1 x 8 ml
4	Substrate B	Buffer containing tetraethylenediamine (TMB). Preservative, 0.1% ProClin™ 300	5 x 8 ml

5	Stop Solution	0.5M Sulfuric acid	1 x 8 ml	5 x 8 ml
6	HCV Negative Control	Normal serum non-reactive for HCV, HBeAg, HbS, and HbA2c	1 x 0.4 ml	5 x 0.4 ml
7	HCV Positive Control	Inactivated serum containing antibodies to HCV and negative for HBeAg, HbS, and HbA2c	1 x 0.4 ml	5 x 0.4 ml
	Plate sealers	Preservative, 0.1% ProClin™ 300	3	15
	Package insert		1	1

Materials Required But Not Provided

- Freshly distilled or deionized water
- Sodium hypochlorite solution for decontamination
- Absorbent paper or paper towel
- Water bath or incubator capable of maintaining 37°C ± 2°C
- Calibrated automatic or manual microwell plate washer capable of aspirating and dispensing 350 µl/well
- Disposable gloves
- Automated processor (optional)
- Calibrated micropipettes with disposable tips capable of dispensing 10, 50, and 100 µl
- Graduated cylinders for wash buffer dilution
- Vortex mixer for specimen mixing (optional)
- Disposable reagent reservoirs
- Calibrated microplate reader capable of reading at 450 nm with a 630-700 nm reference filter, or reading at 450 nm without a reference filter

DIRECTIONS FOR USE

Allow reagents and specimens to reach room temperature (15-30°C) prior to testing. The procedure must be strictly followed. Assay must proceed to completion within time limits. Arrange the controls so that well A1 is the Blank well. From well A1, arrange the controls in a horizontal or vertical configuration. The procedure below assigns specific wells arranged in a vertical configuration. Configuration may depend upon software.

Step	Detailed Procedure	Simplified Procedure
0	• Prepare Working Wash Buffer by diluting the Concentrated Wash Buffer 1:25. Pour the contents of the bottle containing the concentrated wash buffer in a graduated cylinder and fill it with freshly distilled or deionized water to 1250 ml for 96 wells/plate testing. The Working Wash Buffer is stable for 2 weeks at 15-30°C. Note: If crystals are present in the Concentrated Wash Buffer, warm it up at 37°C until all crystals dissolve. • Remove unused strips from the microwell plate, and store in the original resealable pouch at 2-8°C.	• Prepare Working Wash Buffer by diluting the Concentrated Wash Buffer 1:25 • Remove and store unused strips at 2-8°C
1	• Add 100 µl Specimen Diluent in respective wells including Negative Control, Positive Control, Blank and specimen wells. (Green Reagent)	• Add 100 µl Specimen Diluent
2	• Add 10 µl of Negative Control in wells B1 and C1. (Blue Reagent) • Add 10 µl of Positive Control in wells D1 and E1. (Red Reagent) • Add 10 µl of specimen to assigned wells starting at F1. Then a color change from green to blue will occur to verify that the specimen has been added.	• B1 and C1: Add 10 µl Negative Control • D1 and E1: Add 10 µl Positive Control • Starting F1: Add 10 µl specimen
3	• Mix gently by swirling the microwell plate on a flat bench for 30 seconds. • Cover the microwell plate with the Plate Sealer and incubate in a water bath or an incubator at 37°C ± 2°C for 30 minutes ± 2 minutes.	• Mix gently • Cover the microwell plate with the Plate Sealer and incubate at 37°C for 30 min
4	• Remove the Plate Sealer. • Wash each well 5 times by filling each well with 350 µl of Working Wash Buffer, then remove the liquid. • Turn the microwell plate upside down on absorbent tissue for a few seconds. Ensure that all wells have been completely washed and dried. Note: Improper washing may cause false results.	• Remove the Plate Sealer • Wash each well 5 times with 350 µl of Working Wash Buffer • Turn the microwell plate upside down on absorbent tissue for a few seconds
5	• Add 100 µl of Conjugate to each well except for blank wells. (Red Reagent)	• Add 100 µl of Conjugate to each well except for the Blank well



6	• Cover with microtiter plate with the Plate Sealer and incubate in a waterbath or an incubator at 37°C ± 2°C for 30 minutes.	• Repeat Step 4	• Cover the microtiter plate with the Plate Sealer and incubate at 37°C for 30 min
7	• Add 50 µl of Substrate A to each well (Clear Reagent). • Add 50 µl of Substrate B to each well (Clear Reagent). • Substrate A and Substrate B can also be mixed together as a Working Substrate Solution before use. Add a volume of Substrate A to an equal volume of Substrate B in a clean glass or plastic vessel, mix well. Add 100 µl Working Substrate Solution to each well. NOTE: The Working Substrate Solution should be used within 30 minutes.	• Repeat Step 4	• Add 50 µl of Substrate A to each well • Add 50 µl of Substrate B to each well • Or add 100 µl Working Substrate Solution to each well
8	• Mix gently then cover microtiter plate with Plate Sealer and incubate in a water bath or incubator at 37°C ± 2°C for 10 minutes ± 1 minute.	• Mix then cover microtiter plate with Plate Sealer and incubate at 37°C for 10 min	• Mix then cover microtiter plate with Plate Sealer and incubate at 37°C for 10 min
9	• Remove the Plate Sealer. • Add 50 µl of Stop Solution to each well. (Clear Reagent) Then a yellow color should develop in wells containing Positive specimens.	• Remove the Plate Sealer • Add 50 µl of Stop Solution to each well	• Remove the Plate Sealer • Add 50 µl of Stop Solution to each well
10	• Read at 450/630-700 nm within 30 minutes. Note: Microtiter plate can also be read at 450 nm, but it is strongly recommended to read it at 450/630-700 nm for better results.	• Read at 450/630-700 nm within 30 min	• Read at 450/630-700 nm within 30 min

AUTOMATED PROCESSING

Automatic EIA microplate processors may be used to perform the assay after validating the results to ensure they are equivalent to those obtained using the manual method for the same specimens. Incubation times may vary depending on the processors used but do not program less incubation times than the procedure listed above. When automatic EIA microplate processors are used, periodic validation is recommended to ensure proper results.

VALIDATION REQUIREMENTS AND QUALITY CONTROL

1. Calculate the Mean Absorbance of Negative Control and Positive Control by referring to the table below.

Example of Negative Control Calculation

Item	Absorbance
Negative Control Well B1	0.014
Negative Control Well C1	0.012
Total Absorbance of Negative Control	0.014 + 0.012 = 0.026
Mean Absorbance of Negative Control	0.026/2 = 0.013
Blank Absorbance: Well A1	0.006
NCX: Mean Absorbance of Negative Control - Blank Absorbance	0.013 - 0.006 = 0.007

2. Check the validation requirements below to determine if the test results are valid.

Item	Validation Requirements
Blank Well	Blank Absorbance should be < 0.050 if read at 450/630-700 nm Note: It should be < 0.100 if read at 450 nm
Negative Control	Mean Absorbance after subtraction of Blank Absorbance should be < 0.100
Positive Control	Mean Absorbance after subtraction of Blank Absorbance should be > 1.000

NOTE: The test results are considered invalid if the above validation requirements are not met. Repeat the test or contact your local distributor.

3. Calculate the Cut-Off Value using the following formula if the test results are valid.

Example of Cut-Off Value Calculation

Item	Absorbance
NCX	0.007
Cut-Off Value: NCX + 0.145	0.007 + 0.145 = 0.152

INTERPRETATION OF RESULTS

NON-REACTIVE: Specimens with absorbance less than the Cut-Off Value are non-reactive for antibodies to HCV and may be considered negative.

REACTIVE: Specimens with absorbance greater than or equal to the Cut-Off Value are considered initially reactive for antibodies to HCV. The specimen should be retested in duplicate before final interpretation. Specimens that are reactive in at least one of the re-test are presumed to be

repeatedly reactive and should be confirmed using confirmatory testing. Specimens that are non-reactive on both retests should be considered non-reactive.

NOTE: Specimens with values within ±10% of the Cut-Off Value should be retested in duplicates for final interpretation.

LIMITATIONS

- The HCV Antibody EIA Test Kit is used for the detection of antibodies to HCV in human serum or plasma. Diagnosis of an infectious disease should not be established based on a single test result. Further testing, including confirmatory testing, should be performed before a specimen is considered positive. A non-reactive test result does not exclude the possibility of exposure. Specimens containing precipitate may give inconsistent test results.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- As with other sensitive immunoassays, there is the possibility that non-repeatable reaction may occur due to inadequate washing. The results may be affected due to procedural or instrument error.
- The Positive Control in the test kit is not to be used to quantify assay sensitivity. The Positive Control is used to verify that the test kit components are capable of detecting a reactive specimen provided the procedure is followed as defined in the kit and the storage conditions have been strictly adhered to.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The HCV Antibody EIA Test Kit has correctly identified specimens of a seroconversion panel and has been compared to a leading commercial EIA HCV test using clinical specimens. The results show that the clinical sensitivity of the HCV Antibody EIA Test Kit is >99.9%, and the clinical specificity is 99.8%.

HCV Antibody EIA vs. Other EIA

Method	Results	Other EIA	Total Results
HCV Antibody EIA	Positive	Positive	307
	Negative	Negative	5,897
Total Results	297	5,907	6,204

Clinical Sensitivity: >99.9% (97.0-100.0%)*
Clinical Specificity: 99.8% (99.5-100.0%)*
Overall Agreement: 99.8% (99.6-100.0%)*
95% Confidence Interval

Reproducibility

Intra-Assay: Within-run precision has been determined by using 15 replicates of three specimens: a low positive, a medium positive and a high positive.

Inter-Assay: Between-run precision has been determined by 3 independent assays on the same three specimens: a low positive, a medium positive and a high positive. Three different lots of the HCV Antibody EIA Test Kit have been tested using these specimens over a 5-day period.

Specimen	Intra-Assay			Inter-Assay		
	Mean Absorbance/ Cut-Off	Standard Deviation (%)	Coefficient of Variation (%)	Mean Absorbance/ Cut-Off	Standard Deviation (%)	Coefficient of Variation (%)
1	3.259	0.213	6.535	3.731	0.312	8.362
2	6.168	0.404	6.549	7.811	0.630	8.086
3	16.712	0.970	5.804	14.445	0.983	6.805

BIBLIOGRAPHY

- Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science, 1989;244:359.
- Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science, 1988;244:362.
- Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N. Lelie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. Lancet, 1991;337:317.
- Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection. A Review. J. Clinical Immunology, 1993;16:204.

Index of Symbols

Symbol	Consult Instructions for use	Tests per kit	Manufacturer
1	For in vitro diagnostic use only	Use by	Store between 2-8°C
W	HCV	Lot Number	Catalog #
Wash Buffer (25x)	Wash Buffer (25x)	Substrate A	Substrate B
Control	Negative Control	Conjugate	Positive Control
Microtiter Plate	Microtiter Plate	Stop Solution	Package Insert
Specimen Diluent	Specimen Diluent	Plate Sealer	

ACOM

ACOM Laboratories, Inc.
10125 Mesa Rim Road
San Diego, CA 92121, USA



ATLAS C-REACTIVE PROTEIN (CRP) LATEX KIT

For the qualitative and semi-quantitative measurement of C-reactive protein (CRP) in human serum.

[IVD] For in-vitro diagnostic and professional use only

Store at 2-8°C

INTENDED USE

Atlas C-Reactive Protein (CRP) is used to measure the CRP in human serum qualitatively and semi-quantitatively.

INTRODUCTION

C-reactive protein (CRP), the classic acute-phase of human serum, is synthesized by hepatocytes. Normally, it is present only in trace amounts in serum, but it can increase as much as 1,000-fold in response to injury or infection. The clinical measurement of CRP in serum therefore appears to be a valuable screening test for organic disease and a sensitive index of disease activity in inflammatory, infective and ischemic conditions. Macleod and Avery found that antibody produced against purified CRP provided a more sensitive test than the C-polysaccharide assay. Since that time a number of immunological assays have been devised to measure CRP such as capillary precipitation, double immunodiffusion and radical immunodiffusion.

The CRP reagent kit is based on the principle of the latex agglutination assay described by Singer and Plotz. The major advantage of this method is the rapid two (2) minute reaction time.

PRINCIPLE

The CRP reagent kit is based on an immunological reaction between CRP Antisera bound to biologically inert latex particles and CRP in the test specimen. When serum containing greater than 6 mg/L CRP is mixed with the latex reagent, visible agglutination occurs.

MATERIALS

MATERIALS PROVIDED

- CRP Latex Reagent: Latex particles coated with goat IgG anti-human CRP, pH 8.2. MIX WELL BEFORE USE.

- CRP Positive Control Serum: A stabilized pre-diluted human serum containing >20mg/L CRP.
- CRP Negative Control Serum: A stabilized pre-diluted animal serum.
- Glass Slides.
- Stirring Sticks.

MATERIALS REQUIRED BUT NOT PROVIDED

- Mechanical rotator with adjustable speed at 80-100 r.p.m.
- Vortex mixer.
- Pipettes 50 µL.
- Glycine Buffer (20x): add one part to nineteen parts of distilled water before use.

PRECAUTIONS

- Reagents containing sodium azide may be combined with copper and lead plumbing to form highly explosive metal azides. Dispose of reagents by flushing with large amounts of water to prevent azide buildup.
- For In Vitro diagnostic use.
- Positive and negative controls prepared using human serum found negative for hepatitis B surface antigen (HBsAg) by FDA required test; however, handle controls as if potentially infectious.
- Accuracy of the test depends on the drop size of the latex reagent (40µL). Use only the dropper provided with the latex and hold perpendicularly when dispensing.
- Glass slides should be thoroughly rinsed with water and wiped with lint-free tissue after each use.

STORAGE AND STABILITY

- Reagents are stable until specified expiry date on bottle label when stored refrigerated (2 - 8°C). **DO NOT FREEZE.**
- The CRP latex reagent, once shaken must be uniform without visible clumping. When stored refrigerated, a slight sedimentation may occur and should be considered normal.
- Do not use the latex reagent or controls if they become contaminated.

SPECIMEN COLLECTION AND STORAGE

- Use fresh serum collected by centrifuging clotted blood.

- If the test cannot be carried out on the same day, store the specimen for 7 days at 2-8°C and for 3 months at -20°C.
- For longer periods the sample must be frozen.
- As in all serological tests, hemolytic or contaminated serum must not be used.
- Do not use plasma.

PROCEDURE

A. QUALITATIVE TEST:

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 40 µL of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
- Mix the CRP-latex reagent vigorously or on a vortex mixer before using and add one drop (40 µL) next to the samples to be tested.
- Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
- Place the slide on a mechanical rotator at 80-100 r.p.m. for 2 minutes. False positive results could appear if the test is read later than two minutes.

B. SEMI-QUANTITATIVE TEST:

- Make serial two fold dilutions of the sample in 9 g/L saline solution.
- Proceed for each dilution as in the qualitative method.

QUALITY CONTROL

Positive and Negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

All result different from the negative control result, will be considered as a positive.


INTERPRETATION OF RESULTS

A. QUALITATIVE TEST:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the CRP Negative Control.

A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared to the CRP Negative Control (Fig. 11).





Positive
Negative
Figure 1

B. Semi-QUANTITATIVE TEST:

The approximate CRP concentration in the patient sample is calculated as follows:

6xCRP titer = ---- mg/L

INTERFERENCES

NONE INTERFERING SUBSTANCES:

- Hemoglobin (10g/dl)
- Bilirubin(20mg/dl)
- Lipemia(10g/dl)
- Other substances interfere, such as RF (100IU/ml).

NOTE

- High CRP concentration samples may give negative results. Retest the sample again using a drop of 20µl.
- The strength of agglutination is not indicative of the CRP concentration in the samples tested.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

LIMITATIONS

1. Reaction time is critical. If reaction time exceeds two (2) minutes, drying of the reaction mixture may cause false positive results.
2. Freezing the CRP Latex Reagent will result in spontaneous agglutination.
3. Intensity of agglutination is not necessarily indicative of relative CRP concentration; therefore, screening reactions should not be graded.
4. A false negative can be attributed to a prozone phenomenon (antigen excess). It is recommended, therefore, to check all negative sera by retesting at a 1:10 dilution with glycine buffer.

REFERENCE VALUES
Up to 6 mg/L. Each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

- Sensitivity: 6(5-10) mg/L
- Prozone effect: No prozone effect was detected up to 1600 mg/L
- Diagnostic sensitivity: 95.6 %.
- Diagnostic specificity: 96.2 %.

REFERENCES

1. Pepys, M.B., Lancet 1:653 (1981).
2. Werner, M., Clin.Chem. Acta 25:299 (1969).
3. Macleod, C.M., et al., J. Exp. Med 73:191 (1941).
4. Wood, HF, et al., J. Clin. Invest. 30: 616 (1951).
5. Mancini, G., et al. Immunochemistry 2:235 (1965).
6. Singer, J.M., et al., Am. J. Med 21: 888 (1956).
7. Fischer, C.L., Gill, C.W., In Serum Protein Abnormalities. Boston, Little, Brown and Co., (1975).

ATLAS MEDICAL

William James House,
Cowley Road, Cambridge,
CB4 0WX, UK

Tel: +44 (0) 1223 858 910
Fax: +44 (0) 1223 858 524

PP1005A01

Rev H (06.06.2017)

	REF	Catalogue Number		Store at
	IVD	For In-Vitro Diagnostic use		Caution
		Number of tests in the pack		Read product insert before use
	LOT	Lot (batch) number		Manufacturer
		Fragile, handle with care		Expiry date
		Manufacturer fax number		Do not use if package is damaged
		Manufacturer telephone number		



ATLAS RHEUMATOID FACTOR (RF) LATEX KIT

latex slide test for the qualitative and semi-quantitative measurement of RF in human serum.

IVD For In-Vitro diagnostic and professional use only

Store at 2-8°C

INTENDED USE

A latex slide test for the qualitative and semi-quantitative measurement of RF in human serum.

INTRODUCTION

Rheumatoid factors (RF) are antibodies directed against antigenic sites in the Fc fragment of human and animal IgG. Their frequent occurrence in rheumatoid arthritis makes them useful for diagnosis and monitoring of the disease. One method used for rheumatoid factor detection is based on the ability of rheumatoid arthritis sera to agglutinate sensitized sheep red cells, as observed by Waaler and Rose. A more sensitive reagent consisting of biologically inert latex beads coated with human gamma globulin was later described by Singer and Plotz. The RF kit is based on the principle of the latex agglutination assay of Singer and Plotz. The major advantage of this method is rapid performance (2 minute reaction time) and lack of heterophile antibody interference.

PRINCIPLE

The RF reagent is based on an immunological reaction between human IgG bound to biologically inert latex particles and rheumatoid factors in the test specimen. When serum containing rheumatoid factors is mixed with the latex reagent, visible agglutination occurs.

MATERIALS

MATERIALS PROVIDED

- RF Latex Reagent: Latex particles coated with human gamma-globulin, pH, 8.2. Preservative. Contains N, N-dimethylformamide.
- RF Positive Control Serum: Human serum with a RF concentration > 30 IU/mL. Preservative.

- RF Negative Control Serum: Animal serum.
- Reaction Slide
- Stirring sticks

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer
- Test Tubes (for dilution)
- Serological pipettes (for sample addition and for dilution)
- Rotator (optional)
- Glycine Buffer (20x): add one part to nineteen parts of distilled water before use.

PRECAUTIONS

- All reagents contain 0.1 % (w/v) sodium azide as a preservative.
- Reagents containing sodium azide may be combined with copper and lead plumbing to form highly explosive metal azides. Dispose of reagents by flushing with large amounts of water to prevent azide buildup.
- For In Vitro diagnostic use.
- Positive and negative controls prepared using human serum found negative for hepatitis B surface antigen (HBsAg) by FDA required test; however, handle controls as if potentially infectious.
- Accuracy of the test depends on the drop size of the latex reagent (40 µl). Use only the dropper supplied with latex and hold it perpendicularly when dispensing.
- Use a clean pipette tip and stirring stick for each specimen, and glass slides should be thoroughly rinsed with water and wiped with lint-free tissue after each use.
- Check reactivity of the reagent using the controls provided.

STORAGE AND STABILITY

- Reagents are stable until specified expiry date on bottle label when stored refrigerated (2-8°C).
- Do not freeze.
- The RF latex reagent, once shaken must be uniform without visible clumping. When stored refrigerated, a slight sedimentation may occur and should be considered normal.
- Do not use the latex reagent or controls if they become contaminated.

SPECIMEN COLLECTION AND STORAGE

- Use fresh serum collected by centrifuging clotted blood.
- If the test cannot be carried out on the same day, store the specimen for 7 days at 2-8°C and for 3 months at -20°C.
- As in all serological tests, hemolytic or contaminated serum must not be used.
- Do not use PLASMA.

PROCEDURE

Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
 - Place 50 µl of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
 - Mix the RF-latex reagent rigorously or on a vortex mixer before using and add one drop (50 µl) next to the sample to be tested.
 - Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
 - Place the slide on a mechanical rotator at 80-100 r.p.m. for 2 minutes. False positive results could appear if the test is read later than two minutes.
- Semi-quantitative method**
- Make serial two fold dilutions of the sample in 9 g/L saline solution.
 - Proceed for each dilution as in the qualitative method.

READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator. The presence of agglutination indicates a RF concentration equal or greater than 8 IU/mL (Note 1). The titer, in the semi-quantitative method is defined as the highest dilution showing a positive result.

CALCULATIONS

The approximate RF concentration in the patient sample is calculated as follows:

$$8 \times \text{RF Titer} = \text{IU/mL}$$



INTERFERENCES

NON INTERFERING SUBSTANCES:

- Hemoglobin (10g/dl)
 - Bilirubin(20mg/dl)
 - Lipemia(10g/dl)
- Other substances may interfere.

QUALITY CONTROL

1. RF Positive and Negative Control should be included in each test batch.
2. Acceptable performance is indicated when a uniform milky suspension with no agglutination is observed with the RF Negative Control and agglutination with large aggregates is observed with the RF Positive Control.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

8(6-16) IU/ml, under the described assay conditions.

PROZONE EFFECT

No prozone effect was detected up to 1500 IU/ml.

DIAGNOSTIC SENSITIVITY

100%.

DIAGNOSTIC SPECIFICITY

100%.

The diagnostic sensitivity and specificity have been obtained using 118 samples compared with the same method of a computer.

LIMITATIONS

- Reaction time is critical. If reaction time exceeds 2 minutes, drying of the reaction mixture may cause false positive result.
- Freezing the RF Latex Reagent will result in spontaneous agglutination.
- Intensity of agglutination is not necessarily indicative of relative RF concentration; therefore, screening reactions should not be graded.
- Increased levels of RF may be found in some diseases other than rheumatoid arthritis such as infectious mononucleosis, sarcoidosis, lupus erythrematosus, Sjogren's syndrome.
- Certain patients with rheumatoid arthritis will not have the RF present in their serum.

- The incidence of false positive results is about 3-5 %.
- Individuals suffering from infectious mononucleosis, hepatitis, syphilis as well as elderly people may give positive results.
- Diagnosis should not be solely based on the results of latex method but also should be complemented with a Waaler Rose test along with the clinical examination.

REFERENCE VALUES

Up to 8 IU/ml. Each laboratory should establish its own reference range.

NOTES

1. Results obtained with a latex method do not compare with those obtained with Waaler Rose test. Differences in the results between methods do not reflect differences in the ability to detect rheumatoid factors.

REFERENCES

1. Robert W Dorner et al. Clinica Chimica Acta 1987; 167: 1 – 21.
2. Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951-960.
3. Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528 – 534.
4. Adalbert F. Schubart et al. The New England Journal of Medicine 1959; 261: 363 – 368.
5. Charles M. Plotz 1956; American Journal of Medicine; 21:893 – 896.
6. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.














ATLAS MEDICAL

William James Hous Cowley Road,
Cambridge, CB4 0WY, UK.

Tel: +44 (0) 1223 858 910

Fax: +44 (0) 1223 858 524

PP1008A01, Rev H (17.06.2017)

	Catalogue Number		Store at
	For In-Vitro Diagnostic use		Caution
	Number of tests in the pack		Read product insert before use
	Lot (batch) number		Manufacturer
	Fragile, handle with care		Expiry date
	Manufacturer fax number		Do not use if package is damaged
	Manufacturer telephone number		



ANTISTREPTOLYSIN-O (ASO) LATEX SLIDE

TEST

For the qualitative and quantitative measurement of antibodies to Antistreptolysin-O in human serum.

IVD For in-vitro diagnostic and professional use only

Store at 2-8°C

INTENDED USE

ATLAS ANTISTREPTOLYSIN-O (ASO) latex slide Test is used for the qualitative and quantitative measurement of antibodies to Antistreptolysin-O in human serum.

INTRODUCTION

The group A β -hemolytic streptococci produces various toxins that can act as antigens. One of these exotoxins streptolysin-O, was discovered by Todd in 1932.

A person infected with group A -hemolytic streptococci produces specific antibodies against these exotoxins, one of which is antistreptolysin-O. The quantity of this antibody in a patient's serum will establish the degree of infection due to the -hemolytic streptococcal.

The usual procedure for the determination of the antistreptolysin titer is based on the inhibitory effect that the patient's serum produces on the hemolytic power of a pre-titrated and reduced streptolysin-O. However, the antigen-antibody reaction occurs independently of the hemolytic activity of streptolysin-O. This property enables the establishment of a qualitative and quantitative test for the determination of the antistreptolysin-O by agglutination of latex particles on slide.

PRINCIPLE

ASO test method is based on an immunologic reaction between streptococcal exotoxins bound to biologically inert latex particles and streptococcal antibodies in the test sample. Visible agglutination occurs when increased antibody level, are present in the test specimen.

MATERIALS

MATERIALS PROVIDED

- ASO Latex Reagent: latex particles coated with streptolysin O, pH, 8.2. Preservative
- ASO Positive Control (Red cap): Human serum with an ASO concentration > 200 IU/mL. Preservative
- ASO Negative Control (Blue cap) Animal serum. Preservative
- Reaction Slide.
- Stirring Sticks.

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer.
- Test Tubes 12x75mm.
- Test Tube Rack.
- Serological pipettes.
- High intensity light.
- Saline Solution, 0.9% NaCl.

PRECAUTIONS

- All reagents contain 0.1% (w/v) sodium azide as a preservative. Store all reagents at 2-8°C. **DO NOT FREEZE.**

- Reagents containing sodium azide may be combined with copper and lead plumbing to form highly explosive metal azides. Dispose of reagents by flushing with large amounts of water to prevent azide build-up.

- For In Vitro diagnostic use.
- Positive and negative controls prepared using human serum found negative for hepatitis B surface antigen (HBsAg) and HIV-III by FDA required test; however, handle controls as if potentially infectious.

REAGENT STORAGE AND STABILITY

- Reagents are stable until specified expiry date on bottle label when stored refrigerated (2-8°C).
- **DO NOT FREEZE.**
- The ASO Latex Reagent, once shaken must be uniform without visible clumping. When stored refrigerated, a slight sedimentation may occur and should be considered normal.
- Do not use the latex reagent or controls if they become contaminated.

SPECIMEN COLLECTION AND STORAGE

- Use fresh serum collected by centrifuging clotted blood.
- If the test cannot be carried out on the same day, store the specimen for 7 days at 2-8°C and for 3 months at -20°C.

- For longer periods the sample must be frozen.
- As in all serological tests, hemolytic or contaminated serum must not be used.
- **DO NOT USE PLASMA.**

PROCEDURE

Qualitative method

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50 μ L of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
3. Mix the ASO-latex reagent vigorously or on a vortex mixer before using and add one drop (50 μ L) next to the sample to be tested.
4. Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
5. Place the slide on a mechanical rotator at 80-100 r.p.m. for 2 minutes. False positive results could appear if the test is read later than two minutes.

Semi-quantitative method

1. Make serial two fold dilutions of the sample in 9 g/L saline solution.
2. Proceed for each dilution as in the qualitative method.

QUALITY CONTROL

Positive and Negative Controls should be included in each test batch.

Acceptable performance is indicated when a uniform milky suspension with no agglutination is observed with the ASO Negative Control and agglutination with large aggregates is observed with the ASO Positive Control.


RESULTS

A. QUALITATIVE TEST:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the ASO Negative Control.

A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared to the ASO Negative Control.





Positive
Negative
Figure 1

B. QUANTITATIVE TEST

A positive reaction is indicated by any observable agglutination in the reaction mixture. Record the last dilution showing a positive reaction. Concentration of ASO can be determined by multiplying the last positive dilution factor of the sample with the concentration of the positive control (200 IU/ml).

The titer of the serum is the reciprocal of the highest dilution which exhibits a positive reaction.

IU/ml of sample = conc. of positive control (200) x specimen titer

DILUTION	IU/ml
1:1	200
1:2	400
1:4	800
1:8	1600
Etc.	

REFERENCE VALUES

Up to 200 IU/ml (adults) and 100 IU/ml (children < 5 years old)⁶. Each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity:

200 (±50) IU/ml.

PROZONE EFFECT

No prozone effect was detected up to 1500 IU/ml.

SENSITIVITY

98%.

SPECIFICITY

97%.

INTERFERENCES

NON INTERFERING SUBSTANCES:

- Hemoglobin (10g/dl)
 - Bilirubin (20mg/dl)
 - Lipemia (10g/dl)
- Other substances may interfere

REFERENCES

1. Haftejee . Quarterly Journal of Medicine 1992. New series 84; 305: 641-658.
2. Ahmed Samir et al. Pediatric Annals 1992; 21: 835-842.
3. Spaun J et al. Bull Wild Hlth Org 1961; 24: 271-279.
4. The association of Clinical Pathologists 1961. Broadsheet 34.
5. Picard B et al. La Presse Medicale 1983; 23: 2-6.
6. Klein GC. Applied Microbiology 1971; 21: 999-1001.
7. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

ATLAS Medical

William James House,
Cowley Road, Cambridge, CB4 4WX, UK
Tel: ++44 (0) 1223 858 910
Fax: ++44 (0) 1223 858 524

PP1003A01

Rev H (09.09.2017)

	REF	Catalogue Number		Store at
	IND	For In-Vitro Diagnostic use		Caution
	LOT	Number of tests in the pack		Read product insert before use
	LOT	Lot (batch) number		Manufacturer
		Fragile, handle with care		Expiry date
		Manufacturer fax number		Do not use if package is damaged
		Manufacturer telephone		





по применению набора калпаторов для определения концентрации фибриногена набором реагентов «МультиТех-Фибриноген»

ИНСТРУКЦИЯ

Набор предназначен для получения калпаторовых значений времени свертывания при определении концентрации фибриногена в плазме крови модифицированным методом Class без предварительного разведения исследуемой плазмы на автоматических и полуавтоматических коагулометрах. Фибриноген-калпатор не предназначен для калпаторовых и других методов определения концентрации фибриногена, в том числе набора «Тех-Фибриноген-тест».

НАЗНАЧЕНИЕ

Включается в определение времени свертывания цитратной плазмы избыток тромбина (модифицированный метод Class). Время свертывания при этом пропорционально концентрации фибриногена, которую определяют по калпаторовому графику.

Состав набора:

1. Калпатор №1 (иофильно высушенный) - 1 фла.
 2. Калпатор №2 (иофильно высушенный) - 1 фла.
 3. Калпатор №3 (иофильно высушенный) - 1 фла.
 4. Калпатор №4 (иофильно высушенный) - 1 фла.
 5. Калпатор №5 (иофильно высушенный) - 1 фла.
- Концентрация фибриногена для каждого калпатора указана в Листопре к набору.

ХАРАКТЕРИСТИКА НАБОРА

Принцип метода

Время свертывания при этом пропорционально концентрации фибриногена, которую определяют по калпаторовому графику.

АНАЛИТИЧЕСКИЕ

ХАРАКТЕРИСТИКИ НАБОРА

Линейность определения: 0,9-10,0 г/л.
Коэффициент вариации результатов определения концентрации фибриногена при использовании набора калпаторов не превышает 10 %.
Допустимый разброс результатов определения концентрации фибриногена в одной пробе плазмы наборами калпаторов одной серии не превышает 10 %.

МЕРЫ

ПРЕДОСТОРОЖНОСТИ

Потенциальный риск применения набора – класс 2а (ГОСТ Р 51609-2000).
Все реагенты, входящие в набор, используются только для применения *in vitro*.
Все компоненты набора в используемых концентрациях не токсичны.
При работе с набором следует надевать одноразовые резиновые или пластиковые перчатки, так как образцы плазмы крови человека следует рассматривать как потенциально инфицированные, способы длительного время сохранять и передавать ВИЧ, вирус гепатита В или любой другой возбудитель вирусной инфекции.
Все использованные материалы дезинфицировать в соответствии с требованиями МУ-287-113.

ОБОРУДОВАНИЕ, МАТЕРИАЛЫ, РЕАГЕНТЫ

- Коагулометр;
- Дозаторы пипеточные на 1,0 мл;
- вода дистиллированная;
- перчатки резиновые хирургические;
- набор реагентов «МультиТех-Фибриноген» (заказывается дополнительно, кат. № 712 и кат. № 711).

Каталожный номер реагента:

714

ООО фирма «Технология-Стандарт»
456037, Борный, д/п 1351, тел./факс (3652) 22-99-37, 22-99-38, 22-99-39, 27-13-00

ПРИГОТОВЛЕНИЕ РЕАГЕНТОВ И ПРОВЕДЕНИЕ АНАЛИЗА

1. ПОДГОТОВКА РЕАГЕНТОВ К РАБОТЕ

В каждый из пяти флаконов калпаторов фибриногена внести по 1,0 мл дистиллированной воды и растворить со держимое при комнатной температуре (+18...+25°С) и слабом помешивании в течение 15 мин. В результате получают образцы с указанной в Листопре к набору калпаторов концентрации фибриногена.

2. ПРОВЕДЕНИЕ АНАЛИЗА

Используя инструкцию для набора реагентов необходимо определить время свертывания разведенных калпаторов №1, №2, №3, №4 и №5.
Для построения калпаторовой кривой необходим набор для определения фибриногена «МультиТех-Фибриноген» (заказывается дополнительно).
В зависимости от типа коагулометра существуют два варианта набора реагентов «МультиТех-Фибриноген» для автоматических (кат. № 712) и полуавтоматических (кат. № 711) коагулометров.

3. ЧТЕНИЕ РЕЗУЛЬТАТОВ

Время свертывания калпаторового образца плазмы составляет 5-100 с, в зависимости от концентрации фибриногена.

УСЛОВИЯ ХРАНЕНИЯ И ПРИМЕНЕНИЯ


Набор рассчитан на выполнение не менее 10 калпаторовых кривых при расходе по 0,1 мл на одно исследование. Однако в большинстве ситуаций рекомендуется дублирование результатов для построения калпаторовой кривой. Хранение набора должно проводиться при температуре +2...+8°С в течение всего срока годности набора (15 мес). Допускается транспортировка при температуре до +25°С в течение 30 сут.

После разведения раствора калпаторов пригодны для построения калпаторовой кривой в течение 4 часов при комнатной температуре. Разведенные калпаторы не следует замораживать.

ЛИТЕРАТУРА

1. Баркаган З.С., Момот А.П. Диагностика и контроль рутинной терапии нарушений гемостаза. - М.: «НьюДиамед-АО», 2008. - 292 с.
2. Момот А.П. Патология гемостаза. Принципы и алгоритмы клинико-лабораторной диагностики. - СПб.: Формат, 2006. - 208 с.





МультиТех-Финанс

ИНСТРУКЦИЯ
по применению набора реагентов для
определения концентрации фибриногена
(для полуавтоматических коагулометров)

НАЗНАЧЕНИЕ

Набор предназначен для количественного определения содержания фибриногена в плазме крови на полуавтоматическом анализаторе.

ХАРАКТЕРИСТИКА НАБОРА

Принцип метода

Заключается в определении времени свертывания цитрат-ной плазмы изъятком тромбина (модифицированный метод

Class). Время сверявания при этом пропорционально количеству точек, которые необходимо считать. В зависимости от количества точек, которые необходимо считать, время сверявания может быть разным. Например, если необходимо считать 10 точек, то время сверявания будет в 10 раз больше, чем если необходимо считать 1 точку. Таким образом, время сверявания пропорционально количеству точек, которые необходимо считать. Это можно выразить следующей формулой:

Состав набора:

1. Тромбин (лиофилизно высушенный реагент, 500 ед. NIH) - 2 фл.

7. Растворитель для трюмина, 10,5 мл - 2 фл.

АНТИИЗНАШЕНСКИЕ
ХАРАКТЕРИСТИКИ НАБОРА

АНАЛИТИЧЕСКИЕ

КАРАКТЕРИСТИКИ НАБОРА

коэффициент вариации результатов определения кон-
центрации фибриногена не превышает 10 %.

Содержание гепарина в плазме до 1,0 ЕД/мл не влияет на результаты определения.

Допустимый разброс результатов определения концентрации фибриногена в одной пробе плазмы разными наборами

данной серии не превышает 10 %.

МЕРЫ
ПРЕДОСТОРОЖНОСТИ

Потенциальный риск применения набора – класс 2a (ГОСТ

Все реагенты, входящие в набор, используются только для

Все компоненты набора в используемых концентрациях не

При работе с набором следует надевать очки.

...и/или лифтиковые нехватки, так как образцы п/з-

инфицированные, клинические длительные время сохранять и передавать ВИЧ, вирус гепатита В или любой другой возбудитель.

Все использованные материалы дезинфицировать в со-

ответствии с требованиями МУ-287-113.

ОБОРУДОВАНИЕ,
МАТЕРИАЛЫ, РЕАГЕНТЫ

- Центрифуга лабораторная;

- Дозаторы пипеточные на 0,1-0,2, 10,0 мл;

- вода дисциplinированная;
- печатки резиновые хирюпречские;
- набор казачий.

«Фирма «Финно-Калибр» (кат. 714, производитель ООО фирма «Технология-Стандарт»,

...экспрессируется дополнительно).

ИПНТОБЛЕНИЕ
АНАЛИЗИРУЕМЫХ ОБРАЗЦОВ

Кровь для исследования забирают из локтевой вены в

8 % раствор натрия лимоннокислого трехзамещенного (цит-

на - 9:1. Кровь центрифугируют при 3000-4000 об/мин

200 г) в течение 15 мин. В результате получают белую

Каталожный номер набора: 711

ООО фирма "Технология-Стандарт"

22-99-36, 22-99-37, 22-99-38, 22-99-39, 27-13-00