STORZ

KARL STORZ – ENDOSCOPE

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our products, this product is also the result of years of experience and careful manufacturing. You and your organization have decided in favor of a modern high quality product from KARL STORZ. This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the ENDOFLATOR® 40. All of the necessary details and actions are clearly explained. Please read these instructions carefully. Keep this manual in a convenient and conspicuous location and in its protective casing close to the device.

RxOnly Federal (USA) law restricts this device to sale by or on the order of a physician.

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It is recommended to check the suitability of the product for the intended procedure prior to use.

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


Symbols on label and packaging of the device

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4 General information

4.1 Description of the device

The ENDOFLATOR® 40 is an insufflation device for universal application in laparoscopic and thoracoscopic examinations and operations as well as endoscopy of the upper and lower gastrointestinal tract (e.g., TEO®, colonoscopy). In addition, the ENDOFLATOR® 40 is designed to expel ambient air in open and endoscopically assisted cardiac surgery as well as endoscopic vessel harvesting. The very latest technologies for pressure and flow measurement and control allow for various operating modes, including modes tailored to specific situations.

The 'PEDIATRICS' and 'HIGH FLOW' modes enable various applications to be accommodated. Easy handling and clearly arranged controls together with several safety circuits ensure maximum patient safety.

For particularly sensitive applications, the device features a 'PEDIATRICS' mode which is characterized by low insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow ranges (up to max. 15 l/min). This mode also includes tighter control over the flow rate at low flows.

To quickly compensate for the considerable loss of gas occurring in complex laparoscopic operations, the ENDOFLATOR® 40 has been designed for a high flow rate of up to 40 l/min. The ENDOFLATOR® 40 is operated directly via a color screen with a touch-sensitive surface (touch screen) or by an external prime control unit via the SCB interface. The operation displays the information required in the specific context to enable a comprehensive assessment of the situation in a clear manner, and guides and aids the user with all device-related work.

Thus the current patient pressure, current flow and the amount of gas used during the current intervention are displayed simultaneously during operation.

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4.1.1 Property rights

This product is protected in the USA by (at least one of the following) US Patent(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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5 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

5.1 Explanation of warnings and cautions

Please read this instruction manual through carefully and follow its instructions exactly. The words Warning, Caution, and Note convey special meanings. Wherever they are used in this instruction manual, the following text must be carefully reviewed to ensure the safe and effective operation of this device. In addition, to make the signal words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note provides special information regarding the operation of the product, or clarifies important information.

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WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel, or yourself at risk.

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Safety instructions Warnings and cautions

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WARNING: Always unplug the device from the power supply before carrying out any cleaning and maintenance work.

WARNING: Danger of electric shock! Do not open the device. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only. Any opening of the device by unauthorized persons will void warranty claims.

WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

WARNING: Only operate the device under the specified environmental conditions.

WARNING: Always keep a spare device to hand in case the first device should fail.

WARNING: For safety reasons, do not simultaneously touch the device output sockets and the patient.

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WARNING: Always use a hydrophobic bacterial filter to prevent cross-contamination of patients.

WARNING: You must insert a sterile CO₂ filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. Replace the CO₂ filter after every use.

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Before using the device on the patient for the first time it is imperative that you be acquainted with how the device operates and is controlled.

5.2 Intended use

CO₂ insufflators and their accessories are used to create and maintain a cavity during diagnostic or therapeutic endoscopic interventions as well as to expel ambient air during open, endoscopically assisted or endoscopic interventions.

5.2.1 Indication

The insufflator is used to create a cavity in the following diagnostic and therapeutic interventions:

- Laparoscopy
- Thoracoscopy
- Endoscopy of the upper and lower gastrointestinal tract: (e.g., TEO®, colonoscopy)
- Endoscopic vessel harvesting

In addition, the ENDOFLATOR® 40 is also designed to expel ambient air in open and endoscopically assisted cardiac surgery.

WARNING: Only sterile CO₂ gas may be used.

Use of the device in fields other than those indicated above is not allowed for safety reasons. The ENDOFLATOR® 40 may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

WARNING: No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.

WARNING: Unauthorized conversions or modifications to the device are not allowed for safety reasons.

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WÄRME: Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see Clause 16 of IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

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5.2.3 Warnings

5.2.3.1 Warnings

Excessive Absorption of CO₂ from Excessive Pressure/Flow

Excessive absorption of CO₂ results from either excessive flow and/or excessive pressure. The abdomen can be adequately distended by pressure in the range of 15–20 mmHg. It is seldom necessary to use an abdominal pressure greater than 20 mmHg. Little intravasation should occur at these levels. Pressures over 20 mmHg are virtually never needed and will increase the amount and rapidity of intravasation. Adequate respirations help avoid problems related to CO₂.

Metabolic Acidosis and Resultant Cardiac Irregularity

Prolonged intra-abdominal pressures greater than 20 mmHg should be avoided. This can cause any of the following:

- Decreased respiration with compromised diaphragmatic excursion
- Decreased venous return
- Decreased cardiac output
- Azidosis

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Pneumolabium/pneumoskrotum:

PÄDIATRICS MODE

Children are at risk of a pneumolabium or pneumoskrotum.

Increased airway pressure:

PÄDIATRICS MODE

When laparoscopic procedures are performed on children, the increased intraabdominal pressure also increases the risk for higher airway pressures. Always strictly monitor respiration and airway function when performing laparoscopic procedures on children younger than 12 years of age.

Compression of the vena cava:

PÄDIATRICS MODE

When insufflating the abdomen of a child with medical CO₂, an increased risk of compressing the vena cava exists. This risk can be reduced by monitoring the systolic and diastolic blood pressure during the entire surgery.

Recommended work settings:

PÄDIATRICS MODE

The flow values listed* for laparoscopic procedures performed on newborns, infants, and children are only suggested values. The selection of the suitable flow and pressure values is solely the responsibility of the attending physician. However, adhering to the values listed* ensures an optimal performance of the PEDIATRICS MODE of the insufflator.

* see section 6.4.5

Gas flow limit: PEDIATRICS MODE

The gas flow may not exceed 14 l/min when performing a laparoscopy on newborns or patients weighing less than 25 kg approx. 55 US pounds.

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5.2.3.3 Warnings: Endoscopic Vessel Harvesting CO₂ absorption: Endoscopic Vessel Harvesting

Due to the special surgical procedures - start of the heart bypass operation, and the endoscopic removal of the vessel - special care has to be taken as CO₂ is always absorbed through the tissue of the patient during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption.

Instruments: Endoscopic Vessel Harvesting
Before using the insufflator to endoscopic harvest vessels, please check whether the used instrument is intended for CO₂ insufflation.

Metabolic and Cardiac Reactions:
Endoscopic Vessel Harvesting
Due to the special surgical conditions - start of the heart bypass surgery and vessel harvesting - it is especially important to remember the existing risk of metabolic acidosis when insufflating with CO₂. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

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Dehydration: Endoscopic Vessel Harvesting
Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

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Embolism: Endoscopic Vessel Harvesting
Improper placement of the insufflation instrument could cause insufflation of gas into a vessel, resulting in air or CO₂ embolisms. To reduce the risk of air or CO₂ embolism, perform initial insufflation at a low flow rate and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO₂ embolisms can also be caused by a high pressure. Avoid high-pressure settings and close damaged blood vessels at once.

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5.3 User qualification
The ENDOFLATOR® 40 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

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Must have been trained thoroughly at least once in the use of the device.
No physical impairments which limit or prevent perception of acoustic or visual alarm signals.

5.5 Intended conditions of use

5.5.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

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5.5.2 Other intended conditions

- Frequency of use: one or more times a day.
- Length of use: from a few minutes to several hours a day.
- Place of installation: positioning on a level, vibration-free surface.
- Mobility: can be moved if positioned on a cart.
- Combination: can be used on the patient at the same time as other equipment required for the operations.
- Control: can be controlled via the KARL STORZ SCB.
- The device is placed outside the sterile area.
- El equipo está fuera del área estéril.

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5.5.3 User position

- To operate the device, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$, at a distance of approx. 30–70 cm to the front panel.
- For observation of the actual values during the application, a visual distance from the device of 2 m is assumed.

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
5.6 Safety precautions at the installation site

The device may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

WARNING: DO NOT use the device in the presence of flammable anesthetics.

WARNING: The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

WARNING: The device must not be used in oxygenated environments.

This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with the applicable national regulations.

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5.7 Safety precautions when operating the device

It is the user's responsibility to make sure the equipment is safe and operates properly before use.

This device has been tested and found to comply with the EMC limits of the Medical Device Directive 93/42/EEC (EN 55011 Class B and IEC 60601-1-2).

Please observe the Electromagnetic Compatibility information in section 10.

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5.8 Safety features

The ENDOFLATOR® 40 has the following safety features:

5.8.1 Self-test

A self-check procedure is performed each time the device is switched on (for this the device must be supplied with gas). If the test is negative, a corresponding error message is shown and the device can no longer be used as normal. The following components are tested:

- Leakproofness of the system (valves)
- Checking pressure measurement
- Checking flow measurement

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5.8.2 Monitoring during operation

• **Overpressure:** If, during insufflation, the intracorporeal pressure rises above the set setpoint value by a limit dependent on the setpoint value for more than 3 seconds, an acoustic signal sounds (see chapter 8.2.2) and the gas supply is stopped. After a time period of 5 seconds, the overpressure is actively reduced via a blow-off valve or shown by a symbol that it is deactivated.

• **Mode-specific pressure monitoring**

Monitoring of gas supply

If the system recognizes a lack of CO₂ gas when switched on, a warning signal is given. If a missing gas supply is detected and insufflation is stopped, because the gas bottle is empty, a visible and audible alarm is monitored.

Monitoring of gas preheating

The internal gas preheating is monitored by a second circuit, independently of the control. If a fault arises the heating is switched off.

Control of the preheating of the gas

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NOTE: If the monitoring system recognizes an error which renders safe functioning seemingly impossible, the device is transferred into a "safe state". From this state, the device can only be reset by turning off the power supply.

6 Installation and operating instructions

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6.1 Unpacking the equipment

Carefully remove the ENDOFLATOR® 40 and the accessories from their packaging. Check the delivery for missing items and possible damage. Please file any complaints with the manufacturer or supplier immediately.
If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

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- Place the device on a flat surface (Fig. 1).
- ! WARNING:** Ensure that there is sufficient air circulation for the device.
- ! WARNING:** The device must be set up in such a way that the power cord can be simply pulled out at any time.

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
- ! WARNUNG:** Das Gerät nur mit der auf dem Typenschild angegebenen Spannung betreiben.
- Push the power cord into the power socket as far as it will go (Fig. 3).
- ! WARNUNG:** Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.
- ! WARNUNG:** The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.
- ! WARNUNG:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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4

6.3.3 KARL STORZ SCB

NOTA: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device. Therefore, to unplug the SCB cable, pull out the connector – do not pull directly on the cable!

Connect the SCB cable to an (equivalent) SCB socket on the back of the device (Fig. 4).

Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control device or other SCB devices (see KARL STORZ SCB control NEO System Instruction Manual).

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6.3.5 Connecting the CO₂ bottle

Connect the gas outlet of the CO₂ bottle to the device's gas connection using the high-pressure tube supplied (Fig. 6-7).

WARNING: CO₂ bottles connected up to the device must be secured to prevent them from falling over.

WARNING: Follow the safety instructions when changing the CO₂ bottle.

WARNING: The CO₂ bottle must be placed vertically while operating the device, otherwise the functioning of the device cannot be guaranteed.

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6.3.6 CO₂ bottles equipped with German-standard or ISO connection

Place the high-pressure tube (e.g. 20 4002 22; see the section 'Accessories') onto the outlet of the CO₂ bottle and screw tight using the universal wrench provided (Fig. 8).

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6.3.7 CO₂ bottles equipped with PIN-Index connection

Place the connecting piece of the high-pressure tube (e.g., 20 4000 22; see the section 'Accessories') onto the CO₂ bottle and screw this tight.

Use the wrench provided (see the illustration) to open the CO₂ bottle (Fig. 9).

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

6.3.8 Opening the valve of the CO₂ bottle, if necessary

Slowly open the valve stopcock of the CO₂ bottle by turning it approx. one-half turn counterclockwise (Fig. 10-11).

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[12]

6.3.9 Connecting to the central gas supply

Connect the low-pressure tube (e.g., UI001; see the section 'Accessories') to the device's gas connection and connect the other end directly to the wall outlet of the central gas supply (3.3...7 bar) (Fig. 12).

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6.4 Commissioning

NOTE: The gas supply must be connected when the device is switched on, otherwise an error message is emitted.

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NOTE: Changing/adjustment of the settings for language, volume and the gas source are described over the next pages.

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
- Den Netzschalter ① betätigen.
Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo (Fig. 13). After this, a white screen appears.

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- Press the power switch ①.

- Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo (Fig. 13). After this, a white screen appears.

The device is now running a self-test. Once the test has been successfully completed, an acoustic availability signal sounds and a mode selection screen appears (Fig. 14).

- Select the desired mode.

- ‘PEDIATRICS’ (see Fig. 18, page 27)
 - ‘HIGH FLOW’ (see Fig. 19, page 27)
- NOTE:** Selection of the desired mode is confirmed by an audio sound.

WARNING: For safety reasons, the device may only ever be used if, immediately after selecting the mode, the acoustic availability signal sounded!
See section 8.1.5.

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

NOTE: Normal use is only possible after successful conclusion of the tests.

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

If the self-test fails, appears a screen with an information message and an acoustic information signal sounds.

NOTE: If restarting (turning the device off and on again) fails, verify that the device is properly connected and the preparations have been made according to the instructions.

Further information on troubleshooting can be found in the chapter Technical description – Troubleshooting.

In normal cases

The switch-on message is followed by the selected screen PEDIATRICS or HIGH FLOW (Fig. 18-19).

NOTE: On the work screens 'PEDIATRICS' (Fig. 18) and 'HIGH FLOW' (Fig. 19), insufflation can be started directly, and device settings can be changed.

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6.4.2 Normal commissioning

After a successful start (see section 'Commissioning for the first time'), the selected 'Working screen 'PEDIATRICS' or 'HIGH FLOW' (Fig. 27 or 28) appears.

The device saves the operating parameters which were last used and reloads these (with a few exceptions) when the device is restarted (see tables below).

HIGH FLOW mode	Setting	Setting after restart
Pressure	1...20 mmHg	1...20 mmHg
Pressure	21...30 mmHg	20 mmHg
Flow	1...20 l/min	1...20 l/min
Flow	≥ 21 l/min	20 l/min

PEDIATRICS mode	Setting	Setting after restart
Pressure	1...12 mmHg	1...12 mmHg
Pressure	13...15 mmHg	12 mmHg
Flow	0,1...1 l/min	0,1...1 l/min
Flow	≥ 1,1 l/min	1 l/min

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NOTA: On the work screens 'PEDIATRICS' (Fig. 27) and 'HIGH FLOW' (Fig. 28), inflation can be started directly and device settings can be changed.

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- Attach the sterile CO₂ gas filter to the insufflation connection ③ (Fig. 29). **! WARNING:** You must insert a sterile CO₂ filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. The CO₂ filter must be replaced after each use.
- Attach the insufflation tube to the CO₂ gas filter (Fig. 29), and connect the other end of the insufflation tube to the VERESS needle.

NOTE: In the case of the disposable insufflation tube set, the filter is firmly connected to the tube.

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NOTE: The flow display may only show values greater than zero for short periods. If the set pressure is only reached slowly even when the instrument's stopcock is closed or if you note stationary flow rates of 0,2 l/min or more, there is a leak in the system filter – patient tube – instrument.

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6.4.3 Test for proper functioning

WARNING: Test the device prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.

Pressure stability and leakproofness

- Set the mode to 'PEDI/ATRICA'.
- Use 12 mmHg and 1.0 l/min as setpoint values.
- Check the connection of the patient tube to the device.
- Close the stopcock on the instrument.
- Start insufflation.
- Pressure then builds up in the tube which increases above the setpoint value and which, after a few seconds, is released via the blow off valve.
- Observe the instrument's conduct over a 30 second period.

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WARNING: Do not work with a leaky device. There is a risk that the intra-abdominal pressure will increase in an uncontrollable manner.

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6.4.4 Preparing insufflation Enter the operating parameters for the insufflation mode*, pressure and flow.

- If desired, change the selected working screen (»PEDIATRICS«/»HIGH FLOW« – Fig. 30 or 31) via the respective button.
- On the selected working screen (»PEDIATRICS«/»HIGH FLOW« – Fig. 30 or 31) set the setpoint values for pressure and flow using the ± buttons.

*The ENDOFLATOR® 40 is equipped with two different insufflation modes:

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- NOTA:** To compensate for the considerable loss of gas occurring in complex operations, the HIGH FLOW mode has been fitted with a flow rate of up to 40 l/min. For particularly sensitive applications (with small cavities), the device offers the PEDIATRICS' mode. This mode is characterized by a lower insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow range (up to max. 15 l/min).
- WARNING:** Always check the mode set. For safety reasons, use the PEDIATRICS' mode if the required cavity pressure is smaller or equal to 15 mmHg and a maximum flow of 15 l/min is sufficient.

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Installation and operating instructions

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6.4.5 Recommended settings for pediatric subcategories

Group	Weight	Flow range
Children younger than 1 year	Approx. 1...9 kg	0,1...0,5 l/min
Children between 1 and 3 years	10...15 kg	0,5...1,0 l/min
Children between 3 and 4 years	16...19 kg	1,0...2,0 l/min
Children between 4 and 14 years	≥ 20 kg	> 2,0 l/min
All Children	< 25 kg	≤ 14,0 l/min

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6.4.6 Recommended settings for endoscopic vessel harvesting

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Displays

The insufflation pressure, flow and volume of the gas used are continuously displayed in the Working screen (Fig. 32).

If the patient pressure rises above the preselected value, an audible warning signal will be sounded. The pressure is reduced via the overpressure valve after a delay of 5 s.

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6.4.7 Performing CO₂ insufflation

WARNING: Before starting insufflation rinse the entire system with 1 l CO₂.

WARNING: Always use the lowest possible pressure and flow for insufflation.

1. Reset the volume display using the '0' button (Fig. 32).
2. Insert the VERESS needle and open the locking lever. Tap the START/STOP button to initiate the insufflation process (Fig. 32).

NOTE: In order to use the maximum device flow in the HIGH FLOW mode, HiCap® instruments with a suitable connection must be used. The LUER-lock connection fitted to the end of the insufflation tube must be removed in this case.

NOTE: In the case of the disposable tubing set, cut the tube off at the mark behind the LUER-lock.

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
6.4.8 Decommissioning

Close the valve on the CO₂ bottle on ending the intervention (when using a high-pressure gas supply) and remove the insufflation tube (incl. filter) used from the device.

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6.5.1 Device settings

The Device Settings (Fig. 35) are described in the section on commissioning for the first time.

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
6.5.2 Service

The Service screen is protected by a password (Fig. 36).

- NOTE:** The service programs are only available to service personnel. The settings are described in the Service Manual.

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6.5.5 Device information

The following elements are displayed (Fig. 44):

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- 7 Maintenance**
- 7.1 Fuse replacement**
- Turn the device off and disconnect from the power supply (Fig. 45).

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- WARNING:** Only use fuses of the correct rating.

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7.2 Reprocessing

- UI 400 – ENDOFLATOR® 40
- UI 004 Insufflation tube, can be sterilized

General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: Risk of infection: The reusable insufflation tube is not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect the reusable insufflation tube for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the reusable insufflation tube before initial use and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: Risk of infection/product damage: KARL STORZ has validated the procedures given in the instruction manual. Only use these procedures. The reprocessing method must be validated by the user on site.

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7.2.2 Reprocessing of insufflation tube set with gas filter for single use, sterile

WARNING: Insufflation tube set with gas filter for single use, locks and hydrophobic bacteria filters must not be reused. Sterile disposable products are not approved for reprocessing.

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7.2.3 Reprocessing of reusable insufflation tube

WARNING: The reusable insufflation tube must be cleaned, disinfected and sterilized prior to use and subsequent reuse using validated reprocessing procedures.

7.2.3.1 Accessories

Accessories required for carrying out reprocessing:
Brushes: 27652, 27650 A, 27650 B

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Immediately after use of the medical device, remove heavy soiling, corrosive solutions and pharmaceuticals from the medical device. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

7.2.3.3 Manual precleaning Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

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7.2.3.6 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D):

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method. The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

Connecting

In order to ensure effective machine cleaning and disinfection, as well as adequate rinsing out of the instrument, it must be connected up to the washer and disinfector.

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[49]

cleaned and subjected to a full cleaning and disinfection procedure once more.

- Damaged medical devices must be withdrawn from use.
- There should be no holes, cuts, dents or scratches on the tube surface. If the reusable insufflation tube is damaged as described above, it must be replaced. If the reusable insufflation tube is discolored as a result of the reprocessing process, it must be removed immediately from the sterile supply cycle.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

The integrity of the reusable insufflation tube can be tested in the following way:

Attach a syringe to one end of the tube and kink the other end. Place the insufflation tube in completely demineralized water or microbiologically pure, sterile water. Force air into the insufflation tube using the syringe. If any bubbles appear in the water while forcing air through the reusable insufflation tube, the tube must not be reused.

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7.2.3.9 Sterilization

WARNING: Prior to the sterilization of the reusable insufflation tube, particularly in cases of manual reprocessing, it must be ensured that all chemical residue has been rinsed off or removed.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

KARL STORZ has validated and approved steam sterilization using the fractionated prevacuum procedure for this medical device

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132 °C - 137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes
- **NOTE:** This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.
- or 134 °C - 137 °C with a minimum exposure time of 3 minutes and maximum of 18 minutes

7.2.3.10 Limits of reprocessing

7.2.3.10 Limitación de la repreparación
The end of the product's lifespan is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

NOTE: The storage period for the CO₂ supply tubes^{*} is 5 years. After this, the supply tube must be replaced.
* Art. no. Low-pressure tubing:
UJ 001, UJ 002 and UJ 003
Art. no. High-pressure tubing:
20 4000 21, 20 4000 22, 20 4000 27,
20 4000 28 and 20 4002 22)

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7.3 Maintenance and safety check

7.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Information regarding maintenance services can be obtained from your local representative or from the manufacturer.

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7.3.2 Safety checks/repeat inspections as per IEC 62353

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WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that inscriptions relevant to safety are legible.

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Test for proper functioning

A functional inspection must be carried out in accordance with the instruction manual (see section 6.4.3) as well as checking the alarm functions in sections 8.3.1/8.3.2.

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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Servicing and repair


WARNING: Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ, all repair work must employ original KARL STORZ parts only.

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7.5 Disposal

This device has been marked in accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE).

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7.6 Repair program

During the repair period, you will generally be loaned a device, which must then be returned to KARL STORZ as soon as the repaired device is returned to you.

For repairs, please contact your local KARL STORZ subsidiary or authorized dealer. In Germany, you can refer the repairs directly to:
KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbeplatz 83
78579 Neuhausen
Service hotline: +49 7461/708 980
E-mail: technicalsupport@karlstorz.com

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7.7 Infection prevention

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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7.8 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ,
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the device has been used in accordance with its operating instructions at all times.

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7.9 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

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
Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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8 Technical description

8.1 Alarm specification

The device issues two alarms – ‘303: High pressure’ and ‘373: CO₂ empty’ (Fig. 50 and 51). The alarm signals are not self-retaining. This means that the alarm signal is only issued whilst the signal conditions are present. The alarms are displayed for at least 5 s and at least one sequence is played.

All alarm thresholds and delays are permanently programmed. In case of an operating pause or a voltage loss, the most recently set volume is applied.

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8.1.1 Alarm thresholds

Overpressure alarm

Medium priority overpressure alarm (‘303: High pressure’ – see 8.1) indicates that an overly high cavity pressure was measured at the device output. The overpressure alarm is activated 0.5 s after the cavity pressure measurement.

PEDIATRICS mode: The overpressure alarm is triggered if, as a product of the setpoint, the measured cavity pressure is 3 to 4 mmHg above the set setpoint value or if the cavity pressure is the same or larger than 16 mmHg.

HIGH FLOW mode: The overpressure alarm is triggered if, as a product of the setpoint, the measured cavity pressure is 4 to 5 mmHg above the set setpoint value or if the cavity pressure is the same or larger than 31 mmHg.

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	Setting	Reduction from an overpressure of
PEDIATRICS mode	> 12mmHg	4 mmHg
HIGH FLOW mode	> 26 mmHg	5 mmHg

If insufflation is turned off, the pressure is not reduced automatically.

WARNING: Observe the functioning of the device when an overpressure alarm is issued. If the insufflation or the blow-off valve is deactivated, suitable measures to reduce pressure must be initiated.

The following table shows the overpressure thresholds with triggering to one decimal place as the threshold is calculated depending on the setpoint value.

	Pressure set-point value	Overpressure alarm threshold
PEDIATRICS mode	1 mmHg ... 15 mmHg	2,6 mmHg ... 4 mmHg or pressure \geq 16 mmHg
HIGH FLOW mode	1 mmHg ... 30 mmHg	3,6 mmHg ... 5 mmHg or pressure \geq 31 mmHg

CO₂ empty alarm

Low priority CO₂ empty alarm ('373: CO₂ empty' – see 8.1) indicates that the gas supply failed during the application and that insufflation has been deactivated.

The following table shows the activation and deactivation thresholds for the CO₂ empty alarm.

Activation threshold	Input pressure < 1 bar
Deactivation threshold	Input pressure > 1,5 bar

WARNING: Reconnect the gas supply in the event of a CO₂ empty alarm and restart insufflation.

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
8.1.2 Visual signals

The alarms are shown in the title line (e.g., ‘373’ in Fig. 373).⁵²

Medium priority alarms are displayed with blue font on a yellow background; low priority alarms are displayed with blue font on a cyan-colored background.
Medium priority alarms overwrite low priority alarms.


8.1.3 Acoustic signals

The acoustic alarm signals are specified below. If a condition for an acoustic signal is met, at least one complete sequence is played. The volume of the audio signals can be adapted in ‘Settings’. Furthermore, the alarm audio signals can be deactivated for 30 s by tapping on the ‘Bell’ (alarm audio paused). The bell turns yellow during this time (Fig. 52).

At the maximum setting the alarm volume is 74 dBA and at the minimum setting 49 dBA.

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The medium priority alarm signal is a burst of 3 tones (frequency 254 Hz) each lasting 170 ms and with breaks of 230 ms. The burst is repeated every 10 seconds. Four different harmonics are generated.

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This means that the signal is only issued whilst the signal conditions are present. If the signal conditions are only very brief, the minimum duration of the display is likewise only 5 s.

The signals given are ordered according to priority. The alarms have a higher priority than the information signals. There are five priorities for the information signals.


Any low priority signal will be overwritten by an alarm signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If numerous signal conditions with the same priority occur, the most recently detected condition overwrites earlier ones in the title line of the display.

8.2.1 Visual signals

The information signals are also shown in the title line.

Information signals are displayed with blue font on a white background (inversely to normal display).

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
The acoustic information signal is a double tone c-e (263 Hz – 330 Hz) lasting 300 ms.

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
Key click

The pitch is modulated by $+/- 2$ Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

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
Depending on the priority of the message, the sequence is either repeated or played once. The volume of the maximum setting is 69 dBA. In the case of messages which display that the device is no longer operable or if the message ‘371’ (bottle pressure < 10 bar) is coupled with a repetitive audio, the sequence is repeated every 20 s. In all other cases, the sequence is played just once.

Availability signal


On completion of the successful self-test after switching on the device, an availability signal with the following characteristics sounds:

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8.4 Troubleshooting

WARNING: Always unplug the device before all maintenance work.

Description of error:

Complete failure of the device.

Possible causes:

1. Power supply failure.
2. Defective line fuse.

Measure:

1. Check that there is electricity to the wall outlet.
2. Change fuses as described in the instruction manual. Use only fuses of the correct rating.

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Device defective.

Measure:

Switch unit off and back on. If the warning message still appears after restarting, contact Service. Prior to this, open the System Log under Settings and note the ID, Word1, Word2 and information text of the last 5 entries.

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No gas flow.

Measure:

Observe the warnings (device temperature exceeded, temperature of high pressure regulator too low, gas supply not sufficient). If the warning message still appears after restarting, contact Service.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF
- Type of IP protection: IP21

According to IEC 60601-1-2, FCC-B:
Observe the information on electromagnetic compatibility in chapter 10 (pages 65-81).

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According to the Medical Device Directive (MDD):

This medical product belongs to Class IIb.
This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

- **NOTE:** The code number after the CE mark indicates the responsible notified body.

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8.6 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer. Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device. Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

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NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Spare parts, recommended accessories

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9 Spare parts, recommended accessories

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
WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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10 Electromagnetic Compatibility (EMC) information

The ENDOFLATOR® 40 UI 400 is suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

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Accessories and cables used for EMC Compliance

Cable type	Shielded	Length [m]	Ferrite	Used for
PE (potential equalization)	No	> 3	No	Potential equalization
SCB Cable	Yes	≥ 1	No	connecting SCB network, other SCB Units
Power Cord	No	3	No	Power supply

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WARNING: The use of an accessory or cable with the ENDOFLATOR® 40 other than those specified in this manual may result in increased emissions or decreased immunity of the ENDOFLATOR® 40. When using an accessory or cable with the ENDOFLATOR® 40 other than those specified in this manual, it becomes the responsibility of the user of the ENDOFLATOR® 40 to determine compliance with IEC 60601-1-2 when using this item.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ENDOFLATOR® 40 including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for Immunity tests

Guidance and manufacturer's declaration – electromagnetic immunity

The ENDOFLATOR® 40 is intended for use in the electromagnetic environment specified below. The user of ENDOFLATOR® 40 should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100 % for 250/300 cycles	Voltage Dip: Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ENDOFLATOR® 40 require continued operation during power mains interruptions, it is recommended that the ENDOFLATOR® 40 be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the ENDOFLATOR® 40 further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	6 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	

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Table 2
Test levels for Proximity fields from RF wireless communications equipment

Test Frequency MHz	Band MHz	Service	Modulation	Immunity Test Level V/m	Compliance level V/m
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27	27
450	430 – 470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz Sine Wave	28	28
710					
745	704 – 787	LTE Band 13 & 17	Pulse modulation 217 Hz	9	9
780					
810		GSM 800/900 TETRA 800 IDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28	28
870	800 – 960				
930					
1720		GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS	Pulse modulation 217 Hz	28	28
1845	1700 – 1990				
1970					
2450	2400 – 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28	28
5240					
5500	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
5785					

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**Table 3 – Test levels for Radiated and conducted Immunity Tests
Guidance and manufacturer's declaration – electromagnetic immunity**

The ENDOFLATOR® 40 is intended for use in the electromagnetic environment specified below. The user of ENDOFLATOR® 40 should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	6 V _{rms}	<p>Portable and mobile RF communications equipment should be used no closer to any part of the ENDOFLATOR® 40, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter:</p> <p>Recommended separation distance:</p> $d = 1.2 \cdot \sqrt{P}$ <p>where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b.</p>
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \cdot \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$	<p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ENDOFLATOR® 40 is used exceeds the applicable RF compliance level above, the ENDOFLATOR® 40 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ENDOFLATOR® 40.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group

Guidelines and manufacturer's declaration – electromagnetic emissions

The ENDOFLATOR® 40 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 40 should ensure that it is used in such an environment.

RF emissions as per CISPR 11	Group 1	The ENDOFLATOR® 40 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The ENDOFLATOR® 40 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5

Recommended separation distances between portable and mobile RF communications equipment and the ENDOFLATOR® 40

The ENDOFLATOR® 40 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the ENDOFLATOR® 40 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ENDOFLATOR® 40 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7177 Millcreek Drive, Mississauga, Ontario
N1W 4T8 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Total toll free phone: 1-800-268-4880 (Canada only)
Total toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoskop Suisse SA
Talavatte 5
CH-1610 Vantaa, Finland
Phone: +41 905 824774, Fax: +358 (0)968247755
E-Mail: info-ro@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestúcio st. 59 / Lenkovič st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lv@karlstorz.com

KARL STORZ Endoskop Dansk A/S
Skovlyrteren 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162809
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy UK Ltd.
4151 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B. V.
Dipsweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 033 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1815 N.W. 57th Avenue, Suite 480
Ft. Lauderdale, CA 33312-2042, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino America, Inc.
1815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8380, Fax: +1 305 262-8986
E-Mail: info-latam@karlstorz.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántico, Oficina 36
Calle 1 D e 1ra y 3ra
Col. 10400 Verdad, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: ks1.karlstorz@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Laguna Constanza No.326,
Col. 15200 Ciudad de México
Phone: +52 55 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Avda. Presidente Juan Domingo Perón 1000
C1427AAG Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 3526-4680
E-Mail: info-arg@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Ruta Joaquín Flóriano, n° 413, 20º andar - Itaim Bibi,
C.P. 0534-1011 São Paulo, Brasil
Phone: +55 11 5526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Marcha, Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigliere, 3
23135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Endoscopia (South Africa) (Pty) Ltd.
PO. Box 8061, Roggebaai 8012
Capetown, South Africa
Phone: +27 21 417 7977, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
No. 1 West Lintiang Road, Wuhou District,
610014, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 205, 20/F, 20F-N MEDIA International Center,
No. 167 Youth Avenue, Shenyang District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181119, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 1005, 10/F, 20F-N MEDIA International Center,
No. 167 Youth Avenue, Shenyang District,
110014, Shenyang, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: info-asia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: info-sia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungnareo-dong, Gangnam-gu
Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: info@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, Taipei City, Taiwan
Phone: +886 9 933 01 41 60, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Teqipin Ctr 1 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Shibata Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-00071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
15. Onion Road Lane Cove NSW 2095
P.O. Box 501 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 2503-4, 25/F AXA Tower, Landmark East,
No. 100 Hong Ming Street, Kuan Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 1805-1807, Building B, 18/F Beijing IFC
No. 8, Jiangtunmawan Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 56381788, Fax: +86 10 56381799
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shantou Branch
Room 500, Building 5 & Room 501, Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
20123, Shantou, People's Republic of China
Phone: +86 21 60338888, Fax: +86 21 60338908
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
No. 1 West Lintiang Road, Wuhou District,
610014, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
No. 167 Youth Avenue, Shenyang District,
110014, Shenyang, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: info-sia@karlstorz.com

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
No. 1 West Lintiang Road, Wuhou District,
610014, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com